OmniFilter[®]

Main Flow Bacteria Filter

REF

4-006306-00

4-006308-00

English	3
Français	7
Deutsch	11
Italiano	15
Español	19

PURITAN
BENNETT


Pressure Drop Test Setup

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
Main Flow Bacteria Filter

REF# 4-006306-00; 4-006308-00

Indicates flow direction; filter is designed for use only in this orientation.

— "Use by" mark; anticipated filter retirement date.

Directions for Use

Description

Designed as a main flow and/or exhalation filter in a heated exhalation system. The transparent plastic shell allows visual inspection of the filter.

When used as a main flow filter, *OmniFilter* minimizes particles and bacteria delivered to the patient, and reduces the possibility of bacteria reentering the ventilator, where reinfection or cross-infection might result. Used as an exhalation filter in a heated exhalation system, *OmniFilter* minimizes bacteria in exhaled gas from venting to room air. The filter also protects the ventilator's exhalation and spirometry systems, and is designed for repeated autoclavings.

CAUTION

- Using the filter in the exhalation limb of a patient circuit while using a nebulizer will trap unused medication aerosol, and cause premature buildup of resistance. This will shorten the life of the filter. To ensure proper installation, observe the flow direction markings on the filter.

Maintenance

Inspect and check the pressure drop across the filter before every reuse, and after 15 days of continuous use in the exhalation limb. Replace if cracked, if it rattles, or if its case is broken. Maximum allowable pressure drop across the filter is 4 cmH₂O (0.06 psi) at 100 L/min.

WARNING

- **Remove filter from service if pressure drop is greater than 4 cmH₂O (indicating filter medium is occluded). Failure to pressure test the filter, or to perform a ventilator check, may allow a ventilator/filter malfunction to go undetected, resulting in inadequate patient ventilation.**

Puritan Bennett (PB) recommends replacing the filter after a maximum of one year of service, or 100 autoclave cycles, whichever comes first. When putting filter into service, write the retirement date on the filter label.

PB recommends recording the number of autoclave cycles to which the filter has been subjected. The filter medium may become discolored due to repeated sterilizations and/or use of nebulized medications.

WARNING

- **Failure to detect a blocked filter may result (depending on use) in increased inspiratory or expiratory resistance, causing elevated baseline pressure. If undetected, this will block the inspiratory/expiratory path completely and cause patient injury.**

CAUTION

- Abusive handling can damage filter medium. Discard any filter that rattles, or that has bulk liquid or debris visible inside the filter shell (visible bulk liquid may indicate breakdown of the hydrophobic barrier), or visible rips or tears (the seam in the filter medium is not a defect). Disinfect filters before disposal.

Cleaning Instructions

If it is necessary to clean the filter housing, wipe with a damp cloth. Any cleaning solution used for external wipe-down must be compatible with polysulfone (the filter housing material).

WARNING

- **To avoid damage, do not immerse filter in liquid, or scrub or touch the filter medium.**
- **Do not use ethylene oxide (ETO) to disinfect bacteria filters. Always follow your hospital's infection control guidelines for handling infectious material.**

The filter is not sterile, and should be disinfected before use. It is designed for steam autoclaving only, at temperatures not to exceed 135 °C (275 °F). Note that effective disinfection occurs at a minimum cycle of 133 - 135 °C (270 - 275 °F) for 20 minutes for gravity-displacement cycles, or 4 minutes for prevacuum cycles (requires minimum 20 minute drying cycle after sterilization).

Separate couplings from filters, and wrap them in muslin or equivalent paper wrapper before autoclaving. Follow the steam autoclave manufacturer's instructions.

Note: Puritan Bennett recognizes that cleaning, sterilization, sanitation, and disinfection practices vary widely among health care institutions. It is not possible for Puritan Bennett to specify or require specific practices that will meet all needs, or to be responsible for the effectiveness of cleaning, sterilization, and other practices carried out in the patient care setting. Puritan Bennett recommends that users of its products who require cleaning and sterilization/disinfection consider the *Standards and Recommended Practices* for sterilization published by the Association for the Advancement of Medical Instrumentation (AAMI), as well as the following Center for Disease Prevention and Control (CDC) publication: *Guideline for Prevention of Nosocomial Pneumonia*.

Specifications

Connections:	Barbed fittings for patient circuits for 7200® Series, MA-1, MA-2 Series, and other ventilator applications.
Size:	3.5-in. diameter x 6.0-in. length (8.9 cm x 14.6 cm).
Filter medium:	Over 300 square in. of hydrophobic filter medium.
Filter Medium Efficiency:	When new, minimum efficiency of 99.97% retention of 0.3- μ m and larger particles at a flow of 100 Lpm, and pressure drop of approximately 2.5 cmH ₂ O at 100 Lpm. Maximum allowable pressure drop across the filter is 4 cmH ₂ O (0.06 psi) at 100 L/min.

Ordering Information

Item	REF	Description
1	4-006306-00	<i>OmniFilter</i> with Coupling
2	4-006308-00	<i>OmniFilter</i> without Coupling
3	---	0-10 cmH ₂ O gauge (not orderable from Puritan Bennett)

OmniFilter and *7200 Series* are trademarks of Puritan-Bennett Corporation.

© 2006 Puritan-Bennett Incorporated. All rights reserved.

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OmniFilter y la *Serie 7200* son marcas comerciales de Puritan-Bennett Corporation.
Tyco Healthcare UK LTD