



buddy™ 2



BELMONT INSTRUMENT CORPORATION

'buddy™ 2' Blood/Fluid Warmer

Operator's Manual

CE 0843



Belmont Instrument Corp.
780 Boston Road
Billerica, MA 01821



EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands
Tel +31(0) 70 345 8570
Fax +31(0) 70 346 7299

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician
702-00211/Rev B

Table of Contents

CHAPTER 1: SYSTEM OVERVIEW

Introduction.....	1
Indications for Use and Contraindications.....	1
System Overview.....	2
System and Temperature Monitoring.....	3
Alarms and Alarm Messages.....	3
Control Panel: Display and Keys.....	4

CHAPTER 2: OPERATION

Introduction.....	5
Step-by-Step Summary of Operating Procedure.....	7
Mounting the Power Module Unit.....	7
Connecting the Heater Unit.....	8
Installing and Priming the Disposable Set.....	8
Connect to Fluid Bag and Administration Set.....	9
Prime the Disposable Set.....	9
Initiate Infusion.....	10
Maintain Infusion.....	10
End of Procedure.....	11
Flow Characteristics.....	12
Operational Alarms and Troubleshooting.....	13

CHAPTER 3: SERVICE AND MAINTENANCE

Introduction.....	17
Cleaning, Disinfection and Maintenance.....	18
System Operational Checkout.....	19
Periodic Maintenance.....	20
Electromagnetic Compatibility.....	22

CHAPTER 4: TECHNICAL SPECIFICATIONS

Technical Specifications.....	23
-------------------------------	----

CHAPTER 5: ORDERING INFORMATION

Ordering Information.....	28
---------------------------	----

System Overview

INTRODUCTION:

The Belmont® buddy™2 in-line Blood/Fluid Warmer combines advanced microprocessor technology with an efficient heating system to provide a safe and simple system for low flow fluid warming in a miniature package. The system warms blood or replacement IV fluid to physiological temperature at flow rates up to 100 ml/min (6 L/hr) for crystalloids at 20 °C and up to 50 ml/min (3 L/hr) for packed red cells at 10 °C. The heating unit is self-contained in a hard shell case that is placed near the infusion site, minimizing IV line cool down and ensuring that fluids are delivered normothermically. The IV pole-mounted power module, supplies power to the Heater Unit, and displays the output temperature and alarm messages which alert the user to abnormal conditions such as low flow, empty set, no heat, and over and under temperature. The audible alarm may be muted, but visual indications of alarm states will always be displayed.

INDICATIONS FOR USE AND CONTRAINDICATIONS

The Belmont® buddy™2 Blood/Fluid Warmer is to be used to administer blood, blood products and intravenous solutions. It is intended to be used by healthcare professionals in a clinical environment to prevent hypothermia. The system is not meant for the infusion of granulocyte suspensions.

WARNING

The control system can only be used with the Belmont supplied disposable.

Do not alter the disposable, Heater Unit or Power Module.

Use only with crystalloid or anticoagulated blood product.

Not for infusion of drugs

Do not use in the presence of flammable anesthetic.

CAUTION

The system should not be used with pressure infusers. Blood and fluid bags contain air. The set can vent large amounts of air from crystalloid, but with blood only small amounts of air can be vented. The hand squeeze pumps, also known as a ball pump, supplied with some blood sets, may be used to momentarily increase flow, but use caution when doing so.

The system is an efficient blood and fluid warmer designed for flow rates below 100 ml/min. It is neither a resuscitation device nor a substitute for a rapid infusion system.

SYSTEM OVERVIEW

The buddy™2 Blood/Fluid Warmer complete system consists of the **Power Module (with AC/DC power supply)**, which is mounted on an IV pole, the **Heater Unit** that is placed near the patient at the infusion site, and the **disposable set**. See Figure 1.

The Power Module incorporates the power supply, display, alarm circuitry, and a magnet for storing the Heater Unit when not in use (shown in Figure 1). The bright alphanumeric display shows the output temperature of the fluid and alerts the user to alarm conditions.

The Heater Unit contains the heater plates and temperature monitoring and control system.



Figure 1: The buddy™2 Complete System

The disposable set has a sterile fluid path, and is **intended for single use only**.

The disposable set has standard luer connectors and can be connected to any standard catheter.

It contains two thin films, which contact the heater plates for rapid heating of the fluid, and two internal microporous membranes to vent out-gassed air that is generated during heating, see Figure 2.




Figure 2: The buddy™ Disposable set

The set is color-coded and keyed to ensure proper placement in the Heater Unit, blue at the input from a fluid administration set. Both heater plates are heated during operation. There is a built in valve in the disposable set to prevent air entrainment into the patient and a pressure-regulating valve at the input to protect the set and the patient from unintended exposure to high pressure applied to the IV line. This valve will allow increase flow by application of pressure, up to 300 mmHg, but will prevent pressure higher than this from reaching the set or IV line distal to it. There is also a check valve at the output to prevent a back flow.

SYSTEM AND TEMPERATURE MONITORING

Blood or replacement fluid is heated as it passes through the disposable that is in contact with the heater plates in the Heater Unit. Thermistors on the heater boards control and measure the temperature of heater plates that in turn control the temperature of the infusate. The Power Module displays the **infusate** temperature and alerts the user to alarm conditions. An over temperature condition causes the unit to stop heating and issue an audible alarm. The low temperature message (no audible alarm) informs the user that the flow rate is too rapid for the device to maintain the temperature above 31 degrees Centigrade.

ALARMS AND ALARM MESSAGES:

Under any alarm condition, the system sounds an audible alarm, flashes a red LED, and displays an alarm message on the screen. Certain alarms, such as **Over Temperature, No Heat/Check Connection, Over Current, Probe Fault** and **System Reset** indicate system problems and are always active.

The **Empty Set/Check For Air** alarm may be enabled or disabled from the front panel. If this alarm is disabled, the message is still displayed, and the red LED is lit to alert the Operator.

CONTROL PANEL: DISPLAY AND KEYS

The Power Module of the **buddy™2** Blood/Fluid Warmer is shown below:

POWER

The display is illuminated when the unit is plugged into the AC outlet.

The **ON/STANDBY** key controls power from the Power Module to the Heater Unit.

DISPLAY

The display unit shows the infusate temperature.

Under most alarm conditions the temperature display will be replaced by the alarm message. An audible alarm is activated whenever an alarm condition is detected. It can be silenced for 60 seconds by depressing the **ALARM MUTE** key.



DISABLE/ENABLE EMPTY SET/CHECK FOR AIR ALARM

If the ALARM MUTE key is depressed and held for 3 seconds at any time after power up, the **Empty Set/Check For Air** audible alarm will be **disabled** and the red LED will be illuminated. If this alarm condition occurs while the audible alarm is disabled, the red LED will flash and the alarm message will still display, but without the audible tone.

To enable the audible alarm after it has been disabled, depress and hold the ALARM MUTE key for 3 seconds. The red LED lamp within will turn off.

Depressing ALARM MUTE key when the empty set condition occurs can also disable the Empty Set/Check For Air audible alarm. Temperature is displayed on screen indicating that fluid is being heated, the red LED is lit, and the alarm message is displayed.

Operation of the System

INTRODUCTION

This section reviews the set up and operation of the Belmont® buddy™2 Blood/Fluid Warmer to ensure safe and effective operation of the system.

CAUTION:

Use of this device requires detailed familiarity with its operating procedures.

Read this manual thoroughly before operating this device.

WARNINGS:

Use only the Belmont supplied disposable. The disposable set is a single-use, sterile, non-pyrogenic fluid path.

Do not change or alter the disposable. Use extreme caution when handling the disposable; do not expose the disposable to sharp objects.

Do not use in the presence of flammable anesthetics.

Do not open the Heater Unit during or after priming.

Do not open the Heater Unit during operation, as this could damage the disposable.

Use only crystalloid to prime the unit. Do not use blood for the initial prime.

Do not attempt to prime unit or subject the disposable to flow outside of the Heater Unit, as the disposable will be damaged.

The patient line must be completely free of air before fluid is administered.

Do not use for administration of drugs, cryo-precipitate, or granulocyte suspensions.



Use only anticoagulated blood products.

Do not mix lactated Ringer's or other solutions containing calcium with citrated blood products.

Practice standard hospital policy when handling blood products. Treat blood as if it were infected and clean up all spills immediately.

To avoid risk of electric shock, this device must only be connected to a supply mains with protective earth.

STEP-BY-STEP SUMMARY OF OPERATING PROCEDURES

Operating Steps	Detailed Information
<p data-bbox="235 426 321 457">Set-up</p>  <p data-bbox="370 821 488 842">Heater Plates</p>	<p data-bbox="1024 426 1393 495">Make sure that you have all the necessary supplies.</p> <ul data-bbox="1024 537 1442 804" style="list-style-type: none"><li data-bbox="1024 537 1442 716">• Inspect the Heater Unit to ensure that the heater plates are clean and dry. Wipe off any particles or debris on the plate with a soft cloth.<li data-bbox="1024 737 1442 804">• There are no adjustments on the Heater Unit.
<p data-bbox="235 930 472 961">IV Pole Mounting</p>  <p data-bbox="370 1360 472 1413">48" above the ground</p>	<p data-bbox="1024 919 1442 1024">Mount the Power Module on the IV Pole at 48" above the ground for the operator to see and access.</p> <p data-bbox="1024 1098 1442 1266">Open the mounting screw on the rear of the unit by turning it counter clockwise until the IV pole can be placed within the holder.</p> <p data-bbox="1024 1308 1442 1455">Tighten the mounting screw by turning it in a clockwise direction until the unit is secured on the pole.</p> <p data-bbox="1024 1528 1442 1591">Plug the Power Module to the wall outlet.</p>

Operating Steps

Detailed Information

Connect Heater Unit to Power Module




Connect the Heater Unit to the Power Module by aligning the dot on receptacle to the pointed line on the plug.

To disconnect, twist the plug toward the arrow to release.

Install the Disposable Set



Line up orientation hub on the disposable with the notch in the heater unit (RED to RED). The heater unit will not close if they are not aligned.




Carefully remove the disposable from its pouch.

Install it into the Heater Unit, being careful to line up the orientation notches on the disposable with the orientation pins on the Heater Unit.

Take care not to damage the disposable.

Operating Steps	Detailed Information
 <p>Diagram showing the Heater Unit with a blue luer fitting labeled "BLUE (INPUT)".</p>	<p>The Heater Unit will not close if the disposable is not installed properly.</p>
 <p>Diagram showing the Heater Unit connected to a patient and a fluid bag. Labels include "PATIENT", "MALE LUER", "FEMALE LUER", and "CONNECT TO FLUID".</p>	<p>Close and latch the case.</p>
<p>Connect to the Fluid Bag and Administration Set</p>	<ul style="list-style-type: none"> • Hang a fluid bag and administration set. • Connect the fluid administration set to the disposable luer fitting with the blue slide clamp.
<p>Prime the Disposable Set</p> <p>WARNING:</p> <p>The case must remain closed during and after priming.</p>	<p>Open the roller clamps of the administration set and prime the unit making sure there is no air visible in the IV line or patient extension.</p> <p>Tap the Heater Unit while priming to expel any air bubbles that may collect in the disposable.</p>
<p>Establish Power to the Heater Unit</p>	<p>Depress ON/STANDBY button on the Power Module to supply power to the Heater Unit.</p>

Operating Steps	Detailed Information
<p>Initiate Infusion</p>  <p>The diagram illustrates the setup for initiating infusion. A tube with a male luer connector is connected to the patient. The tube passes through a heater unit. A bed sheet clamp is used to secure the heater unit to a flat surface. A connector for fluid is also shown.</p>	<ul style="list-style-type: none"> • Select an appropriate catheter size for the intended flow rate. • Using aseptic technique, make patient connection without entrapping air. • Adjust the roller clamp to the desired flow rate. • Secure the Heater Unit in a horizontal position, on a flat surface near the patient using bed sheet clamp ensuring that there is no strain on the patient line. • If a longer extension is desired, it should be no longer than 6 inches for optimal thermal transfer. Use of longer extensions will result in some line cooling from heater to patient. • The preferred position is at or below the level of the insertion site.
<p>Maintain Infusion</p> <p>There are no adjustments on the Heater Unit.</p> <p>You may see some condensation around the outer-edge of the Heater Unit after running for an extended period of time. This is normal and will not interfere with the system.</p> <p><i>If withdrawing fluid from distal end of disposable set via a syringe, do not exceed a rate of withdrawal of 100 ml/min.</i></p>	<ul style="list-style-type: none"> • Observe the power module display to ensure the unit is sensing flow and the flow is within the normothermic range, as indicated by the infusate temperature reading with no alarm message. • Respond to and correct system alarms.

Operating Steps	Detailed Information
<p>End of Procedure</p>	<ul style="list-style-type: none"> • Clamp roller clamp on the administration set off. • Clamp the input end of disposable set off using the blue slide clamp. • Turn Heater Unit off by pressing ON/STANDBY button on the power module. Unplug from the wall outlet. • Open the heater unit and remove disposable set. • Dispose of disposable set using hospital practices for bio-hazardous materials.

FLOW CHARACTERISTICS

The Belmont® **buddy™ 2** Blood/Fluid Warmer is designed for gravity flow. The maximum flow rates achievable will be determined by the size of the catheter employed.

The disposable has been designed to withstand up to 300 mmHg infusion pressure, provided the disposable set is within the Heater Unit and latched. The system is designed to maintain normothermic temperature at flow rates specified (50 ml/min for undiluted RBC's at 10°C, and 100 ml/min for 20°C crystalloid). However, at higher flows, the fluids administered may not heat to normothermia. **There are no adjustments for the Heater Unit.**

Pressure Drop:

The disposable and heating unit does not restrict flow. If flow appears to be restricted:

- Check the infusion site to ensure catheter has not infiltrated the surrounding tissue
- Check to make sure line is not crimped or clamped off
- Out gassed fluids may be causing air build up in the disposable. If air restriction is suspected:
 - Detach the unit from the patient connection
 - Holding the case with the patient connection port upward, tap the case gently while allowing a few seconds of free flow to remove air. Reattach the unit to the patient with a fluid-to-fluid connection.
 - If flow continues to be restricted, replace the disposable set.

Air Venting

As fluid flows through the disposable set, it passes across microporous membranes, which can vent air from the fluid. The air is then released through the side vents of the set. **When infusing crystalloid, the set can vent large amounts of air. With blood or blood products, the set can vent only small amounts of air. CAUTION: When using blood or blood product, the Operator should be vigilant to avoid infusing air especially if using a hand squeeze pump during the infusion.** If the set were to become filled with air, the system will activate the EMPTY SET/CHECK FOR AIR alarm to alert the Operator to the problem. The fluid path of the set is only about 7 ml, and therefore, only a few ml of air will trigger the alarm. In addition, there is a pressure-regulating valve, incorporated into the disposable set, prevents potential air entrainment into the system by stopping fluid flow if a partial vacuum condition is present.

During **normal operation** – water vapor, along with out-gassed air from warmed fluid, will be vented and will cause condensation to appear within the side vents after prolonged use. This should **not** be mistaken as a seal failure due to materials or workmanship.

OPERATIONAL ALARMS AND TROUBLESHOOTING

The Belmont® buddy™ 2 Blood/Fluid Warmer detects and displays the following alarm conditions:

Empty Set/Check for Air
Low Flow
Low Temperature
No Heat, Check Connection
Over Current
Over Temperature
Probe Fault
System Reset

The chart below details these alarms, the system response and the recommended user actions.

Alarm Message	Condition and System Response	Probable Causes	Actions
Empty Set/Check for Air	<p>System detects an empty set condition. Flash LED, display message, and activate audible alarm. Heater continues to function normally.</p> <p>Audible alarm may be disabled in advance by depressing ALARM MUTE for 3 seconds or depressing ALARM MUTE key when condition occurs. Flashing LED and message will occur at alarm condition.</p> <p>System automatically resets and continues normal operation when air is cleared from the disposables.</p>	<ol style="list-style-type: none"> 1. Fluid source is empty or air in the fluid path. 2. Line is clamped. 	<ol style="list-style-type: none"> 1. If appropriate, replace fluid supply. 2. Open clamps.
Low Flow	<p>System automatically resets and continues normal operation when flow resumes.</p> <p>Display message, no audible alarm. Heating continues.</p>	<ol style="list-style-type: none"> 1. Fluid source is empty. 2. Line is clamped 3. Flow is too low to be detected. 4. Disposable has clogged and is restricting flow. 	<ol style="list-style-type: none"> 1. If appropriate, replace fluid supply. 2. Open clamps. 3. Increase flow rate until flow is detected. 4. Replace disposable.

Alarm Message	Condition and System Response	Probable Causes	Actions
Low Temperature	System detects fluid temperature below 31EC for 10 sec. Display message, no audible alarm . Heating continues.	Flow rate is set too fast.	Slow the flow rate using roller clamp on the administration set.
No Heat Check Connection	System detects that communication is lost between the Heater Unit and power module. Stop heating, flash LED, alarm, and display message. If connection is restored, system will resume normal operation. If not, the system turns power off and displays 'Check Connection, Press ON/STANDBY' message.	1. Cable between the Heater Unit and Power Module is disconnected or loose. 2. Wire is broken.	1. Check and secure connection between power module and Heater Unit. If this fails to rectify the problem obtain a replacement Heater Unit. 2. Replace unit.
**Over Current	System detects short circuit in heater or heater cable. Stop heating, flash LED, alarm, and display message.	1. Heater circuit board defective. 2. Wire is broken.	1. Unplug power module from wall outlet, and re-power. 2. Call service.
**Over Temperature	System detects fluid temperature (i) Moderate to High Flow: >42.5EC for 2 min, or ≥ 48EC for 2.5 sec, or (ii) Low Flow: >42EC for 2.5 sec or ≥ 44EC for 1 sec, or (iii) Output Plate Temperature reaches 50.5EC. Stop heating, flash LED, alarm, and display message.	1. Repeatedly stopping and restarting during high flow conditions. 2. Failure of the microprocessor or temperature sensor.	1. Clamp line, turn unit off and unplug power module from wall outlet, and re-power. 2. If failure reoccurs obtain replacement unit and contact service.
**Probe Fault	Internal malfunction.	Hardware temperature override circuit not working properly.	Unplug power module from wall outlet, and re-power. If not successful, call service.

Alarm Message	Condition and System Response	Probable Causes	Actions
**System Reset	Internal computer malfunction. Alarm continuously and display message.	Loose watchdog function on the Power Module Processor.	Unplug power module from wall outlet, and re-power. If not successful, call service.

**** At these four (4) alarms, system stops heating, flashes a red LED, alarms, and displays message. Unplug Power Module from wall outlet, and re-plug into the wall outlet. Press ON/STANDBY to apply power to the Heater Unit.**

USER MANAGEMENT OF ALARMS

1. To silence an audible alarm, depress the ALARM MUTE key. If the alarm condition persists, the audible alarm will resume after 60 seconds, except the EMPTY SET/CHECK FOR AIR alarm, which will remain disabled until air disappears and flow resumes for at least 1 minute after which it can again be activated automatically.
2. The EMPTY SET/CHECK FOR AIR audible alarm may be disabled by depressing the ALARM MUTE key for 3 seconds or depress ALARM MUTE key when this condition occurs. The red LED will be lit steadily indicating that the alarm tone has been disabled. If the Empty Set condition occurs during this time, the LED will flash and message will appear on the screen.

To re-enable the alarm tone, depress and hold ALARM MUTE key for 3 seconds, the red LED will no longer lit indicating that normal operation with alarm tone enabled have been restored.

3. The LOW TEMPERATURE message is not an alarm condition and does not produce an alarm tone. It usually indicates that the flow rate is set too fast; reduce the flow to get rid of the message.
4. To reset the NO HEAT/CHECK CONNECTION alarm after the system turns the power off, press ON/STANDBY button on the Power Module to turn the heater ON to resume operation.
5. The following alarms indicate a probable system fault: OVER CURRENT, OVER TEMPERATURE, PROBE FAULT, and SYSTEM RESET.

To reset or retry, remove Power Module from the wall outlet then re-connect to the wall outlet. Press ON/STANDBY key to turn the heater on to resume operation.

Service and Maintenance

INTRODUCTION

The Belmont® **buddy™ 2** Blood/Fluid Warmer has no user serviceable parts, and does not require calibration. Routine maintenance consists of cleaning and inspecting the Heater Unit after each use and testing the system operational and routinely test for leakage current to ensure against electrical shock hazard. The circuit diagram, component part lists, descriptions, and calibration instructions will be provided to assist service personnel in device repair.

WARNING:

Test leakage current periodically to insure against electrical shock hazard

CAUTION

Do not open the Power Module. There are no serviceable parts inside the housing. Opening the system will void the warranty.

CAUTION

The heater plates in the heater unit are sealed. Inspect the seals to insure their integrity after each use. If the seal is not intact, do not use. Do not try to open or to service the heater unit.

CAUTION

Check the integrity of the wires between the heater unit and power module. If there is evidence of fraying or damage do not use and return the unit for service.

CLEANING, DISINFECTION AND MAINTENANCE

ROUTINE MAINTAINENCE

Turn heater to **STANDBY**, unplug the power cord from the wall outlet.

1. **Clean and/or Disinfect Exterior:**

A. Power Module

The Power Module may be wiped down with water or a 10% bleach solution using a soft cloth. Take particular care not to scratch the display window.

Note: Avoid the use of acetone or other solvents that may damage the surface.

B. Heater Unit

Wiping down the outside surfaces with a soft cloth moistened with water or a 10% bleach solution.

Note: Avoid the use of acetone or other solvents that may damage the surface.

Open the Heater Unit and inspect the seal around the heater plates. If the seal is intact, using a soft cloth wipe the heater plates with a 10% bleach solution and dry. If the seal around the heater plates is **not** intact, do not use and return the unit for service or replacement.

Note: Do not use abrasive cleaners on the Heater Unit or heater plates.

WARNING

Do not submerge the heater unit.

Practice standard hospital policy when handling blood products. Treat blood as if it were infected and clean up all spills immediately.

2. **Visual Inspection:**

Once the unit is cleaned, inspect all wiring and connections for fraying or worn wiring. If noted, do not use and return the unit for service. If the wear is confined to the Heater Unit, replace the heater. All Heater Units are interchangeable with the Power Module units.

SYSTEM OPERATIONAL CHECKOUT

This system operational checkout including the alarm test should be done periodically. Prior to performing the test, have the following equipment and supplied on hand:

- IV Pole
- Saline or other crystalloid for testing
- Disposable Set
- Administration set
- Graduated cylinder and timing device for measuring flow rate

Setup:

1. Mount the Power Module on the IV Pole at a comfortable height for the operator to see and access using the mounting screw at the rear of the unit.
2. Plug the Power Module to the wall outlet.
3. Connect the Heater Unit to the Power Module by aligning the dot on receptacle to the pointed line on the plug.
4. Install disposable set. The Heater Unit will not close if the disposable set is not installed properly. Close and latch the case.

Connect to the Fluid bag and Administration Set:

1. Hang a fluid bag and administration set. Connect the fluid administration set to the disposable luer fitting with blue slide clamp. Connect administration set with 14 gauge catheter to the output luer fitting.
2. Open the roller clamps and prime the unit. Tap the Heater Unit while priming to expel any air bubbles that may collect in the disposable.
3. Secure the Heater Unit in a horizontal position, making sure that the output of the catheter is at least 5" higher than the Heater Unit to simulate a venous pressure of approximately 10 mmHg.

Operational Test:

1. Depress ON/STANDBY button on the Power Module to supply power to the Heater Unit. Set the flow rate to approximately 80 ml/min.
2. Observe the Power Module display to ensure the unit is sensing flow and the flow is within the normothermic range, $38^{\circ} \pm 2^{\circ}\text{C}$, with no alarm message.
3. Briefly stop the flow by closing off the roller clamp. Verify the 'LOW FLOW' message.
4. Resume flow and infuse until the fluid bag is empty. Squeeze air in the set and verify that when the disposable is filled with air, the system sounds an audible alarm with 'EMPTY SET/CHECK FOR AIR' message display on screen and the red LED flashes.
5. Turn Heater Unit off by pressing ON/STANDBY button on the Power Module. Unplug from the wall outlet and discard the disposable set.

PERIODIC MAINTENANCE

An electrical safety test should be performed at least once a year by a qualified technician. Prior to performing the electrical safety test, have the following equipment and supplied on hand:

- Electrical Safety Analyzer
- Saline or other crystalloid for testing
- Disposable Set
- Administration set

Set-Up: Plug the **buddy™2** Blood/Fluid Warmer into AC outlet on the front of Safety Analyzer.

CAUTION:

Before applying voltage to Safety Analyzer, make sure input line voltage is correct for the **VOLTAGE OF UNIT UNDER TEST**

1. **Ground Wire leakage Current:** This leakage current is done when there is no equipotential connector on the device under test.
 - a. Plug the Safety Analyzer into an appropriate power source. Turn power ON.
 - b. Press ON/STANDBY key on the **buddy™2** to turn the heater on.
 - c. Press the **gnd leakage** key on Safety Analyzer. Do not attach the large red clamp to the device under test.
 - d. The **µA**, **gnd leak**, and **open gnd** LEDs will illuminate. Record the leakage current displayed.
 - e. All measurements should be $< 300 \mu\text{A}$ (for Domestic unit) and $< 500 \mu\text{A}$ (for 230V unit).
 - f. Press ON/STANDBY key on the **buddy™2** to turn the heater off.
 - g. Turn the Safety Analyzer power off.

2. Patient Leakage Current

- a. Install the disposable set.
- b. Attach 12 to 16 gauge stainless steel catheter to the end of the disposable set tubing.
- c. Install saline solution and the administration set. . Open the roller clamp and adjust for maximum flow.
- d. Turn the Safety Analyzer power on.
- e. Press ON/STANDBY key on the **buddy™2** to turn the heater on.
- f. Attach large clamp (red lead) from Safety Analyzer to the catheter.
- g. Press the **chassis leak** key on Safety Analyzer.
- h. The **μA** LED will illuminate. Record the leakage current displayed for each of the following conditions, with Neutral switch in NORM position.

<u>Polarity</u>	<u>Ground</u>	<u>Tolerance Limit (μA)</u>
NORM	NORM	<10 μA
REVERSE	NORM	<10 μA
REVERSE	OPEN	<50 μA
NORM	OPEN	<50 μA

- i. Turn the Safety Analyzer power off.
- j. Press ON/STANDBY key on the **buddy™2** to turn the heater off. Close the roller clamp. Remove the disposable set.
- k. Remove the **buddy™2** from the Safety Analyzer.

ELECTROMAGNETIC COMPATIBILITY

The Belmont® buddy™2 is intended for use in the electromagnetic environment specified below. You should assure that it is used in such environment.


Test	Compliance	Electromagnetic environment - guidance
RF emissions, CISPR 11	Group 1, Class B	The buddy2 may emit electromagnetic energy so that nearby electronic equipment may be affected.
Harmonic emissions IEC 61000-3-2		
Voltage fluctuations/ flicker emissions IEC 61000-3-3		
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated RF IEC 61000-4-3	3 V/m over the frequency range of 80 MHz to 1 GHz	Radiated RF frequency is negligible at any single frequency.
Electric Fast Transient Burst Immunity IEC 61000-4-4	± 2kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Fast Surge Immunity IEC 61000-4-5	± 2kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
RF Common Mode Immunity IEC 61000-4-6	3 Vrms	If interference occurs, it may be necessary to position the buddy2 further from sources of power frequency magnetic field.
Power Frequency Magnetic Field Immunity IEC 61000-4-8	3 A/m	If interference occurs, it may be necessary to position the buddy2 further from sources of power frequency magnetic field.
Voltage Dip and Interrupt Immunity IEC 61000-4-11	> 95% dip for 0.5 cycle > 60% dip for 5 cycles > 30% dip for 25 cycles > 95% dip for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the power main interrupts, the system automatically switches to battery mode.

Technical Specifications

DIMENSIONS	
Size (W x H x D)	Power Module: 4.33" x 7.48" x 2.39" AC/DC Supply: 4.25" x 8.42" x 1.85" Heater Unit: 5.1" x 0.89" x 1.7"
Weight	Power Module: 1.00 lbs AC/DC Supply: 2.65 lbs Heater Unit: 3.5 oz
POWER	
AC Input Voltage, Universal	100 V~ to 240 V~
Operating Frequency	47 to 63Hz
Maximum Power	210 VA
Ground Wire Leakage Current	<300 μ A (For Domestic Unit) < 500 μ A (For 240 V~ Unit)
Patient Leakage Current (Open Ground)	< 50 μ A
Electrical Compliance	EN 60601-1, UL60601-1, CAN/CSA C22.2-NO. 601.1-M90
Power Cord	3/18 AWG, type SJT, 3-wire grounded with hospital grade plug

SAFETY AND MONTORING	
Infusate Temperature	Heater power shut off if fluid temperature exceeds specified limits
Independent Safety Circuit	Shut off the heater if temperature exceeds a limiting value
ALARMS/ALERTS	
System Status (Alert, no audible alarm)	LOW FLOW LOW TEMPERATURE
Operator Setting, User Correctable (Alarm and display message)	EMPTY SET/CHECK FOR AIR NO HEAT, CHECK CONNECTION CHECK CONNECTION, PRESS ON/STANDBY
System Failures (Alarm and display message)	OVER CURRENT OVER TEMPERATURE PROBE FAULT SYSTEM RESET
DISPOSABLE SET	
Disposable Set	Sterile, Non-Pyrogenic Fluid Path, Single-Patient Use Only, Not made with natural rubber latex
Sterilization Method	Ethylene Oxide
Sterility Assurance Level	Greater than or equal to 10 ⁻⁶

CLASSIFICATION	
Type of Protection Against Electrical Shock	Class I
Degree of Protection Against Electric Shock	CF defibrillator-proof
Degree of Protection Against Harmful Ingress of Water	IPX1, Drip proof
Degree of Safety in Presence of Flammable Anesthetics	Not suitable
Mode of Operation	Continuous
<p style="text-align: center;">Medical Equipment</p>   <p style="text-align: center;"> UL 60601-1 CAN/CSA-C22.2 No. 601.1 ANSI/AAMI/ES60601-1 (2005) CAN/CSA-C22.2 No. 60601.1 (2008) </p>	<p>Medical – General Medical Equipment</p> <p>As to electrical shock, fire and mechanical hazards only in accordance with</p> <p>UL 60601-1, CAN/CSA-C22.2 No. 601.1, ANSI/AAMI/ES60601-1 (2005, 3rd ed.), CAN/CSA- C22.2 No. 60601-1 (2008).</p>
Medical Device Directive: Council Directive 93/42/EEC	Hardware: Class IIb
	Disposable Set: Class IIa

SYMBOLS AND DEFINITIONS	
Symbol	Description
	Compliance to Medical Device Directive 93/42/EEC and 2011/65/EU
	Alternating current
 or 	Attention, consult accompanying documents/refer to manual
	Defibrillator-proof type CF equipment
IPX1	Protected against dripping water
SN	Serial Number
	Class I equipment: the protection against electric shock is provided by a combination of basic insulation and the use of the EARTH connection
	Do not immerse into water or any other liquid
	Manufactured by
	Authorized European Representative

Ordering Information

To order buddy System and the Disposable Set call or write the following:

Belmont Instrument Corporation
780 Boston Road
Billerica, MA 10821, USA

T: (866) 663-0212 US/Canada

T: (978) 663-0212 Worldwide

F: (978) 663.0212

Email: sales@belmontinstrument.com

REF: 905-00037 buddy™2 System, 100-240 VAC, English

REF: 905-00010 buddy Family Disposable Set