GE Healthcare

Optima EMS ECG Management System InSite 3.5 Installation and Configuration Guide

Software Version 1.0
2069773-004 Revision C

Optima EMS
English
© 2012-2013 General Electric Company.
All Rights Reserved.

Publication Information

The information in this manual applies only to Optima EMS ECG Management System Version 1.0. It does not apply to earlier product versions. Due to continuing product innovation, specifications in this manual are subject to change without notice.

MUSE and MAC are trademarks owned by GE Medical Systems *Information Technologies*, Inc., a General Electric Company going to market as GE Healthcare. All other trademarks contained herein are the property of their respective owners.

The document part number and revision appear at the bottom of each page. The revision identifies the document's update level. The revision history of this document is summarized in the following table.

Revision	Date	Comments
A	1 November 2012	Internal Release
B	12 December 2012	Initial release for customer
C	22 July 2013	Per SPR HCSDM00221958, added the Install pre-requisites. Updated sections: Original Serial Number and Uninstalling InSite.

To access other GE Healthcare Diagnostic Cardiology manuals, go to the Common Documentation Library (CDL), located at www.gehealthcare.com/documents, and click **Cardiology**.

To access Original Equipment Manufacturer (OEM) manuals, go to the device manufacturer's Web site.

Service Manual Language Information

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Contents

1	Introduction	
	Manual Information	11
	Intended Audience	11
	Manual Purpose	11
	Document Conventions	11
	Responsibility of the Purchaser/Customer	12
	Additional Assistance.....	12
	Installation Pre-requisites.....	13
2	InSite 3.5 Configuration	
	Serial Number/Device Name/CRM Number	15
	Original Serial Number (from serialNo.txt)	15
	Device Name (Auto-Generated).....	16
	CRM Number (Auto-Generated).....	16
	Agent Configuration.....	17
	Device Type.....	17
	Device Prefix.....	17
	Product Name.....	17
	Model Type	17
	Version	18
	File Watcher	18
	File Transfer	18
	Poll Rate	19
	Default Enterprise Server	19
	Required Properties	19
	InSite Service Platform Files Customized/Changed for Product Integration	20
3	Installing InSite 3.5	
	Before You Begin	21
	Configuring InSite 3.5.....	22
	Uninstalling InSite	25
4	InSite BackOffice	
	Remote Desktop Access	27
A	InSite Log Level	

Reducing the InSite Log Level	31
Adding FORCE to Log Level Choices	31
Configuring the Log Level to FORCE	33

1

Introduction

This chapter provides general information required for the proper use of the Optima EMS system (also referred to as "the system" throughout this manual) and this manual. Familiarize yourself with this information before using the system.

InSite 3.5 is the remote service application on the system enabling remote access to the system by a GE Healthcare service representative. Access to the system via InSite 3.5 is fully secured and available only to authorized GE Healthcare service representatives.

Manual Information

This section provides information for the correct use of this manual.

Keep this manual with the equipment at all times and periodically review it. You should request training assistance from GE Healthcare, if needed.

Intended Audience

This manual is intended for use by service personnel and IT personnel.

Manual Purpose

The purpose of this manual is to provide instruction and guidance during the installation and configuration of InSite 3.5.

Document Conventions

This manual uses the following conventions.

Typographical Conventions

Convention	Description
Bold Text	Indicates keys on the keyboard, text to enter, or hardware items such as buttons or switches on the equipment.
<i>Italicized-Bold Text</i>	Indicates software terms that identify menu items, buttons or options in various windows.

Convention	Description
CTRL+ESC	Indicates a keyboard operation. A plus (+) sign between the names of two keys indicates that while holding the first key, you should press and release the second key. For example, Press CTRL+ESC means to press and hold the CTRL key and then press and release the ESC key.
<space>	Indicates that you must press the spacebar. When instructions are given for typing a precise text string with one or more spaces, the point where you must press the spacebar is indicated as <space> . This ensures that the correct number of spaces is inserted in the correct positions within the literal text string. The purpose of the < > brackets is to distinguish the command from the literal text within the string.
Enter	Indicates that you must press the Enter or Return key on the keyboard. Do not type <i>Enter</i> .
>	The greater than symbol, or right angle bracket, is a concise method to indicate a sequence of menu selections. For example, the statement "From the main menu, select System > Setup > Options to open the Option Activation window" replaces the following: <ol style="list-style-type: none"> 1. From the main menu, select System to open the System menu. 2. From the System menu, select Setup to open the Setup menu. 3. From the Setup menu, select Options to open the Option Activation window.

Illustrations

All illustrations in the manual are provided as examples only. Depending on system configuration, screens in the manual may differ from the screens on your system.

All patient names and data are fictitious. Any similarity to actual persons is coincidental.

Notes

Notes provide application tips or additional information that, while useful, are not essential to the correct operation of the system. They are called out from the body text through a flag word and indentation, as follows:

NOTE:

The tip or additional information is indented below the **NOTE** flag word.

Responsibility of the Purchaser/Customer

The customer is responsible for providing appropriate desks, chairs, electrical wall outlets, network connections, analog phone lines, and for locating any of the Optima EMS System components described in this manual in compliance with all local, state, or national codes.

Additional Assistance

GE Healthcare maintains a trained staff of application and technical experts to answer questions and respond to issues and problems that may arise during the installation, maintenance, and use of this system.

Contact your local GE Healthcare representative to request additional assistance.

Installation Pre-requisites

Before installing *InSite Integration*, make sure you have installed the following software:

1. **Net framework 4.0**
You can find the installation file in [OptimaEMS CD]\Microsoft .NET Framework\
2. **VC++ 2010**
You can find the installation file in [OptimaEMS CD]\VC Redistributable Maintenance\
3. **MSXML sp2**
You can find the installation file in [OptimaEMS CD]\InSite Preinstallation\
4. **vcredist_x86_sp1-2.0.50757.4053.exe**
You can find the installation file in [OptimaEMS CD]\InSite Preinstallation\

NOTE:

If you have installed requisites 1 and 2 when installing the Optima EMS system, you may ignore them and continue with 3 and 4.

2

InSite 3.5 Configuration

Before installing InSite 3.5, check and update various configurations in the file *InstallOption.xml*.

Serial Number/Device Name/CRM Number

In order to populate the device name and CRM number used in the BackOffice, the serial number of the device must first be extracted. InSite 3.5 provides a mechanism to read and parse this information from *D:\serialNo.txt*. It is the product team's responsibility to populate this file with the appropriate serial number. Only the single string should be present. No delimiting characters are necessary.

Upon reading the *serialNo.txt*, the device name and CRM number used in the connectivity configuration are assembled as follows.

Original Serial Number (from serialNo.txt)

You need to create a TXT file named *serialNo*. Type the product serial number in to the file and save it to the D drive.

To change the **serialNo** file default (D:**serialNo.txt**), use the following steps:

1. Open the disk where you installed the Optima EMS system server and navigate to the **Insite Installation** folder.
2. Open **GetSerialNumber.bat** with Notepad.
3. Change the value of **SerialNoFilePath** to a new path.

NOTE:

The file name **serialNo.txt** cannot be changed, you can only change the file path.

```
@echo off
rem Returns the serial number from "serialNo.txt" file
rem Returns the serial number up to the manufacturer code if "-1"
set SerialNoFilePath=D:\serialNo.txt
if not exist "%SerialNoFilePath%" (
    echo Failed Finding %SerialNoFilePath% 1>&2
    goto END
)

FOR /F "usebackq tokens=1*" %%a IN ("%SerialNoFilePath%") DO (
    set SerialNumber=%%a
)

if "%SerialNumber%"==" " (
    goto FAILED
)

if "%SerialNumber%"=="-1" (
```


NOTE:

SerialNo is the serial number of the system without the Prefix **OP**.

Device Name (Auto-Generated)

You can find the auto-generated Device Name in "Configuring InSite 3.5" on page 22.


The following graphic shows the generation rule of the **Device Name**.

CRM Number (Auto-Generated)

You can find the auto-generated CRM Number in "Configuring InSite 3.5" on page 22.

The following graphic shows the generation rule of the **CRM Number**.

NOTE:

CRM number is not the same as the system serial number.

The CRM number **MUST** match the system ID registered in the CRM database at the time of sale.

If the device name and CRM number are not auto-generated, make sure that the default names are unique, for example, not **UNKNOWN** but **UNKNOWN+timestamp**. A common occurrence for this is if the **serialNo.txt** is not yet populated. For older releases, you may need to enter **GetDeviceInfo.pl** so that the CRM number is generated.

Agent Configuration

The following items must be negotiated in order to register with the BackOffice.

Device Type

Device Type is comprised of three components: **Business**, **Modality**, and **Product** and you must configure it in the file **InstallOptions.xml**.

For example, **CS_DCR_OPTIMA_EMS_3.5**.

The system uses the following device types:

```
<EnvVar Varname="DeviceType">CS_DCR_OPTIMA_EMS_3.5</EnvVar>
```

Device Prefix

The prefix used for the system is: **OP**.

The following element was modified in **InstallOptions.xml**:

```
<EnvVar varname="Prefix">OP</EnvVar>
```

Product Name

The product name is **Optima**.

The following element was modified in **InstallOptions.xml**:

```
<EnvVar varname="ProductName">Optima</EnvVar>
```

Model Type

The model type is **EMS** and the following element was modified in **InstallOptions.xml**:

```
<text symbol="__SA_ASSET_TYPE_MODEL__">EMS</text>
```

Version

The version is configured as follows in *InstallOptions.xml*:

```
<Version>3.5</Version>
```

File Watcher

The export directory is defined in *InstallOptions.xml* (located in the *Insite Installation* folder where you installed the Optima EMS system server) as follows:

```
<EnvVar varname="ExportDir">D:\export</EnvVar>
```

The file watcher is created in the following directory:

```
<text symbol="__FILE_WATCHER_ALIAS__">Export</text>
```

```
<text symbol="__FILE_WATCHER_DIR__">%ExportDir%</text>
```

The alias represents the string that is visible from the BackOffice. You may alter this if you want.

You can also add additional directories by defining them as virtual. For example:

```
<EnvVar varname="LogDir">D:\log</EnvVar>
```

```
<text symbol="__VIRTUAL_DIR_ALIAS__">Log</text>
```

```
<text symbol="__VIRTUAL_DIR__">%LogDir%</text>
```

Finally, the file filter **.zip* is used.

```
<text symbol="__FILE_WATCHER_FILTER__">*.zip</text>
```

The configuration shown here causes the agent to upload files ending in *.log* from the *D:*. There are no automated dumps of *.log* files into the *Export* directories.

NOTE:

For the system, there is a separate log folder named *Log* in the system server installation folder.

File Transfer

This section explains how to configure multiple file transfer folders (to transfer files between BackOffice and the device), and applies only if you need to configure multiple file transfer folders. only

NOTE:

By default the File Repository directory should be configured to *%INSITE2_DATA_DIR%\etc*. The service team uses this directory at the BackOffice for online support. You should keep this directory as it was configured and create a new file transfer folder as the file repository location.

Do not use Windows System32 Config Folder (*C:\Windows\system32\Config*) directly for file transfer.

When you try to upload the *app event log* file, the process of uploading the file causes the file to be updated and prevents the upload task from being granted access to the file. The task manager locks the file to populate the *Questra* file upload messages into the log.

Poll Rate

The poll rate for a headless server product must be set to 60 seconds.

Default Enterprise Server

The default URL as it appears in **Clearcase** must be **Production**. You can change this to **Development** or **Pilot** as needed but the default must always remain **Production**. It is configured in **InstallOptions.xml** as follows:

```
<text symbol="__ENT_URL__">https://us0-ws.service.gehealthcare.com:443</text>
```

```
<text symbol="__TUN_URL__">https://us0-rd.service.gehealthcare.com:443</text>
```

Required Properties

The following static registration properties are negotiated and defined as part of the product team's **InstallOptions.xml**. These property values are consistent across the device type's install base.

Name	Description	SITEMAP Tag
CS_GSP_MODALITY	DCAR	__SA_MODALITY__
CS_GSP_POLLTYPE	CONSTANT: polls continuously	__SA_POLLTYPE__
CS_GSP_OPSYSTEM	WINOTHER	__SA_OP_SYSTEM__
CS_GSP_OPREV	Revision of device's Operating System	__SA_OP_REVISION__
CS_GSP_APREV	Revision of device's Application Software	__SA_AP_REVISION__
CS_GSP_SRVREV	Revision of device's Service Platform	__SA_SRV_REVISION__
CS_GSP_SERVICE_CENTER	Online Centers across CS poles JROC: Jupiter, FLA	__SA_SERVICE_CENTER__
CS_GSP_HRFS_SUPPORT	FALSE : Not supported	__SA_HRFS_SUPPORT__
CS_GSP_MRFS_SUPPORT	FALSE : Not supported	__SA_MRFS_SUPPORT__
CS_GSP_SWRLOAD_SUPPORT	FALSE : Not supported	__SA_SWRLOAD_SUPPORT__
CS_GSP_SWDLLOAD_SUPPORT	FALSE : Not supported	__SA_SWRLOAD_SUPPORT__

Checklist

- Device Type** defined in **InstallOptions.xml**
- Device Prefix** defined in **InstallOptions.xml**
- Product Name** defined in **InstallOptions.xml**
- Model Type** defined in **InstallOptions.xml**
- Version** defined in **InstallOptions.xml**

Checklist (cont'd.)

- Export Dir** defined in *InstallOptions.xml*. File update workflow is defined
- Poll rate** set to 60 seconds
- Default Enterprise** set to *Prod* in *InstallOptions.xml*
- Required Properties** defined in *InstallOptions.xml*
- Optional Properties** defined in *InstallOptions.xml*

InSite Service Platform Files Customized/ Changed for Product Integration

This section explains the list of **InSite Service Platform** files customized/changed for the system.

Most files were changed by the engineering team. To install InSite 3.5 with no issues, check them first.

Item	File Name	Location	Change(s) Description
1	<i>InstallOption.xml</i>	\	Config changes for SerphyDose devicetype
2	<i>GetSerialNumber.bat</i>	\	To locate SerialNo.txt file in C:\SerphyDose directory
3	<i>GetDeviceInfo.pl</i>	\	Updated so CRM number=serial number
4	<i>hscds.xml</i>	\AgentInstall\templates	Removed the telnet and CSD connecttype

3

Installing InSite 3.5

Before You Begin

NOTE:

Before you install InSite 3.5, make sure you have installed the Optima EMS system server successfully.

1. Create a **serialNo.txt** file on the D drive.
See [“Original Serial Number \(from serialNo.txt\)” on page 15](#) for detailed information.
2. Create an **export** folder on the D drive.
See [“File Watcher” on page 18](#) for detailed information.

3. Run **InsiteWrapper.exe** in the **Insite Installation** folder where you installed the Optima EMS system server.

If the installation is successful, the following window opens and closes automatically.


```

c:\ Command Prompt
Configuring Questra Agent using InstallOption.xml
 1 file(s) copied.
 1 file(s) copied.
 1 file(s) copied.
 1 file(s) copied.
 1 file(s) copied.
 1 file(s) copied.
 1 file(s) copied.
 1 file(s) copied.
 1 file(s) copied.
 1 file(s) copied.
 1 file(s) copied.
Configuring installed components using InstallOption.xml
Configuring components
Updating UNC registry settings
Installed Service Applications Successfully.

Starting Service Platform components that are configured to start automatically.

Successfully installed Service Platform.

*****
See the installation details in "C:\InSite2\Install.log"
Removing the install files and directory

```

4. If the service is stopped, perform the following steps to restart the service manually.
 - a. Right-click on the **InSite utility** icon  in the system tray.
 - b. Select **Start InSite 3.5 and Remote Desktop Services**.

The service is now started.

Configuring InSite 3.5

1. Right-click on the **InSite utility** icon  in the system tray.
The utility menu opens.
2. Select **Configuration Tool...**
The **InSite 3.5 Configuration Tool** window opens.

Agent Configuration

Device Name: OP_SJE12090001W **CRM No.:** SJE12090001WP

Display Name: Description:

Continent: ASIA **Country:** CHINA

Addr Line1:

Addr Line2:

City: WUXI **State (Prov):** JIANGSU **Postal Code:**

Latitude: Longitude:

Institution: GE HEALTHCARE **Department:**

Building: **Floor:** **Room:**

Advanced Configuration

Enterprise Server: PRODUCT **Service Center:** OTHERS **Log Level:** INFO

Enterprise Server URL: https://us0-ws.service.gehealthcare.com:443

Enterprise Tunnel URL: https://us0-rd.service.gehealthcare.com:443

File Repository: C:\InSite2\Questra\GeHealthcare\Agent\etc

File Watcher: Enable **Dir:** D:\export **Filter:** *.totrans

Proxy Configuration

Proxy: Disable

Name: Proxy Disabled

IP Addr: **Port:**

Proxy Authentication: Disable **Scheme:** NONE

Prozy User: **Password:**

If you do not need to change anything in the **InSite 3.5 Configuration Tool** window, review the information and then close the window.

If you do need to change the information in the window, the fields in bold text are required.

- Complete the required fields (in bold) on the **InSite 3.5 Configuration Tool** window using the following table as a guide. **Enterprise Server URL**, **Enterprise Tunnel URL** and **File Repository** are automatically set when you select other fields and you should not change them. The rest of the fields on the window are not required and you can leave them blank.

InSite 3.5 Configuration Tool

Field	Description/Value
Agent Configuration	
Device Name	Unique identifier for the server in the BackOffice. This field should already be filled in with the prefix and serial number that were entered during installation of InSite.
CRM No.	The number the Customer Relationship Management system (Seibel) uses to track service and warranty on the product. This field should already be filled in with the serial number entered during the installation of InSite.



InSite 3.5 Configuration Tool (cont'd.)

Field	Description/Value
Continent	Use the drop-down arrow to select the continent where this system is located.
Country	Use the drop-down arrow to select the country where this system is located.
City	Type the name of the city where this system is located.
State (Prov)	Use the drop-down arrow to select the state (or province) where this system is located.
Institution	Type the name of the institution where this system is located. This field should match the Customer Name in the Seibel Service Request, minus the leading customer number. For Asia and Europe, follow the naming convention specified by your service operations team.
Advanced Configuration	
Enterprise Server	Select Product for a server installed at a customer site. This selection determines with which BackOffice server the product is registered.
Service Center	For Americas, select JROC . For Europe, select EURO . For Asia, select OTHER .
File Watcher	Use the drop-down arrow to select Disable .
Log Level	Use the drop-down arrow to select ERROR . NOTE: For customers concerned about the number of "qsa" error log entries generated by InSite in the Windows Event Viewer (Application log), see Appendix A .
Proxy Configuration	
If the customer site uses a proxy server to connect to the Internet, you need to complete the proxy configuration fields.	
Proxy	Enable proxy settings if you are using a proxy server.
Name	Type the name of the proxy server.
IP Address	Type the IP address of the proxy server.
Port	Enter the port number the proxy server uses.
Proxy Authentication	Enable whether the proxy server requires authentication.
Scheme	Select the appropriate authentication protocol: NTLM, Digest, or Basic.
Proxy User	Type a user name for the proxy server.
Password	Type the password for the user.

4. Click **Submit Changes**.

5. Verify the following services have started:

- **QUESTRA SERVICE AGENT**
- **TightVNC Server**

NOTE:

If the **Startup Type** of these two services is **manual**, make sure you have changed them to **automatic**.

6. Verify with GE Healthcare Service Technical Support team that the device is registered in the InSite 3.5 BackOffice and is now available for remote access.
7. If the system application or HL7 interface will be installed remotely, the local administrator should have the system Application and Support DVD set available, as directed by GE Healthcare Service, so they can be copied onto the system.

Uninstalling InSite

Run **UninstallSvcPlatform.bat** in the **Insite Installation** folder where you installed the Optima EMS system server.

If the uninstallation is successful, the following window opens and closes automatically.


```

Uninstalling BackupRestore if installed
Uninstalling BackupRestore
Uninstalling install directory
Uninstalling Utilities
Deleting C:\InSite2\bin from PATH
Deleting UpdateAgent registry Keys
Uninstalled Service Applications Success .

Uninstalling InSite?
Stopping Questra Agent
Uninstalling Apache/Tomcat Web Server.
Apache/Tomcat Web Server does not exist at installed proper . Apache/Tomcat
Web Server uninstall is skipped.
Uninstalling Questra Agent.
Removing Questra Service Agent
Dockable Device API not installed.
Done removing Questra service agent.
CKM_JAVA_HOME does NOT exist
JRE does not exist or not installed properly . JRE uninstall is skipped
Uninstalling Perl
Testing whether there is an Export dire
Uninstalled InSite2 Successfully.

```

After the uninstallation, check whether you still have a folder named **Insite2** in your C drive. If yes, delete the folder.

If you are using the Windows 7 system, right-click on **My Computer**, select **Properties > Advanced system settings > Advanced > Environment Variables**. The **Environment Variables** window opens. Delete all **INSITE** related variables.

4

InSite BackOffice

The following sections describe how to use the InSite remote desktop.

Remote Desktop Access

The following procedure describes how to access the remote desktop.

Selecting the Device Type


1. Log on to <https://insite2.health.ge.com/qss/gelogin.jsp> with your GE SSO and windows password.
2. Select device type **CS_DCR_OPTIMA_EMS_3.5**.
A device list of the Optima EMS is displayed.
3. Select a device that is polling.

Device Status


Remote Device Status	Indicator
Device is Polling	
First expected polling missed	
Second expected polling missed	

Connecting to the Device

1. Select the **Service Panel**.
2. Select **Connect to device**.


3. Initiate a connection to the device.
 - a. Under the **UL_CSD connection type** click **Connect**.

- b. Wait for the countdown timer to time-out.


Proxy Configuration

1. Click **OK** on the **Proxy Config** window.

Allowing Total Access

1. If the **Total Access: Allow Command** window opens, click **OK**.

NOTE:

If you selected the **Do not ask me about this command again** in a previous window, this box does not display subsequent times.

Opening the Remote Desktop

1. Press **Ctrl+Alt+Del**, to send the command to the device using the **CAD** button in the VNC tool bar.


If the connection was a success, the **InSite utility** icon in the system tray will be displayed in red.

To get a faster connection,

- a. Right-click on the *InSite utility* icon in the system tray
- b. Select **Connect to GE**.

2. Log on to the device with Windows authentication.

Closing the VNC Session

1. Disconnect the **Total Access Session**.

TotalAccess Session	
User	Device
212007043	CV_1212

UL_VNC		On Demand
Status	Connecting to https://ins2d.med.ge.com:9109	<input type="button" value="Disconnect"/> <input type="checkbox"/> Bypass Tunnel
Client Status	In Session	<input type="button" value="Local Settings"/>

InSite Log Level

Reducing the InSite Log Level

Use the following procedure to address customer concerns about the number of **qsa** error log entries generated by **InSite** in the **Windows Event Viewer** (Application log).


This is a two-step process:

1. Add the **FORCE** option to the list of **Log Level** choices in the **InSite 3.5 Configuration Tool**.
2. Configure the **Log Level** to **FORCE** in the **InSite 3.5 Configuration Tool**.

Adding FORCE to Log Level Choices

1. Using **Notepad**, open **C:\InSite2\Questra\AgentConfigMetaData.xml**.
2. Scroll down to the bottom of the file, to line 304.

3. Add a new line: `<Level>FORCE</Level>`

```
AgentConfigMetaData.xml - Notepad
File Edit Format View Help
<St>SOUTH DAKOTA</St>
<St>TENNESSEE</St>
<St>TEXAS</St>
<St>UTAH</St>
<St>VERMONT</St>
<St>VIRGINIA</St>
<St>WASHINGTON</St>
<St>WASHINGTON DC</St>
<St>WEST VIRGINIA</St>
<St>WISCONSIN</St>
<St>WYOMING</St>
</State>
<SvcCenter continent="NORTH AME
  <Svc>JROC</Svc>
  <Svc>RPROC</Svc>
</SvcCenter>
<SvcCenter continent="EUROPE">
  <Svc>EURO</Svc>
</SvcCenter>
<SvcCenter continent="OTHERS">
  <Svc>OTHERS</Svc>
</SvcCenter>
<LogLevels>
  <Level>ERROR</Level>
  <Level>INFO</Level>
  <Level>WARN</Level>
  <Level>TRACE</Level>
  <Level>FORCE</Level>
</LogLevels>
</MetaData>
Ln 304, Col
```

4. Save the file.
5. Close the file.

Configuring the Log Level to FORCE

1. Open the *InSite 3.5 Configuration Tool*.
Right-click the *InSite* icon in the task tray and select *Configuration Tool*.
2. Change the *Log Level* from *WARN* or *ERROR* to *FORCE*.

The screenshot shows the 'InSite ExC Configuration Tool' window. The 'Agent Configuration' section includes fields for Device Name (MC_SHR11200010C), CRM No. (SHR11200010GA), Display Name, Description, Continent (NORTH AMERICA), Country (UNITED STATES), Address lines, City (ROCHESTER), State (MINNESOTA), Postal Code, Latitude, Longitude, Institution (MAYO FOUNDATION), Department, Building, Floor, and Room. The 'Advanced Configuration' section includes Enterprise Server (PRODUCT), Service Center (JROC), Log Level (FORCE), Enterprise Server URL, Enterprise Tunnel URL, File Repository, File Watcher (Enable), Dir (C:\export), and Filter (*.*). The 'Proxy Configuration' section includes Proxy (Disable), Name (Proxy Disabled), IP Addr, Port, Proxy Authentication (Disable), Scheme (NONE), Proxy User, and Password. The 'Submit Changes' button is circled in red.

3. Click *Submit Changes*.
If you are connected remotely you are disconnected after several seconds and need to log back in to complete the next step.
4. At the *Configuration saved* prompt, click *OK*.

GE Medical Systems
Information Technologies, Inc.
8200 West Tower Avenue
Milwaukee, WI 53223 USA
Tel: +1 414 355 5000
+1 800 558 7044 (US Only)
Fax: +1 414 355 3790

Asia Headquarters

GE Medical Systems
Information Technologies, Inc.
Asia; GE (China) Co., Ltd.
1 Huatuo Road
Zhangjiang Hi-tech Park Pudong
Shanghai, People's Republic of China 201203
Tel: +86 21 3877 7888
Fax: +86 21 3877 7451

GE Medical Systems *Information Technologies, Inc.*, a General Electric Company, going to market as GE Healthcare.

www.gehealthcare.com