



# Responder<sup>®</sup> 2000

---

defibrillator/monitor

## **Service Manual**

2025653-048 Revision D  
Applies to Version 2.5X Software

## Revision History

The revision letter identifies the document's update level and changes with every update of the manual.

Part Number and Revision	Date	Comment
2025653-048A	10/2006	Initial Release
2025653-048B	10/2007	Updated LCD cable routing picture; Updated from version 2.12 software to version 2.22; Added safety information for Hellige test lamp; Added Service Sticker orderable part number
2025653-048C	02/2008	Multiple updates on pages 27-45 to update captions and images from Phillips to Torx screws. Updated contact address for CSC from Deerfield to Bothell. Update to software version V2.38.
2025653-048D	02/2012	Update images for software screens; Added printer firmware update instructions; added power cords for China and Brazil. Updated to software version 2.5X.

# TABLE OF CONTENTS

REVISION HISTORY .....	II
<b>SECTION 1: INTRODUCTION .....</b>	<b>8</b>
OVERVIEW .....	8
DESCRIPTION .....	9
CONTROLS AND INDICATORS .....	10
GENERAL SERVICE INFORMATION .....	11
Navigation .....	11
Passwords .....	11
Battery Maintenance .....	11
<b>SECTION 2: MAINTENANCE .....</b>	<b>12</b>
OVERVIEW .....	12
EQUIPMENT SETUP .....	13
Required Tools .....	13
Procedure .....	13
Upgrade Software .....	14
UPDATE SERIAL AND MODEL NUMBERS .....	16
<b>SECTION 3: TROUBLESHOOTING .....</b>	<b>18</b>
OVERVIEW .....	18
SAFETY PRECAUTIONS .....	19
WHO SHOULD PERFORM REPAIRS .....	19
SERVICE/REPLACEMENT PARTS .....	19
DISPOSAL .....	19
REPAIR TOOLS AND EQUIPMENT .....	20
TROUBLESHOOTING GUIDE .....	20
General Troubleshooting .....	20
No Boot .....	21
A. If the AC Power LED is not lit .....	21
B. If the AC Power LED is lit .....	21
Boot, but no Response .....	21
A. If the ECG trace is not moving .....	21
B. The Rotary Selection knob works properly .....	21
C. The Rotary Selection knob does not work properly .....	21
Bad Front Panel LEDs .....	22
Buttons do not Work .....	22
A. The Rotary Selection knob does not work properly .....	22
B. The Rotary Selection knob works properly .....	22
Bad Speaker .....	23
Display is Dim, Dark, Fuzzy, or Unreadable .....	23
A. If the Display is Dim .....	23
B. If the Display is Dark, Fuzzy, or Unreadable .....	23
Printing Problems .....	23
Battery Troubleshooting .....	24
SpO <sub>2</sub> Not Working .....	24
Shock into Paddle Dock Fails .....	24
Shock into Patient Fails .....	25
No ECG from Paddles/Pads/Spoons .....	25

Service Required LED is ON .....	25
<b>SECTION 4: REPAIR .....</b>	<b>26</b>
OVERVIEW .....	26
REQUIRED TOOLS .....	27
DISASSEMBLY OVERVIEW .....	27
HOW TO REPLACE SPECIFIC COMPONENTS .....	28
Capacitors .....	28
Therapy Board .....	28
SpO2 Board .....	29
ECG Board .....	29
Display/Main Board/Printer .....	30
Main Board .....	30
Printer .....	30
Power Supply and Power Control Board .....	30
Switches and Fan .....	30
ASSEMBLY .....	31
Overview .....	31
Paddle Latch Assembly .....	31
Main Board and Printer Installation .....	32
Power Control Board and Power Supply Installation .....	36
ECG Board Installation .....	40
Display Installation .....	42
SPO <sub>2</sub> Board Installation .....	44
Therapy Board Installation .....	45
Capacitor Installation .....	46
Rear Cover Installation .....	47
SECTION 5: PERFORMANCE VERIFICATION AND SAFETY TESTING .....	50
OVERVIEW .....	50
REQUIRED EQUIPMENT .....	51
PREVENTIVE MAINTENANCE .....	51
ANNUAL INSPECTION .....	51
VERIFY OPERATION .....	52
Checkout Procedures .....	52
MAINTENANCE MENU REFERENCE .....	53
Service Tests .....	54
Button Test .....	54
Printer Speed Test .....	55
Battery Readings .....	55
Display Test .....	56
ECG Lead Test .....	56
Temperature Readings .....	57
Fan Test .....	57
Device Configuration .....	58
Select Language .....	58
Stored Data Management .....	58
Clear Event Logs .....	58
Clear Service Required Indicator .....	58
Clear Settings Menu Password .....	58
Clear All Stored Data .....	58
Send Event Logs to PC .....	58

Send/Receive Stored Parameters .....	59
VISUAL INSPECTION.....	60
OPERATIONAL TESTS .....	60
Initial Power Up .....	60
Maintenance Mode Tests .....	61
Verify Serial Communications .....	61
Required Equipment .....	61
Procedure .....	61
Defibrillator Tests.....	62
Energy Timeout Test.....	62
Defibrillator Test Chart.....	62
Paddles Test.....	64
Pads Test.....	64
Spoons Test.....	64
Pacing Test.....	65
ECG Tests .....	67
Required Equipment .....	67
ECG Connections Test.....	67
Lead Detection Test.....	68
ECG Lead Test.....	69
Paddles ECG Test.....	70
Thermal Printer Operation.....	70
Cardioversion Test .....	70
SpO <sub>2</sub> Operation .....	70
SAFETY TESTS.....	71
Hipot Test .....	71
Required Equipment .....	71
ECG to Serial 4.0 KVAC Test .....	72
AC Line/Neutral to AC Ground 1.5 KVAC Test.....	72
ECG to SpO <sub>2</sub> 1.5 KVAC Test.....	72
Apex Test Load to Apex and Sternum 3.0 KVDC Test .....	73
ECG to Apex and Sternum 3.0 KVDC Test.....	73
Apex and Sternum to Service Port 3.0 KVAC Test .....	73
Leakage Current Test.....	74
Required Equipment .....	74
Equipment Setup .....	74
Test Procedure .....	77
ALTERNATE SAFETY TEST PROCEDURES.....	77
Test Equipment .....	77
Leakage Current “Measurement Device” (MD) Characteristics .....	78
Tests .....	78
Visual Inspection.....	78
Leakage Current .....	79
Patient Leakage Current to Ground .....	82
Patient Leakage Current, Mains on Applied Part (All SIP/SOPs Grounded) .....	85
Insulation Resistance.....	87
<b>SECTION 6: PARTS AND ACCESSORIES.....</b>	<b>88</b>
OVERVIEW .....	88
KITS .....	89
ACCESSORIES.....	95
POWER CORDS.....	97

<b>SECTION 7: THEORY OF OPERATION .....</b>	<b>99</b>
OVERVIEW .....	99
SYSTEM OVERVIEW .....	100
Functions .....	100
Features .....	100
SYSTEM INTERCONNECTION BLOCK DIAGRAM .....	101
COMPONENT DESCRIPTIONS .....	102
Main Board .....	102
Therapy Board .....	102
ECG Board .....	102
Power Control Board .....	102
SpO <sub>2</sub> Board .....	102
<b>SECTION 8: SPECIFICATIONS AND SAFETY .....</b>	<b>103</b>
OVERVIEW .....	103
SPECIFICATIONS .....	104
PHYSICAL DIMENSIONS .....	106
ENVIRONMENTAL REQUIREMENTS .....	106
Operating Conditions .....	106
Storage and Shipping Conditions .....	106
RHYTHMx <sup>®</sup> ECG ANALYSIS ALGORITHM .....	107
Detection Rate .....	107
Fine VF .....	107
Asystole .....	107
Noise Detection .....	107
Non-Committed Shock .....	107
Sync Mode .....	108
SVT (Supraventricular Tachycardia) Discriminators .....	108
SVT Rate .....	108
Continuous Monitoring For Shockable Rhythm .....	108
Pacemaker Pulse Information .....	108
STAR <sup>®</sup> BIPHASIC DEFIBRILLATION WAVEFORM .....	109
Energy Levels and Patient Impedance .....	110
ELECTROMAGNETIC COMPATIBILITY REQUIREMENTS .....	111
Emissions .....	111
Immunity .....	111
ENVIRONMENTAL STANDARDS .....	112
Shock and Vibration .....	112
Storage and Shipping .....	112
ELECTROMAGNETIC EMISSIONS TABLE .....	113
ELECTROMAGNETIC IMMUNITY TABLE .....	114
RF COMMUNICATIONS TABLE .....	117
CONNECTORS .....	118
Case .....	118
Power Connector – IEC 320 type – AC Power .....	118
Paddles/Pads Connector – GE Defibrillator Paddles .....	118
Paddle Connector ID Codes .....	119
ECG Connector – GE ECG Connector – Five lead ECG .....	119
ECG Cable Identification Encoding .....	120

SpO <sub>2</sub> Connector – GE SpO <sub>2</sub> .....	121
Serial Service Connector – DB9-Female .....	121
Paddle Cradles (Test Load) .....	121
MAIN CPU PCBA .....	122
J211 - Power Control Interface to Power Control PCBA .....	122
J214 – LCD Panel .....	123
J213 – LCD CCFL Backlight .....	124
J218 – UI Interface to Power Control PCBA .....	125
J216 – SPEAKER .....	127
J219 – PRINTER HEAD .....	127
J220 – PRINTER MOTOR .....	128
J221 – SpO <sub>2</sub> INTERFACE .....	128
J210 – THERAPY PCBA .....	129
J212 – ECG PCBA .....	129
THERAPY PCBA .....	130
J103 - Power Supply .....	130
J110 – Main CPU PCBA interface .....	130
J109 – Paddles Control Connector .....	131
J101, J102, J108 – Paddles HV Connectors .....	131
J104, J105, J106, J107 – Energy Storage Capacitor Connectors .....	132
J126, J127 – Internal Test Load Resistor .....	132
ECG PCBA .....	132
J412 – Interface to main CPU PCBA .....	132
J423 – Patient ECG Connection .....	133
POWER CONTROL PCBA .....	133
J315 – Power input from DC power supply .....	133
J325 – Smart Battery Connector – Data and Power .....	133
J303 – Therapy power .....	134
J311 – Power Control interface from main CPU PCBA .....	134
J318 – UI Features - Interface to main CPU PCBA .....	135
J317 – RS-232 SERIAL SERVICE INTERFACE .....	136
J323 – FAN .....	137
J324 – ENCODER “Trim Knob” .....	137
J2 - SpO <sub>2</sub> HOST INTERFACE .....	137
SpO <sub>2</sub> PATIENT INTERFACE .....	138
POWER SUPPLY MODULE .....	138
TB1 – POWER INPUT .....	138
POWER OUTPUT .....	138
BATTERY .....	139
J325 – Smart Battery Connector – Data and Power .....	139
CONTACT INFORMATION/CUSTOMER SERVICE .....	140

## SECTION 1: Introduction

### Overview

This Service manual provides information needed to service the Responder™ 2000. This manual should only be used by technical personnel trained to service the Responder 2000.

This chapter contains general information for servicing the Responder 2000.

This Service manual assumes familiarity with the controls and basic operation of the Responder 2000. Detailed information regarding controls, operation, set-up, and regular maintenance procedures are found in the Operator's manual. If necessary, review the Operator's Manual before servicing the Responder 2000.




TOPIC	PAGE
Description	9
Controls and Indicators	10
General Service Information	11



## Description

The Responder 2000 is a defibrillator/monitor/pacemaker intended for use by personnel trained in its operation. The device is lightweight, portable, easy to use and reliable. It incorporates a 320 x 240 transmissive color TFT color display for wide viewing angles in all light conditions. The device operates using either an AC power supply or internal rechargeable Li-Ion battery. The device provides continuous ECG monitoring and three types of therapies: defibrillation, cardioversion and external pacing. Defibrillation can be applied manually or semi-automatically. Pacing therapy can be either fixed or demand. The device employs patented RHYTHMx<sup>®</sup> software which provides ECG rhythm analysis. STAR<sup>®</sup> Biphasic waveform delivers impedance-compensated energy ranging from 2-270 Joules. Features and options include external paddles, spoons, disposable pads, 3- and 5-lead ECG, pulse oximetry (SpO<sub>2</sub>), built-in 60 mm thermal printer, internal storage of event history and remote synchronization to bedside monitor.

## Controls and Indicators



## General Service Information

Refer to the Operator's Manual for setup (and configuration options) required before placing the Responder 2000 into service.

### Navigation

Use the Rotary Selection knob to navigate the Responder 2000 user interface. The Rotary Selection knob is used for:

- Scrolling through menus and sub-menus
- Selecting soft keys
- Setting values

### Passwords

The default service password is "HEART1". See the note under Maintenance Mode Tests in Section 5 (on page 54) for detailed password instructions.

### Battery Maintenance

The Responder 2000 has a rechargeable battery requiring periodic calibration. Calibration consists of a full charge, full discharge, and full charge of the battery. The cycle may take up to 20 hours to complete.

To calibrate the battery using the optional external charger:

1. Insert the battery into the external charger.
2. Press the calibration button on the battery. The Mode light turns red to indicate calibration in progress.

To stop calibration, press the calibration button again. The Mode light turns green and the battery begins to charge.

To calibrate the battery without the external charger:

1. Insert the battery into the Responder 2000.
2. Plug the Responder 2000 into an AC outlet.
3. Allow the battery to fully charge.
4. Disconnect the AC power cord and leave on the Responder 2000 until the battery is fully discharged (full discharge takes about four hours).
5. Reconnect the AC power cord and allow the battery to fully charge.

## SECTION 2: Maintenance

### Overview

Use the following procedures to upgrade the Responder 2000 software:

- Upgrade Software to update one or more software files
- Upgrade Printer Software to upgrade printer software and fonts

<b>TOPIC</b>	<b>PAGE</b>
Equipment Setup	13
Upgrade Software	14
Upgrade printer software	16
Update Serial and Model Numbers	17

## Equipment Setup

### Required Tools


- Serial Cable
- PC (Windows 2000 or XP)

### Procedure


1. Plug AC power cord.

**Note:** When performing a software upgrade, the Responder 2000 must be running on AC power.

2. Connect serial cable between the communications port of the PC and the serial port of the Responder 2000.
3. At the desktop on the PC, double-click the Responder 2000 Upgrade icon to open the CodeLink 2000 program.
4. Select the appropriate Com port and click **OK**.




5. Complete the following on-screen instructions and click **OK**.



## Upgrade Software


1. Perform the **Equipment Setup** procedure starting on page 13 until the CodeLink 2000 selection window is displayed:  
**Note:** Close all other application on your PC before you start with the upgrade.



2. Click **Choose folder** to select all the update files in a folder or **Choose file** to select a single update file.
3. Browse to the folder (or file) to be uploaded. Select the folder that contains the files listed in the screen shot below.  
**Note:** Software numbers and revisions may differ slightly from those shown in the screen capture below.
4. Click **OK** (or **Open**).

While the files are updating, the status of each file is displayed.

As the files are updated, the Shock button and then the Manual button flash on the Responder 2000. When all files are updated, a confirmation screen is displayed:



5. Verify Pass is displayed for each file and click **OK**.

**Note:** If a file fails to update, retry updating the file individually.

The files must be updated without reboot. If the Responder 2000 is rebooted during the installation, the printer may not work (the Main Board firmware must then be reinstalled at Cardiac Science or the Main Board has to be replaced by a new Main Board).

#### Process of upgrading files individually

The Upgrade Complete Notification is displayed:



6. Click **Yes** to upgrade more files (and repeat this procedure from step 2) or **No** to exit.
7. Click **Exit**.
8. Verify the file upgrade on the Responder 2000.  
Compare the updated software version numbers on the Responder 2000 to those in the upgrade folder.
  - a. Press and hold the **Power** button for 5 seconds to reboot the Responder 2000.
  - b. Highlight and click the **System** Menu.
  - c. Click **About**.
  - d. Click **Next** to scroll through the version screens. Verify all software versions are updated and match the revisions indicated by the files in the upgrade folder.

## Upgrade Printer Software


When installing the printer software, the files must be upgraded in the following order without reboot:

- printerflash.bin
- printerfont.bin
- printerflash.bin again

### Caution: Procedure Failure.

The printer files must be updated without reboot. If the unit is rebooted during the installation, the printer may not work (the Main Board firmware must then be reinstalled at a repair center).


1. Perform the Equipment Setup procedure starting on page 13 until the CodeLink 2000 selection window is displayed:



2. Click **Upgrade Printer** to select the update folder.
3. Navigate to the **Printer Files** folder. It should contain the files shown in the screen shot of step 5.
4. Click **Open**. While the files are updating, the status of each file is displayed.

As the files are updated, the **Shock** button and then the **Manual** button flash on the Responder 2000.

5. After the files are updated, a confirmation screen is displayed:



**Note:** If a file fails to update, retry updating the file individually in the following order without rebooting:

- printerflash.bin
- printerfont.bin
- printerflash.bin again

The files must be updated without reboot. If the Responder 2000 is rebooted during the installation, the printer may not work (the Main Board firmware must then be reinstalled at Cardiac Science or the Main Board has to be replaced by a new Main Board).



6. **Pass** is displayed for the file and click **OK**.

**Note:** If the file fails to update, try updating the file again.


The Upgrade Complete Notification is displayed:

7. Click **Exit**.
8. Verify the file upgrade on the Responder 2000.
  - a. Press and hold the Power button for 5 seconds to reboot the Responder 2000.
  - b. Highlight and click the System Menu.
  - c. Click About.
  - d. Click Next to scroll through to the printer version screen. Verify the Printer and Font software versions are updated and match the file versions in the printer upgrade folder.

## Update Serial and Model Numbers

When replacing the Main Board, the Serial and Model numbers must be updated.

1. Before replacing the **Main Board**, record the Serial and Model numbers from the **About** screen.  
Note: If the Responder 2000 cannot be booted, record the model and serial number from the back panel label.
2. Replace the **Main Board**.
3. Perform the **Equipment Setup** procedure starting on page 13 until the CodeLink 2000 selection window is displayed:



4. Click **Factory Data** and enter the previously model and serial numbers.  
**Note:** If the **Main Board** is not replaced, the fields are not editable.
5. Click **OK** to accept.

## SECTION 3: Troubleshooting

### Overview

This section describes how to troubleshoot the Responder 2000. These instructions are intended for use only by service providers who are specifically trained to service the Responder 2000.

<b>TOPIC</b>	<b>PAGE</b>
Safety Precautions	19
Who Should Perform Repairs	19
Service/Replacement Parts	19
Disposal	19
Repair Tools and Equipment	20
Troubleshooting Guide	20

## Safety Precautions

### **WARNING: Shock Hazard**

The Responder 2000 is designed to deliver high-voltage therapeutic shock. Before performing any service on the equipment, read and follow all safety precautions and instructions in the Operator's Manual.

### **WARNING: Shock Hazard or Equipment Damage**

Before servicing the Responder 2000, disconnect the AC power cord and remove the battery.

### **WARNING: Shock Hazard or Equipment Damage**

Internal components of the Responder 2000 may still contain high voltages even after the AC power cord and battery are removed. Before working on any internal component, verify high voltages are not present.

### **WARNING: Shock Hazard or Equipment Damage**

Some service activities require the Responder 2000 to be energized with covers removed. Ensure all personnel and equipment is clear while the equipment is energized.

### **WARNING: Biological Contamination**

During normal operation, the Responder 2000 may be contaminated by blood, body fluids, or other biological agents. Always assume the Responder 2000 is contaminated and use appropriate safety procedures until decontamination is performed. Always decontaminate the Responder 2000 in accordance with hospital or facility procedures before servicing or returning to service. Refer to the Responder 2000 Operator's Manual for recommended cleaning agents and instructions

### **Caution: Electrostatic Damage**

Always use a wrist grounding strap and anti-static mat while performing service on any internal components.

### **Caution: Voiding Product Warranty**

Any service performed on the Responder 2000 must be provided by authorized service representatives only. Unauthorized repair voids the product warranty.

## Who Should Perform Repairs

Repair and service of the Responder 2000 must be performed by qualified service technicians trained in safe and proper servicing of the Responder 2000.

## Service/Replacement Parts

For service, please contact your local GE agency. For additional information, please visit our Web site at: <http://www.gehealthcare.com>.

Please have the serial and model numbers available when contacting Customer Service. (The serial and model numbers are located on the back (and on the inside) of the Responder 2000.)

## Disposal

Always dispose of the Responder 2000, any unserviceable parts or accessories, or batteries in accordance with any local disposal regulations for equipment containing electronic parts. Note the following precautions:

### **WARNING: Shock Hazard**

Disposal of the Responder 2000 with the battery inserted presents a potential shock hazard.

### **WARNING: Environmental Contamination**

Disinfect the Responder 2000 appropriately prior to disposal. Also recycle or dispose of the lithium-ion battery in accordance with applicable local regulations.

### **WARNING: Fire or explosion hazard**

Do not burn or incinerate the battery. Recycle or dispose of the lithium-ion battery in accordance with applicable local regulations.

### **Caution: Environmental Contamination**

Dispose of the pads or electrodes in accordance with any local disposal regulations.

## Repair Tools and Equipment

The following tools are needed to perform the procedures in this section.

- 1.5mm Hex Driver
- T8 Torx Driver
- TR15 Security Torx
- TR25 Security Torx
- Serial cable for software upgrades (with Windows 2000 or XP PC)

## Troubleshooting Guide

### General Troubleshooting

Before performing any troubleshooting, check the following:

Ensure all external cables and connections are tight and undamaged.

- If possible, ensure the unit is plugged in to a receptacle with appropriate voltage available (see SECTION 8: Specifications and Safety for voltage requirements). If no receptacle is available, use a known good, charged battery to perform any tests.
- Visually inspect the unit for any obvious external damage, including cracks in the display, case, or connectors.
- Ensure all connector pins and sockets are clean, free of debris, and intact.
- Discuss the issues with the operator:
  - Have the operator explain and demonstrate the problem.
  - Ask about any previous repairs or problems.
  - Has the unit been stressed (including extreme heat/cold, submersion, falls, etc.)?

## No Boot

Use this procedure to troubleshoot the Responder 2000 when plugged into an AC outlet. For battery troubleshooting, including an operation failure while using the battery, see "Battery Troubleshooting".

### A. If the AC Power LED is not lit

1. Verify power at the AC outlet.
2. Verify AC power available from the AC power cord.
3. Replace power supply.
4. Replace the **Power Control Board**.

### B. If the AC Power LED is lit

1. Check **Power** button for mechanical operation.
2. Check cables from **Power Control Board** to **Main Board**.
3. Ensure cables are properly seated, tight, and undamaged.
4. Replace **Power Control Board**.
5. Replace **Main Board**.

## Boot, but no Response

Use this procedure to troubleshoot the Responder 2000 when the unit seems to boot up normally (i.e., the Front Panel LEDs flash and a speaker tone is heard when the Power button is pressed), but does not respond to any subsequent user input.

### A. If the ECG trace is not moving

1. Replace the **Main Board**.

### B. The Rotary Selection knob works properly

1. Check buttons for mechanical operation.
2. Check cables from **Power Control Board** to **Main Board**.  
Ensure cables are properly seated, tight, and undamaged.
3. Replace **Power Control Board**.
4. Replace **Main Board**.

### C. The Rotary Selection knob does not work properly

1. Replace **Rotary Selection knob**.

**Note:** The figure below shows the area where the Rotary Selection knob cable can be easily damaged if flexed several times or otherwise stressed.



2. Check cables from **Power Control Board** to **Main Board**.  
Ensure cables are properly seated, tight, and undamaged.
3. Replace **Power Control Board**.
4. Replace **Main Board**.

## Bad Front Panel LEDs

Use this procedure to troubleshoot the Responder 2000 when the **Front Panel** LEDs are not displaying correctly. Normal operation is indicated by:

At system boot, all LEDs will flash.

Green AC Power LED is lit when the unit is plugged in.

Yellow Battery Charging LED is lit when the battery is lit when the unit is plugged in and a discharged battery is inserted.

Red Service Required LED should be off.

Blue Manual button LED is lit.


1. Ensure the unit is plugged in and receiving power from the wall outlet.
2. If the battery charge LED is not lit when a battery is inserted, replace the battery with a known good battery.
3. Check cables from **Power Control Board** to **Main Board**.
4. Ensure cables are properly seated, tight, and undamaged.
5. Replace **Power Control Board**.
6. Replace **Main Board**.

## Buttons do not Work

### A. The Rotary Selection knob does not work properly

1. Replace Rotary Selection knob.

**Note:** The figure below shows the area where the Rotary Selection knob cable can be easily damaged if flexed several times or otherwise stressed.



2. Check cables from **Power Control Board** to **Main Board**.
3. Ensure cables are properly seated, tight, and undamaged.
4. Replace **Power Control Board**.
5. Replace **Main Board**.

### B. The Rotary Selection knob works properly

Perform Button Tests (refer to the Button Test procedure on page 54).

If *Manual, Charge, Shock* fails:

1. Check cables from **Power Control Board** to **Main Board**.
2. Ensure cables are properly seated, tight, and undamaged.
3. Replace **Power Control Board**.
4. Replace **Main Board**.

If *Paddles, Charge, Shock* fails:

1. Replace paddles.
2. Check cable from rear connector to **Therapy Board**.
3. Check cable between **Therapy Board** and **Main Board**.
4. Replace **Therapy Board**.
5. Replace **Main Board**.

## Bad Speaker

Use this procedure to troubleshoot the Responder 2000 when the **Speaker** is not working correctly.

**Note:** The Speaker cannot be removed from the Front Panel. If the speaker must be replaced, the entire Front Panel must be replaced.

Normal operation is indicated by:

At system boot, the speaker sounds a brief tone (listen carefully because the tone is easy to miss in a noisy environment). When audio indication is enabled in the *Settings* menu, the speaker should emit loud, non-distorted tones.

1. Check speaker plug and cable.
2. Check speaker resistance. It should be between 4 and 10 ohms.
3. Replace **Front Panel**.
4. Replace **Main Board**.

## Display is Dim, Dark, Fuzzy, or Unreadable

Normal operation is indicated by:

At system boot, all LEDs will flash.

The speaker sounds a brief tone (listen carefully because the tone is easy to miss in a noisy environment).

Green **AC Power** LED is lit when the unit is plugged in.

Yellow **Battery Charging** LED is lit when the unit is plugged in and a partially or fully discharged battery is inserted.

Red **Service Required** LED should be off.

Blue **Manual** button LED is lit.

### A. If the Display is Dim

1. Check the software setting for display brightness.
2. Check the backlight cable
3. Run display test from the Maintenance menu to verify the display is working and there are no bad pixels.

If these checks are ok, continue with the following steps.

### B. If the Display is Dark, Fuzzy, or Unreadable

1. Check the display cable.
2. Replace the **Display**.
3. Replace the **Main Board**.

## Printing Problems

**Note:** Check the paper roll to ensure the customer is using approved paper listed in the Accessories list. Non-approved or generic paper may cause printer damage or failure.

1. Ensure paper roll is in good condition (dry and undamaged).
2. Ensure paper is loaded correctly
3. If the *Print* icon is not selectable, ensure the printer door is fully closed and the paper is properly inserted.
4. Replace the printer door.
5. Check printer cables.
6. Replace printer.
7. Replace **Main Board**.

After troubleshooting, always print a test strip (page 55) to ensure proper operation.

## Battery Troubleshooting

Check battery statistics. Good batteries have a capacity of at least 4500 mAH and will be charging when the Responder 2000 is plugged into an AC outlet.

1. Ensure battery is fully charged.
2. Perform Battery Calibration on page 11.
3. Replace the **Battery**.
4. Replace **Power Control Board**.
5. Replace **Main Board**.

## SpO<sub>2</sub> Not Working

Verify the Responder 2000 has a SpO<sub>2</sub> **Main Board** installed.

Check the *About* screen. "SpO<sub>2</sub> - - -" indicates a non-SpO<sub>2</sub> **Main Board** is installed in the unit.

1. Connect the Responder 2000 to a SpO<sub>2</sub> simulator with the appropriate sensor.  
If the Responder 2000 SpO<sub>2</sub> display does not agree with the simulator setting, substitute a known good SpO<sub>2</sub> sensor.
2. Check SpO<sub>2</sub> cable.
3. Check cable from **Front Panel** to **SpO<sub>2</sub> Board**.
4. Replace **SpO<sub>2</sub> Board**.
5. Replace **Main Board**.

## Shock into Paddle Dock Fails

If shocking into the paddle dock does not work or shows high impedance.

Perform a button test for each Paddle.

- If the test fails, replace the Paddle.
- If the test does not fail, connect the Paddles to a simulator and attempt to shock.

If the shock is successful:

1. Check the cables from the **Therapy Board** to the paddle dock (one on each side).
2. Replace the **Therapy Board**.

If the shock is unsuccessful:

1. Ensure the paddles are securely paced in the paddle dock.
2. Replace the Paddles.
3. Check the cables between the **Rear Panel** and the **Therapy Board**.
4. Check the cable between the **Therapy Board** and the **Main Board**.
5. Replace the **Therapy Board**.
6. Replace the **Main Board**.



**Shock into Patient Fails**

1. Perform a button test for each Paddle.  
If the test fails, replace the Paddle.  
If the test does not fail, connect the Paddles/Pads/Spoons to a simulator and attempt to shock.
2. If the shock is successful:  
The unit is working correctly. Check patient preparation.
3. If the shock is unsuccessful:
  1. Replace the Paddles/Pads/Spoons.
  2. Check the cables between the **Rear Panel** and the **Therapy Board**.
  3. Check the cable between the **Therapy Board** and the **Main Board**.
  4. Replace the **Therapy Board**.
  5. Replace the **Main Board**.

**No ECG from Paddles/Pads/Spoons**

1. Replace the paddles/pads/spoons.
2. Check the cables between the **Rear Panel** and the **Therapy Board**.
3. Check the cable between the **Therapy Board** and the **Main Board**.
4. Replace the **Therapy Board**.
5. Replace the **Main Board**.

**Service Required LED is ON**

If the Service Required LED is on and the fault has been corrected, perform this procedure to clear the fault.

1. Go to *Maintenance > Stored Data Management > Clear Service Required Indicator*.
2. Press *Accept*.
3. Cycle power.

If the fault does not clear, go to *History > Event Log* to find the fault log entry.

**Note:** All service required Error codes are displayed in the History with a red exclamation mark.

## SECTION 4: Repair

### Overview

This section describes how to assemble and disassemble the Responder 2000. These instructions are intended for use only by service providers who are specifically trained to service the Responder 2000.

This section is divided into two parts:

How to Replace Specific Components: Provides high-level disassembly steps with references to the assembly steps.

Assembly: Lists production assembly steps applicable to field repair.

<b>TOPIC</b>	<b>PAGE</b>
Required Tools	27
Disassembly Overview	27
How to Replace Specific Components	28
Assembly	31

## Required Tools

The following tools are necessary for Assembly/ Disassembly the Responder 2000.

**Note:** The Security Torx is also known as Tamper Resistant Torx.

T8 Torx driver  
T15 Security Torx driver  
T25 Security Torx driver  
1.5mm hex driver  
Philips Screwdriver  
Tube of Silicone

## Disassembly Overview

**WARNING: Lethal Shock Hazard.**

In the event of equipment failure, the two main capacitors may retain dangerous voltages even if the Responder 2000 is disconnected from AC power and the battery is removed. Normally, the capacitors are discharged when power is shut off; however, it is possible for equipment damage to prevent the capacitors from discharging properly.

Always assume the capacitors are fully charged.

**Caution: Shock Hazard or Equipment Damage.**

Before opening the case, ensure the AC power cord is disconnected and the battery is removed.

**Caution: Equipment Damage.**

Always wear a grounding wrist strap and use an anti-static mat when handling parts.

**Caution: Procedure Failure.**

Even though some assembly steps may not be applicable to a specific replacement procedure, all other steps must be performed in the order listed.

The Responder 2000 has two access points:

**Rear Cover:** Provides access to **Capacitors**, **Therapy Board**, optional **SpO<sub>2</sub> Board**, and **ECG Board**, and **Front Bezel**.

Removing the Bezel allows access to the **Main Board**, **Printer**, and **Display**.

**Base Cover:** Provides access to the **Power Supply**, **Power Control Board**, **Fan**, and **Switches**.

**Note:** The Speaker is part of the Front Body Assembly.

## How to Replace Specific Components

Use the following procedures as a guide to replace specific internal components in the Responder 2000. Not all components are listed. (For example, if the **Front Body** is damaged and must be replaced, the entire unit must be disassembled and then reassembled using the more detailed Assembly instructions).

### Capacitors

(see page 46)

#### **WARNING: Lethal Shock Hazard.**

In the event of equipment failure, the two main capacitors may retain dangerous voltages even if the Responder 2000 is disconnected from AC power and the battery is removed. Normally, the capacitors are discharged when power is shut off; however, it is possible for equipment damage to prevent the capacitors from discharging properly.

Always assume the capacitors are fully charged.

1. Remove **Rear Cover** (see *Figure 67: Rear Body Installed*).

**Note:** Two cover screws are located under the Bed Hook Covers.

2. Cut tie wraps to free capacitor leads.
3. Disconnect **Capacitor** leads (see *Figure 58: Capacitors Installed*) from **Therapy Board**.

#### **Caution: Equipment Damage.**

Carefully note the capacitor lead connection points and double check the lead and jack labels. Connecting the **Capacitors** to the wrong jack could cause equipment damage.

4. Reconnect **Capacitors**. Ensure lead connections are tight.
5. Tie wrap leads and cables. Refer to the assembly instructions for proper cable routing.
6. If necessary replace or reseal rubber tubing along the edge of the **Front Body**.
7. Replace the **Rear Cover**. Ensure the cover seats properly—check for gaps or bulges around the entire edge.

**Note:** One screw is located under the label of the Rear Cover. (See *Figure 67: Rear Body Installed*)

8. Replace **Bed Hook Covers**. Check for proper orientation (see *Figure 68: Bed Hook Cover Improperly Installed*).

### Therapy Board

(see page 45)

1. Follow the **Capacitor** instructions (above) to remove the **Capacitors**.
2. Disconnect all other leads from the **Therapy Board** (see *Figure 56: Therapy Board Cables Installed*).
3. Remove the screws holding the **Therapy Board** (see *Figure 54: Therapy Board Installed*).
4. Exchange the board and replace screws.
5. Reconnect all cables and leads.

#### **Caution: Equipment Damage.**

Carefully note the capacitor lead connection points and double check the lead and jack labels. Connecting the **Capacitors** to the wrong jack could cause equipment damage.

6. Follow the **Capacitor** instructions (above) to complete the reassembly.

## SpO<sub>2</sub> Board

(see page 44)

**Note:** Most SpO<sub>2</sub> failures are the result of bad sensors. Before replacing the **SpO<sub>2</sub> Board**, be sure to complete the troubleshooting procedure on page 24 to verify the board is the cause of the fault.

1. Follow the **Capacitor** and **Therapy Board** instructions (above) to remove the **Capacitors** and **Therapy Board**.
2. Remove the screws holding the **SpO<sub>2</sub> Board** (see *Figure 51: SpO<sub>2</sub> Board Installed*).
3. Lift the **SpO<sub>2</sub> Board** out to disconnect it from the **Main Board**.
4. Disconnect the SpO<sub>2</sub> cable.
5. Exchange the board and reconnect the SpO<sub>2</sub> cable.
6. Reconnect the **SpO<sub>2</sub> Board** to the **Main Board**.
7. Replace the screws.
8. Follow the **Therapy Board** and **Capacitor** instructions (above) to replace the **Therapy Board** and **Capacitors**.

## ECG Board

(see page 40)

**Note:** Most ECG failures are the result of bad cables or sensors. Before replacing the **ECG Board**, be sure to complete the troubleshooting procedure on page 25 to verify the board is the cause of the fault.

1. Follow the **Capacitor** and **Therapy Board** instructions (above) to remove the **Capacitors** and **Therapy Board**.
2. Remove the **ECG Board** EMI shield.
3. Disconnect the ECG cable (see *Figure 41: ECG Cable Installation Detail*).
4. Remove the screws holding the **ECG Board**.
5. Exchange the board and replace the screws.
6. Reconnect the ECG cable.
7. Replace the **ECG Board** EMI shield.
8. Follow the **Therapy Board** and **Capacitor** instructions (above) to replace the **Therapy Board** and **Capacitors**.

## Display/Main Board/Printer

(see page 42)

1. Follow the **Capacitor**, **Therapy Board**, **SpO<sub>2</sub> Board**, and **ECG Board** instructions (above) to remove the **Capacitors Therapy Board**, **SpO<sub>2</sub> Board**, and **ECG Board**.
2. Remove the screws holding the **Front Bezel** (see *Figure 45: Front Bezel Installed*).
3. Remove the screws holding the **Display** (see *Figure 43: LCD Screw Locations*).

At this point the **Display** can be replaced (disconnect **Display** cables).

If the **Main Board** or **Printer** must be replaced, perform the following steps:

### Main Board

- a. Remove the screws (see *Figure 11: Main Board Installed*).
- b. Disconnect all cables
- c. Exchange the Main Board.
- d. Reconnect all cables.
- e. Replace the screws.

### Printer

- a. Remove the screws (see *Figure 14: Printer Screws*).
- b. Disconnect the ribbon and Power cables.
- c. Exchange the **Printer**.
- d. Reconnect ribbon and Power cables.

**Note:** Removing the printer bracket may help when connecting the Printer ribbon cable to the Main Board (see *Figure 10: Printer Bracket Removal*).

- e. Replace the screws.
4. Connect the **Display** cables.
5. Replace the **Display** screws.
6. Follow the **ECG Board**, **SpO<sub>2</sub> Board**, **Therapy Board**, and **Capacitor** instructions (above) to replace the **ECG Board**, **SpO<sub>2</sub> Board**, **Therapy Board**, and **Capacitor**.

## Power Supply and Power Control Board

(see page 36)

Remove the **Base Cover**, disconnect all cables, and remove screws.

**Note:** The Power Supply is held in place only by the two larger Base Cover screws.

## Switches and Fan

(see page 36)

Remove **Base Cover** and **Power Control Board**.

## Assembly

### Overview

The following section details the assembly steps starting from a fully disassembled Responder 2000. When replacing a single board or other component, not all reassembly steps may be applicable.

**WARNING: Lethal Shock Hazard.**

In the event of equipment failure, the two main capacitors may retain dangerous voltages even if the Responder 2000 is disconnected from AC power and the battery is removed. Normally, the capacitors are discharged when power is shut off; however, it is possible for equipment damage to prevent the capacitors from discharging properly.

Always assume the capacitors are fully charged.





**Caution: Shock Hazard or Equipment Damage.**

Before opening the case, ensure the AC power cord is disconnected and the battery is removed.

**Caution: Procedure Failure.**

Even though some assembly steps may not be applicable to a specific replacement procedure, all other steps must be performed in the order listed.

### Paddle Latch Assembly

Assembly Step	Details			
<p>Install <b>Paddle Latch</b>, three springs, and <b>Paddle Latch Retainer</b> (using two Phillips screws) on each side of the <b>Front Body</b>.</p> <p><b>Note:</b> Ensure the <b>Paddle Latches</b> are oriented correctly and the Paddle Latch notch fits into the grove on the Front Body.</p> <p>Apply Silicone sealant to the lower half of the joint between the <b>Front Body</b> and the <b>Paddle Latch Retainer</b>.</p> <p><b>Note:</b> The seal must be complete (no gaps).</p>				
	<p><b>Figure 1: Front Body Assembly (Front View)</b></p>	<p><b>Figure 2: Orientation Detail</b></p>	<p><b>Figure 3: Assembled Latch</b></p>	<p><b>Figure 4: Silicone Application</b></p>

## Main Board and Printer Installation

Assembly Step	Details
---------------	---------

Install Ribbon Cables on **Main Board** at **J210** and **J218**.

**Note:** Ensure the colored strip on cable into **J218** is visible when board is oriented as shown.

**Note:** Ensure cables are installed correctly as shown below.




Figure 5: Main Board

## Cable Installation Notes

**Figure 6** shows the cable correctly inserted into the jack.

**Figure 7** shows the jack latch not fully engaged (latches on both sides must be fully depressed).

**Figure 8** shows the cable not fully inserted into the jack.




Figure 6: Good: Proper Installation



Figure 7: Bad: Jack Latch not Fully Engaged




Figure 8: Bad: Cable not Fully Inserted

**Note:** the following steps detail replacing the entire printer. If the only the printer door is replaced, see Figure 46 for printer door installation.

Remove the **Printer Roller** and connect the printer cables to **J219** and **J220**.

**Note:** Ensure the cable and cable plug are fully inserted into the cable jacks.

**Note:** Removing the printer bracket may help when connecting the Printer ribbon cable to the **Main Board**.




Figure 9: Printer Connections









Figure 10: Printer Bracket Removal



Assembly Step	Details
<p>Place the <b>Main Board</b> and <b>Printer</b> into the <b>Front Body</b>.</p> <p>Tuck the ribbon cables connected to <b>J210</b> and <b>J218</b> as shown.</p> <p>Secure the <b>Main Board</b> with seven Torx screws and install the copper <b>Contact Spring</b> as shown.</p> <p><b>Note:</b> Ensure the printer ribbon cable does not twist during installation.</p>	 <p>Figure 11: Main Board Installed</p>  <p>Figure 12: Cable Detail</p>  <p>Figure 13: Copper Contact Spring Installation Detail</p>
<p>Turn the assembly over (support the <b>Printer</b> so the cable does not twist) and secure the printer with three flat head hex socket cap screws.</p> <p><b>Note:</b> Tighten the center screw first to properly align the <b>Printer</b>.</p>	 <p>Figure 14: Printer Screws</p> <p>Tighten this screw first</p>
<p>Connect the speaker cable to <b>Main Board J216</b>.</p>	 <p>Figure 15: Speaker Connection</p>

## Assembly Step

## Details

Connect the EGG cable (and optional SpO<sub>2</sub> cable) to the cable cover.

Secure to the **Front Body** assembly with four Torx screws

**Note:** Ensure the heart symbol is not upside-down.




Figure 16: ECG Cable Cover Installed



Figure 17: Correct Orientation

Insert the **Rotary Selection** switch as shown.

Turn over the **Front Body** assembly and place the washer over the encoder switch.

Tighten the nut and press on the **Rotary Selection Knob**.

**Note:** When properly installed, the knob rotates freely and clicks when pressed.

**Caution:** The connection between the encoder switch and cable is fragile. Do not excessively bend or twist the cable



Figure 18: Encoder Switch (Rear View)

Apply a 5 mm (3/16 in) bead of silicone to seal the joint between the **Front Body** and the **Lower Frame** (between the two points as shown).

**Note:** The seal must be complete (no gaps).




Figure 19: Lower Frame Silicone Application

Assembly Step	Details
---------------	---------

Connect the **CPU Main Cable** to **Main Board J211**.



**Figure 20: CPU Main Cable**

## Power Control Board and Power Supply Installation

### Assembly Step

### Details

Pass the three cables through the holes in the **Lower Frame** as shown.

Secure the **Lower Frame** to the **Front Body** with six Torx security screws.




Figure 21: Lower Frame (Bottom View)

Route the **ECG Cable** as shown and pass it back up into the main compartment.



Figure 22: ECG Cable Route through Lower Frame

Install the **Fan**.

**Note:** An arrow on the fan indicates air flow direction. Ensure the fan blows upwards into the unit.










Figure 23: Fan Installed

Assembly Step	Details
<p>Snap the two black <b>Brackets</b> into place on the <b>Lower Frame</b>.</p> <p><b>Note:</b> Check the <b>Power Supply</b> green ground wire to ensure it is properly soldered on both ends and is in good condition.</p> <p>Place the <b>Power Supply</b> on the <b>Brackets</b>.</p> <p><b>Note:</b> The <b>Power Supply</b> is not secured until the <b>Base Cover</b> is installed.</p>	 <p>Figure 24: Brackets Installed</p>  <p>Figure 25: Power Supply Installed</p>
<p>Place the <b>AC Cable Connector</b> into the <b>Lower Frame</b> as shown.</p> <p><b>Note:</b> The single ground pin must be at the top.</p>	 <p>Figure 26: AC Cable Connector</p>  <p>Figure 27: Proper Connector Orientation</p>
<p>Insert the <b>Power Button</b> and <b>Control Button</b> assemblies into the <b>Lower Frame</b>.</p>	 <p>Figure 28: Power and Control Button Assemblies Installed</p>
<p>Pass <b>Therapy Power Cable</b> into the <b>Lower Frame</b> as shown. Pull the slack through, leaving enough length to connect to <b>Power Control Board</b> (installed in the next step).</p>	 <p>Figure 29: Therapy Power Cable</p>



## Assembly Step

## Details

Secure the **Power Control Board** with five Torx screws.

Place the **Service Connector** (serial connector) on the **Lower Frame**.

**Note:** Pin 1 must be at the top.

Route and connect all cables as shown.

**Note:** See the Ribbon Cable Installation Notes for proper installation (Figures 6, 7 and 8).

Install the **Battery Release** and extension spring.

**Note:** To help the **Battery Release** slide under the **Power Control board**, back off the screws holding the **Power Control board** a turn or two and then retighten after installation.

If necessary, reinstall rubber tubing along the edge of the **Lower Frame**.

Two lengths are required:

- 705 mm  $\pm$  3 mm (27-3/4 in  $\pm$  1/8 in)
- 148 mm  $\pm$  2 mm (5-15/16 in  $\pm$  1/16 in)

**Note:** Do not stretch the tubing as it is inserted into the groove.

**Note:** The gap between the end of the groove and the start of the tubing must be 2 mm (1/16 in) or less.




Figure 30: Power Control Board Installed



Figure 31: Service Connector (Correct Orientation)



Figure 32: Battery Release Installed



Figure 33: Tubing Properly Installed




Figure 34: Tubing Not Properly Installed

Assembly Step	Details
---------------	---------

Install the **Battery Latch** in the **Base Cover**.

Apply Silicone to the notch for the AC power plug.



Figure 35: Battery Latch



Figure 36: Silicone Applied

Install the **Base Cover** to the **Lower Frame**. Using five Torx security screws (white arrows).

**Note:** Tighten the screws slowly and evenly to ensure the rubber tubing is seated correctly and no wires are pinched.

When all screws are tightened, check for any gaps or bulges in the seal.

Install the two larger Torx security screws (black arrows) to hold the **Power Supply** in place.

Install the four rubber feet if necessary.




Figure 37: Base Cover

If necessary, check the joint between the **Lower Frame** and the **Front Body** for gaps in the silicone. Add more silicone if necessary.



Figure 38: Front Body and Lower Frame Joint (with Correctly Applied Silicone)

## ECG Board Installation

### Assembly Step

### Details

Install the **ECG Insulator** in the **Front Body**.



Figure 39: ECG Insulator Installed

Install the **ECG Board** to the **Front Body** by connecting **J412** on the **ECG Board** to **J212** on the **Main Board**.

Install the seven Torx screws to hold the **ECG Board** in place.



Figure 40: ECG Board Installed

Connect the ECG cable to the **ECG Board**.

- J4231 White
- J4232 Black
- J4233 Red
- J4234 Brown
- J4235 Green
- J4236 Gray

**Notes:** Gently pull on each connector to ensure it is locked in place.

The extra gray connector will be connected later.






Figure 41: ECG Cable Installation Detail



Assembly Step	Details
---------------	---------

**Note:** The next section details the **Front Bezel** and **Display** replacement. That procedure includes steps to route the ECG cable and install the ECG Shield for the **ECG Board**. Even if the **Display** or **Front Bezel** is not replaced, those steps must still be completed (see Figure 47 and Figure 48)

## Display Installation

Assembly Step	Details
<p>If necessary, snap the <b>LCD</b> into the <b>LCD Retainer</b>.</p> <p>Connect the LCD ribbon cable to <b>J214</b> on the <b>Main Board</b>.</p> <p><b>Note:</b> See the Ribbon Cable Installation Notes for proper installation (Figures 6, 7 and 8).</p> <p>Secure the <b>LCD</b> to the <b>Front Body</b> with five Torx screws (two in front, three in back).</p>	  <p><b>Figure 42: LCD Installed in LCD Retainer</b></p> <p><b>Figure 43: LCD Screw Locations</b></p>
<p>Connect the LCD backlight cable to <b>J213</b> on the <b>Main Board</b>.</p> <p>Remove the protective film from the <b>LCD</b>.</p> <p><b>Caution: Equipment Damage.</b> Do not touch the front of the LCD.</p>	 <p><b>Figure 44: LCD Backlight Cable (Properly Routed)</b></p>
	<p><b>Caution: Shock Hazard or Equipment Damage.</b> The Responder 2000 is energized for the following step to verify proper operation. Exposed circuit boards may contain potentially dangerous voltages. Do not touch exposed electronics and keep all tools clear while performing the verification.</p> <p>At this point the assembly may be verified by connecting AC Power and turning on the unit. The display should be clear and bright. When the check is complete, turn the power off and disconnect the AC power cord.</p> <p><b>Note:</b> This test will cause a service code error because the <b>Therapy board</b> is not installed. Delete this error when reassembly is complete.</p>

## Assembly Step

## Details

If necessary, remove protective paper and plastic from inside the **Front Bezel**.

**Caution: Equipment Damage.** Do not touch the inside of the **Front Bezel**.

Clean the LCD and inside of the Front Bezel with a brush or compressed air. Ensure the plastic is clean and not scratched.

Install the Front Bezel with six Torx security screws (white arrows) and one Torx screw (black arrow).

**Note:** Install the **Front Bezel** before assembling the circuit boards on the main compartment of the unit; otherwise you might not be able to access the Front Bezel screws.

Install the **Printer Door** into the **Front Body** (insert the left side and then snap the right side into place).

Note the roller orientation. The left side of the door is inserted first.

Open and close the door several times to ensure proper installation.




Figure 45: Front Bezel Installed



Figure 46: Printer Door Installation

Route the ECG cable and connect the remaining gray ECG cable connector to the clip in the **ECG Shield**.

Install the **ECG Shield** in the **Front Body** with three Torx screws.




Figure 47: ECG Cable (Properly Routed)












Figure 48: ECG Shield Installed

## SpO<sub>2</sub> Board Installation

Assembly Step	Details	
<p>Perform the Following step only if the SpO<sub>2</sub> option is installed.</p> <p>Connect the ECG/ SpO<sub>2</sub> cable to <b>X16</b> on the <b>SpO<sub>2</sub> Board</b>.</p> <p>Connect <b>X14</b> on the <b>SpO<sub>2</sub> Board</b> to <b>J221</b> on the <b>Main Board</b> and then secure with four Torx screws.</p>	 A close-up photograph showing a multi-colored cable being plugged into a connector on a small circuit board.	 A photograph showing the SpO2 board being inserted into a slot within a larger electronic assembly.
	 A photograph of the SpO2 board fully seated in its slot, with four yellow arrows pointing to the Torx screws used for securing it.	<p><b>Figure 51: SpO<sub>2</sub> Board Installed</b></p>

## Therapy Board Installation

Assembly Step	Details
<p>Connect the ribbon cable to <b>J110</b> on the <b>Therapy Board</b>.</p> <p><b>Note:</b> See the Ribbon Cable Installation Notes for proper installation (Figures 6, 7 and 8).</p> <p>Fold the connected cable under the <b>Therapy Board</b>.</p> <p>Install the <b>Therapy Board</b> with four Torx screws.</p>	 <p><b>Figure 52: Therapy Board Ribbon Cable Installation</b></p>  <p><b>Figure 53: Ribbon Cable Routing</b></p>  <p><b>Figure 54: Therapy Board Installed</b></p>
<p>Connect the following cables to the <b>Therapy Board</b>:</p> <ol style="list-style-type: none"> <li>1. J103: Therapy Power cable</li> <li>2. J217 Charge Drain cable (left)</li> <li>3. J216 Charge Drain cable (right)</li> </ol> <p><b>Note:</b> Use one finger to support the <b>Therapy Board</b> when making the connections.</p>	 <p><b>Figure 55: Therapy Board Support</b></p>  <p><b>Figure 56: Therapy Board Cables Installed</b></p>

## Capacitor Installation

### Assembly Step

Insert one tie wrap through each capacitor support bracket. The tie wraps are used to secure the **Therapy Board** power cable, left Charge Drain cable, and the **Capacitor** lead wires after the lead wires are connected.

Place each **Capacitor** in its support bracket.

Connect the lead wires to the **Therapy Board**:

- J105 & J107: White lead wires
- J104 & J106: Red lead wires

Secure the **Therapy Board** power cable, left Charge Drain cable, and the **Capacitor** lead wires with the tie wraps.

**Caution:** Equipment Damage. Misrouting or leaving too much slack in the wires may cause wires to contact hot components or become damaged during reassembly. Ensure all wires are routed and all plugs are oriented exactly as shown in Figure 58: Capacitors Installed.

Specific areas include:


- Red capacitor wires must have no slack between the board connectors and the tie wrap to keep the wires away from the **Therapy Board** components.
- The red/black cable bundle must be routed down the side of the left **capacitor**, not behind it.
- Ensure the lead wires are routed behind the plastic center mount.

**Caution:** Equipment Damage. When installing the capacitors, avoid excess strain on the capacitor leads. Excessive force can break or damage the lead connection.

### Details




**Figure 57: Tie Wrap Installation Detail**



**Figure 58: Capacitors Installed**



**Figure 59: Lead Wire Installation Detail**



**Figure 60: Capacitor Lead Break**



## Rear Cover Installation




Assembly Step	Details
<p>If necessary, insert Rubber Tubing (279 mm <math>\pm</math> 3 mm 11 in <math>\pm</math> 1/8 in) into the groove on the <b>Front Body</b>.</p> <p><b>Note:</b> See Figure 33 for installation note.</p>	
<p>Install <b>Handle</b> to the <b>Rear Body</b> with two Torx security screws and two washers.</p>	
<ol style="list-style-type: none"> <li>1. Tie cables to <b>Rear Body</b>.</li> <li>2. Pull patient cables to a length of 5.5" to 6", measured from the plastic rib shown to the top of the plugs.</li> <li>3. Tuck service loop under patient connector.</li> </ol>	

Figure 63: Rear Cover Cable Installation

## Assembly Step

## Details

1. Connect **P101** (Red Paddle wire) from **Rear Body** to **J101** on the **Therapy Board**. Then connect **P102** to **J102**.
2. Connect gray Paddle Cable to **J109** of **Therapy Board**.

**Notes:** Ensure the wires are routed exactly as shown.

The red patient wires (labeled 1) must go under the white cap wires. Also, both plugs of the red patient connector wires must be pointing left.

**Caution: Equipment Damage.** J122 is not used. Ensure no cable is connected to this jack.

To attach the **Rear Body**, lift the **Rear Body** and twist 180 degrees counterclockwise to assemble into position.

**Note:** Ensure the red and gray wire bundles are to cross as shown.




Figure 64: Paddle Cables Connected



Figure 65: Rear Body Placement






Figure 66: Rear Body Assembly



Assembly Step	Details
<p>Install <b>Rear Body</b> to <b>Front Body</b> using five small Torx security screws and two large Torx security screws.</p> <p><b>Note:</b> Tighten all screws slowly and evenly. Check for bulges, gaps, or pinched wires before tightening completely.</p> <p>Install <b>Bed Hook Covers</b> using two Phillips screws.</p> <p><b>Note:</b> Ensure <b>Bed Hook Covers</b> are properly installed.</p>	 <p>Figure 67 shows the rear body of the device installed on the front body. The rear body is a light-colored plastic component with a central handle and two side-mounted hook covers. A label on the bottom of the rear body reads: "Not for Human Use", "Pending 510(k) submission", "Not available for sale in the United States".</p>  <p>Figure 68 shows a close-up of the bed hook cover on the left side of the device. The cover is improperly installed, appearing to be on the right side instead of the left side.</p> <p><b>Figure 68: Bed Hook Cover Improperly Installed (Right Cover on Left Side)</b></p> <p><b>Figure 67: Rear Body Installed</b></p>

## Section 5: Performance Verification and Safety Testing

### Overview

This section describes how to perform service maintenance and testing of the Responder 2000 after repairs are performed.

<b>TOPIC</b>	<b>PAGE</b>
Required Equipment	51
Preventive Maintenance	51
Annual Inspection	51
Verify Operation	52
Maintenance Menu Reference	53
Visual Inspection	59
Operational Tests	60
Safety Tests	71
Alternate Safety Test Procedures	77

## Required Equipment

- Fluke Model Index 2 SpO2 simulator or (equivalent)
- Fluke Model Impulse 4000 Defib analyzer (or equivalent)
- Metron QA-90 Safety Analyzer (or equivalent)
- UadTech Sentry 30 AC/DC/IR Hipot Tester

**PRECAUTION:** Do not shock into devices with highly non-linear resistance versus the voltage, (e.g. an overvoltage protector). **This may damage the Responder 2000. Do not use** the HELLIGE Test Lamp 301 495 00 for testing. This test lamp includes an overvoltage protector. Shocking into it may damage the Responder 2000.



HELLIGE Test Lamp

## Preventive Maintenance

The following preventative maintenance and operational checks are normally performed only by the user:

- Cleaning the Responder 2000 or accessories
- Calibrating the battery fuel gauge
- Daily maintenance per Defibrillator Checklist
- Other routine maintenance

Refer to the Operator's Manual for more information.

## Annual Inspection

A qualified service technician should check the Responder 2000 once per year for the following:

1. Visual Inspection (page 59)
2. Operational Check (starting on page 60)
3. Printer Test (page 70)

### Verify Operation

Before performing service on the Responder 2000, use the following procedures to configure and test the Responder 2000 to verify proper operation (or locate a fault condition):

- Visual Inspection (page 59)
- Operational Tests (starting on page 60)

### Checkout Procedures

After service, the following operational and safety procedures must be performed to verify proper operation. All tests indicated must be performed before the Responder 2000 can be returned to the customer. Depending on the type of repair, the following procedures must be performed for each of the following cases:

- External Repairs or Replacement

External repairs/replacements do not require opening the case. After repair or replacement, perform a Visual Inspection (page 59).

- Printer Door Replacement

The printer door may be replaced without opening the case. After the printer door is replaced (and the case was not opened) the following tests must be performed:

1. Visual Inspection (page 59)
2. Printer Test (page 70)

- Internal Repairs

If the case is opened, the following tests must be performed after the repair:

1. Visual Inspection (page 59)
2. Operational Check (starting on page 60)
3. Safety Checks (starting on page 71)

## Maintenance Menu Reference

### Caution: Loss of Patient Data


Do not enter Maintenance mode while the Responder 2000 is connected to a patient. Patient data may be lost. Perform any Maintenance mode procedures outside the patient vicinity.

**Note:** The Responder 2000 battery must have at least a 20% charge before performing these procedures.

The following Maintenance Menu items are listed as a reference for quick troubleshooting or verification and are not intended as a comprehensive operational or safety checkout. The entire checkout procedures listed above must be performed before returning the Responder 2000 to the customer.



**Note:** A password is required to access the Maintenance menu. The default password is HEART1, but may have been changed by the facility. To enter a password:




1. Select **Maintenance** from the **System** menu. The password field is displayed, highlighted.
2. Click to highlight the first character field.
3. Click the highlighted field, and then rotate the encoder button to the first letter of the password.
4. Click to accept the letter and then rotate to the next field.
5. After all password characters are entered, rotate one space past the last character field. The border around the password field highlights.
6. Click to highlight the password field and then rotate and click **Accept**.




The Maintenance Menu has the following items:

- Service Tests: These tests, described below, are used to verify operation
- Device Configuration: Used to set the language
- Stored Data Management: Used to clear stored data or settings, transfer stored logs to a PC, and save or load stored settings

## Service Tests




## Button Test




- Press and hold each button to test. **On** while the button is pressed (and **Off** when not pressed) indicates the button is working properly.
- To test the **Rotary Selection** knob, rotate to highlight **Rotary Test** and press the **Rotary Selection** knob. **Rotary** should indicate **On**.

### Printer Speed Test

Prints a test strip at a constant speed of 25 mm/s with a 25 mm period and 25 mm amplitude:




### Battery Readings




Heading	Description
Full Capacity	Capacity of a fully charged battery
Rem. Capacity	Current battery charge
Rel. Capacity	Current charge as a percentage of the full capacity

Normally this screen is used to verify the battery is charging (or fully charged) at the appropriate voltage. Service tests must be performed with the battery at least 20% charged.

- New batteries have a full capacity of at least 6000 mAh and will charge when the Responder 2000 is plugged into an AC outlet.
- If the battery is fully charged and has a capacity of less than 4500 mAh, calibrate the battery (see Battery Calibration on page 11).
- If the battery still has a capacity of less than 4500 mAh after calibration, then replace the battery.

**Display Test**

The following colors fill the entire display in succession: White, Red, Green, Blue, and Black. Verify no pixels are stuck or missing. After the test, the Responder 2000 then automatically cycles power.

**ECG Lead Test**

Displays lead status (on or off). See the ECG Lead Test Procedure on page 69 for proper indications.

**Notes:**

- If the lead RA, LA, LL or V1 is not connected (to either a Patient or ECG Simulator) or the lead is broken, the Lead status displays **OFF**.
- If (Driven) RL lead is not connected or broken, all lead status are displayed **OFF**.



### Temperature Readings



Displays internal temperatures for the **Power Board**, **Therapy Board**, **Printer**, and **Battery**. These readings may vary according to fan speed or external conditions. An abnormally high reading may indicate a fan not working or other fault.

The internal fan should turn on at the following temperatures:

	Fan ON (High)	Fan ON (Low)	Fan OFF
<b>Power Board</b>	65°	60°	55°
<b>Therapy Board</b>	65°	60°	55°
<b>Printer</b>	N/A	N/A	N/A
<b>Battery</b>	N/A	N/A	N/A

### Fan Test



Select **Off**, **Low**, or **High**

**Note:** The Low fan speed may be difficult to hear in a noisy environment.

## Device Configuration

### Select Language



Sets the system language to: English, French, German, Italian, Dutch, Portuguese, Russian, Spanish, or Swedish. Once a new language is selected, the Responder 2000 cycles power.

**Note:** **Maintenance** mode text is always in English, regardless of the language selected.

## Stored Data Management

### Clear Event Logs

Clears all stored event logs. Always clear event logs after the **Main** board is replaced.

### Clear Service Required Indicator

Clears the **Service Required** LED. Once a service fault is corrected, use the **Clear Service Required** menu item to turn off the **Service Required** LED. If the LED does not clear, another fault is present.

### Clear Settings Menu Password

Clears the **Settings Menu** password.

### Clear All Stored Data

Clears all stored data, including:

- Event Logs
- Service Required Indicator
- Settings Menu Password
- All other Settings.


### Send Event Logs to PC

Event logs may be downloaded to a PC and sent to Technical support to help troubleshoot a problem.

1. Go to **Service Mode** (password HEART1).
2. Select **Stored Data Management > Send Event Logs To PC**.
3. On the PC, start the Windows-based ServiceLink2000 Application.
4. Select the appropriate Com Port and click **Start Com**.
5. Click **Receive Events**. A progress bar is displayed.
6. When downloaded, save the events file to the desired location. A progress bar is displayed while the CSV file is saved.
7. The ServiceLink2000 program closes after the file is successfully saved.

### Send/Receive Stored Parameters

Stored parameters should be downloaded before replacing a board, and then uploaded after the board is replaced to return the system to its original state.



To Retrieve Stored Parameters:

1. Go to **Service Mode** (password HEART1).
2. Select **Stored Data Management > Send/Receive Stored Parameters**.
3. On the PC, start the Windows-based ServiceLink2000 Application.
4. Select the appropriate Com Port and click **Start Com**.
5. Click **Receive Parameters**. A progress bar is displayed.
6. When downloaded, save the parameter file to the desired location. A progress bar is displayed while the file is saved.
7. The ServiceLink2000 program closes after the file is successfully saved.

To Send Stored Parameters

1. Go to **Service Mode** (password HEART1).
2. Select **Stored Data Management > Send/Receive Stored Parameters**.
3. On the PC, start the Windows-based ServiceLink2000 Application.
4. Select the appropriate Com Port and click **Start Com**.
5. Click **Send Parameters**.
6. Select the parameters to be uploaded and click **Open**. A progress bar is displayed.
7. Check the Responder 2000 display. If *Params Set* is displayed, the upload was successful.
8. Click No to upload the parameters file again. If the parameters file was not uploaded successfully, retry as necessary.
9. The ServiceLink 2000 program closes after the file is successfully uploaded.

To Retrieve Stored Waveforms:

1. Go to **Service Mode** (password HEART1).
2. Select **Stored Data Management > Send/Receive Stored Parameters**.
3. On the PC, start the Windows-based ServiceLink2000 Application.
4. Select the appropriate Com Port and click **Start Com**.
5. Click **Receive Waveforms**. A progress bar is displayed.
6. When downloaded, save the waveform file to the desired location. A progress bar is displayed while the file is saved.
7. The ServiceLink2000 program closes after the file is successfully saved.

### Visual Inspection

Visually inspect the unit and accessories for foreign substances, damage, cracks, bent or discolored pins, broken parts, or extreme wear, including:

- Cracked, bulging, or otherwise damaged case
- Cracked, loose, or bent connectors or plugs
- Rattling when gently shaken
- Discoloration or signs of overheating
- Obvious scratches or damage to the display
- Loose battery or paddle clips
- Loose printer door

### Operational Tests

#### Initial Power Up

Perform this procedure for the following scenarios:

- With a charged battery (at least 25%) installed and connected to AC power
- With the battery removed and connected to AC power
- With a charged battery (at least 25%) and AC power disconnected

Turn on the Responder 2000 and verify the following during power up:

- No alarms (at any time during power up or operation)
- The speaker emits a short beep, then is silent
- The following indicators momentarily light:
  - Service Required: Red
  - Charge: off
  - Manual: Blue
  - Shock: Red
- The Responder 2000 is ready for use in about 4 seconds.

The AC Power LED and Battery Charging LED are lit when a battery is installed and the Responder 2000 is connected to AC power

- The AC Power LED is not lit and Battery Charging LED is not lit when a battery is installed and the Responder 2000 is not connected to AC power

## Maintenance Mode Tests

Perform each Maintenance Mode Tests (starting on page 54) to ensure proper operation.

### Verify Serial Communications

#### Required Equipment


- Windows Based PC
- Serial Cable

#### Procedure

1. Turn off the Responder 2000.
2. Shut down the PC.
3. Connect the Serial cable from the Serial port on the Responder 2000 to the Serial port on the PC.
4. Turn on the PC.
5. Open Responder 2000 Upgrade program on the PC.
6. Select the appropriate COM port and press **OK** (normally COM 1).
7. Follow the on-screen instructions as shown:



8. Verify the CodeLink 2000 window is displayed as shown:



9. Press **Exit** to close the window.
10. Power down the Responder 2000.
11. Shut down the PC.
12. Disconnect the Serial cable from the Responder 2000 and the PC.

### Defibrillator Tests

#### Warning! Shock Hazard.

Before delivering any shock, ensure all personnel and equipment are clear of the test area.

#### Energy Timeout Test

1. Ensure the Responder 2000 has a charged battery (at least 25%) installed and is connected to AC power
2. Seat each paddle in its side holder. Ensure paddles are securely seated.
3. Press the **Manual** button on the Responder 2000.
4. Select an **Energy Level** of **100J**.
5. Press the **Charge** button on the front panel.
6. The message **Stand Clear** is displayed.
7. Do not press the Shock button. After 30 Seconds the message **Charge Time Out** should be displayed.

#### Defibrillator Test Chart

For each of the defibrillator tests below (Paddles, Pads, and Spoons), perform the test at the energy level indicated in the procedure. Then perform each test at the energy levels indicated in the chart below.

**Note:** Use only a calibrated Impulse 4000 (or equivalent) for pads and spoons. Energy values may vary slightly depending on the simulator used.

## Section 5: Performance Verification and Safety Testing

Selected Energy	Energy Delivered into 50 Ohm (Tolerance +/- 15 % or +/-3 Joule)	Test
2	0-5	Paddles, Pads, and Spoons
3	0-6	Paddles, Pads, and Spoons
5	2-8	Paddles, Pads, and Spoons
7	4-10	Paddles, Pads, and Spoons
10	7-13	Paddles, Pads, and Spoons
15	12-18	Paddles, Pads, and Spoons
20	17-23	Paddles, Pads, and Spoons
30	25.5-34.5	Paddles, Pads, and Spoons
50	42.5-57.5	Paddles, Pads, and Spoons
70	59.5-80.5	Paddles and Pads only
100	85-115	Paddles and Pads only
150	127.5-172.5	Paddles and Pads only
200	170-230	Paddles and Pads only
270	229.5-310.5	Paddles and Pads only

### Paddles Test

1. Ensure the Responder 2000 has a charged battery (at least 25%) installed and is connected to AC power
2. Connect defibrillator paddles connector to the connector on the Responder 2000.
3. Seat each paddle in its side holder. Ensure paddles are securely seated.
4. Press the **Manual** button on the Responder 2000.
5. Select **Energy Level of 270J**.
6. Press **Charge** button on APEX (right) paddle. Verify the Responder 2000 charges to a 270J level and the audible alarm sounds.
7. Simultaneously press and hold the **Shock** buttons on both paddles.
8. Verify that energy delivered is shown on the Responder 2000 display is between 229.5 to 310.5J.
9. Repeat this procedure with input power from:
  - With the battery removed and connected to AC power
  - With a charged battery (at least 25%) and AC power disconnected
10. Repeat this procedure using the energy levels listed in the table above with the following conditions:
  - Do not perform step 9 (this is only required once to verify the energy path)
  - Wait at least 30 seconds between each test shock

### Pads Test

1. Ensure the Responder 2000 has a charged battery (at least 25%) installed and is connected to AC power
2. Connect the pads connector to the Responder 2000.
3. Connect the pads to the Impulse 4000.
4. Press the **Manual** button on the Responder 2000.
5. Select **Energy Level of 270J**.
6. Press **Charge** button on front panel. Verify the Responder 2000 charges to a 270J level and the audible alarm sounds.
7. Press the **Shock** buttons on the front panel.
8. Verify that energy delivered is shown on the Responder 2000 display is between 229.5 to 310.5J.
9. Repeat this procedure using the energy levels listed in the table above with the following conditions:
  - Wait at least 30 seconds between each test shock

### Spoons Test

1. Ensure the Responder 2000 has a charged battery (at least 25%) installed and is connected to AC power
2. Connect the spoons connector to the Responder 2000.
3. Connect the spoons to the Impulse 4000.
4. Press the **Manual** button on the Responder 2000.
5. Select **Energy Level of 50J**.
6. Press **Charge** button on front panel. Verify the Responder 2000 charges to a 50J level and the audible alarm sounds.
7. Press the **Shock** buttons on the front panel.
8. Verify that energy delivered is shown on the Responder 2000 display is between 42.5-57.5J.
9. Repeat this procedure using the energy levels listed in the table above with the following conditions:
  - Do not select an Energy level greater than 50 Joules.
  - Wait at least 30 seconds between each test shock



## Pacing Test

1. Connect pads cable to the back of therapy connector on the back of Responder 2000.
2. Connect pads to a Calibrated Impulse 4000 Defib/Pacer Analyzer or equivalent.
3. Turn on the Responder 2000.
4. Use **Rotary Selector Knob** to move focus to and select **Pacing Info** box in the lower right of the display.  
Note: If pacing is not enabled, the **Pacing Info** box is not selectable.
5. Turning the **Rotary Selector Knob** clockwise, the following items should be selectable::
  - Pacing Type (Fixed, Demand), Default is Demand.
  - PPM ( 30-180 PPM, increments of 5), Default is 60 PPM
  - Exit Pacing Info Box ("x" in corner)
  - Pacing On/Off/Pause (Pause is only available when Pacing is turned On)
  - Pacing Current (0mA-140mA, increments of 5). This option is only selectable when pacing is turned on.
6. Turn the **Rotary Selector Knob** to turn on the pacing.
7. When Pacing is turned on, focus automatically goes to the **Pacing Current** selection. Set desired pacing current.
8. Turn the **Rotary Selector Knob** to select pacing rate and set desired options for pacing rate.
9. Turn the **Rotary Selector Knob** to select pacing mode and set desired options for pacing mode.
10. Check each pacing current value using the chart below.

**Notes:**

- Use only a calibrated Impulse 4000 Defib/Pacer analyzer.
- The measured values may vary slightly depending on the simulator used.
- Wait for at least 30 seconds between each current setting.

## Section 5: Performance Verification and Safety Testing

Pacing Mode = Fixed Pacing Rate = 60 PPM	
Selected Current in mA on Responder 2000	Current Delivered into 50 Ohm setup into Impulse 4000 or equivalent. Tolerance +/- 15% in mA
40	34 – 46
45	38.25 – 51.75
50	42.5 – 57.5
55	46.75 – 63.25
60	51 - 69
65	55.25 - 74.75
70	59.5 - 80.5
75	63.75 - 86.25
80	68 - 92
85	72.25 - 97.75
90	76.5 - 103.5
95	80.75 - 109.25
100	85 - 115
105	89.25 - 120.75
110	93.5 - 126.5
115	97.75 - 132.25
120	102 - 138
125	106.25 - 143.75
130	110.5 - 149.5
135	114.75 - 155.25
140	119 - 161

## ECG Tests


### Required Equipment

- Impulse 4000 (or equivalent)

### ECG Connections Test

1. Verify Responder 2000 is powered on.
2. Connect ECG leads (RL, RA, LA, LL, V1) to the corresponding connectors on the ECG Simulator.
3. Plug in the ECG Cable the ECG Cable Port on the Responder 2000.
4. Change ECG source to **ECGII**.
5. Select **ECG > Norm > 80** on the ECG Simulator and verify the following on the Responder 2000:
  - The ECG signal on screen of Responder is correct amplitude and polarity
  - An audible beep sounds for each ECG beat

## Lead Detection Test



The diagram above indicates where the lead source is displayed.

1. Verify Responder 2000 is powered on.
2. Verify default Lead Source on Responder 2000 is **ECG II**.
3. Connect the Spoons cable between Responder 2000 and Impulse 4000 (Apex and Sternum).
4. Select **Spoons** on the Responder 2000 display as the lead source and ensure the signal is displayed.
5. Disconnect Spoons and connect Paddles cable to the Responder 2000. Ensure paddles are securely stored on left and right side of Responder 2000.
6. Select **Paddles** on the Responder 2000 display as the lead source and ensure the signal is displayed.
7. Disconnect paddles cable from Responder 2000 and connect pads cable to Responder 2000 and Apex/Sternum of Safety Analyzer.
8. Select **Pads** on the Responder 2000 display as the lead source and ensure the signal is displayed.
9. Connect ECG Sync Cable to Responder 2000's ECG connector.
10. Place the Responder 2000 in External Sync mode by performing the following:
  - a. Press **Manual** button on Responder 2000.
  - b. Select **2J** with Rotary Selection knob.
  - c. Turn Rotary Selection knob until **No Sync** is highlighted.
  - d. Press Rotary Selection knob to change to **Sync**.
11. Verify **Manual – Ext Sync** is shown on upper left corner of Responder 2000 display.
12. Disconnect Sync cable from Responder 2000 and connect a 3 Lead ECG cable.
13. Use Rotary Selection knob to exit previous menu and then select the Lead source.
14. Verify the following leads may be selected:
  - o Paddles
  - o ECG I
  - o ECG II
  - o ECG III
15. Disconnect 3 Lead ECG cable from Responder 2000 and connect 5 Lead ECG cable to Responder 2000.
16. Verify the following leads may be selected:
  - o Paddles
  - o ECG I
  - o ECG II
  - o ECG III
  - o aVR
  - o aVF
  - o aVL
  - o V

**ECG Lead Test**

1. Connect the ECG leads to the patient simulator.
2. On the Responder 2000, select **ECG Lead Test** from the **Maintenance** menu.
3. Verify **ECG Lead Status** display as shown:

ECG Lead Status	
RA:	On
RL:	Driven
LA:	On
LL:	On
V1:	On

4. Disconnect RA (White Lead) from Patient Simulator and verify the display as shown:
5. Reconnect RA lead to patient simulator.

ECG Lead Status	
RA:	Off
RL:	Driven
LA:	On
LL:	On
V1:	On

6. Disconnect RL (Green Lead) from Patient Simulator and verify the display as shown:
7. Reconnect RL lead to patient simulator.

ECG Lead Status	
RA:	Off
RL:	Driven
LA:	Off
LL:	Off
V1:	Off

8. Disconnect LA (Black Lead) from Patient Simulator and verify the display as shown:
9. Reconnect LA lead to patient simulator.

ECG Lead Status	
RA:	On
RL:	Driven
LA:	Off
LL:	On
V1:	On

10. Disconnect LL (Red Lead) from Patient Simulator and verify the display as shown:
11. Reconnect LL lead to patient simulator.

ECG Lead Status	
RA:	On
RL:	Driven
LA:	On
LL:	Off
V1:	On

12. Disconnect V1 (Brown Lead) from Patient Simulator and verify the display as shown:
13. Disconnect all leads and return to the **Main** menu.

ECG Lead Status	
RA:	On
RL:	Driven
LA:	On
LL:	On
V1:	Off

### Paddles ECG Test

1. Select **Settings > Channel Settings > Paddles**.
2. Connect the Paddles to the Impulse 4000 and verify both the ECG and the HR value are displayed.
3. Select a Heart Rate of 60 BP/min on the simulator and verify if the display shows this value.

### Thermal Printer Operation

1. Open printer door and verify the **Door Open** message is displayed.
2. Remove paper and close printer door and verify **Printer Paper Out** message is displayed.
3. Re-install paper in printer and close door and verify no printer error messages are displayed.
4. Select **Printer** icon and verify smooth, good quality printing with no stuck pixels.
5. Select **Printer** icon to stop printer and verify the patient/device information page is printed and that printing stops.
6. Tear off the paper strip and verify the paper tears properly.

### Cardioversion Test

1. Connect the Paddles to the Impulse 4000.
2. Set Responder 2000 to **Manual Mode**.
3. Change the status from **No Sync Mode** to **Sync Mode**. Verify the **Manual-Sync** and the trigger is displayed to the QRS complex.
4. Simultaneously press and hold the **Shock** buttons on both paddles. Triggering is affected by the peak of the QRS complex of the beat and the Cardioversion function is automatically reset.
5. If the delay time is greater than 60ms, the test has failed.

### SpO<sub>2</sub> Operation

Perform this test only if the Responder 2000 is equipped with the SpO<sub>2</sub> option.

1. Plug in the SpO<sub>2</sub> sensor connector.
2. Clip the SpO<sub>2</sub> sensor to test subject.
3. Verify a SpO<sub>2</sub> reading is displayed with a tolerance of  $\pm 2$  digits.
4. Verify the SpO<sub>2</sub> waveform is printed on the printout.
5. Verify the SpO<sub>2</sub> alarm operation.

## Safety Tests

**WARNING: Explosion or Shock Hazard.**

Before performing any test, read and follow all safety instructions for the test equipment indicated.

**WARNING: Explosion or Shock Hazard.**

The following instructions are written for the test equipment indicated. Performing any test with equipment other than the specific brand and model indicated may cause explosion, electrical shock, and equipment damage.

All safety tests are required whenever the case is opened, regardless of the type of maintenance or service performed.

Safety tests include:

- Hipot Test
- Leakage Current Test

**Note:** If the Hipot and Leakage current tests in this section cannot be performed for any reason, perform the Alternate Safety Test Procedures starting on page 77.

### Hipot Test

**WARNING: Explosion Hazard.**

Internal damage or improper assembly could cause an internal component to explode during Hipot testing. To minimize risk

- Ensure all Responder 2000 covers are in place and securely attached before testing
- Use appropriate personal safety equipment, including a face shield and gloves
- Ensure the area is clear before starting the test.

**WARNING: Shock Hazard.**

Do not touch the Hipot tester, the Responder 2000, or any leads during the test. Ensure all Hipot leads are electrically isolated before starting the test. If the test is performed on an anti-static mat, ensure the leads do not contact the mat during the test.

### Required Equipment

- UadTech Sentry 30 AC/DC/IR Hipot Tester
- Test fixtures
  - ECG Connector Test Fixture. ECG connector with all leads and ground shield shorted together.
  - AC Plug Test Fixture. AC connector with line and neutral shorted together.
  - SpO<sub>2</sub> Connector Test Fixture. Connection for the ground shield.
  - ECG Paddle Test Fixture. Apex and Sternum pins shorted together.
  - Serial Port Test Fixture. Serial port connector with single lead connect to pin 5.

### ECG to Serial 4.0 KVAC Test

1. Setup the Hipot tester:
  - 4.0 KVAC
  - High current limit 4.0 mA
  - Test 60 Sec.: Ramp 2.0 Sec.
  - 120 VAC ( $\pm 10\%$ ) Input
2. Plug in the ECG Connector Test Fixture into the ECG port.
3. Plug the Serial Port Test Fixture into the Serial port.
4. Connect the red hipot tester lead to the ECG Connector Test Fixture lead.
5. Connect the black hipot tested lead to the Serial Port Test Fixture lead.

#### **WARNING: Shock Hazard.**

Before starting the test, ensure all lead connections are electrically isolated.

6. Press **Start** (green button) on the Hipot tester. The test will take approximately 30 seconds.
  - If the test passes (short beep and green Pass LED is lit), continue with the next test.
  - If the test fails (audible alarm and red Fail LED is lit) discontinue testing. The unit must be returned to the manufacturer for troubleshooting and repair.

### AC Line/Neutral to AC Ground 1.5 KVAC Test

1. Setup the Hipot tester:
  - 1.5 KVAC
  - High current limit 4.0 mA
  - Test 60 Sec.: Ramp 2.0 Sec.
  - 120 VAC ( $\pm 10\%$ ) Input
2. Plug in the AC Plug Test Fixture to the AC input connector.
3. Connect the red hipot tester lead to the AC Plug Test Fixture Line and Neutral lead.
4. Connect the black hipot tested lead to the AC Plug Test Fixture Ground lead.

#### **WARNING: Shock Hazard.**

Before starting the test, ensure all lead connections are electrically isolated.

5. Press **Start** (green button) on the Hipot tester. The test will take approximately two minutes.
  - If the test passes (short beep and green Pass LED is lit), continue with the next test.
  - If the test fails (audible alarm and red Fail LED is lit) discontinue testing. The unit must be returned to the manufacturer for troubleshooting and repair.

### ECG to SpO<sub>2</sub> 1.5 KVAC Test

**Note:** This test is only performed on models with the SpO<sub>2</sub> option installed.

1. Setup the Hipot tester:
  - 1.5 KVAC
  - High current limit 4.0 mA
  - Test 60 Sec.: Ramp 2.0 Sec.
  - 120 VAC ( $\pm 10\%$ ) Input
2. Plug in the ECG Connector Test Fixture into the ECG port.
3. Plug the SpO<sub>2</sub> Test Fixture into the Serial port.
4. Connect the red hipot tester lead to the SpO<sub>2</sub> Test Fixture lead.
5. Connect the black hipot tested lead to Serial Port Test Fixture lead.

#### **WARNING: Shock Hazard.**

Before starting the test, ensure all lead connections are electrically isolated.

6. Press **Start** (green button) on the Hipot tester. The test will take approximately two minutes.
  - If the test passes (short beep and green Pass LED is lit), continue with the next test.
  - If the test fails (audible alarm and red Fail LED is lit) discontinue testing. The unit must be returned to the manufacturer for troubleshooting and repair.



### Apex Test Load to Apex and Sternum 3.0 KVDC Test

**Note:** The previous test settings were testing DC voltage. This test procedure requires setup for AC voltage.

1. Setup the Hipot tester:
  - 3.0 KVAC
  - High current limit 2.0 mA
  - Test 60 Sec.: Ramp 2.0 Sec.
  - 120 VAC ( $\pm 10\%$ ) Input
2. Plug in the Apex and Sternum Test Fixture to the ECG Paddle Connector.
3. Connect the red hipot tester lead to the Apex and Sternum load test lead.
4. Connect the black hipot tester lead to the Apex and Sternum lead.

**WARNING: Shock Hazard.**

Before starting the test, ensure all lead connections are electrically isolated.

5. Press **Start** (green button) on the Hipot tester. The test will take approximately two minutes.
  - If the test passes (short beep and green Pass LED is lit), continue with the next test.
  - If the test fails (audible alarm and red Fail LED is lit) discontinue testing. The unit must be returned to the manufacturer for troubleshooting and repair.

### ECG to Apex and Sternum 3.0 KVDC Test

1. Setup the Hipot tester:
  - 3.0 KVAC
  - High current limit 2.0 mA
  - Test 60 Sec.: Ramp 2.0 Sec.
  - 120 VAC ( $\pm 10\%$ ) Input
2. Plug in the ECG Connector test fixture.
3. Plug in the Apex and Sternum Test Fixture to the ECG Paddle Connector.
4. Connect the red hipot tester lead to ECG Connector Test Fixture lead.
5. Connect the black hipot tester lead to the Apex and Sternum lead.

**WARNING: Shock Hazard.**

Before starting the test, ensure all lead connections are electrically isolated.

6. Press **Start** (green button) on the Hipot tester. The test will take approximately two minutes.
  - If the test passes (short beep and green Pass LED is lit), continue with the next test.
  - If the test fails (audible alarm and red Fail LED is lit) discontinue testing. The unit must be returned to the manufacturer for troubleshooting and repair.

### Apex and Sternum to Service Port 3.0 KVAC Test

1. Setup the Hipot tester:
  - 3.0 KVAC
  - High current limit 2.0 mA
  - Test 60 Sec.: Ramp 2.0 Sec.
  - 120 VAC ( $\pm 10\%$ ) Input
2. Plug in the Apex and Sternum Test Fixture to the ECG Paddle Connector.
3. Plug in the Service Port Connector test fixture.
4. Connect the red hipot tester lead to ECG Connector Test Fixture lead.
5. Connect the black hipot tester lead to the Serial Port Connector Test lead.

**WARNING: Shock Hazard.**

Before starting the test, ensure all lead connections are electrically isolated.

6. Press **Start** (green button) on the Hipot tester. The test will take approximately two minutes.
  - If the test passes (short beep and green Pass LED is lit), continue with the next test.
  - If the test fails (audible alarm and red Fail LED is lit) discontinue testing. The unit must be returned to the manufacturer for troubleshooting and repair.

## Leakage Current Test

### Required Equipment

- Metron QA-90 Safety Analyzer (or equivalent)

### Equipment Setup

**Note:** This setup assumes the use of a Metron QA-90 Safety Analyzer. If an equivalent safety analyzer is used, refer to the Operator's manual for detailed instructions.

- Configure the Safety analyzer for following parameters.

**Note:** The following parameters may be saved for future use.

Leads	Name	Limit
	Earth Leakage Current (OS)	1000 µA
	Earth Leakage Current (NC)	500 µA
	Earth Leakage Current (OSRM)	1000 µA
	Earth Leakage Current (NCRM)	500 µA
	Enclosure Leakage Current (OS)	500 µA
	Enclosure Leakage Current (NC)	100 µA
	Enclosure Leakage Current (OE)	5000 µA
	Enclosure Leakage Current (OSRM)	500 µA
	Enclosure Leakage Current (NCRM)	100 µA
	Enclosure Leakage Current (OERM)	500 µA
1-2	Patient Leakage Current AC (OS)	50 µA
1-2	Patient Leakage Current AC (NC)	10 µA
1-2	Patient Leakage Current AC (OE)	50 µA
1-2	Patient Leakage Current AC (OSRM)	50 µA
1-2	Patient Leakage Current AC (NCRM)	10 µA
1-2	Patient Leakage Current AC (OERM)	50 µA
3-7	Patient Leakage Current AC (OS)	50 µA
3-7	Patient Leakage Current AC (NC)	10 µA
3-7	Patient Leakage Current AC (OE)	50 µA
3-7	Patient Leakage Current AC (OSRM)	50 µA
3-7	Patient Leakage Current AC NCRM	10 µA
3-7	Patient Leakage Current AC OERM	50 µA

## Section 5: Performance Verification and Safety Testing

Leads	Name	Limit
8	Patient Leakage Current AC (OS)	500 $\mu$ A
8	Patient Leakage Current AC (NC)	100 $\mu$ A
8	Patient Leakage Current AC (OE)	500 $\mu$ A
8	Patient Leakage Current AC OSRM	500 $\mu$ A
8	Patient Leakage Current AC NCRM	100 $\mu$ A
8	Patient Leakage Current AC OERM	500 $\mu$ A
1-2	Patient Auxiliary Current AC (OS)	50 $\mu$ A
1-2	Patient Auxiliary Current AC NC	10 $\mu$ A
1-2	Patient Auxiliary Current AC OE	50 $\mu$ A
1-2	Patient Auxiliary Current AC OSRM	50 $\mu$ A
1-2	Patient Auxiliary Current AC NCRM	10 $\mu$ A
1-2	Patient Auxiliary Current AC OERM	50 $\mu$ A
3-7	Patient Auxiliary Current AC (OS)	50 $\mu$ A
3-7	Patient Auxiliary Current AC (NC)	10 $\mu$ A
3-7	Patient Auxiliary Current AC (OE)	50 $\mu$ A
3-7	Patient Auxiliary Current AC OSRM	50 $\mu$ A
3-7	Patient Auxiliary Current AC NCRM	10 $\mu$ A
3-7	Patient Auxiliary Current AC OERM	50 $\mu$ A
1-2	Mains on Applied Part (SFC)	50 $\mu$ A
1-2	Mains on Applied Part (SFCRM)	50 $\mu$ A
3-7	Mains on Applied Part SFC	50 $\mu$ A
3-7	Mains on Applied Part SFCRM	50 $\mu$ A
8	Mains on Applied Part SFC	5000 $\mu$ A
8	Mains on Applied Part SFCRM	5000 $\mu$ A
1-2	Patient Leakage Current DC (OS)	50 $\mu$ A
1-2	Patient Leakage Current DC (NC)	10 $\mu$ A
1-2	Patient Leakage Current DC (OE)	50 $\mu$ A
1-2	Patient Auxiliary Current DC OSRM	50 $\mu$ A
1-2	Patient Auxiliary Current DC NCRM	10 $\mu$ A
1-2	Patient Auxiliary Current DC OERM	50 $\mu$ A

## Section 5: Performance Verification and Safety Testing

Leads	Name	Limit
3-7	Patient Auxiliary Current DC (OS)	50 $\mu$ A
3-7	Patient Auxiliary Current DC (NC)	10 $\mu$ A
3-7	Patient Auxiliary Current DC OE	50 $\mu$ A
3-7	Patient Auxiliary Current DC OSRM	50 $\mu$ A
3-7	Patient Auxiliary Current DC NCRM	10 $\mu$ A
3-7	Patient Auxiliary Current DC OERM	50 $\mu$ A
1-2	Patient Leakage Current DC (OS)	50 $\mu$ A
1-2	Patient Leakage Current DC NC	10 $\mu$ A
1-2	Patient Leakage Current DC OE	50 $\mu$ A
1-2	Patient Leakage Current DC OSRM	50 $\mu$ A
1-2	Patient Leakage Current DC NCRM	10 $\mu$ A
1-2	Patient Leakage Current DC OERM	50 $\mu$ A
3-7	Patient Leakage Current DC (OS)	50 $\mu$ A
3-7	Patient Leakage Current DC (NC)	10 $\mu$ A
3-7	Patient Leakage Current DC OE	50 $\mu$ A
3-7	Patient Leakage Current DC OSRM	50 $\mu$ A
3-7	Patient Leakage Current DC NCRM	10 $\mu$ A
3-7	Patient Leakage Current DC OERM	50 $\mu$ A
8	Patient Leakage Current DC (OS)	50 $\mu$ A
8	Patient Leakage Current DC (NC)	10 $\mu$ A
8	Patient Leakage Current DC OE	50 $\mu$ A
8	Patient Leakage Current DC OSRM	50 $\mu$ A
8	Patient Leakage Current DC NCRM	10 $\mu$ A
8	Patient Leakage Current DC OERM	50 $\mu$ A

2. Plug Cord from 120V, 60hz receptacle of Safety Analyzer into the wall outlet.
3. Plug cord from **Auxiliary Power to Equipment Under Test** receptacle of Safety Analyzer into Step-up transformer output (264VAC).
4. Plug power cord from Responder 2000 into outlet on front of Safety Analyzer.
5. Plug jumper cable from **ENCL** connector to **EARTH** connector.
6. Plug the 5 ECG cable leads from ECG Cable Port into **Patient Lead** connectors **3-7**.
7. Plug cable from Dfib Paddle Port into **Patient Lead** connectors **1** and **2**.
8. (For SpO<sub>2</sub> equipped Responder 2000 only) Plug cable from SpO2 connector to **Patient Lead** connector 8.

**Test Procedure**

**Note:** This procedure assumes the use of a Metron QA-90 Safety Analyzer. If an equivalent safety analyzer is used, refer to the Operator's manual for detailed instructions.

1. Press Responder 2000 **Power** button to turn on.
2. Turn on the Safety Analyzer and wait for the self test to complete.
3. If previously saved, select the stored Responder 2000 parameters.
4. Press **START** on the Safety Analyzer to begin the tests.
5. After test sequence is complete, verify TEST PASSED is displayed on the Safety Analyzer.


**Alternate Safety Test Procedures**

The following procedures are provided in case the recommended test procedures (above) cannot be performed.


**Test Equipment**

Number	Equipment Description	Qty
	Fluke Model Index 2 SpO2 simulator	1
	Fluke Model ECG Simulator Defib analyzer or adequate	1
	Metron QA-90 Safety Analyzer or adequate	1
MarqIII-Kit	Marq-3 EKG Simulator or adequate	1
2025269-003	RESPONDER 2000/PH ECD PACER TESTER	1


## Leakage Current “Measurement Device” (MD) Characteristics



a) Measuring device



b) Frequency characteristics

NOTE The network and voltage measuring instrument above are replaced by the symbol  in the following figures.

a) Non-inductive components

b) Resistance  $\geq 1 \text{ M}\Omega$  and capacitance  $\leq 150 \text{ pF}$

c)  $Z(f)$  is the transfer impedance of the network, i.e.  $V_{out}/I_{in}$ , for a current of frequency  $f$ .

To perform the leakage current measurements, the unit under test (UUT) or device under test (DUT) has to be separated from any interconnection to a system. If the UUT is part of a system, extended tests according to IEC 60601-1-1 have to be performed. The figure above shows the Measuring Device (MD) circuit required for leakage current. The reading in mV corresponds to  $\mu\text{A}$  (leakage current). The Safety Testers generally work with this MD circuit and the displayed values are already converted to leakage current.

## Tests

The “Electrical Safety” tests may be performed under normal ambient conditions of temperature, humidity and pressure and line voltage. Use test “Devices” or “Measuring Devices” as shown or equivalent.

## Visual Inspection

Verify instrument and accessories to ensure that;

- Fuse links have the rating proclaimed by the manufacturer
- Safety labels and inscriptions on the device are clearly legible
- The mechanical condition will allow the device to be put to further use
- Any soiling has no effect on safety of operation.

The defibrillator electrodes as well as handles and holder recesses must be free of any cream residue.

The defibrillator electrodes, pacemaker cable, ECG cable, SPO<sub>2</sub> cable and Responder 2000 power cord should be checked for any visible external damage to the insulation and strain relief.

### Leakage Current

The leakage currents correspond to 110 % of rated voltage for the tested unit. Most Safety Testers take this into account; otherwise the measured values have to be calculated.

- Use of a “Safety Tester” as show or equivalent test setup (i.e. equivalent test configurations need not require/use S4).
- Ensure the “Safety Tester” and UUT are configured (if applicable) for the correct voltage.
- Connect the MD (Measuring Device) to the applicable ports on the “Safety Tester” or equivalent test points per measurement.
- Connect the “Safety Tester” to the mains supply.
- Connect the unit under test to your Safety Tester.

See section 0 for required measuring device circuit and characteristics.




Figure 69: Basic measurement configuration

## Section 5: Performance Verification and Safety Testing

### Earth leakage (AC Line)

Perform the following tests, with the applicable switch settings on the “Safety Tester” and record the results in the applicable test form(s).

- UUT power switch in the “ON” position.


Test Description	S1: DPDT Polarity a) Normal b) Reverse		S2: SPDT Neutral a) Closed b) Open		S3: SPDT Ground a) Closed b) Open		S4: SPDT Test point a) Internal b) TP1		Maximum Limits  (uA)		
	Position		Position		Position		Position				
	a	b	a	b	a	b	a	b	EN/IEC	UL	
Earth leakage current measurements											
TP1 – no connection											
Polarity: Normal Neutral: Normal (Closed) Ground: Open	X		X			X	X		500	300	
Polarity: Normal Neutral: Open Ground: Open	X			X		X	X		1,000	1,000	
Polarity: Reverse Neutral: Open Ground: Open		X		X		X	X		1,000	1,000	
Polarity: Reverse Neutral: Normal (Closed) Ground: Open		X	X			X	X		500	300	



### Enclosure (Chassis) leakage

Perform the following tests, with the applicable switch settings on the "Safety Tester" and record the results in the applicable test form(s).

- UUT power switch in the "ON" position.
- TP1 connected to accessible metal parts (see tables below for individual test points).

Test Description	S1: DPDT		S2: SPDT		S3: SPDT		S4: SPDT		Maximum Limits	
	Polarity		Neutral		Ground		Test point			
	a) Normal b) Reverse		a) Closed b) Open		a) Closed b) Open		a) Internal b) TP1			
	Position		Position		Position		Position			
	a	b	a	b	a	b	a	b	EN/IE C	UL
										
Enclosure (Chassis) leakage current measurements: TP1 connected to “RS-232 DB9 shell and Metal part of Printer.										
Polarity: Normal										
Neutral: Normal (Closed)	X		X		X			X	100	100
Ground: Normal (Closed)										
Polarity: Normal										
Neutral: Open	X			X	X			X	500	300
Ground: Normal (Closed)										
Polarity: Normal										
Neutral: Normal (Closed)	X		X			X		X	500	300
Ground: Open										
Polarity: Reverse										
Neutral: Normal (Closed)		X	X			X		X	500	300
Ground: Open										
Polarity: Reverse										
Neutral: Open		X		X	X			X	500	300
Ground: Normal (Closed)										
Polarity: Reverse										
Neutral: Normal (Closed)		X	X		X			X	100	100
Ground: Normal (Closed)										

## Section 5: Performance Verification and Safety Testing


### Patient Leakage Current to Ground

Perform the following tests, with the applicable switch settings on the “Safety Tester” and record the results in the applicable test form(s).


- UUT power switch in the “ON” position.
- TP1 connected to “Applied Parts” (see tables below for individual test points).

**Note:** The following applied parts have a separate isolated input and have to be tested separately.

- Leads
- Paddles
- SpO<sub>2</sub>

<div>Test Description</div> <div></div>	S1: DPDT		S2: SPDT		S3: SPDT		S4: SPDT		Maximum Limits  (uA)	
	Polarity		Neutral		Ground		Test point			
	a) Normal b) Reverse		a) Closed b) Open		a) Closed b) Open		a) Internal b) TP1			
	Position		Position		Position		Position			
a	b	a	b	a	b	a	b	EN/IEC	UL	
Patient leakage current measurements: LEADS										
TP1 connected to “LEADS” (all LEADS connected together), Type CF defibrillation proof limits										
Polarity: Normal										
Neutral: Normal (Closed)	X		X		X			X	10	10
Ground: Normal (Closed)										
Polarity: Normal										
Neutral: Open	X			X	X			X	50	50
Ground: Normal (Closed)										
Polarity: Normal										
Neutral: Normal (Closed)	X		X			X		X	50	50
Ground: Open										
Polarity: Reverse										
Neutral: Normal (Closed)		X	X			X		X	50	50
Ground: Open										
Polarity: Reverse										
Neutral: Open		X		X	X			X	50	50
Ground: Normal (Closed)										

## Section 5: Performance Verification and Safety Testing

Test Description	S1: DPDT Polarity a) Normal b) Reverse		S2: SPDT Neutral a) Closed b) Open		S3: SPDT Ground a) Closed b) Open		S4: SPDT Test point a) Internal b) TP1		Maximum Limits  (uA)	
	Position		Position		Position		Position			
	a	b	a	b	a	b	a	b	EN/IEC	UL
										
Polarity: Reverse Neutral: Normal (Closed) Ground: Normal (Closed)		X	X		X			X	10	10
Patient leakage current measurements: PADDLES										
TP1 connected to "PADDLES" (all PADDLES connections tied together), Type CF defibrillation proof limits										
Polarity: Normal Neutral: Normal (Closed) Ground: Normal (Closed)	X		X		X			X	10	10
Polarity: Normal Neutral: Open Ground: Normal (Closed)	X			X	X			X	50	50
Polarity: Normal Neutral: Normal (Closed) Ground: Open	X		X			X		X	50	50
Polarity: Reverse Neutral: Normal (Closed) Ground: Open		X	X			X		X	50	50
Polarity: Reverse Neutral: Open Ground: Normal (Closed)		X		X	X			X	50	50
Polarity: Reverse Neutral: Normal (Closed) Ground: Normal (Closed)		X	X		X			X	10	10

## Section 5: Performance Verification and Safety Testing

Test Description	S1: DPDT		S2: SPDT		S3: SPDT		S4: SPDT		Maximum Limits  (uA)	
	Polarity		Neutral		Ground		Test point			
	a) Normal b) Reverse		a) Closed b) Open		a) Closed b) Open		a) Internal b) TP1			
	Position		Position		Position		Position			
	a	b	a	b	a	b	a	b	EN/IEC	UL
<b>Patient leakage current measurements: SpO<sub>2</sub></b>										
<b>TP1 connected to “SpO<sub>2</sub>” (all SpO<sub>2</sub> connections tied together), Type BF defibrillation proof limits</b>										
Polarity: Normal Neutral: Normal (Closed) Ground: Normal (Closed)	X		X		X			X	100	100
Polarity: Normal Neutral: Open Ground: Normal (Closed)	X			X	X			X	500	500
Polarity: Normal Neutral: Normal (Closed) Ground: Open	X		X			X		X	500	500
Polarity: Reverse Neutral: Normal (Closed) Ground: Open		X	X			X		X	500	500
Polarity: Reverse Neutral: Open Ground: Normal (Closed)		X		X	X			X	500	500
Polarity: Reverse Neutral: Normal (Closed) Ground: Normal (Closed)		X	X		X			X	100	100

**Patient Leakage Current, Mains on Applied Part (All SIP/SOPs Grounded)**

Perform the following tests, with the applicable switch settings on the "Safety Tester" and record the results in the applicable test form(s).

1. Disconnect the "Safety Tester" from the "Mains Supply"
2. Ensure the battery is removed
3. Configure the test set up as shown (or equivalent)
  - UUT hot/neutral/ground connected to ground
  - All accessible metal parts (i.e. DB9 shell and exposed metal of printer) connected to ground
4. UUT power switch in the "ON" position.

**Caution:**


Use care when connecting the "Safety Tester" to the "Mains Supply", mains voltage will appear on the applied parts.

**Caution:**


Use care when selecting switch settings (S2 should be open, S3 = safety ground)

Note: The following applied parts have a separate isolated input and have to be tested separately.

- Leads
- Paddles
- SpO<sub>2</sub>



## Section 5: Performance Verification and Safety Testing

Test Description	S1: DPDT		S2: SPDT		S3: SPDT		S4: SPDT		Maximum Limits	
	Polarity		Neutral		Ground		Test point			
	a) Normal b) Reverse		a) Closed b) Open		a) Closed b) Open		a) Internal b) TP1			
	Position		Position		Position		Position			
	a	b	a	b	a	b	a	b	EN/IEC	UL

Patient leakage current, Mains on Applied Part (All SIP/SOPs Grounded): LEADS										
MD connected to "LEADS" (all LEADS connected together), Type CF defibrillation proof limits										
Caution: mains voltage on "LEADS"										
Polarity: Normal										
Neutral: Open	X			X	X			X	50	50
Ground: Normal (Closed)										
Polarity: Reverse										
Neutral: Open		X		X	X			X	50	50
Ground: Normal (Closed)										

Patient leakage current, Mains on Applied Part (All SIP/SOPs Grounded): PADDLES										
MD connected to "PADDLES" (all PADDLES connections tied together), Type CF defibrillation proof limits										
Caution: mains voltage on "PADDLES"										
Polarity: Normal										
Neutral: Open	X			X	X			X	100 (60601-2-4)	50
Ground: Normal (Closed)										
Polarity: Reverse										
Neutral: Open		X		X	X			X	100 (60601-2-4)	50
Ground: Normal (Closed)										

Patient leakage current, Mains on Applied Part (All SIP/SOPs Grounded): SpO <sub>2</sub>										
MD connected to "SpO <sub>2</sub> " (all SpO <sub>2</sub> connections tied together), Type BF defibrillation proof limits										
Caution: mains voltage on "SpO <sub>2</sub> "										
Polarity: Normal										
Neutral: Open	X			X	X			X	5,000	5,000
Ground: Normal (Closed)										
Polarity: Reverse										
Neutral: Open		X		X	X			X	5,000	5,000
Ground: Normal (Closed)										

### Insulation Resistance

Perform the following tests, with the applicable resistance range settings and record the results in the applicable test form(s).


1. Disconnect from the SUPPLY MAINS
2. Remove battery
3. Remove I/O (RS-232) cable
4. UUT power switch in the "ON" position.

The resistance values shall not be less than:

- 2 M $\Omega$  between MAINS PART and all other parts separated by basic insulation, including TYPE B APPLIED PARTS
- 7 M $\Omega$  between MAINS PART and all other parts separated by double insulation, including TYPE B APPLIED PARTS
- 70 M $\Omega$  between MAINS PART to any F-TYPE APPLIED PART


### Accessible Metal Parts Insulation Resistance

- The resistance value shall not be less than: 7 M $\Omega$



### MAINS PART to "F-TYPE APPLIED PARTS" insulation resistance

- The resistance value shall not be less than: 70 M $\Omega$



## SECTION 6: Parts and Accessories

To order service parts for the responder 2000, contact:

GE Healthcare Technologies Service Logistic Center  
79111 Freiburg  
Munzinger Straße 3-5  
Germany


### Overview






This section lists parts and accessories for the Responder 2000.





TOPIC	PAGE
Kits	89
Accessories	95
Power Cords	97






## Kits


<i>Kit Part Number</i>	<i>Kit Description</i>	<i>Part Photograph</i>
<b>2029882-001</b>	<b>Therapy Board Kit</b>	
Therapy Board with Cables	1. <b>Therapy Board</b>	
<b>2029882-002</b>	2. Left Drain Cable	
Therapy Board Cables Only	3. Right Drain Cable	
	4. <b>Therapy Board</b> to Power Control Board Cable	
	5. <b>Therapy Board</b> to Main Board Ribbon Cable	
<b>2029882-003</b>	<b>Main Board Kit</b>	
Main Board basic with Cables (For models 2025653-xxx)	1. Main Board	
<b>2029882-004</b>	2. Main Board to Printer Ribbon Cable	
Main Board pacer with cables (For models 2026109-xxx)	3. Power Control Board to Main Board Ribbon Cable	
<b>2029882-005</b>	4. Power Control Board to Main Board Ribbon Cable	
Main Board Pacer/SpO <sub>2</sub> with cables (For Models 2026114-xxx)	5. Power Control Board to Main Board Cable	
<b>2029882-006</b>		
Main Board Cables Only		

<i>Kit Part Number</i>	<i>Kit Description</i>	<i>Part Photograph</i>
<b>2029882-008</b>	<b>Power Control Board Kit</b>	
Power Control Board and Cables	<ol style="list-style-type: none"> <li>1. Power Control Board</li> <li>2. Service Connector (with Cable)</li> </ol>	
<b>2029882-009</b>		
Wire Set for Power Control Board	<ol style="list-style-type: none"> <li>3. Power Control Board Main Board Cable</li> <li>4. Power Control Board to Therapy Board Cable</li> <li>5. Power Control Board to Main Board Ribbon Cable</li> </ol>	
<b>2029882-007</b>	<b>Power Supply Kit</b>	
Power Supply Kit	Power Supply	
<b>2029882-010</b>	<b>ECG Board Kit</b>	
<b>ECG Board</b> with Isolation Shield	<ol style="list-style-type: none"> <li>1. ECG Board Insulator</li> <li>2. ECG Board EMF Cover</li> <li>3. <b>ECG Board</b></li> </ol>	
<b>2029882-011</b>	<b>SpO<sub>2</sub> Board Kit</b>	
SpO <sub>2</sub> Board Kit	SpO <sub>2</sub> Board	



<i>Kit Part Number</i>	<i>Kit Description</i>	<i>Part Photograph</i>
<b>2029882-012</b> High Voltage Capacitor Set	<b>Main Capacitors Kit</b> Main Capacitors (2)	
<b>2029882-013</b> Rear Cover with Defib Connector	<b>Rear Cover Kit</b> 1. Rear Cover (Including Paddle Cables and Paddle Wires) 2. Handle 3. Bed Rail Hook Covers 4. Rear Cover Label	
<b>2029882-016</b> Device Handle Only		
<b>2029882-017</b> Bed Rail Hook Covers Only		
<b>2029882-015</b> ECG Input Connector, wired	<b>ECG Cable Kit</b> ECG Connector and Cable	
<b>2029882-014</b> ECG/SpO <sub>2</sub> Input Connector wired	<b>ECG/SpO<sub>2</sub> Cable Kit</b> ECG/SpO <sub>2</sub> Connectors and Cables	

<i>Kit Part Number</i>	<i>Kit Description</i>	<i>Part Photograph</i>
<b>2029882-018</b>	<b>Display Kit</b>	
LCD Set with Retainer	1. Display Retainer	
<b>2029882-019</b>	2. Copper Contact Spring	
LCD Retainer only	3. Display Assembly (including Backlight cable)	
	4. Display to Main Board ribbon cable	
<b>2029882-020</b>	<b>Printer Kit</b>	
Printer Assembly	1. Printer Assembly	
	2. Printer Roller (Not Shown)	
<b>2029882-021</b>	<b>Printer Door Kit</b>	
Printer Door Assembly	Printer Door	
	<b>Note:</b> The printer door does not come with the Roller. The roller must be ordered as part of the Printer Kit.	
<b>2029882-025</b>	<b>Base Cover Kit</b>	
Base Cover Kit	1. Extension Spring	
<b>2029882-022</b>	2. Battery Latch	
Battery Release Kit (Includes Spring, Latch, and Release Only)	3. Battery Release	
	4. Base Cover	
	<b>Note:</b> This kit does not include a Front Label.	



<i>Kit Part Number</i>	<i>Kit Description</i>	<i>Part Photograph</i>
<b>2029882-023</b>  Front Bezel Kit	<b>Front Bezel Kit</b>  Front Bezel	
<b>2029882-024</b>  Front Body Kit	<b>Front Body Kit</b>  1. Paddle Latches (with Springs) (2) 2. Paddle Latches (2) 3. Front Body 4. Speaker (not shown)  <b>Note:</b> Speaker is integral to Front Body and not available separately.	
<b>2029882-026</b>  Wear and Tear Kit	<b>Wear and Tear Kit</b>  1. Power Button Assembly 2. Charge, Manual, Shock Button Assembly 3. Fan 4. Encoder Cap 5. Encoder Button and Cable Assembly 6. Front Label	
<b>2029882-027</b>  Hardware Starter Kit	Includes supply of screws,	
<b>2029882-028</b>  Tube of Silicone Sealant	washers, silicone tubing, silicone sealant, rubber grommets, cable ties, Torx driver, screw drivers, and a communications cable.	

<i>Kit Part Number</i>	<i>Kit Description</i>	<i>Part Photograph</i>
<b>2036667-001</b>  Service Stickers	<b>Service Stickers</b>  Includes a roll of 500 tamper proof service stickers	

## Accessories


<i>Accessory Part Number</i>	<i>Part Description</i>	<i>Part Photograph</i>
<b>OXY-RWL</b> Datex-Ohmeda Foam Wrap	<b>Datex-Ohmeda Replacement Large Foam Wrap</b>	
<b>OXY-RTW</b> Datex-Ohmeda Wide Adhesive Tape	<b>Datex-Ohmeda Replacement Wide Adhesive Tape</b>	
<b>2030137-001</b> Paddle Pair, Responder	<b>Paddles with two adult surface plates</b>	
<b>2030134-001</b> Adapter Electrode Paddle Responder	<b>Responder Paddle Electrode Adaptors</b>	
<b>2030247-001</b> Cable Interconnection Responder Electrode/Pad	<b>Responder Electrode/Pad Cable Interconnection</b>	
<b>2030249-001</b> Electrode Pair Responder Internal	<b>Responder Internal Electrode Pair</b>	



<i>Accessory Part Number</i>	<i>Part Description</i>	<i>Part Photograph</i>
<b>2026327-001</b> Thermal Paper	<b>60-mm Thermal Paper (Box of 50)</b>	
<b>2025267-001</b> Rechargeable Battery	<b>Rechargeable Battery</b>	



## Power Cords

<b>Europe</b> <b>2020387-002</b>	<b>Argentina</b> <b>2020387-003</b>	<b>Denmark</b> <b>2020387-004</b>
		
<b>India and South Africa</b> <b>2020387-005</b>	<b>Italy</b> <b>2020387-006</b>	<b>Japan</b> <b>2020387-007</b>
		
<b>UK</b> <b>2020387-009</b>	<b>Israel</b> <b>2020387-010</b>	<b>Australia</b> <b>2020387-011</b>
		

<i>China</i> <i>2020387-012</i>	<i>Brazil</i> <i>2020387-013</i>	
		

## SECTION 7: Theory of Operation

### Overview

This section provides descriptions of the Responder 2000 and its components. These instructions are intended for use only by service providers who are specifically trained to service the Responder 2000.

TOPIC	PAGE
System Overview	100
System Interconnection Block Diagram	101
Component Descriptions	102

## System Overview

The Responder 2000 system has the following functions and features.

**Note:** Refer to the Operator's Manual for specific instructions and intended use.

### Functions

Defibrillation (with dedicated Defibrillator operation buttons)

- 2J to 270J Biphasic waveform (can deliver 270J Biphasic shock into a 50 Ohm load)
- Paddle connector accepts pads, paddles, or spoons

Transcutaneous cardiac pacing

40ms rectilinear current source pace pulse, 0 to 140 mA, 0 to 180 PPM, demand and asynchronous modes

ECG Monitoring through therapy electrodes, paddles, or a separate five wire ECG input

Optional Pulse Oxymetry (GE506)

### Features

Color LCD to display operational, ECG, and SPO<sub>2</sub> information

The display is 5.67", backlit, with a 320 x 240 resolution and over 260k colors

Trim knob for input via the display

Thermal printer for ECG, Pulse Oxymetry, Therapy Actions and Event history

- 384 pixel resolution at 25mm/s on a 60mm paper roll
- Prints in real-time two-trace graphic, text, and raster modes

AC or rechargeable battery operation


- AC: 90 VAC to 264 VAC, 47Hz-63Hz, 90 Watts Cont., 115 Watts intermittent
- Battery: Lithium Ion with built in charge indicator, delivers up to 10A at 11.0V (Nominal)

Internal Defibrillator and Pacer test load

Speaker for audible indication

Serial port for software upgrade

## System Interconnection Block Diagram



## Component Descriptions

### Main Board

The Main Board controls all major subsystems and user interface.

Controls defibrillation (J210), pacing, ECG (J212), and SpO<sub>2</sub> (J221) functions

Controls the display interface and provides 3.3V (J214) and high voltage backlight power (J213) for LCD.

Controls (J219) and provides 5V (J220) for the Printer

Controls the speaker (J216)

Software upgrade path through the serial port (pass-through from the Power Control Board)

### Therapy Board

The **Therapy Board** delivers energy for defibrillation and pacing. It is also a pass-through to the main board for controls, ID features on defibrillator paddles, and ECG signals. Specifics include:

3 kV electrical isolation to separate the patient from the unit and earth ground

Delivers up to 270J biphasic defibrillator shock (J101, J102, & J109)

Delivers up to 140mA 40ms pace pulses

Charges two 100uF Capacitors to 2 kV (J105 & J107 negative; J104 & J106 positive)

Connects to the Main Board (J110)

Receives power from the Power Control Board (J103)

### ECG Board

The optional **ECG Board** receives input from ECG leads connected. The **ECG Board** then amplifies and digitizes the data and passes it on to the Main Board

4 kV electrical isolation to separate the patient from the unit and earth ground

Connects to the Main board (J412)

ECG Leads connect at J23

### Power Control Board

The Power Control Board provides power to the unit (J303, J311), charges the battery (J325), and provides UI pass-through (J318). Specifics include:

Receives power from the Power Supply (J315)

Provides RS-232 service connection (J317)

Buttons connected to the Power Control board include: Power (on/off), Charge, Shock, Manual shock, and Trim knob (J324)

LEDs include: AC Power, Shock, Manual, Charge, Power, Fault, Battery Charge

Powers the Fan (J323).

### SpO<sub>2</sub> Board

The optional SpO<sub>2</sub> board (P221) mounts to the Main board and patient SpO<sub>2</sub> sensor connection is directly to the SpO<sub>2</sub> module (J11).

## SECTION 8: Specifications and Safety

### Overview

This section presents the specifications and safety standards of the Responder 2000.

<b>TOPIC</b>	<b>PAGE</b>
Specifications	104
Physical Dimensions	106
Environmental Requirements	106
RHYTHMx <sup>®</sup> ECG Analysis Algorithm	106
STAR <sup>®</sup> BIPHASIC Defibrillation Waveform	108
Electromagnetic Compatibility Requirements	111
Environmental Standards	111
Electromagnetic Emissions Table	112
Electromagnetic Immunity Table	114
RF Communications Table	115
Connectors	118

## Specifications

### Display

Size:	115.2mm X 86.4mm
Type:	Transmissive Color TFT LCD
Resolution:	320 X 240 pixels
Number of waveform channels:	2

### Defibrillator

Waveform:	Biphasic truncated exponential
Charge Time:	7 seconds nominal
Delivery Method:	Via pads, paddles or spoons
Paddles and Pads Energy Selections (50Ω load):	2, 3, 5, 7, 10, 15, 20, 30, 50, 70, 100, 150, 200, 270
Internal Paddle Energy Selections (50Ω load):	2, 3, 5, 7, 10, 15, 20, 30, 50
Energy Accuracy:	± 15% onto a 50Ω load or ± 3 Joules, whichever is greater

### Non-Invasive Pacing

Output Waveform:	Rectilinear, constant current
Delivery Method	Via Pads
Pulse Width:	40 milliseconds ±4mS
Pacing Modes:	Demand or Asynchronous (Fixed Rate)
Pacing Rate:	Operator adjustable: 30 to 180 ppm, ± 5%
Pacing Current:	0mA to 140mA

### SpO<sub>2</sub>

Display	Plethysmogram, digital value of percent saturation, and upper and lower alarm limits.
Display Update Period	8 seconds or less
Saturation Range	1 to 100%, in 1% increments
Low Saturation Alarm	81-98% in 1% increments
High Saturation Alarm	Off, 95-100% in 1% increments
Pulse Rate Range	20 to 255 BPM, in 1 BPM increments
Audible Alarm Delay	10 seconds or less
Visual Alarm Delay	2 seconds or less
Pulse Rate Accuracy	30 to 250 bpm: ± 2 digits or ± 2%, whichever is greater (without motion) 30 to 250 bpm: ± 3 digits (during low perfusion)
Saturation Accuracy	+/- 2 digits from 70% SpO <sub>2</sub> to 100% SpO <sub>2</sub> with D-O probes, except +/- 3 digits for D-O Ear Probe. OEM board accuracy +/-3 digits from 70%

### ECG Monitoring



**SECTION 8: Specifications and Safety**

Detection Lead:	Through either 3-lead and 5-lead patient cable, pads, or paddles
Lead Selection:	3-Lead Cable: (I, II, III) and 5-Lead (I, II, III, aVL, aVR, aVF, and V), PADS: (Modified Lead II)
Heart Rate Range:	25 to 300 bpm
Display Frequency Response (ECG Leads):	0.5 to 100 Hz ( $\pm 10\%$ )
Display Frequency Response (paddles):	3 Hz to 33 Hz ( $\pm 10\%$ ); (-3db)
ECG Effective Sampling Rate:	1000 samples per second (ECG Leads), 128 samples per second (paddles)
Post Defibrillation Recovery:	8 seconds
Heart Rate Alarms:	Off, 5-300 bpm
Leads Off Sensing	DC current; Sensing Leads $<0.1\text{ }\mu\text{A}$ ; Reference Lead $<1\text{ }\mu\text{A}$ ; when all leads are connected
ECG Leads and paddles leakage current	Less than 10 $\mu\text{A}$ in normal; Less than 50 $\mu\text{A}$ in single fault condition
Dynamic Range: Input ECG signal amplitude:	$\pm 5\text{ mV}$
Dynamic Range: DC Offset voltage:	$\pm 500\text{ mV}$ ECG from ECG cable; $\pm 1000\text{ mV}$ (ECG from pads, paddles and spoons)
Asystole Threshold	0.2mV ( $\pm 0.1\text{ mV}$ )

**Rechargeable Battery**

Battery Voltage:	11.1 V Nominal
Chemistry:	Lithium-ion
Compatibility:	Compatible with Responder 2000
Battery Capacity:	50 shocks, or 240 minutes monitoring time, or 72 minutes monitoring time with pacing
Battery Charge Time:	8 hours in Responder 2000, 4 hours in external charger, 20 hours for calibration cycle in external charger
Battery Standby:	6 months
Battery Life:	2.5 years or 300 Battery charge-discharge cycles, whichever comes first
Battery Weight:	1 lb. 3 oz.; .54 kg

**AC Power Supply**

Universal input:	100V to 240V~ 50Hz-60Hz 200VA
------------------	-------------------------------------

**Printer**

Speed	25 mm/S feed rate
Paper Size	60 mm paper width

## Physical Dimensions

Height: 10.8 inches/27.4 cm  
Width: 11.7 inches/29.7 cm  
Depth: 7.4 inches/18.8 cm  
Weight: Less than 10lbs/4.5kg, excluding battery, paddles, and a full roll of paper.

## Environmental Requirements

### Operating Conditions

Temperature: 0°C to 50°C (32 °F to 122°F)  
Humidity: 10% to 95% RH, non condensing  
Air Pressure altitude: -500 ft (103kPa) to 15000 ft (57kPa)

### Storage and Shipping Conditions

Temperature (Responder 2000): -20°C to 60°C (-4°F to 140°F)  
Temperature (Pads): -12°C to 43°C (10°F to 110°F) for 2 years.

**NOTE:** Do not exceed 38°C (100°F) for periods greater than 6 months in duration. Always store the pads in their pouch to maintain freshness. Do not use if gel has dried out.

Humidity: 10% - 95% RH, non-condensing  
Air Pressure altitude: -500 ft (103kPa) to 15000 ft (57kPa)

#### **WARNING:**

Electrode performance may be adversely affected by pre-attaching and storing with defibrillator cable or exposure to air for long periods of time. These electrodes are not recommended for electrosurgery.

#### **CAUTION: Environment of use**

Responder 2000 is designed for indoor use. Operator must confirm that the environment of use meets the required operating environmental specifications before using.

#### **CAUTION: Cold Environments**

If the Responder 2000 is stored in an environment with a temperature below the operating temperature, the unit should be allowed to warm up to the needed operating temperature before using.

## RHYTHMx® ECG Analysis Algorithm

The RHYTHMx® ECG analysis algorithm provides superior ECG detection capabilities. The features available with the Responder 2000 include the following:

- Detection Rate
- Asystole Threshold
- Noise Detection
- Non-Committed Shock
- Synchronized Shock
- SVT Discriminators
- Supraventricular Tachycardia (SVT) Rate
- Continuous Monitoring

The Responder 2000 rejects all T-waves that are 1 millivolt or less in the conditions specified in ANSI/AAMI EC 13 section 4.1.2.1c.

The Responder 2000 will alarm tachycardia in the conditions specified in ANSI/AAMI EC 13 section 4.1.2.1 g) in less than 10 seconds.

For the alternating ECG complexes specified in ANSI/AAMI EC 13 figure 3, the Responder 2000 will indicate the following heart rates:

Figure 3a 40 bpm

Figure 3b 52 bpm

Figure 3c 59 bpm

Figure 3d 122 bpm

(Refer to the ANSI/AAMI EC 13 for the figures.)

The following sections describe specific RHYTHMx® information.

### Detection Rate

All ventricular fibrillation (VF) and ventricular tachycardia (VT) rhythms at or above this rate will be classified as shockable. All rhythms below this rate will be classified as non-shockable. This rate is configurable between 120 bpm (beats per minute) and 240 bpm. The default Detection Rate for the Responder 2000 is 160 bpm.

### Fine VF

Fine VF is classified by the signal amplitude-less than 0.2 mV peak-to-peak for eight (8) consecutive seconds, preceded by a shockable arrhythmia or the peak-to-peak amplitude-less than 0.9 mV, the amplitude distribution indicator is less than amplitude distribution threshold, the derivative probability density function is satisfied, the RR interval index is not regular, and it is preceded by a shockable rhythm. Fine VF is a shockable arrhythmia.

### Asystole

The signal amplitude is less than 0.2 mV peak-to-peak for 8 consecutive seconds and is not preceded by a shockable rhythm; the rhythm will be classified as Asystole. Asystole is not shockable.

### Noise Detection

The Responder 2000 will detect noise artifact in the ECG. Noise could be introduced by excessive moving of the patient or electronic noise from external sources like cellular and radiotelephones.

### Non-Committed Shock

After the Responder 2000 advises a shock, it continues to monitor the patient ECG rhythm. If the patient's rhythm changes to a non-shockable rhythm before the shock button is pressed, the shock will be cancelled.

## Sync Mode

The Responder 2000 is designed to deliver synchronized shock on the R-wave for Sync Shock. The Responder 2000 will automatically attempt to synchronize a shock to the R-wave. If delivery cannot be synchronized within two seconds it will not deliver the shock. It is recommended to select/adjust to a lead with a unipolar R wave of about 1 mV peak amplitude for safe and reliable synchronization.

IEC 60601-2-4 (2002) and ANSI/AAMI DF 80 (2003), clause 104c states that the maximum delay from the "peak of the QRS" to the peak of the defibrillator output waveform shall be 60ms. Verification testing has shown that the Responder 2000 meets this requirement of the standards.

Some ECG leads may exhibit a bipolar QRS waveform complex, and in these cases, the RHYTHMx software in the Responder 2000 will pick the highest peak of the QRS complex for synchronization. The peak of the defibrillator output will occur in less than 60ms from this peak, and therefore the Responder 2000 meets the synchronization requirement of the standards.

In the case of these leads, certain defibrillator testers may use a different part of the QRS complex (for example the first, smaller peak) to measure synchronization time of defibrillation. This may give a measurement result that exceeds 60ms. Investigation has shown that the Responder 2000 does synchronize appropriately from the peak of the QRS complex, and meets the requirement of the standards in cases where a particular defibrillator tester does not measure the defibrillation delay from the peak of the QRS complex.

## SVT (Supraventricular Tachycardia) Discriminators

The Responder 2000 is supplied with the SVT Discriminator enabled and the default setting is 240 bpm. SVT Discriminators are sophisticated filters that analyze the morphology of the ECG waveforms and distinguish VF/VT from SVT and Normal Sinus Rhythms (NSR). The SVT Discriminator will only be applied to rhythms that fall between the Detection Rate and the SVT Rate.

## SVT Rate

All rhythms with rates between the Detection Rate and SVT Rate will be screened through a number of SVT Discriminators to classify them into VF/VT or SVT. Rhythms classified as SVT between the two set rates are not shockable. All rhythms above the rates will be classified as shockable. The SVT Rate must be greater than the Detection Rate and is selectable between 140 and 240 bpm. The default SVT rate is 240.

## Continuous Monitoring For Shockable Rhythm

The Responder 2000 can monitor the ECG rhythms continuously.

## Pacemaker Pulse Information

Pacemaker pulses without overshoot, in all of the conditions specified in ANSI/AAMI EC13 section 4.1.4.1, in the range of 20mV to 700mV and 0.1 milliseconds to 2 milliseconds wide, will be rejected by the Responder 2000.

Pacemaker pulses with overshoot, in all the conditions specified in ANSI/AAMI EC 13 section 4.1.4.2, in the range of 20mV to 700mV and 0.1 milliseconds to 2 milliseconds wide, will be rejected by the Responder 2000.

The pacer pulse detector will not respond to the waveform of ANSI/AAMI figure 5d, since this waveform is below the threshold of the Responder 2000 pacer pulse detector. The minimum typical slew rate in V/s RTI that will trip the pacer detector is 6.2 V/s for the 3 and 5 lead ECG. The minimum typical slew rate in V/s RTI that will trip the pacer detector is 9.8 V/s for paddles.

**STAR® BIPHASIC Defibrillation Waveform**

The waveform generated by the Responder 2000 is a BIPHASIC TRUNCATED EXPONENTIAL waveform.

STAR® Biphasic Waveform – 270J into Pads or paddles


Table A - 270J Waveform into Different Resistive Loads (Typical Values)							
Patient's Impedance (Ohms)	Phase 1 Delivered Start Volts	Phase 1 Delivered End Volts	Phase 1 Duration (ms)	Phase 2 Delivered Start Volts	Phase 2 Delivered End Volts	Phase 2 Duration (ms)	Total Energy Delivered (J)
25	1692V	990V	3.25 ms	990V	342V	3.2ms	281J
50	1860V	1234V	4.50 ms	1234V	684V	3.2ms	270J
75	1923V	1338V	5.75ms	1338V	887V	3.2ms	254J
100	1957V	1394V	7.00 ms	1394V	1015V	3.2ms	241J
125	1977V	1429V	8.25ms	1429V	1103V	3.2ms	231J
150	1991V	1453V	9.50 ms	1453V	1166V	3.2ms	223J
175	2002V	1469V	10.75ms	1469V	1214V	3.2ms	217J
180	2003V	1472V	11.00ms	1472V	1222V	3.2ms	216J
200	2009V	1520V	11.00ms	1520V	1283V	3.2ms	202J

STAR Biphasic Waveform – 50J into internal spoon

Table B – 50J Waveform into Different Resistive Loads (Typical Values)							
Patient's Impedance (Ohms)	Phase 1 Delivered Start Volts	Phase 1 Delivered End Volts	Phase 1 Duration (ms)	Phase 2 Delivered Start Volts	Phase 2 Delivered End Volts	Phase 2 Duration (ms)	Total Energy Delivered (J)
10	573V	217V	3.00ms	217V	27V	3.2ms	47J
20	697V	386V	3.00ms	386V	108V	3.2ms	51J
25	728V	426V	3.25 ms	426V	147V	3.2ms	52J
50	800V	531V	4.50 ms	531V	294V	3.2ms	50J
75	828V	576V	5.75ms	576V	382V	3.2ms	47J
100	842V	600V	7.00 ms	600V	437V	3.2ms	45J

## Energy Levels and Patient Impedance

The Biphasic Truncated Exponential (BTE) waveform delivers energy that is variant with the patient impedance. The waveform is designed to deliver the selected energy when the patient impedance is 50 Ohms, as shown in the above waveform table.



## Safety Standards and Compliance Requirements

The Responder 2000 is designed to meet all applicable requirements of the standards listed below.

IEC 60601-1, (1988 + A1:1991 + A2:1995), Medical Electrical Equipment Part 1 General Requirements for Safety

EN 60601-1, (1990 + A1:1993 + A2:1995), 2<sup>nd</sup> Edition Medical Electrical Equipment, Part 1: General Requirements for Safety

IEC 60601-1-1, (2000), Medical Electrical Equipment - Part 1: General Requirements for Safety 1: Collateral Standard: Safety Requirements for Medical Electrical Systems

IEC 60601-2-4, (2002), Medical Electrical Equipment – Part 2-4: Particular Requirements for the Safety of Cardiac Defibrillators

IEC 60601-2-49, (2001), Medical Electrical Equipment - Part 2-49: Particular Requirements For The Safety Of Multifunction Patient Monitoring Equipment

IEC 60601-2-27, (1994), Medical electrical equipment, part 2: Particular requirements for the safety of electrocardiographic monitoring equipment

UL 60601-1, (2003), Medical Electrical Equipment Part 1, General Requirements for safety

CAN/CSA-C22.2 No. 601.1-M90, Medical electrical equipment Part 1: General Requirements for Safety

ANSI/AAMI EC-13-2002, Cardiac monitors, heart rate meters, and alarms

ANSI/AAMI DF80-2003, Medical Electrical Equipment – Part 2-4: Particular Requirements for the Safety of Cardiac Defibrillators (including automated external defibrillators)

## Electromagnetic Compatibility Requirements

The Responder 2000 meets the requirements of the following EMC standards, as required by IEC 60601-2-4.:

IEC 60601-1-2 (2001), Medical electrical equipment Part 1: General requirements for safety 2. Collateral standard: electromagnetic compatibility - Requirements and tests.

### Emissions

Electromagnetic Fields: CISPR 11 (2003), Industrial, scientific and medical (ISM) radio-frequency equipment - radio disturbance characteristics - limits and methods of measurement; Group 1, Class B. IEC 60601-2-4 (2002), Section 36.201.1.

Harmonic distortion: IEC 61000-3-2 (2004), Electromagnetic Compatibility (EMC) Part 3-2: Limits - Limits For Harmonic Current Emissions (Equipment Input Current Less Than Or Equal To 16 A Per Phase).

Voltage fluctuations and flicker: IEC 61000-3-3 (2002), Electromagnetic Compatibility (EMC) - Part 3-3: Limits - Limitation Of Voltage Changes, Voltage Fluctuations And Flicker In public Low-Voltage Supply Systems, For Equipment With Rated Current Less Than Or Equal To 16 A Per Phase.

### Immunity

Electromagnetic: IEC 61000-4-3 (2003), Electromagnetic compatibility (EMC) - part 4-3: Testing and measurement techniques - radiated, radio-frequency, electromagnetic field immunity test; Level 3 (10V/m) and X (20V/m). IEC 60601-2-4 (2002) Section 36.202.3.

Magnetic: IEC 61000-4-8 (1994), Electromagnetic compatibility (EMC) - part 4. Testing and measurement techniques - section 8. Power frequency magnetic field immunity test basic EMC publication; Level X (3 A/m). IEC 60601-2-4 (2002), Section 36.202.8.

ESD: IEC 61000-4-2 (2001), Electromagnetic compatibility (EMC) - part 4-2: testing and measurement techniques - electrostatic discharge immunity test; Level 3. IEC 60601-2-4 (2002), Section 36.202.2.

Conducted: IEC 61000-4-6 (2003), Electromagnetic compatibility (EMC) - part 4-6: testing and measurement techniques - immunity to conducted disturbances, induced by radio-frequency fields. IEC 60601-2-4 (2002), Section 36.202.6.

## **SECTION 8: Specifications and Safety**

Fast transients and bursts: IEC 61000-4-4 (2001), Electromagnetic compatibility (EMC) - part 4: Testing and measurement techniques - section 4: Electrical fast transient/burst immunity test. IEC 60601-2-4 (2002), Section 36.202.4.

Surges: IEC 61000-4-5 (2001), Electromagnetic compatibility (EMC) - part 4: Testing and measurement techniques - section 5: Surge immunity test. IEC 60601-2-4 (2002), Section 36.202.5.

Voltage dips, short interruptions and voltage variations on power supply input lines: IEC 60601-4-11 (2004), Electromagnetic Compatibility (EMC) - Part 4-11: Testing And Measurement Techniques - Voltage Dips, Short Interruptions And Voltage Variations Immunity Tests.

## **Environmental Standards**

### **Shock and Vibration**

The Responder 2000 is tested per the following when in the unpackaged condition:

Bump: IEC 60068-2-29 (1987), Test EB: bump; 25g, 6 ms, 0.9 m/s  $\Delta V$ , and 1000 bumps in each direction

Sine Vibration: IEC 60068-2-6 (1995), Environmental testing - part 2: tests - test FC: Vibration (sinusoidal); 0.15mm displacement amplitude, 10-55Hz, 10 sweep cycles in each axis

Random Vibration: IEC 60068-2-64 (1993), Environmental testing - part 2: test methods - test FH: Vibration broadband random (digital control) and guidance: 1-100Hz, 0.01g<sup>2</sup>/Hz 30 minutes.

Free Fall Drop: IEC 60068-2-32 (1975 + A1:1990), Environmental testing - test methods - test ED: free fall; 18 inches

Enclosure Protection: IEC 60529 (2003), Degrees of protection provided by enclosures (IP code); IP22.

### **Storage and Shipping**

When packaged in the shipping container, the Responder 2000 meets the requirements of:

ISTA Preshipment Test 2A (2001), Simulation Performance Test Procedure - Packaged-Products 150lb(68 kg) or Less




## Electromagnetic Emissions Table

Guidance and manufacturer's declaration – electromagnetic emissions		
The Responder 2000 is intended for use in the electromagnetic environment specified below. The customer or the user of the Responder 2000 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Responder 2000 uses RF energy only for its internal function. Therefore its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.  The Responder 2000 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

## Electromagnetic Immunity Table

Guidance and manufacturer's declaration – electromagnetic immunity			
The Responder 2000 is intended for use in the electromagnetic environment specified below. The customer or the user of the Responder 2000 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD)  IEC 61000-4-2	±6 kV contact  ±8 kV air	±6 kV contact  ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst  IEC 61000-4-4	±2 kV for power supply lines  ±1 kV for input/output lines	±2 kV for power supply lines  ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge  IEC 61000-4-5	±1 kV differential mode  ±2 kV common mode	±1 kV differential mode  ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines  61000-4-11	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle  40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles  70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles  <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle  40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles  70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles  <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Responder 2000 requires continued operation during power mains interruptions, it is recommended that the Responder 2000 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field  IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE $U_T$ is the a.c. mains voltage prior to application of the test level.			

## SECTION 8: Specifications and Safety

Guidance and manufacturer's declaration – electromagnetic immunity			
The Responder 2000 is intended for use in the electromagnetic environment specified below. The customer or the user of the Responder 2000 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF  IEC 61000-4-6	3 Vrms  150 kHz to 80 MHz outside ISM bands <sup>a</sup>	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the Responder 2000, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.  <b>Recommended separation distance</b>  $d = 1.2 \sqrt{P}$
	10 Vrms  150 kHz to 80 MHz in ISM bands <sup>a</sup>	10 Vrms	
Radiated RF  IEC 61000-4-3	10 V/m  80 MHz to 2.5 GHz	10 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz  $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz  where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m) <sup>b</sup> .  Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>c</sup> should be less than the compliance level in each frequency range. <sup>d</sup>  Interference may occur in the vicinity of equipment marked with the following symbol:  
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

## SECTION 8: Specifications and Safety

- <sup>a</sup> The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 to 40.70 MHz.
- <sup>b</sup> The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in these frequency ranges.
- <sup>c</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the RESPONDER 2000 is used exceeds the applicable RF compliance level above, the RESPONDER 2000 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the RESPONDER 2000.
- <sup>d</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 1 V/m.

## RF Communications Table

Recommended separation distances between portable and mobile RF communications equipment and the Responder 2000				
The RESPONDER 2000 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the RESPONDER 2000 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the RESPONDER 2000 as recommended below, according to the maximum output power of the communications equipment.				
Rated maximum output power of transmitter  W	Separation distance according to frequency of transmitter m			
	150 kHz to 80 MHz outside ISM bands  $d = 1.2\sqrt{P}$	150 kHz to 80 MHz in ISM bands  $d = 1.2\sqrt{P}$	80 MHz to 800 MHz  $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz  $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.12	0.23
0.1	0.38	0.38	0.38	0.73
1	1.2	1.2	1.2	2.3
10	3.8	3.8	3.8	7.3
100	12	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.				
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.				
NOTE 2 The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 to 40.70 MHz.				
NOTE 3 An additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.				
NOTE 4 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

## Connectors

The following sections describe connectors in the subassemblies within Responder 2000 Defibrillator.


### Case

The case serves both as a housing and mechanical structure for the Responder 2000.

#### Power Connector – IEC 320 type – AC Power

Connector: Heyco 0916		
Pin Number	Name	Description
1	AC1	AC1 – AC2: 90-264VAC, 47-63Hz, <10A. AC1 or AC2 to GND: <264VAC.
2	AC2	
3	GND	Earth ground connection

#### Paddles/Pads Connector – GE Defibrillator Paddles



View looking at socket on the back of the Responder 2000

Connector: GE 43252483		
Pin Number	Name	Description
PC11	SW_CHG/SHCK	Pull Down Input – Charge/Shock Button
PC12	SW_SHCK	Pull Down Input – Shock Button
PC13	N/A	Not Used
PC14	ID_PADS/SPNS	Pull Down Input – ID Bit – Pads

Connector: GE 43252483		
Pin Number	Name	Description
PC15	RET_PC	Return for Buttons and ID bits
PC16	ID_PDL/SPNS	Pull Down Input – ID Bit – Paddles/Spoons
APEX	SHOCK_P	Patient Apex
STERNUM	SHOCK_M	Patient Sternum

**Paddle Connector ID Codes**

Cable P/N	Description	PC14	PC16
22338301	Pacer or Defib Pads	Short to PC15	
22334601	Flat Pads	Short to PC15	
21730801	Internal Spoons	4.7K to PC15	Short to PC15
30344625	Responder Paddles		Short to PC15

**ECG Connector – GE ECG Connector – Five lead ECG**

Connector: GE 401760-1		
Pin Number	Pin Name	Description
1	RA	Input – Right Arm
2	NC	No Connection
3	NC	No Connection
4	LA	Input – Left Arm
5	RL	Input – Right Leg
6	V	Input – V
7	NC	No Connection

**SECTION 8: Specifications and Safety**

8	LL	Input – Left Leg
9	NC	No Connection
10	NC	No Connection
11	SHLD	Shield

**ECG Cable Identification Encoding**

ECG Cables available for use with the Cardiac Science Responder 2000 have identification features to enable identification of the cable type. Cable connections are shown here:

Cable Connections for 5 lead, 3 lead and sync cables.				
ECG Connector		Cable Connections		
Pin Number	Pin Name	5 Lead Cable	3 Lead Cable	Sync Cable
1	RA	RA	RA	- Attenuated Sync Input
2	NC			
3	NC			
4	LA	LA	LA	Center-tap Attenuated Sync Input
5	RL	RL	130M resistor Pin 5 to Pin 6	33M resistor Pin 5 to Pin 6
6	V	V		
7	NC			
8	LL	LL	LL	+ Attenuated Sync Input
9	NC			
10	NC			
11	SHLD	Shield	Shield	



SpO<sub>2</sub> Connector – GE SpO<sub>2</sub>

Connector: GE 401762		
Pin Number	Name	Description
1	SPO2_RCAL	Input – Calibration Data
2	SPO2_LED-	Output – Led Neg
3	SPO2_LED+	Output – Led Pos
5	SPO2_ANODE	Input – Detector Pos
6	GNDF	Isolated Floating Ground
7	SPO2_RCAL_GND	Ground – Calibration Data
9	SPO2_CATHODE	Input – Detector Neg
4, 8, 10, 11	NC	No Connection

## Serial Service Connector – DB9-Female

RS232 - DTE. Hardware handshaking not supported.

Connector: Amp 747905-2		
Pin Number	Name	Description
2	RS232_RxD	Input – RS232 Data
3	RS232_TxD	Output – RS232 Data
5	GND	Ground
1, 4, 6-9	NC	No Connection

## Paddle Cradles (Test Load)

Pin Number	Name	Description
RIGHT	Apex	A 50Ohm internal defibrillator test load is connected across the paddle cradles.
LEFT	Sternum	

**MAIN CPU PCBA**

The main CPU PCBA is the nucleus of the system, with all major subsystems connecting to this board.

**J211 - Power Control Interface to Power Control PCBA**

<b>AMP 1-103638-3</b>		
<b>Pin Number</b>	<b>Name</b>	<b>Description</b>
1	NC	Not Used
2	Fan_/Drive	Pull Down Active Fan Drive
3	Serv_DCE_Rx	RS232 Serial input from Service Connector.
4	Serv_DCE_Tx	RS232 Serial Output from Service Connector.
5	SMBUS_CLK	Output – SMBus Serial Clock
6	SMBUS_DATA	Bi-directional – SMBus Data
7	/PWR_ON	Output – Active low control to power switch on the Power Control PCBA.
8,9,10	+12V_SW	Power Input – Power switched on Power Control PCBA
11,12,13	GND	System Ground
14	+12V_UNSW	Power Input – Power always present from Power Control PCBA

## J214 – LCD Panel

Hirose FH12-40-S-0.5SH		
Pin Number	Name	Description
1-4, 6, 31,36	NC	No Connection
5, 10, 14, 18, 22, 26, 30, 32, 34	GND	System Ground
7	LCD_R0	Output - Red Bit 0
8	LCD_R1	Output - Red Bit 1
9	LCD_R2	Output - Red Bit 2
11	LCD_R3	Output - Red Bit 3
12	LCD_R4	Output - Red Bit 4
13	LCD_R5	Output - Red Bit 5
15	LCD_G0	Output – Green Bit 0
16	LCD_G1	Output – Green Bit 1
17	LCD_G2	Output – Green Bit 2
19	LCD_G3	Output – Green Bit 3
20	LCD_G4	Output – Green Bit 4
21	LCD_G5	Output – Green Bit 5
23	LCD_B0	Output - Blue Bit 0
24	LCD_B1	Output - Blue Bit 1
25	LCD_B2	Output - Blue Bit 2
27	LCD_B3	Output - Blue Bit 3
28	LCD_B4	Output - Blue Bit 4
29	LCD_B5	Output - Blue Bit 5
33	LCD_FPSHFT	Output - LCD Pixel Shift

Hirose FH12-40-S-0.5SH		
Pin Number	Name	Description
35	LCD_DRDY	Output - LCD Frame Shift
37-40	+3.3V	Power Output - LDC Panel Power

**J213 – LCD CCFL Backlight**

JST SM92B-BHSS-1-TB		
Pin Number	Name	Description
1	BACKLIGHT1	Power Output – CCFL High Voltage
2	BACKLIGHT2	Power Output – CCFL High Voltage

**J218 – UI Interface to Power Control PCBA**

(Has UI Functions)

<b>Molex 71226-2425</b>		
<b>Pin Number</b>	<b>Name</b>	<b>Description</b>
1-2	GND	System Ground
3	/Manual_SW	Pull Down Input Pushbutton – Manual
4	/Spare1_SW	Pull Down Input Pushbutton – Spare1
5	/Shock_SW	Pull Down Input Pushbutton – Shock
6	/Charge_SW	Pull Down Input Pushbutton – Charge
7	CPU_ADC5	Analog Input from UI/PC PCBA Temperature Sensor
8	/Select_SW	Pull Down Input Push-knob – Trim knob Select
9	SQRB	Digital Input – Trim knob phase 2
10	SQRA	Digital Input – Trim knob phase 1
11	/Power_SW	Pull Down Input Pushbutton – Power
12	Power_LED	LED Drive – Pull-Down Output Power LED
13	Spare2_LED	LED Drive – Pull-Down Output Spare2 LED
14	Spare1_LED	LED Drive – Pull-Down Output Spare1 LED
15	Charge_LED	LED Drive – Pull-Down Output Charge LED

Molex 71226-2425		
Pin Number	Name	Description
16	Manual_LED	LED Drive – Pull-Down Output Manual LED
17	Shock_LED	LED Drive – Pull-Down Output Shock LED
18	Fault_LED	LED Drive – Pull-Up Output Fault LED
19	CHG_IND	Digital Input – Battery Charging
20	AC_IND	Digital Input – AC Present
21	Buzzer	Pull Down Output – 5V Buzzer
22, 23	+5V	Power Output
24	UI_/Present	Pull Down Input – UI PCBA Connected

**J216 – SPEAKER**

<b>AMP 103638-1</b>		
<b>Pin Number</b>	<b>Name</b>	<b>Description</b>
1	SPEAKER1	Output – Audio
2	SPEAKER2	Output – Audio

**J219 – PRINTER HEAD**

<b>Molex 39-51-3263</b>		
<b>Pin Number</b>	<b>Name</b>	<b>Description</b>
1, 2, 25, 26	+5V_PWR	Power Output – Print head high current supply.
24	+5V_DIG	Power Output – Print head digital supply.
3	CLK	Output – Pixel clock
4	STROBE	Output – Shift Register Latch
5	DIN	Output – Serial Pixel Data
7	OE_0-63	Output – Pixel Burn – Pixels 0-63
8	OE_64-127	Output – Pixel Burn – Pixel 64-127
11	OE_128-191	Output – Pixel Burn – Pixel 128-191
13	THRMSTR	Input – Thermistor (to Ground)
15	OE_192-255	Output – Pixel Burn – Pixel 192-255
17	OE_256-319	Output – Pixel Burn – Pixel 256-319
18	OPTO_COL	Input – Open Collector Active Low paper sense.
20	PAPER_LED	Output – LED Drive – Paper sensor.
22	OE_320-383	Output – Pixel Burn – Pixel 320-383
23	DOUT	Input – Return from end of shift register.

Molex 39-51-3263		
Pin Number	Name	Description
6, 9, 10, 12, 14, 16, 19, 21	GND	System Ground

**J220 – PRINTER MOTOR**

JST S6B-PH-K-S		
Pin Number	Name	Description
1	COIL1A	Output – Motor Coil 1 – Pos
2	COIL1B	Output – Motor Coil 1 – Neg
3	COIL2A	Output – Motor Coil 2 – Pos
4	COIL2B	Output – Motor Coil 2 – Neg
5	DOORSW	Pull Down Input – Door Switch
6	GND	System Ground

**J221 – SpO<sub>2</sub> INTERFACE**

SAMTEC ESW-107-44-S-D		
Pin Number	Name	Description
8	RxD	Output - Isolated Logic Level Serial
9	TxD	Input - Isolated Logic Level Serial
10, 13	+5V	Isolated Power Output
2, 3, 5, 11, 12, 14	GND	Isolated Ground
1, 4, 6, 7	NC	No Connection



**J210 – THERAPY PCBA**

<b>MOLEX 52610-2090</b>		
<b>Pin Number</b>	<b>Name</b>	<b>Description</b>
1	THRP_/PRSNT	Pull Down Input – Therapy PCBA plugged in
2	THRP_HVEN	Output – HV enable
3	THRP_/SHOCK	Output – Shock enable (active low)
4	THRP_ISOPWRON	Output – Therapy isolated power enable
5	THRP_Rx	Input – Logic Level Serial
6	THRP_Tx	Output – Logic Level Serial
7	+3.3V	Power Output
8	GND	System Ground
9	PDL_ID_PC5	Pull Down Input – Paddle ID bit
10	PDL_ID_PC4	Pull Down Input – Paddle ID bit
11	PDL_CHGSHCK	Pull Down Input – Button – Charge/Shock
12	PDL_SHCK	Pull Down Input – Button – Shock
13	PDL_ID_PC6	Pull Down Input – Paddle ID bit
14, 15	PDL_GND	System Ground
16-19	NC	No Connection
20	THP_PRSNT_GND	Grounded output to Therapy – Looped back on THRP_/PRSNT

**J212 – ECG PCBA**

<b>SAMTEC DW-04-11-F-D-500</b>		
<b>Pin Number</b>	<b>Name</b>	<b>Description</b>
2	TxD	Output - Logic Level Serial
4	RxD	Input - Logic Level Serial

SAMTEC DW-04-11-F-D-500		
Pin Number	Name	Description
6	/ECG_EN	Output – Active Low ECG power enable
7, 8	+12V	Power Output
1, 3, 5	GND	System Ground

## THERAPY PCBA

The therapy PCBA provides defibrillation, pacing, and ECG amplification and digitization from the therapy electrodes. All operations are performed per serial communications with the main CPU PCBA. The **Therapy Board** also serves as a pass through board for controls and ID features on defibrillator paddles.

### J103 - Power Supply

Connects to power management PCBA

Molex 39-30-2047		
Pin Number	Name	Description
2, 3	+12V_SW	Power Input – Switched on Power Management PCBA.
1, 4	GND	System Ground

### J110 – Main CPU PCBA interface

Molex 52610-2090		
Pin Number	Name	Description
1	THRP_/PRSNT	Pull Down Input – Therapy PCBA plugged in
2	THRP_HVEN	Output – HV enable
3	THRP_/SHOCK	Output – Shock enable (active low)
4	THRP_ISOPWRON	Output – Therapy isolated power enable
5	THRP_Rx	Input – Logic Level Serial
6	THRP_Tx	Output – Logic Level Serial
7	+3.3V	Power Output

Molex 52610-2090		
Pin Number	Name	Description
8	GND	System Ground
9	PDL_ID_PC5	Pull Down Input – Paddle ID bit
10	PDL_ID_PC4	Pull Down Input – Paddle ID bit
11	PDL_CHGSHCK	Pull Down Input – Button – Charge/Shock
12	PDL_SHCK	Pull Down Input – Button – Shock
13	PDL_ID_PC6	Pull Down Input – Paddle ID bit
14, 15	PDL_GND	System Ground
16-19	NC	No Connection
20	THP_PRSNT_GND	Grounded output to Therapy – Looped back on THRP_/PRSNT

**J109 – Paddles Control Connector**

AMP 103670-7		
Pin Number	Name	Description
8	PDL_CHGSHCK	Pull Down Input – Button – Charge/Shock
7	PDL_SHCK	Pull Down Input – Button – Shock
6	PDL_ID_PC6	Pull Down Input – Paddle ID bit
5	PDL_ID_PC4	Pull Down Input – Paddle ID bit
3	PDL_ID_PC5	Pull Down Input – Paddle ID bit
1, 2, 4	PDL_GND	Paddle Control Ground – Isolated from Patient

**J101, J102, J108 – Paddles HV Connectors**

Keystone 1287-ST		
Pin Number	Name	Description
J101	PDLHV_STERN	Patient Connection – Sternum

**SECTION 8: Specifications and Safety**

J102	PDLHV_APEX	Patient Connection – Apex
------	------------	---------------------------

**J104, J105, J106, J107 – Energy Storage Capacitor Connectors**

Keystone 1287-ST		
Pin Number	Name	Description
J104	HVCAP1_POS	Energy Storage Capacitor 1 – Pos
J105	HVCAP1_NEG	Energy Storage Capacitor 1 – Neg
J106	HVCAP2_POS	Energy Storage Capacitor 2 – Pos
J107	HVCAP2_NEG	Energy Storage Capacitor 2 – Neg

**J126, J127 – Internal Test Load Resistor**

Keystone 1287-ST		
Pin Number	Name	Description
J126	TESTLOAD1	Connects to a paddle cradle
J127	TESTLOAD2	Connects to a paddle cradle

**ECG PCBA**

The ECG PCBA is a monitoring ECG subsystem. All operations are performed per serial communications with the main CPU PCBA.

**J412 – Interface to main CPU PCBA**

Power and serial communications are carried through this connector.

SAMTEC BCS-104-L-D-PE-BE		
Pin Number	Name	Description
1	GND	System Ground
2	SER_CPU2ECG	Logic Level Serial Data Input
3	GND	System Ground
4	SER_ECG2CPU	Logic Level Serial Data Input
5	GND	System Ground

6	/ECG_EN	ECG Isolated power enable
7	+12V	System +12V (nominal)
8	+12V	System +12V (nominal)

**J423 – Patient ECG Connection**

Keystone 4903		
Pin Number	Name	Description
J4231	RA	Electrode - Right Arm
J4232	LA	Electrode - Left Arm
J4233	LL	Electrode - Left Leg
J4234	V1	Electrode - V1
J4235	RL	Electrode - Right Leg
J4236	SHLD	Cable Shield

**POWER CONTROL PCBA**

Responsible for battery charging and DC power switching and distribution. Also includes UI functions

**J315 – Power input from DC power supply**

Molex 39-30-3058		
Pin Number	Name	Description
1	GND	System Ground
2	GND	System Ground
3	GND	System Ground
4	+12V_MAIN	Power Input
5	+12V_MAIN	Power Input

**J325 – Smart Battery Connector – Data and Power**

Suyin 250137MR010G101ZU
-------------------------

**SECTION 8: Specifications and Safety**

Pin Number	Name	Description
1, 2, 3, 4	GND	System Ground
5	BAT_TEMP	Input – 10K to Ground
6	SMBUS_DATA	Bi-directional – SMBus Data
7	SMBUS_CLK	Bi-directional – SMBus Serial Clock
8, 9, 10	+12V_BAT	Battery Positive

**J303 – Therapy power**

Molex 39-30-3047		
Pin Number	Name	Description
1	GND	System Ground
2	+12V_SW	Power Output
3	+12V_SW	Power Output
4	GND	System Ground

**J311 – Power Control interface from main CPU PCBA**

AMP 1-103673-3		
Pin Number	Name	Description
1	+12V_UNSW	Power Output – Power always present.
2, 3, 4	GND	System Ground
5, 6, 7	+12V_SW	Power Output – Power switched on Power Control PCBA
8	/PWR_ON	Input – Active low control to power switch on the Power Control PCBA.
9	SMB_DAT	Bi-directional – SMBus Data
10	SMB_CLK	Input – SMBus Serial Clock
11	SERV_TX	Service Port Tx RS232 Data From Main PCBA

12	SERV_RX	Service Port Rx RS232 Data To Main PCBA
13	/FAN_DRV	Fan Drive – Pulled Down ON Passes through Power Control PCBA to J232
14	N/C	

**J318 – UI Features - Interface to main CPU PCBA**

Molex 52207-2490		
Pin Number	Name	Description
1	UI_/Present	Grounded Output (UI PCBA Connected)
2, 3	+5V	Power Input
4	Buzzer	Pull Down Input – 5V Buzzer
5	AC_IND	Digital Output – AC Present
6	CHG_IND	Digital Output – Battery Charging
7	Fault_LED	LED Drive – Fault LED Pull-Down Input from Main PCBA
8	Shock_LED	LED Drive – Shock LED Pull-Down Input from Main PCBA
9	Manual_LED	LED Drive – Manual LED Pull-Down Input from Main PCBA
10	Charge_LED	LED Drive – Charge LED Pull-Down Input from Main PCBA
11	Reserved	Reserved for: LED Drive – Spare1 LED
12	Reserved	Reserved for: LED Drive – Spare2 LED
13	Power_LED	LED Drive – Power LED Pull-Down Input from Main PCBA

**SECTION 8: Specifications and Safety**

<b>Molex 52207-2490</b>		
<b>Pin Number</b>	<b>Name</b>	<b>Description</b>
14	/Power_SW	Pull Down Output to Main PCBA Pushbutton – Power
15	SQRA	Digital Output to Main PCBA Trim knob phase 1
16	SQRB	Digital Output to Main PCBA Trim knob phase 2
17	/Select_SW	Pull Down Output to Main PCBA Push-knob – Trim knob Select
18	CPU_ADC5	Analog Output to Main PCBA Temperature Sensor
19	/Charge_SW	Pull Down Output to Main PCBA Pushbutton – Charge
20	/Shock_SW	Pull Down Output to Main PCBA Pushbutton – Shock
21	Reserved	Reserved for: Pushbutton – Spare1
22	/Manual_SW	Pull Down Output to Main PCBA Pushbutton – Manual
23, 24	GND	System Ground

**J317 – RS-232 SERIAL SERVICE INTERFACE**

<b>AMP 103673-2</b>		
<b>Pin Number</b>	<b>Name</b>	<b>Description</b>
1	Tx	RS232 Output
2	Rx	RS232 Input
3	GND	System / Serial Ground



## J323 – FAN

AMP 103673-1		
Pin Number	Name	Description
1	+12V	Power Output
2	/FAN_OUT	Pull Down Output – Fan

## J324 – ENCODER “Trim Knob”

AMP 0-215079-6		
Pin Number	Name	Description
1,2	GND	System Ground
3	/Select_SW	Pull-down switch input
4	SQRB	Digital Input from encoder Trim knob phase 2
5	SQRA	Digital Input from encoder Trim knob phase 1
6	+5VDC	Power Output

SpO<sub>2</sub> PCBA – GE506 Module mounts on the main CPU PCBA with interfacing header connectors plugged directly together with no cable.

J2 - SpO<sub>2</sub> HOST INTERFACE

Samtec FTSH-107-01-L-D-RA		
Pin Number	Name	Description
8	RxD	Output – Logic level serial
9	TxD	Input – Logic level serial
11	CTS	Input – Tied low on Main CPU PCBA
10, 13	+5V	80mA
2, 3, 5, 12, 14	GND	Power and Signal Ground Return
1, 4, 6, 7	NC	

SpO<sub>2</sub> PATIENT INTERFACE

Samtec FTSH-105-01-L-D-RA		
Pin Number	Name	Description
1	DET+	Input – Detector Pos
2	DIGICAL_GND	Ground – Calibration Data
3	DIGICAL_DAT	Input – Calibration Data
4	DET-	Input – Detector Neg
5	SHLD_INNER	Inner Shield
6	N/C	No Connection
7	LED-	Output – Led Neg
8	N/C	No Connection
9	LED+	Output – Led Pos
10	SHLD_OUTER	Outer Shield

## POWER SUPPLY MODULE

Provides DC power for the instrument via the Power Control PCBA.

## TB1 – POWER INPUT

Heyco 0916 – Wire connected - Mounts on back of instrument		
Pin Number	Name	Description
1	L	AC Input 90-264VAC 47-63Hz
2	G	
3	N	

## POWER OUTPUT

Wire Connected – Plugs into Power Control PCBA		
Pin Number	Name	Description
1	GND	NOMINALLY 13.5V

2	GND	POWER OUTPUT IS FLOATING WITH RESPECT TO:  EARTH GROUND  POWER INPUTS
3	GND	
4	+12V	
5	+12V	

**BATTERY**

Provides DC power for the instrument via the Power Control PCBA while the instrument is not plugged into an AC power source.

**J325 – Smart Battery Connector – Data and Power**

Suyin 250132FB010G200ZU		
Pin Number	Name	Description
1, 2, 3, 4	GND	Battery Ground
5	BAT_TEMP	Output – 10K to Ground
6	SMBUS_DATA	Bi-directional – SMBus Data
7	SMBUS_CLK	Bi-directional – SMBus Serial Clock
8, 9, 10	+12V_BAT	Battery Positive

## Contact Information/Customer Service

To order supplies or accessories, contact your representative or distributor. For technical support, contact your local GE customer service.

Please have the serial and model numbers available. The serial and model numbers are located on the back of the Responder 2000.

Responder 2000 is manufactured for:

GE Medical Systems Information Technologies, Inc.  
8200 West Tower Avenue, Milwaukee, WI 53223 USA

**Tel.: 800 558 5120 (USA only)**

Fax: 800 421 6841

**Canada Tel: 800 668 0732**

GE Medical Systems Information Technologies GmbH  
Munzinger Str. 3, D-79111 Freiburg, Germany

**Tel.: +49 761 4543 0**

Fax: +49 761 4543 233

Responder 2000 is manufactured by:



Cardiac Science Corporation  
3303 Monte Villa Parkway  
Bothell, WA 98021



MDSS GmbH  
Schiffgraben 41  
D-30175 Hannover  
Germany

---

Responder is a trademark owned by GE Medical Systems Information Technologies, a General Electric Company, going to market as GE Healthcare. STAR, IntelliSense, RHYTHMx, are trademarks and registered trademarks of Cardiac Science Corporation. All other trademarks are property of their respective owners.

© 2012 General Electric Company. All rights reserved.

END OF DOCUMENT – LAST PAGE



Responder<sup>®</sup> 2000

---

defibrillator/monitor