GE Medical Systems

Technical Publications

Direction 2217718-100

Revision 0

MX150 Tube Replacement Hoist Procedure

Copyright© 1998 By General Electric Co.

Operating Documentation

WARNING

- THIS SERVICE MANUAL IS AVAILABLE IN ENGLISH ONLY.
- IF A CUSTOMER'S SERVICE PROVIDER REQUIRES A LANGUAGE OTHER THAN ENGLISH, IT IS THE CUSTOMER'S RESPONSIBILITY TO PROVIDE TRANSLATION SERVICES.
- DO NOT ATTEMPT TO SERVICE THE EQUIPMENT UNLESS THIS SERVICE MANUAL HAS BEEN CONSULTED AND IS UNDERSTOOD.
- FAILURE TO HEED THIS WARNING MAY RESULT IN INJURY TO THE SERVICE PROVIDER, OPERATOR OR PATIENT FROM ELECTRIC SHOCK, MECHANICAL OR OTHER HAZARDS.

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Direction 2217718-100

Revision 0

MX150 Tube Replacement Hoist Procedure

IMPORTANT! . . . X-RAY PROTECTION

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. The General Electric Company, Medical Systems Group, will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide

complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that everyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission

on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that the General Electric Company, Medical Systems Group, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment.

Various protective material and devices are available. It is urged that such materials or devices be used.

CAUTION: United States Federal law restricts this device to use by or on the order of a physician.

THIS PAGE INTENTIONALLY LEFT BLANK.

If you have any comments, suggestions or corrections to the information in this document, please write them down, include the document title and document number, and send them to:

GENERAL ELECTRIC COMPANY MEDICAL SYSTEMS
 MANAGER – INFORMATION INTEGRATION,
 AMERICAS W-622
 P.O. BOX 414
 MILWAUKEE, WI 53201-0414

CERTIFIED ELECTRICAL CONTRACTOR STATEMENT

All electrical installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations, and testing

shall be performed by qualified GE Medical personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE's electrical work on these products will comply with the

requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE's field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation "**damage in shipment**" written on **all** copies of the freight or express bill **before** delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier **immediately**

upon discovery, or in any event, within **14** days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this **14** day period.

Call Traffic and Transportation, Milwaukee, WI (414) 827-3449 /

8*285-3449 **immediately** after damage is found. At this time be ready to supply name of carrier, delivery date, consignee name, freight or express bill number, item damaged and extent of damage.

Complete instructions regarding claim procedure are found in Section "S" of the Policy & Procedure Bulletins.

THIS PAGE INTENTIONALLY LEFT BLANK.

TABLE OF CONTENTS

SECTION	TITLE	PAGE
	REVISION HISTORY	ix
1	GENERAL	1
2	PROCEDURE FOR TC/TCV SYSTEMS	2
3	PROCEDURE FOR LCA/LCV+/LC+ SYSTEMS	3

TABLE OF CONTENTS (cont.)

REVISION HISTORY

REV	DATE	REASON FOR CHANGE
A	May 29, 1998	Preliminary release.
0	July 13, 1998	Initial release.

LIST OF EFFECTIVE PAGES

PAGE NUMBER	REVISION NUMBER	PAGE NUMBER	REVISION NUMBER	PAGE NUMBER	REVISION NUMBER
Title Page	0				
i thru x	0				
1 thru 8	0				
Back Page	-				

THIS PAGE INTENTIONALLY LEFT BLANK.

**SECTION 1
GENERAL**

A tube lifting crane is provided inside the MX150 X-ray tube crate. The procedures described here enable the tube to be mechanically lifted into place on the TC/TCV systems (Section 2), and on the LCA/LCV +/LC+ systems (Section 3). One person can perform these procedures in most instances.

SECTION 2
PROCEDURE FOR TC/TCV SYSTEMS

1. Remove the outer package and the top of the wooden crate.
2. Take out the 4 blocks: block A, block B, and the two blocks D. (See Illustration 1.)
3. Place the blocks D on the floor, positioned such that when the side door of the crate is lowered, it rests on the blocks. (See Illustration 1.)
4. Put the block B (the smaller of the two remaining blocks) in a convenient spot on the floor, near to the gantry but away from the working area and without blocking any access routes.
5. Position the cart and lift the new tube from the crate and place it on the block B. (See Illustrations 2 and 3.)
6. Remove the TC table top using either a safe lifting method or the procedure in the Quick Reference Guide.
7. Position the C arm at 0 degrees, SID to maximum, and the C arm long movement to about the center of travel.
8. Switch off the system, and perform Lock out/Tag out.
9. Disconnect the tube and make it ready to be removed, as explained in the appropriate positioner manual.
10. Remove the collimator and the collimator plate, remove the mini-tank, and the anode HT cable, disconnect the 2 stator cables, the 2 water pipes and remove the 2 fixing bolts from the tube mount and the strengthening plate that prevents the tube being lifted. The collimator can normally be rested on the C arm without having to disconnect it. (See Illustration 6.)
11. Install the lifting handle on the defective tube. (The system was originally supplied with a handle that should be available for this purpose.)
12. Position the crane (See Illustration 6) and lift the tube as high as possible. Withdraw the tube, you may have to tilt the tube in order for the HT cable port to clear the gantry. Tilt safely, at the end of the open part of the crane's hook. Take care that there is no risk of pushing the tube out of the hook.
13. Take the defective tube and place it in the tube crate; the oil pump end should be towards the bottom of the crate.
14. Use the crane to lift the replacement tube off block B and install it in the TC gantry. You may have to tilt the tube in order for the HT cable port to clear the gantry. Tilt safely, at the end of the open part of the crane's hook. Take care that there is no risk of pushing the tube out of the hook.
15. Reconnect the replacement tube and continue with the calibration process. A Quick Reference Guide and a more detailed procedure for the calibration process is available.
16. Repack the defective tube, the crane and the 4 blocks in the crate for return. A packing guide is contained in the tube documentation. Illustration 7 shows the positions of the securing straps.

SECTION 3
PROCEDURE FOR LCA/LCV+/LC+ SYSTEMS

1. Remove the outer package and the top of the wooden crate.
2. Take out the 4 blocks: block A, block B, and the two blocks D. (See Illustration 1.)
3. Place the blocks D on the floor, positioned such that when the side door of the crate is lowered, it rests on the blocks. (See Illustration 1.)
4. Put the block B (the smaller of the two remaining blocks) in a convenient spot on the floor, near to the gantry but away from the working area and without blocking any access routes.
5. Position the cart and lift the new tube from the crate and place it on the block B. (See Illustrations 2 and 3.)
6. Position the L and C to zero degrees. Move the Pivot to between 43 and 45 degrees, until the block A (the remaining block with the wooden base) fits just underneath the tube and almost touching it. (See Illustration 4.)
7. Switch off the system, and perform Lock out/Tag out.
8. Disconnect the tube and make it ready to be removed, as explained in the appropriate positioner manual.
9. Remove the collimator and the collimator plate, remove the mini-tank, and the anode HT cable, disconnect the 2 stator cables, the 2 water pipes and remove the 2 fixing bolts from the tube mount. The collimator can normally be rested on the L arm base without having to disconnect it.
10. Install the lifting handle on the defective tube. (The system was originally supplied with a handle that should be available for this purpose.)
11. Using a safe lifting position, unhook the defective tube from the positioner and lower it such that it rests on block A. (See Illustration 5.)
12. Use the crane to lift the defective tube off block A and place it in the tube crate; the oil pump end should be towards the bottom of the crate.
13. Use the crane to lift the replacement tube from block B and place it on block A. Use the handle and a safe lifting position to hook the replacement tube onto the gantry.
14. Reconnect the replacement tube and continue with the calibration process. A Quick Reference Guide and a more detailed procedure for the calibration process is available.
15. Repack the defective tube, the crane and the 4 blocks in the crate for return. A packing guide is contained in the tube documentation. Illustration 7 shows the positions of the securing straps.

ILLUSTRATION 1
MX 150 X-RAY TUBE CRATE


ILLUSTRATION 2
LIFTING NEW TUBE FROM CRATE


ILLUSTRATION 3
PLACING NEW TUBE ON BLOCK B

ILLUSTRATION 4
LCA/LCV+/LC+
POSITIONING L, C, AND PIVOT





ILLUSTRATION 5
LCA/LCV+/LC+
RESTING DEFECTIVE TUBE ON BLOCK A

**ILLUSTRATION 6
TC/TCV
LIFTING DEFECTIVE TUBE**

**ILLUSTRATION 7
REPACKING DEFECTIVE TUBE**

THIS PAGE INTENTIONALLY LEFT BLANK.

GE Medical Systems

*GE Medical Systems: Telex 3797371
P.O. Box 414, Milwaukee, Wisconsin 53201 U.S.A.
(Asia, Pacific, Latin America, North America)*

*GE Medical Systems — Europe: Telex 698626
283, rue de la Minière, B.P. 34, 78533 Buc Cedex
France*