Instructions for use
Instrucciones de uso
使用説明書

elements™ 8:1

Kerr[™]
ENDODONTICS

Manufactured for:

Kerr Corporation
1717 West Collins Avenue
Orange CA 92867
U.S.
+1-800-KERR-123
www.kerrdental.com

**Manufacturer:**

KaVo Dental GmbH
Bismarckring 39
88400 Biberach
Germany
www.kavo.com

PN: 077-0649
06/2021
Rev. A

Distributed by:

Kerr Australia Pty. Limited
Unit 6, 12 Mars Road
Lane Cove West, New South Wales 2066
Australia
+61-2-8870-3000

**Importer**

Ormco B.V.
Basicweg 20 NL-3821 BR
Amersfoort, The Netherlands

Overview

Instructions for use	3
Instrucciones de uso	91
使用説明書	185

Instructions for use elements™ 8:1

Manufactured for:

Kerr Corporation
1717 West Collins Avenue
Orange CA 92867
U.S.
+1-800-KERR-123
www.kerrdental.com

**Manufacturer:**

KaVo Dental GmbH
Bismarckring 39
88400 Biberach
Germany
www.kavo.com

PN: 077-0649
06/2021
Rev. A

Distributed by:

Kerr Australia Pty. Limited
Unit 6, 12 Mars Road
Lane Cove West, New South Wales 2066
Australia
+61-2-8870-3000

**Importer**

Ormco B.V.
Basicweg 20 NL-3821 BR
Amersfoort, The Netherlands

Table of contents

1	User instructions	8
2	Safety.....	16
2.1	Infection hazard	16
2.2	Improper use	17
2.3	Technical condition	19
2.4	Accessories and combination with other equipment	21
2.5	Qualification of personnel	22
2.6	Service and repair	23
3	Product description	25
3.1	Purpose – Intended use	26
3.2	Technical Specifications.....	27
3.3	Symbols on product and rating plate	29
3.4	Transportation and storage conditions	32
4	Startup and shut down.....	34
5	Operation.....	37

5.1	Attaching the medical device	37
5.2	Removing the medical device	39
5.3	Inserting the dental bur or diamond grinder	40
5.4	Removing the milling tool or diamond grinder	46
6	Checking for malfunctions and troubleshooting.....	48
6.1	Check for malfunctions.....	48
6.2	Troubleshooting	50
6.2.1	Replacing the O-rings on the motor coupling.....	52
7	Reprocessing steps in accordance with ISO 17664.....	53
7.1	Preparations at the site of use.....	53
7.2	Manual Reprocessing	55
7.2.1	Manual external cleaning.....	56
7.2.2	Manual internal cleaning	56
7.2.3	Manual external disinfection	58
7.2.4	Manual internal disinfection	61
7.2.5	Manual drying	63

7.3	Automated reprocessing.....	64
7.3.1	Automated internal and external cleaning and internal and external disinfection...	67
7.3.2	Automated drying	67
7.4	Care products and systems - Servicing	69
7.4.1	Servicing with KaVo Spray	70
7.4.2	Servicing with KaVo QUATTROcare PLUS	71
7.4.3	Care with KaVo SPRAYrotor.....	75
7.4.4	Servicing with KaVo QUATTROcare	76
7.5	Packaging	79
7.6	Sterilisation.....	80
7.7	Storage	83
8	Tools and consumables	84
9	Terms and conditions of warranty	86

1 User instructions

Dear User

Congratulations on purchasing this Kerr Endodontics product. By following the instructions below you will be able to work smoothly, economically and safely.

KaVo and elements 8:1 are either registered trademarks or trademarks of KaVo Dental GmbH.

All other trademarks are property of their respective owners.

Kerr Customer Care

For repairs, please contact your local dealer or Kerr Customer Care directly:

Toll-free: +1-800-KERR-123

Homepage: www.kerrdental.com

Email: KerrCustCare@kavokerr.com

Target group

This document is for dentists and dental office staff.

The section on startup is also intended for the service staff.

General marks and symbols

See Chapter on User Instructions/Hazard Levels

	Important information for users and service technicians
	Action request
	CE mark (European Community). A product bearing this mark meets the requirements of the applicable EC directive.
	Medical device, labelling of medical devices
	Can be steam-sterilised at 134 °C -1 °C / +4 °C (273 °F -1.6 °F / +7.4 °F)
	Thermodisinfectable

Information on the packaging

	Follow instructions for use
	Please note the instructions for use
	Caution
	Temperature range
	Air pressure

	Relative Humidity
Rx only	Caution: Federal (USA) law restricts device to sale by or on the order of a dental professional.
	CE mark (European Conformity mark)
	Medical device, labelling of medical devices
	Importer
	Transport upright with the arrows pointing upwards!
	Fragile - protect against impact!

	Stacking limit by number
	Protect from moisture (Keep dry)!
	Do not dispose this product into the ordinary municipal waste or garbage system
	Properly dispose of non-corrugated fiberboard
	Properly dispose of polyethylene
	Properly dispose of all other plastics such as polycarbonate (PC), polyamide (PA), styrene acrylonitrile (SAN), polymethylmethacrylate (PMMA), polylactides (PLA), etc.

Hazard levels

The warning and safety notes in this document must be observed to prevent personal injury and material damage. The warning notes are designated as shown below:

DANGER

In cases which – if not prevented – directly lead to death or severe injury.

WARNING

In cases which – if not prevented – can lead to death or severe injury.

**CAUTION**

In cases which – if not prevented – can lead to minor or moderate injury.

NOTICE

In cases which – if not prevented – can lead to material damage.

2 Safety

Note

All serious events occurring in relation to the product must be reported to the manufacturer and the competent authority of the member state, in which the user and/or patient resides.

The instructions for use are a component of the product and must be read carefully prior to use and be accessible at all times.

The device may only be used in accordance with the intended use, any other type of use is not permitted.

2.1 Infection hazard

Patients, users or third parties could be infected by contaminated medical devices.

- ▶ Take suitable personal protective measures.
- ▶ Follow the instructions for use of the components.

- ▶ Before initial startup and after each use, reprocess the product and accessories appropriately.
- ▶ Carry out the reprocessing as described in the instructions for use. The procedure has been validated by the manufacturer.
- ▶ If you deviate from this procedure, it is essential to make sure that the reprocessing is effective.
- ▶ Reprocess the product and accessories appropriately before disposal.
- ▶ If there is any injury to soft tissue, do not continue treatment in the oral cavity with compressed air-driven instruments.
- ▶ Use gloves or finger guard whenever you test, insert, and remove the tool.

2.2 Improper use

The improper use of the device could lead to burns or injuries.

- ▶ Check the technical condition before each use.

See also:

 2.3 Technical condition, Page 19

- ▶ Never press the push-button during operation of the device.
- ▶ Never use the instrument to keep the cheek, tongue or lip at a distance.
- ▶ Never touch soft tissue with the handpiece head or instrument cover.
- ▶ Do not use the medical device as a light probe.
- ▶ Use an appropriate light probe for illumination of the oral cavity or site of preparation.
- ▶ After treatment, place the medical device properly in the cradle without the tool.

2.3 Technical condition

A damaged device or components could injure patients, users and third parties.

- ▶ Only operate devices or components if they are undamaged on the outside.
- ▶ Check that the device is working properly and is in satisfactory condition before each use.
- ▶ Have parts with sites of breakage or surface changes checked by the Service.
- ▶ If the following defects occur, stop working and have the service personnel carry out repair work:
 - Malfunctions
 - Damage
 - Irregular running noise
 - Excessive vibration
 - Overheating

- Tool is not seated firmly in the handpiece

Observe the following instructions in order to guarantee optimum functioning and prevent material damage:

- ▶ Service the medical device with care products and systems regularly as described in the instructions for use.
- ▶ The device should be reprocessed and stored in a dry location, according to instructions, if it is not be used for a longer period.

High torque of micromotors can lead to severe burn injuries.

- ▶ Service micromotors regularly.
- ▶ Do not use any damaged motors.
- ▶ Do not use motors for unauthorised purposes.

2.4 Accessories and combination with other equipment

Use of un-authorised accessories or un-authorised modifications of the device could lead to injury.

- ▶ Only use accessories that have been approved for combination with the product by the manufacturer.
- ▶ Only use accessories that are equipped with standardised interfaces.
- ▶ Do not make any modifications to the device unless these have been approved by the manufacturer of the product.

The lack of control equipment for changing the speed range and the direction of rotation can lead to injury.

- ▶ Control facility for changing the speed and the direction of rotation must be present.

- ▶ The medical device may only be combined with a treatment centre / control facility released by KaVo.
- ▶ Comply with the Instructions for Use of the treatment centre / control unit.

2.5 Qualification of personnel

Application of the product by users without the appropriate medical training could injure the patients, the users or third parties.

- ▶ Make sure that the user has read and understood the instructions for use.
- ▶ Only employ the device if the user has the appropriate medical training.
- ▶ Observe national and regional regulations.

2.6 Service and repair

Repairs, servicing and safety checks may only be performed by trained service personnel. The following persons are authorised to do this:

- Service technicians of KaVo branches after the appropriate product training
- Service technicians of KaVo authorised dealers after the appropriate product training

Observe all the following items during servicing work:

- ▶ Have the service and testing tasks carried out according to the Medical Device Operator Ordinance.
- ▶ Following expiry of the warranty, have the tool holding system checked once a year.
- ▶ KaVo recommends specifying in-house service intervals where the medical device is brought to a professional shop for cleaning, servicing and a function check. Define the service interval depending on the frequency of use.

As a result of the use of NON-KaVo original spare parts during the repair, parts such as covers may become undone and injure the patient, user or other people. This may result in aspiration, swallowing of parts and possibly even a risk of suffocation.

- ▶ Only use spare parts that comply with the specification for repair; original KaVo spare parts comply with the specification.

**Note**

If a repair is done with NON-KaVo original spare parts, this may constitute a product modification that leads to the loss of CE conformity. In the event of damage, the responsibility is with the service company or the operator.

The introduction into the market of a modified product, where there is reasonable suspicion that the safety and health of patients or users may be jeopardised, is prohibited by the German medical device law §4, section 1 no. 1 and requires a separate conformity check.

3 Product description

elements™ 8:1, PN: 815-1655

3.1 Purpose – Intended use

Indications for use:

This medical device is

- intended for dental treatment only. All other types of use of or modifications to the product are not permitted and can be hazardous. The medical device is intended for the following uses: cavity preparation and endodontics.
- A medical device according to relevant national statutory regulations.

Proper use:

According to these regulations, this product may only be used for the described application by a properly trained user. You need to comply with the following:

- the applicable health and safety regulations
- the applicable accident prevention regulations

- these Instructions for use

According to these regulations, the user is required:

- to only use equipment that is operating correctly
- adhere to the specified intended use
- to protect him or herself, the patient and third parties from hazards
- to prevent contamination from the product

3.2 Technical Specifications

Drive speed	max. 40,000 rpm
Speed transmission:	8:1
Cooling air flow	5.5 to 9.5 l/min

With push-button chuck.

Contra-angle drills and root canal instruments can be inserted.

The contra-angle handpiece can be mounted on all INTRAmatic (LUX) motors and motors with a connector in accordance with ISO 3964 / DIN 13940.

3.3 Symbols on product and rating plate

The rating plates are affixed to the underside of the unit.

Accompanying documents

	Caution
	Follow instructions for use
	Consult instructions for use
	Please note the electronic instructions for use

Certification

	CE mark (Communauté Européenne)
	Medical device, labelling of medical devices

Product characteristics

	Manufacturer
Type	Device type
SN	Serial number
REF	Material number

	Type B applied part
	Supply voltage
	Operating mode: continuous operation with intermittent load
	Class II equipment
	Do not dispose this product into the ordinary municipal waste or garbage system

3.4 Transportation and storage conditions

NOTICE

Startup after refrigerated storage.

Malfunction.

- Prior to startup, strongly refrigerated products must be allowed to warm up to a temperature of 20 °C to 25 °C (68 °F to 77 °F).

	Temperature: -20 °C to +70 °C (-4 °F to +158 °F)
	Relative humidity: 5% RH to 95% RH absence of condensation

Air pressure: 700 hPa to 1060 hPa (10 psi to 15 psi)

Protect from moisture (Keep dry)

4 Startup and shut down

WARNING

Hazard from non-sterile products.

Infection hazard for dentist and patient.

- ▶ Prior to initial startup and after each use, reprocess the product and accessories.

WARNING

Dispose of the product in appropriate manner.

Infection hazard.

- ▶ Reprocess and sterilise the product and accessories before disposal.

See also:

-  7 Reprocessing steps in accordance with ISO 17664, Page 53

NOTICE**Damage from soiled and moist cooling air.**

Contaminated and moist cooling air can cause malfunctions.

- Make sure that the supply of cooling air is dry, clean, and uncontaminated according to EN ISO 7494-2.

NOTICE**Damage to the medical device caused by spray air and spray water.**

Property damage

- ▶ Un-select spray air and spray water on the supply unit before startup!

5 Operation

5.1 Attaching the medical device

CAUTION

Detachment of the medical device during treatment.

Injury or property damage

A medical device that is not properly locked in place can become disconnected from the motor coupling and fall off.

- ▶ Carefully pull on the medical device before each treatment to ensure that it is securely locked onto the motor coupling.

NOTICE**Removing and attaching the medical device while the drive motor is rotating.**


Damage to the driver.

- ▶ Never attach or remove the medical device while the drive motor is rotating.

NOTICE**Pressing the foot switch while attaching or detaching the medical device.**

Property damage to the medical device.

- ▶ Do not connect or remove the medical device while pressing the foot switch.

- ▶ Lightly spray O-rings on the motor coupling with KaVo Spray.
- ▶ Place the medical device on the motor coupling and lock it into place.
- ▶ Pull on the medical device to make sure that it is securely affixed to the coupling.

5.2 Removing the medical device

- ▶ Unlock the medical device from the motor coupling by twisting it slightly and then pulling it along its axis.

5.3 Inserting the dental bur or diamond grinder

Note

Only use root canal instruments with shafts that comply with ISO 1797-1 type 1 and ISO 3630-1:

- Shaft diameter: Ø 2.334 to Ø 2.350 mm
- Shaft clamping length: at least 12 mm

**WARNING****Use of unauthorised dental burs.**

Risk of injury.

- ▶ Comply with the instructions for use and the intended use of the dental bur.
- ▶ Only use dental bur that do not deviate from the specified data.

**⚠ WARNING****Hazard from rotating dental bur.**

Lacerations.

- ▶ Do not touch the rotating dental bur!
- ▶ Remove the dental burs from the medical device after treatment to avoid injury and infection during storage.

**⚠ CAUTION****Dental bur with worn or damaged shafts.**

Risk of injury, dental bur may fall out during treatment.

- ▶ Never use a dental bur with damaged or worn shafts.

**CAUTION****Hazard from defective chuck system.**

The dental bur can fall out and cause injury.

- ▶ Pull on the dental bur to check if the clamping system works properly and if the dental bur is firmly clamped. Wear gloves or a thimble when you check, insert, or remove the bits to prevent injury and infection.

NOTICE

Dental bur shaft slips inside the chuck due to excessive speed of the dental bur or abrupt engagement of the dental bur.

Material damage to dental bur shaft and chuck system, reduction of the service life of dental bur and chuck system.


- ▶ Do not operate the dental bur at a higher speed than recommended by the manufacturer.

NOTICE

Dental bur with worn or damaged shafts.

Material damage to the chuck system, dental bur is difficult or impossible to remove from the chuck system.

- ▶ Never use a dental bur with damaged or worn shafts.

- ▶ Insert the dental bur into the segment of the head drive by twisting the dental bur slightly, and push it to the end stop.
- ▶ Check if the dental bur is seated securely by pulling on it.

5.4 Removing the milling tool or diamond grinder


WARNING

Hazard from rotating dental bur.

Lacerations and damage to the chuck system.

- ▶ Do not touch the rotating dental bur!
- ▶ Never push the push-button while the dental bur is rotating.
- ▶ Remove the dental burs from the medical device after treatment to avoid injury and infection during storage.

- ▶ After the dental bur has stopped rotating, press the push-button down with your thumb and simultaneously remove the dental bur.

6 Checking for malfunctions and troubleshooting

6.1 Check for malfunctions

CAUTION

Product heats up.

Burn injury or product damage due to over-heating.

- ▶ Do not use the product if it is irregularly heated.

NOTICE

Missing or damaged O-rings.

Malfunction and premature failure.

- ▶ Make sure that all O-rings are on the coupling and are undamaged.
-
- ▶ The medical device is too hot while idling:
Check the amount of cooling air.
 - ▶ The medical device is too hot while working:
Service the medical device.
 - ▶ When the speed drops or is uneven:
Service the medical device.
 - ▶ An O-ring is missing on the motor coupling:
Replace O-ring.

6.2 Troubleshooting

WARNING

Repair WITHOUT using KaVo original spare parts.

Parts such as the cover can come loose and cause injury.
Aspiration, swallowing of parts and danger of suffocation.

- ▶ Only use spare parts that comply with the specification for repair; original KaVo spare parts comply with the specification.

**Note**

If a repair is done with NON-KaVo original spare parts, this may constitute a product modification that leads to the loss of CE conformity. In the event of damage, the responsibility is with the service company or the operator.

The introduction into the market of a modified product, where there is reasonable suspicion that the safety and health of patients or users may be jeopardised, is prohibited by the German medical device law §4, section 1 no. 1 and requires a separate conformity check.

6.2.1 Replacing the O-rings on the motor coupling

NOTICE

Improper care of the O-rings.

Malfunction or complete failure.

- ▶ Do not use Vaseline or other grease or oil.

Note

The O-rings on the motor coupling may only be lubricated with a cotton ball wetted with KaVo spray.

- ▶ Press the O-ring between your fingers to form a loop.
- ▶ Push the O-ring to the front, and remove it.
- ▶ Insert new O-rings into the grooves.

7 Reprocessing steps in accordance with ISO 17664

7.1 Preparations at the site of use

WARNING

Hazard from contaminated products.

Contaminated products are associated with an infection hazard.

- ▶ Take suitable personal protective measures.

WARNING

Sharp tool in the medical device.

Injury hazard from sharp and/or pointed tool.

- ▶ Remove the tool.

- ▶ Reprocess the medical device as soon as possible after treatment.
- ▶ The medical device must be dry when transported to reprocessing.
- ▶ To minimise the risk of infection during reprocessing, always wear protective gloves.
- ▶ Remove the tool from the medical device.
- ▶ Remove all residual cement, composite or blood immediately.
- ▶ Do not place in solutions or similar substances.

7.2 Manual Reprocessing

WARNING

Sharp tool in the medical device.

Injury hazard from sharp and/or pointed tool.

- ▶ Remove the tool.

NOTICE

Never reprocess this medical device in an ultrasonic cleaner.


Malfunction and material damage.

- ▶ Clean manually or in a washer disinfectant only.

7.2.1 Manual external cleaning

Accessories required:

- Tap water $30\text{ }^{\circ}\text{C} \pm 5\text{ }^{\circ}\text{C}$ ($86\text{ }^{\circ}\text{F} \pm 10\text{ }^{\circ}\text{F}$)
- Brush, e.g. medium-hard toothbrush

- ▶ Brush off under flowing tap water.

7.2.2 Manual internal cleaning

Validated internal cleaning (removal of residual protein) can be accomplished with KaVo CLEANspray.

- ▶ Cover the medical device with the KaVo Cleanpac bag, and place it on the corresponding care adapter.

- ▶ Hold the can vertically.
- ▶ Press the spray button three times for 2 seconds each time.
- ▶ Remove the medical device from the spray attachment and let the cleanser act for 1 minute.

See also:

 KaVo CLEANspray Instructions for Use

- ▶ If a manual external and internal disinfection do not follow directly, dry the medical device with KaVo DRYspray.

See also:

 7.2.5 Drying, Page 63

7.2.3 Manual external disinfection

WARNING

Incomplete disinfection.

Infection hazard.

- ▶ Only use disinfection procedures that are verified to be bactericidal, fungicidal and virucidal.
- ▶ If the disinfectants used do not meet these requirements, the process must be concluded by disinfection of the unit(s) without packaging using a steam steriliser.

NOTICE

Never disinfect the medical device with chloride-containing products.

Malfunction and material damage.

- ▶ Only disinfect in a washer disinfector or manually.

KaVo recommends the following products based on the compatibility of the materials. The microbiological efficacy must be ensured by the disinfectant manufacturer and proven by an expert opinion.

Approved disinfectants:

- CaviWipes and CaviCide made by Metrex (intermediate disinfection)
- Mikrozid AF made by Schülke & Mayr (liquid or cloths)
- FD 322 made by Dürre

Consumables required:

Cloths for wiping the medical device.

- ▶ Spray the disinfectant on a cloth, then wipe down the medical device and allow the disinfectant to act according to the instructions of the disinfectant manufacturer.
- ▶ Follow the instructions for use of the disinfectant.

7.2.4 Manual internal disinfection

WARNING

Incomplete disinfection.

Infection hazard.

- ▶ Only use disinfection procedures that are verified to be bactericidal, fungicidal and virucidal.
- ▶ If the disinfectants used do not meet these requirements, the process must be concluded by disinfection of the unit(s) without packaging using a steam steriliser.

NOTICE

Never disinfect the medical device with chloride-containing products.

Malfunction and material damage.

- ▶ Only disinfect in a washer disinfector or manually.

The efficacy of manual internal disinfection must be demonstrated by the manufacturer of the disinfection agent. With KaVo products, use only disinfection agents that have been released by KaVo with respect to the compatibility of materials (e.g. WL-cid / made by ALPRO).

- ▶ Cover the medical device with the KaVo CLEANpac bag, and place it on the corresponding care adapter.
- ▶ Hold the can vertically.

- ▶ Press the spray key for at least 3 seconds.
- ▶ Remove the medical device from the spray attachment and let the disinfectant act for 2 minutes.
- ▶ Follow the instructions for use of the disinfectant.

7.2.5 Manual drying

Use KaVo DRYspray for subsequent drying of the air, water and gear unit ducts.

- ▶ Cover the medical device with the KaVo CLEANpac bag, and place it on the corresponding care adapter.
- ▶ Hold the can vertically.
- ▶ Press the spray key for at least 3 seconds.

See also:

 KaVo DRYspray Instructions for Use

- ▶ Immediately after drying, lubricate the KaVo medical device with care agents from the KaVo care system.

See also:

- 📄 7.4 Care products and systems - Servicing, Page 69

7.3 Automated reprocessing

WARNING

Incomplete disinfection.

Infection hazard.

- ▶ Only use disinfection procedures that are verified to be bactericidal, fungicidal and virucidal.
- ▶ If the disinfectants used do not meet these requirements, the process must be concluded by disinfection of the unit(s) without packaging using a steam steriliser.

**WARNING****Sharp tool in the medical device.**

Injury hazard from sharp and/or pointed tool.

- ▶ Remove the tool.

NOTICE**Never disinfect the medical device with chloride-containing products.**

Malfunction and material damage.

- ▶ Only disinfect in a washer disinfector or manually.

NOTICE

Never reprocess this medical device in an ultrasonic cleaner.

Malfunction and material damage.

- ▶ Clean manually or in a washer disinfectant only.

7.3.1 Automated internal and external cleaning and internal and external disinfection

KaVo recommends washer disinfectors according to EN ISO 15883-1, which are operated using alkaline cleaning agents having a maximum pH value of 10.

The validation was performed in a Miele washer disinfector using the "VARIO-TD" program, the "neodisher mediclean" cleaning agent, the "neodisher Z" neutralizer, and the "neodisher miel-clear" rinsing agent.

- ▶ For programme settings as well as cleansers and disinfectants to be used, please refer to the Instructions for Use of the thermodisinfectant.

7.3.2 Automated drying

The drying procedure is normally part of the cleaning programme of the washer disinfector.

**Note**

Please comply with the instructions for use of the washer disinfectant.

- ▶ In order to prevent impairment of the KaVo medical device, make sure that the inside and outside of the device is dry after the end of the cycle.
- ▶ Remove any residual liquids with KaVo DRYspray.

See also:

7.2.5 Manual drying, Page 63

- ▶ Immediately after drying, lubricate the KaVo medical device with care agents from the KaVo care system.

7.4 Care products and systems - Servicing

WARNING

Sharp tool in the medical device.

Injury hazard from sharp and/or pointed tool.

- ▶ Remove the tool.

CAUTION

Improper service and care.

Risk of injury.


- ▶ Service regularly with suitable agents.

**Note**

KaVo guarantees the proper function of KaVo products only if the care products listed as accessories are used, since these were tested for proper use on our products.

7.4.1 Servicing with KaVo Spray

KaVo recommends servicing the product as part of the reprocessing after each use, i.e. after each cleaning, disinfection, and before each sterilisation, but no later than after 30 minutes of operation and generally after every 6th patient.

- ▶ Remove the tool from the medical device.
- ▶ Cover the medical device with the KaVo CLEANpac bag, and place it on the corresponding care adapter.
- ▶ Press the spray key for 1 to 2 seconds.

Servicing the chucking system

KaVo recommends servicing the chucking system once weekly.

- ▶ Remove the tool from the medical device.
- ▶ Position the tip of the spray nipple in the opening, and apply the spray.
- ▶ Press the spray key for 1 to 2 seconds.

7.4.2 Servicing with KaVo QUATTROcare PLUS

Servicing and cleaning device with expansion pressure for the interior cleaning of inorganic residues and optimum care.

(no validated cleaning of the interior in accordance with German RKI requirements)

KaVo recommends servicing the product as part of the reprocessing after each use, i.e. after each cleaning, disinfection, and before each sterilisation, but no later than after 30 minutes of operation and generally after every 6th patient.

- ▶ Remove the tool from the medical device.
- ▶ Service the product in the QUATTROcare PLUS.

See also:

 Instructions for use KaVo QUATTROcare PLUS

Servicing the chuck


KaVo recommends servicing the chuck system once a week using the chuck servicing program integrated in the device.

Note

Handpieces must be taken off the service couplings before the chuck service can be started and performed.

- ▶ Close the front door and press the chuck service button for at least three seconds until the spray canister control LED flashes three times consecutively.
- ⇒ The device is in chuck service mode.
- ▶ Remove the service coupling of the chuck from the side hatch of the QUATTROcare PLUS and attach it to coupling service point four, on the far right. A MULTIflex adaptor must be mounted there.

- ▶ Press the handpiece together with the guide bush of the chuck to be serviced against the tip of the service coupling.
- ▶ Press the button marked with the chuck service symbol.

**Note****Close the chuck service mode.**

Option 1: Place the dental handpieces in the QUATTROcare PLUS 2124 A, close the front door and start the service procedure.

Option 2: After three minutes with no service procedure running, the device automatically switches back to normal service mode.

See also:

 Servicing with KaVo QUATTROcare PLUS

7.4.3 Care with KaVo SPRAYrotor


**Note**

KaVo SPRAYrotor is no longer included in the current delivery programme.

Follow-up product:

- ▶ QUATTROcare PLUS 2124 A

KaVo recommends servicing the product as part of the reprocessing after each use, i.e. after each cleaning, disinfection, and before each sterilisation, but no later than after 30 minutes of operation and generally after every 6th patient.

- ▶ Cover the medical device with the Cleanpac bag, and place it on the corresponding servicing adapter on the KaVo SPRAYrotor.
- ▶ Service the product.

See also:

 Instructions for use KaVo SPRAYrotor

7.4.4 Servicing with KaVo QUATTROcare

**Note**

QUATTROcare 2104 / 2104 A is no longer included in the current delivery programme.


Follow-up product:

- ▶ QUATTROcare PLUS 2124 A

Servicing and cleaning device with expansion pressure for the interior cleaning of inorganic residues and optimum care.

(no validated cleaning of the interior in accordance with German RKI requirements)

KaVo recommends servicing the product as part of the reprocessing after each use, i.e. after each cleaning, disinfection, and before each sterilisation.

- ▶ Remove the tool from the medical device.
- ▶ Service the product in the QUATTROcare.

See also:

📄 Instructions for use KaVo QUATTROcare 2104 / 2104A / 2124A

Servicing the chuck

KaVo recommends servicing the chucking system once weekly.

- ▶ Remove the tool from the medical device.

- ▶ Plug the spray nipple of the chuck servicing set onto the QUATTROcare plus Spray.
- ▶ Position the tip of the spray nipple in the opening, and apply the spray.
- ▶ Press the spray key for 1 to 2 seconds.

7.5 Packaging

Note

The sterile goods package must be large enough for the product so that the packaging is not stretched.

The quality and use of the packaging of the items to be sterilised must satisfy the applicable standards and be appropriate for the sterilisation process!

- ▶ The medical device must be packed before sterilisation.

7.6 Sterilisation

Sterilisation in a steam steriliser (autoclave) in accordance with EN 13060 / EN ISO 17665-1

CAUTION

Improper service and care.

Risk of injury.

- ▶ Service regularly with suitable agents.

NOTICE

Contact corrosion due to moisture.

Damage to product.

- ▶ Immediately remove the product from the steam steriliser after the sterilisation cycle.

The KaVo medical device has a maximum temperature resistance up to 138 °C (280.4 °F).

Sterilisation parameters:

Select a suitable process from the following sterilisation processes (depending on the available steriliser):

- Steriliser with triple pre-vacuum:
 - at least 3 minutes at 134 °C -1 °C / +4 °C (273 °F -1.6 °F / +7.4 °F)
- Steriliser using the gravity method:
 - at least 10 minutes at 134 °C -1 °C / +4 °C (273 °F -1.6 °F / +7.4 °F)
 - at least 60 minutes at 121 °C -1 °C / +4 °C (250 °F -1.6 °F / +7.4 °F)
- ▶ Remove the medical device immediately after the completion of the sterilisation cycle from the steriliser.
- ▶ Use according to the manufacturer's Instructions for Use.

7.7 Storage

Prepared products must be stored appropriately in a dry, dark, cool room such that they are protected from germs (as far as possible) and dust.

Note

Comply with the expiry date of the sterilised items.

8 Tools and consumables

Available from dental suppliers.

Availability may vary depending on region.

Material summary	Mat. No.
INTRA instrument stand	3.005.5204
Cleanpac 10 units	0.411.9691
Cellulose pad 100 units	0.411.9862
Spray head INTRA (KaVo Spray)	0.411.9911
Service coupling INTRA (QUATTROcare)	1.009.6143

Material summary	Mat. no.
Adaptor INTRAmatic (CLEANspray and DRYspray)	1.007.1776
KaVo CLEANspray 2110 P	1.007.0579
KaVo DRYspray 2117 P	1.007.0580
KaVo Spray 2112 A	0.411.9640
ROTAspray 2 2142 A	0.411.7520
QUATTROcare plus Spray 2140 P	1.005.4525
Chuck servicing set	1.003.1253

9 Terms and conditions of warranty

The following warranty conditions apply to this KaVo medical device:

KaVo provides the end customer with a warranty of proper function and guarantees zero defects in respect of material and workmanship for a period of 12 months from the date of the invoice, subject to the following conditions:

In case of justified complaints, KaVo will honour its warranty with a free replacement or repair. Other claims of any kind whatsoever, in particular with respect to compensation, are excluded. In the event of default, gross negligence or intent, this shall only apply in the absence of mandatory legal regulations to the contrary.

KaVo shall not be liable for defects and their consequences that have arisen or may arise from natural wear, improper handling, cleaning or maintenance, non-compliance with operating, main-

tenance or connection instructions, calcination or corrosion, contaminated air or water supplies or chemical or electrical factors deemed abnormal or impermissible in accordance with KaVo's instructions for use or other manufacturer's instructions. The warranty granted does not usually extend to lamps, optical fibres made of glass and glass fibres, glassware, rubber parts, and the colourfastness of plastic parts.

All liability is excluded if defects or their consequences originate from manipulations or changes to the product made by the customer or a third party that is not authorised by KaVo.

Warranty claims will only be accepted if the product is submitted along with proof of purchase in the form of a copy of the invoice or note of delivery. The dealer, purchase date, type, and serial number must be clearly evident from this document.

Instrucciones de uso elements™ 8:1

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Kerr Corporation
1717 West Collins Avenue
Orange CA 92867
EE.UU.
+1-800-KERR-123
www.kerrdental.com

**Fabricante:**

KaVo Dental GmbH
Bismarckring 39
88400 Biberach
Alemania
www.kavo.com

PN: 077-0649
06/2021
Rev. A

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Kerr Australia Pty. Limited
Unit 6, 12 Mars Road
Lane Cove West, New South Wales 2066
Australia
+61-2-8870-3000

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