

CE



# INSTALLATION INSTRUCTIONS



## **AUT SERIES SURGICAL LIGHTS**

**(STANDARD RADIAL ARM MODELS ONLY)**

Read this manual before starting to work! This information is necessary  
for the safe and efficient operation of the equipment.

## COMPETENCY AND INSTALLATION REQUIREMENTS

This product is a Class 2 medical device that is subject to FDA Part 820 requirements. Installation can only take place by qualified and trained individuals. An Installation Qualification Report is required as proof of system operational validation prior to clinical use. Contact Skytron for installation needs.

### Distributed by:

US - SKYTRON  
5085 Corporate Exchange Blvd. S.E.  
Grand Rapids, MI 49512 (616) 656-2900  
[www.skytron.us](http://www.skytron.us)



### Manufactured by:

DKK Dai-Ichi Shomei Co., LTD  
32-26 Sakashita 1-Chome  
Itabashi-Ku, Tokyo 174-0043  
JAPAN



SKYTRON EUROPE BV  
Floresstraat 52  
8022 AD Zwolle  
THE NETHERLANDS

The base language for this document is ENGLISH. Any translations must be from the base language document.

Printed copies are not controlled documents.

Although current at the time of publication, SKYTRON'S policy of continuous development makes this manual subject to change without notice. If current manuals are required, contact your local SKYTRON representative or contact SKYTRON directly at the distribution addresses listed above.

# TABLE OF CONTENTS

TITLE	PAGE
SECTION 1. SAFETY INFORMATION.....	2
1-1. Special User Attention .....	2
1-2. Safety Precautions .....	2
1-3. Label Symbols .....	3
SECTION 2. EQUIPMENT SPECIFICATIONS / REQUIREMENTS.....	4
2-1. Permissible Environmental Conditions.....	4
2-2. Electrical Requirements .....	4
2-3. ESD Sensitive Devices.....	4
2-4. Connection Means .....	4
2-5. National Electric Code, NFPA & ANSI Requirements.....	4
2-6. Bending Requirements.....	4
2-7. Conduit Requirements.....	5
2-8. Protective Means.....	5
2-9. Final Assembly .....	5
2-10. Fail Safe Compliance .....	5
2-11. Shipping.....	5
SECTION 3. MODEL IDENTIFICATION.....	6
SECTION 4. PRE-INSTALLATION INSTRUCTIONS .....	7
4-1. Pre-Installation Requirements .....	7
4-2. Installation Notes.....	7
4-3. Uncrating .....	7
4-4. Specialty Tools and Equipment .....	7
SECTION 5. INSTALLATION.....	8
5-1. Typical Installation Sequence / Component Identification .....	8
5-2. Install Wall Control.....	9
5-3. Optional Camera Ready Lighthead .....	12
5-4. Install Fixture Mounting Plate .....	12
5-5. Install Radial Arm Assembly and Ceiling Cover .....	12
5-6. Install VST and BOM .....	13
5-7. Install AUT7 Lightheads (If Applicable).....	14
5-8. Install Model AUT5 Lightheads (If Applicable).....	15
5-9. Fixture Adjustments.....	15
SECTION 6. INSTALLATION DRAWINGS.....	16
6-1. Typical Light Fixture and Wiring Requirements (AUT55 Shown) .....	16
6-2. Typical Wiring Diagram.....	17
6-3. Typical Light Fixture and Wiring Requirements (TV Models) .....	18
6-4. Typical Wiring Diagram (TV Models).....	19
6-5. Standard Mounting Structure Guideline .....	20
SECTION 7. REVISION HISTORY.....	21

## SECTION 1. SAFETY INFORMATION

## 1-1. Special User Attention

The procedures described in this manual will be performed by representatives of the owner (staff or contracted service), therefore it is the responsibility of the owner to ensure that all safety precautions are followed. Only qualified and trained individuals should attempt the installation of this product.

## 1-2. Safety Precautions

The following is a summary of DANGERS, WARNINGS, and CAUTIONS denoted in this manual. These precautions are found throughout the manual where they are applicable. Carefully read the manual before proceeding to operate or service the equipment.

**DANGER**

**DANGER with the safety alert symbol, is used to indicate a hazardous situation that, if not avoided, will result in death or serious injury.**

**WARNING**

**WARNING with the safety alert symbol, is used to indicate a hazardous situation that, if not avoided, could result in death or serious injury.**

To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth ground.

DO NOT remove lighthead when support arm is in down position. The BOM will be severely damaged and it may result in bodily injury.

**CAUTION**

**CAUTION with the safety alert symbol, is used to indicate a hazardous situation that, if not avoided, could result in minor or moderate injury.**

Connection of the fixture wires must be made using crimp connectors. Main terminal devices shall be so located or shielded that, should a wire of a stranded conductor escape when the conductors are fitted, there is no risk of accidental contact between live parts and accessible parts. Acceptable shielding methods include UL approved shrink tubing and electrical tape. DO NOT use damaged wire.

**CAUTION**

**CAUTION without the safety alert symbol, is used to address practices not related to personal injury but with a possibility of damage to equipment.**

This equipment is intended for use by healthcare professionals only. This equipment may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the lighthead or shielding the location.

This fixture requires that electrical connections are made by a licensed electrician in accordance with state, local, and national electrical codes using UL (Underwriters Laboratory) recognized materials.

This fixture requires two (2) dedicated conduit raceways at the wall control to separate the 100-240VAC facility supply lines from the DC supply lines to the lighting fixture. Failure to observe this requirement will allow the migration of electrical magnetic interference and will disrupt the operation of the lights.

**CAUTION (CONT'D)**

This fixture requires a properly circuit protected, appropriately sized, dedicated circuit. An isolated power supply circuit must be protected by an appropriately sized double pole, single throw circuit breaker.

SKYTRON surgical lights are packaged in special containers designed to prevent damage from vibration or shock. Always use SKYTRON supplied containers for shipment.

Aurora 3 lightheads operate on DC VOLTAGE. The PC boards are susceptible to static charges even when not powered. Pay close attention to wiring diagrams, wire labeling and color

codes. Wires must remain separate and not touch any other wires or metal parts. Incorrect wiring may result in incorrect polarity being supplied to the lighthead. This WILL DAMAGE internal circuitry and components VOIDING WARRANTY.








The mounting plate must be accurately leveled within 0.1° to prevent lighthead drift.

To prevent support arm damage, the longer screws must be installed in the holes towards the lighthead.

**NOTICE**

**Indicates important information not related to personal injury.**

**1-3. Label Symbols**

Symbol	Description
	With the word DANGER, indicates a hazardous situation that, if not avoided, will result in death or serious injury.
	With the word WARNING, indicates a hazardous situation that, if not avoided, could result in death or serious injury. With the word CAUTION, indicates a hazardous situation that, if not avoided, could result in minor or moderate injury.
	Indicates AC power supply.
	Indicates authorized representative in the european community.
	Indicates Manufacturer.
	Indicates Dangerous Voltage 100-240V ~, 50/60Hz.
	Risk of overbalancing! To avoid injuries or equipment damage, DO NOT push with excessive force, lean on, or rest on the lighting fixture.

## SECTION 2. EQUIPMENT SPECIFICATIONS / REQUIREMENTS

**CAUTION**

This equipment is intended for use by healthcare professionals only. This equipment may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the lighthouse or shielding the location.

**2-1. Permissible Environmental Conditions****a. During Transport and Storage (in Original Packaging Materials)**

- *Ambient Temperature:*  
14° to 140° F (-10 to 60 °C)
- *Relative Humidity:*  
10% to 85% (No Condensation)
- *Atmospheric Pressure:*  
14 to 31 inHg (500 to 1060 hPa)

**b. During Use - For Dry Locations**

- *Ambient Temperature:*  
60° to 85° F (15 to 30 °C)
- *Relative Humidity:*  
30% to 60% (No Condensation)
- *Atmospheric Pressure:*  
20.7 to 31.3 inHg (700 to 1060 hPa)

**2-2. Electrical Requirements****CAUTION**

This fixture requires that electrical connections are made by a licensed electrician in accordance with state, local, and national electrical codes using UL (Underwriters Laboratory) recognized materials.

DO NOT turn on main power to fixture until all lighthoods are installed, connections are complete, and the fixture has been reviewed by a SKYTRON representative.

**ELECTRICAL HAZARDS EXIST!**

Exercise caution when working on this fixture, the installation of this fixture must be made only by qualified and authorized personnel familiar with the essential knowledge and techniques.

**2-3. ESD Sensitive Devices**

When installing devices with electronic circuit boards (e.g., lighthoods, wall control units), appropriate precautions should be taken to prevent damage caused by electrostatic discharge (ESD). These precautions include as a minimum, the use of an ESD wrist strap that is properly connected to an ESD ground.

**2-4. Connection Means****CAUTION**

**Connection of the fixture wires must be made using crimp connectors. Main terminal devices shall be so located or shielded that, should a wire of a stranded conductor escape when the conductors are fitted, there is no risk of accidental contact between live parts and accessible parts. Acceptable shielding methods include UL approved shrink tubing and electrical tape. DO NOT use damaged wire.**

**2-5. National Electric Code, NFPA & ANSI Requirements**

The installation of connecting cords between equipment parts shall meet the requirements of the National Electrical Code, ANSI/NFPA70, IEC 60601-1 and all local codes, as applicable.

**2-6. Bending Requirements**

Connection leads shall be constructed in such a manner that moveable leads in normal use are not bent around a radius of less than five times the outer diameter of the lead concerned. Avoid conditions employing severe bends to ensure the integrity of conductors.

## 2-7. Conduit Requirements

### CAUTION

This fixture requires two (2) dedicated conduit raceways at the wall control to separate the 100-240VAC facility supply lines from the DC supply lines to the lighting fixture. Failure to observe this requirement will allow the migration of electrical magnetic interference and will disrupt the operation of the lights.

Use of approved metal conduit shall be employed throughout the fixture's wiring circuit where applicable.

## 2-8. Protective Means



### WARNING

**To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth ground.**

### CAUTION

This fixture requires a properly circuit protected, appropriately sized, dedicated circuit. An isolated power supply circuit must be protected by an appropriately sized double pole, single throw circuit breaker.

Proper performance and safety of this fixture can only be achieved by an adequate grounding system. Fixture ground must be a dedicated ground point ultimately bonded to the facilities grounding system to prevent the migration of electrical interference generated by other devices.

## 2-9. Final Assembly

All installations of SKYTRON surgical lights should be under the direct supervision of a SKYTRON authorized representative.

Prior to the fixture being placed in service, the SKYTRON authorized representative must initialize the fixture and complete the installation report.

To maintain product warranty and performance, this product requires routine service. Contact your SKYTRON representative for factory service or preventive maintenance contracts.

## 2-10. Fail Safe Compliance

In order for dual or triple lighthead systems to maintain fail safe compliance, a battery back up (UPS) or generator back up power system must be provided in the mains wiring prior to the wall control which will restore power in five (5) seconds or less.

### NOTICE

Fail safe devices are not supplied by SKYTRON.

## 2-11. Shipping

### CAUTION

SKYTRON surgical lights are packaged in special containers designed to prevent damage from vibration or shock. Always use SKYTRON supplied containers for shipment.

### SECTION 3. MODEL IDENTIFICATION (Standard Radial Arm Configurations)



**MODEL AUT5, AUT5TV**  
**1 - 24" (610mm) Diameter Lighthead**



**MODEL AUT7, AUT7TV**  
**1 - 30" (762mm) Diameter Lighthead**



**MODEL AUT55, AUT5TV5, AUT55TV**  
**2 - 24" (610mm) Diameter Lightheads**



**MODEL AUT555, AUT55TV5, AUT555TV**  
**3 - 24" (610mm) Diameter Lightheads**



**AUT5**

**AUT7**

**MODEL AUT75, AUT7TV5, AUT75TV**  
**1 - 30" (762mm) Diameter Lighthead**  
**1 - 24" (610mm) Diameter Lighthead**



**AUT5**

**AUT5**

**AUT7**

**MODEL AUT575, AUT57TV5, AUT575TV**  
**1 - 30" (762mm) Diameter Lighthead**  
**2 - 24" (610mm) Diameter Lightheads**



## SECTION 4. PRE-INSTALLATION INSTRUCTIONS

### 4-1. Pre-Installation Requirements

Appropriate metal conduit and wiring must be installed from wall control mount to ceiling mounting structure. Flexible conduit to extend 18" below finished ceiling.

100-240VAC, circuit protected, dedicated power supply line in separate conduit to be provided at wall control.

Painting and flooring must be complete prior to fixture installation.

Finished ceiling height must be verified.

Optional camera system models require a 1" metal conduit from mounting plate to control unit J-box.

### 4-2. Installation Notes

Follow the installation instructions and utilize the SKYTRON Surgical Light Installation Report to assure proper installation and to meet the installation qualification requirements.

Special adapter plates for mounting SKYTRON surgical lights on existing mounting structures are available. Contact your SKYTRON representative for special application details.

Additional materials required for proper installation include Loctite® compound.

Aurora Series lighting fixtures require a wall mounted control box. 3/4" metal conduit and minimum 12 AWG wire is required between wall control and fixture.

### 4-3. Uncrating

The SKYTRON surgical lighting fixture is normally shipped in two (2) to four (4) crates, depending on the model. A carton containing the vertical support tubes (VSTs), miscellaneous hardware, and various instructional materials is packed separately.

Should any damage to the fixture be noted while uncrating, further unpacking should be stopped and the container with all the wrappings held for inspection. The transportation company should be notified immediately so an inspector can be sent. Consult the Damaged Shipment Claim Procedure sheet for further details.

Personnel uncrating SKYTRON surgical lights should be aware that they are delicate medical equipment and special care in handling should prevail throughout installation.

### a. Uncrating Light Fixture

Open the top of the lighthead box and remove the packing material. Remove the sterilizable positioning handle and lighthead from the crate.

### NOTICE

Details may vary depending upon model and support structure fabrication. All fixtures use METRIC fasteners.

### b. Radial Arm Assembly

When the radial arm assembly is removed from the crate and set on the floor prior to mounting it, be sure the mounting flange is toward the floor. If the arm assembly is set on the floor upside down the hub cover may be damaged (Figure 4-1).




Figure 4-1. Radial Arm Assembly

### c. Lightheads

Use extreme caution when removing the contents from the crates to prevent damage to the lights. Leave the lightheads in their crates until you are ready to install.

If the lighthead must be set down after it is removed from the crate, always lay it on the foam shipping block. DO NOT lay lighthead on the front face.

### 4-4. Specialty Tools and Equipment


The following specialty tools and equipment are required for performing this installation:

- Metric hand tools
- Appropriately rated lifting device (Genie Lift)
- Digital level
- Torque wrench (ft-lbs)
- Metric hex sockets

## SECTION 5. INSTALLATION

### 5-1. Typical Installation Sequence / Component Identification

The Aurora 3 lighting fixture should be installed in the following sequence (Figure 5-1):



**Figure 5-1. Aurora 3 Installation Sequence (AUT75 Shown)**

1. Wall Control
2. Mounting Plate
3. Radial Arm Assembly (RAA)
4. Ceiling Cover
5. Vertical Support Tube (VST)
6. Balance Mechanism (BOM)
7. Lighthead
8. Sterilizable Positioning Handle

**CAUTION**

Aurora 3 lightheads operate on DC VOLTAGE. The PC boards are susceptible to static charges even when not powered. Pay close attention to wiring diagrams, wire labeling and color codes. Wires must remain separate and not touch any other wires or metal parts. Incorrect wiring may result in incorrect polarity being supplied to the lighthead. This WILL DAMAGE internal circuitry and components VOIDING WARRANTY.

**5-2. Install Wall Control****NOTICE**

3/4" metal conduit and minimum 12AWG wire (3 wires per lighthead plus fixture ground) is required between wall control and fixture. Flexible conduit should extend 18" below finished ceiling.

**NOTICE**

Separate dedicated conduit required for 100-240VAC supply lines to wall control.

All wiring to be in accordance with local, state and national electrical codes.

- a. Remove the front panel assembly from the wall control box for ease in wire connection. Remove the (4) screws. Set the front panel assembly aside.
- b. Install the wall control box enclosure as desired for the application (surface or recessed mount) as shown in the wall control illustration (Figure 5-2).
- c. Attach recess mount flange if required for recessed applications (Figure 5-2).

**NOTICE**

Room placement of the wall control will vary by application. Always follow current standards from the NFPA (National Fire Protection Agency), NEC (National Electrical Code) and IEC (International Electrotechnical Commission) for proper compliance.

The selection of anchorage fasteners shall be determined by the engineer of record and will vary by application. The selected fasteners must not interfere with wall control components.





Seismic applications require the use of approved fasteners.

- d. Connect the electrical conductors from front face plate assembly to the wiring from the fixture. Observe wire markings and colors. Avoid undue stress on conductors and internal components.

**CAUTION**

**Connection of the fixture wires must be made using crimp connectors. Main terminal devices shall be so located or shielded that, should a wire of a stranded conductor escape when the conductors are fitted, there is no risk of accidental contact between live parts and accessible parts. Acceptable shielding methods include UL approved shrink tubing and electrical tape.**

- e. Make electrical connections using approved crimp connectors. Observe wire markings and colors (Figure 5-3).
- f. Attach the front panel assembly using the (4) screws removed in Step a. Use care to avoid pinching conductors and creating excessive bends in wiring.

**OPTIONAL BACK-BOX INSTALLATION****SURFACE MOUNT****RECESSED MOUNT****SIDE VIEW****BACK VIEW****DIMENSIONS**

	SINGLE/DUAL	TRIPLE
A	8-5/8" (220)	10-5/8" (270)
B	4" (100)	4" (100)
C	10-5/8" (270)	14" (355)
D	6-7/8" (175)	8-5/8" (220)
E	7-5/8" (195)	11" (280)
F	6" (153)	6-3/8" (162)
G	10" (255)	13-1/2" (343)
H	7-7/8" (200)	10" (253)
I	11-3/4" (298)	13-3/4" (348)
J	13-7/8" (353)	17-1/4" (438)
K	12-3/8" (315)	15-3/4" (400)
L	10-1/4" (260)	12-1/4" (310)

**Figure 5-2. Wall Control Mounting Options**




Figure 5-3 Wall Control Wiring

### 5-3. Optional Camera Ready Lighthead

The optional camera ready lighthead system includes a connector faceplate with a 65 ft. cable that connects to the light fixture junction box.

The faceplate is mounted in a 2" x 4" junction box and the cable is run through a 1" metal conduit to the light fixture junction box. D-sub 9-pin connectors are provided for cable connection at the fixture junction box and at the RAA/VST connection point.

The camera ready lighthead must be installed with the corresponding camera ready BOM and VST. Refer to AUT TV fixture diagrams (Sections 6-3 and 6-4).


### 5-4. Install Fixture Mounting Plate

a. Check the strength and stability of the mounting structure. It should be fabricated of steel and welded or bolted to the structural ceiling. It should be braced in a manner that will allow no twisting or lateral motion. In standard installations a steel stiffener plate should be used to connect the 3/4" diameter "all-thread" support rods and to provide an attachment base for the angle-iron sway bracing. The 3/4" diameter support rods should be mounted in a 9-1/2" (240mm) square pattern and should extend 2-1/4" (60mm) below the finished ceiling.

#### NOTICE

See mounting structure guideline in Section 6.

b. Install the SKYTRON mounting plate on the threaded rods between jam nuts (Figure 5-4). The plate should normally be located 1-1/4" (32mm) off the finished ceiling (measured from the bottom of the plate) and accurately leveled using a digital level. Tighten the jam nuts securely.



**Figure 5-4. Mounting Plate Installation**

### CAUTION

The mounting plate must be accurately leveled within 0.1° to prevent lighthead drift.

### 5-5. Install Radial Arm Assembly and Ceiling Cover

#### NOTICE

An appropriate lifting device is recommended for installing each radial arm assembly. The weights of the various radial arm assemblies are:

Radial Arm Assembly	Weight
Single AUT 5 Radial Arm	58 pounds (26 kg)
Single AUT 7 Radial Arm	62 pounds (28 kg)
Dual Radial Arm	86 pounds (39 kg)
Triple Radial Arm	117 pounds (53 kg)

#### NOTICE

The multiple arm assemblies are easier to handle during installation if the arms are left taped and tied together.


In some cases it may be necessary to connect the electrical wires from the wall control to the radial arm junction box before the arm assembly can be bolted to the mounting plate.

Use two 6" to 8" bolts to hold fixture while connecting wiring.

a. Install the radial arm assembly (RAA) onto the mounting plate using the (6) bolts provided (Figure 5-5). Torque the 1/2"-13 mounting bolts to 38 ft-lbs (52Nm).

#### NOTICE

Radial arm wires are tagged for proper connection to the wall control (top arm #1, next arm #2, bottom arm [triple arm models] #3).



**Figure 5-5. Radial Arm Installation**

- b. Observe wire tags and color codes and connect the electrical wires from the wall control to the radial arm junction box wires.



### CAUTION

**Connection of the fixture wires must be made using crimp connectors. Main terminal devices shall be so located or shielded that, should a wire of a stranded conductor escape when the conductors are fitted, there is no risk of accidental contact between live parts and accessible parts. Acceptable shielding methods include UL approved shrink tubing and electrical tape.**


- c. Install the ceiling cover and secure.

### 5-6. Install VST and BOM

#### NOTICE


Determine correct placement for each BOM/VST on the RAA. The longest VST goes into the top radial arm.

- a. Install the VST on the BOM, apply Blue Loctite® to screw threads and secure VST with the Allen screws provided (Figure 5-6).




**Figure 5-6. Balance Mechanism**

- b. Follow the tightening sequence and torque the screws to 9 ft-lbs (12 Nm) (Figure 5-7).
- c. Insert the VST of VST/BOM assembly into the RAA receptacle (Figure 5-8). Apply Loctite® to screw threads, and secure the BOM/VST assembly with the mounting screws following the tightening sequence (Figure 5-7).




**Figure 5-7. Screw Tightening Sequence**

- d. Remove cap from RAA (Figure 5-8).



**Figure 5-8. VST/BOM Assembly to RAA**

- e. Pull VST and RAA wires (and connectors for camera ready VST only) up through access hole in RAA (Figure 5-9).



**Figure 5-9. Connect VST/RAA Wires**

f. Observe the wire colors and connect the wires from the RAA to the corresponding VST wires using crimp connectors (Figure 5-9). Connect D-sub (9-pin) connectors for camera ready VST (if applicable).

g. Ensure that all wire connections are secure, then place crimped wires and mated connectors inside RAA access hole. Re-install cap removed in Step d.


h. Repeat procedure for any remaining BOM/VST assemblies.

## 5-7. Install AUT7 Lightheads (If Applicable)

### NOTICE

Refer to Model Identification (Section 3) for correct lighthead placement.

a. To make it easier to install the lighthead, locate the support arm of the BOM so that it points inward toward the ceiling cover (Figure 5-10). This will prevent the radial arm from moving when installing the lighthead.




**Figure 5-10. Position for Lighthead Installation**

- b. Remove the six (6) screws from the lighthead support arm.
- c. Install the lighthead mounting collar onto the support arm and secure with the screws removed in Step b (Figure 5-11).

### CAUTION

To prevent support arm damage, the longer screws must be installed in the holes towards the lighthead.




**Figure 5-11. AUT7 Lighthead Installation**

- d. Pull the lighthead down and remove the shipping block from the BOM (Figure 5-12).



### WARNING

**DO NOT** remove lighthead when support arm is in down position. The BOM will be severely damaged and it may result in bodily injury.




**Figure 5-12. Shipping Block**



### 5-8. Install Model AUT5 Lightheads (If Applicable)


a. To make it easier to install the lighthead, locate the support arm of the BOM so that it points inward toward the ceiling cover (Figure 5-13). This will prevent the radial arm from moving when installing the lighthead.



**Figure 5-13. Position of Lighthead Installation**

b. Remove the four (4) screws from the lighthead mounting stub.

c. Install the lighthead mounting stub into the support arm and secure with the screws previously removed (Figure 5-14).




**Figure 5-14. AUT5 Lighthead Installation**

d. Pull the lighthead down and remove the shipping block from the BOM (Figure 5-15).



### WARNING

**DO NOT** remove lighthead when support arm is in down position. The BOM will be severely damaged and it may result in bodily injury.



**Figure 5-15. Shipping Block**


e. Repeat Steps a through d for each additional AUT5 lighthead.

### 5-9. Fixture Adjustments

Check all fixture positioning axis adjustments and adjust as required. Refer to the service manual for instructions on positioning axis adjustments.

## SECTION 6. INSTALLATION DRAWINGS


### 6-1. Typical Light Fixture and Wiring Requirements (AUT55 Shown)




#### **NOTICE:**

- 2 DEDICATED CONDUIT RUNS REQUIRED AT WALL CONTROL TO SEPARATE 100-240VAC INPUT LINES FROM DC OUTPUT LINES TO LIGHT FIXTURE TO PREVENT THE MIGRATION OF ELECTRICAL MAGNETIC INTERFERENCE WHICH WILL DISRUPT THE OPERATION OF THE LIGHT

## 6-2. Typical Wiring Diagram




### 6-3. Typical Light Fixture and Wiring Requirements (TV Models)



#### NOTICE:

- 2 DEDICATED CONDUIT RUNS REQUIRED AT WALL CONTROL TO SEPARATE 100-240VAC INPUT LINES FROM DC OUTPUT LINES TO LIGHT FIXTURE TO PREVENT THE MIGRATION OF ELECTRICAL MAGNETIC INTERFERENCE WHICH WILL DISRUPT THE OPERATION OF THE LIGHT.
- DO NOT USE A 3P-2P CONVERSION ADAPTER WHEN INSERTING THE PLUG OF THE CAMERA CONTROLLER INTO AN OUTLET.

## 6-4. Typical Wiring Diagram (TV Models)





1. 3/4" support rods located for total support of light, all labor and materials for fabrication supplied by General Contractor. 3/4" nuts and washers for support of SKYTRON mounting plate supplied by contractor (8 ea. required).
2. The mounting structure must be attached to structural ceiling and **BRACED TO ALLOW NO TWISTING OR LATERAL MOTION** and shall be designed not to provide a degree of deflection greater than two-tenths of a degree at the mounting plate.
3. 3/4" metal conduit and minimum 12AWG wire size (3 wires per lighthouse plus fixture ground) required between fixture and SKYTRON supplied wall control. All metal conduit, wiring, and other electrical materials as well as installation labor for such materials associated with the installation of the SKYTRON surgical light to be provided by Electrical Contractor. All installations of SKYTRON surgical lights should be under the direct supervision of a SKYTRON representative. All wiring to be in accordance with local codes and by a certified electrician.
4. Optional Camera Ready Systems require a 1" diameter conduit from the fixture junction box to the camera control connector junction box. 65' (20M) camera control cable with faceplate and connector supplied by SKYTRON. Conduit and 2" x 4" junction box supplied by General or Electrical Contractor.
5. CONTRACTOR HAS FINAL RESPONSIBILITY for the strength and stability of the Mounting Structure.

**This is a GENERAL GUIDELINE ONLY.**

## SECTION 7. REVISION HISTORY

Date	Revision	Revision History
05/07/2013	0	Initial release.
05/28/2013	1	Added information to inside cover. Throughout the document, changed 1-1/4" conduit requirement to 1" as a correction.
02/11/2014	2	New logos, deleted "white (neutral) on lightheads" in wiring caution, added callout to wall control wiring figure, deleted separate ground note from installation drawings, corrected junction box wiring diagrams.
08/28/2014	3	Added EC/REP on inside cover and Label Symbols table on Page 3 .



5085 Corporate Exchange Blvd. S.E.  
Grand Rapids, MI 49512 • 616.656.2900 • FAX 616.656.2906