Fabius[®]- CE



Inhalation Anesthesia Machine Software 4.n Instructions for Use

Contents

For Your Safety and that of Your Patients	4		
		Fitting Microbial Filters	
Intended Use		Connecting the Breathing Hoses	
Safety Features		Connecting the O2 Sensor Capsule	
Copyright, Trademark, and Disclaimer Notices	5	Connecting the Pressure Sensor	
		Connecting the Flow Sensor	
Typical Fabius Configurations	6	Connect the Magill Inhalation Device	
		Connecting the O2 Sensor	32
Components	10	Installing Anesthetic Gas Scavenging Hose to the	
Fabius Inhalation Anesthesia Machine		Compact Breathing System	
Configurations	10	Anesthetic Gas Scavenging System AGS	33
Anesthesia Ventilator		Scavenger System Connections for the	
Compact Breathing System		Semi-open Compact Breathing System	
Semi-Open Adapter		A-Cone Switch	
Magill Breathing Circuit	11	Installation of the Semi-Open Adapter	34
A-Cone Switch		Removing the Semi-Open Adaptor and	
Dräger-Vapor® Anesthetic Agent Vaporizer		Installing the CO ₂ Absorber	
Vapor Interlock System		Additional Equipment	36
Dräger-Vapor [®]			
Selectatec®	12	Operation	37
Airway Monitor with Control Units		Nitrogen Wash-out (when Required)	
, and the state of		Adjusting the Fresh Gas Composition	
Ventilation and Monitor Operating Concept	14	Setting the Vapor	
User Interface		Ventilation Modes with the Compact and	
Rotary Control		Semi-open Compact Breathing System	38
Displays		Spontaneous Breathing	
		Manual Ventilation	
Screen Display Configuration	15	Automatic Ventilation	
Mode Selection Buttons	15	Pressure limited ventilation	
Ventilation Parameter Buttons		Using the Magill Breathing System	
		A-Cone Lever Selection	
Function Buttons		O2 Flush	
Setting Alarm Limits		Replacing CO2 Absorbent	
Calibrating the Oxygen Sensor		Replacing Used Absorbent (Compact System)	
Calibrating The Flow Sensor		p.monig bood risociponi (bompast ojoioni)	•
Menu Display Button	20	Shut-down	45
Configuration Menu		Switch Off the Anesthetic Agent Vaporizer	-10
Language Selection Menu		(Dräger-Vapor®)	۸۶
Gas Compensation Setting Menu	21	Close the Metering Valves	
.	~~	Switch off anesthesia ventilator	
Preparation		Disconnect the Fresh-gas Supply	
Caution During Transport		De-pressurizing the Anesthesia Machine	
Gas supply		Switch Off the Control Unit	
Cylinders with Threaded Connectors		Remove the O2 Sensor	
Cylinders with Pin-index Mounting		Nemove the O2 Sensor	71
Electrical Supply		Fault - Cause - Remedy	45
Equipotential Bonding	24	rauit - Cause - Remedy	4 1
Attaching Manual Ventilation Bag		Discoulling Day of	40
Preparing the Ventilator	25	Dismantling Down	45
Fitting the CO ₂ Absorber on the Compact		Dismantling the Compact Breathing System	
Breathing System	26	for Cleaning	
Fitting the Inspiratory Valve	26	Dismantling the Inspiratory Valve	
Fitting the Expiratory Valve		Dismantling the Expiratory Valve	
Fitting the Pressure-limiting Valve		Dismantling the Flow Sensor	
Inserting the Flow Sensor		Dismantling the Pressure-limiting Valve	
Fitting the Waste Gas Outlet Port	27	Dismantling the Absorbent Canister	49

Minimum Flow of Anesthesia	50
Disinfecting/Cleaning Wipe Disinfecting Disinfecting and Cleaning Automatic Cleaning/Disinfecting Machine Sterilizing	51 51 52
Checking Readiness for Operation	52
Maintenance Intervals	53
Disposal	53
What's what Compact Breathing System, Front Compact Breathing System, Back	54
Technical Data	56
Diagrams	62 63
Abbreviations Used	65
Daily Pre-Use Check for Fabius	66
Index	71

For Your Safety and that of Your Patients

Strictly follow the Instructions for Use

Any use of the apparatus requires full understanding and strict observation of these instructions.

The apparatus is only to be used for purposes specified here.

Observe the Instructions for Use of the supplementary equipment used)

Maintenance

The apparatus must be inspected and serviced by trained service personnel at six month intervals.

Repair and general overhaul of the apparatus may only be carried out by trained service personnel.

We recommend that a service contract be obtained with DrägerService and that all repairs also be carried out by them.

Only authentic Dräger spare parts may be used for maintenance.

Observe chapter "Maintenance Intervals".

Accessories

Do not use accessory parts other than those in the order list (86 06 304, revision index 00).

Even reusable accessories (e.g. parts which can be cleaned and/or sterilized) have a limited useful life. Wear and tear may be increased and the service life reduced considerably as a result of various factors when handling and conditioning such parts (e.g. disinfectant residues left behind after autoclaving may corrode the material). Such parts must be replaced immediately if external signs of wear become evident, such as cracks, deformation, discoloration, peeling, etc.

Not for Use in Areas of Explosion Hazard

This apparatus is neither approved nor certified for use in areas where combustible or explosive gas mixtures are likely to occur.

Safe connection with other electrical equipment

Electrical connections to equipment which is not listed in these Instructions for Use should only be made following consultations with the respective manufacturers or an expert.

Liability for proper function or damage

The liability for the proper function of the apparatus is irrevocably transferred to the owner or operator to the extent that the apparatus is serviced or repaired by personnel not employed or authorized by DrägerService or if the apparatus is used in a manner not conforming to its intended use Dräger cannot be held responsible for damage caused by non-compliance with the recommendations given above. The warranty and liability provisions of the terms of sale and delivery of Dräger are likewise not modified by the recommendations given above.

Dräger Medical AG & Co. KG

Intended Use

Fablus is an inhalation anesthesia machine with continuous fresh gas flow for patients with a body weight of at least 5 kg.

Fabius is intended for use in operating, induction and recovery rooms.

It may be used with O₂, N₂O, and AIR supplied by a medical gas pipeline system or by externally mounted gas cylinders. Fabius can be equipped with a compact breathing system, providing fresh gas decoupling, PEEP, and pressure limitation. The following ventilation options are available:

- Automatic Ventilation (IPPV),
- Pressure limited ventilation (IPPV/PLV),
- Manual Ventilation (MAN),
- Spontaneous Breathing (SPONT).

The Fabius can be equipped with an electrically operated and electronically controlled ventilator, which monitors the airway pressure (Paw), the volume (V) and the inspiratory oxygen concentration (FiO2).

According to EN 740 (Anaesthesia Workstations and Their Modules - Particular Requirements) the airway pressure (Paw), the volume (V) and the oxygen concentration (FiO2) must be additionally monitored in the case of a Fabius not equipped with the ventilator.

As per EN 740 (Anesthetic Workstations and their Modules – Special Requirements), additional monitoring of the concentrations of CO2 and anesthetic agent is required.

EN 740 requires that a manual ventilation bag must be available for emergency use.

Fresh gas enrichment is provided by the Dräger-Vapor® anesthetic vaporizer.

Do not use readily flammable anesthetic agents such as ether, cyclopropane, etc.

Safety Features

- Monitoring of P, V, FiO2
- Automatic O₂ shortage alarm and N₂O cut-off
- Integrated S-ORC = Sensitive Oxygen Ratio Controller (control device to ensure minimum O2 concentration of 21 vol.%).

Burns may occur if antistatic or electrically conductive ventilation tubes are used in combination with high-frequency electrical surgery equipment. Therefore, these types of breathing tubes are not recommended.

Caution: Do not use Fabius in the environment of NMR tomography equipment. Malfunctions may result, thereby endangering the patient.

Caution: Do not use mobile phones within a distance of 10 meters from the machine. Mobile phones can cause malfunctions in electrical medical equipment, thereby endangering the patient.

Copyright, Trademark, and Disclaimer Notices

Copyright

Copyright by Dräger Medical, Inc. All rights reserved. No part of this publication may be reproduced, transmitted, transcribed, or stored in a retrieval system in any form or by any means, electronic or mechanical, including photocopying and recording, without written permission of Dräger Medical, Inc.

Registered trade marks

- DrägerService[®]
- Fabius[®]
- Dräger-Vapor[®]
- Drägersorb[®] 800 Plus

are registered trade marks of Dräger Medical AG & Co. KG.

Selectatec[®]

is a registered trade mark of Datex-Ohmeda.

Typical Fabius Configurations

General Description

The Fabius Inhalation Anesthesia Machine is a modular system consisting of a basic gas-delivery module with a variety of components and configuration designs to meet the requirements of various anesthesia delivery applications.

Versions

2-gas version:

O2, N2O

3-gas version:


O2, N2O, Air

Optionally available with

two cylinder connections:


O2, N2O

Wall mounted Fabius with single vaporizer and Magill breathing system (2-gas version), (illustrated with vaporizer holder).




Wall mounted Fabius with dual vaporizers, compact breathing system, ventilator, and control box (3-gas version).

Alternatively: mounted on CE trolley (not illustrated).




Wall mounted Fabius with dual vaporizers, compact breathing system, ventilator, control box, and PM 8050 gas monitor (3-gas version).


Alternatively: mounted on CE trolley (not illustrated).




Pin index trolley mounted Fabius with single vaporizer and Magill breathing system (2-gas version).




Pin index trolley mounted Fabius with dual vaporizers, compact breathing system, ventilator, and control box (3-gas version).




Pin index trolley mounted Fabius with dual vaporizers, compact breathing system, ventilator, control box, and PM 8050 gas monitor (3-gas version).



CE trolley mounted Fabius with dual vaporizers, compact breathing system, ventilator, control box, A-cone switch, and Magill breathing system (2-gas version).




Ceiling mounted Fabius with compact breathing system, dual vaporizers, ventilator, and control box (3-gas version).



Components


Fabius Inhalation Anesthesia Machine Configurations

- anesthesia ventilator
- compact breathing system
- Dräger-Vapor[®] anesthetic agent vaporizers
- airway monitor, complete with control unit




Anesthesia Ventilator

The Fabius can be equipped with a ventilator for controlled ventilation of adults and children during anesthesia. The ventilator is electrically-powered, time-controlled, and volume-constant.




Compact Breathing System

Re-breathing system with CO2 absorption, optimized for a reduced fresh gas flow, with direct switchover from spontaneous breathing mode to manual mode at the pressure limiting valve. Maximum pressure limiting is possible during automatically controlled ventilation and manual ventilation. A PEEP function is available for automatically controlled ventilation. The inspiratory volume is independent of the fresh gas flow since the system is fresh gas decoupled.



Semi-Open Adapter


The semi-open compact breathing system is configured with a semi-open adapter to perform as a semi-open system with no rebreathing. The breathing system is used in the same way as the compact breathing system except that the CO2 absorbent is not used. The fresh gas flow rate must be adjusted higher than the patient's minute volume.






Magill Breathing Circuit

The Magill Circuit is a semi-open breathing system and can be used in combination with a semi-closed breathing system or instead of a semi-closed breathing system. The Magill breathing system makes no allowances for pressure or O2 monitoring.



A-Cone Switch

The A-Cone switch provides a means for switching the fresh gas flow between two separate patient breathing circuits without changing fresh-gas hose connections. Switching the flow is accomplished by moving the lever, located on top of the A-cone switch, to the desired flow position. This lever controls a spring loaded switch that has a positive set position in either the right-hand or the left-hand side of the A-cone switch. This spring loaded feature ensures positive sealing for the fresh gas flow to the breathing circuit selected.




Dräger-Vapor® Anesthetic Agent Vaporizer

The Dräger-Vapor® is used to enrich the fresh gas with a precisely metered quantity of vapor from the liquid anesthetic agent being used, i.e. Isoflurane, Halothane, Enflurane, or Sevoflurane.

When using Desflurane:

230 V Mains	Devapor*, D-Vapor
110 V Mains	D-Tec*, D-Vapor

^{*} Devapor and D-Tec are available through your local Desflurane representative.




Vapor Interlock System

If the mounting for two separate vaporizers is used, an interlock system is provided to ensure that only one vaporizer can be used at a time.

Dräger-Vapor®

Note that the selector lever is shown in the center position. This ensures that both vaporizers are in the locked position. Also, this is the recommended position for the selector lever when moving the Fabius.

1 Moving the selector lever away from the desired vaporizer allows that vaporizer to be utilized and the other to be locked out of use.



Selectatec[®]

The interlock system for the Selectatec[®] is built into the vaporizers. When a vaporizer is selected for use, the interlocking index pins will protrude from the sides of the vaporizer thereby not allowing the neighboring vaporizer to be opened. For more specific information on the Selectatec[®], refer to the Selectatec[®] Vaporizer's Instructions for Use.

Airway Monitor with Control Units

Measuring and monitoring device for ventilation parameters during anesthesia.

The following parameters are displayed:


- Airway pressure (Paw), Peak, PEEP, Pmean
- Expiratory minute volume (MV), tidal volume (VT), breathing frequency (f)
- Inspiratory oxygen concentration (FiO2)

The following parameters are used for monitoring:

- Airway pressure (Paw)
- Inspiratory oxygen concentration (FiO₂)
- Expiratory minute volume (MV)

The monitor comes with a control unit for setting and adjusting all ventilation parameters as outlined in the table below.

Maximum pressure (Pmax)	10 to 70 hPa
Tidal volume (VT)	50 to 1400 mL
Breathing frequency (fIPPV)	6 to 60 L/min
Insp./Exp. time ratio (TI: TE)	1:3 to 2:1
Inspiratory pause (TIP: TI)	5 to 50 %
Positive end expiratory pressure (PEEP)	0 to 15 hPa



Ventilation and Monitor Operating Concept

User Interface

Both the monitor and the ventilator provide a user interface by way of

- buttons
- rotary controls
- displays, and
- warning tones.

Rotary Control

The rotary control is used to select functions by:

... turning

Turning the rotary control moves the cursor (on the screen) or changes the value displayed on the screen.

... pressing

By pressing the rotary control, the selected value is confirmed, or a selection process is started or ended.

Displays

All necessary information is displayed on the screen. Additional LED indicators within buttons display the operating mode or status.

LED lamps in the upper right corner of the control panel indicate the highest degree of urgency of currently active alarms.

- Warning Red Blinking
- Caution Yellow Blinking
- Advisory Yellow Steady

Warning Tones

The warning tones provide an audible alert to the message displays. Each message is assigned a tone or sequence of tones to indicate its degree of urgency.

- 49 --- ---- warning (continuous)
- ← → → → → → ← caution (every 30 seconds)

 Output

 Description

 Output

 Description

 Description

 Output

 Description

 Description
- → advisory (single signal)

Screen Display Configuration

1 Mode:

STANDBY

IPPV

MAN/SPONT SELFTEST

2 Power source:

mains battery

3 Messages:

Warning

Caution

Advisory

4 Specific screen menus

1	2	3
Mode	Ð	Messages
	4	

Mode Selection Buttons


The mode is selected by one of these buttons and confirmed by pressing the rotary control.


- For switching to manual ventilation or spontaneous breathing.
- 6 Button for IPPV.

When the user selects the IPPV mode, by pressing the IPPV button on the front panel, an additional screen appears after the mode change has been confirmed. It questions the user if the APL valve is in the IPPV position. The APL valve must be switched to the IPPV position and this question must be confirmed before the ventilator begins running.

If the user fails to confirm that the APL valve is in the IPPV position, the mode remains as it was before the user pushed the IPPV button on the front panel. This action to select the IPPV mode is treated as a CAUTION, therefore yellow LEDs are momentarily lit and an audible caution is heard.

7 Standby button.






Ventilation Parameter Buttons

Buttons located below the display window:

- 1 Button for setting the maximum pressure for ventilation (IPPV).
- 2 W Button for setting the tidal volume.
- 3 Button for setting the ventilation frequency.
- 4 Button for setting the time ratio between inspiration and expiration phases.
- 5 Button for setting the relative inspiratory pause.
- Button for setting a positive end expiratory pressure (PEEP) during automatic ventilation in the IPPV mode.




Function Buttons

1 Alarm Silence Button


For muting the alarm beep for two minutes. The yellow LED in this button lights up while the audible alarms are suspended.

Any new messages are signalled by the appropriate tone sequence.




2 🙆 Standard Screen Button

Selects the standard screen. The ventilation pressure is graphically indicated in a bar graph format, complete with alarm limits.




3 Additional Screens Button Selects additional screens.

Paw curve and measured values screen



Ventilator settings screen




Currently active alarms screen


IPPU			(1 -)	MV LOW	1
60 40 20	V _{Te} f P _{mean} PEEP	330 12 7 3	ml 1/min hPa hPa	BATTERY POWERED Fi02 LIM < 21%	03037441

Example: Setting Ventilator Parameter Values Using Pmax as an Example.


1 Press 🛌 button.



The current setting in this example is »25«. This unchanged value is displayed in the bottom line.



2 Set the new value by turning the rotary control. In this example, the new value is »30«. The unchanged original value of 25 is still displayed in the bottom line.



 $T_{\mathbf{I}}:T_{\mathbf{E}}$

fippv

 $T_{IP}:T_{I}$

PEEP

[hPa]

0

Confirm the new setting by pressing the rotary control.

The previous value disappears from the bottom line and the new value in the Pmax box immediately becomes valid.

[ml][1/min] [2] [hPa] 10 30 600 12 1:2.0Ventilator Settings

IPPU


 P_{max}

 V_{T}

If the rotary control is neither pressed nor turned within 10 seconds, the device quits the setting mode and the original settings remain unchanged.

Setting Alarm Limits


- 1 Press (**) Alarm Limit button.
- 2 Select the required limit by turning the rotary control.
- 3 Confirm by pressing the rotary control.



Example:

Lower alarm limit for inspiratory O2 concentration.


- 2 To modify the setting, turn the rotary control.
- 3 To confirm the new setting, press the rotary control.



Calibrating the Oxygen Sensor

The oxygen sensor can be calibrated with 21 % O₂ or with 100 % O₂.


- 4 Switch the control unit to standby.
- 5 Press Sensor Calibration Button and select the required calibration.
- 6 Confirm the selected mode by pressing the rotary control.



To calibrate the O_2 sensor, expose it to the appropriate calibration gas.

- Calibrating at 21 vol.% expose the sensor to the ambient air.
- Calibrating at 100 vol.% ~ expose the sensor for at least 2 minutes to 100 vol.% oxygen.

Note: Sensor calibration takes approximately 2 minutes.



Calibrating The Flow Sensor

(only in Standby)

When calibrating the flow sensor, make sure that no flow is present at the sensor. Flow sensor calibration takes 1 second.

Symbols for Status of Calibration:

- 🗸

calibrated

~ (<u>°</u>


calibration running

- 0

calibration failed, no measurement

- no symbol

measurement uses old calibration



Note: If the optional PM 8050 gas monitor (SW 2.05) is attached to the Fabius, the buttons for setting the alarm limits and for calibrating the O2 and flow sensor are no longer used. These functions are now controlled by the PM 8050 gas monitor. For more information on the PM 8050 gas monitor, please refer to the specific instructions for Use.

Menu Display Button

1 Menu Button

Battery Status Indications:

FULL: E

Battery is fully charged.

CHARGING:

Battery is charging from mains power.

POWERING:


Fabius is operating on battery.

TESTING:

Battery test in progress.

LOW:

Battery is discharged.



Configuration


Configures machine (select and confirm by turning and pressing the rotary control).


Language

Select language (select and confirm by turning and pressing the rotary control).

Gas Compensation Setting

Select nearest fresh-gas mixture ranges.





Configuration Menu

Key Code:

For DrägerService only

Display timeout: After this time, the standard screen is automatically set.

Power on setting: LAST (last setting used) or

STRD (standard defaults)

In standard mode (STRD), the following parameter settings are active when the machine is switched on:

Pmax = 25 hPa

VT = 600 mL

fIPPV = 12/min

 $T_1:T_E = 1:2$

TIP:TI = 10 %

PEEP = 0 hPa

FiO2 high = ---

FiO₂ low = 21 %

MV high = 12 L/min

MV low = 3 L/min

Paw high = 40 hPa

Paw low = 8 hPa

Curve IPPV:

The pressure curve may be free running

(FREE) or synchronized (SYNC)

with the breathing cycle.

Display invert:

toggles the screen background between

light and dark (activation in standby mode).

Language Selection Menu

Main language:


Alphabetic languages for displaying

messages and on-screen information.

Add. language:

Messages in non alphabetic

languages (alarm messages).




Configuration Language Language Gas compens. setting

IPPU

₫ FULL \varTheta Last warning No entry

Gas Compensation Setting Menu

- To adapt flow measurement to the gas mixture currently in use, select »Gas compens. setting« and confirm.
- With »Gas mix. O2/N2O«, select the gas composition that comes closest to the fresh gas setting.
- Confirm the selected mixture. If the selected gas composition is not close to the actual gas composition, flow measurement data will be inaccurate.



Ð

Preparation

Caution During Transport

Prior to moving the Fabius, remove monitors and extraneous equipment from the top of the unit that is not permanently fastened to the unit. Take care when rolling the machine over door thresholds and similar obstacles.

Gas supply

Medical gases must be dry and free from dust and oil.

 Rear view of a 3-gas CE trolley version with five gas connections.


Medical Gas Pipeline Supply of O2, N2O, and AIR

- 1 Connect O2 hose to connector on Fabius and
- 4 to wall terminal unit of medical gas pipeline system.
- 2 Connect N2O hose to connector on Fabius and
- 4 to wall terminal unit of medical gas pipeline system.
- 3 Connect AIR hose to connector on the Fabius and
- 4 to wall terminal unit of medical gas pipeline system.

Caution When Handling O2 Cylinders

- Do not oil or grease the O2 cylinder valves and O2 pressure regulator. Do not touch these areas with greasy fingers. There is a risk of explosion.
- If cylinder valves are leaky or difficult to open or close, they
 must be repaired in accordance with the manufacturer's
 specifications.

Even if the gas supply is connected to a medical gas pipeline, the cylinders should remain on the device in reserve.



Cylinders with Threaded Connectors

On the back of the unit:


- Place the full cylinders in the cylinder holders and secure them in position.
- Screw the pressure regulators on to the cylinder valves.
- Screw the compressed gas hoses to the pressure regulators and to the lower connections of the gas inlet block.
- Open the cylinder valves.

The cylinder valves must only be opened or closed manually. Never use tools.

Cylinders with Pin-index Mounting

When attaching a cylinder, ensure that only one washer is installed between the cylinder and the yoke gas inlet. The use of multiple washers will inhibit the pin-index safety system. Be sure to verify the presence of the index pins each time a cylinder is installed. Never attempt to override the pin-index safety system.


- 1 To connect a gas cylinder to its yoke:
- 2 Remove the old washer and install a new washer on the seat of the yoke gas inlet connection.
- 3 Verify that the two index pins below the
- 4 gas inlet are present.
- 5 Insert the head of the gas cylinder into the yoke from below.
- Ensure that the gas outlet and indexing holes on the cylinder head align with the gas inlet
- 6 and index pins of the yoke assembly.
- Engage the indexing holes with the index pins.
- 7 Turn the yoke handle clockwise against the cylinder head, so that the point of the yoke handle bolt is aligned with the indent on the back of the cylinder head.
- Verify that the washer is in place, the index pins are engaged, and the cylinder hangs vertically. Tighten the yoke firmly.
- 8 When required, the cylinder valve
- 9 is opened using the cylinder wrench that is provided.
- 10 When a cylinder is removed, place the yoke plug in the yoke assembly and tighten.



Electrical Supply

Fabius can be operated at mains voltages from 100 V to 127 V. Observe the label beside the mains voltage input socket!


- Push power plug into main supply socket.
- Switch on the machine. The on/off switch is on the rear of the machine.




Equipotential Bonding

For intra-cranial and intra-cardiac operations Trolley-mounted anesthesia machines:

- 2 Use cable ref. 83 01 349 cable.
- Connect the terminal on the back of the trolley to an equipotential bonding point in the operating room.




Pin Index version



Attaching Manual Ventilation Bag

e.g., Resutator 2000 for emergency ventilation

 Hang the fully prepared and tested bag on the rail at the right.




Preparing the Ventilator

Use only disinfected/sterilized components.

- 1 Swing out the ventilator door.
- Unlock the clasps to
- 2 remove the cover.
- 3 Insert the diaphragm.
- 2 Fit the cover and
- 4 lock both locking screws.
- 1 Swing the ventilator unit back into position.


Ventilator Safety Feature

- A Over pressure valve
- B Additional air valve



Fitting the CO₂ Absorber on the Compact Breathing System

- Fill the absorber with fresh CO2 absorbent. Dräger recommends the use of Drägersorb® 800 Plus.
- Tighten the absorber by turning it to the right into the compact breathing system.
- Ensure that no CO2 absorbent dust/particles have been deposited between the gaskets and the sealing surfaces.
 Such dust and particles can cause leaks in the system.




Fitting the Inspiratory Valve

- 1 Place the valve disc in the valve seat.
- 2 Fit the inspection cap (with port).
- 3 Tighten the retaining nut securely.

Fitting the Expiratory Valve


- 4 Place the valve disc in the valve seat.
- 5 Fit the inspection cap.
- 6 Tighten the retaining nut securely.



Fitting the Pressure-limiting Valve

7 Tighten the pressure-limiting valve securely into place with the retaining nut.

Only use a valve marked with "SPONT/IPPV" and "MAN".




Inserting the Flow Sensor

- 1 Unscrew and remove the expiration port.
- Insert the flow sensor, observing the indicated direction of flow.


Fitting the Waste Gas Outlet Port

3 Screw the waste gas port into the compact breathing system from underneath.




Connecting the Compact Breathing System

- The sealing rings on the threaded and conical connectors must be undamaged and clean.
- Only hand-tighten the threaded connectors. Do not use tools.
- Slightly loosen the knurled screw on the compact breathing system.
- 4 Fit the compact breathing system to the pin on the anesthesia machine.
- Tighten the knurled screw.



- 5 Screw the fresh gas hose to the Fabius and to the compact breathing system.
- 6 Screw the ventilation hose to the ventilator and to the compact breathing system.
- 7 Plug the control hose to the connection port on the expiration valve and the connection port on the ventilator.



Fitting Microbial Filters

 Push the microbial filters into the inspiratory port and expiratory port, until they noticeably click into place.


Dismantling instructions:

 Squeeze at the points marked "PRESS" and simultaneously pull the filter off the port.

Condensation may increase the flow resistance of the filters, making breathing more difficult.

Carefully monitor the airway pressure.


Observe the specific Instructions for Use for the microbial filter.




Connecting the Breathing Hoses

Note:

- Take care not to damage the breathing hoses.
- When connecting and disconnecting, always hold the breathing hoses by the end sleeve, not by the spiral reinforcement. Otherwise, the spiral reinforcement may be torn loose, e.g. at the end-sleeve.
- Breathing hoses with damaged spiral reinforcement can kink, thereby interrupting breathing.
- Before each use, check the breathing hoses for damage.




- 1 Push a patient breathing hose on to both the inspiratory cone and the expiratory cone or on to the microbial filters.
- 2 Connect both patient breathing hoses to the Y-piece.
- 3 Connect the manual breathing bag to the bag breathing hose via the double port.
- 4 Connect the bag hose to the elbow port on the compact breathing system.




Connecting the O2 Sensor Capsule

- Push the O2 sensor into the port of the inspiratory valve, and
- 2 plug in the connector on the back of the control unit.




Connecting the Pressure Sensor

- 3 Press the probe of the pressure measuring line into the socket on the underside of the compact breathing system until it engages.
- 4 Connect the hose of the pressure measuring line to the bacterial filter and plug it firmly into the port on the back of the control unit.



Connecting the Flow Sensor


1 Push the measuring hoses on to the connection ports on the flow sensor and at the back of the unit. Please note the connections are different.



Connect the Magill Inhalation Device

For use with mask:


- 1 Screw Magill connector onto the freshgas outlet.
- 2 Connect breathing bag.
- 3 Connect breathing hose E (110 cm).
- 4 Plug ISO mask connector into breathing hose.
- 5 Connect mask to ISO mask connector.
- 6 Connect sampling hose for monitoring to the Magill connector.



For use with tube:

- Remove ISO mask connector.
- Use ISO connector.
- Plug connector for tracheal tube into ISO connector.


Caution: The Magill Circuit shall not be connected to the Fabius ventilator.



Connecting the O2 Sensor


Inserting a new O2 sensor capsule:

- 1 Unscrew the screw cap from the sensor housing.
- Remove the new sensor capsule from its packaging, or use a disinfected sensor capsule.
- 2 Insert the capsule in the housing, with the ring-shaped conductor against the contacts in the housing.
- 1 Screw the screw cap on firmly by hand.



Installing Anesthetic Gas Scavenging Hose to the Compact Breathing System


- Connect the transfer hose to the waste gas port of the Compact Breathing System and to the anesthetic gas scavenging line or an anesthetic agent filter.
- A second transfer hose is required for the Semi-open compact breathing system.



Anesthetic Gas Scavenging System AGS

- 1 Output connection from the scavenging system.
- Connection to scavenger system from Fabius waste gas outlet port.
- 3 Flow indicator. During use, the flow indicator must be between the upper and lower marks on the tube.
- 4 Connection for Semi-open breathing system.

For more detailed information on the scavenging system, refer to the specific Instructions for Use.




Scavenger System Connections for the Semi-open Compact Breathing System

Both exhaust ports


- 5 one on the semi-open adapter and
- 6 the other on the Compact Breathing System housing must be connected to the AGS scavenger.

Remove the existing scavenger plug if necessary.




A-Cone Switch

- 7 The A-Cone switch is connected to the fresh gas outlet via the conical fitting. It must be mounted parallel to the front of the Fabius.
- 8 The horizontal port of the A-Cone switch is connected to the alternate breathing system.
- 9 The lower vertical port of the A-Cone switch is connected to the Compact or Semi-open Compact breathing system on the Fabius.




Installation of the Semi-Open Adapter


- Disconnect the Fabius from the mains power supply.
- Disconnect the Fabius from the central pipeline gas supply.
- Close all gas cylinders (if applicable).
- Remove all hoses, sensors, and control lines from the Compact Breathing System.
- 1 Remove the APL valve.
- 2 Remove the absorbent canister and store properly.
- 3 Loosen the locking screw that secures the Compact Breathing System to the mounting stud.
- Gently lift the breathing system and place upside down on a firm surface. It is recommended that it be placed on a soft surface such as a towel to prevent marring of the unit.



- 4 Remove the three mounting screws (M 5 x 16 mm) and washers that hold the canister mount to the compact breathing system housing.
- 5 Ensure that all "O" rings are removed with the canister
 - Store this mount assembly, hardware, and ${}^{\tt m}O^{\tt m}$ rings with the absorbent canister.




- Prepare the semi-open adapter.
- 1 Ensure that the "O" rings are in the proper position. These "O" rings are provided with the adapter. Do not use the "O" rings from the compact canister mount, as they are not interchangeable.
- Place the adapter onto the compact breathing assembly and
- 2 secure with the three screws provided (M 5 x 80 mm) with the adapter.
- 4 These screws each come with a flat washer and
- 3 and four belleville washers. The belleville washers go on first followed by the flat washer. Please note that the belleville style washer is a curved spring washer and must be installed by opposing each other. Do not over tighten these screws.
- Place the semi-open adapter back on the mounting pole of the Fabius and secure the locking screw.
- Connect all hoses, sensors, and control lines.
- Install the APL valve.
- Connect Fabius to the mains powers supply, and central pipeline gas supply.




Removing the Semi-Open Adaptor and Installing the CO₂ Absorber

- Disconnect the Fabius from the mains power supply.
- Disconnect the Fabius from the central pipeline gas supply.
 Close all gas cylinders (if applicable).
- Remove all hoses, sensors, and control lines from the Compact Breathing System.
- Remove the APL valve.
- Loosen the locking screw that secures the Compact Breathing System to the mounting stud.
- Gently lift the breathing system and place upside down on a firm surface.
 It is recommended that it be placed on a soft surface such
 - as a towel to prevent marring of the unit.


 Remove the three mounting screws (M.5 x 80 mm) and
- 5 Remove the three mounting screws (M 5 x 80 mm) and washers that hold the semi-open breathing system adapter to the compact breathing system housing.



- 1 Ensure that all five "O" rings are in the proper position of the canister mount.
- Install the three mounting screws (M 5 x 16 mm) and washers that hold the canister mount to the compact breathing system housing.
- 3 There are three belleville washers and
- 4 a flat washer for each mounting screw. Please note that the belleville style washer is a curved spring washer and must be installed by opposing each other. The flat washer is installed next. Do not over tighten these screws.



- 5 Gently lift the breathing system and place it on the mounting stud. Tighten the mounting screw.
- Connect all hoses, sensors, and control lines.
- 6 Install the APL valve.
- Connect Fabius to the mains powers supply, and central pipeline gas supply.
- 7 Install the absorbent canister. Ensure it is filled with fresh CO2 absorbent.



Additional Equipment

Prepare additional equipment as specified in the specific Instructions for Use.

Caution: If monitors and other equipment are placed on top

of Fabius, the risk of tipping over the unit is increased, especially when rolling over thresholds etc. Remove all monitors and other equipment from the top of the Fabius before moving the unit.

Caution: Remove all monitors and other equipment from the

top of the Fabius before moving the unit.

Operation


Nitrogen Wash-out (when Required)

During anesthesia induction, air containing about 79 % nitrogen (N2) remains in the compact breathing system (and in the patient's lungs).

When the next intended use of the unit may be for a low-flow anesthesia case, the nitrogen content should first be flushed with a time-limited high flow of fresh gas (O2/N2O) in order to ensure correct denitrogenization.


Adjusting the Fresh Gas Composition

With the metering valves:
 Set fresh gas flow for O2, N2O and AIR.



Setting the Vapor

- 2 Lock the unused Vapor by pushing the lever as far as it will go in the direction of the unused Vapor (in this example left hand Vapor is locked).
- 3 On the Vapor to be used, press and hold down the 0 button and at the same time
- 4 Turn the handwheel counter-clockwise to the desired anesthetic agent concentration.
- Regularly check the filling level on the sight glass. When reaching the minimum mark, fill the Vapor with anesthetic agent.
- Please refer to the specific Instructions for Use for Dräger-Vapor[®].



Ventilation Modes with the Compact and Semi-open Compact Breathing System


There are three ventilation modes available in the Fabius:

- Spontaneous
- Manual and
- Automatic


Spontaneous Breathing

On the control unit

- 1 Press the www button.
- 2 Confirm by pressing the rotary control.



- 3 Set the lever on the pressure-limiting valve to SPONT/IPPV. The valve is now open, regardless of the set pressure limit. Pressure cannot build up in the compact breathing system.
- 4 Administer enough fresh gas to fill the manual ventilation bag on the compact breathing system.



To prevent false alarms, the monitoring function during spontaneous breathing can be reduced to the following parameters:


- Lower alarm limit for O2
- Upper alarm limit for Paw

If other alarm limits are required, see "Setting Alarm Limits" on page 19.


Manual Ventilation

On the control unit

- 1 Press the San button.
- 2 Confirm by pressing the rotary control.



- 3 Set the lever on the pressure-limiting valve to MAN.
- 4 Set the desired pressure limit by rotating the lever until the desired value is indicated on the scale.
- Administer sufficient fresh gas.
- 5 Ventilate manually with the manual ventilation bag,



To prevent false alarms, the monitoring function during manual ventilation can be reduced to the following parameters:

- Lower alarm limit for O2
- Upper alarm limit for Paw

If other alarm limits are required, see "Setting Alarm Limits" on page 19.


Automatic Ventilation


On the control unit set the following controls to the desired values:

1 Pmax; VT; f; Ti:TE; TIP:Ti and PEEP

Then

- 2 Press the button.
- 3 Confirm by pressing the rotary control.
- The IPPV Confirmation Request screen appears.






- 4 Set the pressure-limiting valve to SPONT/IPPV.
- 3 Confirm by pressing the rotary control.

Note:


If the user fails to confirm that the APL valve is in the IPPV position, the mode remains as it was before the user pushed the IPPV button on the front panel. This action to select the IPPV mode is treated as a CAUTION, therefore yellow LEDs are momentarily lit and an audible caution is heard.

Administer sufficient fresh gas.



Pressure limited ventilation (PLV – Pressure Limited Ventilation)

The Fabius CE features an adjustable pressure limit. The adjustable pressure limit Pmax applies for the expiration valve. At Pmax, pressure limited ventilation is possible in IPPV mode. Ventilation is then volume inconstant.





At Pmax, pressure limited ventilation is possible (PLV).

- Set Pmax to the required inspiration pressure, e.g. 15 mbar.
- Set the tidal volume VT to a value higher than the VT corresponding to the patient's weight.

Recommendation: Set the VT to e.g. double the value of the tidal volume actually applied.

As soon as Pmax is reached in the active flow phase (TI – TIP), the pressure remains constant and leaks, e.g. due to an unblocked tube are compensated. At the same time, volume still present in the ventilator is directed into the ventilation bag. When making the settings, ensure that the inspiration flow is not excessively high, resulting in a dynamic overshoot of Pmax. The inspiration flow is influenced by the parameters VT; f; TI:TE; TIP:TI.





Recommended parameter settings

Compliance patient (mL/mbar)	fippv	T::TE	TIP:TI (%)	Set Pmax (mbar)	PEEP (hPa)	Set VT (mL)	Effective Pmax (mbar)	Exp. volume (mL)
50	10	1:2	10	15	5	1000	19	440
50	10	1:2	5	15	5	1000	19	440
50	10	1:1	10	15	5	1000	17	390
50	10	1:1	5	15	5	1000	17	390
5	25	1:2	10	15	5	100 to 150	15	40
5	25	1:2	5	15	5	100 to 150	15	40
5	25	1:1	10	15	5	100 to 150	15	40
5	25	1:1	5	15	5	100 to 150	15	40

Depending on the ventilator settings (VT; f; Ti:TE) and patient compliance, a significant overshoot of the set Pmax limit is possible at the expiration valve.

Set the alarm limit "Paw High" to Pmax + 5 mbar !

Using the Magill Breathing System

Note:

With the Magill system, anesthetic gases escape via

1 the relief valve.

This can be a health hazard!

- Ensure adequate ventilation.
- Set the O2/N2O fresh gas flow on the metering valves on the Fabius.

The fresh gas flow must be equal to at least 2.5 to 3 times the minute volume in order to exclude the risk of rebreathing and the associated accumulation of CO₂ in the inspiratory gas.

A lower fresh gas flow could result in a partial rebreathing or shortage of breathing gas.

- 2 The breathing bag could deflate and suffocate the patient.
- Switch on the vaporizer. Press and hold down the release button and select the desired concentration on the handwheel. Refer to the Instructions for Use of the Dräger Vapor[®].

2

For Spontaneous Breathing

- 1 Turn the relief valve counterclockwise (fully open). This will result in minimum ventilation pressure.
- 2 Observe movement of the breathing bag to check for spontaneous breathing.


For Manual Ventilation

- 1 Set the relief valve.
- Set a correspondingly high fresh gas flow in order to minimize the amount of rebreathed expiratory gas.
- 2 Ventilate via the breathing bag.

Caution: The Magill Breathing System shall not be connected to the Fabius ventilator.


A-Cone Lever Selection

- Facing the A-Cone switch, when the lever (located on the top of the A-Cone switch) is to the right, the fresh-gas flow is towards the lower vertical port.
- When the lever is moved to the left position, the fresh-gas flow is towards the horizontal port.
- The selector lever for the A-Cone is spring loaded and ensures that the switch is either in the right or left position thereby ensuring the flow of fresh gas to only one port at a time.



O₂ Flush

- 1 Press the O2 Flush button.
- About 50 L/min O2 flows into the compact breathing system and to the patient. The flowmeter block and the anesthetic agent vaporizer (Vapor) are bypassed.
- Make sure that the airway pressure does not rise unacceptably high.



Replacing CO₂ Absorbent


- The CO2 absorbent on the compact breathing system must be replaced before two-thirds of the CO2 absorbent has changed color.
- Dräger recommends the use of Drägersorb[®] 800 Plus. The color change indicates that the CO₂ absorbent can no longer absorb CO₂ (Drägersorb[®] 800 Plus changes from white to violet).
- Do not flush CO2 absorbent for long periods with dry gas because the CO2 absorbent will dry out.

When the moisture content falls below a specified minimum level, the following undesirable reactions can occur, regardless of the type of CO2 absorbent and the anesthetic agent used, e.g. Halothane, Enflurane, Isoflurane, Sevoflurane or Desflurane:

- reduced CO2 absorption
- formation of CO
- absorption and/or decomposition of the inhalation anesthetic agent
- increased heat generation in the absorber, leading to higher breathing gas temperatures.


These reactions can result in danger to the patient in the form of CO intoxication, insufficient depth of anesthesia and airway burns.

Please refer to the specific Instructions for Use for Drägersorb[®] 800 Plus.



Replacing Used Absorbent (Compact System)


- Empty the spent CO₂ absorbent in the absorber. Dispose
 of spent soda lime (CO₂) in accordance with the local
 waste disposal regulations.
- Fill the absorber with fresh CO2 absorbent.
- Ensure that no CO2 absorbent dust/particles have been deposited between the gaskets and sealing surfaces. Such dust and particles can cause leaks in the system.
- Dräger recommends the use of Drägersorb[®] 800 Plus.



Shut-down

Switch Off the Anesthetic Agent Vaporizer (Dräger-Vapor®)

Turn the handwheel to 0.
 The button engages.




Close the Metering Valves

- Close the N2O valve.
- Close the O2 valve and AIR valve.


N2O Pollutes the Environment

Dry fresh gas reduces the moisture content of the CO2 absorbent. Therefore, avoid unnecessary fresh gas flow.




Switch off anesthesia ventilator

- Switch the anesthesia ventilator to standby by pressing the
 (b) « button.
- 2 Confirm by pressing the rotary control. Fabius is now in standby mode.



Disconnect the Fresh-gas Supply

- Remove all plug-in couplings from the wall terminal units.
- Close gas cylinders.




De-pressurizing the Anesthesia Machine

Open the metering valves in the order listed below and hold them open until the floats in the flowmeter tubes have dropped to zero.


Proceed as follows:

- open N2O valve
- open O2 valve
- open AIR valve
- close all the metering valves.



Switch Off the Control Unit

 Switch off the unit using the switch at the back and disconnect the power plug.



Remove the O2 Sensor

 Remove the O2 sensor from the inspiratory valve and leave exposed to air. This precaution prolongs the service life of the O2 sensor.

Fault - Cause - Remedy

Display	Cause	Remedy
APNEA PRESSURE	Inadequate fresh-gas supply.	Ensure adequate fresh-gas supply.
	Breathing/ventilation stops.	Check ventilator.
	Leak in hose system.	Check hose system.
	Permanent high pressure.	Push lever on APL-valve to release the
		gas.
APNEA VOLUME	Breathing/ventilation stops – no expired volume for 15 seconds.	Check ventilator.
	Inadequate fresh-gas supply.	Ensure adequate fresh-gas supply.
	Blocked/kinked hose.	Check hose system.
	Leak in hose system.	Check hose system.
APNEA PRESS OFF	Lower limit Paw off in IPPV.	Set adequate lower limit.
BATTERY LOW	Battery is flat.	Connect mains.
BATTERY POWERED	Mains not connected.	Connect mains.
CHECK BATTERY	Backup Battery not functional.	Replace fuse.
		Inform Dräger Service.
CHECK VENTILATOR	Ventilator not assembled correctly.	Check diaphragm and close cover.
		Select Standby Mode and switch back to IPPV.
CONT HIGH PRESS	APL-Valve in Position MAN during IPPV.	Set APL-valve to IPPV/SPONT.
DEVICE TEMP HIGH	Temperature of electronic circuit too high.	Check for external heat sources. Check fan.
		Inform Dräger Service.
EX-VALVE LEAKAGE	Expiratory valve does not correctly close.	Check tubing of expiration control line.
		Calibrate Flow sensor.
		Inform Dräger Service.
FiO ₂ HIGH	O ₂ flush in use.	Check O2 setting on flowmeter.
	Inspiratory O ₂ concentration exceeds the upper alarm limit.	
FiO ₂ INOP	O2 sensor has not been correctly calibrated.	Calibrate O2 sensor.
	O ₂ sensor replaced and/or not calibrated.	Calibrate O2 sensor.
	O2 sensor used up.	Replace sensor capsule and calibrate.
	O2 sensor disconnected.	Replace O2 sensor housing assembly.
	Faulty sensor cable.	Replace sensor cable.
FiO ₂ LIMIT < 21%	Lower limit FiO2 between 18 vol.% and 20 vol.%.	Set adequate lower limit.
FiO2 LOW	Inspiratory O2 concentration is below lower alarm limit.	Check O2 supply. Check setting on flowmeter.
FLOW INOP	Sensor has not been calibrated.	Calibrate sensor.
	Sensor faulty.	Replace sensor and calibrate.
	-	Inform Dräger Service.
MEMORY ERROR	Internal data loss.	Switch machine off and on. Inform DrägerService.
		<u> </u>

Display	Cause	Remedy
MV HIGH	Upper alarm limit for minute volume has been exceeded.	
	Flow sensor not calibrated or faulty.	Calibrate flow sensor.
		Replace if necessary.
MV LOW	Minute volume has fallen below lower alarm limit.	
	Blocked/kinked hose.	Check hose system.
	Leak in breathing system.	Make breathing system leak proof.
	Reduced volume due to pressure limitation.	Correct the breathing pattern.
	Reduced lung compliance.	
	Flow sensor not calibrated or faulty.	Calibrate flow sensor.
		Replace if necessary.
Paw HIGH	Upper alarm limit for airway pressure has been exceeded, ventilation hose is kinked, stenosis.	Check hose system on anesthesia machine.
	Pressure limit has been set too high.	Correct Pmax or alarm limit value.
Paw NEGATIVE	Insufficient fresh gas supply.	Increase fresh gas flow.
POWER SUPPLY INOP	Machine fault.	Inform DrägerService.
PRESSURE INOP	Faulty sensor.	Inform DrägerService.
PRESSURE LIMIT	Set tidal volume is not completely applied.	Increase Pmax.
SET APL TO IPPV	IPPV mode has been selected.	Set APL-valve to SPONT/IPPV.
SETTINGS LOST	Setting values and calibration values lost.	Check sensor calibrations and ventilator settings.
VENTILATOR INOP	Internal machine fault.	Switch machine off and on.
	IPPV mode not operational.	Inform DrägerService.
VOL ALARM OFF	Lower limit MV off in IPPV.	Set adequate lower limit.
Only continual audible alarm without	Mains not connected.	Connect mains.
display	Battery flat.	Recharge Battery.
	Machine fault.	Inform DrägerService.

Dismantling Down

Dismantling the Compact Breathing System for Cleaning

- Leave the Dräger-Vapor[®] on the machine.
- Remove all breathing hoses.
- Remove the breathing bag and hose.
- Remove both microbial filters and prepare in accordance with the specific Instructions for Use.
- Remove the ventilation hose.
- Remove the fresh gas hose from the breathing system.
 The fresh gas hose can remain on the machine.
- · Remove the anesthetic scavenging hose.
- Detach the pressure measuring hose from the breathing system. Pull back the coupling sleeve and remove the hose and filter from the back of the control unit.
- Remove the flow measuring hoses from the back of the control unit. The measuring hoses remain on the flow sensor.
- Remove the O2 sensor lead.
- · Remove the compact breathing system.

Dismantling the Inspiratory Valve

- Unscrew the retaining nut.
- Remove the inspection cap.
- Extract the valve disc.

Dismantling the Expiratory Valve

- Unscrew the retaining nut.
- Remove the inspection cap.
- Extract the valve disc.

Dismantling the Flow Sensor

- Loosen fitting on the expiration port.
- Extract the flow sensor.

Dismantling the Pressure-limiting Valve

- Unscrew the retaining nut.
- Remove the pressure-limiting valve.
- Unscrew the waste gas outlet port.

Dismantling the Absorbent Canister


- Turn the absorber counter-clockwise and remove by pulling down.
- Remove spent CO2 absorbent. Dispose of spent soda lime (CO2) in accordance with the local waste disposal regulations.

Minimum Flow of Anesthesia

When long-term flow of anesthesia is below 0.5 L/min, increased humidity in the freshgas hose is a natural occurrence. Disconnect the freshgas hose and clean before and after long term procedures.


Magill Breathing Device

- 1 Stripping Down the Magill Inhalation Device.
- 2 Remove mask from ISO mask connector.
- 3 Disconnect ISO mask connector from the breathing hose.
- 4 Disconnect breathing hose from ISO connector.
- 5 Remove sampling hose.
- 6 Unscrew ISO connector from fresh gas outlet of the anesthesia machine.



Dismantling Parts of the Ventilator

- 7 Swing out the ventilator door.
- 8 Unlock the clasps to remove the cover.
- 9 Remove the diaphragm.



Disinfecting/Cleaning

Disinfect and clean the unit daily. If using microbial filters:

- Daily disinfecting/cleaning is required for the following components: breathing hoses, breathing bag, Y-piece.
- Weekly disinfecting/cleaning is required for the following components: compact breathing system and its components, and ventilator components.

Wipe off impurities with a disposable cloth.

Wipe Disinfecting

The following components should be wipe-disinfected:

- Fabius basic unit and trolley
- Dräger-Vapor[®]
- Pressure hoses
- Fresh gas hose

Disinfecting and Cleaning

Use surface disinfectants for disinfecting. For reasons of material compatibility, preparations based on the following active agents are suitable:

- aldehydes
- alcohols
- quaternary ammonium compounds.

Due to possible material damage, preparations based on the following active agents are not suitable:

- compounds containing alkylamine phenols
- compounds that release halogen
- strong organic acids
- compounds that release oxygen.

For users in the Federal Republic of Germany we recommend that only disinfectants on the current DGHM list are used (DGHM: German Society for Hygiene and Microbiology). The DGHM list (published by mph-Verlag GmbH, Wiesbaden) also classifies each disinfectant by its active agents. Note the manufacturers' directions for use.

For countries in which the DGHM list is not available, we

recommend the types of disinfectant given above.

All removable parts can also be disinfected with steam (93 °C/10 minutes).

- Do not disinfect microbial filters.
- Do not disinfect or sterilize the O2 sensor.

Automatic Cleaning/Disinfecting Machine

Only with detergent (93 °C/10 minutes).

The following components may be disinfected with moist heat:

- ventilation hose
- all parts of the compact, semi-open compact breathing systems, except for the O2 sensor
- ventilator cover
- diaphragm
- flow sensor with measuring hoses
- fresh gas hoses
- Magill System Components (Pressure gauge must be removed)

Do not add disinfectant, because these parts may corrode.

Sterilizing

The following components can be sterilized at 134 °C in hot steam:

- ventilation hose
- all parts of the compact breathing systems except for the O2 sensor
- ventilator cover
- diaphragm
- flow sensor with measuring hoses.

Checking Readiness for Operation

Refer to "Preparation" on page 22 and "Daily Pre-Use Check for Fabius" on page 66 for reassemble and pre-use checkout prior to operating the Fabius.

Maintenance Intervals

Clean and disinfect the machine and components before each service (and also when returning for repair).

When Required:

- Replace the O2 sensor as soon as calibration is no longer possible.
- Replace the pressure-measuring line (silicone rubber hose and sleeve), flow sensor or flow sensor plug-in probe if damaged.

Every 6 Months:

Inspection and service by trained service personnel. Dräger recommends DrägerService.

- Fabius
- Breathing systems
- Dräger-Vapor[®]
- Sensors

Dräger-Vapor® with report, in connection with Fabius.

Annually:

- Replace the bacterial filter on the pressure-measuring line.
- Replace the measuring line for pressure and flow.
- Replace the diaphragm in the ventilator.

Inspect and verify proper working order:

- Magill Circuit

After 3 Years:

By trained service personnel:

- Replace the lead gel rechargeable battery for the back-up power supply.
- Replace the O-rings and diaphragm of the ventilator.

After 6 Years:

Fabius anesthesia machine:

 General overhaul of the pressure regulator (pressure reducer) by trained service personnel.

Disposal

This device is subject to EU Directive 2002/96/EC (WEEE).It is not registered for use in private households, and may not be disposed of at municipal collection points for waste electricaland electronic equipment.

Dräger Medical has authorized a firm to dispose of this device in the proper manner: for more detailed information, please contact your local Dräger Medical organization.

Batteries and O2 sensors


- Batteries must be disposed of in conformity with the local waste disposal regulations.
- Spent O2 sensors can be returned to Dräger Medical AG & Co. KG Moislinger Allee 53–55 23542 Lübeck
- Do not open forcibly: danger of chemical burns.
- Do not incinerate: danger of explosion.

Bacterial filters

 Must be disposed of as infectious special waste. Can be incinerated at temperatures above 800 °C with minimal environmental pollution.


What's what

Compact Breathing System, Front



- 1 O2 sensor on inspiratory valve
- 2 Inspiratory valve
- 3 Switching lever MAN and SPONT/IPPV on pressure limiting valve
- 4 Connector for breathing bag
- 5 Fresh gas outlet
- 6 Dräger-Vapor® anesthetic agent vaporizer
- 7 Oxygen flush
- 8 Flowmeter tube block (2-gas or 3-gas)
- 9 Anesthesia ventilator
- 10 Control box (settings for ventilation parameters and airway monitoring)
- 11 Ventilator hose

Compact Breathing System, Back



- 1 Power cable
- 2 On/off switch
- 3 Internal interface
- 4 Connectors for medical gas pipeline supply
- 5 Connectors for gas cylinder supply
- 6 Oxygen outlet port (for secretion aspirator)
- 7 Flow sensor
- 8 Ventilator hose
- 9 Connector for pressure-measuring hose
- 10 Pressure-measuring hose
- 11 Control line PEEP/Pmax
- 12 COM Medibus interface/RS 232
- 13 Connection for PEEP/Pmax control line
- 14 Expiratory valve

Technical Data

Ambient Conditions

During operation

15 to 35 °C

Temperature

(machine should be at ambient temperature)

Atmospheric pressure

700 to 1060 hPa

Relative humidity

20 to 80 %

During storage

Temperature

-10 to 60 °C

Atmospheric pressure

500 to 1060 hPa

Relative humidity

0 to 80 %

Machine data

Gas inlet from medical gas pipeline system

Pressure at machine connector

O2, N2O, Air

270 to 550 kPa

Gas supply connectors

NIST or DISS (where required)

Gas quality:

Dew point

>5 °C below ambient temperature

Oil content

<0.1 mg/m³

Particle content

Dust-free air (filtered with pore size <1 µm)

Gas supply from supplementary O2 and N2O cylinders (with threaded connections)

Pressure at machine connector

O2, N2O

500 kPa

Each inlet is fitted with a non-return valve

Gas supply from supplementary O2 and N2O cylinders (with pin-index connections)

Cylinder Connections

Pin-indexed hanger yokes (CGA V-1-1994)

Cylinder Gas Pressure

O₂, AlR 1900 PSI (13100 kPa)

(typical full loads at 21 °C)

N2O 745 PSI (5130 kPa)

Fresh-gas outlet

Male cone 22 ISO, Female cone 15 ISO,

(with thread to secure)

Equipment Class

↑ Class 1, Type B, IXPO

Dimensions and weight

Weight (without vaporizers)

app. 90 kg (depending on trolley version)

Dimensions

app. 960 x 1320 x 670 mm (WxHxD)(depending on trolley version)

Power supply

230 to 240 V (±10 %), 50/60 Hz, 0.7 A or

Rating Non-configurable

100 to 127 V (+10 %, -15 %), 50/60 Hz, 1.4 A

Rechargeable battery

24 V; 3.5 Ah

Operation time with fully charged battery

1.5 hours

Recharging time at least 4 hours on the mains (power switch ON) for full operation time.

Never allow the battery to completely discharge. If the battery does discharge completely, recharge immediately.


In case of power failure, the machine will continue without interruption.

A "BATTERY POWERED" caution message is signalled every five minutes.

Fuses


1 Mains fuse at the rear:

For 230 to 240 V supply voltage: 2 x T1AL 250 V IEC 127/III For 100 to 127 V supply voltage: 2 x T2AL 250 V IEC 127/III




Fuses located on circuit board:

- 1 x T4AL 250 V IEC 127/III
- 1 x T3.15AL 250 V IEC 127/III
- 1 x T2.5AL 250 V IEC 127/III
- 1 x T1.6AL 250 V IEC 127/III



Internal battery fuse:

1 x T3.15AL 250 V IEC 127/III



Electromagnetic Compatibility (EMC)

Conforming to EN 60601-1-2 (1993)

The operation of this anesthetic workstation or module may be adversely affected by electromagnetic interference exceeding the levels specified in EN 60601-1-2 (1993).

f

Technical Data of the Ventilator

(electronically controlled, electrically driven, without drive gas consumption)

Pmax Pressure limiting 10 to 70 hPa* ±5 hPa

VT Tidal volume 50 to 1400 mL \pm max {5 %; 20 mL}

No compliance correction

Breathing frequency 6 to 60 breath/min (±1 breath/min)

TI: TE Inspiration/expiration ratio $1:3 \text{ to } 2:1 \pm 100 \text{ ms}$

TIP/TI Inspiration pause 5 % to 50 % ±100 ms

PEEP End-expiratory pressure 0.2 to 15 hPa ±2 hPa or ±20 %

Pressure-limiting valve MAN mode 5 to 70 hPa ±15 %

SPONT/IPPV mode 1.5 ±1 hPa ±15 %

Performance Characteristics

Fresh gas metering Low-flow flowmeters (calibrated for 20 °C, 1013 hPa)

O2: 0.02 to 0.5 L/min ±10 %
O2: 0.55 to 10.0 L/min ±10 %
N2O: 0.02 to 0.5 L/min ±10 %

N2O: 0.55 to 10.0 L/min ±10 %

Air: 0.2 to 12 L/min ±10 %

O2 flush (bypass) Volume supplied

at 500 kPa: max. 70 L/min at 270 kPa: min. 28 L/min

Pressure limiting 80 kPa ±5 kPa at the gas outlet

Oz outlet: for secretion aspirator

300 to 500 kPa, self-closing

max. 20 L/min

Anesthetic Agent Flow Control

Dräger-Vapor® quick-change plug-in system for two anesthetic agent vaporizers.

The connections are automatically closed and sealed against the atmosphere as soon as the Vapor is removed.

Halothane vapor

Enflurane vapor

Isoflurane vapor

Sevoflurane vapor

Devapor, D-Vapor or D-Tec for

Desflurane

Technical data of anesthetic agent vaporizers - see specific Instructions for Use.

Note conditions for use of supplementary equipment. This may limit the range of use of the overall system.

^{*} hPa = mbar

Monitori	ng and Measurement	Range	in steps of	Accuracy	Condition
Paw	Airway pressure	–9 to 99 hPa	1 hPa	±4 % of the measured value or ±2 hPa, whichever is greater	
Ve	Expiratory minute volume	0 to 99.9 L/min	0.1 L/min	±20 %	ATPS
	Expiratory tidal volume	40 to 9995 mL	5 mL	±20 %	ATPS
f	Breathing frequency	0 to 60 breath/min	±1 breath/min	±1 breath/min	
FiO ₂	O2 measurement in the main gas flow	5 to 100 vol.%	1 vol.%	±3 vol.%	with reference to ambient pressure during calibration
	Reaction time	Less than 15 s			
	Service life of sensor cell	300.00 vol.% h			

Use adequate monitoring for CO₂ and anesthetic agent. Sampling rate max. 250 mL/min.

Connect the sampling line to the Y-piece with an adapter or insert a T-piece in the patient connection of the Y-piece.

Breathing System

	Semi-Ope Compact	en Breathing	System	Compact Breathing System			
	Complian	0.2 L + bag ce: 0.22 ml tic mode (II	•	Volume: 2.8 L + bag Compliance: 0.22 mL/hPa in automatic mode (IPPV) Absorber volume: 1500 mL			
		Resistance reathing Sy		Resistance of Breathing System			
	5 L/min	30 L/min	60 L/min	5 L/min	30 L/min	60 L/min	
Inspiratory Resistance	0.4 hPa	1.1 hPa	3.0 hPa	0.5 hPa	1.3 hPa	2.8 hPa	
Expiratory Resistance	1.3 hPa	3.9 hPa	7.1 hPa	0.7 hPa	2.4 hPa	5.3 hPa	

Classification IIb

In accordance Directive 93/42/EEC Annex IX

UMDNS Code 10-134

Universal Medical Device Nomenclature System

Interface

Serial interface

COM

Protocol

MEDIBUS

Connector

9-pin sub-D (socket/female), galvanic isolation 1.5 kV

Pin assignment

1	Not used
2	TX, transmit
3	RX, receive
4	Not used
5	GND, ground
6	Not used
7	Not used
8	Not used
9	Not used
Settings	9600 Baud

Even parity 8 data bits 1 start bit

1 stop bit

Device identification

ID Number Name

8085 Fabius

Refer to Instructions for Use of Medibus for Fabius CE 90 37 975

Safety Features

The Fabius inhalation anesthesia machine is equipped with the following safety features:

- Positive ventilation pressure limiting 75 ±5 hPa
- Negative pressure safety valve (additional air valve)
 -4 ±2 hPa

O2 Shortage Signal/N2O Cut-off

Alarm limit Warning signal as soon as the

pressure drops below 1.8 bar

Alarm signal High priority alarm, duration at least

7 seconds.

Cannot be deactivated.

During O2 failure N2O is cut off.

S-ORC (Sensitive Oxygen Ratio Controller; control device to ensure a minimum oxygen concentration)

S-ORC is a control element which guarantees a minimum O2 concentration in the fresh gas flow.

As from a flow rate of approx. 300 mL/min, the N2O concentration in the fresh gas can be freely set between 0 and 75 %.

During O2 shortage S-ORC limits the N2O concentration in the fresh gas, so that the O2 concentration does not drop below 21 vol.%.

N2O metering valve open and O2 metering valve closed or O2 flow 0.2 L/min: S-ORC prevents N2O flow.

During N2O failure O2 may still be administered. No alarm. The float in the N2O measuring tube drops to zero.


Reactivate the O2 Shortage Alarm: after each O2 supply failure. After the O2 supply has been restored, a supply pressure of at least 2.7 bar must be applied for at least 20 seconds before another O2 shortage signal can be emitted. During this period, do not activate any devices that consume O2 (e.g. O2 flush, O2 fresh gas flow or secretion aspiration).

However, S-ORC is not an oxygen-specific monitoring device and therefore cannot provide any protection against the effects of accidental use of the wrong gas.


Therefore: always monitor the O2 concentration.

Diagrams


Compact Breathing System



Semi-Open Compact Breathing System



Gas Delivery Unit (3-Gas version)



Abbreviations Used

Abbreviation	Description	Unit
fIPPV	Ventilation frequency	1/min
FiO ₂	Inspiratory oxygen concentration	vol.%
FGE	Fresh Gas Decoupling	
MV	Expiratory minute volume	L/min
Paw	Airway pressure	hPa
Peak	Peak pressure	hPa
PEEP	Positive ventilation pressure in the end-expiratory phase	hPa
Pmean	Mean pressure	hPa
Vī	Tidal volume	mL

Fabius

Carial Mumber

Daily Pre-Use Check for Fabius

Please note that this Daily Pre-use check list takes into consideration all possible configurations of the Fabius. The clinician need only use those areas that apply to their specific Fabius configuration.

All checks must be carried out daily before equipment is used. The person who carries out the checks must be fully conversant with the Instruction for Use. Checks marked with a P must be carried out before each patient use.

These pages should be removed and copied to establish a daily record of machines checks. Mark each function when checks have been satisfactorily completed.

26	i ai i	
Dra	-co	nditions
		Inspection intervals for machine and accessories are current.
P		Machine fully assembled and connected
		Monitors (O2, P, V, CO2 anesthetic agent) switched on, self test carried out satisfactorily.
P		Sampling line for gas monitoring connected to Luer lock on the Y-piece, correct anesthetic agent selected
P		Devapor (when being used) switched on
Ch	ecki	ing Reserve Power
Þ		Verify that battery is fully charged. (If the battery does not show full a charge, the battery operation time is not guaranteed to be 1.5 hours.)
Ch	ecki	ing the Medical Gas Connections
		Visually inspect all supply hoses from the medical gas pipeline system to make sure that they fit tightly and perfectly and have no leaks.
		Verify that all medical gas pipeline supplies are within acceptable pressure ranges.
		Open reserve gas cylinders (when present).
		O2 pressure more than 70 bar.
		N2O pressure greater than 43 bar.
		Close reserve gas cylinders.
O2	Flu	sh Function
		Press O ₂ flush: A strong flow of gas should be emitted from the patient connection.

A-Cone Switch Function (If present)

Toggle the A-cone switch to alternate positions
 Press O2 flush. A strong flow of gas should emit from the alternate patient connection.

Checking the Flow Control/Metering System

- ☐ Fully open the O₂ metering valve. O₂ flow of at least 10 L/min present.
- ☐ Fully open the N2O metering valve. N2O flow of at least 10 L/min present.
- ☐ Cut off the O2 supply. Remove the O2 connector or close the O2 cylinder valve. The O2 shortage alarm signal must be activated N2O must not flow.
- Restore the O2 supply. N2O must flow again.

Checking the Flow Control/Metering System

- ☐ Set O2 metering valve to 1.5 L/min. N2O flow = 3 L/min to 5 L/min.
- ☐ Close the O₂ metering valve: No N₂O flow.
- Open the AIR flow control valve. Air flow of at least 10 L/min present.
- ☐ Close all metering valves.

Sensor Calibration

Before calibrating the O2 and flow sensors, wait at least 5 minutes after switching on for the control unit to warm up to its normal operating temperature.

- ☐ Calibrate O2 Sensor.
- ☐ Calibrate flow sensor.

Checking the Gas Type

- ☐ Unscrew the fresh gas hose from the Fabius.
- ☐ Set the O₂ metering valve to approx. 3 L/min.
- □ Remove the O₂ sensor from the inspiratory valve; hold it under the fresh gas outlet and allow fresh gas to flow over it. The control unit must indicate an O₂ concentration of approx. 100 vol.%.
- ☐ Set the N2O metering valve to approx. 3 L/min. The displayed O2 concentration must now be approx. 50 vol.% O2.
- ☐ Reclose the metering valves for O₂ and N₂O.
- ☐ Replace the O₂ sensor on the inspiratory valve.
- ☐ Reconnect the fresh gas hose to the Fabius.

^{*} The O2 shortage alarm is ready to operate again after an operating pressure of at least 2.7 bar has been present for 20 seconds, provided that no O2 has been used during this time.

Va P P P P	por	19.n, Vapor 2000 (Tec 5) Fastening; Latched down firmly and set vertically. Handwheel; In zero position and engaged. Filling level between min. and max. Interlock; Locking function OK (when present). Key-indexed filling system; Sealing key or pin inserted and closed tight (when present). Quik Fil or Funnel filling system; Locking screw tight (when present).	Leak Testing the Compact (Semi-Open) Breathing Systems Test once without Dräger-Vapor® – after removing the vaporizer and once with each Dräger-Vapor® with the handwheel set to zero. In the following check list, please note that there are three boxes for each item to be checked. Box A = Fabius, No Vaporizers Box B = Fabius, Single Vaporizer Box C = Fabius with Dual Vaporizers					
De	vap	or, D-Vapor (if fitted)		A	В	С		
P -		Fastening; Latched down firmly and set vertically.	P				Press the key on the control unit.	
P -		Handwheel; In zero position and engaged.	P				Remove the breathing bag and connector.	
P		Filling level between min. and max.	Р				Fit the free end of the breathing hose to the	
Р	Ц	Mains cable plugged in, LED »Operational« on.		_	_		Y-piece.	
Κe	v.inr	dexed filling system	P				Set the pressure-limiting valve to MAN and set to 70 hPa.	
P		Filling opening locked shut and pointing downward.	P				With O ₂ flush, build up a pressure of 25 to 30 hPa, while monitoring the airway pressure.	
Se	lect	atec [®]	P				•	
P		Fastening; Latched down firmly and set vertically.	•	_		_	The airway pressure must rise or remain	
Р		Handwheel; In zero position and engaged.	_	_	_	_	constant.	
Ρ		Filling level between min. and max.	Р	Ц		Ц	Switch pressure-limiting valve to SPONT. Pressure drops.	
Р	Ц	Interlock; Locking function OK (when present).	Р				Switch pressure-limiting valve to MAN.	
Ch	ecki	ing the Condition of CO2 Absorbent	Р				Set Vapor to lowest concentration shown on	
(C	omp	act Breathing Systems)	Р			П	handwheel.	
P		Color change is no more than half the column of CO2 absorbent.	r	Ц	ш		Press O ₂ flush; let pressure build to 30 mbar; pressure remains constant or rises	
			P				Reset Vapor to zero.	
Ch	ecki	ing Non-Rebreathing of Semi-Open	P				Release pressure.	
		act Breathing System)	Р				Push lever on pressure-limiting valve.	
		Disconnect the fresh gas hose.	P				Reopen the patient connection and	
		Put the test adapter on the fresh gas hose and connect to the patient Y-piece.					reconnect the manual ventilation bag.	
		Occlude the semi-open exhaust port.					not perfectly leak-proof	
		Set the APL to 70 mbar.	(i.e				reduced):	
		Set the O ₂ metering valve to 0.25 L/min.	 Check that all plug-in, push-fit and screw connectors tightly. Replace any missing or damaged seals. 					
		With the O2 flush, build up a pressure of 60 mbar.						
		Airway pressure must rise or remain constant.	•			_	, call DrägerService.	
ре	rforn	ystem fails, locate and correct any external leaks by ning the leak test of the compact breathing system peat this test. If necessary, call DrägerService.			= - -	· ,	,	

		e stin orese		agiil Circuit						
Tes	st on	ce w	itho	ut Dräger-Vapor [®] – after removing the	Ch	ecki	ng Ventilator Operation			
				nce with each Dräger-Vapor® with the	P		Press the key on the control unit.			
	handwheel set to zero. In the following check list, please note that there are three						Connect another breathing bag to the Y-piece to act as test lung.			
box	kes f	or ea	ach i	tem to be checked.	P		Check that ventilation pressure is available.			
				ius, No Vaporizers	Р		Check that the ventilator piston is working correctly.			
	Вох	B=	Fab	ius, Single Vaporizer	Р		Monitor the functions of the inspiratory and			
	Вох	C =	Fab	ius with Dual Vaporizers			expiratory valves only on Compact and semi-open Compact breathing systems.			
	Α	В	С		Ρ		Check that the breathing bag (test lung) on the			
P				Close off relief valve and seal patient connection.			Y-piece is ventilated.			
Р				Gently pinch off the sampling hose	Air	way	Monitor			
Р				(if present). Open O2 metering valve on anesthesia unit,			rm function can be tested by setting alarm limits to hat are certain to trigger an alarm.			
				fill breathing bag, and then close the O2						
_	_	_	-	metering valve.			Test the flow monitor and alarm module.			
Р	Ш	Ц		The breathing bag must not deflate. If the bag deflates, determine where the leak is coming from and take corrective action.			Test the airway monitor and alarm module.			
Ρ		☐ ☐ Open relief valve and sampling hose.				Additional Measuring Units (when present)				
							Check the CO2 monitor and alarm module.			
							Check the anesthetic agent monitor and alarm			
				Expiratory Valves thing System)			module.			
Se	t pre	ssur	e-lin	niting valve to 30 mbar in manual mode.	Anesthetic Gas Scavenging Line					
Pre	ess (D2 flu	ush;	let pressure build up to 30 mbar.	P		Float in AGS flowmeter between minimum and			
P		Bre	athir	ng bag for manual ventilation fills.			maximum.			
Р		O2	on n	nonitor greater than 80 vol.%.	P		Check the hose connections.			
Р				sc in inspiratory valve moves when breathing queezed during manual ventilation.	Р		Press and hold O ₂ Flush: airway pressure <10 cm/H ₂ O with Y piece occluded.			
P				sc in expiratory valve moves when breathing eleased.						
P				g on patient connection fills and empties cally during manual ventilation.						
P				stance during manual ventilation which might stenosis.						
Pre	essu	re-Li	imiti	ng Valve (Compact Breathing System)						
P		Set	APL	valve to MAN and 30 mbar.						
				flow-control valve to 5 L/min. on monitor; 26 to 35 mbar.						

		l Ve r	ntilation Bag for Emergency Ventilation ent)						
			Check that the bag is functioning correctly by pumping manually.						
		tanç	When the bag is squeezed, air must audibly and tangibly flow out of the mask cone; when the bag is released, it must rapidly recover its original shape.						
		you	ck off the mask connector (cone) with the ball of r thumb: you should only be able to squeeze the a little.						
P		Bef	ore Connecting to Patient						
		•	Pre-set tidal volume and frequency to suit patient's needs.						
	 Set alarm limits for O2, P, V, CO2, and anesthetic agent to suit patient's needs. 								
		 Adjust the pressure-limiting valve to meet the patient's requirements (Compact) or adjust the pressure-relief valve if the Magill breathing circuit is being used. 							
		•	Set monitor to anesthetic agent being used.						
		•	When there is electrical contact to the patient, establish potential equalization between machines.						
		•	For monitors with sidestream gas measurement, empty water trap.						
ma	chìr	ne m	k can not be carried out satisfactorily, the ust not be used.						
		e ch	neck carried out						
Nai									
Dat	e								

Index

2-gas version	Disinfecting/Cleaning the Fabius 51
3-gas version	Dismantling Parts of the Ventilator 50
	Dismantling the Absorbent Canister
Abbreviations Used	Dismantling the Expiratory Valve
A-Cone Lever Selection	Dismantling the Flow Sensor
	Dismantling the Inspiratory Valve
A-Cone Switch	Dismantling the Pressure-limiting Valve
Additional Screens Button	Displays
Adjusting the Fresh Gas Composition	Disposal of Bacterial Filter 53
Airway Monitor with Control Unit	Disposal of Used Batteries and O2 Sensors 53
Airway pressure (Paw) 5	Dräger Vapor Anesthetic Agent Vaporizer
Alarm limit button	Dräger Vapor®
Alarm Limits	Diagot raporo 1,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,
Alarm silence button	
Ambient Conditions	Electrical Supply
Anesthesia Ventilator	Electromagnetic Compatibility (EMC)
Anesthetic Agent Flow Control	EN 740 (Anaesthesia Workstations and Their Modules - Parti-
Anesthetic Agent Vaporizer	cular Requirements) 5
APL Valve 15, 36	Enflurane 12
Automatic Cleaning/Disinfecting Machine 52	Equipment Class 56
Automatic Ventilation 40	Equipotential Bonding24
	Expiratory port
Battery 56	Expiratory Valve
Belleville Washers	
	False Alarms
Buttons	
	Fault - Cause - Remedy
Calibrating The Flow Sensor	Fitting Microbial Filters
Calibrating the Oxygen Sensor	Fitting the Pressure-limiting Valve
Caution During Transport	Fitting the Waste Gas Outlet Port
Caution When Handling O2 Cylinders	Fresh Gas Composition
Ceiling mounted Fabius 9	Fresh Gas Decoupling 5
Central Pipeline Gas Supply	Fresh Gas Flow
Checking Readiness for Operation	Fresh gas hose
Classification	Fresh-gas outlet 56
COM 55, 60	Function Buttons
Compact Breathing System	Fuses 57
Compact Breathing System, Back	
Compact Breathing System, Front	Gas Compensation Setting
Configuration Menu	Gas Compensation Setting Menu
Connect a Gas Cylinder to Its Yoke	Gas Delivery Unit
Connecting the Breathing Hoses	Gas inlet from medical gas pipeline system
Connecting the Compact Breathing System	Gas supply
Connecting the Flow Sensor	Gas supply from supplementary O2 and N2O cylinders 56
	The state of the s
Connecting the O2 Sensor Capsule	General Description 6
Connecting the Pressure Sensor	
Control hose	Halothane 12
Currently active alarms screen	Hand-tighten 27
Curve IPPV	
Cylinders with Pin-index Mounting	Inserting the Flow Sensor
Cylinders with Threaded Connectors	Inserting the Flow Sensor
	Inspiratory Volve
Desflurane	Inspiratory Valve
Diagrams	Installation of the Semi-Open Compact Breathing Adapter
Dimensions and weight	34
Disinfecting and Cleaning the Magill Components 51	Intended Use
Distinctioning and Oleaning the magnit Compension 11111. OT	Interface 55, 60

IPPV 5,		Replacing Used Absorbent (Compact System)	
SO Mask		Rotary Control	14
soflurane	12	RS 232	55
Key Code	21	Safety	. 4
Knurled Screw	27	Safety Features	. 5
·		Scavenging System for Fabius	
Longuago	00	Screen Display Configuration	
Language		Sealing Rings	
Language Selection Menu		Selectatec	
Liability for proper function or damage		Selector Lever	
Lower alarm limit for O2	39	Semi-open Compact Breathing System	
		Sensor Calibration Button	
Machine data	56	Setting the Vapor	
Magill Breathing Circuit	11	Sevoflurane	
Magill Breathing System	. 9	Shut-down	
Magill-Breathing Device		S-ORC (Sensitive Oxygen Ratio Controller)	
Maintenance Intervals		SPONT	
MAN		Spontaneous Breathing	
Manual Ventilation		Standard Mode	
Manual Ventilation Bag		Standard Screen Button	
Maximum pressure (Pmax)	18		
	60	Standby button	
Medical Gas Pipeline Supply		Sterilizing	
Menu Button		Stripping Down Fabius Compact Breathing System	
		Stripping Down the Fabius	
Messages		Switching the Flow	11
Metering valve			
Minimum Flow of Anesthesia		Technical Data	56
Mode		Technical data of the interface	
Mode Selection Buttons		Technical Data of the Ventilator	
Monitoring and Measurement	59	Transfer hose	
Never Use Tools		UMDNS Code	59
Nitrogen	37	Upper alarm limit for Paw	
		User Interface	
02 Flush	43	Using the Magill Breathing System	
O2 Monitoring		Osing the Magin Breathing System	72
O2 Sensor Capsule			
O2 Shortage Signal/N2O Cut-off		Vapor Interlock System	12
Operation		Vaporizers	. 6
Oxygen concentration (FiO2)		Ventilation and Monitor Operating Concept	14
ongon conconsument (1102)		Ventilation hose	27
		Ventilation Modes	38
PEEP 10,		Ventilation Parameter Buttons	16
Performance Characteristics	58	Ventilator Safety Feature	25
Pin Index Trolley		Ventilator settings screen	
Pin-index	23	Volume (V)	
PM 8050	. 7	(,, ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	
Pmax	40		
Power On Setting	21	Warning Tones	
Power Source		Washer	
		Wipe Disinfecting	51
Preparation			
Preparing the Ventilator		Yoke Handle	9.2
Pressure limited ventilation		TORO Haridio Triffic T	20
	71		
Decation to the CO Charles At	0.4		
Reactivate the O2 Shortage Alarm			
SEDIACION LAUZ ADSORBÊNT	44		



If no Serial No, has been filled in by Dräger these Instructions for Use are provided for general information only and are not intended for use with any specific machine or device.



Directive 93/42/EEG Concerning Medical Devices

Dräger Medical AG & Co. KG

Germany

☆ Moislinger Allee 53 – 55D-23542 Lübeck

Ø +49 451 8 82- 0

FAX +49 451 8 82- 20 80 http://www.draeger.com

90 37 441 - GA 5330.400 en
© Dräger Medical AG & Co. KG
5th edition - January 2006
Subject to alteration