CONSTRUCT









A guide for using the Construct Pretreated Reinforcing Fiber System

Thank you for your decision to purchase another quality Kerr product. For over 100 years, Kerr has built a solid reputation for leading-edge products that represent a great value to the dental community. Please take a moment to read these instructions before using the material.










CONSTRUCT

Reinforced Crowns











CONSTRUCT

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
CONSTRUCT

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
CONSTRUCT

Lingual Splints - Direct Technique

Lingual Splints - Indirect Technique

CONSTRUCT

CONSTRUCT™ is an ultra-high strength reinforcing polyethylene braid that can be used to provide internal matrices for all Polymer-Ceramic metal-less crowns bridges and splints. The surface of the CONSTRUCT™ fibers are cold gas plasma treated, silanated, and impregnated with resin to maximize its adhesive qualities.

CONSTRUCT™ can be handled by fingers prior to applying the tinted resins without fear of contamination. To apply the resin, slightly spread the braided weave by pushing the ends inward and impregnate it with the tinted resin.

CONSTRUCT™ is packaged in a convenient spool. Simply measure and cut the amount needed. There is no wastage of material. A special pair of scissors is provided in the kit to cut the CONSTRUCT™.

NOTE: The Light Resin is to be used in the lighter shades such as A-1 to A-3. The Dark Resin is to be used with the darker shades such as A-3 to A-4 and will work on all shades from the A range to the D range. You can mix the shades to custom tint the CONSTRUCT™ fibers. The Neutral Resin is used on light and extra light shades.

CONSTRUCT™ is available in 1mm, 2mm and 3mm widths. Our unique braid utilized in the CONSTRUCT™ can be expanded or contracted to the exact width of the preparation.

CONSTRUCT™ INTRODUCTORY KIT (PN# 30869)

1mm, 2mm and 3 mm CONSTRUCT™ Braid
Neutral, Light and Dark CONSTRUCT™ Resin
Scissors
Tweezers

CONSTRUCT™ REFILL PACKAGES:

PN#	30870	1mm CONSTRUCT™ Braid
PN#	30871	2mm CONSTRUCT™ Braid
PN#	30872	3mm CONSTRUCT™ Braid
PN#	30873	Neutral CONSTRUCT™ Resin
PN#	30874	Light CONSTRUCT™ Resin
PN#	30875	Dark CONSTRUCT™ Resin

APPLICATIONS OF CONSTRUCT™ CONSTRUCT™ REINFORCED CROWNS

1. Lightly mark with red (or blue) margin liner, and block out all undercuts. Apply a separator (Rubbersep or belleGlass HP Separator A/B) following the recommended instructions. Let separator completely dry.
2. Dispense 1 to 2mm of the desired shade of Opaqueous Dentin onto the Mixing Palette.
3. Place the appropriate belleGlass HP Opaqueous Dentin shade over the die. Start at the top and use a brushing stroke to shape the Dentin material over the entire prep to the margins. This should result in a “tent” formation. The coping should be approximately 0.2mm thick when properly shaped. Do not cure at this point.
4. Cut CONSTRUCT™ to desired length. Using tweezers, lift the appropriate length of CONSTRUCT™

and cover the entire circumference of the crown allowing for a slight overlap at the joint.

Hint: Any dead soft material, i.e. dental floss, can be used to measure the length of CONSTRUCT™ required for the procedure.

Expand the width of the CONSTRUCT™ and impregnate it with the appropriate shade of CONSTRUCT™ Resin.

- Using tweezers, place the fiber on the coping and wrap the fibers around the die approximately 1/3 up from the margin.
- Work the fibers in with a spatula to remove all air. Continue adding a thin layer of Opaque Dentin to mask out the CONSTRUCT™. Cure all surfaces for 40 seconds by exposure to a light-curing unit. NOTE: If you feel it necessary to darken the fibers on the coping, paint on a thin layer of the dark resin and light cure for 40 seconds. After curing, lightly grasp the coping and loosen it from the die. Immediately replace it onto the die again. This will help in the removal of the crown after processing. Build up the appropriate Opaque/Translucent Dentin shade shy of the occlusion and check for proper occlusal clearances. Light cure for 40 seconds.
- Place belleGlass HP Enamel composite as required. Light cure 40 seconds on all surfaces. Cure in the belleGlass™ HP Unit 10 minutes.
- Finish and polish using provided green stones, silicone rubber wheels, carbide burs, and Composite Polishing Kit.
- Use MICROETCHER to etch inner surface of crown(s) to enhance adhesive qualities for cementation.
- Rinse/Dry. Silanate the bonding surface with Silane Primer. Thoroughly dry.

CONSTRUCT™ REINFORCED LINGUAL SPLINTS (INDIRECT TECHNIQUE)

- Lightly mark with red (or blue) margin liner, and block out all undercuts. Apply a separator (Rubbersep or belleGlass HP Separator A/B) following the recommended instructions. Let separator completely dry or set.
- Select correct shade/opacity of belleGlass HP Opaque Dentin . Lay down a base of Opaque Dentin 0.2mm thick in the outline on stone cast.
- Cut the length of CONSTRUCT™ to be used. Expand the width of the CONSTRUCT™ and impregnate it with the appropriate shade of CONSTRUCT™ Resin.
- Place the CONSTRUCT™ over the Opaque Dentin base. Work it in with a spatula to remove all air and light cure for 40 seconds.
- Add a second layer of belleGlass™ 0.3mm thick to laminate fibers. Light cure for a minimum of 40 seconds per surface. Cure in the belleGlass™ HP Unit 20 minutes. Grind to contour and rubber wheel for smoothness. Remove splint from cast and finish to outline.
- Polish and re-fit to check margins. Sandblast the tooth contact side with 50 micron aluminum oxide using the Microetcher. Clean and disinfect the splint with alcohol or acetone. Silanate the bonding surface. Protect case during shipment.

CONSTRUCT™ REINFORCED ANTERIOR BRIDGE

- Lightly mark with red (or blue) margin liner, and block out all undercuts. Apply a separator (Rubbersep or belleGlass HP Separator A/B) following the recommended instructions. Let separator completely dry or set.
- Cut a length of CONSTRUCT™ fibers long enough to wrap around the entire circumference of the abutments and pontic to make the triple beam with one continuous piece.
- Dispense the desired shade of Opaque Dentin onto the Mixing Palette.
- Lay down a base of belleGlass Opaque Dentin 0.2mm thick over the abutments forming coping(s). With a brushing motion shape composite to the margins.
- Bridge the copings together with a bar of Opaque Dentin about 1 mm in width. DO NOT LIGHT CURE! NOTE: THE MAXIMUM ALLOWABLE SPAN IS 15 mm.

6. Expand the width of desired length CONSTRUCT™ and saturate it with the appropriate shade of CONSTRUCT™ tinted Resin. Holding the CONSTRUCT™ with tweezers, layer the fiber from the mesial lingual of the distal abutment. Hold in a vertical position and press into the soft coping.
7. Extend the CONSTRUCT™ across the bridge and wrap it around the mesial abutment pressing it into the soft coping. Return it across the pontic and continue wrapping it around the distal abutment back to the start position of the layering. Pinch the contacts together and adapt the beam to the desired shape.
8. Press the CONSTRUCT™ along the beam making a third layer of CONSTRUCT™ ending it slightly beyond the contact area. Pinch the contacts together and adjust the shape of the substructure to clear the opposing bite. Light cure all surfaces for 40 seconds.
9. Add a thin layer of belleGlass HP Opaque Dentin over the cured CONSTRUCT™ to cover the braid and to add bulk to the pontic. The substructure is now ready for the Translucent Dentin build-up.
Note: When building up the pontic, select a lighter shade of belleGlass HP Opaque Dentin. Due to the thickness of the material used to form the pontic, the shade after curing will appear darker than the abutment if the same shade is used.
10. Build-up Dentin/Body of both abutments and cure for 40 seconds (allow for incisal). Build up pontic(s) as needed.
11. Apply desired belleGlass HP Enamel over tooth. Light cure 40 seconds on all surfaces. Cure in the belleGlass™ HP Unit for 20 minutes.
12. Finish and polish using provided green stones, rubber wheels, separating disks and carbide burs.

CONSTRUCT™ REINFORCED POSTERIOR BRIDGE

NOTE: FOR SHORT SPANS USE ANTERIOR BRIDGE TECHNIQUE.

1. Lightly mark with red or blue margin liner, and block out all undercuts. Apply a separator (Rubbersep or belleGlass HP Separator A/B) following the recommended instructions. Let separator completely dry or set.
2. Cut a length of CONSTRUCT™ fiber long enough to wrap around the entire circumference of each abutment.
3. Lay down a base of belleGlass Opaque Dentin 0.2mm thick over the abutments forming coping(s). With a brushing motion shape composite to the margins. **DO NOT LIGHT CURE!**
4. Expand the width of desired length CONSTRUCT™ and saturate it with the appropriate shade of CONSTRUCT™ tinted Resin.
5. Adapt the saturated CONSTRUCT™ braid around the coping keeping the fibers as close to the cervical margins as possible, as in the CONSTRUCT™ Crown.
6. Repeat this procedure for any other abutments.
7. Adapt the third piece of the saturated CONSTRUCT™ braid starting on the occlusal surface of the mesial abutment. Continue adapting the CONSTRUCT™ down the distal wall. Extend the fibers across the pontic area and carry it up the distal wall to the occlusal surface of the corresponding abutment, forming a strap. The CONSTRUCT™ braid terminates on the occlusal area.
DO NOT LIGHT CURE AT THIS TIME.

NOTE: THE MAXIMUM ALLOWABLE SPAN IS 15mm.

8. Laminate and pre-cure a support beam using the following technique: Measure and cut two lengths of CONSTRUCT™ braid. Lay a 1mm thick layer of belleGlass HP Opaque Dentin onto the first piece. Place the second piece on top of the Opaque Dentin forming a sandwich. Light cure for 40 seconds. Trim and fit the beam into the pontic area with the layers in a vertical orientation.
9. Adapt the strap to the bottom of the beam, and lute in place with belleGlass HP Opaque Dentin. Make sure to adjust the substructure for proper occlusal clearance
LIGHT CURE ALL SURFACES FOR 40 SECONDS.
Place a final layer of CONSTRUCT™ over the laminated beam. If occlusal space permits, this horizontal oriented beam may extend onto the occlusal surface of the abutments. Light cure with for 40

seconds. Cover all CONSTRUCT™ with a thin layer of belleGlass HP Opaque Dentin to mask out the braid, and finally light cure for 40 seconds. Continue the build-up as in the anterior bridge technique.

CONSTRUCTING PERIODONTAL SPLINTS (DIRECT TECHNIQUE)

1. Measure the teeth to be splinted by laying a narrow strip of dead-soft material (i.e. Dental Floss) on the appropriate teeth to determine the length of CONSTRUCT™ needed.
2. Cut and place the CONSTRUCT™ on a clean pad or slab, impregnate it with any flowable composite, such as KERR'S REVOLUTION™ and cover it until ready to use.
Note: Mobile teeth can be stabilized in their desired position by tack-curing any hybrid resin such as KERR'S PRODIGY or XRV to the interproximals on the facial surfaces of these teeth. The composite resin can be removed after the splint is placed.
3. Prepare teeth for bonding. Clean the teeth with a sandblaster or use a bur to roughen the enamel prior to cleaning. Prepare the teeth for bonding, apply bonding agent such as KERR'S OPTIBOND FL. Apply a 0.2 mm hybrid resin base (PRODIGY or XRV). The Unidose delivery system makes application easier. DO NOT LIGHT CURE.
4. Adapt CONSTRUCT™ to the teeth using washed, gloved fingers to compress the CONSTRUCT™ through the composite layer so that it touches the teeth and is compressed. Use an instrument to tuck it interproximally so that it adapts to the contours of the teeth. To avoid displacing the already adapted CONSTRUCT™, hold it down with your gloved finger and continue to advance one tooth at a time as you tuck the CONSTRUCT™ interproximally with an instrument.
Once the CONSTRUCT™ is adapted, remove the excess hybrid resin. Feather the hybrid resin at the gingival and incisal margins of the splint.
Light-cure the splint from both the lingual and the labial. Cure one tooth at a time 40 seconds per surface.
Add a final layer of hybrid resin (PRODIGY or XRV) to sandwich the CONSTRUCT™ between resin layers.
Remove the excess hybrid resin that has flowed into the interproximal spaces and light cure for 60 seconds per surface.
5. Check the occlusion and finish with 12 and 30 fluted burs. Polish with KERR Composite Finishing System.

Supplemental KERR products:

29085 KERR Composite Finishing System

014-531 Red Margin liner

014-532 Blue Margin liner

CAUTION: Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact.

LIMITED WARRANTY - LIMITATION OF KERRLAB'S LIABILITY

KerrLab's technical advice, whether verbal or in writing, is designed to assist dentists/laboratories in using KerrLab's product. Such advice does not expand KerrLab's limited warranty or relieve the dentist/laboratory of testing KerrLab's products to determine their suitability for the intended uses and procedures. The dentist/laboratory assumes all risk and liability for damages arising out of the improper use of KerrLab's product. In the event of a defect in material or workmanship, KerrLab's liability is limited, at KerrLab's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to KerrLab. In no event shall KerrLab be liable for any indirect, incidental, or consequential damages.

EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES, BY KERRLAB, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.

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Kerr Italia, SpA
Kerr U.K. Ltd.
Kerr Australia Pty. Ltd.
Sybron Dental Specialties Japan, Inc.