ROTO SILENTA 630 RS

ⓁⒹ	Bedienungsanleitung.....	10
ⓁⒺ	Operating Instructions.....	34
ⓁⒻ	Mode d'emploi	57
ⓁⒼ	Istruzioni per l'uso	81


Fig. 1


Fig. 2

EG-Konformitätserklärung
EC Declaration of conformity
Déclaration de conformité CE
Dichiarazione di conformità CE

des Herstellers / of the manufacturer / du fabricant / del costruttore
Andreas Hettich GmbH & Co. KG • Föhrenstraße 12 • D-78532 Tuttlingen • Germany

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We hereby declare under our sole responsibility that the designated device and its accessories, which are listed in the technical documentation for this device and whose conformity has been assessed together with the device, conform to the Medical Device Directive – Council Directive 93/42/EEC.

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The conformity evaluation process was performed in accordance with appendices V and VII of Directive 93/42/EEC involving the following notified bodies:

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Further partly applicable European directives:

- Machinery Directive 2006/42/EC
- EMC directive 2004/108/EC
- Low voltage directive 2006/95/EC
- RoHS II Directive 2011/65/EC (without involvement of a notified body)

Applied standards:

According to the list of applied standards, which is part of the product file.

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Last validity of this declaration:

In accordance with the last date of validity of the EC certificate of conformity as issued by the above-mentioned notified body in accordance with appendix V of Directive 93/42/EEC which was valid at the time of the issuance of this declaration: 2017-06-09

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- **USA:** QSR, 21CFR 820 "CFR Title 21 - Food and Drugs: TITLE 21- FOOD AND DRUGS, CHAPTER I - FOOD AND DRUG ADMINISTRATION DEPARTMENT OF HEALTH AND HUMAN SERVICES, SUBCHAPTER H - MEDICAL DEVICES, Part 820 QUALITY SYSTEM REGULATIONS"
- **Kanada:** CMDR, SOR/98-282 "Medical Devices Regulations"

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Standards and regulations which apply to this device

The device is a high-end technical product. It is subject to extensive testing and certification procedures according to the following standards and regulations in their respectively valid version:

Electrical and mechanical safety for design and final testing:

Standard series: IEC 61010 (conform to standards of DIN EN 61010)

- IEC 61010-1 "Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 1: General requirements" (Pollution Degree 2, Installation Category II)
- IEC 61010-2-010 "Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-010: Particular requirements for laboratory equipment for the heating of materials" (applied to heated centrifuges only)
- IEC 61010-2-020 "Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-020: Particular requirements for laboratory centrifuges"
- IEC 61010-2-101 "Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-101: Particular requirements for in vitro diagnostic (IVD) medical equipment"

Electromagnetic Compatibility:

- EN 61326-1 "Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements"

Restriction of Hazardous Substances (RoHS II):

- EN 50581 "Technical documentation for assessing electric and electronic devices with regard to the restriction of hazardous substances"

European directives applied for conformity assessment procedures:

- Medical device directive 93/42/EEC (class 2a medical products)
EC conformity assessment procedure according to annex V
"EC DECLARATION OF CONFORMITY (production quality assurance)" and annex VII
"EC DECLARATION OF CONFORMITY (self-declaration by the manufacturer)".
The annex V certificate is issued by following notified body:

MEDCERT GmbH – Notified Body CE 0482
tel: +49 (0)40 - 22633250
fax: +49 (0)40 - 226332515
e-mail: info@medcert.de
website: www.med-cert.com
mail: Pilatuspool 2, D-20355 Hamburg; Germany

- Directive 2011/65/EU for the restriction of use of certain hazardous substances in electric and electronic devices. Carrying out the EC conformity assessment process is the sole responsibility of the manufacturer, without the involvement of a notified body.

Applied medical device regulations outside Europe:

- **USA:** QSR, 21CFR 820 "CFR Title 21 - Food and Drugs: TITLE 21- FOOD AND DRUGS, CHAPTER I - FOOD AND DRUG ADMINISTRATION DEPARTMENT OF HEALTH AND HUMAN SERVICES, SUBCHAPTER H - MEDICAL DEVICES, Part 820 QUALITY SYSTEM REGULATIONS"
- **Canada:** CMDR, SOR/98-282 "Medical Devices Regulations"

Certified quality management system according to

- ISO 9001 "Quality management systems – Requirements"
- ISO13485 "Medical devices - Quality management systems - Requirements for regulatory purposes"

Environmental management system according to

- ISO 14001 "Environmental management systems - Requirements with guidance for use"

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Restrictions relatives aux substances dangereuses (RoHS II):

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- **USA:** QSR, 21CFR 820 "CFR Title 21 - Food and Drugs : TITLE 21- FOOD AND DRUGS, CHAPTER I - FOOD AND DRUG ADMINISTRATION DEPARTMENT OF HEALTH AND HUMAN SERVICES, SUBCHAPTER H - MEDICAL DEVICES, Part 820 QUALITY SYSTEM REGULATIONS"
- **Canada:** CMDR, SOR/98-282 "Medical Devices Regulations"

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Contents

1	Intended use	36
2	Remaining risks	36
3	Technical specifications	37
4	Notes on safety	39
5	Symbol meanings	41
6	Delivery checklist	41
7	Unpacking the centrifuge and installation	42
8	Initial operation	42
9	Opening and closing the lid	42
9.1	Opening the lid	42
9.2	Closing the lid	42
10	Installation and removal of the rotor	43
11	Loading the rotor	43
12	Closing biosafety systems	44
13	Control and display devices	44
13.1	Status symbols	44
13.2	Control knob	44
13.3	Keys and set-up options	44
14	Entering centrifugation parameters	45
14.1	Running time	45
14.1.1	Time selection	45
14.1.2	Continuous running	45
14.2	Speed (RPM)	45
14.3	Relative centrifugal force (RCF/RZB)	45
14.4	Run-up and run-down parameters	46
14.4.1	Run-up step	46
14.4.2	Run-up time	46
14.4.3	Braking step	46
14.4.4	Run-down time	46
14.4.5	Brake de-energisation speed	46
14.5	Radius/Temperature	46
14.5.1	Radius	46
14.5.2	Temperature (for centrifuges with cooling only)	46
14.6	Automatic buffer	46
15	Programming	47
15.1	Entering / modifying programs	47
15.2	Calling up programs	47
16	Centrifugation	47
16.1	Centrifugation with time selection	47
16.2	Continuous running	47
17	Changing settings during the centrifugation run	48
18	Integral RCF (jRCF)	48
19	Displaying the maximum speed of the rotor	48

20	Displaying the maximum RCF of the rotor.....	48
21	Emergency stop	48
22	Audible signal	48
23	Interrogating the operating hours	48
24	Setting the date and time	49
25	Centrifugation data displayed after switching on.	49
26	Immediate display of the centrifugation data after switching on	49
27	Key-operated switch (for centrifuges with key-operated switch only).	49
28	Program linking (for centrifuges with program linking only).	49
28.1	Linking programs	49
28.2	Changing program links	50
28.3	Centrifugation run with program linking	50
28.4	Deleting program linking	50
29	Cooling (only in centrifuges with cooling)	50
29.1	Standby-cooling	50
29.2	Pre-cooling the rotor.....	50
30	Heating (only on centrifuges with heating/cooling option)	50
31	Relative centrifugal force (RCF)	51
32	Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm ³	51
33	Rotor recognition	51
34	Emergency release	52
35	Maintenance and servicing	52
35.1	Centrifuge (housing, lid and centrifuging chamber).....	52
35.1.1	Surface cleaning and care	52
35.1.2	Surface disinfection.....	52
35.1.3	Removal of radioactive contaminants	53
35.2	Rotors and Attachments.....	53
35.2.1	Cleaning and care	53
35.2.2	Disinfection	53
35.2.3	Removal of radioactive contaminants	54
35.2.4	Trunnions.....	54
35.2.5	Rotors and accessories with limited service lives	54
35.3	Autoclaving	54
35.4	Centrifuge containers	54
36	Faults	55
37	Returning Devices	56
38	Disposal	56
39	Anhang / Appendix	105
39.1	Rotoren und Zubehör / Rotors and accessories.....	105

1 Intended use

This device is a laboratory centrifuge which is suitable for medical and in-vitro diagnostic applications.

The therapeutic purpose with the corresponding centrifuge vessels is the separation of blood by means of centrifugation for the purpose of transfusions or autotransfusion of the won blood components.

The centrifuge is also suitable for preparing human bodily fluids by means of centrifugation for in-vitro diagnostic purposes.

Furthermore, the centrifuge is generally also suitable for separating substances or substance mixtures of non-human origin having a density of max. 1.2 kg/dm³.

The centrifuge may only be used by skilled personnel in closed laboratories.

The centrifuge is only meant for the above-named purposes.

Another use or one which goes beyond this, is considered to be non-intended. The company Andreas Hettich GmbH & Co. KG is not liable for damage resulting from this.

Observing all information in the operating instructions and complying with the measures described therein is also a part of the intended use.

2 Remaining risks

The device is built according to the state-of-the-art and the recognized safety regulations.

If used and handled improperly, there could be life-threatening danger to the user or third parties, or the device could be impaired or there could be other property damage.

The device is only to be used for its intended purpose and only when it is in safe working condition.

Malfunctions which could affect safety must be corrected immediately.

3 Technical specifications

Manufacturer	Andreas Hettich GmbH & Co. KG D-78532 Tuttlingen			
Model	ROTO SILENTA 630 RS			
Type	5005, 5005-20, 5005-50, SK 48.06-30	5005-80	5005-90	5005-60
Mains voltage (± 10%)	400 V 3~ +N			
Mains frequency	50 - 60 Hz			
Connected load	9700 VA	9400 VA	6600 VA	
Current consumption	14 A	13.5 A	9.5 A	
Cooling medium	R 404A			----
Max. capacity	12000 ml			
Allowed density	1.2 kg/dm ³			
Speed (RPM)	6000			
Force (RCF)	6498			
Kinetic energy	215000 Nm			
Obligatory inspection (BGR 500)	Yes			
Ambient conditions (EN / IEC 61010-1) – Set-up site – Altitude – Ambient temperature – Humidity – Excess-voltage category (IEC 60364-4-443) – Pollution degree	Indoors only Up to 2000 m above sea level 5°C to 40°C Maximum relative humidity 80% for temperatures up to 31°C, linearly decreasing to 50% relative humidity at 40°C. II 2			
Device protection class	I			
Not suitable for use in explosion-endangered areas.				
EMC – Emitted interference, Interference immunity	EN / IEC 61326-1, Class B			
Noise level (dependent on rotor)	≤ 62 dB(A)		≤ 56 dB(A)	
Dimensions – Width – Depth – Height	813 mm 1015 mm 1050 mm 973 mm			
Weight	approx. 355 kg	approx. 367 kg	approx. 306 kg	approx. 303 kg

Manufacturer	Andreas Hettich GmbH & Co. KG D-78532 Tuttlingen		
Model	ROTO SILENTA 630 RS		
Type	5005-08	5005-88	5005-98
Mains voltage	208–220 V +6/-10% 3~ (+N) +PE		
Mains frequency	50 - 60 Hz		
Connected load	9000 VA	8700 VA	5200 VA
Current consumption	25 A	24 A	14 A
Cooling medium	R 404A		
Max. capacity	12000 ml		
Allowed density	1.2 kg/dm ³		
Speed (RPM)	6000		
Force (RCF)	6498		
Kinetic energy	215000 Nm		
Obligatory inspection (BGR 500)	yes		
Ambient conditions (EN / IEC 61010-1) – Set-up site – Altitude – Ambient temperature – Humidity – Excess-voltage category (IEC 60364-4-443) – Pollution degree	Indoors only Up to 2000 m above sea level 5°C to 30°C 5°C to 40°C Maximum relative humidity 80% for temperatures up to 31°C, linearly decreasing to 50% relative humidity at 40°C. II 2		
Device protection class	I		
Not suitable for use in explosion-endangered areas.			
EMC – Emitted interference, Interference immunity	EN / IEC 61326-1, Class B		
Noise level (dependent on rotor)	≤ 62 dB(A)		
Dimensions – Width – Depth – Height	813 mm 1015 mm 973 mm		
Weight	approx. 401 kg		

4 Notes on safety

No claim of warranty will be considered by the manufacturer unless ALL instructions in this manual have been followed.

- The centrifuge should be installed on a good, stable base.
- Before using the centrifuge absolutely check the rotor for firm placement.
- When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.
- Rotors, suspensions and accessories that possess traces of corrosion or mechanical damage or if their term of use has expired may not be used any longer.
- The centrifuge may no longer be put into operation when the centrifuging chamber has safety-related damages.
- With swing-out rotors the trunnions must be regularly lubricated (Hettich Lubricating Grease No. 4051) in order to ensure consistent swinging out of the hangers.
- For centrifuges without temperature control, when the room temperature is increased and/or if the device is frequently used, the centrifuging chamber could be heated up. Therefore, it can't be ruled out that the sample material might be changed due to the temperature.

- Before the initial operation of your centrifuge you should read and pay attention to the operating instructions. Only personnel that has read and understood the operating instructions are allowed to operate the device.
- Along with the operating instructions and the legal regulations on accident prevention, you should also follow the recognised professional regulations for working in a safe and professional manner. These operating instructions should be read in conjunction with any other instructions concerning accident prevention and environmental protection based on the national regulations of the country where the device is to be used.
- This centrifuge is a state-of-the-art piece of equipment which is extremely safe to operate. However, it can lead to danger for users or others if used by untrained staff, in an inappropriate way or for a purpose other than that it was designed for.
- The centrifuge must not be moved or knocked during operation.
- In case of fault or emergency release, never touch the rotor before it has stopped turning.
- To avoid damage due to condensate, when changing from a cold to a warm room the centrifuge must either heat up for at least 3 hours in the warm room before being connected to the mains, or run hot for 30 minutes in the cold room.
- Only the rotors and accessories approved by the manufacturer for this device may be used (see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories"). Before centrifuge vessels are used which are not listed in the chapter "Appendix, Rotors and accessories", the user must make sure they can be used by asking the manufacturer.
- The centrifuge rotor may only be loaded in accordance with the chapter "Loading the rotor".
- When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.
- The centrifuge may only be operated when the balance is within the bounds of acceptability.
- The centrifuge may not be operated in explosion-endangered areas.
- The centrifuge must not be used with:
 - inflammable or explosive materials
 - materials that react with one another producing a lot of energy.

- If users have to centrifuge hazardous materials or compounds contaminated with toxic, radioactive or pathogenic micro-organisms, they must take appropriate measures.
For hazardous substances centrifuge containers with special screw caps must strictly be used. In addition to the screw cap centrifuge containers, for materials in hazard category 3 and 4 a biosafety system must be used (see the World Health Organisation's "Laboratory Biosafety Manual").
In a biosafety system, droplets and aerosols are prevented from escaping by a bioseal (packing ring).
If the hanger of a biosafety system is used without the lid, the packing ring must be removed from the hanger in order to prevent the packing ring from being damaged during the centrifugation run. Damaged packing rings must not be used to seal the biosafety system.
Without the use of a biosafety system the centrifuge is not microbiologically sealed in the sense of the EN / IEC 610101-2-020 standard.
When closing a biosafety system, please observe the instructions provided in the chapter, "Closing Biosafety Systems".
For further details of available biosafety systems see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". If in doubt, you should obtain relevant information from the manufacturer.
- The centrifuge must not be operated with highly corrosive substances which could impair the mechanical integrity of rotors, hangers and accessories.
- Repairs must only be carried out by personnel authorised to do so by the manufacturer.
- Only original spare parts and original accessories licensed by the Andreas Hettich GmbH & Co. KG company are allowed to be utilised.
- The following safety regulations apply:
EN / IEC 61010-1 and EN / IEC 61010-2-020 as well as their national deviations.
- The safe operation and reliability of the centrifuge can only be guaranteed if:
 - the centrifuge is operated in accordance with the operating instructions,
 - the electrical installation on the site where the centrifuge is installed conforms to the demands of EN / IEC stipulations,
 - the tests for device safety required in the respective countries, e.g. in Germany in acc. with BGV A1 and BGR 500, are carried out by an expert.

5 Symbol meanings

Symbol on the device:
Attention, general hazard area.
Before using the device, make sure you read the operating instructions and observe the safety information!

Symbol in this document:
Attention, general hazard area.
This symbol refers to safety relevant warnings and indicates possibly dangerous situations.
The non-adherence to these warnings can lead to material damage and injury to personal.

Symbol on the device and in this document:
Beware of biohazard.

Symbol on the device and in this document:
Beware of hot surface.
Nonobservance of this warning can lead to material damage and personal injury.

Symbol on the device:
Equipotential: Connector (PE connector) for potential equalization (only for centrifuge with PE connector).

Symbol on the device:
Switch positions of the key-operated switch.

Symbol in this document:
This symbol refers to important circumstances.

Symbol on the device and in this document:
Symbol for the separate collection of electric and electronic devices according to the guideline 2002/96/EG (WEEE). The device belongs to Group 8 (medical devices).
Applies in the countries of the European Union, as well as in Norway and Switzerland.

6 Delivery checklist

- 4 Vibration absorber
- 1 Double open-end wrench 17/19 mm
- 1 Hexagonal pin-type spanner
- 1 Lubricating grease for trunnions
- 1 Operating instructions
- 3 Program-data sheet
- 1 Information on how to remove the packing
- 1 Installation information sheet

The rotor(s) and associated accessories are included in the delivery in the quantity ordered.

7 Unpacking the centrifuge and installation

Unpack and install the centrifuge in accordance with the instruction leaflets AH5005XX and AH5005-01XX provided.

8 Initial operation

- According to the laboratory instrument standards EN / IEC 61010-2-020 an emergency switch to separate power supply in the event of a failure must be installed in the building electrical system. This switch has to be placed remote from the centrifuge, preferred outside of the room in which the centrifuge is installed or near by the exit of this room.
- According to the laboratory instrument standards IEC 61010-1 a switch to separate the power supply agreeable to IEC 60947-5-3 must be installed in the building electrical system. This switch has to be placed close to the centrifuge, must be easily reachable for the user and has to be characterized as a separator for this unit. The switch must also be able to be secured against being switched on again.
- Types 5005-90, 5005-08, 5005-88 and 5005-98 are permanently connected devices.
- Types 5005-80, 5005-90, 5005-08, 5005-88 and 5005-98 may only be installed by authorized skilled personnel.

These devices must be installed according to the enclosed set-up and installation instructions, AH5005-02EN.

- Position the centrifuge in a stable and level manner in a suitable place. During set-up, the required safety margin of 300 mm around the centrifuge is to be kept according to EN / IEC 61010-2-020.**

When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.

- Do not place any object in front of the ventiduct.
Keep a ventilation area of 300 mm around the ventiduct.
- Centrifuge with PE connector:
If required, connect the PE connector at the back of the device with an additional medical potential equalization system.
- Check whether the mains voltage tallies with the statement on the type plate.
- For types 5005, 5005-50, 5005-80 and SK 48.06-30:
Connect the centrifuge with the connection cable to a standard mains socket. For connection ratings refer to Chapter "Technical specification".
- Switch on the mains switch in the building installation.
- Switch on the switch on the front panel. Switch position "I".
The following display follows:
1. Type of centrifuge, 2. The last rotor code detected by the rotor identification and the maximum speed of the rotor (n-max-Rotor), 3. The programme version, 4. **OPEN** **OEFFNEN**.
- Open the lid.
The centrifuge data of the last used program or program 1 will be displayed.

9 Opening and closing the lid

9.1 Opening the lid

The lid can only be opened when the centrifuge is switched on and the rotor is at rest. If it cannot be opened under these circumstances, see the section on "Emergency release".

- Press the button  on the front panel. The lid unlocks via the motor, the illuminated pushbutton  extinguishes and in the display the symbol  changes to the symbol .


9.2 Closing the lid

Do not put your fingers between lid and housing.
Do not bang the lid shut.

- Place the lid and slightly press it down at the tubular handle. The locking action is effected by motor. The button  lights up and the symbol  in the display turns to the symbol .

10 Installation and removal of the rotor


- Clean the motor shaft (C) and the rotor drilling (A), and lightly grease the motor shaft afterwards. Dirt particles between the motor shaft and the rotor hinder a perfect seating of the rotor and cause an irregular operation.
- Place the rotor vertically on the motor shaft. The motor shaft dog (D) has to fit in the rotor slot (B). The alignment of the groove is labelled on the rotor.
- Tighten the rotor tension nut with the supplied wrench by turning in a clockwise direction.
- Check the rotor for firm seating.
- Loosening the rotor: Loosen the tension nut by turning in a counter clockwise direction, and turning until the working point for lifting. After passing the working point for lifting the rotor is loosened from the motor shaft cone. Turn the tension nut until the rotor is able to be lifted from the motor shaft.

11 Loading the rotor


Standard centrifuge containers of glass will not stand RCF values exceeding 4000 (DIN 58970, pg. 2).

- Check the rotor for firm seating.
 - With swing-out rotors all rotor positions must be lined with **identical** hangers. Certain hangers are marked with the number of the rotor position. These hangers may only be used in the respective rotor position. Hangers that are marked with a set number (e.g. S001/4) may only be used in the set.
 - The rotors and hangers may only be loaded symmetrically. The centrifuge containers have to be distributed evenly on all rotor positions. For authorised combinations see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".
- In the case of angle rotors all possible rotor positions must be loaded, see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".


Rotor is evenly loaded

Not permitted!
Rotor is not evenly loaded

- On certain hangers, the weight of the maximum load or the weight of the maximum load and the maximum weight of the completely loaded hanger is specified. This weight may not be exceeded. In case of exception, see chapter "Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm^3 ". The weight specified for the maximum loading includes the total weight of adapter, frame, centrifuging container and content.
- In containers with rubber inserts, the same number of rubber inserts must always be among the centrifuge containers.
- The centrifuge containers may only be filled outside of the centrifuge.
- The maximum filling quantity for the centrifuge containers specified by the manufacturer must not be exceeded.

In the case of angle rotors, the centrifuging vessels may only be filled so far that no fluid can be expelled from them while the centrifuge is running.

- When loading the angle rotors, no liquid may get into the angle rotors or centrifuging chamber.
- When loading the hangers and swing-out rotors and when swivelling out the hangers while the centrifuge is running, no liquid may enter the hangers or the centrifuging chamber.
- In order to maintain the weight differences within the centrifuge container as marginal as possible, a consistent fill level in the containers is to be heeded.


12 Closing biosafety systems

In order to guarantee that it is sealed, the lid to the biosafety system has to be firmly closed.

In order to prevent the packing ring from twisting when opening and closing the cover, the packing ring must be lightly rubbed with talcum powder or a rubber care product.

For further details of available biosafety systems see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". If in doubt, you should obtain relevant information from the manufacturer.

Lid with clip lock:

- Put the lid on the hanger and push down both clamps until they are underneath the loops on the hanger.
- Fold down both clamps as far as they can go. The clamps have to lie flat against the hanger.

13 Control and display devices

See figure on page 2.

Fig. 2: Display and operating panel

13.1 Status symbols

Lid open.

Lid closed.

Rotation indicator. The indication appears while the centrifuge is running as long as the rotor is turning.

STOP

Centrifugation run stopped or finished. Indication after completion of the centrifugation run as long as the rotor is turning. The indicator flashes following an emergency stop.

**LOCK 1,
LOCK 2**

Keyswitch position (for centrifuges with keyswitch only).

**LOCK 4,
LOCK 5**

Program interlocking when serial communication provided (for centrifuges with serial communication only).

PC, PE

Serial communication (for centrifuges with serial communication only).

Operating errors or faults that occur are shown symbolically in the display (see Chapter "Faults").

13.2 Control knob

For setting the individual parameters.

Turning anticlockwise reduces the value. Turning clockwise increases the value.

13.3 Keys and set-up options

Running time, parameter t/min:sec

1. Parameter t/min: Adjustable from 1 - 999 min, in 1 min steps.
2. Parameter t/ :sec Adjustable from 1 - 59 s, in 1 second steps.
3. Continuous running "----:--"

Interrogation of the integral RCF, parameter f_{RCF}.

1. Speed, parameter RPM

A numerical value from 50 RPM up to the maximum speed of the rotor (n-max-Rotor) can be set. Maximum speed of the rotor, see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". Adjustable in steps of 10.

2. Interrogation of the maximum speed of the rotor, parameter n-max-Rotor.

1. Relative centrifugal force, parameter RCF/RZB

A numerical value can be set, which gives a speed between 50 RPM and the maximum speed of the rotor (n-max-Rotor). Adjustable in steps of 1.

2. Interrogation of the maximum RCF of the rotor, parameter RCF-max-Rotor.

Run-up parameters

1. Run-up steps, parameter $\sqrt{1-9}$. Step 9 = shortest run-up time, ... Step 1 = longest run-up time.
2. Run-up time, parameter $\sqrt{\text{min:sec}}$. The time range that can be set is dependent upon the set speed.

Run-down parameters

1. Braking steps, parameter $\sqrt{0-9}$. R = linear braking curve, B = similar to an exponential braking curve. Step R9, B9 = short run-down time, ... Step R1, B1 = long run-down time, Step R0 = unbraked run-down.
2. Run-down time, parameter $\sqrt{\text{min:sec}}$. The time range that can be set is dependent upon the set speed.
3. Brake de-energisation speed, parameter $n^{(r)}/\text{RPM}$
After reaching this speed, unbraked run-down occurs.

1. Temperature setpoint, parameter $T/^{\circ}\text{C}$ (for centrifuges with cooling only). Adjustable from -20°C to $+40^{\circ}\text{C}$, in 1°C steps (heating/cooling option settable between -20°C and $+60^{\circ}\text{C}$ / $+90^{\circ}\text{C}$). The lowest temperature that can be achieved is dependent upon the rotor (see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").
2. Centrifuging radius, parameter r/mm . Entry in mm. Centrifuging radius, see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

START

1. Start centrifugation run. The rotation indicator  appears.
2. Acceptance of changes during the centrifugation run.

STOP

Stop centrifugation run.

The rotor decelerates with the preset run-down parameters. Pressing the key twice will initiate an emergency stop.

Select program location, parameter PROG-Nr.

Calling up of programs.

Saving of programs.

89 programs can be saved (program locations 1 to 89). Note: The program locations "----" and 90 to 99 are used as automatic buffers (see Chapter "Automatic buffer"). No programs can be saved in these program locations.

14 Entering centrifugation parameters

It is only possible to enter a parameter when the input field is displayed inversely (dark background). An inverse input field will go off automatically after 10 seconds.

14.1 Running time

14.1.1 Time selection

- Press the  key repeatedly until the parameter $t/\text{min:}$ or $t/\text{:sec}$ is displayed. The input field will be displayed inversely.
- Set the desired value by means of the control knob .

14.1.2 Continuous running

- Select the parameters $t/\text{min:}$ and $t/\text{:sec}$ one after the other (see Chapter "Time selection"), and set both to zero using the control knob . "---:--" will be displayed in the input field.

14.2 Speed (RPM)

- Press the  key repeatedly until the parameter RPM is displayed and the input field is shown inversely.
- Set the desired value by means of the control knob .

14.3 Relative centrifugal force (RCF/RZB)

- Press the  key repeatedly until the parameter RCF/RZB is displayed and the input field is shown inversely.
- Set the desired value by means of the control knob .

The centrifuging radius indicator will flash while the adjustment is being made.

14.4 Run-up and run-down parameters

The run-up and run-down parameters that have been set are displayed.

x: 1-9 = run-up step, t = run-up time

y: R1-R9, B1-B9 = braking step, R0 = unbraked run-down, t = run-down time, $n^{(c)}$ = brake de-energisation speed

14.4.1 Run-up step

- Press the  key repeatedly until the parameter ₁₋₉ or  is displayed. The input field will be displayed inversely.
- Set the desired step by means of the control knob .

14.4.2 Run-up time

- Press the  key repeatedly until the parameter _{min:sec} is displayed. The input field will be displayed inversely.
- Set the desired value by means of the control knob .

If a run-up time is set that is longer than the running time, the centrifugation run will end before the set speed is reached.

14.4.3 Braking step

- Press the  key repeatedly until the parameter ₀₋₉ or  is displayed. The input field will be displayed inversely.
- Set the desired step by means of the control knob .

B-braking steps can only be set for special rotors.

14.4.4 Run-down time

If a brake de-energisation speed has been set, it is not possible to set a run-down time.

- Press the  key repeatedly until the parameter _{min:sec} is displayed. The input field will be displayed inversely.
- Set the desired value by means of the control knob .

14.4.5 Brake de-energisation speed

- Press the  key repeatedly until the parameter $n^{(c)}/\text{RPM}$ is displayed. The input field will be displayed inversely.
- Set the desired value by means of the control knob .

14.5 Radius/Temperature

14.5.1 Radius

- Press the  key repeatedly until the parameter r/mm is displayed and the input field is shown inversely.
- Set the desired value by means of the control knob .

The value of the RCF/RZB adapts automatically when the radius is changed; this is shown by the indication flashing.

14.5.2 Temperature (for centrifuges with cooling only)

- Press the  key repeatedly until the parameter T/°C is displayed and the input field is shown inversely.
- Set the desired value by means of the control knob .

14.6 Automatic buffer

The buffer comprises the program locations "----" and 90 to 99. Every time a centrifugation run is started, modified centrifugation data are automatically stored in program location "----". The modified centrifugation data for the last 11 centrifugation runs are stored in the buffer and can be called up (see Chapter "Calling up programs").

15 Programming

15.1 Entering / modifying programs

- Set the desired parameters (see Chapter "Entering centrifugation parameters").
- Press the **PROG** key in order to select the parameter PROG-Nr. The input field will be displayed inversely.
- Set the desired program location by means of the control knob $\odot$.
If the program location indicator flashes, then this program location has already had centrifugation data assigned to it. In this case, set a free program location or overwrite the centrifugation data by continuing.
- Press the **STO** key once to save the settings in the desired program location.
Press the **STO** key twice to overwrite centrifugation data already stored.

15.2 Calling up programs

- Press the **PROG** key in order to select the parameter PROG-Nr. The input field will be displayed inversely.
- Set the desired program location by means of the control knob $\odot$.
- Press the **RCL** key. The centrifugation data for the selected program location will be displayed.

16 Centrifugation

When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.

If the permissible weight difference within the rotor loading has been exceeded, the drive shuts down during the start-up, and **IMBALANCE / UNWUCHT** is displayed.

If the speed in the selected program is higher than the maximum speed of the rotor (n-max-Rotor), it will not be possible to start a centrifugation run. The message **N > ROTOR MAX 96** will be displayed (see Chapter "Faults").

A centrifugation run can be stopped at any time by pushing the key **STOP**.

It is possible to switch over between RPM and RCF display at any time using the **n** and **RCF** keys. When working with the RCF display, it is necessary to enter the centrifuging radius.

If **≤ OPEN ≤ OEFFNEN** is displayed, a further operation of the centrifuge is only possible after opening the lid once.

If **R xx n-max xxxxx** is displayed, then no centrifugation run has taken place as the rotor was changed beforehand, refer to Chapter "Rotor Identification".

- Switch on the mains supply switch. Switch position I.
- Load the rotor and close the centrifuge lid.

16.1 Centrifugation with time selection

- Set the time or call up a program with a pre-set time (see Chapter "Entering centrifugation parameters" or "Calling up programs").
- Press the **START** key. The rotation indicator $\odot$ will appear as long as the rotor is turning.
- When the time has expired or if the centrifugation run is interrupted by pressing the **STOP** key, run-down will take place with the selected run-down parameters. When the rotor comes to a stop, the $\square$ symbol will flash in the display until the lid is opened.

During the centrifugation run, the speed of the rotor or the RCF value resulting from this, the sample temperature (for centrifuges with cooling only) and the remaining time will be displayed.

16.2 Continuous running

- Select the ---:-- symbol or call up a continuous running program (see Chapter "Entering centrifugation parameters" or "Calling up programs").
- Press the **START** key. The rotation indicator $\odot$ will appear as long as the rotor is turning. The time count starts from 00:00.
- Press the **STOP** key to stop the centrifugation run. Run-down will take place with the selected run-down parameters. When the rotor comes to a stop, the $\square$ symbol will flash in the display until the lid is opened.

During the centrifugation run, the speed of the rotor or the RCF value resulting from this, the sample temperature (for centrifuges with cooling only) and the time run will be displayed.

17 Changing settings during the centrifugation run.

The running time, the speed, the relative centrifugal force (RCF/RZB), the run-up and run-down parameters as well as the temperature (for centrifuges with cooling only) can be changed while the centrifuge is running.

The parameters can only be changed one at a time and in succession.

- Select the desired parameter and change the value using the control knob $\odot$ (see Chapter "Entering centrifugation parameters").
- Press the **START** key. The modified setting will be stored in program location "----" (see Chapter "Automatic buffer"). The original program will not be overwritten.

18 Integral RCF (JRCF)

The integral RCF (JRCF) is a measure for the sedimentation effect ($\int n^2 dt$). This numerical value is used to compare centrifugation runs. To interrogate the integral RCF, press and hold the **JRCF** key.

19 Displaying the maximum speed of the rotor

- Press the **n** key repeatedly until the parameter RPM is displayed and the input field is shown inversely.
- Press the **n** key once more and hold it pressed; the maximum speed of the rotor (n-max-Rotor) will be displayed.

20 Displaying the maximum RCF of the rotor

- Press the **RCF** key repeatedly until the parameter RCF/RZB is displayed and the input field is shown inversely.
- Press the **RCF** key once more and hold it pressed; the maximum RCF of the rotor (RCF-max-Rotor) will be displayed.

21 Emergency stop

- Press the **STOP** key twice. The **STOP** symbol will flash in the display.

In the case of an emergency stop, run-down takes place with braking step R9 (shortest run-down time). If braking step R0 was selected, the run-down time will be longer than with braking step R9 for technical reasons.

22 Audible signal

The audible signal sounds in accordance with the following philosophy:

OFF	• at 2 s intervals on the occurrence of a fault
ON1	• at 2 s intervals on the occurrence of a fault • at 30 s intervals on completion of a centrifugation run and the rotor coming to rest.
ON2	• at 2 s intervals on the occurrence of a fault • at 30 s intervals on completion of a centrifugation run and the rotor coming to rest. • every time a key is pressed

The audible signal is silenced by opening the lid or pressing any key.

The signal can be set as follows when the rotor is at standstill.

- Open the lid.
- Press and hold the **⏏** key for 8 seconds.
After 8 seconds **SOUND / BELL XXX** will be displayed.
- Set the desired function (OFF, ON1, ON2) by means of the control knob $\odot$.
- Press the **START** key to save the setting.
***** ok ***** will be displayed briefly as a confirmation.

23 Interrogating the operating hours

It is only possible to interrogate the operating hours when the rotor is at standstill.

- Open the lid.
- Press and hold the **⏏** key for 8 seconds.
After 8 seconds **SOUND / BELL XXX** will be displayed.
- Press the **⏏** key once more.
The operating hours (**CONTROL:**) for the centrifuge will be displayed.
The operating hours display will go off automatically after 10 seconds.

24 Setting the date and time

It is only possible to set the date and time when the rotor is at standstill.

- Open the lid.
- Press and hold the **[]** key for 8 seconds.
After 8 seconds **SOUND / BELL XXX** will be displayed.
- Press the **[]** key twice more.
The date and time will be displayed (a: year, mon: month, d: day, h: hours, min: minutes).
- Press the **[]** key repeatedly until the input field for the desired parameter is shown inversely.
- Set the desired value by means of the control knob **○**.
- Press the **[START]** key to save the setting.
***** ok ***** will be displayed briefly as a confirmation.
- To exit the date and time set-up, press any key apart from the **[]**, **[]** and **[START]** keys.

25 Centrifugation data displayed after switching on.

After switching on, the centrifugation data for program 1 or for the program last used will be displayed. This can be set as follows:

- Open the lid.
- Switch the mains supply switch off and on again. Switch position I.
- Press the **[STOP]** key when the first visual change appears in the display (inverse display).
PROGRAM 1, LAST PROGRAM will be displayed.
- Set the desired function by means of the control knob **○**.
- Press the **[START]** key to save the setting.
***** ok ***** will be displayed briefly as a confirmation.

26 Immediate display of the centrifugation data after switching on

- Switch on the mains supply switch. Switch position I.
- Press any key apart from the **[STOP]** key when the first visual change appears in the display (inverse display). The centrifugation data will be displayed immediately.

27 Key-operated switch (for centrifuges with key-operated switch only).

The following program interlocks can be set by the key-operated switch:

- Left key position: **LOCK 1** will be displayed.
Programs can be called up but not changed.
- Right key position: **LOCK 2** will be displayed.
No programs can be called up or changed.
- Middle key position: no status display.
No program interlock. Programs can be called up and changed.

28 Program linking (for centrifuges with program linking only).

Several centrifugation runs can be linked together with the help of the program linking function.

28.1 Linking programs

Program linking is only possible with programs in which run-up and braking steps have been set.

Before linking, the programs must be stored in the desired sequence either by entering the program or by calling up the program (see Chapter "Programming").

The program locations must be consecutive (e.g. program locations 10+11+12).

1. Press the **[PROG]** key in order to select the parameter PROG-Nr. The input field will be displayed inversely.
2. Set the program location for the first program (XX+) by means of the control knob **○**.
3. Press the **[RCL]** key. The centrifugation data for the selected program location will be displayed.
4. Press the **[PROG]** key twice in order to select the parameter PR-PART. The input field will be displayed inversely.
5. Press the **[STO]** key twice. The program will be linked and the program number for the next program location (+XX+) will be displayed.
6. Press the **[RCL]** key twice. The centrifugation data for the selected program location will be displayed.
7. Press the **[STO]** key twice. The program will be linked and the program number for the next program location (+XX+) will be displayed.
8. Keep repeating steps 6 and 7 until all programs have been linked.
9. Press the **[PROG]** key to exit. The program number for the last program (+XX) will be displayed.


28.2 Changing program links

- Call up the desired program (see Chapter "Calling up programs"), change the desired parameters (see Chapter "Entering centrifugation parameters"), and save the modified centrifugation data to the same program location (see Chapter "Entering / modifying programs"). As a result of saving, the program linking will be removed.
- Re-link the programs (see Chapter "Linking programs").

28.3 Centrifugation run with program linking

- Press the **PROG** key twice in order to select the parameter PR-PART. The input field will be displayed inversely.
- Set the program location for the first program (XX+) by means of the control knob $\odot$.
- Press the **RCL** key. The centrifugation data for the selected program location will be displayed.
- Press the **START** key. The rotation indicator $\odot$ will appear as long as the rotor is turning.

The run-up and braking steps for the program linking will be displayed:

- When the time in the last program has expired, run-down will take place with the braking step of the last program. If the centrifugation run is interrupted by pressing the **STOP** key, run-down will take place with the braking step of the program running at the time.

28.4 Deleting program linking

- Press the **PROG** key in order to display the parameter PROG-Nr. The input field will be displayed inversely.
- Set the program location for the first program (XX+) by means of the control knob $\odot$.
- Press the **RCL** key. The centrifugation data for the selected program location will be displayed.
- Press the **PROG** key twice in order to display the parameter PR-PART. The input field will be displayed inversely.
- Press the **STO** key twice.
- Press the **PROG** key.

29 Cooling (only in centrifuges with cooling)

The temperature set-point can be adjusted from -20°C to +40°C. On centrifuges with a heating/cooling option, the temperature set-point is settable from -20°C to +90°C. The lowest obtainable temperature is dependent on the rotor (see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories"). If the actual temperature and the set temperature deviate more than 5 °C, the temperature value in the display will start to flash on and off.

29.1 Standby-cooling

With rotor standstill and closed lid the centrifugal chamber is cooled to the pre-selected temperature. The temperature set-point is shown in the display.

29.2 Pre-cooling the rotor

For rapid pre-cooling of the unloaded rotor and the accessories, it is recommended that the centrifuge is run with the continuous running settings and a speed of ca. 20% of the maximum rotor speed.

30 Heating (only on centrifuges with heating/cooling option)

During operation of the centrifuge the centrifuge chamber can be heated as required to the preselected temperature. The heating is switched off when the rotor is stationary.

Risk of burning! The surface temperature of the heating element in the centrifugal chamber of the centrifuge can reach 500°C / 932°F. Do not touch the heating element.

During a centrifugation run at very high temperature (e.g. +90°C), the inside of the device lid heats up. In this case, do not touch the inside of the lid.

31 Relative centrifugal force (RCF)

The relative centrifugal force (RCF) is given as a multiple of the acceleration of gravity (g). It is a unit-free value and serves to compare the separation and sedimentation performance.

These values are calculated using the formula below:

$$\text{RCF} = \left(\frac{\text{RPM}}{1000} \right)^2 \times r \times 1,118 \quad \Rightarrow \quad \text{RPM} = \sqrt{\frac{\text{RCF}}{r \times 1,118}} \times 1000$$

RCF = relative centrifugal force

RPM = rotational speed (revolutions per minute)

r = centrifugal radius in mm = distance from the centre of the turning axis to the bottom of the centrifuge.
For more on the centrifugal radius see the chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

The relative centrifugal force (RCF) stands in relation to the revolutions per minute and the centrifugal radius.

32 Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm³

When centrifuging with maximum revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.

The speed must be reduced for materials or mixtures of materials with a higher density.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed (n}_{\text{red}}) = \sqrt{\frac{1.2}{\text{Greater density [kg/dm}^3]}} \times \text{maximum speed [RPM]}$$

e.g.: maximum speed RPM 4000, density 1.6 kg/dm³

$$n_{\text{red}} = \sqrt{\frac{1.2 \text{ kg/dm}^3}{1.6 \text{ kg/dm}^3}} \times 4000 \text{ RPM} = 3464 \text{ RPM}$$

In the exceptional case that the maximum loading indicated on the hanger is exceeded, the speed must also be reduced.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed (n}_{\text{red}}) = \sqrt{\frac{\text{maximum load [g]}}{\text{actual load [g]}}} \times \text{maximum speed [RPM]}$$

e.g.: maximum speed RPM 4000, maximum load 300 g, actual load 350 g

$$n_{\text{red}} = \sqrt{\frac{300 \text{ g}}{350 \text{ g}}} \times 4000 \text{ RPM} = 3703 \text{ RPM}$$

If in doubt you should obtain clarification from the manufacturer.

33 Rotor recognition

Rotor recognition is carried out after each start of the centrifugation run.

If the rotor has been changed, the centrifugation run is canceled after rotor recognition. The rotor code (R xx) and the maximum speed (n-max=xxxxx) of the rotor are displayed.

A further operation of the centrifuge is only possible after a single opening of the lid.

If the maximum speed of the rotor being used is less than the set speed, the speed is limited to the rotor's maximum speed.

34 Emergency release

During a power failure the lid cannot be unlocked by motor. An emergency release has to be executed by hand.

For emergency release disconnect the centrifuge from the mains. In the case of permanently connected devices, switch off the mains switch to disconnect the mains supply of the device in the building installation and secure it against being switched on again, e.g. by locking the switch.
Open the lid only during rotor standstill.

See figure on page 2.

- Switch off the switch on the front panel (switch position "0").
- Look through the window in the lid to be sure that the rotor has come to a standstill.
- Insert the hexagon socket head wrench horizontally into the drilling (Fig. 1, A) and carefully rotate by half a turn in clockwise direction until the lid can be opened.
- Pull the hexagon socket head wrench out of the drilling again.
- After turning the centrifuge on again, press the button  so that the motor-driven lid locking once again assumes the normal position (opened).

35 Maintenance and servicing

The device can be contaminated.

Pull the mains plug before cleaning. In the case of permanently connected devices, switch off the mains switch to disconnect the mains supply of the device in the building installation and secure it against being switched on again, e.g. by locking the switch.

Before any other cleaning or decontamination process other than that recommended by the manufacturer is applied, the user has to check with the manufacturer that the planned process does not damage the device.

- Centrifuges, rotors and accessories must not be cleaned in rinsing machines.
- They may only be cleaned by hand and disinfected with liquids.
- The water temperature must be between 20 – 25°C.
- Only detergents/disinfectants may be used which:
 - have a pH between 5 - 8
 - do not contain caustic alkalis, peroxides, chlorine compounds, acids and alkaline solutions
- In order to prevent appearances of corrosion through cleaning agents or disinfectants, the application guide from the manufacturer of the cleaning agent or disinfectant are absolutely to be heeded.

35.1 Centrifuge (housing, lid and centrifuging chamber)

35.1.1 Surface cleaning and care

- Clean the centrifuge housing and the centrifuging chamber regularly, using soap or a mild detergent and a damp cloth if required. For one thing, this services purposes of hygiene, and it also prevents corrosion through adhering impurities.
- Ingredients of suitable detergents: soap, anionic tensides, non-ionic tensides.
- After using detergents, remove the detergent residue by wiping with a damp cloth.
- The surfaces must be dried immediately after cleaning.
- In the event of condensation water formation, dry the centrifugal chamber by wiping out with an absorbent cloth.
- Lightly rub the rubber seal of the centrifuge chamber with talcum powder or a rubber care product after each cleaning.
- The centrifuging chamber is to be checked for damage once a year.

If damage is found which is relevant to safety, the centrifuge may no longer be put into operation. In this case, notify Customer Service.

35.1.2 Surface disinfection

- If infectious materials penetrates into the centrifugal chamber this is to be disinfected immediately.
- Ingredients of suitable disinfectants: ethanol, n-propanol, ethyl hexanol, anionic tensides, corrosion inhibitors.
- After using disinfectants, remove the disinfectant residue by wiping with a damp cloth.
- The surfaces must be dried immediately after disinfecting.

35.1.3 Removal of radioactive contaminants


- The agent must be specifically labelled as being an agent for removing radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants: anionic tensides, non-ionic tensides, polyhydrated ethanol.
- After removing the radioactive contaminants, remove the agent residue by wiping with a damp cloth.
- The surfaces must be dried directly after removing the radioactive contaminants.

35.2 Rotors and Attachments

35.2.1 Cleaning and care

- In order to avoid corrosion and changes in materials, the rotors and accessories have to be cleaned regularly with soap or with a mild cleaning agent and a moist cloth. Cleaning is recommended at least once a week. Contaminants must be removed immediately.
- Ingredients of suitable detergents: soap, anionic tensides, non-ionic tensides.
- After using detergents, remove detergent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after cleaning.
- Angle rotors, container and hanger made of aluminium are to be lightly greased after drying using acid-free grease, e.g. vaseline.
- In the case of biosafety systems (for further details of available biosafety systems see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories"), the packing rings must be checked and cleaned regularly (weekly). The sealing ring is to be replaced immediately upon indication of crack formation, embrittlement or abrasive wear. Lightly rub the packing rings with talcum powder or a rubber care product after each cleaning.
- In order to prevent corrosion as a result of moisture between the rotor and the motor shaft, the rotor should be disassembled and cleaned at least once a month, and the motor shaft should be lightly greased.
- The rotors and accessories have to be checked weekly for wear and corrosion. For swing-out rotors, it is important to check the area of the lifting lugs, for hangers, the grooves and the base should be checked for cracks.

Example: Crack in the groove area:

Rotors and attachments may no longer be utilised upon indication of wear and tear or corrosion.

- Check the firm seating of the rotor on a weekly basis.

35.2.2 Disinfection

- If infectious material should get on the rotors or accessories, they must be appropriately disinfected.
- Ingredients of suitable disinfectants: ethanol, n-propanol, ethyl hexanol, anionic tensides, corrosion inhibitors.
- After using disinfectants, remove disinfectant residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after disinfection.

35.2.3 Removal of radioactive contaminants

- The agent must be specifically labelled as being an agent for the removal of radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants: anionic tensides, non-ionic tensides, polyhydrated ethanol.
- After removing the radioactive contaminants, remove agent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after removing the radioactive contaminants.

35.2.4 Trunnions

With swing-out rotors the trunnions must be regularly lubricated (Hettich Lubricating Grease No. 4051) in order to ensure consistent swinging out of the hangers.

35.2.5 Rotors and accessories with limited service lives

The use of certain rotors, hangers and accessory parts is limited by time.

These are marked with the maximum permitted number of operating cycles or with an expiration date and the maximum permitted number of operating cycles or just with the expiration date; e.g.:

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For safety reasons, rotors, hangers and accessory parts may no longer be used if either the indicated maximum number of operating cycles or the indicated expiration date has been reached.

35.3 Autoclaving

The following accessory can be autoclaved at 121°C / 250°F (20 min):

- Swing-out rotors
- Angle rotors made of aluminum
- Hanger made of metal
- Lid with biocontainment
- Stands
- Reductions

Otherwise you must ask the manufacturer.

No statement can be made about the degree of sterility.

The lids of the rotors and containers must be removed prior to autoclaving.

Autoclaving accelerates the ageing process of plastics. In addition, autoclaving may discolour plastics.

We recommend that the packing rings of the bio-safety system be replaced after autoclaving.

35.4 Centrifuge containers

- With leakiness or after the breakage of centrifuging containers broken container parts and leaked centrifugation material are to be completely removed.
- The rubber inserts as well as the plastic sleeves of the rotors are to be replaced after a glass breakage.

Remaining glass splitters cause further glass breakage!

- If this concerns infectious material, a disinfection process is to be executed immediately.

36 Faults

If the fault cannot be eliminated with the help of the fault table, please inform Customer Service.

Please specify the type of centrifuge and the serial number. Both numbers can be found on the name plate of the centrifuge.

Perform a MAINS RESET:

- Switch off the mains switch (switch position "0").
- Wait at least 10 seconds and then switch on the mains switch again (switch position "1").

Message / fault		Cause	Remedy
No display	---	No voltage. Overvoltage protection tripped out (only for types 5005-08).	– Check supply voltage. – Mains switch ON.
TACHO - ERROR	01	Faulty speedometer.	– Open the cover. – Switch off the mains switch (switch position "0"). – Wait at least 10 seconds. – Turn the rotor vigorously by hand. – Switch on the mains switch again (switch position "1"). The rotor must turn during switch-on.
	02	No rotor installed. Defective motor, frequency converter or drive.	
IMBALANCE / UNWUCHT	---	The rotor is unevenly loaded.	– Open lid. – Check the loading of the rotor, see chapter "Loading the rotor". – Repeat the centrifugation run.
CONTROL - ERROR	04, 06 - 09	Error in lid locking or lid closure.	– Perform a MAINS RESET.
N > MAX	05	Rotation too fast	
N < MIN	13	Rotation too slow	
ROTORCODE	10	Incorrect rotor coding	
MAINS INTERRUPT	---	Power failure during the centrifugation run. (The centrifugation run was not finished.)	– Open lid. – Push START button. – Repeat the centrifugation run if necessary.
VERSIONS-ERROR	12	Mismatch between electronic components	– Perform a MAINS RESET.
SER I/O - ERROR	30 - 38	Error / defective interface	
° C * - ERROR	50 - 56, 58	Error / defective cooling	
LOCK - ERROR	57	Error / defective program locking	
FU / CCI - ERROR	60 - 83	Error / defective motor control	
CONTROL-ERROR	26, 90 - 95, 97 - 99	Error / defective control unit	
N > ROTOR-MAX	96	Speed in the selected program greater than the maximum speed of the rotor (n-max-Rotor).	– Check the set speed. – Reduce the set speed
		The rotor has been changed. The maximum speed of the installed rotor is higher than that of the rotor that was previously used and this has not yet been detected by rotor recognition.	– Set a speed up to the maximum speed of the rotor that was previously used. Press the START key to run rotor recognition. See chapter, "Rotor Recognition".

37 Returning Devices

If the device or its accessories are returned to Andreas Hettich GmbH & Co. KG, in order to provide protection for people, the environment and materials, it has to be decontaminated and cleaned before being shipped.

We reserve the right to refuse contaminated devices or accessories.

Costs incurred for cleaning and disinfection are to be charged to the customer.

We ask for your understanding in this matter.

38 Disposal

Before disposal, the device must be decontaminated and cleaned to protect people, the environment and property.

When you are disposing of the device, the respective statutory rules must be observed.

Pursuant to guideline 2002/96/EC (WEEE), all devices supplied after August 13, 2005 may not be disposed as part of domestic waste. The device belongs to group 8 (medical devices) and is categorized in the business-to-business field.

The icon of the crossed-out trash can shows that the device may not be disposed as part of domestic waste.

The waste disposal guidelines of the individual EC countries might vary. If necessary, contact your supplier.

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




- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 9) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".






4174	5052 + 5057								
Ausschwingrotor 4-fach / Swing out rotor 4-times  ∠ 90°									
	mit Bioabdichtung / with bio-containment ⁹⁾								
	4226	4232	4233		4238		4241	4245	
									
	Microtainer	0509	4234	0551	5127	---		0513	
			 ⁸⁾	 ⁸⁾	 ⁸⁾				
Kapazität / capacity	ml	0,8	15	750	600	250	290	25	50
Maße / dimensions Ø x L	mm	8 x 45	17 x 120	93 x 137	93 x 134	62 x 122	62 x 137	25 x 90	29 x 115
Anzahl p. Rotor / number p. rotor		312	92	4	4	4	4	32	32
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	²⁾	5003	5184	5184	5184	5003	5003	5003	5184
Radius / radius	mm	221	229	229	229	221	221	221	269
 9 (97%)	sec	125							
 9	sec	197							
Temperatur / temperature	°C ¹⁾	-1							

4174	5052 + 5057				5052			
Ausschwingrotor 4-fach / Swing out rotor 4-times  								
	mit Bioabdichtung / with bio-containment ⁹⁾							
	4249	4258	4258 + 4449	SK 09.00	6322			
								
		0512	Corning		Corning			
	⁸⁾							
Kapazität / capacity	ml	50	750	500	14	250		
Maße / dimensions Ø x L	mm	29 x 115	97 x 152	96 x 147	16,5 x 106	60 x 162		
Anzahl p. Rotor / number p. rotor		24	4	4	92	6		
Drehzahl / speed	RPM	4500	4500	4500	4500	4500		
RZB / RCF	²⁾	5094	5184	5184	5117	5003		
Radius / radius	mm	225	229	229	226	221		
 9 (97%)	sec	125						
 9	sec	197						
Temperatur / temperature	°C ¹⁾	-1						

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- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 3) tube will not stand RCF values exceeding 700
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.
- 9) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".

4174	4522-A							
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$								
	4213		4213-93		4214			
								
	0501	0578		0553				
								
Kapazität / capacity	ml	6	7	4,5 - 5	5	2,7 - 3	7,5 - 8,5	8,5 - 10
Maße / dimensions $\varnothing \times L$	mm	12 x 82	12 x 100	11 x 92	12 x 75	11 x 66	15 x 92	16 x 100
Anzahl p. Rotor / number p. rotor		192	192	192	192	192	120	120
Drehzahl / speed	RPM	4500						
RZB / RCF	²⁾	5003						
Radius / radius	mm	221						
 9 (97%)	sec	125						
 9	sec	197						
Temperatur / temperature	°C ¹⁾	-1						

4174	4522-A							
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$								
	4214		4214-93		4215		4216	4218
								
		0518			0519	---	0521	0523
								
Kapazität / capacity	ml	10	15	4,5 - 5	4 - 7	25	30	50
Maße / dimensions $\varnothing \times L$	mm	15 x 102	17 x 100	15 x 75	16 x 75	24 x 100	25 x 110	34 x 100
Anzahl p. Rotor / number p. rotor		120	120	120	120	44	44	24
Drehzahl / speed	RPM	4500						
RZB / RCF	²⁾	5003						
Radius / radius	mm	221						
 9 (97%)	sec	125						
 9	sec	197						
Temperatur / temperature	°C ¹⁾	-1						

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



1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature













2) Centrifuge containers of glass will not stand RCF values exceeding 4000




3) tube will not stand RCF values exceeding 700


8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.

10) May only be used in bucket 4524-A and with both sides loaded. If the holding pins are used the max. RCF of 1000 may not be exceeded.

	Position of the holding pins	Distance of the holding pins from the bottom of the inserts (inside)
4529-O	upper section	199.5 mm
4529-M	middle section	182.0 mm
4529-U	lower section	164.5 mm

4174	4524-A				4546-A	4591-A	4595-A		
Ausschwingrotor 4-fach / Swing out rotor 4-times  ∠ 90°									
	Baxter/ 4525	Biopack/ 4525	Fenwal 4525 6)	4592-A	4559-H/ 4559-N	4592-A		4596	
									
	2 x 4526	2 x 4527	2 x 4528	4-fach/ 4-times	1-fach/ 1-times	4-fach/4-times	4-fach/ 4-times	0550	4597
								0550 8)	
Kapazität / capacity	ml	500	500	500	500 750	450	500 450	2000	100
Maße / dimensions ∅ x L	mm	-----	-----	-----	-----	-----	-----	180x150x100	-----
Anzahl p. Rotor / number p. rotor		8	8	8	8	8	8	4	16
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	3500	1000
RZB / RCF		5683	5683	5683	5683	5479	5705	5569	272
Radius / radius	mm	251	251	251	251	242	252	246	243
 9 (97%)	sec	125							
 9	sec	197							
Temperatur / temperature	°C 1)	10							

4174	4579-A + 4255 / 4255-P 4)									
Ausschwingrotor 4-fach / Swing out rotor 4-times   90°										
	4432		4433				4434			
										
	2078	0536		0578			0500	2079	0518	
										
Kapazität / capacity	ml	1,5	2,0	5	7	2,7 - 3	4,5 - 5	9	10	15
Maße / dimensions Ø x L	mm	11 x 38	12 x 75	12 x 100	11 x 66	11 x 92	14 x 100	17 x 70	17 x 100	
Anzahl p. Rotor / number p. rotor		168	120	120	120	120	120	76	76	76
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	²⁾	5003	4935	4935	4935	4935	4935	5094	5094	5094
Radius / radius	mm	221	218	218	218	218	218	225	225	225
 9 (97%)	sec	125								
 9	sec	197								
Temperatur / temperature	°C ¹⁾	-11								

4174	4579-A + 4255 / 4255-P 4)							
Ausschwingrotor 4-fach / Swing out rotor 4-times   90°								
	4434							
								
								
Kapazität / capacity	ml	10	10	8	4 - 5,5	9 - 10	4 - 7	8,5 - 10
Maße / dimensions $\varnothing$ x L	mm	15 x 102	16 x 80	16 x 81	15 x 75	16 x 92	16 x 75	16 x 100
Anzahl p. Rotor / number p. rotor		76	76	76	76	76	76	
Drehzahl / speed	RPM	4500						
RZB / RCF	²⁾	5094						
Radius / radius	mm	225						
 9 (97%)	sec	125						
 9	sec	197						
Temperatur / temperature	°C ¹⁾	-11						

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

1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)

2) Centrifuge containers of glass will not stand RCF values exceeding 4000

4) 4255-P: special surface treatment for highest hygienic requirements

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





- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 3) tube will not stand RCF values exceeding 700
- 4) 4255-P: special surface treatment for highest hygienic requirements
- 5) 4255 cannot be closed with the lid
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.

4174		4579-A + 4255 / 4255-P 4)						
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$								
		4466	SK 39.93	SK 61.98				
								
		0551	0530					
		 8)						
Kapazität / capacity	ml	600	250	50				
Maße / dimensions $\varnothing \times L$	mm	93 x 134	65 x 115	29 x 115				
Anzahl p. Rotor / number p. rotor		4	4	20				
Drehzahl / speed	RPM	4500	4500	4500				
RZB / RCF	2)	5501	5049	5003				
Radius / radius	mm	243	223	221				
 9 (97%)	sec	125						
 9	sec	197						
Temperatur / temperature	$^\circ\text{C}$ 1)	-11						

4174		4547-A + 5621							
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$									
				4449	4430			4432	
									
		0512	4239	Corning	Corning	Nalgene	Nunc	2078	0536
		 8)	 8)						
Kapazität / capacity	ml	750	1000	500	250	175	200	1,5	2,0
Maße / dimensions $\varnothing \times L$	mm	97 x 152	96 x 176	96 x 147	60 x 162	61,5 x 144,3	60 x 130	11 x 38	11 x 38
Anzahl p. Rotor / number p. rotor		4	4	4	4	4	4	168	252
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	2)	5592	5592	5592	5592	5592	5592	5094	5094
Radius / radius	mm	247	247	247	247	247	247	225	225
 9 (97%)	sec	125							
 9	sec	197							
Temperatur / temperature	$^\circ\text{C}$ 1)	-9							

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



1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature

2) Centrifuge containers of glass will not stand RCF values exceeding 4000

4) 4255-P: special surface treatment for highest hygienic requirements

8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.

4174	4547-A + 5621							
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$								
	4433				4434			
								
	0553	0578			0500	2079		
								
Kapazität / capacity	ml	5	7	2,7 - 3	4,5 - 5	9	10	10
Maße / dimensions $\varnothing \times L$	mm	12 x 75	12 x 100	11 x 66	11 x 92	14 x 100	17 x 70	16 x 80
Anzahl p. Rotor / number p. rotor		120	120	120	120	76	76	76
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	²⁾	5026	5026	5026	5026	5184	5184	5184
Radius / radius	mm	222	222	222	222	229	229	229
 9 (97%)	sec	125						
 9	sec	197						
Temperatur / temperature	$^\circ\text{C}$ ¹⁾	-9						

4174	4547-A + 5621							
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$								
	4434							
								






								
Kapazität / capacity	ml	8	4 - 5,5	7,5 - 8,2	9 - 10	4 - 7	8	8,5 - 10
Maße / dimensions $\varnothing \times L$	mm	16 x 81	15 x 75	15 x 92	16 x 92	16 x 75	16 x 125	16 x 100
Anzahl p. Rotor / number p. rotor		76	76	76	76	76	76	76
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	²⁾	5184	5184	5184	5184	5184	5184	5184
Radius / radius	mm	229	229	229	229	229	229	229
 9 (97%)	sec	125						
 9	sec	197						
Temperatur / temperature	$^\circ\text{C}$ ¹⁾	-9						

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











1) Lowest possible temperature during maximum speed, 1 h running time and 20 °C ambient temperature

2) Centrifuge containers of glass will not stand RCF values exceeding 4000

8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.

4174	4547-A + 5621				4523-A			
Ausschwingrotor 4-fach / Swing out rotor 4-times 								
	4443	4466	SK 39.93	SK 61.98	4516			
								
		0551	0530		4-fach / 4-times	1-fach / 1-times	1-fach / 1-times	
Kapazität / capacity	ml	290	600	250	50	500	750	1000
Maße / dimensions $\varnothing \times L$	mm	62 x 137	93 x 134	65 x 115	29 x 115	---	---	---
Anzahl p. Rotor / number p. rotor		4	4	4	20	4	4	4
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	
RZB / RCF	²⁾	5411	5592	5139	5094	5637	5637	5705
Radius / radius	mm	239	247	227	225	249	249	252
 9 (97%)	sec	125						
 9	sec	197						
Temperatur / temperature	°C ¹⁾	-9				3		

4174	SK 06.07 + 5629						
Ausschwingrotor 4-fach / Swing out rotor 4-times  ∠ 90°							
	4626						
							
	MTP	MTP	MS	CP	DWP	Microtest- platten / plate Terasaki	QP
							
Maße / dimensions ∅ x L mm	86x128x15	86x128x17,5	86x128x46	86x128x22	86x128x44,5	59x84x11	86x128x 83
Anzahl p. Rotor / number p. rotor	24	20	4	16	4	8	4
Drehzahl / speed RPM	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	4324	4324	4324	4324	4324	4324	4324
Radius / radius mm	191	191	191	191	191	191	191
 9 (97%) sec	125						
 9 sec	197						
Temperatur / temperature °C ¹⁾							

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










CP Kulturplatte /
Culture plate

DWP Deep Well Platte /
Deep well plate

MS Micronic System /
Micronic system

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

1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature

2) Centrifuge containers of glass will not stand RCF values exceeding 4000

7) only executable in centrifuges type SK 48.06-30

4176	5052 + 5057								
Ausschwingrotor 6-fach / Swing out rotor 6-times  									
	mit Bioabdichtung / with bio-containment ⁹⁾								
	4213			4213-93		4214			
									
	0501	0578				0507			
									
Kapazität / capacity	ml	6	7	4,5 - 5	5	2,7 - 3	15	7,5 – 8,5	8,5 -10
Maße / dimensions $\varnothing \times L$	mm	12 x 82	12 x 100	11 x 92	12 x 75	11 x 66	17 x 100	15 x 92	16 x 100
Anzahl p. Rotor / number p. rotor		288	288	288	288	288	180	180	180
Drehzahl / speed	RPM	4500							
RZB / RCF	²⁾	5818							
Radius / radius	mm	257							
 9 (97%)	sec	125							
 9	sec	197							
Temperatur / temperature	°C ¹⁾	14							

4176	5052 + 5057								
Ausschwingrotor 6-fach / Swing out rotor 6-times  $\angle 90^\circ$	 								
	mit Bioabdichtung / with bio-containment ⁹⁾								
	4214		4214-93		4215		4216	4218	
									
		0518			0519	---	0521	0523	
									
Kapazität / capacity	ml	10	15	4,5 - 5	4 - 7	25	30	50	100
Maße / dimensions $\varnothing \times L$	mm	15 x 102	17 x 100	15 x 75	16 x 75	24 x 100	25 x 110	34 x 100	44 x 115
Anzahl p. Rotor / number p. rotor		180	180	180	180	66	66	36	24
Drehzahl / speed	RPM	4500							
RZB / RCF	²⁾	5818							
Radius / radius	mm	257							
 9 (97%)	sec	125							
 9	sec	197							
Temperatur / temperature	°C ¹⁾	14							

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- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 3) tube will not stand RCF values exceeding 700
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.
- 10) May only be used in bucket 4524-A and with both sides loaded. If the holding pins are used the max. RCF of 1000 may not be exceeded.

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All rights reserved. No part of this publication may be reproduced without the prior written permission of the copyright owner.