**D4001 Procedure Cart
Operator's Manual**

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
Table of Contents		Page
1.0	Preface	1
1.1	Proprietary Information	1
1.2	Indications for Use / Intended Use	1
1.3	Contraindications	1
1.4	Operating Principle	1
1.5	Terminology	1
1.6	Conventions Used in this Guide	1
1.7	Environmental Directive	1
1.8	Warranty and Repair	1
1.9	Hotline / Technical Support	2
2.0	Introducing the D4001 Procedure Cart	3
2.1	Device Description	3
2.2	General Safety	3
2.3	Unpacking Cart Components and Accessories	4
3.0	General Assembly and Use	5
3.1	Locking the Wheels	5
3.2	Installing the IV Pole	5
3.3	Removing the IV Pole	5
3.4	Adjusting Shelf	5
3.5	Assembling Drawer	6
3.6	Installing Console-to-Shelf	6
3.7	Installing Handle with Cord Hooks	7
4.0	Maintenance	8
4.1	Cleaning	8
4.2	Preventative Maintenance	8
5.0	Technical Specifications	8
6.0	Environmental Limits	9
7.0	Symbols used on Labeling	9

1.0 Preface

1.1 Proprietary Information

This manual contains information deemed proprietary to CONMED Corporation. The information contained herein, including all of the designs and related materials, is the sole property of CONMED and/or its licensors. CONMED and/or its licensors reserve all patent, copyright and other proprietary rights to this document, including all design, manufacturing methodology and reproduction.

This document, and any related materials, is confidential and is protected by copyright laws and shall not be duplicated, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form or by any means, electronic, mechanical, magnetic, manual or otherwise, or disclosed to third parties, in whole or in part, without the prior express written consent of CONMED.

CONMED reserves the right to revise this publication and to make changes from time to time in the contents hereof without obligation to notify any person of such revision or changes, unless otherwise required by law.

© CONMED Corporation 2022. All Rights Reserved

1.2 Indications for Use / Intended Use

To support storage and use of arthroscopic equipment.

1.3 Contraindications

None known.

1.4 Operating Principle

The D4001 is used to store arthroscopic equipment for single or multiple surgical specialties and resides outside of the sterile field. They include adjustable shelves for equipment placement and storage and incorporate an IV pole, cord wraps, and cable management clips. Wheels that are built into the base of the cart allow ease of transport, three (green lever) wheels are locking type and one (red lever) wheel is directional type. See technical specification for more detail.

1.5 Terminology

Term	Definition
D4001	Procedure Cart

1.6 Conventions Used in this Guide

This section contains warnings and precautions that must be thoroughly understood before operating any of the equipment. Lack of understanding or adherence to these warnings and precautions may result in injury.

The words Warning, Precaution, and Note carry special meanings and they must be read carefully.

WARNING: A warning contains critical information regarding serious adverse reactions and potential safety hazards that can occur in proper use or misuse of the equipment. Failure to observe the information or procedures presented in a Warning may result in injury to the patient and / or surgical staff.

Precaution: A precaution contains instructions for any special care to be exercised by the practitioners for the safe and effective use of the equipment. Failure to observe the information or procedures presented in a Precaution may result in damage to equipment.

NOTE: A note is added to provide additional focused information. The information has no critical effect on the patient, user or equipment.

1.7 Environmental Directive

It's the facility's responsibility to dispose of waste equipment by handing it over to a designated collection point for the recycling of Waste Electrical and Electronic Equipment. The separate collection and recycling of your waste equipment at the time of disposal will help conserve natural resources and ensure that it is recycled in a manner that protects human health and the environment. For more information on the disposal of waste electrical equipment at the end of its useful life for recycling, please contact CONMED.

1.8 Warranty and Repair

The following capital components are warranted for a period of two (2) years from the date of delivery. This policy covers the cost of repair or replacement of the system, labor, parts, loaner equipment (if requested) and shipping costs to your facility. ConMed only accepts products or devices that are cleaned and disinfected. Warranty repairs will not extend the warranty period under this policy.

- D4001 Procedure Cart

The warranty becomes void and the manufacturer is not liable for direct or resulting damage if:

- Negligence, such as, abuse or misuse of the systems
- Non-authorized persons perform repairs, adjustments or alterations to the device or accessories.

Note: Please refer to CONMED's Terms and Conditions for all warranty provisions that apply. For more information, please contact:

Tech Services Hotline
1-866-4CONMED (1-866-426-6633)
TechProductSupport@conmed.com

Responsibility of the Manufacturer

CONMED is responsible for the safety, reliability, and performance of the D4001 Procedure Cart under the following circumstance:

- The user has followed the Installation and Setup Procedures in this Instruction Manual.
- Only persons authorized by CONMED performed assembly operation, readjustments, modifications, or repairs.
- The electrical installation of the relevant room complies with local codes and regulatory requirements, such as IEC, UL, Intertek, etc.
- Equipment use is in accordance with the CONMED instructions for use.

Returning the D4001 Procedure Cart for Service

If you need technical assistance regarding the use or application of the D4001 Procedure Cart, or you encounter a problem that requires servicing or repair, contact CONMED Customer Service at 1-866-426-6633 or your CONMED Sales Representative. Outside the U.S. contact your local CONMED Representative or your local CONMED Service Center (refer to www.CONMED.com for the proper contact information).

Report any complaints to CONMED Customer Service.

Products returned for repair must have an authorized Service Request (SR) number prominently displayed on the box and included on all paperwork. Refer to this number if making inquiries about the repair status. Please call CONMED Customer Service and provide the following information to obtain an SR number prior to returning any product for repair:

- Product Number
- Serial/Lot Number
- Reason for Return
- Original Invoice Number
- Date of Purchase
- Detailed description of the problem.

1.9 Hotline / Technical Support

Contact CONMED Customer Service Department or the CONMED Sales Representative if a technical assistance is required regarding the use or application of this product, or the product requires servicing or repair. Outside the U.S. contact the local CONMED Representative.

Attn.: Customer Service Dept.
11311 Concept Boulevard Largo,
Florida 33773-4908 USA

Customer Service

(within U.S.)

Phone: 1-866-426-6633

FAX: (727) 399-5256

(outside U.S.)

Phone: 1-866-426-6633

FAX: + 1 (727)-397-4540

2.0 Introducing the D4001 Procedure Cart

Catalog Number	Description
D4001	Procedure Cart

2.1 Device Description

The D4001 Procedure Cart is used to house arthroscopic equipment for single or multiple surgical specialties.

The cart is used to move/store arthroscopic equipment and other accessories required during a procedure. It is comprised of shelves, storage drawer, IV pole, and cord wraps. The cart is equipped with three locking castors and 1 directional caster to prevent unwanted movement of the cart, and a handle for moving the cart. Refer to the D4001 Procedure Cart Operator's Manual (P000030963) for information related to installation, compatibility, and safety notices.

2.2 General Safety

WARNINGS

1. This equipment is designed for use by medical professionals completely familiar with the required techniques and instructions for use of the equipment. Read and follow all warning and precautionary notices and instructions marked on the product and included in this manual.
2. Use care when moving the cart.
3. Secure all equipment using included hardware.
4. Remove saline bags from hanger when moving cart on inclined surfaces, over bumps, power cords, elevator doorways or thresholds.
5. Ensure the wheel locks are released and path is clear of obstacles.
6. When moving the cart, push, don't pull. Use handle only to move cart.
7. Lock the front wheel brakes when not in transport to prevent unexpected cart movement.
8. Use caution when moving cart on inclined surfaces, over bumps, power cords, elevator doorways or thresholds.
9. Lock wheel brakes when not in transport. Release wheel brakes prior to transportation.
10. When installing the cart, tighten and securely connect each cart component (shelves, drawer, etc.) before placing equipment or accessories to avoid risk of injury.
11. Only one console may be stacked on any shelf at all times.
12. Ensure IV pole is secure in bracket prior to use.
13. Maximum rated loads.
 - Saline bag hanger load not to exceed 30 lbs (13.6 kg); 6 lbs (2.72 kg). per hook
 - Shelf load not to exceed 50 lbs (22.7 kg)
 - Drawer load not to exceed 10 lbs (4.5 kg)
 - Total load not to exceed 300 lbs (136 kg).
14. No modification of this equipment is allowed.
15. D4001 Cart is recommended to have preventative maintenance done periodically. Failure to do routine preventative maintenance may affect device performance.

PRECAUTIONS

1. Inspect the cart prior to each use to ensure proper function.
2. Do not use a damaged cart.
3. Check the wheels daily for possible buildup of materials gathered from the operating room floor.
4. Remove equipment from shelves prior to repositioning the shelves.
5. Distribute load evenly from bottom to top of cart for optimal stability.
6. Prescription Only - Federal law restricts this device to sale by or on the order of a physician
7. The D4001 Procedure Cart is not intended to be used with intravenous infusion equipment.
8. Do not directly spray the cart surfaces with cleaning agent. Fluid can cause corrosion.

2.3 Unpacking Cart Components and Accessories

Carefully unpack and inspect all components shipped with the D4001 Procedure Cart. If any parts are missing or damaged, contact your representative.

Table 1 – Procedure Cart Components

D4001 Procedure Cart Components	Quantity
4 mm Hex Key	1
Cart Instructions for Use	1
Adjustable Shelves	3
Drawer	1
Cord Wrap	1 Set
Cable Management Clips	2
I.V. Pole with Secondary IV Bag Hanger	1
D4001 Procedure Cart	1


Figure 1 – Assembled Procedure Cart

3.0 General Assembly and Use

3.1 Locking the Wheels

1. To lock the front wheel(s), place your foot on the wheel locking lever and push down firmly.
2. To unlock the front wheel(s), place your foot under the wheel locking lever and lift up.

3.2 Installing the IV Pole

1. Lock the rear wheels of the cart.
2. Insert the IV Pole into the mounting post of the cart.
3. Lift the IV pole up slightly to screw the rubber foot onto bottom of IV pole.


Figure 2 – IV Pole Installation and Removal

3.3 Removing the IV Pole

1. Lock the rear wheels of the cart and remove all fluids from the saline bag hangers prior to removing the IV pole.
2. Unscrew the rubber foot from bottom of IV pole.
3. Lift the pole up and out to free it from the mounting post in the cart.

3.4 Adjusting Shelf

1. The shelves of the D4001 Procedure Cart are adjustable in 1" (25.4 mm) increments.
2. To adjust the position of a shelf, use the supplied 4mm hex key to remove the 4 supporting screws.
3. Move the shelf to the desired position.
4. Insert the 4 supporting screws and tighten.
5. To install the accessory shelf, move shelf to desired position, insert the 4 supporting screws with lock washers attached and tighten using the 4 mm hex key.


Figure 3 – Shelf Adjustment

3.5 Assembling Drawer

1. Install drawer to the bottom of shelf using four (4) 6-32 x 1/4" screws and tighten using the 4 mm hex key.
2. Push in drawer.


Figure 4 – Drawer Assembly

3.6 Installing Console-to-Shelf

NOTE: Standoffs can be removed from the shelf if they interfere with placement of equipment.

1. Place each piece of equipment on an individual shelf.
2. Use the shelf standoffs to position and attach the equipment to the shelf, where applicable.
3. Remove thumb screws from standoffs.
4. Place the holes in the bottom edge of the back of the equipment over the standoffs.
5. Replace the thumb screws.


Figure 5 – Console-To-Shelf Installation

3.7 Installing Handle with Cord Hooks

1. Slide metal brackets on handle into back extrusion of the Cart and firmly secure at the desired height.
 - a. Ensure handle is level.
 - b. Tighten all 4 screws on the cart handle with the provided 4mm Allen wrench. See Figure 6 below.


Figure 6 – Handle Assembly

2. Insert extrusion caps into top of extrusion rails by pressing into place. See Figure 7 below.


Figure 7 – Extrusion Cap Placement

4.0 Maintenance

4.1 Cleaning

For routine cleaning, wipe the surfaces of the cart with a damp cloth. Do not directly spray the cart surfaces. Fluid can cause corrosion.

The following cleaning agents can be used to clean the exterior of the cart:

- Mild detergent such as Windex® or Formula 409®
- 70% isopropyl alcohol
- 5% sodium hypochlorite (Clorox®)
- CaviWipes

4.2 Preventative Maintenance





1. Periodically check all hardware for tightness and security.
2. Tighten fasteners as required.
3. Check wheels for damage.
4. Make sure that IV Pole is secure in bracket and raises and lowers properly.

5.0 Technical Specifications

Table 2 – Technical Specification

Part	Specification
Cart Classification	Type I
Shelves	15.2" W x 18.2" D (38.6 cm W x 46.2 cm D)
Wheels	5" (125 mm) Single Wheel Swivel Wheels; 3 with Total Lock Brake; 1 Direction Lock
Handles	1.25" (31.8 mm) Diameter; Aluminum
Paint	Carbide Black textured hybrid & cornflower gray semi-gloss polyester
Drawer	11.75" W x 17.63" D x 5.00" H (28.8 cm W x 44.7 cm D x 12.7 cm H)
Cable Management Clips	Nylon 6/6 Arrowhead Mount, 0.44" x 1.00" wire capacity
Cord Wraps	PC/PET E500 material, swivel, 1.625" x 1.875" wire capacity
Dimension without IV Pole	43.4" H x 20.7" W x 29.5" D (110.5 cm H x 52.6 cm W x 75.0 cm D)
IV Pole	Detachable; 4 IV Bag Hangers; Up to 78" (198 cm) Height
Secondary IV Bag Hanger	Adjustable; 1 IV Bag Hanger
Compatible Devices within a Patient Environment	CV4000 24K D4000 D4000A 10K AES-1 (When Plugged into Separate Wall Outlet)

6.0 Environmental Limits

Environmental Conditions	Operating	Storage and Transport
Temperature:		
Relative Humidity:		

7.0 Symbols used on Labeling

Table 3 – Symbols used on Labeling

	Federal law (USA) restricts this device to sale by, or on order of a licensed healthcare practitioner or physician		Catalog Number
	Refer to Instruction Manual/Booklet (for critical safety instruction). The color of the logo is blue.		Manufacturer of the device
	Authorized Representative in the European Community		Consult Instructions for Use
	Caution		Serial Number
	Non Sterile		Date of Manufacture
	Listing Mark		European Union Conformity (CE Mark Conformity)
	Keep Dry		Identifies the item is considered as a Medical Device
	Maximum Stacking Height	Max Shelf Load 22.7 kg (50 lbs.)	Maximum Shelf Load
	Max Drawer Weight 4.5 kg (10 lbs.)	Max Cart Load 136 kg (300 lbs.)	Maximum Cart Load
		Max Hook Load 2.72 kg (6 lbs.) Max IV Pole Load 13.6 kg (30 lbs.)	Maximum Hook Load Maximum IV Pole Load

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©2022 CONMED Corporation. All Rights Reserved.
MDSS GmbH