
Nellcor™

Adult Respiratory Sensor SpO₂, RR


en	Instructions For Use
fr	Mode d'emploi
de	Gebrauchsanleitung
nl	Gebruiksaanwijzing
it	Istruzioni per l'uso
es	Instrucciones de uso
sv	Bruksanvisning
da	Brugervejledning
no	Bruksanvisning
fi	Käyttöohjeet
pt	Instruções de uso
ru	Инструкция по применению
zh	使用说明
pl	Instrukcja stosowania
cs	Návod k použití
sk	Návod na použitie
sl	Navodila za uporabo
hu	Használati utasítás
el	Οδηγίες χρήσης
tr	Kullanma Talimatı
bg	Указания за употреба
ro	Instrucțiuni de utilizare
et	Kasutusjuhend
ja	使用方法
lt	Naudojimo instrukcija
lv	Lietošanas instrukcijas
sr	Uputstva za upotrebu
hr	Upute za uporabu

REF 10068119


OxiMax™
Technology


1

2

3

Adult Respiratory Sensor

SpO₂, RR

Identification of a substance that is contained or present within the product or packaging

Do not use if package is opened or damaged

Not made with natural rubber latex

Not made with DEHP

Identification of a substance that is not contained or present within the product or packaging

Do Not Use with MRI

IP22

Protected against access to hazardous parts with a finger and protected against vertically falling water drops when sensor tilted up to 15°

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in bio-incompatibility, infection or product failure risks to the patient.

This device is not made with natural rubber latex or DEHP.

Directions for Use

Indications/Contraindications

The Nellcor™ Adult Respiratory Sensor, when used in conjunction with a Nellcor pulse oximeter and the Nellcor Respiration Rate Software, is indicated for single patient use when continuous noninvasive arterial oxygen saturation, pulse rate and respiration rate monitoring are required for adult patients weighing more than 30 kg.

The Nellcor Adult Respiratory Sensor is contraindicated for use on patients who exhibit allergic reactions to the adhesive tape.

Instructions for Use

1. Remove plastic backing from the sensor and locate transparent windows (a) on the adhesive side. Windows cover the sensor's optical components. **(1)**

An index finger is the preferred sensor location. Alternatively, apply the sensor to a small thumb, small finger, or great toe.

Note: When selecting a sensor site priority should be given to an extremity free of an arterial catheter, blood pressure cuff, or intravascular infusion line.

2. Orient the sensor so the dashed line in the middle of the sensor is centered on the tip of the digit. Wrap adhesive flaps on non-cable end around the digit. Note that the cable must be positioned on the top of the hand. **(2)**
3. Fold cable end over top of digit so that windows are directly opposite each other. Wrap adhesive securely around sides of digit. **(3)**
4. Plug the sensor into the oximeter and verify proper operation as described in the oximeter operator's manual.

Note: If the sensor does not track the pulse reliably, it may be incorrectly positioned, or the sensor site may be too thick, thin, deeply pigmented, or otherwise deeply colored (for example, as a result of externally applied coloring such as nail polish, dye, or pigmented cream) to permit appropriate light transmission. If any of these situations occurs, reposition the sensor or choose an alternate Nellcor sensor for use on a different site.

WARNINGS

1. **Do not use the Nellcor Adult Respiratory Sensor or other oximetry sensors during MRI scanning. Conducted current may cause burns. Also, the sensor may affect the MRI image, and the MRI unit may affect the accuracy of oximetry measurements.**

CAUTIONS

1. Federal law restricts this device to sale by or on the order of a physician.
2. In the event of damage to the sterile packaging, do NOT re-sterilize. Follow local governing ordinances and recycling instructions regarding disposal or recycling of sensors.
3. Failure to apply the sensor properly may cause incorrect measurements.
4. While the sensor is designed to reduce the effects of ambient light, excessive light may cause inaccurate measurements. In such cases, cover the sensor with opaque material.
5. Circulation distal to the sensor site should be checked routinely. The site must be inspected every 8 hours to ensure adhesion, application pressure, skin integrity, and correct optical alignment. If skin integrity changes, move the sensor to another site. If the sensor is misapplied with excessive pressure, a pressure injury can occur.
6. Intravascular dyes or externally applied coloring such as nail polish, dye, or pigmented cream may lead to inaccurate measurements.
7. Excessive motion may compromise performance. In such cases, try to keep the patient still, or change the sensor site to one with less motion.
8. Do not immerse in water or cleaning solutions. Do not resterilize. Immersion or resterilization could damage the sensor.
9. If the sensor is wrapped too tightly or supplemental tape is applied, venous pulsations may lead to inaccurate saturation measurements.
10. Do not alter or modify the sensor. Alterations or modifications may affect performance or accuracy.
11. For additional warnings, cautions or contraindications when using this

sensor with Nellcor-compatible instruments, refer to the instrument operator's manual or contact the manufacturer of the instrument.

Note: High oxygen levels may predispose a premature infant to develop retinopathy. Therefore, the upper alarm limit for oxygen saturation must be carefully selected in accordance with accepted clinical standards and considering the accuracy range of the oximeter being used.

Accuracy Specifications

For the accuracy specification range when used with Nellcor monitors, refer to information provided with the monitor, or contact Covidien's Technical Services Department. Outside the U.S., contact your local Covidien representative.

For the accuracy specification range of this sensor when used with Nellcor-compatible instruments, refer to the instrument operator's manual or contact the instrument manufacturer.

Use this sensor only with Nellcor OxiMax™ instruments and instruments containing Nellcor oximetry, or with instruments licensed to use Nellcor OxiMax sensors (Nellcor-compatible instruments). This sensor integrates Nellcor OxiMax technology into its design. When connected to an OxiMax-enabled instrument, this sensor uses OxiMax technology to provide additional advanced sensor performance features. Consult individual manufacturers for features and compatibility of particular instruments and sensor models.

Each Nellcor OxiMax-compatible instrument manufacturer is responsible for determining optimal compatibility conditions and settings for its instruments to provide safe and effective use of each Nellcor OxiMax sensor. This includes specifications and/or warnings, cautions, or contraindications. Refer to each instrument operator's manual or consult manufacturer for complete instructions, warnings, cautions, or contraindications on the use of this sensor with their Nellcor OxiMax compatible instrument.

Additional Copies of Instructions

Additional copies of these Instructions are available at no charge by calling Covidien or its authorized distributors. Also, permission is hereby granted under Covidien copyrights to purchasers of products obtained from Covidien or its authorized distributors to make additional copies of these instructions for use by such purchasers.

Any other use of this sensor is not authorized by Covidien under any patents.

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COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.
EC REP Covidien Ireland Limited,