Clarimat 300

Operator and service manual

- Operator and Service Manual
- Manuel de l'Opérateur et de Service
- Operador y Manual de Servicio
- Manuale dell' Operatore e Servizio

X-ray Systems • Digital Systems • Intraoral Camera Systems

GENDEX

CLARIMAT 300
ENGLISH

Clarimat 300

Contents

<i>1. Introduction and table of symbols</i>	<i>2</i>
<i>2. External components of the machine</i>	<i>3</i>
<i>3. Installation</i>	<i>5</i>
<i>4. Preparing to use the machine</i>	<i>7</i>
<i>5. Operating the CLARIMAT 300</i>	<i>9</i>
<i>6. Routine maintenance</i>	<i>12</i>
<i>7. Fault finding</i>	<i>13</i>
<i>8. Technical data</i>	<i>16</i>

1. Introduction and table of symbols

This manual is intended to assist the user in the safe and efficient operation of the equipment described. The equipment must be used in accordance with the procedures contained in the manual and must not be used for purposes other than those which are described herein. The equipment should only be used by persons having recognised qualifications and, if relevant, having adequate training on that particular equipment. It is the responsibility of the user to ensure that existing legal regulations regarding installation of the equipment are observed. Moreover, the legal regulations regarding operation of the equipment must be observed. Incorrect operation, or failure of the user to maintain the equipment in accordance with the maintenance schedule, relieves the manufacturer or his agent from all responsibilities for consequent non-compliance, damage, injury, defect and/or other malfunction.

Symbol	Description
	Alternating current
	Off (power: disconnection from the mains)
	On (power: connection to the mains)
	Type b equipment
	Film transport switch

Note (in this Operator's Manual)


Warning (in this Operator's Manual)


The CE symbol ensures that the product here in specified, and manufactured from 1996 onwards, needs the provisions of the European Council Directive 89/336/EEC relating to electromagnetic compatibility.

2. External components of the machine

3. View showing internal tanks and transport system

3. Installation

Siting your machine for daylight operation

Although the Gendex is designed for use in daylight, care should be taken to avoid sources of intense light. For example do not sit directly under a window, fluorescent tube or floodlight.

Siting your machine for darkroom operation

It is important to ensure that the darkroom has good ventilation and is adequate size. It should be completely light tight and of sufficient size to allow for loading and unloading of cassettes.

A site near a sink is ideal as this is a useful position invaluable when changing chemicals and cleaning machine. avoid carpeted areas. If your machine has not been supplied with the designed stand then ensure the processor is sited on a firm and a flat surface of correct height (and size, see diagram, page 5. It is important that when the machine is filled with solution that the stand does not rock or shake as this could cause spillage. Never move the machine when chemicals are in the tanks.

Services electrical

The machine should be positioned within approximately 3ft of the power supply. a fused plug should be fitted with a 7a fuse (uk only). In most instances the machine will have been supplied with a moulded plug mains lead.

Plumbing

Study the plumbing diagram, page 5.

Waste

As you will note from the diagram a standpipe waste should be installed similar to that required on a domestic washing machine. The wastes should not be less than 1 1/2" diameter and should rise between 22" and 24" above floor level immediately behind or within 2ft of the left hand side of the machine. The outlet hoses for water and chemical wastes are put into the stand pipe as shown.

Water supply

The machine requires only cold water for film rinsing. this cold water supply should be taken from the buildings water tanks rather than directly off the mains. The water supply should be brought within 3 ft of the stand and fitted with a ball-o-flx isolating valve and terminated with a bibcock tap. this can be seen clearly illustrated in diagram, page 5. the tap should be situated in a convenient position so that it can be turned on and off each day with ease.

Connecting the machine

The machine can then be easily connected with the flexible hose supplied by attaching the hose to the water inlet on the back of the machine and its other end to the bibcock tap. see diagram, page 5.

Water supply must be flushed out before connecting to the machine to remove any dirt/dust, etc.. failure to do so could result in a damaged solenoid valve.

The flexible connecting pipe has a rubber sealing washer one end and a filter at the bibcock tap end ensure that these are in place, see diagram, page 5. When screwing on the fitting ensure they are firmly tightened by hand. Do not overtighten with mole grips or pliers as this could cause damage.

Water flow

The machine is fitted with a solenoid valve which cuts off water supply when the machine is turned off.

The machine requires a water flow of approximately 1 litre per minute. The flow rate is factory set, water flow noise can be reduced by adjustment of the ball-o-fix valve. To set the water flow via this valve proceed as follows: remove processor lid and dryer module from the machine. Locate solenoid at bottom of dryer and using a screwdriver turn screw 1/2 clockwise to vertical position. (see diagram below).

Close ball-o-fix valve by turning it with a screwdriver to the off position. Fully open bibcock tap (on-off tap), then carefully open ball-o-fix valve while looking into the water tank of the machine. Water will be seen to flow out of the water inlet pipe.

When required flow rate is achieved turn solenoid screw back to its original horizontal position (anti-clockwise).

The machine must not be switched on when setting the water flow.

4. Preparing to use the machine

Preparing the machine for daylight operation attaching the daylight loader

See diagram, page 3. extra oral films are extremely light sensitive so it is important that the daylight loader is fitted correctly and sealed around the film entry slot and process switch.

To open the loader release catches and hinge lid backwards.

Proceed as follows:

- 1) Remove white protective strip from sealing material around film entry recess and process switch hole
- 2) Remove the four black knobs from the front of the Gendex and slide daylight loader over projecting studs. couple the loader to the main body of the machine by replacing the black knobs on the studs from the inside of the loader. This should be done with one hand starting at the top first, whilst the other hand should be firmly pressing the panel around the entry slot until all knobs are tight and in position.
- 3) Now check catches on daylight loader.

Practice using the daylight loader with a dummy film in order to familiarise yourself with the location of the film entry slot and film process switch. Also practice loading and unloading cassettes. When processing film do not remove hands from loading section until the film has fully entered the machine. Use this time to reload the cassette. Keep the daylight loader clear of old film wrappers to avoid these entering the machine and to maximise the space available.

Ensure the viewing cover is in place when loading extra-oral film, failure to do this will cause the film to fog.

Filling the machine with chemicals

Unlock and open the lid by means of the protruding grip, see diagram page 3, item 11. Bring it to the upright position where it will rest. The lid may be removed if this is more convenient.

Remove all four transport modules, starting with the dryer module as follows. Lift the module slightly, slide it to the left (to disengage from the drive dog in the side of the machine) and then lift it straight out. The position of each transport is clearly marked including the transport itself.

Before filling the tanks with chemicals check the drain tubes are in the proper place. When using Gendex ready to use chemicals, fill the tanks to the top of the upper level rib (3650 ml). When using Gendex one plus one chemicals, fill to the top of the lower rib (1825 ml) and top up with clean cold water to the top of the upper level rib.

Do not allow chemicals to splash. Cross contamination could occur, drastically shortening the life of either chemical. Any spillage should be immediately wiped clean.

Replacing the transport system

First replace the developer transport by lowering gently into the tank of developer solution. Too rapid immersion will cause excessive displacement of chemical leading to spillage. Engage the drive pin into drive dog and lower the back end of the module behind the module limit stop. Repeat this procedure for all stages of the machine, fixer, wash and dry see diagram, page 4.

Re-check the level of the liquids and top up to the overflow hole in the standpipe of each tank, finally, replace the lid, and lock.

With all modules in place it is important to check that they are engaged and running smoothly. (see page 9, starting of the machine).

This is best achieved by passing a clean up film through the machine, when all the modules should be seen to function smoothly. If any module fails to turn re-check the engagement of the drive pin in the drive dogs, by lifting slightly and re-engaging.

5. Operating the CLARIMAT 300

Control panel

Each day

Turn on water supply.

Check liquid level in coolant and top up if necessary if required.

For cleaning and test film through process as the machine can be initiated by pressing the process button.

Before processing film check that ready light has gone out.

After loading a series of film always to press for stop/start button before taking machine.

Turn off mains and water supply to the processor at night.

Never use the machine without the presence of liquid.

When changing chemicals

1. Drain each tank. In turn, drain together. Tanks are drained by removal of drain tubes.
2. Remove transport system which should be immersed and thoroughly cleaned using hot water (not boiling). Drums and rollers should be cleaned with the brush provided then returned to their respective tanks in the machine.
3. Fill developer and fixer tanks with clean water and run for approximately 5 minutes. Drain chemicals whilst run (chemicals foaming great care not to splash the chemicals on contamination of either developer or fixer solution could occur).
4. Always keep machine clean and dry.

DO NOT OPERATE CLARIMAT 300 IN PRESENCE OF FLAMMABLE GASES OR VAPOURS

DO NOT OPERATE CLARIMAT 300 IN PRESENCE OF FLAMMABLE GASES OR VAPOURS

DO NOT OPERATE CLARIMAT 300 IN PRESENCE OF FLAMMABLE GASES OR VAPOURS

WARNING

THIS EQUIPMENT IS NOT DESIGNED FOR USE IN PRESENCE OF FLAMMABLE GASES OR VAPOURS

When the mains switch is turned on the automatic temperature control will maintain the correct operating temperature for processing. Correct temperature is indicated when the red light is extinguished. The machine can be safely left in standby throughout the day ready for immediate use.

The time taken to achieve the correct temperature for processing depends on the room temperature. Having achieved the required temperature the machine is now ready to accept films (average warm up time is 10-15 mins).

The velopex is equipped with automatic standby mode. The machine is then initiated by pressing the process switch (8) on the front of the machine.

Starting the machine

The machine having been filled with chemicals is now ready for operation, proceed as follows:

Turn on water supply and mains

Turn processor on by pressing the switch on the control panel to the indicate "on" position. The "ready light" will illuminate if the chemicals have reached the correct temperature for processing.

The time taken for the indicator to extinguish will be dependant on the room temperature. From time to time this light will illuminate automatically for short intervals as the machine requires heat throughout the day.

Having come to the required temperature the machine is now ready to accept film. The machine can now be initiated by pressing the process switch (8) on the front of the machine.

When first turned on the machine will operate for approximately 8 minutes after which time it will automatically stop. This is part of the machine's automatic standby operation and is an economy measure. If at this point the "ready light" has extinguished the machine is ready to process films. The machine should be left in the standby mode all day. This enables you to process films instantly at any point throughout the day. The machine need only be turned off at night.

Processing films

Before processing, intra oral films should be wiped clean of all mouth contaminants. Open daylight loader and place in film to be processed, close and lock. put hands through the loading gloves, press the process switch, strip the wrapping from the film and post into film entry slot.

As the velopex processor has a continuous processing mechanism, as soon as the first film has fully entered the mechanism so the next film can be processed. Intra-oral film can be processed via the film guides across the width of the film entry slot for ease of patient identification.

Always remove old film wrappers from inside the loader.

Use velopex micro-dots to mount your films for professional filing of your patients records.

Processing extra oral films

Ensure viewers cover lid is correctly in place before attempting to load cassette or process film. Remember extra-oral films are extremely light sensitive. Feed films squarely into the film entry slot. Do not remove hands until the film has left the loading area. Use this opportunity to reload your cassette to minimise handling time.

Remember

Do not remove hands from loading section until the films have completely entered the machine.

It is recommended that the "process" switch be pressed again after the last film has been posted, this avoids the possibility of the machine going into standby while the films are still in the machine.

Film collector

All films collect in the film collector at the rear of the machine. The small films collected in sectionalised pockets which are aligned with the entry slots on the front of the machine. These pockets separate the films for organisational purposes.

Each days

Before use check level indicators. If illuminated, open machine and top up chemicals where appropriate, before processing run through "clean-up film" for clean up film use a spare intra oval film which may be kept for this purpose for one week. After this discard the old film and use a fresh one. This system helps to clean the transport system by removing trapped crystals and ensures that the transport mechanism is running properly. At the end of the day turn off water and mains switch.

6. Routine maintenance

Cleaning processor

A clean processor is the key to consistent good results and the following cleaning routine should be carried out at each chemical change or approximately once a fortnight, depending on the use of your machine and the number of x-rays passed through it and the sizes of those films.

When the Gendex is fitted with automatic replenishment, the change cycle will be substantially extended, with complete chemical changes being required only every 4-6 weeks, according to use.

Unlock the lid of the machine.

Open the lid by means of the protruding grip, bring it to the upright position where it will rest, remember that the transport modules will now be wet with chemicals, so take care to let them drain a little over their tanks before removing them to a sink.

- a) Remove all four transport modules, starting with the dryer module, as follows: lift the module slightly, slide it to the left and lift it straight out. Hold it above the tank allowing it to drain off before completely removing.

Processing chemicals should be handled with care, if spill they may cause staining or corrosion of surfaces, they can also cause skin and eye irritation, therefore wipe out any drips or spillage immediately. If chemicals are splashed onto the skin or eyes, flush immediately with copious amounts of water. More detailed advice is printed on the containers of developer and fixer.

- b) Particular attention should be paid to the transport modules, they should be cleaned one at a time and always returned to the tank from which they were removed eg developer to developer. Please remember that the fixer wash and dryer modules also require regular cleaning procedure use fresh water for each module to avoid cross-contamination.
- c) The best method of cleaning a module is to immerse it totally in a sink full of hot water (not boiling), using the brush provided, scrub around the gears and the roller ends. Rinse the belts thoroughly, using the special tool provided, turn the modules by hand in order to rotate the gears and belts to assist in thorough cleaning. Do not use detergents, just hot water (not boiling), do not use boiling water straight from a kettle as this will cause damage to the module. Check if the drain plugs are securely in place before filling the machine with fresh chemicals. Then replace modules in correct order as in illustration in diagram, page 4. Keep a note of the date when chemicals are changed. Create a chart so you have a record of this operation.

CLARIMAT 300 system cleaner

Improved cleaning can be achieved by using velopex system cleaner. Instructions for use are enclosed with the product. Do not begrudge time spent cleaning the machine as you will be amply rewarded by quality results and extended chemical and machine life.

Disinfection

Intra oral films should be wiped clean before processing to remove mouth contaminants. Having removed old film wrappers from the daylight loader, clean interior with cold sterilising solution and wipe dry after appropriate period.

7. Fault finding

The processor does not operate

- 1) Check that the mains supply is switched on at the socket and on the machine.
- 2) The machine is fitted with a safety switch, check that the lid has been lifted properly. If the lid is not correctly closed the safety switch will prevent the machine operating. Check by opening and closing again.
- 3) Check the fuses in the machine. These are to be found at the back of the machine (page, item 6).
- 4) Check the fuses in the mains plug
- 5) Liquid levels low, machine will not operate.

Film lost in the machine

- 1) Check that transport modules are correctly located in their drive dogs and running correctly. If not relocate.
- 2) Check that all transport modules have been put through the correct cleaning procedure. A dirty machine can cause loss of film, chemical crystals can gather on the moving parts, causing them to wear or break, and on the webbing, making it sticky enough for the film to slip.
- 3) Check the spring tensioning of the modules and ensure all springs are in place. Use the module tuning tool to rotate the belts to inspect that it is running correctly. Run a dummy film through the module by this method. Should you be unable to rectify the functioning of the module contact the manufacturer for further advice.
- 4) Should there be any damage to the belts such as tearing or twisting order a replacement module.
- 5) If films end up lying across the tops of the transport modules they may have been affected by static from the dryer module. Some surgery environments are more prone to this than other. This problem may easily be overcome by dunking the dryer module (after it has been cleaned in the normal way into a sink of cold water, in which has been dissolved half a cupful of softener. Lift the module out of the liquid, shake it to remove the water and replace it in the machine.
Fabric softener is available in most supermarkets for domestic washing machines and is quite suitable.
- 6) Post small films accurately and squarely in their respective lanes this will ensure that they arrive in their correct position on the collector. Incorrect positioning could cause the film to slide off the collector and be lost.

Film will not enter film entry slot

- A) Press film process switch as the machine may be in standby mode.
- B) Incorrect engagement of the developer transport units. Re-engage where necessary.
- C) Liquid levels low, machine will not operate

Films come out too dark

Test for light fog by feeding an unexposed film through the machine (either in the darkroom or via the daylight loading section). It should process as a transparent piece of film base and there should be no shadows or blackness on the film. If there is :

- A) Check that the top lid is firmly in place

- B) Check that the darkroom is light-tight or that the daylight loading section is firmly on, refer to the previous instructions in this manual. Also check very closely for other possible light leaks. Ensure that, when the daylight loader is used, the machine has not been sited in direct sunlight or under intense lighting conditions.
- C) Is the viewer cover on correctly? Extra-oral films are extremely light sensitive and will fog if the cover is not used.
- D) Did you inadvertently take your hands out of the hand entry gloves before the film had fed completely into the machine?
- E) Could light have got into the box of film? In both cases D and E above, the fog tends to be to one end of the film or in patches across the film.
- F) Check that the Gendax is sited on a flat and level surface. An uneven surface- perhaps caused by an electrical lead being trapped under the machine - could cause distortion of the cabinet and lead to light leaks, and erratic running of the machine.
- G) Check that there has been no chemical mix up leading to cross contamination.
- H) Check temperature of developer and fixer tanks. This is generally set at 27°C if there is any substantial difference this could be a reason for the dark film. Turn off heater switch and contact supplier or manufacturer.
- I) Is the developer temperature light going out? If not, turn off heater switch and contact supplier or manufacturer.
- J) Check at hand entry ports ensuring there is a good light seal around the wrists if not, contact the manufacturer.
- K) Check if film are not out of date, excessive heat can cause premature ageing of the film. Check expiry date on box. Keep films in a cool dry place.

Films come out too light

- A) Chemicals may need replacing.
- B) The chemicals may be contaminated.
- C) The level of chemicals in the tanks may be too low.
- D) The temperature of the chemicals may be too low. If the temperature indicator light does not go out, check with a thermometer. If the temperature is much less than 27°C, consult supplier or manufacturer.
- E) Check that films are correct for the type of intensifying screen used with the x-ray machine.
- F) Check films have been exposed correctly.

Films come out grubby or marked

- A) The water tank may not be cleaned.
- B) Check that the transport are being cleaned correctly.
- C) The transport modules may have been replaced in the wrong position, leading to cross contamination.
- D) Stray light may be getting into the machine, check for leaks. Watch out for sources of strong light such as sunlight, spot lights and strip light sited too close the machine. Check darkroom for stray light particularly doors/vents.
- E) Check films are not out of date, excessive heat can cause premature ageing of the film. Check expiry date on box. Keep films in a cool dry place.

Contamination

When changing chemicals, make sure that the tanks are drained down fully and rinsed out as fresh chemical could be spoiled by contamination leading to poor results.

Should your tanks fail to drain down fully when the drain tubes are removed, it will probably be because the waste outlet pipe is not lying flat but rising before entering the waste pipe. The tank will then take its level from this point. Correct this by ensuring that the waste exit pipe does not rise above bench level. If the unit should emit abnormal odour, become overheated in some parts or produce unusual noises, immediately switch the machine off at machine and mains and contact your supplier or manufacturer.

A clean, well maintained machine makes for good clear x-rays.

A dirty machine just makes problems.

Do not struggle with film processing problems. We, through experience, probably have the answer. If you are having problems, give us a ring or write, enclosing sample of the film and we will do our best to resolve them for you.

8. Technical data

Please refer to TECHNICAL SECTION of this manual for more detailed informations.

Power requirements: 220 - 240 V
 Frequency: 50 Hz
 Power consumption: 10 A max

Tank capacity:
 - Developer 4 l
 - Fixer 4 l
 - Water 1,4 l
 Water feed: 3/4" BSP

Film format:
 - Infraoral all available formats (from 2 x 3 cm. up to 4 x 5 cm)
 - Panoramic and all available formats (from 15 x 30 cm. up to 24 x 30 cm)
 - Cephalometric
 - Max film width 260 mm

Timing Characteristic:
 - Warming-up time 10-15 minutes approx.
 - Processing time 4 minutes
 (fixed and controlled) 2 minutes (Endo)
 - Film feed speed 0,83 cm per s (3,6 s per 15 x 30 cm or 24 x 30 cm)
 - Automatic stand-by after 8 minutes an automatic stand-by facility stops the processor (film processing is started via an operating switch)

Chemicals temperature: Developer 25°C
 (fixed and automatically controlled) Fixer 27°C
 Environment temperature: < 23°C

Weight: approx. 27 kg (with liquid), approx. 21 kg (empty)
 Dimensions:

	Clarimat 300 (with day-light loader)	Clarimat 300 (without day-light loader)
Width	500 mm	500 mm
Depth	600 mm	420 mm
Height	340 mm	340 mm

Clarimat 300 is designed and produced to meet the following standards:


- General (electrical/mechanical) safety: IEC 601-1 (1988)
- Electro-Magnetic interference: IEC 601-1-2 (1993)

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Each day

Run on water supply

Check that water is correct and feed up chemicals if required

For cleaning and test film through processor

The machine can be initiated by pressing the process switch

Before processing film check that water light has gone out

After feeding a meter of film always re-press the transport switch before working machine

Turn off mains and water supply to the processor at night

Never use the machine without the presence of fluid

When changing chemicals

1. Drain each tank. In turn, never together. Tanks are drained by removal of drain tubes
2. Remove a transport system which should be immersed and thoroughly cleaned using hot water (not boiling). Drain and repeat this on the second with the brush provided then return it further respective tanks in the machine
3. Fill developer and fixer tanks with cold water and run for approximately 5 minutes. Drain and refill with fresh chemicals (using great care not to splash the chemicals or contamination of other developer or fixer solutions could occur)
4. Always keep machine clean and dry

Always use chemicals and solutions in the correct order

Always use chemicals and solutions in the correct order

Always use chemicals and solutions in the correct order

WARNING:

THIS EQUIPMENT IS NOT SUITABLE FOR USE
IN PRESENCE OF FLAMMABLE GASES
OR VAPOURS

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Each day

Run on water supply

Check fluid level indicator and top up whenever it is low

Free cleaning and feed film through processor

The machine can be initiated by pressing the power button

Before processing film check that room light has gone out

After loading is over of film always require for finished which leaves leaving machine

Run off water and water supply to the processor at night

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When changing chemicals

1. Open each tank. In turn, immerse together. Tanks are drained by removal of clean tubes
2. Remove tankstand is then which should be immersed and thoroughly cleaned using hot water and brush. Clean and replace should be cleaned with the brush provided then returned to their respective lines in the machine
3. Fill developer tank four times with cold water and run for approximately 5 minutes. Drain and refill with fresh chemical being used. Do not let touch the chemicals or contamination of other chemicals or their solutions should result
4. Always keep machine clean and dry

WARNING: NEVER USE THE MACHINE WITHOUT THE PRESENCE OF FLUID

WARNING: NEVER USE THE MACHINE WITHOUT THE PRESENCE OF FLUID

WARNING: NEVER USE THE MACHINE WITHOUT THE PRESENCE OF FLUID

WARNING:

THIS EQUIPMENT IS NOT SUITABLE FOR USE IN PRESENCE OF FLAMMABLE GASES OR VAPORS

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Clarimat 300

Contents

<i>Information to assist in the removal of parts</i>	<i>3</i>
1) Fuse	
2) Access to internal components	
3) Motor mounting	
4) Replacing components on general drive strip	
5) Replacing drive dogs	
6) Fan and heater assembly	
7) To replace heater element	
8) Temperature and time controls	
9) Tanks installation	
10) Temperature sensor assembly	
11) Solenoid valve assembly	
<i>External layout of the machine</i>	<i>7</i>
<i>Internal layout of the machine</i>	<i>9</i>
<i>Motor board assembly</i>	<i>10</i>
<i>Transport module assembly</i>	<i>11</i>
<i>Wiring diagram</i>	<i>12</i>
<i>Pump assembly and tank installation</i>	<i>13</i>
<i>Tank installation</i>	<i>14</i>
<i>Replenisher pump mounting instructions</i>	<i>15</i>
<i>Replenisher pump assembly</i>	<i>16</i>
<i>Plumbing diagram</i>	<i>17</i>
<i>Component part numbers</i>	<i>19</i>

Information to assist in the removal of parts

NOTE: Always Switch off Mains Power and Remove Electricity Plug before beginning any work or inspection procedure.

1. Fuse

Refer to diagram, page 9. The power cord socket is located on the back of the machine. This contains a drawer section, which, when slid out, reveals both fuses on UK-Continental machines.

USA: no fuse drawer, but circuit breakers are fitted above socket.

2. Access to Internal Components

Refer to diagram, page 9. To access internal workings of the machine unscrew control panel retaining screw (item 15). The control panel may then be hinged upwards. The side panel may then be removed by releasing catch as illustrated and described in diagram.

3. Motor Mounting

Refer to diagrams, pages 9 and 10. Disconnect wires from control PCB (page 9, item 9) and terminal block (page 9, item 5). Release the motor by removing screws (page 10, item 15) when the motor may be withdrawn from gear strip and drive gear. To replace, reverse procedure.

4. Replacing Components on the Gear Strip

Refer to the diagram on page 10 taking special note of the gear strip assembly detail before stripping it down.

The assembly consists of three strips (Items 2, 3 and 4) and the motor board (Item 1) which is clamped between gear strips 2 and 3. It is essential to maintain this assembly order.

To replace the main drive gear (Item 7) on the motor shaft, do not dismantle the gear strip assembly. Simply remove the motor (Item 14) by referring to section 3. Withdraw the main drive gear upwards from the gear strip assembly and replace with the new gear. Finally refer to section 3 and replace the motor.

To replace the other gears (Items 6 and 5) follow the procedure as for replacing the drive dogs (Item 10) in section 5.

5. Replacing the Drive Dogs

Refer to the diagram on page 10 and remove the four clamping nuts (Item 12). This allows the removal of the motor support strip (Item 4) and the gear support strip (Item 3) along with the drive dog shaft (Item 13), the gears (Items 5 and 6), the drive dog springs (Item 9) and the thrust washers (Item 8). Now dismantle the assembly for cleaning. Be careful not to lose any of the components.

Wipe away the old grease from the springs, thrust washers and the drive dog shafts. Assemble the module drive gears (Item 5) onto the drive dog shafts (Item 13) feed the shafts through the gear support strip (Item 3) and apply a little silicon grease to the shafts before fitting the thrust washers (Item 8) and the drive dog springs (Item 9). The grease will hold the springs in place during reassembly. Fit the idler gears (Item 6) onto their spigots on the gear support strip (Item 3) and assemble the motor support strip (Item 4) into place.

Remove the old drive dogs (Item 10) from the motor side of the motor board (Item 1) and wipe clean the holes in the drive dog cover strip (Item 2). Apply a little silicon grease to the outside diameter of the new drive dogs (Item 10) and at their "D" entry holes and insert them into the drive dog cover plate (Item 2).

Offer up the above assembly to its position on the motor board (Item 1) taking care that the motor mounting holes are at the bottom of the motor support strip (Item 4). Starting at one end, align and centre the drive dog shaft (Item 13) into the drive dog (Item 10). Hold the assembly in position and fit the end clamping nut (Item 12) loosely to its gear strip assembly screw (Item 11).

Work along the other three drive dog shafts (Item 13) aligning and entering them into their drive dogs (Item 10) and fitting the clamping nut (Item 12) loosely to each gear strip assembly screw (Item 11) as you go. Finally tighten all four clamping nuts (Item 12). Check that all the gears turn freely and the drive dogs return freely to their outer position after being compressed. Slip the main drive gear (Item 7) into place in the centre of the gear train and refit the motor as described in section 3.

6. Fan and Heater Assembly

Refer to diagrams, pages 9 and 10. To remove fan unit (page 9, item 4) disconnect wires from terminal strip and PCB (page 9, items 5, 6, 9) and dryer terminal block (page 9, item 6). Remove grille (page 10, item 19) retained by four screws (page 10, item 20). Remove heater assembly. Fan unit may then be released by removing two screws (page 10, item 21) located at the back of the dryer element housing (page 10 item 18). To replace reverse procedure.

7. To Replace Heater Element

Refer to diagram, page 10. To replace heater element (item 17) disconnect wires from terminal block (item 24), remove grille (item 19) retained by four screws (item 20). Slide out heater element, withdrawing wires through the holes (item 22) at the back of the element housing (item 10). To replace, reverse the above procedure.

8. Temperature and Replenisher

Temperature and replenisher are located on the control PCB (page 9, item 9). These are pre-set at time of manufacture and are not customer adjustable. Normal factory settings are: Developer 77 degrees (F) (25 degrees C); Fixer 81 degrees (F) (27 degrees C) and 2 minutes (replenisher). If special temperature settings are found to be necessary a qualified engineer can make the adjustments on site after referring to the PCB Operating Instructions I64.

9. Tank Removal and Installation

Refer to diagrams, pages 13 and 14. To replace tank, drain down fully all three sections. Remove pump and manifold assemblies from both developer and fixer sections by undoing pump coupling nuts (page 13, item 10).

Beneath the machine, locate the 'tank retaining screw' (item 37 page 7) found on the base of the machine adjacent to the right hand machine foot and next to the motor board locating slot. This is a 'dome-headed' phillips screw (not flat-headed or counter-sunk like the screws retaining the foot tray). Failure to remove this tank retaining screw will cause serious damage to the tank.

Release lower flexible drain tubes (page 14, item 6) located under foot of machine (page 14, item 19) which is retained by two (counter-sunk) screws (page 14, item 20). Unscrew connectors (page 14, item 10) from locking nuts (page 14, item 9). Connectors can then be pulled clear of back panel, allowing the removal of all the connecting pipework including the water waste and overflow drain pipes (page 14, items 4 and 5) and solenoid assembly (page 141, item 7). Be sure to disconnect Temperature Sensor leads from control pcb. The tank may then be lifted clear of the machine. To install simply reverse the above procedure. The water connector has the flow controller inserted and should be replaced in the same position on left rear of machine (water inlet, item 10(c)).

After installation check thoroughly for leaks by filling the tanks with water and inspecting. Only when satisfied there are no leaks fill with chemicals.

10. Temperature Sensor Assemblies

Before starting to replace the Developer or Filter temperature sensor, be sure to drain the relevant tank of all liquid.

Refer to diagram, page 13. To replace the temperature sensor assembly (Item 7) disconnect sensor leads from the control PCB (page 9, Item 9) and unscrew the sensor assembly nut. Pass the sensor leads through the nut, one at a time, and remove the sensor from the tank. To replace, reverse this procedure.

11. Solenoid Valve Assembly (refer to diagram, page 14)

The solenoid valve assembly consists of components 7, 8 and 15. The solenoid connector (Item 8) is screwed into the main valve body (Item 7) using PTFE tape to effect water sealing. The male stud coupling (Item 15) is screwed into the other end of the main valve body. The stud coupling has its own special seal and does not require the use of PTFE tape.

When installing the solenoid ensure that the valve override screw faces out into the cabinet of the machine for access.

Current Models are fitted with an onboard Flow Controller situated within Component 10 (c). This device controls the water flow to a rate of 1 litre per minute.


12. Dryer Element Cut-Out Re-set Operation (see page 10, Item 17)


Dryer Elements have a small re-set button on the rear of the small black switch inside the Element. If at any time the fan should slow or stop, the Dryer Element will switch off and stay off until re-set button has been operated.

Refer to the PCB operating instructions at rear of manual for details of a further over-temperature cutout now fitted to the PCB.

External componets of the machine

Internal layout of the machine

Access to internal workings of the machine

Release control panel retaining screw 15, then control panel may be hinged back. To remove side panel release one side of the panel at a time by pushing out front and back panel to disengage the bolts while lifting the side panel upwards as illustrated in diagram.

Fuse

A main lead socket 7 is located on the back of the machine. This contains a draw section which when slid out reveals both the fuses.

Temperature and time controls

Temperature and time controls are located on the control PC, item 9. Normal settings are 27°C and 8 mins.




Diagram of side panel lock

- ➊ Motor compartment
- ➋ Motor
- ➌ Gear drive strip
- ➍ Dryer fan
- ➎ Terminal strip
- ➏ Dryer terminal block
- ➐ Main lead socket
- ➑ Safety switch
- ➒ Control PC
- ➓ Main switch
- ➔ -
- ➕ Process switch
- ➖ Circulation pump developer
- ➗ Circulation pump filter
- ➘ Control panel retaining screw
- ➙ Motorboard brace

Motor board assembly


Detail incorporating • Motor mounting • General drive strip • Drive dog and spring assembly • Fan and heater assembly • Terminal blocks

- ① Motor board
- ② Drive dog cover strip A
- ③ Gear support strip B
- ④ Motor support strip C
- ⑤ Module drive gears
- ⑥ Idler gears
- ⑦ Main drive gears
- ⑧ Thrust washer
- ⑨ Drive dog spring
- ⑩ Drive dog
- ⑪ Screw - gear strip
- ⑫ Nut - gear strip
- ⑬ Drive dog shaft
- ⑭ Motor & gearbox
- ⑮ Motor mounting screws
- ⑯ Fan unit
- ⑰ Dryer elements
- ⑱ Dryer element housing
- ⑲ Dryer grille
- ⑳ Screws - dryer grille
- ㉑ Screws - fan unit
- ㉒ Holes for heater wires
- ㉓ Terminal strip
- ㉔ Terminal block


Transport module assembly

Incorporating • Gear Assembly • Positional layout of Y bars and belt


- ① Main drive gear
- ② D shaped center hole gears (Total 4)
- ③ Small idler round hole gear (Total 2)
- ④ Large idler gear
- ⑤ Retaining screw
- ⑥ Gear cover plate

EXTRA-X M & N CONNECTIONS


Pump assembly and tank installation

incorporating pump assembly • Temperature sensor assembly
 • Liquid level sensor

Tank installation

- Detail incorporating • Solenoid valve assembly • Overflow and drain assemblies
- Flexible tank drain down couplings

Replenisher pump mounting instructions

Refer to diagram Replenisher pump assembly:

- 1) With pump motor (item 1) inside, and inner rotor housing (item 2) outside the machine, screw together using 3 M5 screws (item 9).
- 2) Put a washer (item 4) onto the motor shaft.
- 3) With the open end of the spring clip facing outwards, push a rotor Assy. (item 3) onto the motor shaft along the longest flat of the shaft.
- 4) With the rotor assembly in a vertical position, push a hose assembly (item 8) leaving approx. 65mm (2 1/2") hanging outside, into the RH slot of the rotor housing. Feed the remainder of the hose assembly into the rotor housing by turning the rotor assembly, by hand, in an anti clockwise direction. Finally, press the hose assembly into the LH slot of the rotor housing where it will be held captive. The two ends should now be of equal length, if not re-position.
- 5) Place a cover plate (item 5) in position, locating on bosses on the inner rotor housing (item 2).
- 6) Place outer rotor housing (item 6) in position, locating on bosses on the inner rotor housing (item 2).
- 7) Put the remaining washer (item 4) onto the motor shaft.
- 8) With the open end of the spring clip facing outwards push the remaining rotor assembly (item 3) onto the motor shaft along the shortest flat of the shaft ie ensure that the two rotor assemblies are at right angles to one another.
- 9) Fit the remaining hose assembly (item 8) following instructions at note 4
- 10) Place cover plate (with brass bush) (item 7) in position
- 11) Finally screw completed assembly together using 4 no. 8 screws (item 10)

Hose connections:

- 1) Push the end of the red hose assembly onto the shortest replenisher tube protruding through the back panel of the machine, and fix with the 2 remaining hose clips provided. It will be found easier to put the lower of the 2 clips on first.
- 2) Push the end of the black hose assembly onto the longest replenisher tube protruding through the back panel of the machine, and fix with the two remaining hose clips provided. Again it will be found easier to put the lower of the 2 clips on first.

Electrical connections:

Refer to the wiring diagram.


- 1) Connect earth wire (green/yellow) on pump motor to terminal block (item 23)
- 2) Connect live wire (violet) on pump motor to pc (page 14, item 9) position RE.
- 3) Connect neutral wire (blue) on pump motor in terminal block (item 23)

Finally


- 1) Replace all panels.
- 2) Position machine to its final working position, re connect water and drain tubes.
- 3) Remove the fixer bottle cap and replace with the dip tube with the black bung and cap.

Replenisher pump assembly


- 1 Pump motor
- 2 Rotor housing (inner)
- 3 Rotor assembly (total 2)
- 4 Washer (total 2)
- 5 Cover plate
- 6 Rotor housing (outer)
- 7 Cover plate (with brass bush)
- 8 Hose assembly (total 2)
- 9 M5 screw (total 3)
- 10 No. 8 unc screw (total 4)
- 11 Back panel



Plumbing diagram

- 1 Water inlet
- 2 Water outlet
- 3 Water overflow
- 4 To 2nd developer waste
- 5 To drain waste for disposal of chemicals
- 6 All tubes vented safely
- 7 Cold water to developer
- 8 Hot water to developer
- 9 To 2nd developer waste
- 10 To 2nd developer waste
- 11 To 2nd developer waste
- 12 To 2nd developer waste
- 13 To 2nd developer waste
- 14 To 2nd developer waste
- 15 To 2nd developer waste
- 16 To 2nd developer waste
- 17 To 2nd developer waste
- 18 To 2nd developer waste
- 19 To 2nd developer waste
- 20 To 2nd developer waste
- 21 To 2nd developer waste
- 22 To 2nd developer waste
- 23 To 2nd developer waste
- 24 To 2nd developer waste
- 25 To 2nd developer waste
- 26 To 2nd developer waste
- 27 To 2nd developer waste
- 28 To 2nd developer waste
- 29 To 2nd developer waste
- 30 To 2nd developer waste
- 31 To 2nd developer waste
- 32 To 2nd developer waste
- 33 To 2nd developer waste
- 34 To 2nd developer waste
- 35 To 2nd developer waste
- 36 To 2nd developer waste
- 37 To 2nd developer waste
- 38 To 2nd developer waste
- 39 To 2nd developer waste
- 40 To 2nd developer waste
- 41 To 2nd developer waste
- 42 To 2nd developer waste
- 43 To 2nd developer waste
- 44 To 2nd developer waste
- 45 To 2nd developer waste
- 46 To 2nd developer waste
- 47 To 2nd developer waste
- 48 To 2nd developer waste
- 49 To 2nd developer waste
- 50 To 2nd developer waste
- 51 To 2nd developer waste
- 52 To 2nd developer waste
- 53 To 2nd developer waste
- 54 To 2nd developer waste
- 55 To 2nd developer waste
- 56 To 2nd developer waste
- 57 To 2nd developer waste
- 58 To 2nd developer waste
- 59 To 2nd developer waste
- 60 To 2nd developer waste
- 61 To 2nd developer waste
- 62 To 2nd developer waste
- 63 To 2nd developer waste
- 64 To 2nd developer waste
- 65 To 2nd developer waste
- 66 To 2nd developer waste
- 67 To 2nd developer waste
- 68 To 2nd developer waste
- 69 To 2nd developer waste
- 70 To 2nd developer waste
- 71 To 2nd developer waste
- 72 To 2nd developer waste
- 73 To 2nd developer waste
- 74 To 2nd developer waste
- 75 To 2nd developer waste
- 76 To 2nd developer waste
- 77 To 2nd developer waste
- 78 To 2nd developer waste
- 79 To 2nd developer waste
- 80 To 2nd developer waste
- 81 To 2nd developer waste
- 82 To 2nd developer waste
- 83 To 2nd developer waste
- 84 To 2nd developer waste
- 85 To 2nd developer waste
- 86 To 2nd developer waste
- 87 To 2nd developer waste
- 88 To 2nd developer waste
- 89 To 2nd developer waste
- 90 To 2nd developer waste
- 91 To 2nd developer waste
- 92 To 2nd developer waste
- 93 To 2nd developer waste
- 94 To 2nd developer waste
- 95 To 2nd developer waste
- 96 To 2nd developer waste
- 97 To 2nd developer waste
- 98 To 2nd developer waste
- 99 To 2nd developer waste
- 100 To 2nd developer waste

PLEASE NOTE:
Safety valve number 12
is found on all units
from all parts 1, 2 and 3

In areas where it is unacceptable to
discharge chemicals directly into the
drain, waste, developer and fix can be
refined to developer and fix pipes to
their original bottles.

IMPORTANT NOTE

When using the
replenishment system if the
developer and fix waste pipes
are not permanently plumbed in
then their pipes should be
placed into empty containers to
receive waste, all services need
to be within 3-4 ft of the stand.

36

ISS = IS Components	LF = Littlefuse	GDB or GDC = Buismann	Schurter = Circuit breaker	Issue 13
---------------------	-----------------	-----------------------	----------------------------	----------

Page 13

Component part number

PAGE	ITEM	DESCRIPTION	CODE
7	1	Control Panel	VMOG2034P
7	2	Main switch	VELC2030P
7	6	Process switch	VELC2005P
7	7	5A Fused socket	VELC2071P
7	9	Machine lid	VMOG2040P
7	10	Machine lid catch (2 parts)	VMOG2155P
7	12	Film collector	VMOG2152P
7	13	Film collector backplate	VFIT5015P
7	14	Film collector spring	SPR0001P
7	15	Daylight loader coupling knob	VFIT2035P
7	16	Daylight loader	VMAC9300A
7	17	Viewer cover	VFIT5010M
7	18	Hand loader gloves	VFIT3002A
7	19	Daylight loader catch (2 parts)	VMOG2155P
7	20	Loader lid	VMOG3009P
7	21	Regensifter pump	VELC4000P
7	22	Regensifter inlet tubes	VFIT4042P
7	26	Developer waste hose, black	VFIT2031P
7	27	Fixer waste hose, red	VFIT2030P
7	28	Machine foot	VMOG2047P
7	30	Film entry slit (2 parts)	VMOG2060P
7	32	Side panel	VMOG2045P
7	33	Back panel	VMOG2030P
7	34	Front panel	VMOG2025P
7	36	Static strip	VFIT2010P
8	1	Module turning tool	VMOG5145P
8	3	Drive shaft gear	VMOG2030P
8	4	Synthetic film transport belt (set of 2)	WV50012A
8	5	Drive dog (Supplied as set with shaft)	SA351500A
8	7	Developer module assembly	VMOD0007A
8	8	Fixer module assembly	VMOD0008A
8	9	Wash module assembly	VMOD0012A
8	10	Dryer module assembly	VMOD0017A
8	12A	Drain tube developer, black	SA355025A
8	12B	Drain tube fix, red	SA355020A
8	12C	Drain tube water, blue	SA355035A
9	2	Motor	VELC2124P
9	3	Gear drive strip	VMOG2137P
9	4	Dryer fan	VELC2090P
9	5	Terminal strip section	VELC2161P
9	6	Terminal strip end cap	VELC2162P
9	6	Dryer terminal block	VELC2160P
9	7	Power cord socket	VELC2071P
9	8	Safety switch	VELC2037A
9	9	Control pcb	VELC2400P
9	10	Main switch	VELC2000P

Component part number

PAGE	ITEM	DESCRIPTION	CODE
8	12	Processor switch	VELC2008P
8	13	Circulation pump, developer	VELC2150P
9	14	Circulation pump, fixer	VELC2150P
9	15	Control panel retaining screw (2BA) x 3/8	IFDX0085P
9	15	Control panel retaining screw (M5) x 10	IFDX1046P
9	16	Mainboard brace	IFIT2060P
10	1	Main board	VMDG2030M
10	2	Drive dog cover plate A	VMDG2135P
10	3	Main gear support strip B	VMDG2136P
10	4	Gear cover strip C	VMDG2137P
10	5	Module drive gears	VMDG2120P
10	6	Idle gears	VMDG2115P
10	7	Main drive gear	VMDG2120P
10	8-10	Supplied as set with shaft (*)	WASBY 5206A
10	8*	Thrust washer	
10	9*	Drive dog spring	
10	10*	Drive dog	
10	11	Gear strip assembly screws	IFDX0032P
10	12	Clamping nut	IFDX0009P
10	13	Gear axle rod and drive dog shaft	VMDG2130P
10	14	Motor	VELC2124P
10	15	Motor mounting screws	IFDX0075P
10	16	Dryer fan	VELC2090P
10	17	Dryer element assembly	VELC2117A
10	19	Heater grille plate assembly & thermistor	WKS00050A
10	19	Heater grille plate only	IFIT0060P
10	20	Grille screw	IFDX0011P
10	21	Fan heater retaining screw	IFDX0067A
10	23	Terminal strip section	VELC2161P
10	23	Terminal strip end cap	VELC2162P
10	24	Dryer terminal block	VELC2160P
11	6	Gear cover plate	VMDG0077P
11		Complete transport module	
		Developer	VMOD0003A
		Fixer	VMOD0008A
		Washer	VMOD 0013A
		Dryer	VMOD0017A
11	1-4	Module drive gear set (*sold as set)	VMOD0100A
11	1*	Main drive gear	
11	2*	D shaped centre drive gear	
11	3*	Small idle round hole gear	
11	4*	Large idle gear	

PAGE	ITEM	DESCRIPTION	CODE
11	6	Gear cover plate	IMDG00177
13	1	Chemical tank complete	IA3500011
13	2	Circulation pump, developer	IFLC2150P
13	3	Circulation pump, fix	IFLC2150P
13	4	Manifold assembly with pump	IA3501010
13	4	Manifold assembly only	IA3501111
13	5	Thermal cutout	IFLC2050P
13	7	Temperature sensor assembly (Dev)	IA350117A
13	7	Temperature sensor assembly (Fix)	IA350119A
13	8	Liquid level sensor screws	IFIX2076P
13	8	Liquid level sensor nuts	IFIX2100P
13	8	Liquid level sensor washers, flat	IFIX2115P
13	8	Liquid level sensor shockproof washers	IFIX2130P
13	8	Liquid level sensor sealing washers	IFIX2086P
13	9	Tank/manifold seals	IFIT2055P
13	10	Pump coupling nuts	IMDG2075P
13	11	Heater element seal	IFIT2050P
13	12	Heater element	IFLC2096P (A35
13	13	Clamping nut	IFIX2104P
13	14	Tank wall stiffener	IFIT2045P
14	1	Chemical tank	IMDG2148P
14	2	Drain tube, developer, black	IA350025A
14	3	Drain tube fix, red (Europe)	IA350030A
14	3	Drain tube water, blue (Europe)	IA350035A
14	3	Water inlet pipe	IFIT2145P
14	4	Water overflow pipe	IMDG2094P
14	5	Water drain pipe	IMDG2095P
14	6	Developer drain hose and seal, black	IFIT2031P
14	6	Fixer drain hose and seal, red	IFIT2030P
14	7	Solenoid valve	IFLC2106P
14	8-14-10	Also supplied as 10 Valve body (*)	IA352106A
14	8"	Solenoid connector	
14	14"	Connector locking nut	
14	10"	Pipe connector	
14	9	Connector locking nut (Items 4 & 5)	IMDG3075P
14	11	Replenisher tube, developer	IFIT8045P
14	12	Replenisher tube, fix	IFIT8050P
14	13	Washer rubber (solenoid only)	IFIT2047P
14	18	Male stud coupling	IFIT0004P
14	16	Bush (tank water pipe 1/2)	IMDG6035P
14	19	Foot tray	IMDG2047P
14	20	Foot tray screws	IFIX2080P
14	4	"O" ring 11.3 solenoid pipes	IFIT2077P
14	8	"O" ring per seal drain housing	IFIT2094P

Component part number

PAGE	ITEM	DESCRIPTION	CODE
16	8	hose assembly (spare)	WHT2140F
17	10	filter	WHT2028P
17	11	flexible hose - Grey	WHT2027P
		- White	WHT2026P
		- Blue	WHT2025P
17	14	Replenisher pump	WLC4000P
17	15	replacement replenishment hose	
		Assembly Dev. fix USA	WAS2141P
17	18	Deep loader extension	WHT3015P

Italy

Gendex Dental Systems
Via A. Mazoni, 44
20095 Cusano Milanino - MI
Tel +39 02 618008.1
Fax +39 02 618008.09

Germany

Kavo Dental GmbH
Bismarckring 39
D89 400
Biberach
Germany
Tel +49 408 996 8820

Kavo Gendex - Germany - Oct 07