
# BW 485/L

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- user friendly
- microprocessor controlled
- high transfusion rates
- minimized energy consumption

## Service Manual

EG  
(DR) 5

## 1. CONTINUOUS TEST

### *1.1 General*

- Start-up
- After approx. 15 min. measure the temperature at the heat exchanger
- After 24 h measure the temperature again

During this test the initial and final values must be within the tolerance range and the measurements must be taken under the same ambient conditions.

The final temperature may deviate from the initial value by  $\pm 0.5^{\circ}\text{C}$ .


Fig. 1.1.1 Measurement set-up

1    IR thermometer  
2    Heat exchanger  
(BW ring)

Measuring instrument used:

HORIBA infrared thermometer IT 330

Accuracy:  $\pm 1\%$ ,  $\pm 1$  digit at full scale

Emission: black body

Date	Created by	Release	File	Revision no.	Page
February, 1997	Riedl		BW 485 Anweisung Englisch	01	

## Temperature Adjustment

- Connect adapter box to computer and start terminal program (e.g. Procomm)
- Connect adapter box to BW 485 CON2 (on PCB)
- Set switch to „TEST“ position. Device then emits an audio signal at regular intervals
- Preset spindle trimmer:

P3	Gain 1	right-hand end position
P6	Gain 2	left-hand end position
P5	Alarm limit 41°C	left-hand end position
P4	Set temperature	right-hand end position

- Switch on BW 485
- Set temperature controller (P4) to required temperature (display on PC: „Set“). For fast presetting, press „SET“ key on adapter box in order to display set temperature on bar chart.
- Set current externally measured heater temperature with spindle trimmer P6 (display on PC: „SENS2“)
- Turn spindle trimmer P3 slowly to left until device starts to heat up. When heating is on L13 is illuminated. Slowly increase heater temperature to set value.

The two temperatures displayed on the PC as „Sens1“ and „Sens2“ must be as close to each other as possible ( $\pm 0.3$ ) otherwise an alarm is triggered.

The difference between the two temperatures is displayed on the PC under „Diff1“. Because measurements are taken at different points, the internally and externally measured temperatures may deviate slightly from one another.

### Setting the high temperature alarm limit:

1. Set switch in adapter box to „Test“ position
2. Set temperature to 41,5°C with P4 and wait until temperature is constant
3. Turn P5 slowly to right until relays drop out and alarm is triggered
4. Let device cool down. Heat it up again and check cut-off point again
5. Increase set temperature (39°C) slowly to 41°C and control alarm shut off
6. Reset at previous temperature

## 4. Temperature Adjustment

### Service interface:

The BW 485 is equipped with a service interface. This serial interface outputs a range of data which can be very useful for trimming and trouble shooting.  
Connect the BW485 to your PC or terminal with the adapter box. The data can be displayed using any common terminal program e.g. Procomm.  
Data format: 4800,n,8,1

### Key to serial output data:

Set	Set value
Sens2	Temperature measured by Sensor2 (temperature signal for display and 41°C safety cut-off)
Sens1	Temperature measured by Sensor 1 (control temperature)
Diff	Difference between set and actual (only when actual < set, otherwise 00.0)
P%	Current heating output in % of maximum power-rating
Diff 1	Difference between Sens1 and Sens2
Err	Error counter
t	Time until low temperature alarm is activated (0=on)


### Error codes:


The BW485 indicates various errors and displays them on the bar chart as follows:

-  Low temperature alarm (1<sup>st</sup> LED flashes)
-  Temperature difference between Sens1 and Sens2>0.5°C
-  Sens1>41°C
-  Sens2>41°C

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
## EC Certificate of Conformity

No. 0001 V1 / 08.96

Authorized agent: E. BIEGLER G.m.b.H.  
Medical Electronics

Address: Allhangstrasse 18 a  
3001 Mauerbach  
Austria

Manufacturer: E. BIEGLER G.m.b.H.  
Medical Electronics  
3001 Mauerbach  
Austria


Product description: BW 485  
from serial no. V7 2001


BW 485L  
from serial no. V7 7001

The product described above is a medical product (as per 93/42/EEC Section 1, Para. 2a) and conforms to the provisions of the following European Directives:

93/42/EEC Council Directive on the harmonization of statutory provisions for medical products.

E. Biegler G.m.b.H.  
Mauerbach, Austria, 28.05.97  
(Date)

F. Netauschek/Development Manager .....   
(Name, Position)


I. Biegler/Managing Director .....   
(Name, Position)

This document certifies conformity with the above directives but does not constitute warranty of a quality as defined by the Austrian product liability law.

The safety regulations in the product documentation supplied must be observed.

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
*W. WJS*

# QUALITY ASSURANCE

## CERTIFICATE OF COMPLIANCE OF THE QUALITY SYSTEM WITH STANDARDS MENTIONED BELOW

Certificate No.: Q1 95 03 14553 005

The Certification Body for Quality Systems  
of TÜV PRODUCT SERVICE GMBH certifies that

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for the following area:

Ventilators and respiratory training apparatus,  
Blood and infusion warmers,  
Fibrin sealant application systems,  
Cooling systems for medical solutions,  
Olfactometers, Electrical stimulation therapy devices  
Laboratory and diagnosis equipments

has established and is operating a quality system  
which meets the requirement(s) of

ISO 9001 : 1994 / EN 46 001 : 10/93

as documented in the audit report no. 01/338-5-003/02

This certificate is valid until March 1998.  
provided that a periodical surveillance is conducted.

Munich, March 24, 1995

TÜV PRODUCT SERVICE GMBH  
ACCREDITED CERTIFICATION BODY  
for Quality Systems

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Biegler GmbH