


Anaesthetic Ventilators

Ventilog



Anaesthetic Ventilators

Ventilog

Technical training worksheets

Anaesthetic Ventilators

Ventilog

Technical Data

Operating principle	Bellows-in-bottle ventilator with primary/secondary system
Control principle	pneumatic, time-cycled, volume constant
Ventilation frequency	Ventilog: 6 to 40/min $\pm$ 20 % Ventilog 2: 6 to 60/min $\pm$ 15 %
I:E ratio	Ventilog: 1:2 $\pm$ 20 %, permanently set Ventilog 2: 1:1, 1:2, 1:3, $\pm$ 20 % adjustable
Switch «Q»/«I»	Main switch for Ventilog
Switch «Q»/«I:E»	Main and selector switch for Ventilog 2
Inspiratory flow	Ventilog: 20 to 80 L/min $\pm$ 15 %, infinitely adjustable Ventilog 2: 20 to 80 L/min $\pm$ 15 % at 20 mbar counter pressure, infinitely adjustable
Tidal volume	50 to 150 mL with bellows K for children 150 to 1600 mL with bellows E for adults
Internal compliance	2.3 mL/mbar with bellows K for children, corrugated hose 1 m 3.0 mL/mbar with bellows E for adults, corrugated hose 1 m 6.4 mL/mbar with bellows E, corrugated hose 1 m, switching valve, circle system 7 a/8 ISO with 2 absorbers
Minute volume (as per ISO)	up to 25 L/min with a fresh-gas flow of 4 L/min
Working pressure	30 mbar $\pm$ 15 %, permanently set
PEEP	up to 19 mbar $\pm$ 20 %, infinitely adjustable
Compressed-gas supply	Oxygen or oil-free compressed air from a central supply system or from compressed-gas cylinders. Connecting thread M 15 x 1, male
Compressed-gas consumption	Ventilog: $\frac{1}{3}$ of the set inspiratory flow Ventilog 2: 16 L/min (at I:E = 1:2)
Patient system	autoclavable bellows, removable after releasing slide-in unit
Anaesthetic-gas disposal	Excess exhaled gas is routed to the anaesthetic gas exhaust via an exhaust socket at the rear of the unit
Dimensions	W x H x D = 212 x 266 x 300 mm
Weight	12 kg

Deviations of settings and operating ranges are indicated pertinent to desired values.
The frequency/control elements are calibrated with oxygen.
If compressed air is used as drive gas, this results in an increase of frequency of approx. 10 %.



Anaesthetic Ventilators

Ventilog

*Cleaning and disinfection or steam sterilisation
of Ventilog 12 and accessories*

	cleaning			disinfection				sterilisation		
	running water	wiping	bath	bath	wiping	aseptor	purfactor	gas	120°C	steam 734g
Ventilog - Basic unit		X			X	X				
Patient system	X					X		X	X	X
Switching valve	X					X		X	X	X
Exhaust gas valve	X					X		X	X	X
Corrugated bellows	X			X		X	X	X	X	X
Connecting tube to circuit system	X		X			X		X	X	X
Barolog		X			X	X				
Pre-om - pressure gauge		X			X	X		X	X	
Pre-om - alarm unit		X			X	X		X		




Note!

Aluminium parts such as the switching valve and patient system may not be cleaned in the Dräger Purfactor or similar machines. Mechanical cleaning of such parts with chlorinated cleansers may cause corrosion and is therefore not permissible.

Anaesthetic Ventilators

Ventilog

Functional diagram



Key:

- 1 Connection M 15 x 1
- 2 Filter
- 3 Pressure regulator 8402745
- 4 Pressure regulator 8402745
- 5 3/2-way valve 8403006
- 6 Flow control valve 8404680
- 7 Switch "Automatic-manual, spontaneous" 8402441
- 8 Plug-in coupling with check valve 8404948
- 9 Bellows 2M 8138
- 10 Pressure chamber
- 11 Silencer)
- 12 Vent valve) Vent valve 8404660
- 13 Silencer)
- 14 Excess gas discharge valve
- 15 Frequency control valve 8403845
- 16 Vessel 8404675
- 17 Vessel 8404673
- 18 PEEP valve 8404730
- 19 Metering unit 8403795
- 20 Metering unit 8403795
- 21 Plateau valve 8404404 (single)
- 22 Metering unit 8403795
- 23 Anaesthetic exhaust socket 8404690


A, E, B, F Amplifier relays M 25716
C, D, G, H, I, K Amplifier relays M 22260

Anaesthetic Ventilators

Ventilog

The oldest one
Functional diagram, units up to serial No. 3218

Unit ~~off~~ gas on



V_1 = inspiratory time

V_2 = expiratory time

A, E, B, F Amplifier relays M 25716
C, D, G, H, I, K Amplifier relays M 22260

Anaesthetic Ventilators


Ventilog

Modifications as of device No. 3219 *include also*


1. Changes have been made to the material used for the diaphragm 8403377 installed in the pressure-chamber vent valve (12) and in the surplus-gas discharge valve (14). The new diaphragms can be recognized from their colour (old version red, new version translucent).

Due to these modifications a pressure divider (24) had to be additionally fitted to ensure that the diaphragm in the valve (14) is actuated with a higher pressure than the diaphragm in the valve (12). This measure was necessary in order to prevent vibration of the diaphragm in the valve (12).

Metering pipe M 27297 Pipe 3 x 1 DIN 1754 SF-CU F22



Reason for modification: The intention was to prevent a PEEP of > 2 mbar in the system in the "0" setting of the PEEP valve (18).




Anaesthetic Ventilators

Ventilog

2. A metering unit M 27297 is to be placed on the plateau valve 8406266. This metering unit is not to be set to a fixed value, but will rather be manufactured as required. The copper pipe will be pressed together until the vibration noise of the diaphragm in the plateau valve can no longer be heard.

Reason for modification: Noise reduction



3. The metering units (19) and (20) are being converted to fixed metering units.

Reason for modification: No possibility of self-adjustment.


4. A filter 8405241 is being located ahead of the metering unit (22).

Reason for modification: Protection of metering unit (22) and ratio valve (15) against dirt.

Anaesthetic Ventilators

Ventilog

Functional diagram as of device 3219



Key	1	Connection M 15 x 1	A, E, B, F, D	Amplifier relays M 25716
	2	Filter	C, G, H, I, K	Amplifier relays M 22260
	3	Pressure regulator 8402745		
	4	Pressure regulator 8402745		
	5	3/2 -way valve 8403006		
	6	Flow control valve 8404680		
	7	Switch "0 - 1" 8402441		
	8	Plug-in coupling with check valve 8404948		
	9	Bellows 2M 8138		
	10	Pressure chamber		
	11	Silencer)		
	12	Vent valve) vent valve 8404660		
	13	Silencer)		
	14	Excess gas discharge valve		
	15	Frequency control valve 8403845		
	16	Vessel 8404675		
	17	Vessel 8404673		
	18	PEEP valve 8404730		
	19	Metering unit 8406292		
	20	Metering unit 8406314		
	21	Plateau valve 8406266 (single)		
	22	Metering unit 8403795		
	23	Anaesth. exhaust socket 8404690		
	24	Metering unit M 27297		
	25	Metering unit M 27297		
	26	Filter		

Anaesthetic Ventilators


Ventilog

Functional diagram as of device No. 3219

middle age

8403377

428475



we be Ventilog!
disassemble this device
and take the clipper
 428475 and 8403377

Unit off, gas on


Anaesthetic Ventilators

Ventilog

Drive gas

Oxygen or compressed air passes via the connection M 15 x 1 (1) through a filter (2) to the pressure reducer (3), which constantly reduces the pressure of the gas to 1.8 bar. From the pressure reducer the gas is passed to the 3/2-way valve (5). At the same time the gas coming from the connection (1) passes into the pressure reducer (4), which constantly reduces the pressure of the gas to 1.5 bar. From there the gas is passed to the switch "Automatic-manual, spontaneous" (7) and the relay D.

After moving the switch (7) to "Automatic" the gas passes to the distributor "—0" and the plug-in coupling (8). Via a distributor this pressure is applied to all points marked with "0—".



Expiration
Control group

The expiration signal is applied as a result of the relay D being driven.

PEEP setting


The expiration signal is supplied to the relays A, E, H, I, K. As a result the inspiration and frequency signals in relays H and A respectively are blocked. Relay I is driven such that a connection is made between the PEEP valve (18) and the excess gas discharge valve (14). (In the PEEP valve the diaphragm is pretensioned against the seat by way of spring. Via a metering unit (19) the system is filled up to the set spring pressure with a constant flow (0.5 l/min); the remainder is discharged into the atmosphere).

Anaesthetic Ventilators


Ventilog

At the same time the gas coming from the metering unit (20) (the plateau valve (21) has the task of limiting the pressure in the pressure chamber and patient system to 80 mbar) is blocked off in relays I and K.

The 3/2-way valve (5) is also blocked off.




The gas coming via the metering unit (22) into the frequency valve (15) and from there into relay E is routed into the vessel V_2 (17). A time-delayed build-up of the switching pressure for relay F is effected via the vessel (17). Following switchover the output of relay F is relieved of pressure and the 1.5 bar in relay G are passed to relay H which conducts.



Anaesthetic Ventilators

Ventilog


With the commencement of inspiration the relays I, K, D, E, A and the 3/2-way valve (5) are driven via relay H. In relay I the PEEP signal is blocked and as a result the gas passes from the plateau valve (21) to the surplus gas discharge valve (14) which is closed. Via relay K gas also passes from the plateau valve (21) to the vent valve (12) which is closed. In relays D and E the expiration/frequency signals are blocked.



The 1.8 bar at the 3/2-way valve (5) pass via the flow control valve (6) and the silencer (11) into the pressure chamber (10). With the valve (12) closed the bellows (9) are compressed and the gas from the bellows passes via the switching valve into the circle system and from there to the patient.

Anaesthetic Ventilators


Ventilog



The gas coming from the frequency valve (15) is routed into the vessel V_1 (16) via relay A.


When the switching pressure is attained the output of relay B is relieved of pressure (as with expiration) and the 1.5 bar from relay C are passed to relay D as a result of which the expiration signal is again switched through.

Expiration



The 3/2-way valve (5) is reversed thereby ending inspiration.

The pressure on the vent valve (12) is relieved thus enabling the gas in the pressure chamber (10) to escape into the atmosphere via the silencer (13). As a result of their own weight the bellows (9) then unfold and fill up.




Anaesthetic Ventilators

Ventilog

Pneumatic switching valve M 27235 / ISO M 27240



Depending on the setting "automatic" or "manual, spontaneous", the switching valve routes the gas coming from the Ventilog or manual ventilation bellows to the circle system.



Anaesthetic Ventilators

Ventilog

Location of components:



*adj.
working pressure*


Anaesthetic Ventilators

Ventilog

Sub-assemblies

Pressure regulator 8402745 (3, 4)

Given a supply pressure of 2 - 6 bar the pressure reducer keeps the downstream pressure at a constant 1.5 bar or 1.8 bar.




The pressure reducer (3) (set to $1.8 \pm 0.1 \text{ bar}$) generates the drive gas for the bellows.

The pressure reducer (4) (set to $1.5 - 0.1 \text{ bar}$) generates the control pressure for the circuit.

After loosening the locking screw the downstream pressure (P_H) can be changed using the cap. The cap is then to be secured in position again using the screw.

3/2-way valve 8403006 (5)

The 3/2-way valve operates in the same manner as the amplifier relays, however with a greater capacity.




Anaesthetic Ventilators

Ventilog

Plateau valve, single 8404404 (21)


The plateau valve 8404404 is used to limit the working pressure. This is permanently set in the range 80 - 95 mbar.

(corresponding to
min. - max. flow)




3/2-way valve 8402441 (7)

The 3/2-way valve is used as a changeover switch.



Amplifier relays M 22260 (C, D, G, H, I, K)

The identification code for M 22260 is UO B. The double diaphragm relays are subjected to a functional check (dynamic intake) and a leakage test.



Amplifier relays M 25716 (A, B, E, F)


The identification code for M 25716 is UO BS. These double diaphragm relays are additionally subjected to checking of the switching characteristic.

Anaesthetic Ventilators

Ventilog


Precision metering valve 8403845 (15)

The precision metering valve 8403845 is used for frequency adjustment. The frequency is altered by changing the annulus.




PEEP valve 8404730 (18)

During expiration the PEEP valve limits the ventilation pressure (0 - 15 mbar).




Flow valve 8404680 (6)

The flow valve 8404680 is used to regulate the working flow. The flow is altered by changing the annulus.




Anaesthetic Ventilators

Ventilog



connections for venting




Patient system 8404700 (9, 14)


The patient system is designed to separate the patient air from the drive air of the device. During inspiration the valve 14 shuts off the path to the anaesthetic exhaust socket and limits the pressure to 80 - 95 mbar by way of the plateau valve 21. During expiration the PEEP pressure is controlled via the valve 14.

surplus gas discharge
valve 14

to anaesthetic exhaust socket



to circle system




Anaesthetic Ventilators

Ventilog

Anaesthetic exhaust socket 8404690 (23)

The anaesthetic exhaust socket is used for removing the excess fresh gas. The built-in check valve prevents the intake of "extraneous air".




8403377

gettyimages

- Gas in, Luft out

Handwritten sketches of various mechanical components, labeled with letters a through j. The sketches include:

- a**: A gear or pulley with a central shaft.
- b**: A shaft with a central hole and a small circular feature on the side.
- c**: A shaft with a central hole and a small circular feature on the side.
- d**: A shaft with a central hole and a small circular feature on the side.
- e**: A shaft with a central hole and a small circular feature on the side.
- f**: A shaft with a central hole and a small circular feature on the side.
- g**: A shaft with a central hole and a small circular feature on the side.
- h**: A shaft with a central hole and a small circular feature on the side.
- i**: A shaft with a central hole and a small circular feature on the side.
- j**: A shaft with a central hole and a small circular feature on the side.



no pulse during
resp.

Anaesthetic Ventilators

Ventilog

Explanation of symbols:

OK	: —	Check condition	C
Defect	:	Check function	O
Spare parts used	: ○	Check for leakage	L
Report	: /	Enter test value	V

Location of
apparatus:


Serial No. :
Date of delivery :
Startup :

Invoice No. or
delivery No. :



Others :

The anaesthetic apparatus and circle system are to be tested in accordance with the Test Card before checking the Ventilolog.

1. Pneumatic switching valve 8404950
Pneumatic switching valve ISO 8405276 C
- 1.1. Open valve
- * 1.1.1. Diaphragm 8401650 2 x
Replace every 2 years C
- 1.1.2. Thrust piece 8404945 C
- 1.1.3. O-ring 8402234 C
- 1.1.4. O-ring 2M 8777 2 x C
- 1.2. Assemble valve
Leak test





Test set-up:
Seal connection C with cork
Connections B and D open.
Build up pressure of 100 mbar via
connection A.



Pneumatic switching valve


- 1.2.1. Test value
Gas flow ≤ 0.2 L/min 0
- 1.2.1.1. Test set-up:
Remove cork from connection C,
otherwise as item 1.2.
- 1.2.1.1.1. Test value:
Pressure drops to 0 bar 0
- 1.2.2. Test set-up:
Seal connection C with cork.
Connection A open.
Build up pressure of 1.2 bar at
connection D.
Build up pressure of 100 mbar via
connection B.
- 1.2.2.1. Test value:
Gas flow ≤ 0.2 L/min 0
- 1.2.2.2. Test set-up:
Remove cork from connection C,
otherwise as item 1.2.2.
- 1.2.2.2.1. Test value:
Pressure drops to 0 bar 0
- 1.2.3. Test set-up:
Build up pressure of 1.2 bar at
connection D. All other
connections open.
Seal actuator 1.2 bar.
- 1.2.3.1. Test value:
Pressure must not drop in 5 s. L
- 1.2.4. Build up alternating pressure of
1.2 bar and 0 bar at connection D,
otherwise as Item 1.2.3.
- 1.2.4.1. Test value:
Visual inspection through connection
C. Tappet must move from stop to
stop in each case without visible
delay when switching. 0
- 1.3. Manual switching valve 8405305
Manual switching valve ISO 8405295 C
- 1.3.1. Open valve
- 1.3.1.1. Tappet and spring C
- * 1.3.1.2. Sealing rings 8405281
Replace every 2 years C
- 1.4. Assemble valve
Leak test
Test set-up:
Seal connection C with cork
Connection B open.
Build up pressure of 100 mbar
via connection A.
Switch setting man./spont.



Manual switching valve

39 754

- 1.4.1. Test value:
Gas flow ≤ 0.2 L/min
- 1.4.1.1. Test set-up:
Remove cork from connection C,
otherwise as Item 1.4.
- 1.4.1.1.1. Test value:
Pressure drops to 0 bar
- 1.4.2. Test set-up:
Seal connection C with cork
Connection A open.
Build up pressure of 100 mbar via
connection B.
Switch setting automatic
- 1.4.2.1. Test value:
Gas flow ≤ 0.2 L/min
- 1.4.2.1.1. Test set-up:
Remove cork from connection C,
otherwise as Item 1.4.2.
- 1.4.2.1.1.1. Test value:
Pressure drops to 0 bar
- 1.4.3. Switch lever alternately
- 1.4.3.1. Test value:
Visual inspection through connection C.
Tappet must switch.
Tappet pressure must be noticeable in
all lever positions.
- 1.5. Pneumatic switching valve M 27235
Pneumatic switching valve ISO M 27240
- 1.5.1. Open valve:
* Diaphragm 8403352
* O-ring CH 13972 (replace annually)
- 1.5.2. O-ring 8402234
- 1.5.3. O-ring 2M 8777 2x
- 1.6. Assemble valve - see Item No. 11,
File 5324 - Leak test
Test set-up:



Anaesthetic Ventilators


Ventilog

02.85 Steinwarder No. 11

Pneumatic switching valve M 27235
Pneumatic switching valve ISO M 27240

Replacement of diaphragm 8403352 and O-ring CH 13972 (required once a year in accordance with Test card). The following procedure is to be adopted when replacing the diaphragm:

- Loosen screw 1328255
- Unscrew cover M 27242
- Loosen clamping piece M 27244 (fitted with silicone rubber Elastosil E 41)
- Remove insert M 27229
- Replace O-ring CH 13972
- Thoroughly grease tappet and O-ring
- Fit insert M 27229. Observe preferred position!
- Replace diaphragm 8403352
- Re-install clamping piece M 27244 using silicone rubber Elastosil E 41
- Grease O-ring 2M 10296 and cover M 27242
- Screw in cover such that O-ring projects halfway into valve housing.



- Seal connection C, apply a flow of 0.5 L/min. to connection A
Run a hose 2 x 1 into a water bath from connection D.
- Now screw in cover until air bubbles no longer emerge from the hose
(float of flowmeter goes to "0").
- Secure cover with screw (fitted with silicone rubber Elastosil E 41).

Technical Customer Service

B. J. Regel

Test set-up:
Seal connection C with cork.
Connections B and D open.
Build up pressure of 100 mbar via connection A.

1.6.1. Test value:
Gas flow ≤ 0.2 L/min 0

1.6.1.1. Test set-up:
Remove cork from connection C,
otherwise as Item 1.6.

1.6.1.1.1. Test value:
Pressure drops to 0 bar 0

1.6.2. Test set-up:
Seal connection C with cork.
Connection A open.
Build up pressure of 1.2 bar at connection D.
Build up pressure of 100 mbar via connection B.

1.6.2.1. Test value:
Gas flow ≤ 0.2 L/min 0

1.6.2.2. Test set-up:
Remove cork from connection C,
otherwise as Item 1.6.2.

1.6.2.2.1. Test value:
Pressure drops to 0 bar 0

1.6.3. Test set-up:
Build up pressure of 1.2 bar at connection D.
All other connections open.
Seal actuator 1.2 bar.

1.6.3.1. Test value:
Pressure must not drop in 5 s. L

1.6.4. Build up alternating pressure of 1.2 bar and 0 bar at connection D.
Otherwise as item 1.6.3.

1.6.4.1. Test value:
Visual inspection through connection C.
Tappet must move from stop to stop in each case without visible delay when switching. 0


1.7. Gas supply "pneumatic switching valve"

1.7.1. Plug-in coupling 8404948 C 0

1.7.2. Plug-in nipple 8404951 C 0

1.7.3. Hose 2 m, silicone hose 2 x 1.5 1190520 C

1.7.4. Gas supply
Test set-up
Connect plug-in nipple to pressure gauge.
Switch setting "automatic" or "I"



1.7.4.1. Test value:
 $P = 1.5 \pm 0.1 \text{ bar}$

0

1.7.4.2. Test set-up:
 Switch setting "Man/Spont" or "O",
 otherwise as Item 1.7.4.

1.7.4.2.1. Test value:
 $P = 0 \text{ bar}$

0

2. Remove Ventillog from anaesthetic apparatus. Swivel waste-gas socket 8404690 through 45° and pull out. Detach locking screws on slide-in unit, pull out slide-in unit.

2.1. Gas supply

2.1.1. Connection M 15 x 1.

C


* 2.1.1.1. Filter element D 11521
 Replace every 2 years

C

2.1.2. Pressure regulator
 "Working pressure"

* Diaphragm 8404049 (replace every 6 years)



* Inlet valve 8404048 (replace every 6 years)



2.1.2.1. Pressure regulator "working pressure"
 Downstream pressure
 Remove hose screw connection
 Connect valve 2 M 5223 with test hose 7900261.

2.1.2.1.1. Test value
 $P_{\text{stat}} \leq 1.9 \text{ bar}$
 (Flow must be guaranteed)



0 L



*Test point, in older units
 not connected. For modification
 use a.u. parts.*

3. Waste-gas socket and patient system
 Insert slide-in unit and secure screws.
 Remove cover plate.


- 3.1. Waste-gas socket 8404690 C
- 3.1.1. O-ring E 20278 C
- 3. Check valve M 23543 C
- 3.2. Patient system 8204700 C
- 3.2.1. Connection socket C
- 3.2.1.1. O-ring 2M 8777 C
- 3.2.2. Rotary-knob locking
Pull out rotary knob and turn in an
antil-clockwise direction. C
- 3.2.2.1. Locking mechanism released, pull
out patient system O
- 3.2.3. As of device no. 2154 insert
slide-in unit without patient
system. Locking es effected. O
- 3.2.3. Seal (pressure chamber) C
- * Lip seal 8404064 C
- * Seal 8404065 C
- Repalce every 2 years
(fit using silicone adhesive)
- 3.2.6. Leak test: Patient system -
pressure chamber
Test set-up:
Remove bellows.
Seal pressure-chamber vent with
cork. Seal waste-gas socket.



- 2.6.1. Test value:
Patient system locking O
- Gas flow ≤ 1 L/min L

Anaesthetic Ventilators


Ventilog




- 3.3. Leak test, actuation of fresh-gas discharge valve
Test set-up:

Build up pressure of 100 mbar via socket  Pressure gauge

Ventilog



Flowmeter



Detach actuation hose of fresh-gas discharge valve

- 3.3.1. Test value:
Gas flow ≤ 0.05 L/min.

0

- 3.4. Servo pressure test (PEEP control pressure) for discharge valve

Test set-up:

Switch setting "automatic" or "I"

PEEP max.

Connect pressure gauge between DMRI and discharge valve 14 via T-piece

- 3.4.1. Test value:
P = 12 + 1 mbar

0 V

4. Mount apparatus (without bellows and cover plate) and assemble ready for operation.

- 4.1. O₂/Compressed-air connecting
M 17670

C L

- 4.2. Connecting hose M 25050

C L

- 4.3. Hose 1/0.07 m 8404733

C L

- 4.4. Hose 1/1 m 8404731

C L

- 4.5. Hose 2/1 m 8404748

C L

- 4.6. Hose 2/1.5 m 8404732

C L

- 4.7. Ventilator - anaesthetic apparatus attachment


C O

- 4.8. Base plate, latching of Ventilator

C O

5. Functional check
Drive pressure 2 - 5 bar O₂

Anaesthetic Ventilators



adjusted to
Peep-max.

5.1. Working flow
 Test set-up:
 Switch setting "automatic" or "I"
 Ratio min.
 PEEP 0
 Remove bellows, *remove cork!*
 Seal waste-gas socket
 Connect patient system outlet,
 Ventillog to flowmeter 10 - 100 L/min.

5.1.1. Test value:
 Flow min. = 20 ± 2 L/min. 0
 Flow max. = 80 ± 8 L/min. 0

5.2. Working pressure
 Test set-up:
 As Item 5.1.
 However connect patient system
 outlet, Ventillog to pressure gauge.

5.2.1. Test value:
 Flow min. = 86 ± 9 mbar 0
 Flow max. = 86 ± 9 mbar 0

adjustable with plateau valve 21

5.2.2. Safety pressure
 Test set-up:
 As Item 5.2.
 Actuate pressure relief valve with
 150 mbar, remove waste-gas socket.
 Caution! Remove cork.
close control hose with plug!

*Anders als
bei Ventillog II*

5.2.2.1. Test values:
 Flow min. = in each case min. 5 mbar
 in excess of values 0
 Flow max. = from Item 5.2.1. 0
 However ≤ 100 mbar 0

5.2.3. Test set-up:
 As Item 5.2.,
 however fit bellows
 Volume 500 mL, waste-gas socket
 open

5.2.3.1. Test value:
 Flow min. = 80 ± 9 mbar 0
 Flow max. = 80 ± 9 mbar 0

5.3. Ratio setting
 Test set-up:
 Switch setting "automatic" or
 "I"
 Flow max.
 PEEP 0
 Volume 500 mL
 Patient system outlet, Ventillog
 open.


5.3.1. Test value:
 $6 = 6 \pm 1 \text{ min}^{-1}$
 $15 = 15 \pm 2.5 \text{ min}^{-1}$
 $30 = 30 \pm 4.5 \text{ min}^{-1}$ } *EDM 13* 0 V
 0 V
 0 V

Anaesthetic Ventilators

Ventilog


pt. 5.2.2.

Attention: for this procedure the unit must be in the condition than unit with serial no > 3218!



Ventilog

Anaesthetic Ventilators




Pl. 5.3

Anaesthetic Ventilators

Ventilog


4 plateau valve 21



Anaesthetic Ventilators

Ventilog

pt 5:3



5.4. I : E - ratio
 Test set-up:
 as Item 5.3.
 Ratio min

5.4.1. Test value:

Reference value = $\frac{\text{Inspiration time}}{\text{Expiration time}}$

Reference value = 0.43 - 0.58


0

5.5. Check valve function - waste-gas socket
 Test set-up:
 As Item 5.4.
 During inspiration, in plateau phase, seal
 patient system outlet; switch apparatus
 to "man./spontan." or "0".

5.5.1. Test value:
 Bellows must not move downwards

0

5.6. PEEP setting
 Test set-up:



Slip anaesthetic gas filter with tubing
 over waste-gas socket.

Setting "automatic" or "I"
 Flow max
 Ratio min
 Volume 500 mL

5.6.1. Test value:
 PEEP setting on right stop
 Fresh-gas flow 2 L/min
 PEEP $\geq$ 15 mbar

0 V

~~Fresh-gas flow 4 L/min~~
~~PEEP $\geq$ 19 ± 3.5 mbar~~

0 V

Fresh-gas flow 8 L/min
 PEEP $\leq$ 23 mbar

0 V

PEEP setting on left stop
 Fresh-gas flow 8 L/min
 PEEP $\leq$ 2 mbar

0 V

Apparatus must not hum or whistle in
 any setting.

104 6

*Attention!
 Look for control pressure
 12 ± 1 mbar (pt. 3.4)
 out of the tolerance
 change diaphragm from valve 14*

6. Make device available to the user
in a ready-to-operate condition.
Expenditure on such work is classed
as repair services and not included in
the Inspection Service price.

7. Confirmation of test

Name:

Date:

8. Report:

colored tubes, if available:

blue = supply

green = inspiration

red = expiration

brown = control line zoombar

orange = relief

yellow } = free
(natural) clear }

Anaesthetic Ventilators

Ventilog

Modification of Ventilator

1. To get safe starting conditions and no interruptions of the pneumatic time control system during work, we assembled 2 "Non return" valves into the tubes to MR 1 and MR 6.

Ord.no 8406909

In general we plan not to fix this non return valves automatically in older units. Please do this only in case of a.u. problems.



2. After fixing non return valve control and adjust the frequency system when the unit is switched off, but the gas supply is connected it is possible to assemble one more Diaphragm valve. In this case MR 14.

Ord.no 8407324

A additional volume CSS

Ord.no 840389 (fixed with clamp G50 236) with metering unit R6

Ord.no. M27297

will give the safety that the unit is starting after switching "ON" with the exhalation phase.

For old units only in case of problems, not in general.

Anaesthetic Ventilators

Ventilog

3. Loudness of pri valve (21)

To avoid giving noise during the exhalation phase

from the pri valve (21), we change the connections of

MR5 and MR10. So the pri valve is not supplied with control gas during the exhalation phase.


The metering unit R4 is only connected, if it is necessary.

For old we plan not in general the modification.

All units since approx. July 1986 are in the new condition.

Anaesthetic Ventilators

Ventilog



Van 4/09
with modification from July 1986