PADPRO®

Rx ONLY

English

PadPro® IFU

MINI INFANT, <3 kg (2602) and INFANT,
<10 kg (2603)
Non-Sterile

Instructions for Use:

Defibrillation/Pacing/Cardioversion/Monitoring Electrodes

Read and become familiar with all instructions, warnings and cautions before using this product.

Indications for Use:

The CONMED PadPro® radiotranslucent external multifunction electrodes (MFEs) are indicated for use by trained medical professionals in medical facilities to deliver energy for defibrillation, cardioversion, external pacing, and ECG monitoring applications. The MFE is a non-sterile, disposable device for single patient use only. The electrodes provide the conductive interface between the defibrillator and/or the external transcutaneous (noninvasive) cardiac pacemaker and the patient's skin.

Not for use in AED mode.

Manual Use:




When used in manual mode, i.e. under direction of a qualified health care professional, the electrode is intended for use on infant patients (3-10 kg) and mini-infant patients (<3 kg). Follow American Heart Association (AHA) guidelines for administration of energy levels, which recommends a first dose of 2 J/kg, and subsequent doses of 4 J/kg. During Refractory ventricular fibrillation, do not exceed a maximum energy level 10 J/kg.

WARNING: DO NOT use PadPro® MFEs with any model of defibrillator or therapy cable that are not listed on the following table. Defibrillator models that are not listed have not been qualified for use with CONMED MFEs and MFE performance may be impacted.

Compatibility:

- All multifunction electrodes are compatible with biphasic defibrillators.
- PadPro® Infant and Mini Infant multifunction electrodes have been tested with the following defibrillators:

PadPro® MFE Connector Styles			Compatible Equipment Connector Styles		
Connector	Label Color	Connector Style	Defibrillator Model	Therapy Cable	PadPro® MFE Adapters
	Green	Physio-Control QUIK COMBO (M)	LIFEPAK 15	11113-000004	ZMMF
			LIFEPAK 20/20E	11110-000040	
	Orange	Zoll (Z)	R Series R Series Plus	8009-0749* 8009-0750*	MMZF
			X Series Propaq MD	8300-0783	
	Yellow	Anderson	Compatible with the following PadPro® MFE Adapters • MMUF (Physio-Control Male Universal Female) • ZMUF (Zoll Male Universal Female)		

*Note: For Defibrillation/Pacing/Cardioversion/Monitoring Only

PadPro® MFEs and MFE adapters comply with CISPR 11 Class B, Group 1 Emissions and Home Healthcare Immunity Test Level requirements per ANSI AAMI IEC 60601-1-2, 4th Ed. when used with specified compatible devices in the above table.

PadPro® 2603 (Infant) series MFEs include a 60" lead wire. PadPro® 2602 (Mini-infant) series MFEs include a 24" lead wire.

Warnings:

- Not for use in AED mode.
- If the patient is less than 8 years old or weighs less than 55 lb (25 kg), use CONMED Infant or Mini-Infant series in manual mode. Do not delay therapy to determine the patient's exact age or weight.
- Pediatric energy levels must be manually set to site-specific clinical protocols.
- Misuse (including reuse) of multifunction electrodes, use of compromised or altered product, and/or failure to follow product instructions may result in patient burns, inadequate delivery of therapy, and/or loss of ECG trace quality.
- Defibrillation current can cause operator or bystander injury. Do not touch patient, or equipment connected to the patient, during defibrillation.
- Do not conduct chest compressions through the pads. Doing so may cause damage to the pads that could lead to the possibility of arcing and skin burns.
- After patient movement due to muscle contraction or patient repositioning, press pads to skin to ensure good coupling between pads and skin.
- Do not discharge hand-held paddles through these multifunction electrodes.
- Do not allow pads to touch each other or other ECG monitoring electrodes, lead wires, dressings, etc.
- Apply a new pair of electrodes if the electrodes do not adhere properly to the patient.

- Do not use if packaging appears to be compromised. Damaged packaging may cause the conductive polymer gel to dry out.
- Do not use multifunction electrodes in the presence of flammable agents or in an oxygen enriched atmosphere as this presents an explosion and fire hazard.

Storage: Do not fold, trim, crush or store under heavy objects. Avoid storage of electrodes where they may become subject to excessive heat or cold. Multifunction electrodes should be stored in unopened pouches at room temperature.

Suggested Placement: The packaging diagrams for electrode placement are suggestive of Advanced Cardiac Life Support (ACLS) and American Heart Association (AHA) guidelines. As with any multifunction electrodes, follow hospital protocols.

Pre-Connect Option (Models XXXX-PC):

- 1.) Plug the connector into the therapy cable.
- 2.) Do not open electrode compartment until ready to use.
- 3.) Follow instructions below for use.

Instructions for Use: (all models)

WARNING: Electrodes must be used before date indicated on package. Do not open package until immediately prior to use. Dried conductive polymer gel can lead to skin irritation or burns, and can potentially impact MFE performance.

Preparation: Ensure that the patient's skin is clean and dry. Remove any substances from the skin (e.g. lotions, oils, etc.).

- 1.) Tear open the pouch.
- 2.) Using the extended tab of the release liner, remove the release liner from the multifunction electrodes to expose conductive and adhesive areas. (Arrow indicates peel location.)

WARNING: Visually examine the pad before placement. Check that adhesive is intact and undamaged. Do not use any damaged pads, replace pads if necessary.

- 3.) Apply the multifunction electrodes to patient, in accordance with your institution's protocol, by firmly rolling the electrode from top to bottom. Ensure that the total surface area of the electrode is in contact with the prepared skin.

WARNING: Always apply electrodes to flat areas of skin. If possible, avoid folds of skin such as those underneath the breast or those visible on obese individuals. Inadequate pad adhesion can lead to arcing or skin burns. Poor electrode pad-to-patient contact may result in defibrillator alarm, prompt or other indication. Check all electrical and patient connections from the device to the skin. If the indication continues, refer to the defibrillator's instruction for use or consult the defibrillator manufacturer.

- 4.) Connect electrodes to the defibrillator and follow defibrillator instructions or hospital protocol.

Apply a new pair of multifunction electrodes after every 24 hours on the skin, 25 defibrillations or 4 hours of pacing. 2603 Infant series MFEs are rated for maximum output of 100 J. 2602 Mini-Infant series MFEs are rated for a maximum of 30 J. Do not reposition electrodes. Replace with new electrodes if necessary. For best results during pacing, apply both standard ECG monitoring electrodes and PadPro® electrodes to the patient. Pacing therapy electrodes cannot be used to monitor ECG rhythm and deliver pacing current at the same time. During prolonged pacing greater than 30 mins, periodically examine patient's skin for irritation.

Disposal: After use, this product may be a potential biohazard. Handle and dispose of in accordance with hospital protocol and local, state and federal laws.

Notice: Any serious incident that has occurred using this device should be reported to CONMED Corporation and the competent authority of the Member State in which the user and/or patient is established.

Cautions:

- Refer to the defibrillator user's guide for overall operating procedures.
- Maintain at least 2.5 cm (1 in) separation between electrodes. Use electrodes on pediatric patients only if they fit completely on the torso and have at least 2.5 cm (1 in) of space between electrodes.
- Avoid electrode placement near the generator of an internal pacemaker. Defibrillation could cause pacemaker malfunction.
- The ability to effectively clean or sterilize this single use device has not been established and subsequent re-use may adversely affect the performance and safety of the device.
- Materials incorporated in this device could result in sensitization or allergic reaction to patient or end user.
- Federal law (USA) restricts this device to sale by or on the order of a physician.
- Discontinue noninvasive pacing if skin burn develops and another method of pacing is available.
- On cessation of pacing, immediately remove or replace electrodes.
- Do not use in the presence of a Magnetic Resonance Imaging (MRI) device.
- Do not open the multifunction electrode portion of the package until just prior to using the electrodes (pre-connect package only).
- Do not use if conductive polymer gel is dried out.
- Inspect multifunction electrode for damage, wear and loose connections prior to use. Discard if any damage or wear is found.
- Do not apply multifunction electrodes over broken skin.
- Do not apply any substance to the skin surface that will leave a residue (e.g. lotions, oils, etc.).
- Do not reapply multifunction electrodes if removed. Apply a new pair.
- Do not cut, bend, fold, crush, puncture, or modify the multifunction electrode pads. Do not cut or alter pads' size or shape.
- After application, keep the multifunction electrodes and surrounding area free of fluids.
- This product may be used by authorized medical personnel only.

	Operating Environmental Limitations		Defibrillation-Proof Type CF Applied Part		Not for Use in AED Mode
	Storage Environmental Limitations		Transportation Environmental Limitations		

Made in USA for:

CONMED

CONMED CORPORATION

525 French Road

Utica, New York 13502-5994 USA

(315)797-8375 • FAX 1-800-438-3051

Customer Service (USA) 1-866-426-6633

e-mail: CustomerExperience@conmed.com

EC REP

MDSS GmbH

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D-30175 Hannover, Germany

For International orders or inquiries, please contact CONMED International Sales

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© CONMED Corporation