
Nellcor™

SpO₂ Sensor

Multisite Reusable


en	Instructions For Use
fr	Mode d'emploi
de	Gebrauchsanweisung
nl	Gebruiksaanwijzing
it	Istruzioni per l'uso
es	Instrucciones de uso
sv	Bruksanvisning
da	Brugsanvisning
pt	Instruções de uso
ru	Инструкция по применению
zh	使用说明
pl	Instrukcja użytkowania
cs	Návod k použití
hu	Használati útmutató
el	Οδηγίες χρήσης
tr	Kullanma Talimatı


1


2

3

4

SpO₂ Sensor

Multisite Reusable

Identification of a substance that is contained or present within the product or packaging.

Identification of a substance that is not contained or present within the product or packaging.

IP22

Protected against access to hazardous parts with a finger and protected against vertically falling water drops when sensor tilted up to 15°.

This device is not made with natural rubber latex or DEHP.

Directions for Use

Indications/Contraindications

The Nellcor™ SpO₂ Sensor Multisite Reusable, model D-YS, is indicated for use when continuous noninvasive arterial oxygen saturation and pulse rate monitoring are required for patients weighing greater than 1 kg.

The SpO₂ Sensor is applied to the sensor site using Nellcor adhesive wraps, models ADH-A/N (adult/neonatal) and ADH-P/I (pediatric/infant), or Nellcor disposable foam wraps, models FOAM-A/N and FOAM-P/I. Nellcor adhesive and disposable foam wraps are intended for single patient use only. Additionally, the SpO₂ Sensor can be used with D-YSE (SpO₂ Ear Sensor) or D-YSPD (Pediatric SpO₂ Sensor Reusable).

	>1 – <3 kg	3 – 40 kg	>40 kg
ADH-A/N	◆		◆
FOAM-A/N	◆		◆
ADH-P/I		◆	
FOAM-P/I		◆	
D-YSE			◆
D-YSPD		◆	

The SpO₂ Sensor is contraindicated for use on patients who exhibit allergic reactions to the pressure-sensitive adhesive on the wraps.

Instructions for Use

Reusable sensors may be used on the same site for a maximum of 4 hours, provided the site is inspected routinely to ensure skin integrity, correct positioning, and adhesion of the sensor wrap. Because individual skin condition affects the ability of the skin to tolerate sensor placement, it may be necessary to change the sensor site more frequently with some patients.

To apply the SpO₂ Sensor:

1. Remove the small piece of paper backing from the sensor wrap. **(1)**
2. To secure the wrap to the sensor, push the “button” on the black sensor pad through the hole, as shown. **(2)**
3. Remove the remainder of the paper backing. Push the “button” on the white sensor pad through the remaining hole. **(3)**
4. The sensor and wrap are now ready to be applied to the patient. Select an appropriate site as described below. **(4)**

Note: When placing the SpO₂ Sensor on a finger or toe, the WHITE sensor pad must be placed on the NAIL SIDE.

Note: When selecting a sensor site, priority should be given to an extremity free of an arterial catheter, blood pressure cuff, or intravascular infusion line.

Neonate (1–3 kg): The preferred site is a foot, proximal to the toes with the cable positioned along the side of the foot. Alternatively, place the sensor around the palm of a hand, with the cable positioned along the side of the hand.

Infant (3–15 kg): The preferred site is around a great toe with the cable positioned along the side of the foot.

Pediatric (15–40 kg): The preferred site is around an index finger with the cable positioned along the medial side of the finger. Alternative sites are around a thumb, finger, or a great toe.

Adult (> 40 kg): The preferred site is around an index finger with the cable positioned along the medial side of the finger. Alternatively, use a thumb, a smaller finger, or a great toe.

5. Position the sensor so that the electronic components are directly opposite each other and with the sensor aligned on the side of the foot, toe, or finger so that the cable is positioned along the appropriate surface (as described in #4 above).
6. Wrap the remaining portion of the adhesive wrap around the site, loosely enough to ensure good circulation.
7. Plug the sensor into the instrument as described in the instrument operator’s manual, and verify proper operation.

Note: If the sensor does not track the pulse reliably, it may be incorrectly positioned—or the sensor site may be too thick, thin, or deeply pigmented, or otherwise deeply colored (for example, as a result of externally applied coloring such as nail polish, dye, or pigmented cream) to permit appropriate light transmission. If any of these situations occurs, reposition the sensor or choose an alternate Nellcor sensor for use on a different site.

Removal and Cleaning

Remove the disposable wrap from the D-YS by peeling it away from the sensor. The sensor may be surface-cleaned with a solution such as 70% isopropyl alcohol. If low-level disinfection is required, use a 1:10 bleach solution. Do not use undiluted bleach (5%~5.25% sodium hypochlorite) or any cleaning solution other than those recommended here because permanent damage to the sensor could occur.

CAUTION: Do not expose connector pins to cleaning solution as this may damage sensor.

To clean or disinfect the sensor using the recommended wipe method:

1. Saturate a clean, dry gauze pad with the cleaning solution. Wipe all surfaces of the sensor and the cable with this gauze pad.
2. Saturate another clean, dry gauze pad with sterile or distilled water. Wipe all surfaces of the sensor and cable with this gauze pad.
3. Dry the sensor and cable by wiping all surfaces with a clean, dry gauze pad.

To clean or disinfect the sensor using the recommended soak method:

1. Place the sensor in the cleaning solution, such that the sensor head(s) and desired length of cable are completely immersed.

CAUTION: Do not immerse the connector end of the cable as this may damage the sensor.

2. Dislodge air bubbles by gently shaking the sensor and cable.
3. Soak the sensor and the cable for 10 minutes.
4. Remove from cleaning solution.
5. Place the sensor and the cable in room temperature sterile or distilled water for 10 minutes.
6. Remove from the water.
7. Dry the sensor and cable by wiping all surfaces with a clean, dry gauze pad.

WARNINGS

1. **Do not use the SpO₂ Sensor or other oximetry sensors during MRI scanning. Conducted current may cause burns. Also, the sensor may affect the MRI image, and the MRI unit may affect the accuracy of oximetry measurements.**
2. **If the sensor is misapplied with excessive pressure for prolonged periods, a pressure injury can occur.**

CAUTIONS

1. Failure to apply the SpO₂ Sensor properly may cause incorrect measurements.
2. Using the SpO₂ Sensor in the presence of bright lights may result in inaccurate measurements. In such cases, cover the sensor with an opaque material.
3. Reusable sensors must be moved to a new site at least every 4 hours. Because individual skin condition affects the ability of the skin to tolerate sensor placement, it may be necessary to change the sensor site more frequently with some patients. If skin integrity changes, move the sensor to another site.

4. Use only Nellcor wraps designed for use with this sensor. Do not use tape. Use of additional tape or various other types of wraps can cause skin damage.
5. Intravascular dyes or externally applied coloring such as nail polish, dye, or pigmented cream may lead to inaccurate measurements.
6. Excessive motion may compromise performance. In such cases, try to keep the patient still, or change the sensor site to one with less motion.
7. When the sensor is wrapped too tightly or supplemental tape is applied, venous pulsations may lead to inaccurate saturation measurements.
8. As with all medical equipment, carefully route cables to reduce the possibility of patient entanglement or strangulation.
9. Using excessive force when removing the disposable wrap may damage the sensor.
10. Do not sterilize by irradiation, steam, or ethylene oxide. Such sterilization could damage the sensor.
11. Do not alter or modify the SpO₂ Sensor. Alterations or modifications may affect performance or accuracy.
12. For additional warnings, cautions or contraindications when using this sensor with Nellcor-compatible instruments, refer to the instrument operator's manual or contact the manufacturer of the instrument.

Note: High oxygen levels may predispose a premature infant to develop retinopathy. Therefore, the upper alarm limit for oxygen saturation must be carefully selected in accord with accepted clinical standards and considering the accuracy range of the oximeter being used.

If you have questions regarding any of this information, contact Covidien Technical Services Department, or your local Covidien representative.

Warranty

To obtain information about a warranty, if any, for this product, contact Covidien Technical Services or your local Covidien representative.

Accuracy Specifications

For the accuracy specification range when used with Nellcor™ monitors, refer to information provided with the monitor, or (in the U.S.), contact Covidien Technical Services Department. Outside the U.S., contact your local Covidien representative.

For the accuracy specification range of this sensor when used with Nellcor-compatible instruments, refer to the instrument operator's manual or contact the instrument manufacturer.

Use this sensor only with Nellcor instruments and instruments containing Nellcor oximetry, or with instruments licensed to use Nellcor sensors (Nellcor-compatible instruments). This sensor integrates Nellcor OxiMax™ technology into its design. When connected to an OxiMax-enabled instrument, this sensor uses OxiMax technology to provide additional

advanced sensor performance features. Consult individual manufacturers for features and compatibility of particular instruments and sensor models.

Each manufacturer of Nellcor™-compatible instruments is responsible for determining whether and under what conditions its instruments are compatible for safe and effective use with each Nellcor sensor model. This may include different specifications and/or warnings, cautions, or contraindications. Refer to the instrument operator's manual or consult manufacturer for complete instructions for use of this sensor with their Nellcor-compatible instrument.

Additional Copies of Instructions

Additional copies of these Instructions are available at no charge by calling Covidien or its authorized distributors. Also, permission is hereby granted under Covidien copyrights to Purchasers of products obtained from Covidien or its authorized distributors to make additional copies of these instructions for use by such purchasers.

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COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.
Covidien Inc, 15 Hampshire Street,
Covidien Ireland Limited,