

INSTRUCTION MANUAL

OPERATION & SERVICE

326M SERIES

**SUCTION APPARATUS, SURGICAL, AND
GASTROINTESTINAL ABDOMINAL DRAINAGE**

PORTABLE, AC/DC/RECHARGEABLE BATTERY

CONTRACT NO. SPO200-98-D-8703

TABLE OF CONTENTS

| SUBJECT | PAGE |
|--|------|
| LIST OF ILLUSTRATIONS | iii |
| SHIPPING CONTENTS | v |
| ACCESSORIES LIST | v |
| LIMITED COPYRIGHT RELEASE | v |
| CALIBRATION NOTICE | vi |
| UNPACKING | vi |
| LOCATION OF USE | vi |
| WARNINGS REGARDING USE | vi |
| ASSEMBLY, INTERCONNECTIONS AND INITIAL ADJUSTMENTS | vi |
| ASSEMBLY | |
| INTERCONNECTIONS | |
| INITIAL ADJUSTMENTS | |
| USAGE APPLICATIONS | viii |
| SECTION I. OPERATION | 1-1 |
| INTRODUCTION | 1-1 |
| OPERATION | 1-1 |
| DESCRIPTION OF CONTROLS, CIRCUIT BREAKER, CONNECTORS AND INDICATORS | |
| OPERATING POWER SELECTION & STOPPING | |
| SUCTIONING | |
| VACUUM REGULATOR (LIMITER) | |
| COLLECTION CANISTER | |
| OPERATOR PERFORMANCE CHECKS | |
| BATTERY CARE | 2-1 |
| ROUTINE CARE AND MAINTENANCE | 3-1 |
| CLEANING | |
| Component Removal | |
| Exterior Case | |
| Collection Canister | |
| Hydrophobic/Bacterial/Overflow Filter | |
| MAINTENANCE | |
| IN CASE OF DIFFICULTY | 4-1 |
| OPERATOR CORRECTIBLE PROBLEMS | |
| OPERATOR PROBLEMS REQUIRING SERVICE | |
| STORAGE INFORMATION | 5-1 |

TABLE OF CONTENTS (CONT'D)

| SUBJECT | PAGE |
|--|------|
| LIMITED WARRANTY | 5-1 |
| SPECIFICATIONS | 6-1 |
| SECTION II SERVICE | 7-1 |
| INTRODUCTION | 7-1 |
| CAUTIONARY NOTES | 7-1 |
| OPERATING VOLTAGES PRECAUTION | |
| INTERNAL RECHARGEABLE BATTERY | |
| HELPFUL HINTS | 7-1 |
| DISASSEMBLY/REASSEMBLY | 8-1 |
| BOTTOM COVER | |
| BATTERY COMPARTMENT, CLEAR FRONT COVER & TOP COVER | |
| SWITCHER PRINTED CIRCUIT BOARD | |
| REGULATOR PRINTED CIRCUIT BOARD | |
| MAIN PRINTED CIRCUIT BOARD | |
| CONNECTOR PANEL | |
| PUMP/MANIFOLD | |
| FRONT PANEL | |
| CALIBRATION PROCEDURE | 9-1 |
| REQUIRED EQUIPMENT | |
| PROCEDURES | |
| CIRCUIT DESCRIPTIONS | 10-1 |
| EXTERNAL AC POWER SUPPLY | |
| INTERNAL POWER SUPPLIES | |
| DC TO DC CONVERTER | |
| 6-VOLT POWER SUPPLY | |
| 5-VOLT & 25.6-VOLT POWER SUPPLIES | |
| ISOLATED 6-VOLT POWER SUPPLY | |
| PUMP CIRCUIT | |
| SOLENOID CIRCUITS | |
| TIMING CIRCUITS | |
| BATTERY CHARGER CIRCUIT | |
| BATTERY LEVEL MONITORING CIRCUIT | |
| PREVENTATIVE MAINTENANCE INSPECTIONS | 11-1 |
| VISUAL CHECKS | |
| PERFORMANCE CHECKS | |
| CLEANING | |
| TROUBLESHOOTING GUIDE | 12-1 |
| TECHNICAL DOCUMENTATION | 13-1 |
| BLOCK DIAGRAM | |
| BILLS OF MATERIAL | |
| WIRING DIAGRAM & ELECTRICAL SCHEMATICS | |
| WIRE LIST | |

LIST OF ILLUSTRATIONS

| FIGURE # | DESCRIPTION | PAGE |
|----------|--|------------------|
| 1. | Model 326/326M Main Features | iv |
| 2. | Interconnection Diagrams | vii |
| | Suction Apparatus shown with hydrophobic/bacterial/overflow filter | |
| 3. | Panel Controls, Circuit Breaker, Connectors and Indicators | 1-1 |
| 4. | Composite Illustration Depicting Major Sub-Assemblies | 13-1 |
| 5. | Bottom Cover Assembly | 13-2 |
| 6. | Top Cover Assembly | 13-3 |
| 7. | Connector Panel & Battery Pack Assemblies | 13-4 |
| 8. | Pump & Manifold Assemblies | 13-5 |
| 9. | Front Panel Assembly | 13-6 |
| 10. | AC Charger Rectifier Assembly | 13-7 |
| 11. | Collection Canister & Auto Power Cable Assemblies | 13-8 |
| 12. | Wire Harness Assembly | 13-9 |
| 13. | Accessory Kit-Air Force Assembly | 13-10 |
| 14. | Accessory Kit-Army Assembly | 13-11 |
| 15. | Switcher Printed Circuit Board Assembly | 13-12 |
| 16. | Regulator Printed Circuit Board Assembly | 13-13 |
| 17. | Main Printed Circuit Board Assembly | 13-14 |
| 18. | Block Diagram | 13-15 |
| | Wiring Diagram & Electrical Schematics (5 pages) | 13-37 thru 13-41 |




FIGURE 1 MODEL 326/326M MAIN FEATURES

SHIPPING CONTENTS

Each Model 326 is shipped with the following contents:

- | | |
|---|---|
| 1 ea. Apparatus, Suction, Portable | 1 ea. Hose, Clear, PVC, 1' Long |
| 1 ea. Assembly, Auto Power Cable | 1 ea. Hose, Clear, PVC, 18" Long |
| 1 ea. Suction Hose, Sterile, Clear, 6' Long | 1 ea. Hose, Clear, PVC, 2' Long |
| 1 ea. AC/DC Power Supply (115/230 VAC, 50/60/400 HZ input; nominal 12 VDC output) | 1 ea. Universal Canister Attachment Bracket |
| | 1 ea. Fuse, Spare |
| | 1 ea. Instruction Manual, Operation |

Each Model 326M is shipped with the following contents:

- | | |
|---|--|
| 1 ea. Apparatus, Suction, Portable | The following are provided with NSN 6515-01-435-0050 |
| 1 ea. Assembly, Auto Power Cable | |
| 1 ea. Suction Hose, Sterile, Clear, 6' Long | 2 ea. Collection Canister Assembly, Autoclavable |
| 2 ea. Fuse, Spare | 2 ea. Filter, Disposable, Hydrophobic/Bacterial/Overflow (for use with reusable canisters only) |
| 1 ea. AC/DC Power Supply (115/230 VAC, 50/60/400 HZ input; nominal 12 VDC output) | 1 ea. Carry-all Case with Label |
| 1 ea. Hose, Clear, PVC, 1' Long | The following are provided with Air Force Accessory Kit |
| 2 ea. Hose, Clear, PVC, 18" Long | |
| 2 ea. Hose, Clear, PVC, 2' Long | 2 ea. Disposable Collection Canister & Lid |
| 2 ea. Universal Canister Attachment Bracket | 1 ea. Catheter, 14 French |
| 2 ea. Instruction Manual, Operation & Service | 1 ea. Catheter, 18 French |
| 1 ea. Case, Padded, Aspirator | |
| 2 ea. Strap, Velcro® | |

ACCESSORIES LIST

The Accessories List contains common items, required from time to time. Each item is preceded by its part number. Accessories may be ordered direct from Impact. When ordering, please include the part number, description and quantity required.

Send written purchase orders to:

Telephone orders: 973/882-1212
FAX orders: 973/882-4993

Impact Instrumentation, Inc.
P.O. Box 508
27 Fairfield Place
West Caldwell, New Jersey 07006

| <u>PART NUMBER</u> | <u>DESCRIPTION</u> |
|--------------------|---|
| 465-0005-00 | Filter, Disposable, Hydrophobic/Bacterial/Overflow |
| 540-0068-00 | Hose, Clear, PVC, 1' Long |
| 540-0055-00 | Hose, Clear, PVC, 18" Long |
| 540-0051-00 | Hose, Clear, PVC, 2' Long |
| 703-0326-07 | Assembly, Collection Jar |
| 703-0326-06 | AC/DC Power Supply |
| 820-0018-00 | Tubing, Suction, Sterile, 9/32" I.D. X 6' |
| 906-0326-03 | Instruction Manual, Operation & Service, Model 326M |
| 906-0326-01 | Instruction Manual, Operation, Model 326 |
| 906-0326-04 | Instruction Manual, Operation & Service, Model 326 |

LIMITED COPYRIGHT RELEASE

Permission is hereby granted to the Department of Defense to reproduce all material furnished under this contract for use in a military service training program and other technical training programs.

CALIBRATION NOTICE

This device should be incorporated into a regular preventative maintenance program to insure compliance with operating specifications. Calibration measurements should be made on a biannual basis unless significant usage warrants a shorter period between preventative maintenance inspections. A calibration check should be made following each cumulative period of 300 hours of operation. Recommended maintenance checks can be found in the **SERVICE** section of this Manual.

UNPACKING

Check the contents of the shipping case(s) against the enclosed packing list. Examine the instrument for any obvious signs of shipping damage. If there is no apparent sign of mechanical damage, read the instructions contained within this manual before attempting to operate the instrument.

LOCATION OF USE

The Model 326/326M is a transportable device, therefore, its physical area of use will vary. When operated in a wet environment, user's should take precautions and protect this device by covering it with a protective barrier (small tarp, plastic sheet, etc.).

WARNINGS REGARDING USE

Caution: Federal law restricts this device to sale by or on the order of a physician.

This equipment is intended for use by qualified medical personnel or person(s) under the guidance and instruction of certified medical therapists.

Danger - Possible explosion hazard if used in the presence of flammable anesthetics.

Caution - Electric shock hazard, do not remove inside cover. Refer servicing to qualified biomedical equipment technicians only. See section entitled **SERVICE**.

Do not operate this instrument prior to reading the instructions contained within this manual.

Disposable hydrophobic/bacterial/overflow filter is for use only with reusable, autoclavable, collection canisters.

Do not clean collection canister with abrasive cleansers or (See **ROUTINE CARE AND MAINTENANCE** section **CLEANING**).

ASSEMBLY, INTERCONNECTIONS AND INITIAL ADJUSTMENTS

ASSEMBLY: No assembly is required before placing this device into operation.

INTERCONNECTIONS: Tubing interconnections are required before placing this device into operation. Figure 2 depicts the interconnections required to connect the Model 326/326M to the reusable collection canisters (in series) and disposable hydrophobic/bacterial/overflow filter.

INITIAL ADJUSTMENTS: Before placing this device into operation read the section entitled **OPERATION, DESCRIPTION OF CONTROLS, CIRCUIT BREAKER, CONNECTORS AND INDICATORS**. Make control settings and verify device performance prior to interfacing with patient.

ASSEMBLY, INTERCONNECTIONS AND INITIAL ADJUSTMENTS (cont'd)




FIGURE 2 INTERCONNECTION DIAGRAM

Suction Apparatus shown with disposable hydrophobic/bacterial/overflow filter.

USAGE APPLICATIONS - DEEP DRAINAGE/GASTRO-INTESTINAL SUCTION

IMPACT MODELS 306/306M AND 326/326M PORTABLE ASPIRATORS

Deep drainage/gastro-intestinal suction, frequently referred to as Wangensteen suction (after Dr. Owen H. Wangensteen*), requires a thorough understanding of "mild suction" procedures including collection bottle placement and tubing lengths. The following data is for reference only and specific applications may warrant deviations and/or additional considerations.

The Impact Models 306/306M and 326/326M are capable of maintaining precise vacuum and airflow levels. ON/OFF cycling times may be varied from patient to patient effectively optimizing drainage intervals. Cycling intervals are somewhat arbitrary. It is recommended that the user select intervals consistent with past practices and protocols.

The gastro-intestinal tract represents a pliable area which when subjected to continuous negative pressure will collapse around an implanted drainage tube. Such a collapse can cause continuous suction in a specific tissue area, thus preventing additional drainage while fluids accumulate and potential ulcerations can lead to perforations and hemorrhage at the collapsed site.

The intended action of intermittent suction is to provide a hydraulic effect whereby drainage fluids are pulled and pushed during ON and OFF cycles. The "pull" effect occurs as suction is applied to the gastro-intestinal tract, thus removing accumulated fluids. The "push" effect occurs when fluid in the elevated connecting tubing (18" maximum height) flows back into the drainage site, "washing" the implanted catheter away from the intestinal wall, thus relocating it in time for the next ON cycle.

It is important to note that the height of the collection canister and placement "dressing" of the connecting tubing each has a critical effect on overall performance. Under normal conditions, vacuum levels of 90 mmHg and 120 mmHg would raise a column of water approximately 4' and 5' respectively. It is possible that mild suction, involving particulate matter (solids), blood clots, and viscous secretions could present a situation whereby the delivered suction strength can only raise the drained contents to a height of roughly 3' to 4'. To compensate for this effect, the user is cautioned not to exceed an 18" height between collection bottle entry point and the distal end of the implanted suction catheter. The user should also be discouraged from using excessive lengths of connecting tubing. Since mild suction consists of minimal airflow, excessive tubing lengths could prevent fluid from actually entering the collection canister prior to the beginning of the next OFF cycle. The actual length of the ON cycle can be extended to compensate for this, however, the clinician must always recognize that an increased possibility of tissue occlusion exists.

The hydraulic action previously discussed is only possible when the collection canister is at a height greater than the distal end of the implanted catheter. During OFF cycles, the Models 306/306M and 326/326M vent to atmospheric pressure, enabling a gravitational backflow to occur. If the collecting jar is positioned below the height of the distal end of the implanted catheter, a siphoning effect can occur. Siphoning will act as a continuous suction presence between patient and collection canister which is undesirable. Cycling between the Models 306/306M and 326/326M and collection canister will continue despite siphoning. Clinicians may have visual objections to the appearance of drained fluids and seek options regarding collection canister placement. The factors involving the time cycles (ON and OFF), collection canister height and connecting tubing must be considered in order to have safe, effective intermittent suction.

*"INTESTINAL OBSTRUCTIONS" by Owen H. Wangensteen, M.D., C.C. Thomas, Publisher, Springfield, Illinois: Chapter VI.

SECTION I. OPERATION

INTRODUCTION

The Impact Model 326/326M is a self-contained, multi-purpose, suction apparatus designed for removing secretions from the upper airway during oropharyngeal, nasopharyngeal and tracheal suctioning procedures; programmable gastrointestinal and abdominal wound drainage; or surgical site debris in field hospitals. The apparatus includes EMI/RFI suppression circuitry and is suitable for use in desert, tropic, arctic or aeromedical environments. Simultaneous operation and battery recharge is permitted from either 115/230 VAC, 50-400 hz using the provided AC-DC Power Supply or 12 VDC using the provided Auto Power Cable accessories. Simultaneous operation and/or recharge from aircraft electrical systems (24 to 28 VDC) is permissible with an appropriate cable. The Model 326/326M is designed to accept DC power sources ranging from 11 to 30 volts, positive or negative ground.

The Model 326/326M is a completely self-contained suction source complete with vacuum pump, controls and DC power supply circuits; a dual-scale vacuum gauge with adjustable vacuum regulator (limiter); circuit breaker, collection canister mounting brackets, accessories and attachment tubing is provided. The complete device is housed within an injection molded case which includes a locking battery compartment door, a locking control panel door (clear), carrying handle and mounting interface for wall, bulkhead, pole or cart mounting.

This device is intended for use in non-explosive atmospheres. Read the instructions contained within this manual before attempting to operate this instrument.

OPERATION

DESCRIPTION OF CONTROLS, CIRCUIT BREAKER, CONNECTORS AND INDICATORS

Refer to the reference pictorial below. Numbers contained within this text (in parenthesis) correspond to the numbers indicated in the pictorial.

- (1) Circuit Breaker
- (2) Vacuum Inlet Fitting
- (3) External Power Input Jack
- (4) Dual-scale Vacuum Gauge
- (5) "ON TIME" Suction Interval Control
- (6) "OFF TIME" Suction Interval Control
- (7) Battery Charge Indicator
- (8) Charge Lamp
- (9) Power & Mode Selector Switch
- (10) Vacuum Regulator Control
- (11) External Power Lamp




Figure 3. Panel Controls, Circuit Breaker, Connectors and Indicators

OPERATION (cont'd)

DESCRIPTION OF CONTROLS, CIRCUIT BREAKER, CONNECTORS AND INDICATORS (cont'd)

- (1) Circuit Breaker - Protects pump motor from drawing excessive current.
- (2) Vacuum Inlet Fitting - Vacuum source connection to filter and collection canister(s).
- (3) External Power Input Jack - Connection for external operating and battery recharging power. Input accommodates 11 to 30 VDC input, positive or negative ground.
- (4) Dual-scale Vacuum Gauge - Displays vacuum developed within the patient circuit.
- (5) "ON TIME" Suction Interval Control - Setting determines how long intermittent suction will last.
- (6) "OFF TIME" Suction Interval Control - Setting determines how long before the next intermittent cycle begins.
- (7) Battery Charge Indicator - Displays battery charge status. Green area portrays the charged zone, the red area portrays the discharged zone.
- (8) Charge Lamp - Indicates the presence of battery charging current from an external power source.
- (9) Power & Mode Selector Switch - Turns operating power Off or On when the continuous "CONT" or intermittent "INT" suction mode is selected.

CONT - Selects continuous suctioning, vacuum adjustable from 0 to 550 mmHg

INT - Selects intermittent suctioning, vacuum adjustable from 0 to 200 mmHg.

Time interval combinations are selectable between 5 and 40 seconds ON, and 5 and 40 seconds OFF.

- (10) Vacuum Regulator (Limiter) Control - Limits the maximum deliverable vacuum level.
- (11) External Power Lamp - Illuminates when connected to a "live" AC mains.

OPERATING POWER SELECTION & STOPPING

The Model 326/326M is designed to operate from external power or internal rechargeable batteries. The Power & Mode Selector Switch (9) acts as a master power switch to start and stop operation.

SUCTIONING

1. Insure that all suction tubing is properly secured to respective fittings. Verify that collection canister lids are properly secured, fittings in place as shown (Figure 2), and no kinks in connecting tubing.
2. Charge Lamp (8) will illuminate during operation if external operating power is applied.
3. Turn Power & Mode Selector Switch (9) to the appropriate position: CONTinuous or INTermittent. The internal vacuum pump should begin operation.
4. Adjust Vacuum Regulator (Limiter) Control (10) to the maximum desired vacuum level by "pinching" and holding the vacuum tubing going to the collection canister. Deliverable vacuum levels will not exceed the preset maximum level. Adjusted and delivered vacuum levels are displayed on the Dual-scale Vacuum Gauge (4).

OPERATION, (cont'd)

SUCTIONING, (cont'd)

5. A disposable hydrophobic/bacterial/overflow filter is provided for use with Impact's standard autoclavable collection canisters. This filter connects between the Vacuum Inlet Fitting (2) and collection canister (Figure 2). This filter should be replaced when discoloration of its membrane occurs, the membrane contacts aspirate, or following 150 cumulative hours of use. This filter is designed to retain bacteria which would otherwise be exhausted into the immediate vicinity. The filter's retention of aerosolized aspirate will create more resistance to airflow resistance. It is important that the operator be aware of this and replace the filter accordingly.

6. Disposable collection canisters, traditionally include a built-in filter and may be used with your Model 326/326M. It is not necessary to use the disposable filter in conjunction with disposable canisters.

7. Reusable canisters should be used with the disposable filter to prevent inadvertent overflow, damage or contamination of the pump mechanism.

VACUUM REGULATOR (LIMITER)

The Vacuum Regulator (Limiter) (10) works in conjunction with the Vacuum Gauge (4). Vacuum levels may be selected by rotating the Vacuum Regulator (Limiter) (10): clockwise to increase vacuum; counterclockwise to decrease vacuum. Vacuum regulator adjustments should be made to the maximum desired vacuum level by "pinching" and holding the vacuum tubing going to the patient. Deliverable vacuum levels will not exceed this preset maximum. Adjusted and delivered vacuum levels will continuously display during operation on the dual-scale Vacuum Gauge (4).

COLLECTION CANISTER

Protect the suction mechanism from overflows which may permanently damage the vacuum pump if not properly and promptly cleaned. Use the disposable filter in conjunction with reusable collection canisters.

Collection canisters may be repeatedly sterilized via autoclave, ethylene dioxide gas or cold liquid sterilants (diluted in accordance with their respective instructions). **DO NOT** use abrasive cleansing agents. To determine compatibility with commercially available cleanser/disinfectants please note the following material content:

Collection Canister Bottle: Polysulfone/Polycarbonate
Collection Canister Cap: EPDM
Hose Fittings: Stainless Steel

OPERATOR PERFORMANCE CHECKS

Before placing this device into service, the operator can perform various checks to insure proper device performance.

1. Verify operating power selections from external power sources or internal rechargeable batteries.
2. Verify continuous operation, verify intermittent operation.
3. Test the Vacuum Regulator Limiter (10) for correct operation at various vacuum settings.
4. Insure that all hoses and fittings are properly connected.

BATTERY CARE

The Model 326/326M utilizes sealed GEL cell batteries which offer excellent charge retention characteristics, particularly during long periods of storage. This ensures an ample amount of power during emergencies and transitory procedures. The battery pack in this device is not intended for use as the primary power source, therefore, it should be used with discretion and its design understood. To provide long life and maximum performance capabilities, the Model 326/326M requires 16-hours to fully recharge its fully discharged batteries. Of course, the batteries are rarely discharged this much so the subsequent recharge time is usually less. GEL Cell batteries require little user care to provide optimum performance and life expectancy. Because their self-discharge rate is extremely low (approximately 1 1/2% per month), lengthy periods of disuse without replenishment charging is possible. If long-term disuse is common, it would be advisable to recharge the unit once every two months. This will insure that battery charge is maintained at 80% capacity or better. Continuous battery recharging from the AC mains is permissible but not required. Recharging from 115 or 230 VAC is accomplished using the AC/DC Power Supply accessory supplied with the Model 326/326M. Batteries can be recharged from an external 12VDC power source too. For external 12VDC operation or recharge, a connecting cable with automotive plug is furnished. The External Power Input Jack (3) is located on the connector panel adjacent to the Vacuum Inlet Fitting (2).

The External Power Input Jack (3) will accept voltages ranging between 11 and 30 VDC for device operation and/or battery recharging. An appropriate cable, to interface between the External Power Input Jack (3) and the power source is required.

The life of these batteries depends, to a great extent, upon the care they receive. Following these simple guidelines will prevent premature charge depletion and reduction of battery life.

1. **DO NOT** operate this instrument where the temperature range exceeds -60°C to 60°C (-76°F to 140°F).
2. **DO NOT** charge this instrument where the temperature range exceeds -20°C to 50°C (-4°F to 122°F).
3. **DO NOT** store this instrument with the batteries discharged. Always store in a charged condition.
4. For long-term storage, the optimum storage temperature range is 10°C to 30°C (50°F to 80°F).

ROUTINE CARE AND MAINTENANCE

CLEANING

Note: Routine decontaminations which do not involve the removal of aspirate can be effected using a spray disinfectant. With the device operating, simply spray a small amount of disinfectant directly into the collection canister and shortly thereafter into the Vacuum Inlet Fitting (2). This should be performed after each use to avoid risk of bacterial growth. In the event of an aspirate overflow, remove the pump head assembly (SECTION II. SERVICE, DISASSEMBLY/REASSEMBLY). All pump head components may be sterilized using a liquid disinfectant, a mild spray disinfectant or ethylene oxide gas. **DISCONNECT AC/DC POWER SUPPLY OR AUTO POWER CABLE FROM MODEL 326/326M PRIOR TO CLEANING.**

- **Component Removal**

The collection canister(s) and filter should be detached to facilitate cleaning.

- **Exterior Case**

Periodically or when applicable, clean the exterior case using a mild, non- abrasive cleanser. Disconnect Remove collection canister(s) and filter. **DO NOT** immerse or allow liquids to enter the case. A damp cloth will suffice in most instances. Disinfectant spraying is recommended at regular intervals. Allow to dry.

- **Collection Canister**

Impact's reusable collection canisters may be repeatedly sterilized via autoclave, ethylene oxide gas or cold liquid sterilants (diluted in accordance with their respective instructions). **Do not** use abrasive cleansing agents. To determine compatibility with commercially available cleanser/disinfectants please note the following material content:

Collection Canister Bottle: Polysulfone/Polycarbonate
Collection Canister Cap: EPDM
Hose Fittings: Stainless steel

1. Thoroughly clean collection canister(s), and fittings after each use. Handy hint: Before emptying collection canister, cap both hose fittings with a short length of tubing. This will prevent accidental spillage of aspirate.
2. Tubing is considered disposable and should be discarded following each use.
3. Insure that all parts are securely fastened and properly connected after cleaning.
4. Orient collection canister and route tubing as shown in Figure 2 following cleaning.
5. To prevent risk of cross contamination and the spread of airborne particulate matter, the use of bacterial filters is recommended. If filters are not used, this device should be disinfected and cleaned, following each use, as described earlier within this section.

Disposable filters connect between the collection canister and Vacuum Inlet Fitting (2). The filter contains a hose barb at each end to facilitate quick connections to the collection canister and Vacuum Inlet Fitting (2).
6. Disposable filters may be obtained from Impact. When ordering, specify Impact part number 465-0005-00. Filters may be used repeatedly until discolored or contact with aspirate and/or fluids occurs. As filters become occluded with particulate matter during repeated usages, a reduction in device airflow will become evident. Filter replacement will restore the device to its original airflow levels.

ROUTINE CARE AND MAINTENANCE (cont'd)

- **Bacterial/Overflow Filter**

Do not attempt to clean disposable bacterial/overflow filters. This item is disposable and should be replaced whenever it becomes discolored or contacts aspirate, airflow is impeded, or following 150 cumulative hours of use, whichever comes first.

Do not bypass this filter. Its intended use is to retain bacteria which would be expelled through the exhaust port or allowed to accumulate in the pump head. As filters become occluded with particulate matter during repeated usages, a reduction in device airflow will become evident. Filter replacement will restore device performance to its normal airflow levels.

MAINTENANCE

Routine maintenance should be performed on this apparatus at regular intervals and prior to its being placed into service. Routine maintenance should consist of the following:

1. Cleaning checks - as described above.
2. Filter checks - replace when discolored, contact with aspirate occurs, airflow performance diminishes considerably or following 150 hours of cumulative use.
3. Operational checks - as described in **OPERATOR PERFORMANCE CHECKS**.
4. Tubing checks - replace crimped, cracked or worn tubing as required.

IN CASE OF DIFFICULTY

Authorization to service this instrument by other than factory-trained or certified personnel will not be given, nor does Impact Instrumentation, Inc. assume any responsibility and/or liability resulting from such unauthorized servicing.

Impact will, upon request, provide competent biomedical engineering departments with service data and schematics. Such departments are encouraged to contact the factory for assistance when needed and it is recommended that staff members attend a factory training course. Details may be obtained by contacting the Impact Customer Service Department.

OPERATOR CORRECTIBLE PROBLEMS

Common problems may be quickly rectified by users. Should the Model 326/326M fail to operate properly, verify the integrity of all tube connections, tubing, fittings, and control settings. One can quickly isolate problems to an accessory item or the suction apparatus by testing for vacuum at various locations.

To isolate a problem, check for vacuum at the inlet of each item, tracing backwards through the system, i.e.: vacuum from the collection canister to the Vacuum Inlet Fitting; or if a filter is used: vacuum from the collection canister to the filter, then vacuum from the filter to the Vacuum Inlet Fitting.

OPERATOR PROBLEMS REQUIRING SERVICE

If the tests described above do not resolve an operating problem, service is required. Should servicing be necessary, contact your nearest Impact representative or the Impact Customer Service Department (973) 882-1212.

Please have the Model and Serial Numbers ready and any other pertinent data you wish to include in your service request. The Model 326/326M Serial Number is located on the outer case identification label.

STORAGE INFORMATION

For prolonged storage periods, the Model 326/326M should be stored indoors. The environment should be clean, and out of direct sunlight. Storage temperatures should range between 5°F and 104°F (-15°C to 40°C), humidity should be low.

When batteries are in extended storage, it is recommended that they receive a refresh charge at recommended intervals:

| <u>STORAGE AMBIENT</u> | <u>RECHARGE INTERVAL</u> |
|----------------------------|--------------------------|
| Below 68°F (20°C) | 18 months |
| 68° to 86°F (20° to 30°C) | 12 months |
| 86° to 104°F (30° to 40°C) | 6 months |

Following periods of extended storage in non-controlled environments, allow the Model 326/326M sufficient time to stabilize to a temperature within its specified operating range (see **BATTERY CARE**).

LIMITED WARRANTY

Impact Instrumentation, Inc. warrants this instrument to be free from all defects in materials and workmanship for a period of one (1) year. Batteries, which by their nature are consumable and subjected to environmental extremes, will be warranted only for defects of manufacturing origin for a period of ninety (90) days. Disposable accessories, consumable in usage, will be warranted only for defects of manufacturing origin prior to their initial use. This warranty is neither assignable nor transferable, nor does it apply if this instrument is tampered with, misused or serviced by unauthorized personnel. All warranty repairs shall be subject to return postage billing.

SPECIFICATIONS

| | |
|------------------------------|--|
| VACUUM RANGE (Continuous): | 0-550 mm/Hg (0-22 inches/mercury) |
| (Intermittent): | 0-200 mmHg (0-8 inches/mercury) |
| FREE AIRFLOW: | 30 Liters Per Minute (LPM) |
| TEMPERATURE OPERATING RANGE: | -20°C to 49°C (-4°F to 120°F) |
| CONTROLS: | Power & Mode Selector Switch (Power OFF/CONTinuous/INTermittent Vacuum Regulator (Limiter) "ON TIME" Suction Interval Control (5-40 seconds, continuously variable) "OFF TIME" Suction Interval Control (5-40 seconds, continuously variable) Power Source Select Switch (115/230 VAC, located on AC/DC Power Supply) |
| DISPLAYS: | Charge Indicator Lamp External Power Lamp Battery Charge Indicator (Meter) Vacuum Gauge, Dual-Scale (Metric/English), 2 1/2" diameter. |
| CONNECTORS: | Vacuum Inlet Fitting External Power Input Jack |
| PUMP PROTECTION: | Circuit Breaker |
| POWER: | |
| External: | |
| 115/230 VAC, 50-400 Hz: | Continuous operation and/or recharge with supplied AC/DC Power Supply |
| 12 VDC: | Continuous operation and/or recharge with supplied Auto Power Cable |
| 11-30 VDC: | Continuous operation and/or recharge (cable not supplied) |
| Internal: | |
| Battery Pack, 12 VDC: | Operating Time, 2-hours (minimum, when cycled as specified) |
| Recharge Time: | 16-hours (maximum) |
| COLLECTION CAPACITY: | 1200 ml X 2 (Air Force) 1100 ml X 2 (Model 326M NSN: 6515-01-435-0050 only) |
| WARRANTY: | One (1) year, limited, as specified |
| CASE: | |
| Material: | Polycarbonate, color-through, injection-molded |
| Size: | 24.1 cm W x 29.2 cm H x 12.4cm D (9.5" W x 11.5" H x 4.87" D) |
| Weight: | 5.5 kg (12 lbs) |

SPECIFICATIONS, (cont'd)

MODEL 326 (Included Accessories):

- | | |
|---|---|
| 1 ea. Apparatus, Suction, Portable | 1 ea. Assembly, Auto Power Cable, 6' Long |
| 1 ea. Suction Hose, Sterile, Clear, 6' Long | 1 ea. Hose, Clear, PVC, 1' Long |
| 1 ea. Hose, Clear, PVC, 18" Long | 1 ea. Hose, Clear, PVC, 2' Long |
| 1 ea. Universal Canister Attachment Bracket | 1 ea. Fuse, Spare |
| 1 ea. AC/DC Power Supply | 1 ea. Padded Case, Aspirator |
| 1 ea. Instruction Manual, Operation | |

MODEL 326M (Included Accessories):

- | | |
|---|---|
| 1 ea. Apparatus, Suction, Portable | 1 ea. Assembly, Auto Power Cable, 6' Long |
| 1 ea. Suction Hose, Sterile, Clear, 6' Long | 2 ea. Fuse, Spare |
| 1 ea. AC/DC Power Supply | 1 ea. Hose, Clear, PVC, 1' Long |
| 2 ea. Hose, Clear, PVC, 18" Long | 2 ea. Hose, Clear, PVC, 2' Long |
| 2 ea. Universal Canister Attachment Bracket | 1 ea. Padded Case, Aspirator |
| 2 ea. Strap, Velcro® | 2 ea. Instruction Manual, Operation & Service |

The following are provided with NSN 6515-01-435-0050

- 2 ea. Collection Canister Assembly, Autoclavable
- 2 ea. Filter, Disposable, Hydrophobic/Bacterial/Overflow
- 1 ea. Carry-all Case with Label

The following are provided with Air Force Accessory Kit

- 2 ea. Disposable Collection Canister and Lid
- 1 ea. Catheter, 14 French
- 1 ea. Catheter, 18 French

SPECIFICATIONS CONTAINED HEREIN REPRESENT TYPICAL DEVICE PERFORMANCE

SECTION II. SERVICE

INTRODUCTION

The information contained herein is intended only for use by factory-trained, and certified personnel or military personnel trained in the care and servicing of this product. The manufacturer does not authorize or assume any obligations resulting from unauthorized servicing nor will it be held liable for any injuries or damages incurred therefrom.

Impact Instrumentation will provide service training at the manufacturing site at no schooling charge to users; however, travel and meal costs resulting therefrom shall be borne by the user. Impact will perform training at the customer's site if requested. The customer will be responsible for costs incurred by Impact personnel for travel, meal, and time - at prevailing rates.

The Impact service facility encourages dialogue from user service personnel towards rectifying any service related matter. All service requests may be addressed to the Service Manager, Impact Instrumentation, Inc., 27 Fairfield Place, West Caldwell, New Jersey 07006; via telephone 973/882-1212; or fax 973/882-4993.

CAUTIONARY NOTES

OPERATING VOLTAGES PRECAUTION

Prior to servicing this device, be aware of the presence of potentially dangerous operating voltages.

INTERNAL RECHARGEABLE BATTERY

Military contracts may require this product to be shipped with its battery pack contained within a separate carton. Prior to placing this device into operation, insure that its battery pack is installed and recharged. DO NOT operate this unit until its battery pack is installed.

HELPFUL HINTS

Before attempting to repair/calibrate this instrument, please take a few moments to insure that the problem is not accessory related.

Check the integrity of all vacuum hoses and tubing. Verify that the tubing has no crimps or cuts in it.

Insure that collection canisters seal properly, overflow shut-off valves do not stick, and bacterial filters (if used) are not clogged.

Refer to the enclosed schematic, assembly drawing, and pictorials when electrically troubleshooting. Isolate each problem to a functional segment of the circuitry. Verify the integrity of circuit ground and the presence of correct power supply voltages.

ALWAYS safeguard your personal well being when troubleshooting electronic circuitry. Remove jewelry, such as rings and bracelets, and keep liquids away from the vicinity of live circuitry.

DISASSEMBLY/REASSEMBLY

REQUIRED TOOLS

Screwdriver, slotted, 6" to 8" long, small tip
Screwdriver, phillips head, small #1
Screwdriver, phillips head, medium #2
10mm socket wrench
1/4" socket with drive handle
5/16" socket with drive handle
3/8" socket with drive handle
Pliers, needle nosed
Open-end wrench, adjustable, 10"
Open-end wrench set, 1/4" to 5/8"
Bench vise with smooth jaws
Diagonal Cutter
Cable Tie, miniature
3-32 Allen Key
Retaining Ring Extractor
Teflon Tape
Ruler, 6"
Krytox® Fluorinated Grease

DISSASSEMBLY/REASSEMBLY

BOTTOM COVER

Verify that Mode Selector Switch is in OFF position. Disconnect all hoses and cables. Open Battery Compartment Door, disconnect Battery Pack connector and release Velcro\ holddown strap. Remove Battery Pack.

Unscrew two (2) 6-32 Phillips Bind Head screws and washers (located inside Battery Compartment) securing Bottom and Top Covers. Turn unit over so that its Bottom Cover is facing upward. Unscrew five (5) Phillips Bind Head screws securing Bottom and Top Covers. **DO NOT** unscrew the one (1) 6-32 Phillips Bind Head screw that is located on bottom of Bottom Cover adjacent to the Condensed Operating Instructions label.

Connector Panel secures to Bottom Cover with four (4) 4-40 Phillips Bind Head screws. Remove screws. Bottom Cover can now be lifted from Top Cover Assembly.

To reassemble, reverse above process. Insure that Battery Pack connector wires are not crimped between components.

BATTERY COMPARTMENT, CLEAR FRONT COVER AND TOP COVER

Perform following instructions only if servicing Battery Compartment Separator, Battery Compartment Door, Clear Front Cover or Top Cover. Otherwise proceed to next step.

With Diagonal Cutter, cut miniature cable tie that secures battery connector wires to Battery Compartment Separator. **DO NOT** cut any other cable tie. Using 1/4" socket with drive handle, remove three (3) 4-40 keps nuts and three (3) #6 flat washers that secure Battery Compartment Separator, Battery Compartment Door and Clear Front Cover. Remove Battery Compartment Separator, Battery Compartment Door and Clear Front Cover. Remove remaining three (3) #6 flat washers and three (3) 4-40 keps nuts that secure Front Panel to Top Cover.

Pull bottom of Top Cover upwards until it is inbetween the Regulator and Switcher Printed Circuit Board (**NO HIGHER**). Pull top of Top Cover above Manifold Assembly. Remove Top Cover.

To reassemble, reverse this process exactly. Replace cable tie. Insure that Battery Pack connector wires are not crimped by Top Cover and Front Panel.

SWITCHER PRINTED CIRCUIT BOARD

Unscrew four (4) 4-40 screws and four (4) #4 lock washers securing Switcher Printed Circuit Board to Regulator Printed Circuit Board. Unplug 2-Pin connector. Lift Switcher Printed Circuit Board to separate its header from Regulator Printed Circuit Board.

Follow these instructions in reverse to reassemble.

REGULATOR PRINTED CIRCUIT BOARD

Unscrew four (4) 1/4" hex, male-to-female 4-40 spacers securing Regulator Printed Circuit Board. Unplug 5-pin connector from Regulator Printed Circuit Board. Lift Regulator Printed Circuit Board to separate its header from Main Printed Circuit Board.

Follow these instructions in reverse to reassemble. Make sure that header pins and their respective sockets align.

MAIN PRINTED CIRCUIT BOARD

Unscrew four (4) hex, female-to-female 4-40 spacers and one (1) 4-40 keps nut that secures Main Printed Circuit Board to rear of Front Panel. Carefully turn unit over, with its clear door opened. Remove two (2) small collet knobs by removing their caps and loosening their locking screws. Under each collet knob skirt is a nut. Carefully remove each nut. Separate cable connectors from Main Printed Circuit Board. Remove Mode Selector Switch knob using 10mm Socket Wrench. Carefully remove nut and washer with an open-end wrench. Unscrew the Charge Lamp and External Power Lamp. Carefully remove each nut and washer with an open-end wrench. Main Printed Circuit Board can now be removed.

To reassemble, reverse above steps. Be careful not to scratch Front Panel. When reinstalling collet knobs, insure that each control is in its full counterclockwise position - this will permit proper knob pointer alignment. Make sure that ribbon cable connector is properly aligned and attached.

CONNECTOR PANEL

Connector Panel is attached to Bottom Cover. See Bottom Cover disassembly for initial disassembly instructions.

Remove nut securing Circuit Breaker to Connector Panel. Remove Circuit Breaker. Remove two (2) 4-40 screws, two (2) 4-40 keps nuts and two (2) #4 fiber washers securing the External Power Jack and flange to Connector Panel. Remove External Power Jack and flange. Remove tubing to Suction Port.

Reverse above steps to reassemble.

PUMP/MANIFOLD

The Pump/Manifold Assembly should not be disassembled unless it is absolutely necessary. Its alignment and correct positioning of the Vacuum Adjust Knob are important. Under normal conditions treat the Pump/Manifold and Front Panel as one assembly.

PUMP/MANIFOLD/PUMP HEAD

IMPORTANT: Note the orientation of all components prior to disassembly/reassembly. After reassembly, test this device to insure proper operation before returning it for patient use.

Remove 1/8" ID PVC tubing from black swivel barb at end of manifold, and 3/8" ID corrugated hose from 3/8 ID hosebarb attached to manifold. Hold 3/8" hose by cuff when removing, pulling on corrugated section may damage hose.

Remove retaining ring from Plug Valve in manifold using a retaining ring extractor. Remove Plug Valve by carefully pulling and twisting from front side of panel.

PUMP/MANIFOLD/PUMP HEAD (cont'd)

Using nut driver, remove four (4) 6-32 keps nuts securing pump brackets. Lift pump slightly for access to both solenoid valves screws. Using the 3/32 allen key, remove both sets of screws (two per valve). Pull solenoids from manifold. The square-shaped valve has a rubber gasket. The rectangular valve has a metal rod, a circular rubber seal, and a metal disk under the seal. Save these parts for reassembly with the same solenoid and manifold; they are a matched set.

With pump still soldered to wiring, carefully lift pump off of 6-32 stud. Do not pinch wiring attached to pump or solenoid valves. Rotate pump to have access to bottom of pump head. Remove four (4) 8-32 x 3/4 phillips head screws and lockwashers securing pump head to pump body. Pull pump head and manifold from the pump body. Pump head may stick, slightly, to pump body, this is normal

Note: The distance between the manifold and the beginning of the pump head should be 0.150 +/- .020 inches. This will insure proper fit during reassembly.

To remove manifold from pump head, use hands to rotate pump head counterclockwise while holding manifold steady. Use of wrench(es) or a vise should not be necessary. Remove two (2) small phillips flat head screws securing valve plate to head plate. Separate plates to access die cut valve.

To reassemble: Follow the disassembly instructions in reverse. Observe the following notes and cautions.

Make sure die cut valve is seated properly before securing between valve plate and head plate.

Use teflon tape on 1/8" NPT close nipple (no more than 2 1/2 wraps). Use hands when rotating pump head onto manifold. Tighten slowly until you reach the 0.150 +/- .020 inches between pump head and manifold, and bottom of pump head is even with manifold bottom. **CAUTION:** Do not over tighten, pump head could crack. Do not tighten with wrench or secure in vise.

Make sure pump diaphragm is top, dead center, and that it is laying properly on pump before reattaching pump head.

After reattaching solenoids, position pump in its proper setting and slide pump with brackets onto four (4) 6-32 studs. Align stop pin in manifold with notch in vacuum adjust hole on front panel. Manifold and pump head should lay flat on underside of front panel. Check o-rings on plug valve for damage, replace if necessary. Lubricate plug valve o-rings with Krytox®. Carefully insert plug valve, making sure it is in completely. Check for proper 90° rotation.

Note: you may have to loosen the 6-32 X 5/16 phillips head screws securing brackets to pump to align pump/manifold. Tighten screws and four (4) 6-32 keps nuts holding brackets to front panel.

Reattach retaining ring and hoses to complete reassembly.

FRONT PANEL

Unscrew six (6) 4-40 keps nuts. The Front Panel can now be removed from the Top Cover.

Reverse above steps to reassemble.

CALIBRATION PROCEDURE

REQUIRED EQUIPMENT

- A. Oscilloscope, DC, Triggered, with a minimum 5-second horizontal sweep, storage capability desirable.
- B. Small, slotted screwdriver.
- C. Stopwatch, minimum 0.1 second resolution.
- D. Regulated DC-Power Supply (5-Ampere minimum rated output)

PROCEDURES

A. MAXIMUM LOW VACUUM LEVEL LIMIT

- 1. Set controls for either AC or DC operation (insure that batteries have been fully charged if calibrating from battery power).
- 2. Select INTERmittent operation.
- 3. Turn VACUUM ADJUST control to MAXimum (fully clockwise).
- 4. Occlude VACUUM INLET fitting and adjust R20 until Front Panel vacuum gauge reads 200 mmHg.
- 5. This reading may be verified using a calibrated vacuum gauge applied to the vacuum inlet.

B. INTERMITTENT ON/OFF TIMING CIRCUITS

NOTE: An oscilloscope may be used to calibrate ON/OFF Timing Circuits as described below. However, technicians may find it easier to calibrate this unit using a stopwatch. When using a stopwatch, utilize pump turn on and turn off as timing reference points.

- 1. Set controls for either AC or DC operation (insure that batteries have been fully charged if calibrating from battery power).
- 2. Select INTERmittent operation.
- 3. Set ON/OFF times for 5-seconds ON, 5-seconds OFF.
- 4. Turn VACUUM ADJUST control to MAXimum (fully clockwise).
- 5. Trigger oscilloscope sweep to begin when motor turns on and adjust R7 to set ON Time circuit for a 5-second sweep. Close verification can be made using a storage type oscilloscope. Adjustments should be to within +/- 0.5 seconds.
- 6. Trigger oscilloscope to begin when motor turns off and adjust R12 to set OFF Time circuit for a 5-second sweep. Adjustments should be to within +/- 0.5 seconds. Again, the use of a storage type oscilloscope will simplify measurement.

NOTE: Steps 5 and 6 above may be monitored at various points. For simplicity and convenience, the positive voltage motor input should be used and the oscilloscope triggering slope set for the ON circuit, then reset for OFF circuit.

C. LOW BATTERY INDICATOR ADJUSTMENT

- 1. Disconnect and remove battery pack from Aspirator.
- 2. Preset external power supply for 12 VDC and attach to Aspirator battery connector.
- 3. Turn external power supply on, set Aspirator Mode Select Switch to CONTinuous.
- 4. Adjust external power supply for 11.1-volts, and adjust R36 until the Low Battery Indicator needle is centered directly between the red and green divisions.
- 5. Disconnect the external power supply, reconnect battery pack.

CIRCUIT DESCRIPTIONS

Refer to enclosed schematic diagrams.

EXTERNAL AC POWER SUPPLY - Components P1, F1, VS1, T1, BR1, and C20 represent a full wave bridge rectifier circuit which enables simultaneous instrument operation and battery recharging. All components except C20 are located in an external enclosure, while C20 is mounted on the Connector Panel inside the instrument. VS1 allows selection of either nominal 115 VAC or 230 VAC input power.

INTERNAL POWER SUPPLIES - The instruments' internal power supply consists of five separate sections: a DC to DC Converter section which provides pre-regulation for externally applied power, plus separate 6-volt, 5-volt, 25.6-volt and isolated 6-volt sections. Connector J1 incorporates a switch that switches from battery power to external power whenever an external power plug is inserted into the jack. If no power is available at the plug when it is inserted into the power jack, the unit will not function even if a charged battery is installed. For the instrument to operate, the external power supply must be energized with the plug inserted, or the plug must be removed to allow battery operation.

DC TO DC CONVERTER - The DC to DC Converter consists of FLT1, U7, Rtrim, D4, D5 and LP2. Its function is to allow a wide external input voltage range of 11 to 30 VDC, while maintaining a stable 13.5 VDC output. This section of the power supply is energized whenever a valid external DC power source is connected. The output of this section also has two blocking diodes, D4 and D5, that prevents reverse voltage from being applied to the converter during battery operation.

6-VOLT POWER SUPPLY - The 6-volt Power Supply consists of C8 thru C11, and U4. It powers the timing circuitry, and is only operational for the INTERmittent mode. Switch SW1C turns on this section when in the INTERmittent position.

5-VOLT AND 25.6-VOLT POWER SUPPLIES - The 25.6-volt Power Supply powers the battery charge monitoring circuit. The higher voltage is required to prevent the amplifiers from operating at the power rail. The 5-volt Supply provides a stable reference for the comparator circuit that controls the charge indicating lamp. The 5-volt Power Supply consists of C21, C22 and U8, and the 25.6-volt supply of U9, C23, C24, D6 and D7. The 5-volt supply is a simple fixed three terminal regulator. The 25.6-volt supply consists of a voltage doubling circuit built around the ICL7662 integrated circuit. Both power supplies are active whenever external power is present.

ISOLATED 6-VOLT POWER SUPPLY - The isolated 6-volt supply consists of C33, C34, and U13. It is utilized in the battery monitoring circuit, and is active in the CONTInuous and INTERmittent operating modes.

PUMP CIRCUIT - The pump circuit consists of R17, Q4, Q5, C16 thru C19, U6, R19, R20, SW1B, CB1, FLT2, D3 and M1. When set to CONTInuous mode, switch SW1B supplies full system power to the pump. Circuit breaker CB1 provides resettable overload protection, while D3 suppresses reverse EMF when the pump's motor is turned off. Filter FLT2 is required for EMI suppression. In the INTERmittent mode, R17, Q4 and Q5 turn on power to U6 whenever the timing circuit is in an "ON" cycle. U6 is a 3-terminal adjustable regulator. R20 allows adjustment of the voltage applied to the pump, and when properly calibrated, will limit the maximum pump output vacuum to 200mmHg.

SOLENOID CIRCUITS - There are two solenoids with respective circuitry. Both are active only in the INTERmittent mode. The bleed solenoid circuit consists of S1, Q2 and D1. The dump solenoid circuit consists of S2, Q3, C12 thru C14, R15, R16, R18, U5 and D2. When the timing circuit is in the "ON" cycle, the bleed solenoid is activated to prevent the pump from stalling when flow is restricted. The same signal that turns on the pump also turns on Q2 which drives solenoid S1. D1 is a suppression diode that suppresses reverse EMF when the solenoid turns off. The dump solenoid activates at the beginning of an "OFF" cycle to bring the vacuum collection system to atmospheric pressure. It turns off after approximately 3-seconds to conserve power during battery operation. U5 is a 555-timer circuit that is connected to operate as a monostable multivibrator, which means that it will provide a single output pulse when it is triggered. R15 and C12 provide a trigger input pulse to the

device. The 555 is negative edge triggered, so it will trigger only when Q1 turns off. Normally, the output of the timer is off (0-volts). When the device is triggered, its output goes high (to Vcc = 6-volts), and stays high for the programmed time duration. R16 and C14 set the time duration that the output stays on. D2, like D1, is a suppression diode.

TIMING CIRCUITS - The timing circuitry consists of a master clock, two timers and interconnecting circuitry. The clock is a 50 Hz, self-starting square wave generator consisting of U1C, U1D, R1, R2, and C1. The "ON" timer controls how long the pump stays on during an "ON" cycle, and consists of U2, R3 thru R7, C3, C4, and C35. U2 is a 555-timer circuit that is connected to operate as a monostable multivibrator, which means that it will provide a single output pulse when it is triggered. R3 and C2 are driven by buffer U1A, and provide a trigger input pulse to the 555. C4, R4, and R5 set the time period that the timer output is on after it is triggered. R5 is brought out on the instrument front panel to allow user adjustment of the "ON" time. C3, R6, and R7 are used to adjust the delay multiplier in the time period equation $T = RC$. The output of the "ON" timer triggers the "OFF" timer by means of R8 and C5. The "OFF" timer functions essentially the same as the "ON" timer. It consists of U3, R8 thru R12, C6, C7 and C36. The output of the "OFF" timer drives buffer U1B, and gates on buffer U1A whenever its output is low (off). The "ON" and "OFF" timers are slaved to each other, so that the "ON" timer triggers the "OFF" timer whenever the former times out, and the "OFF" timer gates the master clock to only trigger the "ON" timer when the "OFF" timer has timed out. U1B and Q1 act as an output buffer for the timing circuitry, while R13 is a pulldown resistor for the CMOS transistors Q2 and Q4.

BATTERY CHARGER CIRCUIT - The battery charging circuit consists of R21 thru R26, U10, U11A, C25, Q6, and LP1. R21 is a 1-ohm resistor that provides current limiting and also acts as the current sense resistor for the battery charge monitor. The battery charger is a simple charger that provides higher initial charging current when a battery is depleted, and later tapers downward to about 50-75ma as the battery replenishes. U10 is a differential amplifier that amplifies and ground references the small voltage that appears across R21 as the battery is charging. R22 is the gain resistor for U10. U11A is connected as a comparator, and drives charging lamp LP1 through R26 and Q6. R24, R25 and C25 provide a stable reference voltage for the comparator. Since the signal being measured across R21 is essentially at system supply potential, a higher supply (25.6V) powers U10 and U11 to prevent the amplifiers from operating at the power supply rail. The battery charging and monitoring circuit is active whenever an external power source is present.

BATTERY LEVEL MONITORING CIRCUIT - The battery level monitoring circuit consists of SW1A, R27 thru R36, U12, and Meter1. It is powered by the actual battery voltage. Switch SW1A powers this circuitry in both INTERmittent and CONTInuous mode. The nominal battery voltage being monitored is 12-volts. To allow the circuit to operate with good sensitivity at an operating point away from the power rail, the isolated 6-volt supply is used to supply a virtual ground for this circuit. R27 and R29 divide the battery voltage in half. R31 and R32 provide an offset signal. Both signals are fed into summing amplifier U12A via R28 and R30. R33 provides gain for this stage in conjunction with R28 and R30. U12B is an inverting amplifier with a gain of 1. R34 and R35 are part of this stage. R36 can be trimmed to center the meter with a battery voltage of 11.1V being present.

***** NOTE *****

Waveforms and voltage measurements have been noted at various locations on each schematic drawing. In most cases, considerable leeway has been given as to what constitutes an acceptable voltage value in order to maintain device performance over a broad range of conditions.

PREVENTATIVE MAINTENANCE INSPECTIONS

Preventative maintenance inspections should be incorporated on a routine basis to insure proper device performance. These inspections should consist of both visual and performance checks, and cleaning when warranted.

Preventative maintenance inspections (PMI) should be made as follows:

If monthly usage is less than 50 hours - PMI bimonthly.

If monthly usage is greater than 100 hours - PMI monthly.

VISUAL CHECKS: Visual checks should include, but not be limited to:

1. Exterior case housing. Examine exterior of unit for cleanliness and general physical condition. Be sure that housing is intact, that all assembly hardware is tight, and that there are no signs of spilled liquids or other abuse.
2. AC Power Supply Plugs. Examine incoming AC plug for visible signs of damage. Verify that plug blades are secure. Examine the output power plug for visible signs of damage.
3. Input and output line cords. Inspect cords for signs of damage. If evidence of damage exists, replace entire cord.
NOTE: This device may contain shielded line cords to meet EMI/RFI criteria. Contact Impact Customer Service Department before replacing any defective line cord.
4. Strain reliefs. Examine strain reliefs at both ends of each line cord. Each line cord should be held securely by its respective strain relief.
5. Circuit breaker. This device has a switch type circuit breaker. Verify the free movement of its "red" stem.
6. Inspect all tubing to insure that they are not cracked, crimped or dirty. Replace defective tubing immediately.
7. Inspect all fittings and connectors. Replace fittings and connectors that exhibit evidence of wear or suspected leaks.
8. Controls and switches. Examine all controls and switches. Check their positions to verify correct alignment. Record the settings of all controls and switches that must be returned to their original positions following the inspection. Check for physical damage, loose mounting and inappropriate motion. If the control has fixed-limit stops, check for proper alignment and positive stopping. During this inspection, insure that each control and switch performs its proper operation.

Verify correct operation of the vacuum adjustment control and that it is usable over entire suction range indicated in device specifications. Perform this test with vacuum port occluded.
9. Vacuum pump/motor. Confirm proper operation of vacuum pump/motor. Verify deliverability of maximum vacuum as indicated in device specifications.
10. Battery Pack and AC Power Supply. Inspect physical condition of batteries and battery connectors. Check operation on battery power for several minutes to check that battery is charged and can hold a charge. Note battery charge indicator reading. Check condition of AC Power Supply and confirm that it charges the battery. When it is necessary to replace a battery, label it with the date of replacement.
11. Indicators and displays. Confirm the operation of all indicators, meters, gauges and visual displays. Inspect the vacuum gauge for cracks and scale visibility. Make sure the indicator rests on zero without vacuum applied.
12. Labeling. Verify that all labeling is intact and legible.

13. Accessories. Verify that clean canisters, suction catheters, and tubing are available. Inspect canisters for cracks, chips, leaks and general wear. Replace when necessary.

14. Mechanical overflow protection. To verify operation of mechanical overflow protection, liquids must be aspirated into collection canister until protective device is activated. (Observe when performing this test so that liquid is not aspirated into pump if overflow mechanism fails.)

Place a bucket of water on floor adjacent to unit. Connect a short length of hose to patient fitting on collection canister. Turn unit on and suction water into collection canister. Verify activation of overflow protection mechanism.

15. Disposable bacterial and hydrophobic filters. Replace when discolored or contacted with aspirate.

PERFORMANCE CHECKS: Performance checks should include, but not be limited to:

1. Grounding resistance. Perform this test with AC Power Supply connected. Using an ohmmeter, electrical safety analyzer, or multimeter with good resolution of fractional ohms, measure and record resistance between grounding pin of power cord and exposed (*unpainted and not anodized*) metal surfaces. Grounding resistance should not exceed 0.5 ohms.

2. Leakage current. Measure leakage current to ground with the grounding conductor of plug-connected equipment temporarily opened. Record the maximum leakage current with the unit off and on. Leakage current should not exceed 100 microamperes.

3. Maximum flow. Measure the maximum free airflow with a flowmeter. Perform this test with collection canister in place (*do not connect filter when performing this test*). **DO NOT** perform this test with a disposable canister having its own filter. Filters cause restrictions resulting in reduced free airflow. Connect flowmeter to "patient" fitting on collection canister using a small length of connecting hose. Use connecting hose having an inside diameter properly sized to collection canister fitting. Vacuum adjust valve must be set to its maximum vacuum position. This device should produce about 30 liters per minute free airflow.

4. Rate of vacuum rise (pumpdown). With collection canister in place and using a stopwatch, record the amount of time it takes to reach 300 mmHg when the "patient" fitting is occluded. Make sure that vacuum adjust valve is set to the maximum vacuum position. This device should reach 300 mmHg within 4- seconds when connected to a 1000 ml collection canister. **DO NOT** perform this test with a filter or with a collection canister that has its own integral filter.

5. Maximum vacuum. With collection canister in place, turn unit on, occlude "patient" fitting and note the maximum attainable vacuum. Make sure that vacuum adjust valve is set to the maximum vacuum position. This unit should produce about 550 mmHg. If lower vacuums are noted, check for air leaks in collection canister, hoses, and around fittings.

6. Vacuum gauge accuracy. Using a test gauge, compare vacuum readings at various vacuum settings. Readings between the unit's vacuum gauge and the test gauge should be within 10%.

CLEANING: Refer to the ROUTINE CARE AND MAINTENANCE "CLEANING" section in the OPERATION portion of this Manual.

1. Clean the exterior case and reusable collection canister as needed.

2. The vacuum pump/motor requires no lubrication.

3. Replace filters, hoses, tubing, fittings and connectors as required.

TROUBLESHOOTING

Continuous Suction

Symptom: No Vacuum or Weak Vacuum.

Controls: Mode Selector Switch - CONTinuous
Vacuum Adjust - MAXimum, fully clockwise.

Check: Hose and hose connections for cracks or crimps.
Verify that pump runs freely.
Check interconnections with collection canister(s). Verify that canister lid(s) is/are sealed tight.
Momentarily select INTermittent Mode and verify voltage at M1.
If no vacuum is present, check CB1 (3A), reset if necessary, and S1B for voltage presence.
If voltage is not present, verify that EXTERNAL POWER Jack is not damaged.

Symptom: No Internal Battery Power.

Controls: Mode Selector Switch - OFF/CHARGE

Check: Check for charging current going through R16 into BATT. See that L1 illuminates.
Remove AC Charger Rectifier. Select Internal 12 VDC power (CONTinuous Mode) and test output of BATT through S1A (allow adequate recharge time before testing for battery output).

Intermittent Suction

Symptom: Poor Intermittent Operation.

Controls: Mode Selector Switch - INTermittent
Vacuum ADJUST - MAXimum, fully clockwise.
ON Time - 5 seconds
OFF Time - 5 seconds

Check: Check 6-volt regulator (U4) for proper output.
Check U1D, pin 11, for proper clock signal (approximately 50Hz).
Check output of U1D, pin 4, for proper switching.
Check input to U6 for proper switching.
Check output of U6, and try to adjust R20 until motor operates. Adjust R20 for 200 mmHg.
Check overall appearance of PC Boards for shorts, bad solder connections and/or loose connector cables when erratic operation occurs.

Symptom: Motor does not turn on.

Check: Check circuit breaker CB1, reset if necessary.
Perform same checks as for poor intermittent suction above.
Check integrity of switch SW1B.

Symptom: ON or OFF Cycle does not end.

Check: Perform same checks as for poor intermittent suction above.
Check that clock signal appears at U1A, pin 3; verify that when U1A, pin 3, goes low, a trigger signal appears at the trigger input of U2, pin 2.
Check output of ON timer (U2, pin 3); verify that when U2, pin 3, goes low, a trigger pulse appears at U3, pin 4.
Check connections of timing resistors R4 and R5, and R9 and R10.

TECHNICAL DOCUMENTATION

NOTE: Refer to applicable Bill of Material for Part Number Description






FIGURE 4, COMPOSITE ILLUSTRATION DEPICTING MAJOR SUB-ASSEMBLIES



| ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N |
|------|-------------|------|-------------|------|-------------|------|-------------|
| 1 | 312-0018-00 | 6 | 338-0007-00 | 11 | 362-0440-06 | 16 | 416-0326-11 |
| 2 | 325-0326-01 | 7 | 346-0440-01 | 12 | 376-0007-00 | 17 | 450-0008-00 |
| 3 | 325-0326-02 | 8 | 346-0632-01 | 13 | 376-0008-00 | | |
| 4 | 334-0053-00 | 9 | 358-0632-06 | 14 | 392-0750-11 | | |
| 5 | 334-0054-00 | 10 | 358-0632-08 | 15 | 404-0326-11 | | |

SEE BILL OF MATERIAL 703-0326-01 (PAGE 13 - 21) FOR COMPLETE PART NUMBER DESCRIPTION


FIGURE 5, ASSEMBLY, BOTTOM COVER



| ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N |
|------|-------------|------|-------------|------|-------------|------|-------------|
| 1 | 312-0041-00 | 4 | 376-0007-00 | 7 | 392-0750-11 | 10 | 418-0750-31 |
| 2 | 338-0007-00 | 5 | 376-0008-00 | 8 | 404-0750-91 | 11 | 418-0750-41 |
| 3 | 358-0632-06 | 6 | 390-0002-00 | 9 | 416-0326-21 | 12 | 703-0326-02 |


SEE BILL OF MATERIAL 703-0326-03 (PAGE 13 - 23) FOR COMPLETE PART NUMBER DESCRIPTION

FIGURE 6, ASSEMBLY, TOP COVER



704-0750-03 Assembly, Battery Pack

| ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N |
|------|-------------|------|-------------|------|-------------|
| 1 | 016-0029-00 | 3 | 099-0009-02 | 5 | 600-0003-00 |
| 2 | 021-0019-00 | 4 | 310-0016-00 | 6 | 700-0326-17 |



703-0326-15 Assembly, Connector Panel


Subassemblies &
Sealants Not Shown:

12 13 14 15

| ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N |
|------|-------------|------|-------------|------|-------------|------|-------------|
| 1 | 016-0062-00 | 5 | 335-0003-00 | 9 | 404-0326-61 | 13 | 700-0326-01 |
| 2 | 081-0017-00 | 6 | 346-0440-01 | 10 | 422-0326-21 | 14 | 700-0326-24 |
| 3 | 089-0011-00 | 7 | 358-0440-08 | 11 | 480-0220-00 | 15 | 704-0326-01 |
| 4 | 253-1096-41 | 8 | 376-0035-00 | 12 | 602-0001-00 | 16 | 820-0037-00 |


SEE BILLS OF MATERIAL 703-0326-15 (PAGE 13 - 28) and 704-0750-03 (PAGE 13 - 31)
FOR COMPLETE PART NUMBER DESCRIPTION

FIGURE 7, ASSEMBLY, CONNECTOR PANEL
& ASSEMBLY, BATTERY PACK




Items Not
Shown

13 16 17 18



| ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N |
|------|-------------|------|-------------|------|-------------|------|-------------|
| 1 | 031-0002-00 | 6 | 358-0632-05 | 11 | 404-0326-41 | 16 | 700-0326-14 |
| 2 | 031-0010-00 | 7 | 374-0017-00 | 12 | 404-0326-51 | 17 | 700-0326-16 |
| 3 | 041-0013-00 | 8 | 374-0026-00 | 13 | 602-0006-00 | 18 | 703-0326-09 |
| 4 | 047-4935-00 | 9 | 376-0004-00 | 14 | 700-0326-03 | | |
| 5 | 358-0632-04 | 10 | 404-0326-31 | 15 | 700-0326-07 | | |


Not shown: 12



| ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N |
|------|-------------|------|-------------|------|-------------|------|-------------|
| 1 | 016-0004-00 | 6 | 305-0001-00 | 11 | 480-0229-00 | 16 | 700-0326-20 |
| 2 | 016-0016-00 | 7 | 348-0002-00 | 12 | 490-0045-01 | 17 | 700-0326-21 |
| 3 | 047-4935-00 | 8 | 480-0068-00 | 13 | 700-0326-08 | 18 | 704-0326-03 |
| 4 | 092-0006-00 | 9 | 480-0132-00 | 14 | 700-0326-12 | | |
| 5 | 100-0002-00 | 10 | 480-0210-00 | 15 | 700-0326-19 | | |


SEE BILLS OF MATERIAL 703-0326-05 (PAGE 13 - 24); and 703-0326-09
(PAGE 13 - 27) FOR COMPLETE PART NUMBER DESCRIPTION

FIGURE 8 ASSEMBLY, PUMP & ASSEMBLY, MANIFOLD



Hardware Items Not Shown
(for connecting Printed Circuit
Board and Pump Assemblies):

07 08 09 10 12 13




Sub-Assemblies Not Shown:

702-0326-01 703-0326-05
702-0326-02
702-0326-03

| ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N |
|------|-------------|------|-------------|------|-------------|------|-------------|
| 1 | 035-0002-00 | 8 | 346-0632-01 | 15 | 392-0026-00 | 22 | 540-0083-00 |
| 2 | 068-0004-00 | 9 | 358-0440-06 | 16 | 392-0029-00 | 23 | 602-0001-00 |
| 3 | 068-0005-00 | 10 | 368-0010-00 | 17 | 392-0032-00 | 24 | 700-0326-22 |
| 4 | 092-0006-00 | 11 | 368-0025-00 | 18 | 392-0033-00 | 25 | 700-0326-23 |
| 5 | 100-0010-00 | 12 | 368-0028-00 | 19 | 422-0326-31 | | |
| 6 | 315-0008-00 | 13 | 376-0019-00 | 20 | 478-0002-00 | | |
| 7 | 346-0440-01 | 14 | 392-0025-00 | 21 | 480-0219-00 | | |

SEE BILL OF MATERIAL 703-0326-02 (PAGE 13 - 22)
FOR COMPLETE PART NUMBER DESCRIPTION

FIGURE 9 ASSEMBLY, FRONT PANEL




Not Shown: 25


| ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N |
|------|-------------|------|-------------|------|-------------|------|-------------|
| 1 | 016-0004-00 | 8 | 117-0013-00 | 15 | 346-0440-01 | 22 | 374-0016-00 |
| 2 | 016-0034-00 | 9 | 312-0034-00 | 16 | 346-0632-01 | 23 | 402-0326-11 |
| 3 | 016-0043-00 | 10 | 325-0326-03 | 17 | 346-0832-01 | 24 | 414-0326-11 |
| 4 | 023-0020-00 | 11 | 334-0034-00 | 18 | 358-0440-04 | 25 | 606-0001-00 |
| 5 | 031-0009-00 | 12 | 334-0067-00 | 19 | 358-0440-06 | 26 | 708-0014-00 |
| 6 | 047-0005-00 | 13 | 340-0030-00 | 20 | 358-0632-10 | 27 | 708-0750-01 |
| 7 | 081-0020-00 | 14 | 340-0048-00 | 21 | 368-0045-00 | | |

SEE BILL OF MATERIAL 703-0326-06 (PAGE 13 - 25)
FOR COMPLETE PART NUMBER DESCRIPTION

FIGURE 10 ASSEMBLY, AC CHARGER RECTIFIER




| ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N |
|------|-------------|------|-------------|------|-------------|
| 1 | 410-0010-04 | 2 | 416-0016-00 | 3 | 540-0066-00 |



| ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N |
|------|-------------|------|-------------|------|-------------|
| 1 | 016-0068-00 | 2 | 099-0004-02 | 3 | 708-0750-01 |

SEE BILLS OF MATERIAL 703-0326-07 (PAGE 13 - 26); and 708-0326-01
(PAGE 13 - 32) FOR COMPLETE PART NUMBER DESCRIPTION


FIGURE 11 ASSEMBLY, COLLECTION CANISTER & AUTO POWER CABLE



| ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N |
|------|-------------|------|-------------|------|-------------|------|-------------|
| 1 | 016-0004-00 | 6 | 099-0010-02 | 11 | 700-0326-02 | 16 | 700-0326-10 |
| 2 | 016-0044-00 | 7 | 100-0008-00 | 12 | 700-0326-04 | 17 | 700-0326-11 |
| 3 | 031-0002-00 | 8 | 100-0009-00 | 13 | 700-0326-05 | 18 | 700-0326-13 |
| 4 | 092-0006-00 | 9 | 374-0017-00 | 14 | 700-0326-06 | 19 | 700-0326-18 |
| 5 | 092-0007-00 | 10 | 374-0018-00 | 15 | 700-0326-09 | | |

SEE BILL OF MATERIAL 704-0326-01 (PAGE 13 - 30)
FOR COMPLETE PART NUMBER DESCRIPTION


FIGURE 12 ASSEMBLY, WIRE HARNESS



| ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N |
|------|-------------|------|-------------|------|-------------|------|-------------|
| 1 | 081-0020-00 | 5 | 410-0004-00 | 9 | 703-0326-06 | 13 | 820-0018-00 |
| 2 | 334-0020-00 | 6 | 540-0051-00 | 10 | 708-0326-01 | 14 | 906-0326-03 |
| 3 | 334-0030-00 | 7 | 540-0055-00 | 11 | 820-0004-00 | | |
| 4 | 402-0017-00 | 8 | 540-0068-00 | 12 | 820-0005-00 | | |

SEE BILL OF MATERIAL 802-0326-02 (PAGE 13 - 33)
FOR COMPLETE PART NUMBER DESCRIPTION

FIGURE 13 ASSEMBLY, ACCESSORY KIT - AIR FORCE




| ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N | ITEM | IMPACT P/N |
|------|-------------|------|-------------|------|-------------|------|-------------|
| 1 | 081-0020-00 | 5 | 465-0005-00 | 9 | 703-0326-06 | 13 | 820-0018-00 |
| 2 | 334-0020-00 | 6 | 540-0051-00 | 10 | 703-0326-07 | 14 | 906-0326-03 |
| 3 | 334-0030-00 | 7 | 540-0055-00 | 11 | 703-0326-16 | | |
| 4 | 402-0017-00 | 8 | 540-0068-00 | 12 | 708-0326-01 | | |

SEE BILL OF MATERIAL 703-0326-16 (PAGE 13 - 29) and 802-0326-03 (PAGE 13 - 34)
FOR COMPLETE PART NUMBER DESCRIPTION

FIGURE 14 ASSEMBLY, ACCESSORY KIT, ARMY (NSN 6515-01-435-0050)

FIGURE 15



See Bill of Material 702-0326-03, Page 13-20, for complete Part Number Description

df:70203263.DWG

THIS DOCUMENT AND THE
INFORMATION IT CONTAINS ARE
THE SOLE PROPERTY OF IMPACT
INSTRUMENTATION, INC. THE USE
AND OR REPRODUCTION OF THIS
DOCUMENT FOR ANY PURPOSE
WITHOUT THE EXPRESS WRITTEN
PERMISSION OF IMPACT
INSTRUMENTATION, INC. IS
STRICTLY PROHIBITED.

TOLERANCES UNLESS
OTHERWISE SPECIFIED
FRACTION DEC ANGLES
+/- +/- +/-


| | |
|------------|---------|
| APPROVALS | DATE |
| DRAWN AG | 2-14-96 |
| CHECKED LS | 2-14-96 |
| | |
| | |

IMPACT

ASSEMBLY SWITCHER
PRINTED CIRCUIT BOARD

| | | |
|--------|------|--------------|
| SCALE | SIZE | DRAWING NO. |
| NONE | B | 702-0326-03 |
| Rev. D | | SHEET 1 OF 1 |

FIGURE 16



See Bill of Material 702-0326-02, Page 13-19, for complete Part Number Description

df:70203262.DWG

THIS DOCUMENT AND THE INFORMATION IT CONTAINS ARE THE SOLE PROPERTY OF IMPACT INSTRUMENTATION, INC. THE USE AND OR REPRODUCTION OF THIS DOCUMENT FOR ANY PURPOSE WITHOUT THE EXPRESS WRITTEN PERMISSION OF IMPACT INSTRUMENTATION, INC. IS STRICTLY PROHIBITED.

TOLERANCES UNLESS OTHERWISE SPECIFIED
FRACTION DEC ANGLES
+/- +/- +/-


| | |
|-------------------|---------|
| APPROVALS | DATE |
| DRAWN <i>HA</i> | 10/4/00 |
| CHECKED <i>HA</i> | 12/4/00 |

IMPACT

ASSEMBLY, REGULATOR
PRINTED CIRCUIT BOARD

| | | |
|---------------|--------------|----------------------------|
| SCALE NONE | SIZE B | DRAWING NO. 702-0326-02 |
| Rev. E | SHEET 1 OF 1 | |

FIGURE 17



See Bill of Material 702-0326-01, Page 13-18, for complete Part Number Description

df:70203261.DWG

THIS DOCUMENT AND THE INFORMATION IT CONTAINS ARE THE SOLE PROPERTY OF IMPACT INSTRUMENTATION, INC. THE USE AND OR REPRODUCTION OF THIS DOCUMENT FOR ANY PURPOSE WITHOUT THE EXPRESS WRITTEN PERMISSION OF IMPACT INSTRUMENTATION, INC., IS STRICTLY PROHIBITED.

TOLERANCES UNLESS OTHERWISE SPECIFIED
FRACTION DEC ANGLES
+/- +/- +/-

| | |
|------------|---------|
| APPROVALS | DATE |
| DRAWN AG | 2-14-96 |
| CHECKED LS | 2-14-96 |
| | |
| | |

IMPACT

ASSEMBLY, MAIN
PRINTED CIRCUIT BOARD

| | | |
|--------|--------------|-------------|
| SCALE | SIZE | DRAWING NO. |
| NONE | B | 702-0326-01 |
| Rev. E | SHEET 1 OF 1 | |




FIGURE 18 BLOCK DIAGRAM

IMPACT Instrumentation, Inc.

MASTER BOM Model 326M

Approved By: AG
Date: 12-4-00

| ITEM NUMBER | PART NUMBER | QTY | REV | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-----|-----|---|-------------|-------------------|
| 0001 | 701-0326-01 | 1 | B | Assembly, Final Mechanical | N/A | PL701-0326-01 |
| 0002 | 702-0326-01 | 1 | C | Assembly, Main Printed Circuit Board | N/A | PL702-0326-01 |
| 0003 | 702-0326-02 | 1 | E | Assembly, Regulator Printed Circuit Board | N/A | PL702-0326-02 |
| 0004 | 702-0326-03 | 1 | B | Assembly, Switcher Printed Circuit Board | N/A | PL702-0326-03 |
| 0005 | 703-0326-01 | 1 | A | Assembly, Bottom Cover | N/A | PL703-0326-01 |
| 0006 | 703-0326-02 | 1 | C | Assembly, Front Panel | N/A | PL703-0326-02 |
| 0007 | 703-0326-03 | 1 | A | Assembly, Top Cover | N/A | PL703-0326-03 |
| 0008 | 703-0326-05 | 1 | A | Assembly, Pump | N/A | PL703-0326-05 |
| 0009 | 703-0326-06 | 1 | B | Assembly, AC Charger Rectifier | N/A | PL703-0326-06 |
| 0010 | 703-0326-07 | 1 | 0 | Assembly, Collection Canister | N/A | PL703-0326-07 |
| 0011 | 703-0326-09 | 1 | A | Assembly, Manifold | N/A | PL703-0326-09 |
| 0012 | 703-0326-15 | 1 | 0 | Assembly, Connector Panel | N/A | PL703-0326-15 |
| 0013 | 703-0326-16 | 1 | 0 | Assembly, Carry Case, Hard | N/A | PL703-0326-16 |
| 0014 | 704-0326-01 | 1 | 0 | Assembly, Wire Harness | N/A | PL704-0326-01 |
| 0015 | 704-0750-03 | 1 | 0 | Assembly, Battery Pack, 4AH or 5AH | N/A | PL704-0750-03 |
| 0016 | 708-0326-01 | 1 | A | Assembly, Auto Power Cable | N/A | PL708-0326-01 |
| 0017 | 802-0326-02 | 1 | B | Assembly, Accessory Kit, Air Force | N/A | PL802-0326-02 |
| 0018 | 802-0326-03 | 1 | B | Assembly, Accessory Kit, Army | N/A | PL802-0326-03 |
| 0019 | 802-0326-05 | 1 | 0 | Assembly, Accessory Kit, DEPMEDS | N/A | PL802-0326-05 |

| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-------|---|-------------|-------------------|
| 0001 | 016-0004-00 | 2 | Tubing, Shrink, 3/16" I.D. X 1" Long | N/A | |
| 0002 | 305-0001-00 | 2 | Tie Wrap, Miniature, 4" Long | N/A | |
| 0003 | 312-0014-00 | 1 | Tape, Foam, White, P/S, 2" Long | N/A | |
| 0004 | 346-0440-01 | 2 | Nut, Keps, 4-40 | N/A | |
| 0005 | 358-0632-08 | 7 | Screw, Phillips, Pan Head, SS, 6-32 X 1/2 | N/A | |
| 0006 | 376-0007-00 | 8 | Washer, Flat, SS, #6 | N/A | |
| 0007 | 602-0006-00 | ASREQ | Sealant, Vibratite | N/A | |
| 0008 | 703-0326-01 | 1 | Assembly, Bottom Cover | N/A | PL703-0326-01 |
| 0009 | 703-0326-03 | 1 | Assembly, Top Cover | N/A | PL703-0326-03 |
| 0010 | 703-0326-15 | 1 | Assembly, Connector Panel | N/A | |
| 0011 | 704-0750-03 | 1 | Assembly, Battery Pack, 4AH or 5AH | N/A | PL704-0750-03 |

| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|--------------------------------------|---|---------------------|-------------------|
| 0001 | 012-0003-00 | ASREQ | Busswire, 18 AWG, Solid | N/A | |
| 0002 | 051-3643-00 | 1 | Transistor, 2N3643 | Q1 | |
| 0003 | 054-4001-00 | 1 | IC, DIP, CMOS, CD4001AE; Quad, 2-Input, NOR | U1 | |
| 0004 | 055-0358-00 | 1 | IC, DIP, LM358N | U12 | |
| 0005 | 055-0555-00 | 3 | IC, Timer, DIP, MOS, LM555CN | U,2,3,5 | |
| 0006 | 055-4310-00 | 3 | IC, ZVN4310A N-Channel, J-FET | Q2,3,4 | |
| 0007 | 055-7806-00 | 2 | IC, Case TO-220, 7806, +6V Voltage Regulator | U4,13 | |
| 0008 | 068-0005-00 | 1 | Lampholder w/lens, Green | LP1 | |
| 0009 | 093-0013-05 | 1 | Header, Right Angle, Friction Lock, 5-Pin, .100 Centers | BP1 | |
| 0010 | 093-0029-04 | 1 | Header, Right Angle, Friction Lock, 4-Pin, .100 Centers | BP2 | |
| 0011 | 093-0039-07 | 1 | Header, Straight, Dual Row, Male, 7 X 2 Pins, .1 Centers | BP3A | |
| 0012 | 109-0007-00 | 1 | Socket, IC, DIP, 14-Pin, Hi-Rel | N/A | |
| 0013 | 109-0008-00 | 4 | Socket, IC, DIP, 8-Pin, Hi-Rel | N/A | |
| 0014 | 117-0018-00 | 1 | Switch, Rotary, 3-Pole, 3-Position, Non-Shorting, PCB | SW1 | |
| 0015 | 200-0103-02 | 3 | Resistor, Fixed, Carbon Film, 1/4W, 5%, 10K | R3,8,17 | |
| 0016 | 200-0104-02 | 1 | Resistor, Fixed, Carbon Film, 1/4W, 5%, 100K | R2 | |
| 0017 | 200-0203-02 | 2 | Resistor, Fixed, Carbon Film, 1/4W, 5%, 20K | R6,11 | |
| 0018 | 200-0223-02 | 1 | Resistor, Fixed, Carbon Film, 1/4W, 5%, 22K | R1 | |
| 0019 | 200-0473-02 | 2 | Resistor, Fixed, Carbon Film, 1/4W, 5%, 47K | R13, 18 | |
| 0020 | 200-0474-02 | 1 | Resistor, Fixed, Carbon Film, 1/4W, 5%, 470K | R15 | |
| 0021 | 200-0564-02 | 1 | Resistor, Fixed, Carbon Film, 1/4W, 5%, 560K | R16 | |
| 0022 | 213-1002-00 | 1 | Resistor, Fixed, Precision, Metal Film, 1/4W, 1%, 10K | R31 | |
| 0023 | 213-1003-00 | 4 | Resistor, Fixed, Precision, Metal Film, 1/4W, 1%, 100K | R28,30,33,35 | |
| 0024 | 213-1212-00 | 1 | Resistor, Fixed, Precision, Metal Film, 1/4W, 1%, 12.1K | R34 | |
| 0025 | 213-3482-00 | 2 | Resistor, Fixed, Precision, Metal Film, 1/4W, 1%, 34.8K | R4,9 | |
| 0026 | 213-4701-00 | 2 | Resistor, Fixed, Precision, Metal Film, 1/4W, 1%, 4.7K | R27, 29 | |
| 0027 | 214-0254-01 | 2 | Resistor, Variable, PM, Linear Taper, 1/2W, 250K | R5, 10 | |
| 0028 | 216-0202-00 | 1 | Resistor, Vertical Mounting, Trim, 2K | R36 | |
| 0029 | 216-0203-00 | 1 | Resistor, Vertical Mounting, Trim, 20K | R32 | |
| 0030 | 216-0204-10 | 2 | Resistor, Vertical Mounting, Trim, 200K | R7, 12 | |
| 0031 | 240-0326-01 | 1 | Printed Circuit Board, Main | N/A | |
| 0032 | 253-1076-21 | 2 | Capacitor, Aluminum, Electrolytic, Radial, 100 ufd, 16V | C35, 36 | |
| 0033 | 253-2276-11 | 1 | Capacitor, Aluminum, Electrolytic, Radial, 220 ufd, 10V | C11 | |
| 0034 | 253-2276-21 | 1 | Capacitor, Aluminum, Electrolytic, Radial, SU, 220 ufd, 16V | C9 | |
| 0035 | 253-3366-21 | 2 | Capacitor, Aluminum, Electrolytic, Radial, SU, 33 ufd, 16V | C4, 7 | |
| 0036 | 255-4756-11 | 1 | Capacitor, Tantalum, Radial, 4.7 ufd, 10V | C14 | |
| 0037 | 257-1034-12 | 6 | Capacitor, Mylar, 10%, .01 ufd, 100V | C2, 3, 5, 6, 12, 13 | |
| 0038 | 259-1044-51 | 10 | Capacitor, Metallized Film, 5%, .1 ufd, 50V | C1,8,10,26-29,32-34 | |
| 0039 | 346-0440-01 | 2 | Nut, Keps, 4-40 | N/A | |
| 0040 | 358-0440-04 | 2 | Screw, Phillips, Pan Head, Zinc Plated, 4-40 X 1/4 | N/A | |
| 702-0326-01 | | Assembly, Main Printed Circuit Board | | MODEL 326M | |
| PL | | TITLE | | WHERE USED | |
| | | | | C REV | 10/4/00 DATE |

| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-----|---|-----------------------|-------------------|
| 0001 | 012-0008-00 | 1 | Buswire, 18 AWG, 1/2" Long | JP | |
| 0002 | 047-4004-00 | 3 | Diode, 1N4004, (1A, 400PIV) | D5, 6, 7 | |
| 0003 | 047-5825-00 | 2 | Diode, Schottkey, (5A, 40V) | D4, D8 | |
| 0004 | 051-3643-00 | 1 | Transistor, 2N3643 | Q6 | |
| 0005 | 055-0003-00 | 1 | IC, DIP, AMP02EP, Instrumentation Amplifier | U10 | |
| 0006 | 055-0338-00 | 1 | IC, Regulator, 5A, LM338K | U6 | |
| 0007 | 055-0340-00 | 1 | IC, TO-42, LM340LAZ-5.0, Voltage Regulator | U8 | |
| 0008 | 055-0358-00 | 1 | IC, DIP, LM358N | U11 | |
| 0009 | 055-7662-00 | 1 | IC, DIP, ICL7662, 8V Voltage Converter | U9 | |
| 0010 | 055-9540-00 | 1 | IC, Case 221A-03, IRF9540, P-Channel J-Fet | Q5 | |
| 0011 | 093-0036-03 | 1 | Header, Straight, Dual Row, Male, 3 X 2 Pins, .1 Centers | BP6A | |
| 0012 | 093-0037-05 | 1 | Header, Right Angle, 5 Pin, .156 Centers | BP4 | |
| 0013 | 093-0038-07 | 1 | Header, Straight, Dual Row, Female, 7 X 2 Pins, .1 Centers | BP3B | |
| 0014 | 109-0008-00 | 3 | Socket, IC, DIP, 8-Pin, Hi-Rel | N/A | |
| 0015 | 200-0102-02 | 1 | Resistor, Fixed, Carbon Film, 1/4W, 5%, 1K | R23 | |
| 0016 | 200-0181-02 | 1 | Resistor, Fixed, Carbon Film, 1/4W, 5%, 180 | R26 | |
| 0017 | 200-0241-02 | 1 | Resistor, Fixed, Carbon Film, 1/4W, 5%, 240 | R19 | |
| 0018 | 204-0010-02 | 1 | Resistor, Fixed, Metal Oxide Film, 1W, 5%, 1.0 | R21 | |
| 0019 | 213-1003-00 | 3 | Resistor, Fixed, Precision, Metal Film, 1/4W, 1%, 100K | R25, 37, 38 | |
| 0020 | 213-1823-00 | 1 | Resistor, Fixed, Precision, Metal Film, 1/4W, 1%, 182K | R24 | |
| 0021 | 213-5600-00 | 1 | Resistor, Fixed, Precision, Metal Film, 1/4W, 1%, 560 | R22 | |
| 0022 | 216-0502-00 | 1 | Resistor, Vertical Mounting, Trim, 5K | R20 | |
| 0023 | 240-0326-02 | 1 | Printed Circuit Board, Regulator | N/A | |
| 0024 | 252-1066-41 | 2 | Capacitor, Aluminum, Electrolytic, Axial, 10 ufd, 35V | C23, 24 | |
| 0025 | 253-2276-11-1 | 1 | Capacitor, Aluminum, Electrolytic, Radial, KA, 220 ufd, 10V | C19 | |
| 0026 | 255-1066-11 | 1 | Capacitor, Tantalum, Radial, 10 ufd, 16V | C17 | |
| 0027 | 257-1034-12 | 2 | Capacitor, Mylar, 10%, .01 ufd, 100V | C37, 38 | |
| 0028 | 259-1044-51 | 7 | Capacitor, Metallized Film, 5%, .1 ufd, 50V | C16,18,21,22,25,30,31 | |
| 0029 | 310-0024-00 | 1 | Insulator, Heat Sink Pad, TO-220 | N/A | |
| 0030 | 310-0025-00 | 1 | Insulator, Heat Sink Pad, TO-3 | N/A | |
| 0031 | 310-0326-00 | 1 | Heat Sink 6051B, For TO-3 Pkg | N/A | |
| 0032 | 310-0326-10 | 1 | Heat Sink 6072B, For 220A Pkg | N/A | |
| 0033 | 344-0003-00 | 1 | Lug, Solder, #10 | N/A | |
| 0034 | 346-0632-01 | 4 | Nut, Keps, 6-32 | N/A | |
| 0035 | 358-0632-06 | 4 | Screw, Phillips, Pan Head, SS, 6-32 X 3/8 | N/A | |
| 0036 | 700-0326-15 | 1 | Jumper, Red, 18AWG, 1" Long. | N/A | |

| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-------|---|-------------|-------------------|
| 0001 | 026-0001-00 | 1 | Module, Switcher, 11-30 V-in, 13.8 VDC-out | U7 | |
| 0002 | 031-0011-00 | 1 | Filter, Input Module | FLT1 | |
| 0003 | 093-0034-02 | 1 | Header, Right Angle, Friction Lock, .1 Centers, 2-pin | BP7 | |
| 0004 | 093-0035-03 | 1 | Header, Straight, Dual Row, Female, 3 x 2 Pin, .100 Centers | BP6B | |
| 0005 | 200-0513-02 | 1 | Resistor, Fixed, Carbon Film, 1/4W, 5%, 51K | RTRIM | |
| 0006 | 240-0326-03 | 1 | Printed Circuit Board, Switcher | N/A | |
| 0007 | 312-0058-00 | ASREQ | Tape, Copper, EMI/RFI, 1" Wide, 50' Long | N/A | |
| 0008 | 602-0001-00 | ASREQ | Sealant, Silicone Rubber, RTV, Translucent | N/A | |

| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-----|---|-------------|-------------------|
| 0001 | 312-0018-00 | 1 | Tape, Foam, White, P/S, 4" Long | N/A | |
| 0002 | 325-0326-01 | 1 | Label, P/S, Mylar, Condensed Operating Instructions & Warning | N/A | |
| 0003 | 325-0326-02 | 1 | Label, I.D. | N/A | |
| 0004 | 334-0053-00 | 1 | Strap, Velcro, 9" Long, w/D-Ring | N/A | |
| 0005 | 334-0054-00 | 1 | Strap, Velcro, 10" Long, w/ Hook | N/A | |
| 0006 | 338-0007-00 | 2 | Rivet, Blind, 1/8" Dia, Aluminum, Dome Head, 19/32" Long | N/A | |
| 0007 | 346-0440-01 | 4 | Nut, Keps, 4-40 | N/A | |
| 0008 | 346-0632-01 | 2 | Nut, Keps, 6-32 | N/A | |
| 0009 | 358-0632-06 | 2 | Screw, Phillips, Pan Head, SS, 6-32 X 3/8 | N/A | |
| 0010 | 358-0632-08 | 2 | Screw, Phillips, Pan Head, SS, 6-32 X 1/2 | N/A | |
| 0011 | 362-0440-06 | 4 | Screw, Truss Head, Phillips, SS, 4-40 X 3/8 | N/A | |
| 0012 | 376-0007-00 | 2 | Washer, Flat, SS, #6 | N/A | |
| 0013 | 376-0008-00 | 2 | Washer, Flat, SS, #4 | N/A | |
| 0014 | 392-0750-11 | 2 | Knob, Molded, Door | N/A | |
| 0015 | 404-0326-11 | 2 | Bracket, Slide, Molded, 3.20" | N/A | |
| 0016 | 416-0326-11 | 1 | Cover, Bottom, EMI/RFI, (Nickel Coated) | N/A | |
| 0017 | 450-0008-00 | 8 | Bumper, Rubber, Foot, P/S, Round, 1/2" Dia. X 1/8", Blk | N/A | |

| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-------|---|-------------|-------------------|
| 0001 | 035-0002-00 | 1 | Meter, Low Battery Indicator | N/A | |
| 0002 | 068-0004-00 | 2 | Bulb, Incandescent, 14V, 80mA | N/A | |
| 0003 | 068-0005-00 | 1 | Lampholder w/lens, Green | N/A | |
| 0004 | 092-0006-00 | 4 | Crimp Terminal, .1 Spacing | N/A | |
| 0005 | 100-0010-00 | 1 | Connector, 4-Pin, .100 Centers | N/A | |
| 0006 | 315-0008-00 | 1 | Gauge, Vacuum, 2" O.D. CBM, Dual Scale, 1/4" NPTM | N/A | |
| 0007 | 346-0440-01 | 1 | Nut, Keps, 4-40 | N/A | |
| 0008 | 346-0632-01 | 4 | Nut, Keps, 6-32 | N/A | |
| 0009 | 358-0440-06 | 4 | Screw, Phillips, Pan Head, SS, 4-40 X 3/8 | N/A | |
| 0010 | 368-0010-00 | 4 | Spacer, Hex, Threaded, Female, 4-40 X 5/8" Long, Brass/Zinc | N/A | |
| 0011 | 368-0025-00 | 5 | Spacer, Hex, Threaded, Male-to-Female, 4-40 X 3/4" Long, Brass/Nickel | N/A | |
| 0012 | 368-0028-00 | 4 | Spacer, Hex, Threaded, Male to Female, 4-40 X 1/2 Long, Alum. | N/A | |
| 0013 | 376-0019-00 | 4 | Washer, Lock, Internal Tooth, SS, #4 | N/A | |
| 0014 | 392-0025-00 | 2 | Knob, Collet, Gray w/Indicator, 1/8" Shaft | N/A | |
| 0015 | 392-0026-00 | 2 | Knob, Cap, Black w/Indicator, 1/8" Shaft | N/A | |
| 0016 | 392-0029-00 | 2 | Knob, Nut Cover, Black | N/A | |
| 0017 | 392-0032-00 | 1 | Knob, Collet, Black Matte w/Indicator, 21mm 1/4" Shaft | N/A | |
| 0018 | 392-0033-00 | 1 | Knob, Cap, Gray w/Indicator, 1/4" Shaft | N/A | |
| 0019 | 422-0326-31 | 1 | Panel, Front, Square Hole | N/A | |
| 0020 | 478-0002-00 | 1 | Orifice Restrictor, .010" Orifice, 10/32 THD X 1/8 ID Hose Barbs | N/A | |
| 0021 | 480-0219-00 | 1 | Cap, Brass, 1/4" NPTF | N/A | |
| 0022 | 540-0083-00 | 1 | Hose, 8" Long | N/A | |
| 0023 | 602-0001-00 | ASREQ | Sealant, Silicone Rubber, RTV, Translucent | N/A | |
| 0024 | 700-0326-22 | 1 | Jumper, Red, 22AWG, 5" Long | N/A | |
| 0025 | 700-0326-23 | 1 | Jumper, Green, 22 AWG, 5" Long | N/A | |
| 0026 | 702-0326-01 | 1 | Assembly, Main Printed Circuit Board | N/A | |
| 0027 | 702-0326-02 | 1 | Assembly, Regulator Printed Circuit Board | N/A | |
| 0028 | 702-0326-03 | 1 | Assembly, Switcher Printed Circuit Board | N/A | |
| 0029 | 703-0326-05 | 1 | Assembly, Pump | N/A | |

| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-----|--|-------------|-------------------|
| 0001 | 312-0041-00 | 2 | Tape, Foam, White, P/S, 1 1/2" Lg | N/A | |
| 0002 | 338-0007-00 | 1 | Rivet, Blind, 1/8" Dia, Aluminum, Dome Head, 19/32" Long | N/A | |
| 0003 | 358-0632-06 | 2 | Screw, Phillips, Pan Head, SS, 6-32 X 3/8 | N/A | |
| 0004 | 376-0007-00 | 2 | Washer, Flat, SS, #6 | N/A | |
| 0005 | 376-0008-00 | 1 | Washer, Flat, SS, #4 | N/A | |
| 0006 | 390-0002-00 | 1 | Handle (With Mounting Flanges) | N/A | |
| 0007 | 392-0750-11 | 1 | Knob, Molded, Door | N/A | |
| 0008 | 404-0750-91 | 1 | Bracket, Battery Compartment Divider EMI/RFI (Nickel Coated) | N/A | |
| 0009 | 416-0326-21 | 1 | Cover, Top, EMI/RFI, (Nickel Coated) | N/A | |
| 0010 | 418-0750-31 | 1 | Door, Battery Compartment Divider, (Nickel Coated) | N/A | |
| 0011 | 418-0750-41 | 1 | Door, Top Cover | N/A | |
| 0012 | 703-0326-02 | 1 | Assembly, Front Panel | N/A | |

| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-------|---|-------------|-------------------|
| 0001 | 031-0002-00 | 1 | Ferrite | N/A | |
| 0002 | 031-0010-00 | 1 | Filter, Low Pass, EMI | N/A | |
| 0003 | 041-0013-00 | 1 | Pump, Vacuum, Diaphragm, In-Line Ports, 12V | N/A | |
| 0004 | 047-4935-00 | 1 | Diode, 1N4935, Fast Recovery | N/A | |
| 0005 | 358-0632-04 | 2 | Screw, Phillips, Pan Head, SS, 6-32 X 1/4 | N/A | |
| 0006 | 358-0632-05 | 4 | Screw, Phillips, Pan Head, SS, 6-32 X 5/16 | N/A | |
| 0007 | 374-0017-00 | 2 | Terminal, Insulated, Ring, #6 Stud, 16-14 AWG, Blue | N/A | |
| 0008 | 374-0026-00 | 1 | Terminal, Insulated, Ring, 1/4" Stud, 16-14 AWG, Blue | N/A | |
| 0009 | 376-0004-00 | 6 | Washer, Lock, Internal Tooth, SS, #6 | N/A | |
| 0010 | 404-0326-31 | 1 | Bracket, Left | N/A | |
| 0011 | 404-0326-41 | 1 | Bracket, Right | N/A | |
| 0012 | 404-0326-51 | 1 | Bracket, Filter Mounting | N/A | |
| 0013 | 602-0006-00 | ASREQ | Sealant, Vibratite | N/A | |
| 0014 | 700-0326-03 | 1 | Jumper, Black, 18AWG, 4" Long. | N/A | |
| 0015 | 700-0326-07 | 1 | Jumper, Black, 18AWG, 8" Long. | N/A | |
| 0016 | 700-0326-14 | 1 | Jumper, Red, 18AWG, 11" Long. | N/A | |
| 0017 | 700-0326-16 | 1 | Jumper, Red, 18AWG, 1 1/2" Long. | N/A | |
| 0018 | 703-0326-09 | 1 | Assembly, Manifold | N/A | |

| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-------|---|-------------|-------------------|
| 0001 | 016-0004-00 | 5 | Tubing, Shrink, 3/16" I.D. X 1" Long | N/A | |
| 0002 | 016-0034-00 | 2 | Tubing, Shrink, 1/8" I.D. X 1" Long | N/A | |
| 0003 | 016-0043-00 | 1 | Tubing, Shrink, 3/16 I.D. X 4" Long | N/A | |
| 0004 | 023-0020-00 | 1 | Transformer, Power, 115/230VAC Primary to 14VAC, 3A Secondary | N/A | |
| 0005 | 031-0009-00 | 1 | Attenuator, Split Channel, Assembly | N/A | |
| 0006 | 047-0005-00 | 1 | Diode, Bridge Rectifier, 12A, 50PIV | N/A | |
| 0007 | 081-0020-00 | 1 | Fuse, Midget, 2A, 250V | N/A | |
| 0008 | 117-0013-00 | 1 | Switch, Rotary, Voltage Select | N/A | |
| 0009 | 312-0034-00 | 1 | Tape, Foam, White, 3" Lg | N/A | |
| 0010 | 325-0326-03 | 1 | Label, Serial Number, & Heat Sink, 326/326M Power Supply | N/A | |
| 0011 | 334-0034-00 | 1 | Fuseholder, Universal, In-line, Red Leads | N/A | |
| 0012 | 334-0067-00 | 1 | Collar, Strain Relief, Cable Mounting | N/A | |
| 0013 | 340-0030-00 | 1 | Bushing, Strain Relief, Nylon, Lg. w/Nut, LTF 11, Black | N/A | |
| 0014 | 340-0048-00 | 1 | Bushing, Strain Relief, Nylon, Lg. W/Nut, LFT9, Black | N/A | |
| 0015 | 346-0440-01 | 4 | Nut, Keps, 4-40 | N/A | |
| 0016 | 346-0632-01 | 1 | Nut, Keps, 6-32 | N/A | |
| 0017 | 346-0832-01 | 2 | Nut, Keps, 8-32 | N/A | |
| 0018 | 358-0440-04 | 4 | Screw, Phillips, Pan Head, Zinc Plated, 4-40 X 1/4 | N/A | |
| 0019 | 358-0440-06 | 4 | Screw, Phillips, Pan Head, SS, 4-40 X 3/8 | N/A | |
| 0020 | 358-0632-10 | 1 | Screw, Phillips, Pan Head, SS, 6-32 X 5/8 | N/A | |
| 0021 | 368-0045-00 | 1 | Spacer, Nylon, .625" O.D. X .342 I.D. X .375" Long | N/A | |
| 0022 | 374-0016-00 | 2 | Terminal, Insulated, Ring, #10 Stud, 16-14 AWG, Blue | N/A | |
| 0023 | 402-0326-11 | 1 | Case, AC Charger Rectifier | N/A | |
| 0024 | 414-0326-11 | 1 | Chassis, Charger/Rectifier | N/A | |
| 0025 | 606-0001-00 | ASREQ | Silicone Grease Or Heat Sink Compound | N/A | |
| 0026 | 708-0014-00 | 1 | Line Cord, SJT, 18 AWG, 3 Cord, 120" Long | N/A | |
| 0027 | 708-0750-01 | 1 | Assembly, Cable, Shielded, 2-Con-w/molded 2.5 MM Plug, 5ft | N/A | |

| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-----|---|-------------|-------------------|
| 0001 | 410-0010-04 | 1 | Jar, Collection, Polysulfone/Polycarbonate, | N/A | |
| 0002 | 416-0016-00 | 1 | Cap, Molded, 83mm | N/A | |
| 0003 | 540-0066-00 | 2 | Tube, Elbow, SS, .319" I.D. X .375" O.D. | N/A | |

703-0326-07

PL

Assembly, Collection Canister

TITLE

MODEL 326M

WHERE USED

0

REV

10/4/00

DATE

Page 13-26

| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-----|--|-------------|-------------------|
| 0001 | 016-0004-00 | 2 | Tubing, Shrink, 3/16" I.D. X 1" Long | N/A | |
| 0002 | 016-0016-00 | 1 | Tubing, Shrink, 3/8" I.D. X 1 1/2" Long | N/A | |
| 0003 | 047-4935-00 | 2 | Diode, 1N4935, Fast Recovery | N/A | |
| 0004 | 092-0006-00 | 5 | Crimp Terminal, .1 Spacing | N/A | |
| 0005 | 100-0002-00 | 1 | Connector, 5-Pin | N/A | |
| 0006 | 305-0001-00 | 3 | Tie Wrap, Miniature, 4" Long | N/A | |
| 0007 | 348-0002-00 | 1 | Pin, Stop, Knob | N/A | |
| 0008 | 480-0068-00 | 1 | Hose Barb, Chrome, 1/8" NPTM to 3/8" I.D. Tube | N/A | |
| 0009 | 480-0132-00 | 1 | Nipple, Close, Chrome, 1/8" NPTM X 1/8" NPTM , .500 Long | N/A | |
| 0010 | 480-0210-00 | 1 | Elbow, 45 Deg., Street, 1/8"NPTM x 1/8"NPTF | N/A | |
| 0011 | 480-0229-00 | 1 | Elbow, Swivel, 10-32 x 1/8" Tube, Black | N/A | |
| 0012 | 490-0045-01 | 1 | Valve, Plug | N/A | |
| 0013 | 700-0326-08 | 1 | Jumper, Black, 22AWG, 9" Long. | N/A | |
| 0014 | 700-0326-12 | 1 | Jumper, Red, 22AWG, 9" Long. | N/A | |
| 0015 | 700-0326-19 | 1 | Jumper, Orange, 22AWG, 9" Long. | N/A | |
| 0016 | 700-0326-20 | 1 | Jumper, Orange, 22AWG, 11" Long. | N/A | |
| 0017 | 700-0326-21 | 1 | Jumper, Grey, 22AWG, 11" Long. | N/A | |
| 0018 | 704-0326-03 | 1 | Assembly, Manifold/Valve | N/A | |

| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-------|--|-------------|-------------------|
| 0001 | 016-0062-00 | 1 | Tubing Cable Wrap, Blk, 1/2" X 5" | N/A | |
| 0002 | 081-0017-00 | 1 | Circuit Breaker, 5A, 32VDC, 250VAC | N/A | |
| 0003 | 089-0011-00 | 1 | Connector, PMF, 2.5mm PIN, Moulded PC Mount w/Switch | N/A | |
| 0004 | 253-1096-41 | 1 | Capacitor, Aluminum, Electrolytic, Radial, 10,000 ufd, 35V | N/A | |
| 0005 | 335-0003-00 | 1 | Flange, Panel Mounting | N/A | |
| 0006 | 346-0440-01 | 6 | Nut, Keps, 4-40 | N/A | |
| 0007 | 358-0440-08 | 6 | Screw, Phillips, Pan Head, SS, 4-40 X 1/2 | N/A | |
| 0008 | 376-0035-00 | 2 | Washer, Flat, Fiber, #4 | N/A | |
| 0009 | 404-0326-61 | 2 | Bracket, Jack, Power-Type | N/A | B404-0326-61 |
| 0010 | 422-0326-21 | 1 | Panel, Connector | N/A | B422-0326-21 |
| 0011 | 480-0220-00 | 1 | Hose Barb, Union Bulkhead, PM, 3/8" X 3/8 | N/A | |
| 0012 | 602-0001-00 | ASREQ | Sealant, Silicone Rubber, RTV, Translucent | N/A | |
| 0013 | 700-0326-01 | 1 | Jumper, Black, 18AWG, 2 1/2" Long. | N/A | |
| 0014 | 700-0326-24 | 1 | Jumper, Red, 18 AWG, 3" Long | N/A | |
| 0015 | 704-0326-01 | 1 | Assembly, Wire Harness | N/A | |
| 0016 | 820-0037-00 | 1 | Hose, Spiral, 8" Long | N/A | |

| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-----|---|-------------|-------------------|
| 0001 | 313-0009-00 | 1 | Foam, Die Cut, Pad, Lid, 19 3/8" x 12 3/8" x 1" | N/A | |
| 0002 | 313-0010-00 | 1 | Foam, Die Cut, Pad, Base, 18 1/4" x 14 1/4" x 1 | N/A | |
| 0003 | 313-0013-00 | 1 | Foam, Die Cut, Pad, Insert | N/A | |
| 0004 | 325-0300-07 | 1 | Label, P/S, Anodized Mylar, Name & Model | N/A | |
| 0005 | 402-0005-00 | 1 | Case, HDPE, Black, 23" X 16" X 7" | N/A | |

| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-----|---|-------------|-------------------|
| 0001 | 016-0004-00 | 2 | Tubing, Shrink, 3/16" I.D. X 1" Long | N/A | |
| 0002 | 016-0044-00 | 1 | Tubing, Shrink 3/4 I.D. X 3" Long | N/A | |
| 0003 | 031-0002-00 | 4 | Ferrite | N/A | |
| 0004 | 092-0006-00 | 2 | Crimp Terminal, .1 Spacing | N/A | |
| 0005 | 092-0007-00 | 5 | Crimp Terminal, .156 Spacing | N/A | |
| 0006 | 099-0010-02 | 1 | Plug, Black Wire Insulated, | N/A | |
| 0007 | 100-0008-00 | 1 | Connector, 5-Pin, .156 Centers | N/A | |
| 0008 | 100-0009-00 | 1 | Connector, 2-Pin, .100 Centers | N/A | |
| 0009 | 374-0017-00 | 1 | Terminal, Insulated, Ring, #6 Stud, 16-14 AWG, Blue | N/A | |
| 0010 | 374-0018-00 | 2 | Terminal, Insulated, Ring, #4 Stud, 22-18 AWG, Red | N/A | |
| 0011 | 700-0326-02 | 1 | Jumper, Black, 18AWG, 3" Long. | N/A | |
| 0012 | 700-0326-04 | 1 | Jumper, Black, 18AWG, 6" Long. | N/A | |
| 0013 | 700-0326-05 | 1 | Jumper, Black, 18 AWG, 7" Long | N/A | |
| 0014 | 700-0326-06 | 1 | Jumper, Black, 18AWG, 9" Long. | N/A | |
| 0015 | 700-0326-09 | 1 | Jumper, Black, 22AWG, 11" Long. | N/A | |
| 0016 | 700-0326-10 | 1 | Jumper, Red, 18AWG, 9" Long. | N/A | |
| 0017 | 700-0326-11 | 1 | Jumper, Red, 22AWG, 7" Long. | N/A | |
| 0018 | 700-0326-13 | 1 | Jumper, Red, 22AWG, 11" Long. | N/A | |
| 0019 | 700-0326-18 | 1 | Jumper, Orange, 22AWG, 7" Long. | N/A | |

| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-------|--|-------------|-------------------|
| 0001 | 016-0029-00 | 1 | Tubing, Shrink, 123mm Flat X 9 3/4" Long | N/A | |
| 0002 | 021-0019-00 | 2 | Battery, Gel-Cel, 6VDC, 5.0AH, 6VDC, 5.6AH | N/A | |
| 0003 | 099-0009-02 | 1 | Plug, Red Wire Insulated, | N/A | |
| 0004 | 310-0016-00 | 1 | Insulator, Battery | N/A | |
| 0005 | 600-0003-00 | ASREQ | Glue, Stick | N/A | |
| 0006 | 700-0326-17 | 1 | Jumper, Orange, 18AWG, 3" Long. | N/A | |

| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-----|--|-------------|-------------------|
| 0001 | 016-0068-00 | 1 | Tubing, Shrink, 3/8" I.D. X 5/8" Long, | N/A | |
| 0002 | 099-0004-02 | 1 | Plug, CMM, 2-Cond, Auto Type | N/A | |
| 0003 | 708-0750-01 | 1 | Assembly, Cable, Shielded, 2-Con-w/molded 2.5 MM Plug, 5ft | N/A | |


| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-----|--|-------------|-------------------|
| 0001 | 081-0020-00 | 2 | Fuse, Midget, 2A, 250V | N/A | |
| 0002 | 334-0020-00 | 2 | Strap, Velcro/Webbing 1"W X 24"L | N/A | |
| 0003 | 334-0030-00 | 2 | Holder, Collection Bottle, Chrome Plated | N/A | |
| 0004 | 402-0017-00 | 1 | Case, Padded, Aspirator, W/O Accessories | N/A | |
| 0005 | 410-0004-00 | 2 | Jar, Disposable, w/Lid, 1200cc, w/ Impact Name | N/A | |
| 0006 | 540-0051-00 | 2 | Hose, 2' Long | N/A | |
| 0007 | 540-0055-00 | 2 | Hose 18" Long | N/A | |
| 0008 | 540-0068-00 | 1 | Hose, 12" Long | N/A | |
| 0009 | 703-0326-06 | 1 | Assembly, AC Charger Rectifier | N/A | |
| 0010 | 708-0326-01 | 1 | Assembly, Auto Power Cable | N/A | |
| 0011 | 820-0004-00 | 1 | Catheter, 14 French | N/A | |
| 0012 | 820-0005-00 | 1 | Catheter, 18 French | N/A | |
| 0013 | 820-0018-00 | 1 | Tubing,Suction, Sterile, 9/32" I.D. X 6' Long | N/A | |
| 0014 | 906-0326-03 | 2 | Manual, Operation/Service, Model 326M | N/A | |

| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-----|--|-------------|-------------------|
| 0001 | 081-0020-00 | 2 | Fuse, Midget, 2A, 250V | N/A | |
| 0002 | 334-0020-00 | 2 | Strap, Velcro/Webbing 1"W X 24"L | N/A | |
| 0003 | 334-0030-00 | 2 | Holder, Collection Bottle, Chrome Plated | N/A | |
| 0004 | 402-0017-00 | 1 | Case, Padded, Aspirator, W/O Accessories | N/A | |
| 0005 | 465-0005-00 | 2 | Filter, Hydro/Bact, Disposable, w/Hose Barbs | N/A | |
| 0006 | 540-0051-00 | 2 | Hose, 2' Long | N/A | |
| 0007 | 540-0055-00 | 2 | Hose 18" Long | N/A | |
| 0008 | 540-0068-00 | 1 | Hose, 12" Long | N/A | |
| 0009 | 703-0326-06 | 1 | Assembly, AC Charger Rectifier | N/A | |
| 0010 | 703-0326-07 | 2 | Assembly, Collection Canister | N/A | |
| 0011 | 703-0326-16 | 1 | Assembly, Carry Case, Hard | N/A | |
| 0012 | 708-0326-01 | 1 | Assembly, Auto Power Cable | N/A | |
| 0013 | 820-0018-00 | 1 | Tubing, Suction, Sterile, 9/32" I.D. X 6' Long | N/A | |
| 0014 | 906-0326-03 | 2 | Manual, Operation/Service, Model 326M | N/A | |


| ITEM NUMBER | PART NUMBER | QTY | ITEM NAME / DESCRIPTION | DESIGNATION | DRAWING NUMBER |
|----------------|----------------|-----|--|-------------|-------------------|
| 0001 | 081-0020-00 | 2 | Fuse, Midget, 2A, 250V | N/A | |
| 0002 | 334-0020-00 | 2 | Strap, Velcro/Webbing 1"W X 24"L | N/A | |
| 0003 | 334-0030-00 | 2 | Holder, Collection Bottle, Chrome Plated | N/A | |
| 0004 | 402-0017-00 | 1 | Case, Padded, Aspirator, W/O Accessories | N/A | |
| 0005 | 402-0326-02 | 1 | Case, Reusable Transit, 326M | N/A | |
| 0006 | 465-0005-00 | 2 | Filter, Hydro/Bact, Disposable, w/Hose Barbs | N/A | |
| 0007 | 540-0051-00 | 2 | Hose, 2' Long | N/A | |
| 0008 | 540-0055-00 | 2 | Hose 18" Long | N/A | |
| 0009 | 540-0068-00 | 1 | Hose, 12" Long | N/A | |
| 0010 | 703-0326-06 | 1 | Assembly, AC Charger Rectifier | N/A | |
| 0011 | 703-0326-07 | 2 | Assembly, Collection Canister | N/A | |
| 0012 | 708-0326-01 | 1 | Assembly, Auto Power Cable | N/A | |
| 0013 | 820-0018-00 | 1 | Tubing, Suction, Sterile, 9/32" I.D. X 6' Long | N/A | |
| 0014 | 906-0326-03 | 2 | Manual, Operation/Service, Model 326M | N/A | |

Model 326M Wire List

| Part # | Sub Assy # | Description | Color | Gauge | Length | Strip | Tin | Origin | Strip | Tin | Destination |
|-------------|-------------|------------------------------------|-------|-------|--------|-------|-----|----------------|-------|-----|--------------------|
| 700-0326-01 | 703-0326-15 | Jumper, Black, 18AWG, 2 1/2" Long. | BLK | 18 | 2 1/2 | 3/8 | Y | J1-1 | 3/8 | Y | CASE C20(-) |
| 700-0326-02 | 704-0326-01 | Jumper, Black, 18AWG, 3" Long. | BLK | 18 | 3 | 3/8 | Y | J1-1 | 3/8 | | CASE CONN.PANEL(-) |
| 700-0326-03 | 703-0326-05 | Jumper, Black, 18AWG, 4" Long. | BLK | 18 | 4 | 3/8 | | LOW PASS(-) | 3/8 | | CASE (-) |
| 700-0326-04 | 704-0326-01 | Jumper, Black, 18AWG, 6" Long. | BLK | 18 | 6 | 3/8 | Y | CONN PANEL (-) | 3/8 | | CASE (-) |
| 700-0326-05 | 704-0326-01 | Jumper, Black, 18 AWG, 7" Long | BLK | 18 | 7 | 3/8 | | BP4-5 | 3/8 | | CASE (-) |
| 700-0326-06 | 704-0326-01 | Jumper, Black, 18AWG, 9" Long. | BLK | 18 | 9 | 3/8 | | BP4-4 | 3/8 | Y | BATT (-) |
| 700-0326-07 | 703-0326-05 | Jumper, Black, 18AWG, 8" Long. | BLK | 18 | 8 | 3/8 | Y | PUMP(-) | 3/8 | | CASE (-) |
| 700-0326-08 | 703-0326-09 | Jumper, Black, 22AWG, 9" Long. | BLK | 22 | 9 | 3/8 | | BP1-4 | 3/8 | Y | DUMP (-) |
| 700-0326-09 | 704-0326-01 | Jumper, Black, 22AWG, 11" Long. | BLK | 22 | 11 | 3/8 | Y | J1-1 | 3/8 | | SWITCHER (-) |
| 700-0326-10 | 704-0326-01 | Jumper, Red, 18AWG, 9" Long. | RED | 18 | 9 | 3/8 | Y | BATT (+) | 3/8 | | BP4-3 |
| 700-0326-11 | 704-0326-01 | Jumper, Red, 22AWG, 7" Long. | RED | 22 | 7 | 3/8 | Y | J1-4 | 3/8 | | BP4-1 |
| 700-0326-12 | 703-0326-09 | Jumper, Red, 22AWG, 9" Long. | RED | 22 | 9 | 3/8 | Y | DUMP (+) | 3/8 | | BP1-3 |
| 700-0326-13 | 704-0326-01 | Jumper, Red, 22AWG, 11" Long. | RED | 22 | 11 | 3/8 | Y | J1-2 | 3/8 | | SWITCHER (+) |
| 700-0326-14 | 703-0326-05 | Jumper, Red, 18AWG, 11" Long. | RED | 18 | 11 | 3/8 | Y | CB | 3/8 | Y | LOW PASS B(+) |
| 700-0326-15 | 702-0326-02 | Jumper, Red, 18AWG, 1" Long. | RED | 18 | 1 | 3/8 | Y | U6-3 | 3/8 | Y | C19 (+) |
| 700-0326-16 | 703-0326-05 | Jumper, Red, 18AWG, 1 1/2" Long. | RED | 18 | 1 1/2 | 3/8 | Y | MOTOR (+) | 3/8 | Y | LOW PASS A(+) |
| 700-0326-17 | 704-0750-03 | Jumper, Orange, 18AWG, 3" Long. | ORG | 18 | 3 | 3/8 | Y | B1 (-) | 3/8 | Y | B2 (+) |
| 700-0326-18 | 704-0326-01 | Jumper, Orange, 22AWG, 7" Long. | ORG | 22 | 7 | 3/8 | Y | J1-3 | 3/8 | | BP4-2 |
| 700-0326-19 | 703-0326-09 | Jumper, Orange, 22AWG, 9" Long. | ORG | 22 | 9 | 3/8 | Y | CB | 3/8 | | BP1-5 |
| 700-0326-20 | 703-0326-09 | Jumper, Orange, 22AWG, 11" Long. | ORG | 22 | 11 | 3/8 | Y | BLEED (+) | 3/8 | | BP1-1 |
| 700-0326-21 | 703-0326-09 | Jumper, Grey, 22AWG, 11" Long. | GRY | 22 | 11 | 3/8 | Y | BLEED (-) | 3/8 | | BP1-2 |
| 700-0326-22 | 703-0326-02 | Jumper, Red, 22AWG, 5" Long | RED | 22 | 5 | 3/8 | Y | EXTPWR+ | 3/8 | | BP2-3 |
| 700-0326-23 | 703-0326-02 | Jumper, Green, 22 AWG, 5" Long | GREEN | 22 | 5 | 3/8 | Y | EXTPWR- | 3/8 | | BP2-4 |
| 700-0326-24 | 703-0326-15 | Jumper, Red, 18 AWG, 3" Long | RED | 18 | 3 | 3/8 | Y | J1-2 | 3/8 | Y | C20+ |



| | | | |
|---|---------|------------------------------|------------------------------|
| Approval | Date | Impact Instrumentation, Inc. | |
| Drawn | 12/4/00 | West Caldwell, NJ 07008 | |
| Checked | 12/4/02 | 326/326M Aspirator | |
| Rev. | C | SIZE | DRAWING NO. 240-0326-00-4 |
| | | Sheet 1 of 4 | |
| This Document and the information it contains are the sole property of Impact Instrumentation Inc. The use or reproduction of this property for any purpose without the written permission of Impact Instrumentation Inc. is strictly prohibited. | | | |




| | | | |
|----------|------|------------------------------|---------------|
| Approval | Date | Impact Instrumentation, Inc. | |
| Drawn | HK | West Caldwell, NJ 07008 | |
| Checked | FW | 326/326M Main PCB | |
| Rev. | C | SIZE | DRAWING NO. |
| | | | 240-0326-01-4 |
| | | | Sheet 2 of 4 |


This Document and the information it contains are the sole property of Impact Instrumentation Inc. The use or reproduction of this property for any purpose without the written permission of Impact Instrumentation Inc. is strictly prohibited.

| | | | |
|----------------|----------|---|------------------------------|
| Approval | Date | Impact Instrumentation, Inc. West Caldwell, NJ 07006 | |
| Drawn YJA | 12/14/00 | 326/326M Regulator PCB | |
| Checked PAC | 12/19/00 | SIZE | DRAWING NO. 240-0326-02-4 |
| Rev. D | | Sheet 3 of 4 | |

This Document and the information it contains are the sole property of Impact Instrumentation Inc. The use or reproduction of this property for any purpose without the written permission of Impact Instrumentation Inc. is strictly prohibited.



| | | | |
|---|------|------------------------------|---------------|
| Approval | Date | Impact Instrumentation, Inc. | |
| Drawn | HK | West Caldwell, NJ 07006 | |
| Checked | FW | 326/326M Switcher PCB | |
| Rev. | B | SIZE | DRAWING NO. |
| | | | 240-0326-03-4 |
| | | Sheet 4 of 4 | |
| This Document and the information it contains are the sole property of Impact Instrumentation Inc. The use or reproduction of this property for any purpose without the written permission of Impact Instrumentation Inc. is strictly prohibited. | | | |



| | | | |
|---|---------|---|------------------------------|
| Approval | Date | Impact Instrumentation, Inc. West Caldwell, NJ 07008 | |
| Drawn HK | 8/16/99 | 326/326M AC/DC Power Supply | |
| Checked FW | 8/16/99 | SIZE | DRAWING NO. 703-0326-06-4 |
| Rev. A | | Sheet 1 of 1 | |
| This Document and the information it contains are the sole property of Impact Instrumentation Inc. The use or reproduction of this property for any purpose without the written permission of Impact Instrumentation Inc. is strictly prohibited. | | | |