
Nellcor™


Adult SpO₂ Sensor

en	Instructions For Use
fr	Mode d'emploi
de	Gebrauchsanweisung
nl	Gebruiksaanwijzing
it	Istruzioni per l'uso
es	Instrucciones de uso
sv	Bruksanvisning
da	Betjeningsvejledning
no	Bruksanvisning
fi	Käyttöohjeet
pt	Instruções de uso
ru	Руководство по применению
zh	使用说明
pl	Instrukcja użycia
cs	Návod k použití
sk	Návod na použitie
hu	Használati utasítás
el	Οδηγίες χρήσης
tr	Kullanma Talimatları


REF MAXA
REF MAXAL
REF MAXAI
REF MAXALI


1

2

3

MAXA, MAXAL

Adult SpO₂ Sensor

Identification of a substance that is contained or present within the product or packaging.

Identification of a substance that is not contained or present within the product or packaging.

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in bio-incompatibility, infection or product failure risks to the patient.

This device is not made with natural rubber latex or DEHP.

Directions for Use

Indications/Contraindications

The Nellcor™ Adult SpO₂ Sensor, model MAXA(L) is indicated for single patient use when continuous noninvasive arterial oxygen saturation and pulse rate monitoring are required for patients weighing more than 30 kg.

The MAXA(L) is contraindicated for use on patients who exhibit allergic reactions to the adhesive tape.

Note: Adhesive sensor consists of sensor, cable and wrap adhesive flaps.

Instructions for Use

1. Remove plastic backing from the MAXA(L) and locate transparent windows (a) on the adhesive side. Windows cover optical components. **1**

An index finger is the preferred MAXA(L) location. Alternatively, apply the sensor to a small thumb, smaller finger, or great toe.

Note: When selecting a sensor site, priority should be given to an extremity free of an arterial catheter, blood pressure cuff, or intravascular infusion line.

2. Orient the MAXA(L) so the dashed line in the middle of the sensor is centered on the tip of the digit. Wrap adhesive flaps on non-cable end around the digit. Note that the cable must be positioned on the top of the hand. **2**
3. Fold cable end over top of digit so that windows are directly opposite each other. Wrap adhesive securely around sides of digit. **3**
4. Plug the MAXA(L) into the oximeter and verify proper operation as described in the oximeter operator's manual.

Note: If the sensor does not track the pulse reliably, it may be incorrectly positioned—or the sensor site may be too thick, thin, or deeply pigmented, or otherwise deeply colored (for example, as a result of externally applied coloring such as nail polish, dye, or pigmented cream) to permit appropriate light transmission. If any of these situations occurs, reposition the sensor or choose an alternate Nellcor sensor for use on a different site.

To remove the Adhesive Sensor:

1. Peel the tape wrap to gently remove the sensor.
2. Disconnect the sensor from the monitor by unplugging the sensor from the extension cable or monitor.

WARNINGS

1. **Use this sensor only with Nellcor OxiMax™ instruments and instruments containing Nellcor oximetry, or with instruments licensed to use Nellcor OxiMax sensors (Nellcor-compatible instruments).**
2. **If the sensor is misapplied with excessive pressure for prolonged periods, a pressure injury can occur.**
3. **Do not use the MAXA(L) or other oximetry sensors during MRI scanning. Conducted current may cause burns. Also, the MAXA(L) may affect the MRI image, and the MRI unit may affect the accuracy of oximetry measurements.**
4. **As with all medical equipment, carefully route cables to reduce the possibility of patient entanglement or strangulation.**

CAUTIONS

1. In the event of damage to the sterile packaging, do NOT resterilize. Follow local governing ordinances and recycling instructions regarding disposal or recycling of sensors.
2. Do not use a damaged sensor or pulse oximetry cable. Do not use a sensor with exposed optical components.
3. Do not immerse in water or cleaning solutions. Do not resterilize. Such sterilization could damage the sensor and may result in sensor malfunction and/or erroneous oximetry measurements.
4. Failure to apply the MAXA(L) properly may cause incorrect measurements.
5. If the sensor is wrapped too tightly or supplemental tape is applied, venous pulsations may lead to inaccurate saturation measurements.
6. While the MAXA(L) is designed to reduce the effects of ambient light, excessive light may cause inaccurate measurements. In such cases, cover the sensor with opaque material.
7. Circulation distal to the sensor site should be checked routinely. The site must be inspected every 8 hours to ensure adhesion, application pressure, skin integrity, and correct optical alignment. If skin integrity changes, move the sensor to another site. If the sensor is misapplied with excessive pressure, a pressure injury can occur.
8. Deeply pigmented skin, intravascular dyes or externally applied coloring (such as nail polish, dye or pigmented cream) may lead to inaccurate measurements.
9. Excessive motion may compromise performance. In such cases, try to keep the patient still, or change the sensor site to one with less motion.
10. Do not alter or modify the MAXA(L). Alterations or modifications may affect performance or accuracy.
11. For additional warnings, cautions or contraindications when using this sensor with Nellcor-compatible instruments, refer to the instrument operator's manual or contact the manufacturer of the instrument.

Accuracy Specifications

This sensor integrates Nellcor OxiMax technology into its design. When connected to an OxiMax-enabled instrument, this sensor uses OxiMax technology to provide additional advanced sensor performance features. Consult individual manufacturers for features and compatibility of particular instruments and sensor models.

Each Nellcor OxiMax-compatible instrument manufacturer is responsible for determining optimal compatibility conditions and settings for its instruments to provide safe and effective use of each Nellcor OxiMax sensor. This includes specifications and/or warnings, cautions, or contraindications. Refer to each instrument operator's manual or consult manufacturer for complete instructions, warnings, cautions, or contraindications on the use of this sensor with their Nellcor OxiMax compatible instrument.

Sensor Saturation Accuracy

Low Perfusion	70 to 100% $\pm$ 2 digits
Adult	70 to 100% $\pm$ 2 digits
Adult and Neonate with Motion	70 to 100% $\pm$ 3 digits
*as tested with PM1000N monitor	

Refer to each instrument operator's manual for clinical testing summary and clinical study Bland-Altman plots.

Optical Specifications

LED wavelength (Approximate)	Output power
Red: 650 to 670 nm	3.0 mW
Infrared: 880 to 910 nm	4.0 mW

Wavelength range can be especially useful to clinicians.

Environmental Specifications

Operational temperature range	0°C to 50°C
Storage temperature range	-40°C to 70°C
Operational and Storage Relative humidity range	15% to 95% non condensing

Additional Copies of Instructions

Additional copies of these Instructions are available at no charge by calling Covidien or its authorized distributors. Also, permission is hereby granted under Covidien copyrights to purchasers of products obtained from Covidien or its authorized distributors to make additional copies of these instructions for use by such purchasers.

Any other use of this sensor is not authorized by Covidien under any patents.

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registered trademarks of Covidien AG.
Other brands are trademarks of a
EC REP Covidien Ireland Limited,