PadPro[®] Multi-Function Electrode Adapter IFU

ZMUF

Instructions for Use:

Multi-Function Electrode (MFE) Adapters for Defibrillation/Pacing/Cardioversion/Monitoring

Read and become familiar with all instructions, warnings and cautions before using this product.


Indications: The CONMED PadPro[®] are indicated for use by trained medical professionals in medical facilities or medical transport environments to adapt connection systems associated with a specific defibrillator/therapy cable to a different style connection system. The PadPro[®] MFE adapters are intended for delivery of energy for defibrillation, cardioversion, external pacing, and ECG monitoring applications. The PadPro[®] MFE adapter is a non-sterile, reusable device, providing conductive interface between the defibrillator and/or therapy cable and MFE electrode. This device is intended for use on defibrillators whose output is classified as low power (up to 360 joule maximum).

WARNING: DO NOT use PadPro[®] MFE adapters with any model of defibrillator or therapy cable that are not listed on the following table. Defibrillator models that are not listed have not been qualified for use with CONMED MFE adapters and MFE performance may be impacted.

Compatibility:


- All PadPro® MFE Adapters are compatible with biphasic defibrillators.
- These PadPro® MFE Adapters have been qualified with the following defibrillators and PadPro® MFEs Pads:

Figure 1 - Assembly Diagram

PadPro® MFE Adapter catalog number consists of xM (Male Connector Style) + xF (Female Connector Style)

Table 1 - Standardization Solution Options

<u>MFE Therapy Cable Connection (Male)</u>			<u>MFE Pad Connection (Female)</u>	
<u>Compatible Defibrillator</u>	<u>Compatible Therapy Cable</u>	<u>MFE Adapter Male Connector Style</u>	<u>MFE Adapter Female Connector Style</u>	<u>Compatible PadPro® MFE Pads</u>
Zoll R Series R Series Plus	8009-0749 8009-0750	<u>ZM - Zoll</u> 	<u>UF Universal/ Anderson</u> 	2001, 2001-C, 2001-EPS, 2001-PC, 2502, 2516, 2516-PC, 2603, 2603-PC, 2602
Zoll X Series Propaq MD	8300-0783			

Warnings:

- Misuse of PadPro® MFE adapters, use of compromised or altered product, and/or failure to follow product instructions may result in patient burns, inadequate delivery of therapy, and/or loss of ECG trace quality.
- Defibrillation current can cause operator or bystander injury. Do not touch patient, or equipment connected to the patient, during defibrillation.
- Do not use PadPro® MFE adapters in the presence of flammable agents or in an oxygen enriched atmosphere as this presents an explosion and fire hazard.
- No modification of this device is allowed.
- Do not use beyond 200 defibrillations or 2 years' service life after initial use.

Storage: Store PadPro® MFE Adapter in dry, room temperature environment.

Cleaning Information:

- PadPro® MFE Adapters may be cleaned by wiping with disinfectant wipe or cloth soaked in isopropyl alcohol.

WARNING: Do not submerge, soak, autoclave, or sterilize PadPro® MFE Adapters. Device is intended to be provided non-sterile.

Instructions for Use: (all models)

WARNING: Before each use, inspect the condition of the adapter for cracks, bent pins, or loose parts. If the adapter appears to be damaged, discard and replace with a new adapter.

- 1.) Install the PadPro® MFE Adapter to compatible therapy cable by first connecting the tether, and then inserting the male end of the MFE Adapter into the therapy cable until a click or other tactile feel (i.e. end stop) is detected.
- 2.) Connect the compatible PadPro® series electrode to the female end of the PadPro® MFE Adapter until a click or other tactile feel (i.e. end stop) is detected. After connection of electrodes, follow defibrillator instructions or hospital protocol.
- 3.) Disconnect the single use PadPro® electrode after use. When disconnecting the electrodes, be sure to leave the adapter properly inserted into the therapy cable or future use (and/or testing).

WARNING: PadPro® MFE Adapters are reusable devices. Do not dispose of PadPro® MFE Adapter during disposal of single use PadPro® series MFE pad.

- 4.) Clean PadPro® MFE Adapter after each use as instructed in “Cleaning Instructions” section of this Instructions for Use.
- 5.) Inspection and Test:
 - When testing the defibrillator through the therapy cable, make sure the PadPro® MFE adapter is connected and testing into the original manufacturer’s or CONMED provided test load/ shorting block.
 - Refer to the defibrillator Operator’s Manual for additional information regarding the test procedure.
 - If the adapter fails during the routine test, discard and replace with a new adapter.

Disposal: After use, this product may be a potential biohazard. Handle and dispose of in accordance with institution protocol and local, state and federal laws.

Notice: Any serious incident that has occurred using this device should be reported to CONMED Corporation and the competent authority of the Member State in which the user and/or patient is established.

Cautions:

- Refer to the defibrillator user's guide for overall operating procedures.
- Federal law (USA) restricts this device to sale by or on the order of a physician.
- Do not use in the presence of a Magnetic Resonance Imaging (MRI) device.
- Materials incorporated in the device could result in sensitization or allergic reaction to patient or end user.
- This product may be used by authorized medical personnel only.

	Operating Environmental Limitations		Defibrillation-Proof Type CF Applied Part
	Storage Environmental Limitations		Transportation Environmental Limitations

Made in USA for:

CONMED CORPORATION

525 French Road

Utica, New York 13502-5994 USA

(315)797-8375 • FAX 1-800-438-3051

Customer Service (USA) 1-866-426-6633

e-mail: CustomerExperience@conmed.com

© CONMED Corporation

MDSS GmbH

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D-30175 Hannover, Germany

For International orders or inquiries, please contact

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