GE Medical Systems

Technical Publications

sm 2119703–100

Revision 9

MP Phase 4 + Vascular Option Generator Schematics

do not duplicate

Copyright © 1996, 1997, 1998 1999, 2000 by General Electric Co.

ATTENTION

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WARNING

**X-RAY EQUIPMENT IS DANGEROUS TO BOTH PATIENT AND OPERATOR
UNLESS MEASURES OF PROTECTION ARE STRICTLY OBSERVED**

Though this equipment is built to the highest standards of electrical and mechanical safety, the useful x-ray beam becomes a source of danger in the hands of the unauthorized or unqualified operator. Excessive exposure to x-radiation causes damage to human tissue.

Therefore, adequate precautions must be taken to prevent unauthorized or unqualified persons from operating this equipment or exposing themselves or others to its radiation.

Before operation, persons qualified and authorized to operate this equipment should be familiar with the Recommendations of the International Commission on Radiological Protection, contained in Annals Number 26 of the ICRP, and with applicable national standards.

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VOLTAGE PRESENT

Before any intervention:

1. Switch **OFF** CB2 on the AC Power Distribution (VPE1 A2 A25) in Positioner Cabinet.
2. Switch **OFF** Positioner Cabinet Circuit Breaker in Generator Power Cabinet.

LOCK OUT/ TAG OUT

When necessary, perform the following lock out/ tag out procedure:

- 1.. Cut the customer power distribution cabinet from its energy source.
- 2.. Inform all concerned parties of the lockout.
- 3.. Lock and put a lockout label on the cabinet.
- 4.. Verify the absence of voltage in the cabinet.
- 5.. Proceed with service to the equipment.
- 6.. Unlock the cabinet and remove the label from the cabinet.
- 7.. Power-up the system.

Direction 2119703–100

Revision 9

MP PHASE 4 Generator
Installation

IMPORTANT! ... X–RAY PROTECTION

X–RAY equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood before you attempt to place this equipment in operation. The General Electric Company, Medical Systems Group, will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x–radiation other than the useful beam. No practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take

adequate precautions to prevent the possibility of any persons carelessly, unwisely, or unknowingly exposing themselves or others to radiation.

It is important that everyone having anything to do with x–radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814 and of the International Commission on Radiation Protection, and take adequate steps to insure protection against injury.

All persons authorized to use the equipment must be cognizant of the danger of excessive exposure to x–radiation and the equipment is sold with the understanding that the General Electric Company, Medical Systems Group, its agents, and representatives have no responsibility for injury or damage which may result from exposure to x–radiation.

Various protective material and devices are available. It is urged that such materials or devices be used.

Blank page

WARNING

- THIS SERVICE MANUAL IS AVAILABLE IN ENGLISH ONLY.
- IF A CUSTOMER’S SERVICE PROVIDER REQUIRES A LANGUAGE OTHER THAN ENGLISH, IT IS THE CUSTOMER’S RESPONSIBILITY TO PROVIDE TRANSLATION SERVICES.
- DO NOT ATTEMPT TO SERVICE THE EQUIPMENT UNLESS THIS SERVICE MANUAL HAS BEEN CONSULTED AND IS UNDERSTOOD.
- FAILURE TO HEED THIS WARNING MAY RESULT IN INJURY TO THE SERVICE PROVIDER, OPERATOR OR PATIENT FROM ELECTRIC SHOCK, MECHANICAL OR OTHER HAZARDS.

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REVISION HISTORY

REV	DATE	REASON FOR CHANGE
0	Jul. 17, 1991	Schematics for production.
1	Oct. 31, 1991	Review of rev.0.
2	Mar. 04, 1992	Updated.
3	Jun. 04, 1992	
0	May 31, 1993	MP Phase 4 (Ph. 2). New Part No. 36004389 for sm.
0	July 15, 1995	MP Phase 4 (Ph. 4) + Vascular Option Generator. New Part No. 2119703–100 for Schematics (old Part No.: 36004385).
1	31–Jan–96	Updated (Technical Release 2) New schematics for: Measure and Commutation Board, Mini Bridge Board, Regulation Board, Control Board, Interface 2 Board, Grid/Bias Interface Board, AC Distribution Panel and Filter Board.
2	14–Jun–96	Updated (SPR BUCge17973) New schematic for CPU Board.
3	23–Oct–96	Schematic Update.
4	3–Apr–97	Schematic Update: Interface I Board and Measurement and Commutation
5	2–Dec–97	Update.
6	06–Mar–98	Updates for Butterfly: Modification to Schematic 2115419. Update Schematic 2127280 to revB.
7	09–Apr–99	Updates for NF2: Modification to Schematic 2140213–2. Update Schematic 2140213–2 to revA.
8	09–Feb.–2000	New Security Board. Updated Inverters 1 & 2 Block Diagrams.
9	25–Sept–2001	External ECM Power Supply updated on CPU Board schematic. (page 205)

LIST OF EFFECTIVE PAGES

PAGE NUMBER	REVISION NUMBER	PAGE NUMBER	REVISION NUMBER	PAGE NUMBER	REVISION NUMBER	PAGE NUMBER	REVISION NUMBER	PAGE NUMBER	REVISION NUMBER
Title page	9								
Safety Instructions	9								
i thru x	9								
1 thru 250	9								

Blank page

WARNING

**LINE VOLTAGE IS PRESENT ON MP PHASE 42 A1 A4
BOARD
WHEN MAIN CIRCUIT BREAKER IS "OFF"**

REMARK

There is a s/w protection circuit which prevents frequent "ON–OFF" cycling of the MP. The MP can be cycled "ON" and "OFF" only 7 times in a 21 minutes time period. Every 3 minutes one of the switch cycle counts is removed. (e.g. you could cycle indefinitely if you wait 3 minutes between each "ON–OFF" cycle).

This delay is required to protect certain components in the AC Supply board (MP Phase 41 A4 A1), from the effects of capacitor charge and discharge.

Blank page

CONTENTS


REFERENCE DESIGNATOR	TITLE	REF Nb.	PAGE
MP PHASE 41	MP PHASE 4 <i>OVERVIEW</i>	2131108ADW A	1
	MP PHASE 4 <i>GENERAL BLOCK DIAGRAM</i>		3
	MP PHASE 4 <i>CABINETS INTERCONNECTION</i>	2131241WIR B	5
A1	HIGH VOLTAGE TANK <i>DOOR DIAGRAM</i>	2129739PSP B	7
	SWITCHING TANK <i>SCHEMATIC</i>	45312098:030 C	9
	HV TANK <i>BLOCK DIAGRAM</i>		11
	” <i>SCHEMATIC</i>	45204019:030 D	13
A1	MEASURE AND COMMUTATION BOARD <i>MNEMONICS</i>		15
	” ” <i>BLOCK DIAGRAM</i>		17
	” ” <i>SCHEMATIC</i>	2164152SCH B	19
A2	HV DIVIDER BOARD <i>SCHEMATIC</i>	36004682:030 A	27
A3	QUENCH BOARD <i>SCHEMATIC</i>	2128181SCH A	29
	” <i>DOOR DIAGRAM</i>	SEE MP PHASE 41 A1	
A2	INVERTER 1 <i>PANEL DOOR DIAGRAM</i>	2131245PSP B	31
	” <i>BLOCK DIAGRAM</i>		33
	” <i>SCHEMATIC</i>	2131097SCH B	35
A1	GATE 1 BOARD <i>SCHEMATIC</i>	45203847:030 B	37
A2	CURRENT TRANSFORMER 1 BOARD	SEE MP PHASE 41 A2	
A3	MINI BRIDGE BOARD <i>SCHEMATIC</i>	2140213SCH A	39
	” <i>DOOR DIAGRAM</i>	SEE MP PHASE 41 A1	
A3	INVERTER 2 <i>BLOCK DIAGRAM</i>		41
	” <i>DOOR DIAGRAM</i>	2131246PSP B	43
	” <i>SCHEMATIC</i>	45312335:030 C	45
A1	GATE 2 BOARD	SEE MP PHASE 41 A2	
A2	CURRENT TRANSFORMER 2 BOARD	SEE MP PHASE 41 A2	
A4	AC SUPPLY PANEL <i>DOOR DIAGRAM</i>	36004532:006	47
	” <i>BLOCK DIAGRAM</i>		49
	” <i>SCHEMATIC</i>	45203920:030 E	51
A1	AC SUPPLY BOARD <i>SCHEMATIC</i>	45203921:030 H	53
A2	ON/OFF SERVICE BOARD <i>SCHEMATIC</i>	45561065:030 B	59
A5	CONTROL CARD RACK <i>DOOR DIAGRAM</i>	2131243PSP C	63
	” ” <i>SCHEMATIC</i>	2127663SCH A	65
	” ” <i>SCHEMATIC</i>	2131111SCH C	85
A1	REGULATION BOARD <i>MNEMONICS</i>		89
	” ” <i>BLOCK DIAGRAM</i>		91
	” ” <i>SCHEMATIC</i>	2138360SCH A	93
	” ” <i>SCHEMATICS SYMBOLS</i>		105


CONTENTS

REFERENCE DESIGNATOR	TITLE	REF Nb.	PAGE
MP PHASE 41			
A5			
A2	CONTROL BOARD <i>MNEMONICS</i>		107
	” ” <i>BLOCK DIAGRAM</i>		109
	” ” <i>SCHEMATIC</i>	2139533SCH A	111
	” ” <i>SCHEMATICS SYMBOLS</i>		119
A3	SECURITY BOARD <i>MNEMONICS</i>		121
	” ” <i>BLOCK DIAGRAM</i>		123
	” ” <i>SCHEMATIC</i>	45475108–2SCH R0	125
	” ” <i>SCHEMATICS SYMBOLS</i>		137
A5	INTERFACE 1 BOARD <i>MNEMONICS</i>		139
	” ” <i>BLOCK DIAGRAM</i>		141
	” ” <i>SCHEMATIC</i>	36004541:030 F	143
A6	INTERFACE 2 BOARD <i>MNEMONICS</i>		167
	” ” <i>BLOCK DIAGRAM</i>		169
	” ” <i>SCHEMATIC</i>	2138744SCH A	171
A9	CPU BOARD <i>BLOCK DIAGRAM</i>		193
	” <i>SCHEMATIC</i>	1 0235 S P E	195
A10	HEATER SUPPLY BOARD <i>MNEMONICS</i>		207
	” ” <i>BLOCK DIAGRAM</i>		209
	” ” <i>SCHEMATIC</i>	45562521:030 B	211
A11	HEATER SUPPLY BOARD	SEE MP PHASE 41 A5 A10	
A12	AUXILIARY VOLTAGE BOARD <i>BLOCK DIAGRAM</i>		217
	” ” <i>SCHEMATIC</i>	45563250:030 A	219
A13	GRID/BIAS INTERFACE BOARD <i>MNEMONICS</i>		221
	” ” <i>BLOCK DIAGRAM</i>		223
	” ” <i>SCHEMATIC</i>	2137893SCH A	225
A20	LV POWER SUPPLY	SEE MP PHASE 41 A5	
A6	STABILIZER PANEL	SEE MP PHASE 41 A5	
A7	CABINET FANS	SEE MP PHASE 41	


CONTENTS


REFERENCE DESIGNATOR	TITLE	REF Nb.	PAGE
MP PHASE 42	MP PHASE 4 POWER UNIT MAINS CABINET	SEE MP PHASE 41	
A1	AC DISTRIBUTION PANEL <i>DOOR DIAGRAM</i>	2115966PSP B	235
	” ” <i>SCHEMATIC</i>	2115419SCH J	237
A1	VOLTAGE SELECTION PANEL	SEE MP PHASE 42 A1	
A2	SYSTEM ON/OFF SWITCH	SEE MP PHASE 42 A1	
A3	CIRCUIT BREAKER PANEL	SEE MP PHASE 42 A1	
A4	MP PHASE 4 POWER BOARD <i>SCHEMATIC</i>	2101740SCH A	239
A2	TIRC <i>BLOCK DIAGRAM</i>		241
A3	CABINET FANS	SEE MP PHASE 41	
A4	TUBE CONNECTOR PANEL	SEE MP PHASE 41	
A5	FILTER BOARD <i>SCHEMATIC</i>	2138132SCH A	243
MP PHASE 43	3rd TUBE OPTION <i>SCHEMATIC</i>	36004448:030 A	245
MP PHASE 44	3rd/4th TUBE OPTION <i>SCHEMATIC</i>	36004449:030 A	247
MP PHASE 45	GRID/BIAS TANK <i>SCHEMATIC</i>	2127280SCH B	249



Blank page

Blank page

Blank page



Blank page


MNEMONICS		3-0
MP PHASE 41 A1 A1 – MEASURE & COMMUTATION BOARD		
XS	Small Focus	
XL	Large Focus	
kV+	High voltage anode	
kV-	High voltage cathode	
Test kV-mA		
SW1FDBK	Switch 1 Feedback	
SW2FDBK	Switch 2 Feedback	
SW3FDBK	Switch 3 Feedback	
SW1/2/3 T	High voltage switch 1/2/3 Feedback to Interface 1 Board	
SW1/2/3 R	High voltage switch 1/2/3 Feedback to Interface 1 Board	


Blank page

3-1




H	BLOCK PATHNAME										\$PROJECT_VAMP/mesplu_v4_ASM/mesplu_v4										BLOCK SHT										3										OF										8									
	REV 1		PCN		XXXXXXXX		DESCRIPTION				XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX																MADE		ASSE		DATE		25-06-96		GE Medical Systems										DESIGN TITLE				MEASURE AND COMMUTATION											
	MADE BY		NAME				DATE		DD/MM/YY		APPR BY		NAME		DATE		DD/MM/YY		CHECK		DURVILLE		DATE		25-06-96		FIRST MADE FOR				VAMP																													
	REV 2		PCN		XXXXXXXX		DESCRIPTION				XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX																APPR		ENCELLAZ		DATE		25-06-96		BUC - FRANCE										DRAWING NO.				REV		SIZE		SHT/CONT. ON							
	MADE BY		NAME				DATE		DD/MM/YY		APPR BY		NAME		DATE		DD/MM/YY		OTHER				DATE				2164152SCH				B		3												3 / 8															


BLOCK PATHNAME				\$PROJECT VAMP/mesplu_v4 ASM/mesplu_v4				BLOCK SHT 4 OF 8			
REV 1	PCN	XXXXXXXX	DESCRIPTION	XXXXXXXXXXXXXXXXXXXXXXXXXXXX	MADE	ASSE	DATE	25-06-96	GE Medical Systems		
MADE BY	NAME		DATE DD/MM/YY	APPR BY NAME	CHECK	DURVILLE	DATE	25-06-96			
REV 2	PCN	XXXXXXXX	DESCRIPTION	XXXXXXXXXXXXXXXXXXXXXXXXXXXX	APPR	ENCELLAZ	DATE	25-06-96	BUC - FRANCE		
MADE BY	NAME		DATE DD/MM/YY	APPR BY NAME	OTHER		DATE				
									DESIGN TITLE MEASURE AND COMMUTATION		
									FIRST MADE FOR VAMP		
									DRAWING NO. 2164152SCH		
									REV B	SIZE 3	SHT/CONT. ON 4 / 8


BLOCK PATHNAME										\$PROJECT_VAMP/mesplu_v4 ASM/mesplu_v4										BLOCK SHT				5	OF	8															
REV 1		PCN		XXXXXXXX		DESCRIPTION				XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX										MADE		ASSE		DATE		25-06-96		GE Medical Systems				DESIGN TITLE				MEASURE AND COMMUTATION					
MADE BY		NAME				DATE		DD/MM/YY		APPR BY		NAME		DATE		DD/MM/YY		CHECK		DURVILLE		DATE		25-06-96		FIRST MADE FOR						VAMP									
REV 2		PCN		XXXXXXXX		DESCRIPTION				XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX										APPR		ENCELLAZ		DATE		25-06-96		BUC - FRANCE				DRAWING NO.				REV		SIZE		SHT/CONT. ON	
MADE BY		NAME				DATE		DD/MM/YY		APPR BY		NAME		DATE		DD/MM/YY		OTHER				DATE				2164152SCH						B		3		5 / 8					


	1	2	3	4	5	6	7	8	9		10	SHT.	8	CONT.	8																																																																										
	2164152SCH																																																																																								
A	COMPOSANTS A SUPPRIMER :															A																																																																									
	<table><tr><th>Page</th><th>Repere</th><th>DESIGNATION</th><th>Quantite</th></tr><tr><td rowspan="9">3/8</td><td rowspan="2">X11,X9</td><td>RELAIS</td><td rowspan="2">2</td></tr><tr><td>SP4-P-DC-12V</td></tr><tr><td rowspan="2">X13</td><td>RELAIS</td><td rowspan="2">1</td></tr><tr><td>RP-820-012</td></tr><tr><td>D20,D28,D34</td><td>DIODES IN4448</td><td>3</td></tr><tr><td>R105,R73</td><td>100 OHMS</td><td>2</td></tr><tr><td>R98,R72</td><td>1 KOHMS</td><td>2</td></tr><tr><td>R106</td><td>66,5 OHMS</td><td>1</td></tr><tr><td>DS9,DS10</td><td>HLMP 0503</td><td>2</td></tr><tr><td rowspan="9">6/8</td><td colspan="2">XJ15,XJ16</td><td>Prises BURNDY</td><td>2</td></tr><tr><td colspan="2">LS9,LS10</td><td>V250 LA20 A</td><td>2</td></tr><tr><td colspan="2">XJ11,XJ12</td><td>Prises BURNDY</td><td>2</td></tr><tr><td rowspan="2">7/8</td><td>XJY1,ZJZ1</td><td>Prises BURNDY</td><td>2</td></tr><tr><td>XJZ1</td><td>BOUCHON</td><td>1</td></tr><tr><td rowspan="2"></td><td></td><td>36 004 536</td><td rowspan="2"></td></tr><tr><td>XJY1</td><td>BOUCHON</td></tr><tr><td rowspan="2"></td><td></td><td>36 004 537</td><td></td></tr><tr><td>R126,R125 R124,R123</td><td>Resistances 10 OHMS</td><td>4</td></tr><tr><td rowspan="2">F</td><td rowspan="2"></td><td>C66,C67 C68,C69</td><td>Ccondensateurs 10 NF 63V</td><td>4</td></tr><tr><td>D39,D40 D41,D42</td><td>DIODES IN5908</td><td>4</td></tr></table>															Page	Repere	DESIGNATION	Quantite	3/8	X11,X9	RELAIS	2	SP4-P-DC-12V	X13	RELAIS	1	RP-820-012	D20,D28,D34	DIODES IN4448	3	R105,R73	100 OHMS	2	R98,R72	1 KOHMS	2	R106	66,5 OHMS	1	DS9,DS10	HLMP 0503	2	6/8	XJ15,XJ16		Prises BURNDY	2	LS9,LS10		V250 LA20 A	2	XJ11,XJ12		Prises BURNDY	2	7/8	XJY1,ZJZ1	Prises BURNDY	2	XJZ1	BOUCHON	1			36 004 536		XJY1	BOUCHON			36 004 537		R126,R125 R124,R123	Resistances 10 OHMS	4	F		C66,C67 C68,C69	Ccondensateurs 10 NF 63V	4	D39,D40 D41,D42	DIODES IN5908	4					
Page	Repere	DESIGNATION	Quantite																																																																																						
3/8	X11,X9	RELAIS	2																																																																																						
		SP4-P-DC-12V																																																																																							
	X13	RELAIS	1																																																																																						
		RP-820-012																																																																																							
	D20,D28,D34	DIODES IN4448	3																																																																																						
	R105,R73	100 OHMS	2																																																																																						
	R98,R72	1 KOHMS	2																																																																																						
	R106	66,5 OHMS	1																																																																																						
	DS9,DS10	HLMP 0503	2																																																																																						
6/8	XJ15,XJ16		Prises BURNDY	2																																																																																					
	LS9,LS10		V250 LA20 A	2																																																																																					
	XJ11,XJ12		Prises BURNDY	2																																																																																					
	7/8	XJY1,ZJZ1	Prises BURNDY	2																																																																																					
		XJZ1	BOUCHON	1																																																																																					
			36 004 536																																																																																						
		XJY1	BOUCHON																																																																																						
			36 004 537																																																																																						
		R126,R125 R124,R123	Resistances 10 OHMS	4																																																																																					
F		C66,C67 C68,C69	Ccondensateurs 10 NF 63V	4																																																																																					
		D39,D40 D41,D42	DIODES IN5908	4																																																																																					
	FAIRE UN STRAP ENTRE LES CONTACTS 5 ET 6 DE X9 FAIRE UN STRAP ENTRE LES CONTACTS 11 ET 12 DE X9 FAIRE UN STRAP ENTRE LA BROCHE 34 DE XJ7 ET LE OV FAIRE UN STRAP ENTRE LA BROCHE 36 DE XJ7 ET LE OV																																																																																								
H	<table><tr><td colspan="5">BLOCK PATHNAME \$PROJECT VAMP/mesplu_v4 ASM/mesplu_v4</td><td colspan="5">BLOCK SHT 900 OF 8</td><td colspan="5"></td></tr><tr><td>REV 1</td><td>PCN</td><td>XXXXXXXX</td><td>DESCRIPTION</td><td>XXXXXXXXXXXXXXXXXXXXXXXXXXXX</td><td>MADE</td><td>ASSE</td><td>DATE</td><td>25-06-96</td><td colspan="2" rowspan="2">GE Medical Systems</td><td colspan="4">DESIGN TITLE MEASURE AND COMMUTATION</td></tr><tr><td>MADE BY</td><td>NAME</td><td></td><td>DATE DD/MM/YY</td><td>APPR BY NAME</td><td>CHECK</td><td>DURVILLE</td><td>DATE</td><td>25-06-96</td><td colspan="4">FIRST MADE FOR VAMP</td></tr><tr><td>REV 2</td><td>PCN</td><td>XXXXXXXX</td><td>DESCRIPTION</td><td>XXXXXXXXXXXXXXXXXXXXXXXXXXXX</td><td>APPR</td><td>ENCCELLAZ</td><td>DATE</td><td>25-06-96</td><td colspan="2" rowspan="2">BUC - FRANCE</td><td colspan="2">DRAWING NO.</td><td>REV</td><td>SIZE</td><td>SHT/CONT. ON</td></tr><tr><td>MADE BY</td><td>NAME</td><td></td><td>DATE DD/MM/YY</td><td>APPR BY NAME</td><td>OTHER</td><td>....</td><td>DATE</td><td>.....</td><td colspan="2">2164152SCH</td><td>B</td><td>3</td><td>8 / 8</td></tr></table>															BLOCK PATHNAME \$PROJECT VAMP/mesplu_v4 ASM/mesplu_v4					BLOCK SHT 900 OF 8										REV 1	PCN	XXXXXXXX	DESCRIPTION	XXXXXXXXXXXXXXXXXXXXXXXXXXXX	MADE	ASSE	DATE	25-06-96	GE Medical Systems		DESIGN TITLE MEASURE AND COMMUTATION				MADE BY	NAME		DATE DD/MM/YY	APPR BY NAME	CHECK	DURVILLE	DATE	25-06-96	FIRST MADE FOR VAMP				REV 2	PCN	XXXXXXXX	DESCRIPTION	XXXXXXXXXXXXXXXXXXXXXXXXXXXX	APPR	ENCCELLAZ	DATE	25-06-96	BUC - FRANCE		DRAWING NO.		REV	SIZE	SHT/CONT. ON	MADE BY	NAME		DATE DD/MM/YY	APPR BY NAME	OTHER		DATE		2164152SCH		B	3	8 / 8	H
BLOCK PATHNAME \$PROJECT VAMP/mesplu_v4 ASM/mesplu_v4					BLOCK SHT 900 OF 8																																																																																				
REV 1	PCN	XXXXXXXX	DESCRIPTION	XXXXXXXXXXXXXXXXXXXXXXXXXXXX	MADE	ASSE	DATE	25-06-96	GE Medical Systems		DESIGN TITLE MEASURE AND COMMUTATION																																																																														
MADE BY	NAME		DATE DD/MM/YY	APPR BY NAME	CHECK	DURVILLE	DATE	25-06-96			FIRST MADE FOR VAMP																																																																														
REV 2	PCN	XXXXXXXX	DESCRIPTION	XXXXXXXXXXXXXXXXXXXXXXXXXXXX	APPR	ENCCELLAZ	DATE	25-06-96	BUC - FRANCE		DRAWING NO.		REV	SIZE	SHT/CONT. ON																																																																										
MADE BY	NAME		DATE DD/MM/YY	APPR BY NAME	OTHER		DATE				2164152SCH		B	3	8 / 8																																																																										
	1	2	3	4	5	6	7	8	9	10	11																																																																														

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
5-1

MP Phase 41 A1 A3


BLOCK PATHNAME \$PROJECT VAMP/quench_ASM/quench				BLOCK SHT 01 OF 01			
REV 1	PCN XXXXXXXX	DESCRIPTION XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX	DATE DD/MM/YY	MADE CHECK APPR OTHER	ASSE MENARD JP	DATE 19-01-95	DESIGN TITLE QUENCH
MADE BY NAME	DATE DD/MM/YY	APPR BY NAME	DATE DD/MM/YY	NAME	NAME	DATE DD/MM/YY	FIRST MADE FOR VAMP
REV 2	PCN XXXXXXXX	DESCRIPTION XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX	DATE DD/MM/YY	NAME	NAME	DATE DD/MM/YY	DRAWING NO. 2128181SCH
MADE BY NAME	DATE DD/MM/YY	APPR BY NAME	DATE DD/MM/YY	NAME	NAME	DATE DD/MM/YY	REV A SIZE 3 SHT/CONT 01 / 01 ON

Blank page

MP PHASE 41 A3 – INVERTER 2 BLOCK DIAGRAM


6-2



Blank page

Blank page



MP PHASE 41 A2 – INVERTER 1 BLOCK DIAGRAM


9-2

Blank page

Blank page


Blank page

Blank page

MP PHASE 41 A4 – AC SUPPLY PANEL *BLOCK DIAGRAM*


10-2

Blank page


Blank page


Blank page


Blank page


Blank page

61


Blank page


Blank page

Blank page


MP Phase 41 A5


Blank page

Blank page

Blank page


Blank page

2127664A.P07


Blank page


Blank page

Blank page

Blank page

MNEMONICS

13-0


MP PHASE 41 A5 A1 – REGULATION BOARD

VERR	Voltage Error
VF	Input Voltage of Voltage/Frequency Converter
STSADC	Analog/Digital converter Status
FREQ	V/F Frequency Output
kV < kV CORRECTED	kV value measured < kV corrected value
MB ref.	Mini bridge kV ref.
kV > 0.75 kV REF.	kV value measured > 0.75 kV ref. value
kV > kV THRESHOLD	kV value measured > kV threshold value
kV DOWN	Spit detection
kV MEAS.	kV value measurement
kV-ACQ-EN	kV measure reset
AUTO HB	Half bridge enable


Blank page

MP PHASE 41 A5 A1 – REGULATION BOARD *BLOCK DIAGRAM*


13-1

Blank page

Blank page

Blank page

Blank page

Blank page

Blank page


Blank page


Blank page


MNEMONICS		14-0
MP PHASE 41 A5 A2 – CONTROL BOARD		
STQ	Sense TQ (checks operating enable with TQ min interlock) TQ min actuated.	
TQ	Time between thyristor current cancel and voltage applied to its terminals.	
HB CMD	Half Bridge Command	
AUTO HB	Half-bridge enable – enables half-bridge connection when KV threshold is reactivated	
P INV1	Pulse Inverter 1	
P INV2	Pulse Inverter 2	
RFB1 RFB2	Return Full Bridge Inverter 1 or 2	
RHB1 RHB2	Return Half Bridge Inverter 1 or 2	
SYNC	Signal enabled with the first inverter pulse and remains active during inverter running	
MB CDE	Mini bridge command	
MB EN	Mini bridge enable (actived in grided mode)	
FMAX EN	Frequency limitation change	
MB CTRL	Mini bridge activated	
QUENCH CTRL	Quench mode activated	

Blank page

Blank page

Blank page

Blank page

Blank page

Blank page

MNEMONICS


15–0


MP PHASE 41 A5 A3 – SECURITY BOARD

TQ EN	Enable given by TQ min interlock
AUTO HB	Half bridge enable – enables half–bridge connection when KV threshold is reactivated
INV1 CMD	Inverter 1 Command
INV2 CMD	Inverter 2 Command
HB CMD	Half bridge command
SECURITY	One secutity activated
RFB1	Feedback Full bridge 1
RFB2	Feedback Full bridge 2
RHB1	Feedback Half bridge 1
RHB2	Feedback Half bridge 2
I INV 1	Current into inverter 1
I INV 2	Current into inverter 2
E > 380 V	Inverter dc voltage > 380 V
E > 440 V	Inverter dc voltage > 440 V
E > 500 V	Inverter dc voltage > 500 V
E > 560 V	Inverter dc voltage > 560 V
STQ	Thyristor tq limitation

MP PHASE 41 A5 A3 – SECURITY BOARD BLOCK DIAGRAM


15-1


475108-2.PL1


MPPU1-A5-A3


Blank page

Blank page


Blank page

Blank page


Blank page


Blank page


MNEMONICS		16-0
MP PHASE 41 A5 A5 – INTERFACE 1 BOARD		
SW1 R	High voltage switch 1 Feedback	
SW1 T	High voltage switch 1 Feedback	
SW2 R	High voltage switch 2 Feedback	
SW2 T	High voltage switch 2 Feedback	
SW3 R	High voltage switch 3 Feedback	
SW3 T	High voltage switch 3 Feedback	

MP PHASE 41 A5 A5 – INTERFACE 1 BOARD BLOCK
DIAGRAM

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









17-0


MP PHASE 41 A5 A6 – INTERFACE 2 BOARD

MB FDBKA	Mini bridge operation Feedback
GRID FUNCT OK	Grid function OK
GRID 0VLD	Grid current overload
GRID SWITCH OK	36 V grid switch OK
GRID PULSE OK	Feedback to Int2 board
BIAS FUNCT OK	Bias voltage OK
BIAS ON	Bias voltage Feedback to Int 2 Board

MP PHASE 41 A6 – INTERFACE 2 BOARD BLOCK DIAGRAM


17-1


BLOCK PATHNAME \$PROJECT_VAMP/int2_v8_ASM/int2_v8					BLOCK SHT 1 OF 12						
REV 1	PCN XXXXXXXX	DESCRIPTION XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX	DATE DD/MM/YY	MADE STEF	DATE 10/10/96	GE Medical Systems BUC - FRANCE	DESIGN TITLE INTERFACE 2				
MADE BY	NAME	DATE DD/MM/YY	APPR BY NAME	CHECK CLOATRE	DATE 10/10/96		FIRST MADE FOR VAMP				
REV 2	PCN XXXXXXXX	DESCRIPTION XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX	DATE DD/MM/YY	APPR ENCELLAZ	DATE 10/10/96		DRAWING NO.		REV A	SIZE 3	SHT/CONT. ON 1 / 12
MADE BY	NAME	DATE DD/MM/YY	APPR BY NAME	OTHER	DATE		2138744-3SCH				

Blank page


BLOCK PATHNAME										\$PROJECT VAMP/int2_v8_ASM/int2_v8										BLOCK SHT										2 OF 12																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																
REV 1	PCN	XXXXXXX	DESCRIPTION	XXXXXXXXXXXXXXXXXXXXXXXXXXXX	MADE	STEF	DATE	10/10/96	GE Medical Systems BUC - FRANCE										DESIGN TITLE										INTERFACE 2																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																	
MADE BY	NAME		DATE	DD/MM/YY	APPR BY	NAME	DATE	DD/MM/YY											CHECK	CLOATRE	DATE	10/10/96	FIRST MADE FOR										VAMP																																																																																																																																																																																																																																																																																																																																																																																																																																																																																													
REV 2	PCN	XXXXXXX	DESCRIPTION	XXXXXXXXXXXXXXXXXXXXXXXXXXXX	APPR	ENCCELLAZ	DATE	10/10/96											DRAWING NO.																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																											
MADE BY	NAME		DATE	DD/MM/YY	APPR BY	NAME	DATE	DD/MM/YY											OTHER		DATE		2138744-3SCH										REV A										SIZE 3										SHT/CONT. ON 2 / 12																																																																																																																																																																																																																																																																																																																																																																																																																																																																									
1																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																														

Blank page


BLOCK PATHNAME \$PROJECT VAMP/int2_v8 ASM/int2_v8										BLOCK SHT 3 OF 12											
REV 1	PCN	xxxxxxx	DESCRIPTION xxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxx			MADE	STEF	DATE	10/10/96	GE Medical Systems		DESIGN TITLE INTERFACE 2									
MADE BY	NAME		DATE	DD/MM/YY	APPR BY	NAME	DATE	DD/MM/YY	DATE			10/10/96	FIRST MADE FOR VAMP								
REV 2	PCN	xxxxxxx	DESCRIPTION xxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxx			APPR	ENCCELLAZ	DATE	10/10/96	BUC - FRANCE		DRAWING NO.		REV	SIZE	SHT/CONT. ON					
MADE BY	NAME		DATE	DD/MM/YY	APPR BY	NAME	DATE	DD/MM/YY	DATE			10/10/96	2138744-3SCH		A	3	3 / 12				
		2		3		4		5		6		7		8		9		10		11	

Blank page


Blank page


BLOCK PATHNAME				BLOCK SHT				5 OF 12			
REV 1	PCN	XXXXXXXX	DESCRIPTION	XXXXXXXXXXXXXXXXXXXXXXXXXXXX	MADE	STEF	DATE	10/10/96	GE Medical Systems BUC - FRANCE		
MADE BY	NAME		DATE	DD/MM/YY	CHECK	CLOATRE	DATE	10/10/96			
REV 2	PCN	XXXXXXXX	DESCRIPTION	XXXXXXXXXXXXXXXXXXXXXXXXXXXX	APPR	ENCILLAZ	DATE	10/10/96	DESIGN TITLE INTERFACE 2 FIRST MADE FOR VAMP		
MADE BY	NAME		DATE	DD/MM/YY	OTHER		DATE				
									DRAWING NO.	REV	SIZE
									2138744-3SCH	A	3
									SHT/CONT. ON		
									5 / 12		

Blank page


BLOCK PATHNAME										\$PROJECT VAMP/int2_v8 ASM/int2_v8										BLOCK SHT										6		OF		12											
REV 1		PCN		XXXXXXXX		DESCRIPTION		XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX										MADE		STEF		DATE		10/10/96		GE Medical Systems BUC - FRANCE										DESIGN TITLE				INTERFACE 2					
MADE BY		NAME				DATE		DD/MM/YY		APPR BY		NAME				CHECK		CLOATRE		DATE		10/10/96		FIRST MADE FOR												VAMP									
REV 2		PCN		XXXXXXXX		DESCRIPTION		XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX										APPR		ENCCELLAZ		DATE		10/10/96												DRAWING NO.				REV		SIZE		SHT/CONT. ON	
MADE BY		NAME				DATE		DD/MM/YY		APPR BY		NAME				OTHER				DATE				2138744-3SCH												A		3		6 / 12					
1						2				3				4				5				6				7				8				9				10				11			

Blank page


REV 1 PCN XXXXXXXX DESCRIPTION XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX				MADE STEF		DATE 10/10/96		GE Medical Systems BUC - FRANCE		DESIGN TITLE INTERFACE 2											
MADE BY NAME DATE DD/MM/YY APPR BY NAME DATE DD/MM/YY				CHECK CLOATRE		DATE 10/10/96				FIRST MADE FOR VAMP											
REV 2 PCN XXXXXXXX DESCRIPTION XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX				APPR ENCELLAZ		DATE 10/10/96		DRAWING NO. 2138744-3SCH		REV A	SIZE 3	SHT/CONT. ON 7 / 12									
MADE BY NAME DATE DD/MM/YY APPR BY NAME DATE DD/MM/YY				OTHER		DATE															
1		2		3		4		5		6		7		8		9		10		11	

Blank page


Blank page


BLOCK PATHNAME \$PROJECT VAMP/int2_v8 ASM/int2_v8											BLOCK SHT 9 OF 12																					
REV 1	PCN	xxxxxxx	DESCRIPTION xxxxxxxxxxxxxxxxxxxxxxxxxxxxxxx			MADE	STEF	DATE	10/10/96	GE Medical Systems			DESIGN TITLE			INTERFACE 2																
MADE BY	NAME		DATE	DD/MM/YY	APPR BY	CHECK	CLOATRE	DATE	10/10/96				FIRST MADE FOR			VAMP																
REV 2	PCN	xxxxxxx	DESCRIPTION xxxxxxxxxxxxxxxxxxxxxxxxxxxxxxx			APPR	ENCCELLAZ	DATE	10/10/96	BUC - FRANCE			DRAWING NO.			REV	SIZE	SHT/CONT. ON														
MADE BY	NAME		DATE	DD/MM/YY	APPR BY	OTHER		DATE					2138744-3SCH			A	3	9 / 12														
		1			2			3			4			5			6			7			8			9			10			11

Blank page

BLOCK PATHNAME										\$PROJECT VAMP/int2_v8_ASM/int2_v8										BLOCK SHT										10 OF 12																
REV 1	PCN	XXXXXXXX	DESCRIPTION										XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX										MADE	STEF	DATE		10/10/96		GE Medical Systems BUC - FRANCE										DESIGN TITLE				INTERFACE 2			
MADE BY	NAME		DATE		DD/MM/YY		APPR BY		NAME		DATE		DD/MM/YY		CHECK	CLOATRE	DATE		10/10/96		FIRST MADE FOR				VAMP																					
REV 2	PCN	XXXXXXXX	DESCRIPTION										XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX										APPR	ENCILLAZ	DATE		10/10/96		DRAWING NO. 2138744-3SCH										REV		SIZE		SHT/CONT. ON			
MADE BY	NAME		DATE		DD/MM/YY		APPR BY		NAME		DATE		DD/MM/YY		OTHER		DATE				A		3		10 / 12																					

Blank page



BLOCK PATHNAME				\$PROJECT VAMP/int2_v8 ASM/int2_v8				BLOCK SHT 11 OF 12			
REV 1	PCN	XXXXXXXX	DESCRIPTION	XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX				MADE	STEF	DATE	10/10/96
MADE BY	NAME		DATE	DD/MM/YY	APPR BY	NAME		CHECK	CLOATRE	DATE	10/10/96
REV 2	PCN	XXXXXXXX	DESCRIPTION	XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX				APPR	ENCELLAZ	DATE	10/10/96
MADE BY	NAME		DATE	DD/MM/YY	APPR BY	NAME		OTHER		DATE	
								GE Medical Systems		DESIGN TITLE	INTERFACE 2
								BUC - FRANCE		FIRST MADE FOR	VAMP
										DRAWING NO.	2138744-3SCH
										REV	A
										SIZE	3
										SHT/CONT. ON	11 / 12

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





Blank page


Blank page

MNEMONICS		19-0
MP PHASE 41 A5 A10 – HEATER SUPPLY BOARD		
XS	Small Focus	
XL	Large Focus	
XS/XL DEMAND	RMS Filament Current Requested	


Blank page


Blank page


Blank page


Blank page

Blank page

Blank page


Blank page

MNEMONICS		21-0
MP PHASE 41 A5 A13 – GRID/BIAS INTERFACE BOARD		
GRID INV FDBK	Micro grid inverter Feeback (grid operation)	
BIAS DMD A	Bias voltage value	
BIAS DMD B	Bias voltage value	
BIAS ENABLE	Bias enable	
BIAS INV FDBK	Micro grid inverter Feedback (bias operation)	
GRID FUNCT OK	36 V power supply OK	
GRID SWITCH OK	36 V grid switch OK	
GRID INV SUPPLY	Gird power supply	
GRID 0VLD	Grid current overload	
GRID PULSE OK	Feedback to Int 2 Board	
BIAS FUNCT OK	Bias voltage OK	
BIAS INV SUPPLY	Bias power voltage	
BIAS FDBK A	Bias voltage Feedback	
BIAS FDBK B	Bias voltage Feedback	
BIAS ON	Bias voltage Feedback to Int 2 Board	


Blank page

MP PHASE 41 A5 A13 – GRID/BIAS INTERFACE BOARD BLOCK DIAGRAM

21-1


Blank page

MP Phase 41 A5 A13


BLOCK PATHNAME				BLOCK SHT				1 OF 5			
REV 1	PCN	xxxxxxx	DESCRIPTION	xxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxx	MADE	NGO	DATE	07-06-95	GE Medical Systems		
MADE BY	NAME		DATE	DD/MM/YY	CHECK	ENCELLAZ	DATE	07-06-95			
REV 2	PCN	xxxxxxx	DESCRIPTION	xxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxx	APPR	NAME	DATE	DD/MM/YY	BUC - FRANCE		
MADE BY	NAME		DATE	DD/MM/YY	OTHER		DATE				
									DESIGN TITLE		
									G/B INTERFACE		
									FIRST MADE FOR		
									VAMP		
									DRAWING NO.		
									2137893SCH		
									REV	SIZE	SHT/CONT ON
									A	3	1 / 5

MP Phase 41 A5 A13


BLOCK PATHNAME				BLOCK SHT				2 OF 5			
REV 1	PCN	xxxxxxx	DESCRIPTION	xxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxx	MADE	NGO	DATE	07-06-95	GE Medical Systems		
MADE BY	NAME		DATE	DD/MM/YY	CHECK	ENCELLAZ	DATE	07-06-95			
REV 2	PCN	xxxxxxx	DESCRIPTION	xxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxx	APPR	NAME	DATE	DD/MM/YY	BUC - FRANCE		
MADE BY	NAME		DATE	DD/MM/YY	OTHER		DATE		DESIGN TITLE		
									G/B INTERFACE		
									FIRST MADE FOR		
									VAMP		
									DRAWING NO.		
									2137893SCH		
									REV	SIZE	SHT/CONT
									A	3	2 / 5

Blank page

MP Phase 41 A5 A13


BLOCK PATHNAME				BLOCK SHT				3 OF 5			
REV 1	PCN	xxxxxxx	DESCRIPTION	xxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxx	MADE	NGO	DATE	07-06-95	GE Medical Systems		
MADE BY	NAME		DATE	DD/MM/YY	CHECK	ENCELLAZ	DATE	07-06-95			
REV 2	PCN	xxxxxxx	DESCRIPTION	xxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxx	APPR	NAME	DATE	DD/MM/YY	BUC - FRANCE		
MADE BY	NAME		DATE	DD/MM/YY	OTHER		DATE				
									DESIGN TITLE		
									G/B INTERFACE		
									FIRST MADE FOR		
									VAMP		
									DRAWING NO.		
									2137893SCH		
									REV	SIZE	SHT/CONT ON
									A	3	3 / 5

MP Phase 41 A5 A13

BLOCK PATHNAME				BLOCK SHT				4 OF 5			
REV 1	PCN	xxxxxxx	DESCRIPTION	xxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxx	MADE	NGO	DATE	07-06-95	GE Medical Systems		
MADE BY	NAME		DATE	DD/MM/YY	CHECK	ENCELLAZ	DATE	07-06-95			
REV 2	PCN	xxxxxxx	DESCRIPTION	xxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxx	APPR	NAME	DATE	DD/MM/YY	BUC - FRANCE		
MADE BY	NAME		DATE	DD/MM/YY	OTHER		DATE				
									DESIGN TITLE		
									G/B INTERFACE		
									FIRST MADE FOR		
									VAMP		
									DRAWING NO.		
									2137893SCH		
									REV	SIZE	SHT/CONT ON
									A	3	4 / 5

Blank page


MP Phase 41 A5 A13


BLOCK PATHNAME										sPROJECT VAMP/gb_int_v4_ASM/gb_int_v4										BLOCK SHT										5	OF	5	H																																																																											
REV	1	PCN	xxxxxxx			DESCRIPTION			xxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxx			MADE	NGO	DATE	07-06-95			GE Medical Systems										DESIGN TITLE										G/B INTERFACE																																																																						
MADE BY		NAME				DATE			DD/MM/YY			APPR	BY	NAME	DATE													07-06-95			FIRST MADE FOR										VAMP																																																																			
REV	2	PCN	xxxxxxx			DESCRIPTION			xxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxx			MADE	NGO	DATE	DD/MM/YY			BUC - FRANCE										DRAWING NO.										REV										SIZE										SHT/CONT										ON																																								
MADE BY		NAME				DATE			DD/MM/YY			APPR	BY	NAME	DATE													DD/MM/YY			2137893SCH										A										3										5 /										5																																					
1						2						3			4						5													6													7													8													9													10													11									

Blank page


Blank page


C2101741.PL1
NON PCB


Blank page


Blank page


Blank page

Blank page

Blank page