# Optima EMS ECG Management System XML Import Guideline for Third-Party Device

Software Version 1.0  
2069773-003 Revision C

Optiam EMS  
English  
© 2012, 2013 General Electric Company.  
All Rights Reserved.

#### **Publication Information**

The information in this manual applies only to Optima EMS ECG Management System Version 1.0. It does not apply to earlier product versions. Due to continuing product innovation, specifications in this manual are subject to change without notice.

MUSE and MAC are trademarks owned by GE Medical Systems *Information Technologies*, Inc., a General Electric Company going to market as GE Healthcare. All other trademarks contained herein are the property of their respective owners.

The document part number and revision appear at the bottom of each page. The revision identifies the document's update level. The revision history of this document is summarized in the following table.

Revision	Date	Comments
A	1 November 2012	Internal Release
B	10 December 2012	Initial release for customer
C	17 July 2013	Per SPR HCSDM00221960 a new section titled: <i>Setting Up the Connection with NK1550P</i> was added.

To access other GE Healthcare manuals, go to the Common Documentation Library (CDL), located at [http://www.gehealthcare.com/usen/service/biomed\\_tech\\_selfservice/services\\_user\\_doc/products/support.html](http://www.gehealthcare.com/usen/service/biomed_tech_selfservice/services_user_doc/products/support.html) and click **Cardiology**.

To access Original Equipment Manufacturer (OEM) manuals, go to the device manufacturer's Web site.

#### **Service Manual Language Information**

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<b>Responsibility of the Purchaser/Customer .....</b>	<b>12</b>
---	-----------

<b>Additional Assistance.....</b>	<b>12</b>
-----------------------------------	-----------

## 2 Third-Party ECG Data Solution Overview

## 3 Third-Party ECG Data Solution User Guidelines

<b>Setting Up the Connection with NK1550P .....</b>	<b>16</b>
Before Installation.....	16
NK1550P Network Setup.....	17
Communication Setup.....	21
Installing and Setting Up the ECG Viewer II.....	24

<b>Setting Up the Transmission Layer .....</b>	<b>26</b>
FTP Transmission Layer Example.....	26

<b>The Parser Layer .....</b>	<b>39</b>
-------------------------------	-----------

<b>Optima EMS Folder for Each Device .....</b>	<b>45</b>
--	-----------

## 4 Third-Party ECG Data Solution Detailed Design

<b>The System Compatible XML Format .....</b>	<b>47</b>
---	-----------

<b>Third-Party Solution Architecture .....</b>	<b>53</b>
The Top Layer .....	53
The Core Layer .....	53
The Adapters Layer .....	56

## 5 Third-Party ECG Data Extensibility

# 1

## Introduction

This chapter provides general information required for the proper use of the system and this manual. Familiarize yourself with this information before using the system.

## Manual Information

This section provides information for the correct use of this manual.

Keep this manual with the equipment at all times and periodically review it. You should request training assistance from GE Healthcare, if needed.

## Intended Audience

This manual is intended for use by anyone responsible for importing XML files from third-party device.

## Manual Purpose

This manual provides a detailed description of the third-party ECG data solution of the Optima EMS system (also referred to as "the system" throughout this manual).

## Document Conventions

This manual uses the following conventions.

### Typographical Conventions

Convention	Description
<b>Bold</b> Text	Indicates keys on the keyboard, text to enter, or hardware items such as buttons or switches on the equipment.
<b><i>Italicized-Bold</i></b> Text	Indicates software terms that identify menu items, buttons or options in various windows.
<b>CTRL+ESC</b>	Indicates a keyboard operation. A plus (+) sign between the names of two keys indicates that while holding the first key, you should press and release the second key. For example, Press <b>CTRL+ESC</b> means to press and hold the <b>CTRL</b> key and then press and release the <b>ESC</b> key.

Convention	Description
<space>	Indicates that you must press the spacebar. When instructions are given for typing a precise text string with one or more spaces, the point where you must press the spacebar is indicated as <space>. This ensures that the correct number of spaces is inserted in the correct positions within the literal text string. The purpose of the < > brackets is to distinguish the command from the literal text within the string.
Enter	Indicates that you must press the <b>Enter</b> or <b>Return</b> key on the keyboard. Do not type <b>Enter</b> .
>	The greater than symbol, or right angle bracket, is a concise method to indicate a sequence of menu selections.  For example, the statement "From the main menu, select <b>System</b> > <b>Setup</b> > <b>Options</b> to open the <b>Option Activation</b> window" replaces the following: <ol style="list-style-type: none"> <li>1. From the main menu, select <b>System</b> to open the <b>System</b> menu.</li> <li>2. From the <b>System</b> menu, select <b>Setup</b> to open the <b>Setup</b> menu.</li> <li>3. From the <b>Setup</b> menu, select <b>Options</b> to open the <b>Option Activation</b> window.</li> </ol>

## Illustrations

All illustrations in the manual are provided as examples only. Depending on system configuration, screens in the manual may differ from the screens on your system.

All patient names and data are fictitious. Any similarity to actual persons is coincidental.

## Notes

Notes provide application tips or additional information that, while useful, are not essential to the correct operation of the system. They are called out from the body text through a flag word and indentation, as follows:

### NOTE:

The tip or additional information is indented below the **NOTE** flag word.

## Responsibility of the Purchaser/Customer

The customer is responsible for providing appropriate desks, chairs, electrical wall outlets, network connections, analog phone lines, and for locating any of the Optima EMS system components described in this manual in compliance with all local, state or national codes.

## Additional Assistance


GE Healthcare maintains a trained staff of application and technical experts to answer questions and respond to issues and problems that may arise during the installation, maintenance, and use of this system.

Contact your local GE Healthcare representative to request additional assistance.

# 2

## Third-Party ECG Data Solution Overview

Third-party ECG data solution provides the capability to receive data from third-party devices (such as NK devices, Fukuda devices, and so on).

### Solution Overview

As illustrated in the figure Solution Overview, there are two steps when sending third-party data to the system.

1. Receive data from third-party devices.  
The most preferred transmission method is FTP and shared folder. These two methods need no extra proprietary software. Free software could handle the transmission equally well. Serial transmission, in most cases, needs software from device vendors.
2. Convert the received data to the system compatible format and send it to the system.  
This step converts data to the system compatible format. The solution includes three parsers, one is dedicate for NK devices, one is dedicate for Fukuda devices, and the last is dedicated to the SCP open format. As will be illustrated in the document, other newly developed parsers could also easily be added to this solution.

# 3

## Third-Party ECG Data Solution User Guidelines

Third-party ECG data solution provides the capability to receive data from third-party devices (such as NK devices, Fukuda devices, and so on).

# Setting Up the Connection with NK1550P

## Before Installation

ECG Viewer II is required when connecting to NK1550P.

**NOTE:**

ECG Viewer II is compatible with Windows 2000 and Windows XP. If you are using Windows 7, compatible mode is needed during installing.


## NK1550P Network Setup


- From the main window, click **Menu**.

- Click **System Testing**.


3. Click **Network Setup**.


4. Double-click the network adapter.


5. Setup the **IP Address**, **Subnet Mask**, and **Default Gateway** according to the network environment.

6. Exit the *Network Setup* by selecting *File > Close*.

7. Close the *System Testing* by clicking *End* in the *System Testing* window.


## Communication Setup


- From the main window, click **Menu**.


- Click **System Setup**.

3. Click **Communication Setup**.

4. Click **Ethernet Setup**.


5. Setup the Ethernet according to the following illustrations.


**NOTE:**

Enter the **Directive Server Address** and **File Server Address** with the IP address of the computer where **ECG Viewer II** is installed.

- a. Enter the highlighted information.

- b. Enter the highlighted information.

6. Click **Check in** to save and exit.


## Installing and Setting Up the ECG Viewer II

Install the **ECG Viewer II** with the installation CD.


**NOTE:**

Some programs of the **ECG Viewer** can be run only when the installation CD remains in the drive.


1. Select the installation directory and run **QB903D**.


2. Create the Database accordingly.

3. Access the **FTP Setup** by selecting **Setup > FTP Setup** in the **ECG File Select** window.

4. Make the following settings as shown in the following illustration.

5. The **FTP** window will open, click **OK** to restart **QB903D**.

## Setting Up the Transmission Layer

This section demonstrates how to setup the transmission layer used in the third-party ECG Data solution.

### FTP Transmission Layer Example

#### Demonstration Environment

FTP Server	Quick Easy FTP Server
Third-Party Device	Fukuda CardiMax FX-7202

## Installing the FTP Server

1. Double-click the FTP server installation.

The **Welcome** window opens.

2. Click **Next**.


The **License Agreement** window opens.

3. Select **I accept the agreement**, and click **Next**.

The **Select Destination Location** window opens.

4. Click **Browse** to select the location where you want to install the FTP server.

5. Click **Next**.

The **Select Start Menu Folder** window opens.

6. Click **Browse** to select a different folder.

7. Click **Next**.

The **Select Additional Tasks** window opens.

If you want to create a desktop icon for this FTP server, select the **Create a desktop icon** check box.

8. Click **Next**.


The **Ready to Install** window opens.

Review the information to make sure it is correct. If you want to make a change, click **Back** and you will go back to the previous windows.

9. Click **Install**.

The installation begins.

10. After the installation is finished, click **OK** to finish and exit the installation.


#### Creating a New FTP Server Account

The first time you run the FTP server, the configuration guide will display.


1. Run the FTP server by double-clicking the FTP server shortcut icon on the desktop.

The **Create a new account** window opens.

2. Type the **User Name**.

3. Click **Next**.  
The **Enter the password** window opens.
4. Type the **Password** in the **Password** filed.
5. Type the password again to confirm it in the **Confirm Password** field.


6. Click **Next**.  
The **Root Directory for the account** window opens.
7. Click **Browse** to choose the location that you want to make as the Root Directory.

8. Click **Next**.  
The **Permission settings** window opens.

9. Enable all the functions listed in the window by selecting the check boxes next to each function..

10. Click **Next**.

The **Your account has been created successfully** window opens.

11. Click **Finish** to finish the configuration.

## Starting the Server

After creating a new FTP server account, the FTP server configuration utility opens automatically.

1. Click **Server Configuration** to open the **Server Configuration** window.

2. Enter the following information:

- **Server IP**
- **Server Port**


3. Start the server by clicking the **Start** button.

The server is now started.

## Configuring the CardiMax FX-7202

1. Power on the CardiMax FX-7202.
2. In the Main window, click **Menu**.

The **Mode Select** window opens.

3. Select **Setup**.

The **Setup (Menu)** window opens

4. Select **File**.


The **Setup (File)** window opens.

5. Enter **ECG Device IP Address** and **HOST IP Address** by clicking the buttons next to the field names.

**NOTE:**


The CardiMax FX-7202 does not support DHCP.

The HOST IP address is the IP address of the FTP server.

## Sending Data to the FTP Server

1. Acquire ECG data on the ECG device.
2. In the Main window, click **Menu**.


The **Mode Select** window opens.

3. Select **File**.

The **File (Menu)** window opens.

4. Click **Send**.

The **Sending ECG Data** window opens.

5. Click **Media**.

The ECG records on the media are listed in the window.

6. Select the records that you want to send, and click **Send**.

NO.	ID NO.	DATE/TIME	SEX AGE
1	123456	2012/ 6/14 14:34:03	
2	12345678	2012/ 6/13 15:38:22	
3	88888888	2012/ 6/13 15:28:30	M 65
4	12345678	2012/ 6/13 15:17:36	
5	12345678	2012/ 6/13 15:12:11	
6	369	2012/ 6/ 8 13:34:12	F 56
7		2012/ 6/ 8 13:33:00	
8		2012/ 6/ 5 10:45:46	
9	1	2012/ 6/ 5 10:20:57	F 21
10	1234567890123456	2012/ 6/ 5 10:13:39	M 26


7. A message displays indicating the progress.

Wait until all the records have been sent.

NO.	ID NO.	DATE/TIME	SEX AGE
1	123456	2012/ 6/14 14:34:03	
2	12345678	2012/ 6/13 15:38:22	
3	88888888	2012/ 6/13 15:28:30	M 65
4	12345678	2012/ 6/13 15:17:36	
5	12345678	2012/ 6/13 15:12:11	
6	369	2012/ 6/ 8 13:34:12	F 56
7		2012/ 6/ 8 13:33:00	
8		2012/ 6/ 5 10:45:46	
9	1	2012/ 6/ 5 10:20:57	F 21
10	1234567890123456	2012/ 6/ 5 10:13:39	M 26

### Checking Transmitted Data

Check the FTP root folder to make sure an ECG file exists.

# The Parser Layer

The parser layer parses the data from a third-party device, converts it to the system compatible format, and send it to the system server. The parser layer is the core of this solution. It consists of a framework, a configurator, and several parsers.

The configurator modifies the configuration file, which defines how the framework works.


The parsers are each dedicated to a specific device, parsing data from the device.

The framework reads the configuration file, which must be named **ConverterConfig.xml**, loads the parser, and allows everything to work together.

## Configuring the Parser

1. Run the configurator.


The **3rd Party Convert Configurator** window opens.


The configurator has two major functions.

- Defining parsers, which is done in the **Adapters** tab.
- Defining monitors, which is down in the **Monitors** tab.


2. At the **Adapters** tab, click **Add**.

3. Type a new parser name in the **Name** field and click **Load**.


4. Select an assembly where there is a parser, such as **FUKUDA2XMLAdapter.dll**.


The parser/adapter in the assembly will be listed.


For example, as indicated in the following figure the **FUKUDA2XMLAdapter.SCAdapter** and **FUKUDA2XMLAdapter.FDAdapter** are listed.


5. Select the appropriate adapter and click **Add**.

6. Click the **Monitors** tab, and click **Add**.

7. Select an Adapter/Parser.


Add the information to the other fields as follows:

Field	Description	Required
Source	This is the path where the third-party data is located.	Yes
Dest	This is the path where the converted data will be located.	Yes
Vender	Vender of the device.	Yes
Module	Device name.	Yes
Serial	Device serial number.	Yes

Each added item will start a thread for converting data.

8. Save the configuration by selecting **File > Save**. The saved file name must be **ConverterConfig.xml**.

## Running the Parser

To run the parser, click the **ThirdPartyDataConverter.exe**.

The following window opens.

## Transmitting Files

1. Transmit a file from the FUKUDA device.
  2. Open the **D:\Converted** folder.

If the file transmitted successfully, you will find an XML file similar to the following, in the **D:\Converted** folder:

# Optima EMS Folder for Each Device

A device needs to be configured before any acquisition work is done.

The system server will create a folder for each device under the current site's folder.

For example, a folder on the server named **D:\somefolder** with the shared name of **MACIT ECG Data** on the **Minhang** site, should have a folder named **D:\somefolder\Minhang** on the **Minhang** site.

If the device **p1025** is added to the **Minhang** site, there will be a folder named **D:\somefolder\Minhang\p1025**. This folder is where the data from **p1025** will be stored and all data going into that folder will be coming from the device **p1025**.

It is the configurator's responsibility to route the converted xml file to the correct folder.

# 4

## Third-Party ECG Data Solution Detailed Design

### The System Compatible XML Format

This XML format is the data convention between the system and any third-party system. The system can process data from any system as long as it supports the Compatible XML Format.

Following is the XML format schema.

```
<?xml version="1.0" encoding="utf-8"?>
<xsschema attributeFormDefault="unqualified" elementFormDefault="qualified"
xmlns:xs="http://www.w3.org/2001/XMLSchema">
  <xs:element name="RestingECG">
    <xs:complexType>
      <xs:sequence>
        <xs:element name="Patient" minOccurs="1" maxOccurs="1">
          <xs:complexType>
            <xs:attribute name="id" type="xs:string" use="optional" />
            <xs:attribute name="name" type="xs:string" use="optional" />
            <xs:attribute name="gender" type="xs:string" use="optional" />
            <xs:attribute name="ageUnit" type="xs:string" use="optional" />
            <xs:attribute name="age" type="xs:unsignedInt" use="optional" />
            <xs:attribute name="height" type="xs:unsignedInt" use="optional" />
            <xs:attribute name="weight" type="xs:unsignedInt" use="optional" />
          </xs:complexType>
        </xs:element>
        <xs:element name="Device" minOccurs="1" maxOccurs="1">
          <xs:complexType>
```

```
<xs:attribute name="vender" type="xs:string" use="required" />
<xs:attribute name="module" type="xs:string" use="required" />
<xs:attribute name="serial" type="xs:string" use="required" />
</xs:complexType>
</xs:element>
<xs:element name="DiagnosisResult" minOccurs="0" maxOccurs="1">
<xs:complexType>
<xs:attribute name="text" type="xs:string" use="required" />
</xs:complexType>
</xs:element>
<xs:element name="Measurements" minOccurs="0" maxOccurs="1">
<xs:complexType>
<xs:attribute name="atrialRate" type="xs:unsignedShort" use="optional" />
<xs:attribute name="ventRate" type="xs:unsignedShort" use="optional" />
<xs:attribute name="pOn" type="xs:unsignedShort" use="optional" />
<xs:attribute name="pOff" type="xs:unsignedShort" use="optional" />
<xs:attribute name="qOn" type="xs:unsignedShort" use="optional" />
<xs:attribute name="qOff" type="xs:unsignedShort" use="optional" />
<xs:attribute name="tOff" type="xs:unsignedShort" use="optional" />
<xs:attribute name="sysBp" type="xs:unsignedShort" use="optional" />
<xs:attribute name="diaBp" type="xs:unsignedShort" use="optional" />
<xs:attribute name="prInterval" type="xs:unsignedShort" use="optional" />
<xs:attribute name="qrsDuration" type="xs:unsignedShort" use="optional" />
<xs:attribute name="qtInterval" type="xs:unsignedShort" use="optional" />
<xs:attribute name="qtc" type="xs:unsignedShort" use="optional" />
<xs:attribute name="qtcFred" type="xs:unsignedShort" use="optional" />
<xs:attribute name="qtcFram" type="xs:unsignedShort" use="optional" />
<xs:attribute name="pAxis" type="xs:short" use="optional" />
<xs:attribute name="rAxis" type="xs:short" use="optional" />
<xs:attribute name="tAxis" type="xs:short" use="optional" />
<xs:attribute name="qrsCount" type="xs:unsignedShort" use="optional" />
</xs:complexType>
</xs:element>
<xs:element name="WaveData" minOccurs="1" maxOccurs="1">
<xs:complexType>
<xs:sequence>
<xs:element name="Summary" minOccurs="1" maxOccurs="1">
<xs:complexType>
```

```
<xs:attribute name="precision" type="xs:decimal" use="required" />
<xs:attribute name="sampleRate" type="xs:unsignedShort" use="required" />
<xs:attribute name="rhythmsDuration" type="xs:unsignedShort" use="required" />
<xs:attribute name="lowPassFilter" type="xs:double" use="optional" />
<xs:attribute name="highPassFilter" type="xs:double" use="optional" />
</xs:complexType>
</xs:element>
<xs:element name="Medians" minOccurs="0" maxOccurs="1">
<xs:complexType>
<xs:sequence>
<xs:element minOccurs="8" maxOccurs="12" name="Lead">
<xs:complexType>
<xs:attribute name="id" type="xs:string" use="required" />
<xs:attribute name="samples" type="xs:unsignedShort" use="required" />
<xs:attribute name="value" type="xs:string" use="required" />
</xs:complexType>
</xs:element>
</xs:sequence>
</xs:complexType>
</xs:element>
<xs:element name="Rhythms" minOccurs="1" maxOccurs="1">
<xs:complexType>
<xs:sequence>
<xs:element minOccurs="8" maxOccurs="12" name="Lead">
<xs:complexType>
<xs:attribute name="id" type="xs:string" use="required" />
<xs:attribute name="samples" type="xs:unsignedShort" use="required" />
<xs:attribute name="value" type="xs:string" use="required" />
</xs:complexType>
</xs:element>
</xs:sequence>
</xs:complexType>
</xs:element>
</xs:sequence>
</xs:complexType>
```

```

    </xs:element>
  </xs:sequence>
  <xs:attribute name="acquiredDateTime" type="xs:dateTime" use="optional"/>
</xs:complexType>
</xs:element>
</xs:schema>

```

## Schema Definitions

Element	Attribute	Definition
<i>Patient</i>	<i>ID</i>	Patient ID, which should be unique in site scope.
	<i>Name</i>	Patient name/family name/given name. The maximum length is 20 bytes.
	<i>Gender</i>	Enumeration values: Male, Female.
	<i>AgeUnit</i>	Enumeration values: Year, Month, Week Day.
	<i>Age</i>	Patient age. The attribute works with the <i>AgeUnit</i> attribute to describe patient age.
	<i>Height</i>	Patient height. Unit is centimeter.
<i>Device</i>	<i>Vender</i>	Device Vender. Maximum length is 50 bytes.
	<i>Module</i>	Device Type. Maximum length is 50 bytes.
	<i>Serial</i>	Device serial number. Maximum length is 50 bytes.
<i>DiagnosisResult</i>	<i>Text</i>	Diagnosis text. Maximum length is 1024 bytes.
<i>Measurements</i>	<i>atrialRate</i>	Atrial rate. Unit is beats/minute.
	<i>ventRate</i>	Ventricular rate. Unit is beats/minute.
	<i>pOn</i>	P onset time. Unit is ms.
	<i>pOff</i>	P offset time. Unit is ms.
	<i>qOn</i>	Q onset time. Unit is ms.
	<i>qOff</i>	Q offset time. Unit is ms.
	<i>tOff</i>	T offset time. Unit is ms.

## Schema Definitions (cont'd.)

		<b>sysBp</b>	Systolic blood pressure. Unit is mmHg.
		<b>diaBp</b>	Diastolic blood pressure. Unit is mmHg.
		<b>prInterval</b>	PR interval. Unit is ms.
		<b>qrsDuration</b>	QRS duration. Unit is ms.
		<b>qtInterval</b>	QT interval. Unit is ms.
		<b>qtc</b>	QTC. Unit is ms.
		<b>qtcFred</b>	QTC Fridericia. Unit is ms.
		<b>qtcFram</b>	QTC Fram. Unit is ms.
		<b>pAxis</b>	P Axis. Unit is degree.
		<b>rAxis</b>	R Axis. Unit is degree.
		<b>tAxis</b>	T Axis. Unit is degree.
		<b>qrsCount</b>	QRS count. Unit is heartbeat.
<b>WaveData</b>	<b>Summary</b>	<b>Precession</b>	Must be 4.88uV.
		<b>SampleRate</b>	Must be 500Hz, if not, resample original data.
		<b>RhythmsDuration</b>	Must be 10 second. The system only supports 10s ECG.
		<b>LowPassFilter</b>	Low pass filter when acquiring the wave data. Unit is Hz.
		<b>HightPassFilter</b>	High pass filter when acquiring the wave data. Unit is Hz.
	<b>Medians</b>	<b>ID</b>	ECG Lead ID definition, see the Lead ID Definition table.
		<b>Samples</b>	Must be 600 (1.2 s)
		<b>Value</b>	The lead sample values separated with white spaces.
	<b>Rhythms</b>	<b>ID</b>	ECG lead ID definition, see the Lead ID Definition table.
		<b>Samples</b>	Must be 5000 for 10 seconds ECG data.
		<b>Value</b>	The lead sample values separated with white spaces.
		<b>acquiredDateTime</b>	The date and time when this data is acquired.

Available lead IDs in the system are: (All lead IDs are case-sensitive).

### Lead ID Definition


Lead ID	Comments
I	I lead in waveform.
II	II lead in waveform.
III	III lead in waveform.
aVL	aVL lead in waveform.
aVR	aVR lead in waveform.
aVF	aVF lead in waveform.
V1	V1 lead in waveform.
V2	V2 lead in waveform.
V3	V3 lead in waveform.
V4	V4 lead in waveform.
V5	V5 lead in waveform.
V6	V6 lead in waveform.
V7	V7 lead in waveform.
V8	V8 lead in waveform.
V9	V9 lead in waveform.
V2R	V2R lead in waveform.
V3R	V3R lead in waveform.
V4R	V4R lead in waveform.
V5R	V5R lead in waveform.
V6R	V6R lead in waveform.
V7R	V7R lead in waveform.
V8R	V8R lead in waveform.
V9R	V9R lead in waveform.
X	X lead in waveform.
Y	Y lead in waveform.
Z	Z lead in waveform.
mX	mX lead in waveform.
mY	mY lead in waveform.
mZ	mZ lead in waveform.
Xd	Xd lead in waveform.
Yd	Yd lead in waveform.
Zd	Zd lead in waveform.
mXd	mXd lead in waveform.

### Lead ID Definition (cont'd.)

<i>mYd</i>	mYd lead in waveform.
<i>mZd</i>	mZd lead in waveform.
<i>NEHB_D</i>	NEHB_D lead in waveform.
<i>NEHB_A</i>	NEHB_A lead in waveform.
<i>NEHB_J</i>	NEHB_J lead in waveform.
<i>ML</i>	ML lead in waveform.
<i>CH</i>	CH lead in waveform.
<i>A1</i>	A1 lead in waveform.
<i>A2</i>	A2 lead in waveform.
<i>A3</i>	A3 lead in waveform.
<i>A4</i>	A4 lead in waveform.

## Third-Party Solution Architecture

The following figure is the architecture of the third-party ECG data solution.

### The Top Layer

The top layer is the transmission layer responsible for receiving data from third-party devices. The methods for receiving data in this solution are general FTP, shared folder, and external storage (SD Card, USB Disk, and so on). The serial line approach needs specialized software support, which is generally bundled with devices.

### The Core Layer

The middle layer is the core part of this solution, it contains three parts: the configuration, the configurator, and the converters.

The configuration reads and writes the configuration file. Configurator modifies the configuration files. Converters are organized by the rules in configuration files and converts the third-party data.

The configuration file schema is as follows:

```
<?xml version="1.0" encoding="utf-8"?>
<xss:schema xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
  xmlns:xsd="http://www.w3.org/2001/XMLSchema" xmlns:xs="http://www.w3.org/2001/
  XMLSchema" attributeFormDefault="unqualified" elementFormDefault="qualified"
  targetNamespace="http://parser.3rdParty.com">
  <xss:element name="Configuration">
    <xss:complexType>
      <xss:sequence>
        <xss:element name="Adapters">
          <xss:complexType>
            <xss:sequence>
              <xss:element name="Adapter">
                <xss:complexType>
                  <xss:attribute name="Name" type="xsd:string" use="required" />
                  <xss:attribute name="Assembly" type="xsd:string" use="required" />
                  <xss:attribute name="Type" type="xsd:string" use="required" />
                </xss:complexType>
              </xss:element>
            </xss:sequence>
          </xss:complexType>
        </xss:element>
      </xss:sequence>
    </xss:complexType>
  </xss:element>
  <xss:element name="Monitors">
    <xss:complexType>
      <xss:sequence>
        <xss:element name="Monitor">
          <xss:complexType>
            <xss:attribute name="Adapter" type="xsd:string" use="required" />
            <xss:attribute name="Source" type="xsd:string" use="required" />
            <xss:attribute name="Dest" type="xsd:string" use="required" />
            <xss:attribute name="Vender" type="xsd:string" use="required" />
            <xss:attribute name="Module" type="xsd:string" use="required" />
            <xss:attribute name="Serial" type="xsd:unsignedByte" use="required" />
          </xss:complexType>
        </xss:element>
      </xss:sequence>
    </xss:complexType>
  </xss:element>
</xss:schema>
```

```

    </xsd:sequence>
    </xsd:complexType>
    </xsd:element>
</xs:schema>
```

There are two elements in the configuration file: **Adapters** and **Monitors**.

**Adapters** describe the adapters which will be created. **Monitors** describe the monitoring threads that will be created.

**Adapters** and **Monitors** are related to the **Adapter.Name** and **Monitor.Adapter** attributes.

Following is a sample of a configuration file. This configuration indicates there is a monitoring thread, which makes use of the NKType2 adapter, and monitoring the **D:\Source** folder. If there are any files in the **D:\Source** folder, the monitor will try to convert it to the system compatible XML file and then save the XML file to **D:\dest** folder.

```

<?xml version="1.0" encoding="utf-8"?>
<Configuration xmlns="http://parser.3rdParty.com">

    <Adapters>
        <Adapter Name="NKType2" Assembly="NK2XMLAdapter.dll" Type=
            "NK2XMLAdapter.Adapter" />
    </Adapters>
    <Monitors>
        <Monitor Adapter="NKType2" Source="D:\Source" Dest="d:\Dest" Vender="NK"
            Module="P1025" Serial="12345" />
    </Monitors>
</Configuration>
```

### Configuration File Schema Descriptions

Element	Attribute	Description
<b>Adapter</b>	<b>Name</b>	Adapter Name
	<b>Assembly</b>	The assembly file, where this adapter is located.
	<b>Type</b>	The full name, including namespace, of the adapter.

### Configuration File Schema Descriptions (cont'd.)

<b>Monitor</b>	<b>Adapter</b>	Name of the adapter used by this monitor.
	<b>Source</b>	The folder being monitored.
	<b>Dest</b>	The destination folder of the converted system compatible XML file.
	<b>Vender</b>	Device vender. The value of this field will be copied to the converted XML file.
	<b>Module</b>	Device module. The value of this field will be copied to the converted XML file.
	<b>Serial</b>	Device serial. The value of this field will be copied to the converted XML file.

## The Adapters Layer

The Adapters layer provides flexibility that it could be extended by arbitrary adapters. To achieve this flexibility a configuration file, a configurator, and dependency injection technology are employed.

All adapters will need to implement the following interface.

```
public interface IAdapter
{
    string ToXML(string rawFileName);
}
```

This interface converts third-party data to the system compatible XML format.

# 5

## Third-Party ECG Data Extensibility

To add an additional adapter, perform the following procedure:

1. Develop an adapter which implements the IAdapter interface. (For developer only).
2. Copy the assembly containing the adapter to the install path.
3. Run the configurator, and open the existing configuration file.
4. Use the adapter editor to load the newly added assembly.
5. Add an adaptername, the adapter, and set the assembly and type.
6. Use the monitor editor to add a monitor. Set the adapter filed to the newly added adapter, and set the source and destination folder according to the requirement.
7. Save the configuration file.
8. Restart the *ThirdPartyDataConverter*.

GE Medical Systems  
*Information Technologies, Inc.*  
8200 West Tower Avenue  
Milwaukee, WI 53223 USA  
Tel: +1 414 355 5000  
+1 800 558 7044 (US Only)  
Fax: +1 414 355 3790

#### Asia Headquarters

GE Medical Systems  
*Information Technologies, Inc.*  
Asia; GE (China) Co., Ltd.  
1 Huatuo Road  
Zhangjiang Hi-tech Park Pudong  
Shanghai, People's Republic of China 201203  
Tel: +86 21 3877 7888  
Fax: +86 21 3877 7451

GE Medical Systems *Information Technologies, Inc.*, a General Electric Company, going to market as GE Healthcare.

[www.gehealthcare.com](http://www.gehealthcare.com)