


INSTRUCTIONS FOR USE AND OPERATION OF THE

intelect™ model 700

39 Lbs
Ship weight



CHATTANOOGA CORPORATION

foreword

This manual has been prepared for the owners and operators of Intellect Model 700™. It contains general instructions on operation, safety practices, maintenance and parts information. In order to obtain maximum life and efficiency from your Model 700™ and to aid in the safe operation of the unit, read and understand this manual thoroughly and become totally familiar with the controls on the panel and the various electrodes that come with the unit before operating it. The specifications put forth in this manual were in effect at the time of publication. However, owing to Chattanooga Corporation's policy of continuous improvement, changes to these specifications may be made at any time without obligation on the part of Chattanooga Corporation.

safety instruction

1. Read, understand and practice the safety and operating instructions. Know the limitations and hazards associated with Ultrasound. Observe the safety and operational decals placed on the unit.
2. Grounding — Make certain that the unit is electrically grounded by plugging into an electrical outlet with a ground terminal receptacle (U-ground outlet). Follow the National Electrical Code.
3. Intellect 700™ should not be connected to any other electrical device when in use.
4. **CAUTION:** Federal law restricts this device to sale by, or on the order of, a physician or licensed practitioner.
5. The generator should be routinely checked before each use to determine that all controls function normally; especially that the INTENSITY control does properly adjust the intensity of ultrasonic power output in a stable manner. Also determine that the TREATMENT TIME control does actually terminate ultrasonic output power when the control is turned to zero time (off).
6. **"CAUTION"** — Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous exposure to ultrasonic energy."

full one year warranty

Chattanooga Corporation ("Company") warrants that Intellect Model 700™ ("Product") is free of defects in material and workmanship.

This warranty shall remain in effect for one (1) year from the date of the original consumer purchase of this Product and extends to any owner of the Product during the warranty period. If this Product fails to function during the one year warranty period because of a defect in material and workmanship, Company or the selling dealer will replace or repair the Product without charge within a period of 30 days from the date on which the defective product is returned to the Company or the dealer. Company or the dealer will ship the replacement or the repaired product to the consumer's residence.

THIS WARRANTY DOES NOT COVER:

1. Replacement parts or labor furnished by anyone other than Company, the dealer or an approved Company service agent.
2. Defects or damage caused by labor furnished by someone other than Company, the dealer or an approved Company service agent.
3. Any malfunction or failure in the Product while it is in the possession of the owner during the warranty period if the malfunction or failure is not caused by a defect in material and workmanship or if the malfunction or failure is caused by the unreasonable use, including the failure to provide reasonable and necessary maintenance.

COMPANY SHALL NOT BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES TO PROPERTY OR BUSINESS

Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

TO OBTAIN SERVICE from Company or the selling dealer under this warranty, the owner must do or abide by the following:

1. A written claim must be made within the warranty period to Company or the selling dealer. If the claim is made to Company, the written claim should be sent to P.O. Box 4287, 101 Memorial Dr., Chattanooga, Tennessee, 37405.
2. The Product must be returned to Company or the selling dealer by the owner.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

Company does not authorize any person or representative to create for it any other obligation or liability in connection with the sale of this Product. Any representative or agreement not contained in the warranty shall be void and of no effect.


installation

Remove the Intellect Model 700™ and any additional items ordered with the unit from the carton and inspect for damage that may have occurred during shipment.


The following is a list of accessories that should be included with the unit as standard accessories:

Qty.	Description	Catalog No.
1	Hand-Held Probe	79001
1	Rectangular Applicator for Probe	79008
1	Sponge for Rectangular Electrode	79062
1	Spot Applicator for Probe	79005
1	Sponge for Spot Electrode	79059
1	Active Electrode Pad, 4" Red Round	72851
1	Active Electrode Pad, 4" Black Round	72850
1	Dispersive Electrode Pad 8" x 10"	72854
1	Active Lead, Red and Black (72")	72951
1	Dispersive Lead (72")	72955
1	Active Lead, Red Bifurcated (18")	72849
1	Active Lead, Black Bifurcated (18")	72855
2	Active Electrode Pads, 3" Red Round	72853
2	Active Electrode Pads, 3" Black Round	72852
1	Active Lead, Red (18")	72953
1	Active Lead, Black (18")	72954
2	Short Nylatex Straps 2.5" x 24"	10648
1	Long Nylatex Strap 2.5" x 48"	10832
2	Weight Bags	79036
1	Instruction Booklet	73297
1	Ultrasound Gel	72201

Check the voltage rating on the serial plate located on the back of the unit. Plug the unit into a 120 Volt or 220 Volt A.C. outlet as required. DO NOT attempt to use direct current. Follow the procedures indicated in the safety instructions. DO NOT attempt to use the unit if it is not properly grounded.





120 VOLT HOSPITAL GRADE



220 VOLT

Plug in 2-lead cable into the 3-terminal connector receptacle. Both of these receptacles are located at the front of the unit. Also, plug in the hand-held probe into the 4-terminal receptacle. These connectors will not mate incorrectly.



Plug the ultrasound applicator cable into the cable connector on the right front of the cabinet. Make sure the threaded connector is securely tightened. Slide the rubber covering (boot) over the metal connector.

By activating the timer, the functions of the stimulator and ultrasound can be checked out as per operating instructions in the following sections.

control pan

PULSES/SECOND CONTROL - This knob selects the pulse rate, from one pulse per second (1 pps) to a maximum of 120 pulses per second (120 pps).

INTENSITY METER - This meter indicates the intensity of either voltage or current, as selected. Voltage is read on the upper scale from 0 to 500 volts. Peak current is read on the lower scale from 0 to 2500 milliamps.

METER SELECTOR BUTTON - Depressing this button selects whether the meter reads in volts, as indicated by the voltage readout light (green); or milliamps, as indicated by the current readout light (red).

INTENSITY/RESET CONTROL - Any time the Function Switch is operated, e.g. switching from Probe to Continuous, the Intensity Control must be reset by rotating counterclockwise until a click is heard (or felt) at the end of rotation. Output will remain at zero unless this control is reset properly.

INTENSITY BALANCE - This knob controls relative individual electrode pad intensity IN THE RECIPROCATATE MODE ONLY. If the knob is turned to black, the black active electrode becomes relatively stronger by lessening the intensity of the red electrodes and vice versa. At the beginning of each treatment the knob should be centered to equalize pad intensities.


POLARITY SWITCH - This switch selects positive (+) or negative (-) polarity of the active electrode pads or probe.

SURGE MODE - The new independently timed SURGE ON and SURGE OFF functions allow for greater flexibility in application of high voltage stimulation. By setting the Surge On time to a longer or shorter period than the Surge Off time, the user can achieve various degrees of muscle fatigue. NOTE: The Surge On and Surge Off control knobs work only when the FUNCTION SELECTOR SWITCH is in the SURGE position. The Surge controls are inoperative when the Function Selector Switch is in the Probe, Continuous or Reciprocate 2.5, 5, or 10 second positions.

SURGE ON - Both electrode pads or sets of pads are activated for the time indicated around the knob from 0.5 seconds to 60 seconds. An "ON SURGE RAMP" up is included in this interval and the duration of the ramp is a 1 to 3 ratio of Surge On time selected. EXAMPLE: A Surge On time of 10 seconds creates a 3.3 second ramp up from 0 intensity to maximum intensity.

SURGE OFF - The rest period between surges is independently selected by this control. Surge Off intervals range from 0.5 seconds to 60 seconds.

TREATMENT TIMER - This control is both POWER ON/OFF switch and TIMER (i.e., a Timer).



FUNCTION SELECTOR SWITCH - This control knob selects the means of treatment—either with the hand-held probe or with the active electrode pads.

When the knob is in the "Probe" position, the output is through the hand-held probe. Output is on continuously and intensity is controlled remotely at the probe only.

In the Continuous and Reciprocate positions the output is through the electrode pads and intensity is controlled by the Intensity knob on the panel.


In the Continuous position both pads or sets of pads are on continuously.

The Reciprocate mode includes three selections: 2.5, 5, or 10 seconds. These three positions determine the amount of time that one of the two active electrodes is

el operation

wer Switch). This control functions to apply and .rove line input power to the generator circuitry. It electrically connects both sides of the incoming power line.

OUTPUT METER - This control is an analog-indicating meter which provides the user with accurate information pertaining to the level of ultrasound power available at the transducer at any given time during the treatment interval. The meter has one arc with dual calibrations of Power in watts and Intensity in W/cm^2 .



DUTY CYCLE SELECTOR - This control determines the type of ultrasound output waveform. The three (3) types of output available are: 100% duty cycle, 50% duty cycle, and 20% duty cycle. In the "100%" Mode, the ultrasound output is a CONTINUOUS sinusoidal waveform at a frequency of 1 MHz nominal.

In the "50%" Mode, the ultrasound output is pulsed at 100 PULSES/SECOND, with the MARK/SPACE ratio being 1. This produces an output pulse of 1 MHz ultrasound energy having a duration of 5 milliseconds, with an off-time of 5 milliseconds between pulses. The pulse shape is true-rectangular in nature, to allow precise measurement of ultrasound intensity.

In the "20%" Mode, the ultrasound output is again pulsed at 100 PPS, with the MARK/SPACE ratio being 4, to produce rectangular pulses of 2 milliseconds duration with an off-time of 8 milliseconds between pulses.

INTENSITY/RESET CONTROL - This control is continuously variable to produce the desired level of ultrasound power from 1 to 20 watts, corresponding to INTENSITIES from 0.1 WATT/CM² to 2.4 W/CM² as indicated on the OUTPUT METER.

HEAD MAX. TEMP. CONTROL - This indicator is included as an additional measure of patient safety.

While the head temperature remains at 140° F. or above, the generator output remains completely off and the HEAD MAX. TEMP. indicator light remains on. The OUTPUT METER indication returns to zero also, to clearly indicate the condition of zero output energy.

on while the other is off. At each 2.5, 5, or 10 second interval, as selected, the electrode pad that is on will switch off and the off pad will switch on. This alternate switching continues throughout the total treatment period as set on the timer.

In the Surge mode both pads or sets of pads are cycled on and off together at rates selected by the Surge On and Surge Off controls.

In the Ultrasound mode you may select pulsed modes (20% or 50% duty cycle), or the continuous (100% duty cycle) mode of operation.

In the Combination mode both ultrasound and high voltage stimulation are transmitted by the applicator head.

A. INDICATIONS FOR THERAPY

The Intellect 700™ is a combination ultrasound and electrical muscle stimulator designed to deliver therapeutic deep heat and muscle stimulation. These treatment modalities can be delivered simultaneously in the combination mode or separately in their respective modes.

B. CONTRAINDICATIONS OF ULTRASOUND THERAPY ❖

Ultrasound should not be used over the eyes or the reproductive organs.

Ultrasound should not be used over a pregnant uterus.

Other contraindications include acute infection or sepsis, deep vein thrombosis, or arterial disease, and over anesthetized areas or conditions that cause impairment of sensations, such as chemotherapy.

Ultrasound is not to be used over cancerous lesions.

C. CONTRAINDICATIONS OF HIGH VOLTAGE THERAPY

This device should not be used on patients with cardiac pacemakers.

This device should not be used over the carotid sinus area.

This device is not to be used transcranially.

This device should not be used to relieve pain syndromes until etiology has been established.

This device should not be used over a pregnant uterus.

This device should not be used over or near cancerous malignancy.

❖ REF: Lehmann, J.F., Therapeutic Heat and Cold; 13: 367 - 378; 1972.

D. OPERATING PROCEDURE

This section is divided into three sections: high voltage, ultrasound and combination.

a. HIGH VOLTAGE OPERATION

The controls for this mode are located on the left side of the front panel.

1. Dial the Function Selector Switch in the lower left of the panel to the desired setting: either Probe, Continuous, Surge, Reciprocate 2.5 secs., Reciprocate 5 secs., or Reciprocate 10 secs.
2. Set the polarity switch to Positive (+) or Negative (-).
3. Dial the Pulses/Sec to the determined rate.
4. Attach the moistened dispersive pad to the patient.
5. Attach the moistened active pad or pads to the patient unless you are going to use the probe with the roller electrodes. If you are using the roller electrodes use a conductive gel on the area of the patient where high voltage is going to be used.
6. Turn Treatment Time Knob to the determined total time of the treatment. Intensity will remain at zero if treatment timer is not turned on.
7. Turn Intensity Balance control to center position.
8. Turn Intensity control counterclockwise to below Min to the Reset position where a click will take place. After this, slowly turn clockwise until patient feels sensation, and continue to desired level of tolerance.

9. If patient has more sensory perception of electrical stimulation in one active electrode pad or pads over the others, or if a stronger muscle contraction occurs under one electrode pad or pads and not the other, then the Intensity Balance control can be used to readjust relative strengths with the mode selector switch in Reciprocate 2.5, 5 or 10 only. The active pads are color-coded for easy identification. The Intensity Balance control will adjust the relative intensity of each pad individually, by turning the Intensity Balance control in the direction corresponding to the weaker color pad. Example: If the red pad feels weaker than the black pad, turn the Intensity Balance control counterclockwise toward the red a slight amount, then readjust the Intensity control as permitted by patient tolerance.
10. When the treatment is complete, remove pads or pad, return probe (if used) to its holder and return intensity control to reset.

b. ULTRASOUND OPERATION

The controls for this mode are generally located on the right side of the control panel.

1. Set the function selector switch to the ultrasound position.
2. Set INTENSITY control fully counterclockwise to the RESET position then the MIN position. A normal, low output of approximately one watt is present, as indicated by the OUTPUT meter.
CAUTION: Do not attempt to adjust the OUTPUT meter mechanical zero when the unit is on.
3. Select desired mode of ultrasound output. PULSED ultrasound in the 20% or 50% DUTY CYCLE position, or CONTINUOUS ultrasound in the 100% DUTY CYCLE position of the selector switch.
4. Turn TREATMENT TIME control to the desired preset time. For settings of 5 minutes or less, turn knob past 5 minutes and then return to the desired time.
5. Begin normal treatment by applying INTELECT ULTRASOUND GEL to the treatment area of patient's body. Contact the applicator to the patient's body with firm, uniform pressure. Adjust the desired treatment INTENSITY while continuously moving the applicator over the affected area. Do not maintain the applicator stationary over any given area for extended time periods. This may result in hazardous exposure to ultrasonic energy.
6. Continue procedure described in Step 5 above for duration of treatment time. Insure adequate coupling by adding sufficient gel as required. Inadequate coupling is apparent by noting variations in the intensity meter indication while moving the applicator.
7. When treatment is complete, return setting of INTENSITY control to RESET, then store applicator in the holder.

c. COMBINATION OPERATION

1. Set the function selector switch to the combination position.
2. Attach the moistened dispersive pad to the patient.
3. Set the polarity switch to Positive (+) or Negative (-).
4. Dial the Pulses/Sec to the determined rate.
5. Turn Intensity Balance control to center position.
6. Turn both the high voltage and the ultrasound intensity switches to the RESET position.

7. Select desired mode of ultrasound output. PULSED ultrasound in the 20% or 50% DUTY CYCLE position, or CONTINUOUS ultrasound in the 100% DUTY CYCLE position of the selector switch.
8. Turn TREATMENT TIME control to the desired preset time. For setting of 5 minutes or less, turn knob past 5 minutes and then return to the desired time.
9. Begin normal treatment by applying INTELECT ULTRASOUND GEL to the treatment area of patient's body. Contact the applicator to the patient's body with firm, uniform pressure. Adjust the desired treatment INTENSITIES, ✕ first the Ultrasound and then the High Voltage, while continuously moving the applicator over the affected area. Do not maintain the applicator stationary over any given area for extended periods. This may result in hazardous exposure to ultrasonic energy.
- ✕ It is possible that high voltage stimulation may have significant pain reduction effect, to the point where the patient may have no appreciable response to an over-dosage condition of the ultrasound energy. Therefore, it is suggested to apply the ultrasound first, then the stimulation.
10. Continue procedure described in Step 9 above for duration of treatment time. Insure adequate coupling by adding sufficient gel as required. Inadequate coupling is apparent by noting variations in the intensity meter indication while moving the applicator.
11. When treatment is complete, return settings of INTENSITY controls to RESET, then store applicator in the holder.

E. MAINTENANCE PROCEDURES

The following equipment must be available to perform service and calibration adjustments outlined below in Section E.

1. RADIOMETER; 0-20 WATTS minimum (ULTRASONIC WATTMETER, for use with applicator coupled through water).
2. WIDE-BAND OSCILLOSCOPE; 10 MHz minimum frequency response.
3. LO-CAPACITANCE OSCILLOSCOPE probe.
4. FREQUENCY COUNTER; 1% accuracy or better, required. Range at least 10 MHz.
5. VOLTMETER; AC and DC ranges, at least 20 Kilohms/Volt. DVM preferred.
6. LEAKAGE CURRENT TESTER - 100 μ A. AC-RMS range, 50/60 Hz.
7. HIGH POTENTIAL TESTER - 2500 VOLTS RMS.
8. Additional requirement - source of de-oxygenated water (< 5ppm) for use in the wattmeter, as required. Distilled, degassed water is to be used when measuring ultrasound power output.

F. SERVICE INSTRUCTIONS, ULTRASOUND

To fully maintain compliance with Federal Regulation Title 21 (21CFR) this unit must be recalibrated annually. It is also recommended that all Chattanooga Corporation ultrasonic products be returned to the factory or authorized servicing dealer for repairs or recalibration (see SECTION D. for equipment requirements.)

Measure the ultrasound power output utilizing a suitable wattmeter of known accuracy. If required, recalibration can be attained in the following sequence.

1. Couple the transducer through water to the wattmeter. Place front panel controls as follows: TREATMENT TIME to 30 minutes; DUTY CYCLE selector to 100%; INTENSITY control to MAX.
2. Adjust TEMP. CAL. - P8 fully clockwise (located on the CONTROL BOARD).

3. Adjust INTENSITY CAL. - P1 for maximum voltage at C4, negative side. Reduce this maximum setting by 2 VOLTS.
4. Set DUTY CYCLE selector to 50%. Connect oscilloscope probe (10:1) to J3-pin 3. Adjust scope sweep for one full cycle of the 100 PPS nominal pulse repetition rate. Adjust 50% DUTY CYCLE calibration P7 for exactly 50% DUTY CYCLE of the pulse.
5. Set DUTY CYCLE selector to 20%. Adjust 20% DUTY CYCLE calibration P6 to exactly 20% duty cycle.
6. Adjust OSCILLATOR TUNING adjustment for maximum output as indicated by the external wattmeter. Set DUTY CYCLE selector to 100%. Set indicated wattmeter power to 20 watts using the front panel INTENSITY control. Calibrate the front panel OUTPUT meter to full scale using the 100 M. CALIBRATION control P3.
7. Set the front panel INTENSITY control fully clockwise and the incoming line voltage at 120 VAC. Adjust the INTENSITY CAL. - P1 for a full scale indication on the OUTPUT meter.
8. Set DUTY CYCLE selector to 50%. Adjust indicated wattmeter power to 10 watts using the INTENSITY control. Calibrate OUTPUT meter to full scale using the 50 M. CAL - P4.
9. Set DUTY CYCLE to 20%. Adjust indicated wattmeter power to 4 watts using INTENSITY control. Calibrate OUTPUT meter to full scale using 20 M. CAL - P5.
10. Place the applicator in a 140 degrees Fahrenheit water bath. Adjust the TEMP. CAL - P8 until the HEAD MAX TEMP. light just comes on. This calbrates the temperature at which ultrasound power is terminated should the transducer overheat for any reason. The ultrasound generator is now fully calibrated, and ready for normal use.

"CAUTION: Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous exposure to ultrasonic energy."

G. SPECIFICATIONS, ULTRASOUND

Frequency - 1.0 MHz \pm 5% *950,000 - 1,050,000*

Duty Cycle - 100% (continuous mode)

50% \pm 10% (pulse mode)

20% \pm 10% (pulse mode)

Pulse Repetition Rate - 100 Hz \pm 20%

Ultrasonic Power - Variable from 1 watt to 20 watts.

Output Meter Accuracy - \pm 20% (for any output above 10% of maximum)

Temporal Peak/Average Intensity Ratio - 2:1 \pm 20% for 50% D.C.

5:1 \pm 20% for 20% D.C.

Output:

1. Continuous - 1 MHz signal that is on as long as the timer is running.
2. Pulse - 1 MHz signal modulated 100% by the 100 Hz rectangular wave with the selected Duty Cycle.

Timer Accuracy:

1. Less than 0.5 seconds for settings less than 5 minutes
2. 10% for settings from 5 minutes to 10 minutes
3. 1 minute for settings greater than 10 minutes

Applicator:

1. Effective radiating area - $8.5 \text{ cm}^2 \pm 1.5 \text{ cm}^2$
2. Maximum beam non-uniformity ratio - 6.0:1
3. Beam type - Collimating

✱ Input Power Requirements:

(Domestic) 120V/60 Hz $\pm 10\%$, 1.25 Amps

(Export) 220V/50 Hz $\pm 10\%$, .8 Amps


✱ Voltages in parentheses in the procedure are for 220VAC units—e.g., 108(198)VAC indicates 108VAC for 120V unit, and 198VAC for a 220V unit.

H. DESCRIPTION OF ULTRASONIC FIELD

The spatial distribution of the radiated field is essentially a collimated beam of ultrasonic energy having a cross-sectional area of 9.5 CM^2 when measured at a point 5 millimeters from the transducer face.

The energy distribution within the radiated field is 2.0 W/CM^2 maximum, and takes a generally conic shape having decreasing intensity at progressively increasing distance from the face of the transducer.

This field distribution applies for the radiation emitted into the equivalent of an infinite medium of distilled degassed water at 30° C. and with line voltage variations in the range of ± 10 percent of 120 VOLTS RMS.



I. HIGH VOLTAGE SYSTEM DESCRIPTION

Pulse Charge - approximately 14 microcoulombs per twin peak.

Pulse Frequency - One pulse per second to 120 pulses per second.

Output - 0 to 500 volts.

Output Current - 0 to 2500 milliamps peak.


Average Current - Approximately 1.5 milliamps at 120 pulses per second.

Area of Conducting Surface of Electrodes

4-inch active electrodes: 3.75" dia. = 11.04 in² (7.125mm²)


3-inch active electrodes: 3.12" dia. = 7.64 in² (4,932mm²)

Dispersive electrodes: 7.5" x 9.5" = 71.25 in² (46,000mm²)



SECTION J.

EXPLODED VIEW




SECTION J.

REF. NO.	PART NO.	DESCRIPTION	QTY.
1	70882	Hinge Intelect 700	1
2	70881	Hinge Intelect 700	1
3	72098	Cover Dust Intelect 700	1
4	70879	Bracket Power Supply	1
5	60157	Cord Set 183 SJT	1
6	79102	Holder, Cord	2
7	70430	Ultrasound Applicator	1
8	90678	Probe Holder	1
9	79113	Holder Probe Intelect 500	1
10	79088	Cover End Right Side	1
11	79089	Cover End Left Side	1
12	70919	Cabinet Intelect 700	1
13	73247	Insulator Black Rubber	1

SECTION K.

EXPLODED VIEW



SECTION K.

REF. NO.	PART NO.	DESCRIPTION	QTY.
1	70335	Knob, Timer Rhodes	1
2	70334	Knob Alco	4
3	70436	Meter, Ultrasound	1
4	72139	Timer Rhodes, 30 min.	1
5	72355	Meter, Intellect 500	1
6	72740	Switch Rotary CTS I-700	1
7	70127	Board PC Intellect 200-B1 (See page 15 for Board Layout)	1
8	72687	Switch Rotary Surge Time	2
9	72518	Board PC Intellect 500	1
10	72806	Oscillator Sub-Assembly	1
11	72437	Transformer 241-5-20	1
12	70417	Transformer 241-5-16	1
13	70415	Transformer 241-8-1757	1
14	71766	Fuse .25 Amp ABC	2
15	70843	Fuse MDL 1¼	2
16	72097	Knob Alco	5

CONTROL BOARD LAYOUT

200-B1



CONTENTS

	Page
Ultrasound Specifications	17
Ultrasound Calibration Procedure	17 - 21
Ultrasound Calibration Checklist	22
Stimulator Calibration Procedure	22 - 25
Stimulator Calibration Checklist	26
Logic Truth Tables	27
Timing Diagram for Surge Circuit	28
Timing Diagram for Pulse Generation.....	29
Surge On and Off Indicator	30
Current Transformer Adapter.....	31
Intelect 200 & 700 Calibrating Adapter.....	32

ULTRASOUND SPECIFICATIONS

Frequency - 1.0 MHz \pm 5%

Duty Cycle - 100% (continuous mode)
50% \pm 10% (pulse mode)
20% \pm 10% (pulse mode)

Pulse Repetition Rate - 100 Hz \pm 20%

Ultrasonic Power - Variable from 1 watt to 20 watts

Ultrasound Output Meter Accuracy - \pm 20% (for any output above 10% of maximum)

Temporal Peak/Average Intensity Ratio - 2:1, \pm 20% for 50% D.C.
5:1, \pm 20% for 20% D.C.

Output:

1. Continuous - 1MHz signal that is on as long as the timer is running
2. Pulse - 1MHz signal modulated 100% by the 100Hz rectangular wave with the selected Duty Cycle

Timer Accuracy:

1. Less than 0.5 minutes for settings less than 5 minutes
2. 10% for settings from 5 minutes to 10 minutes
3. 1 minute for settings greater than 10 minutes

Applicator:

1. Effective radiating area - 8.5 cm² \pm 1.5 cm²
2. Maximum beam non-uniformity ratio - 6.0:1
3. Beam type - collimating

*Input Power Requirements:

(Domestic) 120V/60Hz \pm 10%, 1.25 amps
(Export) 220V/50Hz \pm .8 amps

*Voltages in parentheses in the procedure are for 220VAC units; e.g., 108(198)VAC indicates 108VAC for 120V unit and 198VAC for a 220VAC unit.

ULTRASOUND CALIBRATION

TEST EQUIPMENT REQUIRED

1. Power line monitor (expanded scale voltmeter for rated line voltage \pm 10%), VIZ model WV-120B or equivalent for 120VAC line.
2. Autotransformer, adjustable from 90% to 110% of rated line voltage, 150 watts or greater.

3. Ultrasound Power Meter, Ohmic Instruments Model UPM-30 or equivalent.
4. Oscilloscope, Tektronix T922 or equivalent.
5. Probe, voltage, X10, Scope, low capacitance.
6. Probe, current, Hewlett Packard Model 1110 or equivalent
7. Voltmeter, Digital, 3-1/2 digits, Simpson Model 461 or equivalent
8. Probe, temperature, Fluke Model 80T-150 or equivalent.
9. Source of approximately 1/2 gallon of distilled de-oxygenated (<5 PPM) water at 30 degrees Celsius for use in UPM-30 power meter (Item 3).
10. Counter, frequency, 10MHz, Triplet 7000 or equivalent.
11. Stopwatch, Siliconix Model 705 or equivalent.
12. Applicator current transformer adapter.

PROCEDURE:

I. Instrument Preparation:

- A. Make certain that all power is removed from the cabinet.
- B. Disconnect the applicator cable from the front cabinet connector.
- C. Remove the five #8 truss head screws that attach the front panel to the cabinet.
- D. Tilt the top of the front panel toward the front of the cabinet.
- E. Remove the following connectors which attach the front panel to the circuitry inside the cabinet: 16 pin connector with ribbon cable, 20 pin connector with ribbon cable, 9 pin power connector, 2 each 5 pin flat connectors and 9 pin flat connector.
- F. Remove front panel from cabinet.
- G. Disconnect the following connectors which attach power supply bracket to cabinet: 2 pin flat connector, 12 pin flat connector, coaxial jumper with UHF connector on each end and the strain relief on the line cord.
- H. Remove the eight #8 flat head phillips screws which attach the power supply bracket to the cabinet, in the bottom of the cabinet.
- J. Remove the power supply bracket from the cabinet.
- K. Reattach the connectors which are coming from the front panel to the circuitry on the power supply bracket.
- L. Connect applicator current transformer adapter to the UHF connector on the oscillator box and connect the applicator cable to the applicator current transformer adapter.

II. Power Supply and Duty Cycle Adjustment:

- A. Set the front panel controls:
 ULTRASOUND INTENSITY - fully CW
 DUTY CYCLE - 100% (continuous)
 TREATMENT TIME - 0

- B. Pre-set the internal adjustments:
 INTENSITY CAL. (P1) - fully CW
 50% D.C. CAL. (P7) - center of travel
 20% D.C. CAL. (P7) - center of travel
 100M CAL. (P3) - fully CCW
 50M CAL. (P4) - fully CCW
 20M CAL. (P5) - fully CCW
 TEMP. CAL. (P8) - fully CW

- 50% INT. CAL. (P9) - fully CCW
 20% INT. CAL. (P10) - fully CCW

C. Connect the test set-up in Figure 1:

1. Set AC input voltage with the autotransformer to 120(220)VAC line monitor.
2. See OHMIC INSTRUMENTS clinical engineering notes AN-330 for operation of the UPM-30 U.S. Power Meter.
3. Connect the digital voltmeter between common, (+) terminal of BR1, and either side of resistor R6 to measure 0 - 70VDC.




FIGURE 1

D. Set the TREATMENT TIME to either 30 minutes or lock down the dial knob so that the unit is on continuously.

1. Observe:
 - (a) Ultrasound Power Meter on Intellect 700 moves up-scale from zero.
 - (b) Digital voltmeter indicates between -5VDC and -70VDC.

E. Adjust INTENSITY CAL. (P1 on the control board) in the CCW direction until the digital voltmeter indicates a maximum voltage (approx. -70VDC) and starts to decline. Turn P1 in the CW direction past maximum voltage to the position where the voltmeter indicates 2 volts less than the maximum voltage reading.

F. Connect the voltage probe from the oscilloscope to pin #3 of U3 on the control board.

2. Observe: Voltage displayed on the scope is less than +1VDC.

G. Set the DUTY CYCLE selector switch on the front panel to 50%.

1. Observe: A rectangular wave 100 Hz \pm 20% signal operating between 1V and +12V approx.

H. Set the time base on the oscilloscope for 10.0 divisions horizontal display of one cycle of the 100 Hz signal.

1. Adjust 50% D.C. CAL. on the control board for 5.0 divisions horizontal display of the signal at the +12V level.

I. Set the DUTY CYCLE selector switch to 20%.

1. Observe: Rectangular wave 100 Hz \pm 20% signal operating between 1V and +12V.

J. Repeat Step H above; preset time base only.

K. Adjust 20% D.C. CAL. on the control board for 8.0 divisions horizontal display of the signal at the +12V level.

L. Remove oscilloscope probe and voltmeter.

III. Oscillator Peaking and Output Power Adjustment:

A. Set the front panel controls:

ULTRASOUND INTENSITY - fully CW

DUTY CYCLE - 100%

TREATMENT TIME - fully CW (ON)

- B. Connect the voltage probe from the oscilloscope to the oscillator output at the applicator current transformer adapter.
- C. Connect the current probe from the scope around the center conductor of the applicator current transformer adapter.
 1. Observe: With the scope set for dual channel display, the voltage and current waveforms should not be more than 90 degrees out of phase. If the phase displacement is more than 90 degrees, reverse the current probe connection or invert Channel 2 display on the scope. The displacement should now be less than 90 degrees.
- D. Align the U.S. applicator in the UPM-30 per the power meter manufacturer's notes, AN-330.
- E. Adjust OSC TUNING ADJUSTMENT (tuning slug in transformer T4) through access in bottom of the oscillator enclosure for a peak power indication on the UPM-30.
 1. Observe:
 - (a) The voltage and current waveforms on the scope should now be in phase.
 - (b) The peak power reading should be 23 watts or greater.
- F. Connect the frequency counter to the oscillator output (applicator current transformer adapter) and measure the oscillator frequency.
 1. Observe: The frequency should be $1.0 \text{ MHz} \pm 5\%$.
- G. Set the DUTY CYCLE selector on the front panel to 50% and vary the INTENSITY control on the front panel from fully CCW to fully CW.
 1. Observe:
 - (a) The voltage and current waveforms on the scope remain in phase over the full intensity range.
 - (b) A modulation envelope of the 1 MHz output modulated 100% by the 100 Hz rectangular wave signal with 50% duty cycle.
- H. Set the DUTY CYCLE selector on the front panel to 20% and vary the INTENSITY control on the front panel from fully CCW to fully CW.
 1. Observe:
 - (a) The voltage and current waveforms on the scope remain in phase over the full intensity range.
 - (b) A modulation envelope of the 1 MHz output modulated 100% by the 100 Hz rectangular wave signal with 20% duty cycle.
- I. Disconnect the voltage probe, current probe and frequency counter from the ultrasound output.
- J. Set the DUTY CYCLE selector to 100% and the Ultrasound INTENSITY control fully CW.
- K. Check the peak power on the UPM-30 and readjust the OSC TUNING ADJUSTMENT (tuning slug in transformer T4) if necessary for a maximum power indication.
- L. Set the Ultrasound INTENSITY control on the front panel for a power indication of 20 watts on the UPM-30 power meter.
 1. Adjust 100M CAL. (P3 on the I-200 B1 Board) for a 20W indication on the I-200 output meter on the front panel.
- M. Rotate the INTENSITY control fully CW.
 1. Adjust INTENSITY CAL. (P1 on the I-200 B1) for an indication of 23 watts on the I-700 ultrasound output power meter.
- N. Increase the AC input voltage to the INTELECT 700 from 120(220)VAC to 132(242)VAC.
 1. Observe: The power indicated on the UPM-30 power meter should be no more than 28.5 watts.

- O. Check the accuracy of the I-700 ultrasound output meter by comparing its indication to the power indicated on the UPM-30 with the line voltage at 108(198)VAC, 120(220)VAC, and 132(242)VAC.
1. Observe: The following readings at all line voltages:

I-700 Meter Reading Peak Power	UPM-30 Indication Average Power	
	Minimum	Maximum
2.0 W	1.72 W	2.28 W
5.0 W	4.30 W	5.70 W
10.0 W	8.60 W	11.40 W
15.0 W	12.90 W	17.10 W
20.0 W	17.20 W	22.80 W

- P. Return the line voltage to 120(220)VAC and set the DUTY CYCLE selector on the front panel to 50%.

- Q. Set the Ultrasound INTENSITY control on the front panel for a power indication of 10 watts on the UPM-30 power meter. *P4*

1. Adjust 50M CAL. (*P5* on the control board) for a 20 watt reading on the I-200 output power meter.
NOTE: UPM-30 meter indication = DUTY CYCLE (50%) x I-200 reading. Average power is indicated on the UPM-30 and peak power is indicated on the I-200 meter in either of the pulse modes (50% or 20% D.C.).

- R. Rotate the Ultrasound INTENSITY control fully CW.

1. Adjust the 50% INT. CAL. (P9 on the pulse int. cal. board, located on the back of the meter) for a reading of 23 watts on the I-200 output meter.

- S. Check the accuracy of the I-700 output meter by comparing its indication to the power indicated on the UPM-30 with the line voltage at 108(198)VAC, 120(220)VAC, and 132(242)VAC.

1. Observe: The following readings at all line voltages:

I-700 Meter Reading	UPM-30 Indication	
	Minimum	Maximum
2.0 W	.86 W	1.14 W
5.0 W	2.15 W	2.85 W
10.0 W	4.30 W	5.70 W
15.0 W	6.50 W	8.50 W
20.0 W	8.60 W	11.40 W

- T. Return the line voltage to 120(220)VAC and set the DUTY CYCLE on the front panel to 20%.

- U. Set the Ultrasound INTENSITY control on the front panel for a power indication of 4 watts on the UPM-30 power meter.

1. Adjust the 20M CAL. (P5 on the I-200 B1) for a 20W reading on the I-700 ultrasound output meter.
NOTE: UPM-30 indication = DUTY CYCLE (20%) x I-700 reading.

- V. Rotate the ultrasound INTENSITY control on the front panel fully CW.

1. Adjust the 20% INT. CAL. (P10 on the PULSE INT. CAL. board, located on the back of the meter) for a reading of 23 watts on the I-700 ultrasound output meter.

- W. Check the accuracy of the I-700 ultrasound output meter by comparing its indication to the power indicated on the UPM-30 with the line voltage at 108(198)VAC, 120(220)VAC and 132(242)VAC.
1. Observe: The following readings at all line voltages:

I-700 Meter Reading Peak Power	UPM-30 Indication Average Power	
	Minimum	Maximum
2.0 W	.34 W	.46 W
5.0 W	.86 W	1.14 W
10.0 W	1.72 W	2.28 W
15.0 W	2.58 W	3.42 W
20.0 W	3.44 W	4.56 W

- X. Disconnect the test setup.

IV. Adjustment of Head Maximum Temperature Trip.

- A. Connect the I-700 line cord to a 120(220)VAC power source.
- B. Place the applicator in a 140° F. water bath and monitor the temperature with the temperature probe or a thermometer that will accurately indicate 140° F.
- C. Set the front panel controls on the I-700 in the following positions:
TREATMENT TIME - 30 minutes
DUTY CYCLE - 100%
ULTRASOUND INTENSITY - 10 W indication on the output meter
- D. Allow at least 5 minutes warm-up time after the front panel controls are set and power is applied.
- E. Adjust the TEMP. CAL. (P8 on the 200 B1) in the CCW direction until the HEAD MAX. TEMP. light on the front panel just turns on.
- F. Remove the applicator from the bath and allow the head to cool until the HEAD MAX. TEMP. light turns off.
- G. Replace the applicator in the water bath at 140° F. for 5 minutes. If the HEAD MAX. TEMP. light does not turn off, repeat Step E.
- H. Remove applicator from water bath and disconnect all equipment.

V. Timer Accuracy Check:

- A. Connect the Intellect 700 to a 120(220)VAC power source.
- B. Set the timer to 2.5 minutes and start the test stopwatch as the timer knob is released. (Turn past 5 and set to 2.5.)
 1. Observe: Time indicated on the stopwatch when the timer interrupts power to the unit is between 2.0 minutes and 3.0 minutes.
- C. Set the timer to 5 minutes and start the stopwatch as the timer knob is released.
 1. Observe: Time indicated on the stopwatch when the timer interrupts power to the unit is between 4.5 and 5.5 minutes.
- D. Set the timer to 10 minutes and start the stopwatch as the timer knob is released.
 1. Observe: Time indicated on the stopwatch when the timer interrupts power to the unit is between 9.0 and 11.0 minutes.
- E. Set the timer to 30.0 minutes and start the stopwatch as the timer knob is released.
 1. Observe: Time indicated on the stopwatch when the timer interrupts power to the unit is between 29.0 and 31.0 minutes.

ULTRASOUND CALIBRATION AND TESTING CHECKLIST

	BY	DATE
1. Adjust power supply output.	_____	_____
2. Set duty cycles.	_____	_____
3. Tune oscillator.	_____	_____
4. At 120 volts adjust power as follows:		
100%	_____	_____
50%	_____	_____
20%	_____	_____
5. Adjust intensity calibration at 108 volts (line).	_____	_____
6. Check and record maximum intensity at 132 volts line:		
Setting Reading		
100%	_____	_____
50%	_____	_____
20%	_____	_____
7. Check and record meter accuracy at the following (use 120V line):		
Setting Reading		
1 W	_____	_____
5 W	_____	_____
10 W	_____	_____
15 W	_____	_____
20 W	_____	_____
8. Adjust head temperature to 140° F.	_____	_____
9. Check and record timer accuracy at the following times:		
Time Set Actual Time		
15 min	_____	_____
10 min	_____	_____
5 min	_____	_____
1 min	_____	_____

HIGH VOLTAGE STIMULATION CALIBRATION PROCEDURE

EQUIPMENT NEEDED:

1. Oscilloscope - 10 MHz or greater bandwidth with 10:1 probe.
2. Frequency Counter with period function (1 sec.) and time interval function (10 sec.); i.e., Simpson 2725.
3. Digital Voltmeter (DVM).
4. Stop Watch.
5. Surge On and Off Indicator.

PRECALIBRATION SET-UP:

Adjust the following controls to the designated positions before applying power to the unit under test:

Designation	Description	Position
I-500 B1 Board:		
R90	Voltage Regulator Adjust	CCW
R2	1 Hz Adjust	Center
R1	120 Hz Adjust	Center
R18	Surge On Time Adjust	Center
R20	Surge Off Time Adjust	Center
R27	Surge Generator Zero Adjust	Center
R32	Surge Modulator Drive Adjust	N/A
R46	819 Hz Reciprocate Clock Adjust	Center
R14	Output Voltage Adjust	CCW
R68	Surge Modulator Zero Adjust	Center
R70	Meter Full Scale Adjust - Voltage	Center
R104	Meter Full Scale Adjust - Current	Center
R97	Meter Zero Adjust - Current	Center
Front Panel Controls:		
R118	Pulse Rate Control	60 pps
S7	Surge On Control	5 sec.
S8	Surge Off Control	5 sec.
R121	Intensity Control	Reset
S6	Treatment Timer	Zero
S5	Power Switch	Off
S2	Function Switch	Continuous
S3	Meter Range Selector	N/A
S4	Polarity Switch	Positive
R122	Intensity Balance Control	Center

CALIBRATION:

1. Power Supply Voltage Adjustment -

Set DVM on 20 VDC range. Attach negative lead of meter to circuit common (negative side of C34) and positive lead to +15 volt supply (anode of D13). Plug Intellect 700 in, turn Power Switch on, and turn Treatment Timer on. The DVM should read less than 15 volts. Adjust R90 clockwise until DVM reads 15.00 volts \pm .02 volts. Move positive lead of DVM to pin 14 of IC26 and check reference voltage. It should be 10 volts \pm .4 volts.

2. Meter Voltage Adjustment -

Set Meter Range Selector in "Volts" position. Adjust mechanical zero of meter with intensity at zero. Turn Intensity Control full clockwise. Adjust R70 until meter reads 500 volts. Also check collector of Q16 (Heat Sink) or emitter of Q17 with DVM. As intensity varies from 0 to 500 volts, this point should vary from 0 to 15 volts \pm 1 volt. Return Intensity Control to reset position. NOTE: When Intensity Control is CCW, knob should point to reset.

3. Reciprocate Time Adjustment -

Set frequency counter on 1 kHz range. Attach common lead to circuit common (negative side of C34) and attach other lead to TP1 (pin 3 of IC13). Adjust R46 until frequency counter reads 819 Hz \pm 20 Hz. NOTE: This signal is a 15 volt pk-pk rectangular waveform with approximately 50% duty cycle.

4. **Pulse Rate Adjustment -**

Set knob on shaft so that when control is CCW the pointer on the knob will point to the dot below 1 pps. Connect frequency counter to pin 3 of IC2 (or cathode of D1). Turn Pulse Rate Control to 120 pps. (Set pointer exactly on 120 pps). Adjust R1 until frequency counter reads 120 Hz. Put frequency counter in period mode. Turn Pulse Rate control CCW gradually. The period will increase as the control is rotated CCW until you reach the end of the electrical rotation, where the control will rotate further but the period will not increase. Then rotate the control slowly CW until the period just starts to decrease. Without changing the position of the shaft, set the pointer on the knob at exactly 1 pps on the dial and tighten the knob on the shaft. Set Pulse Rate control again in 120 pps position and adjust R1 for 120 Hz on counter. Return counter to period function and set Pulse Rate control to 1 pps and adjust R2 for a period of 1.00 seconds. Repeat 120 pps and 1 pps adjustments one more time, as there is a slight interaction between these adjustments. Also check that the Pulse LED is operating.

5. **Output Voltage Adjustment -**

Set Pulse Rate control at 80 pps. Attach oscilloscope ground to Dispersive Pad and 10:1 probe to Active Pads. Turn Intensity control fully CW (500 volts). Adjust R14 for 500 volts peak output. Check the waveshape of the double pulse and the 75 μ sec spacing between pulses. NOTE: Make sure that the oscilloscope probes are calibrated before this step. Then place a 200 ohm load across the output at 500 volts to check for single pulsing. If there is evidence of single pulsing, put oscilloscope into differential mode and attach inverted channel to Dispersive Pad and normal channel to Active Pads and check for single pulsing again with the 200 ohm load.

6. **Current Meter Calibration -**

Disconnect oscilloscope probe. Push Meter Selector Switch until meter is in peak current readout mode. Adjust R97 (Current Zero Adjust) for zero reading on meter. Make sure that Intensity is at zero volts when adjusting zero. Turn Intensity control up to 500 volts and place 500 ohm 1% resistor across Intellect output. Then adjust R104 for a 1000ma reading on meter. Turn intensity to zero.

7. **Surge Calibration -**

The surge circuit consists of two parts. The first part is a function generator which generates the complex surge waveform and will be abbreviated **Surge Function Generator (SFG)**. The second part is a modulator which causes the D.C. intensity level to change at the same rate as the surge function. This modulator will be abbreviated as the **Surge Modulator (SM)**. The SFG consists of Q3, 4, 5, 6, 8; IC7, 8, 9, 10, 11, 12; and Q24. The SM consists of Q7, IC18 and IC32.

The SFG will be calibrated first. Place unit in SURGE mode. R18 calibrates the SURGE ON time. Set SURGE ON time to 10 seconds and SURGE OFF time to 1 second. Attach time interval counter to TP2 and with counter on 10 second range, adjust R18 until counter reads 10.0 ± 0.1 seconds. Check 1 second time with the counter $\pm 5\%$. Then set SURGE OFF time to 10 seconds and SURGE ON time to 1 second. Attach time interval counter to TP3 and with counter on 10 second range, adjust R20 to calibrate the SURGE OFF time to 10.0 ± 0.1 seconds. Check 1 second time with the counter $\pm 5\%$. Then attach the SURGE ON and OFF Indicator to the unit under test and check the 60 second SURGE ON and SURGE OFF times with the stopwatch.

The surge function zero is adjusted by R27. Attach positive lead of DVM to the base of Q7 or pin 8 of IC7 and the negative lead to the negative side of C34. Set SURGE OFF time to 10 seconds and SURGE ON control to 1 second. Then during the time when the SURGE OFF LED is on, adjust R27 until DVM just reads a null (less than 0.1 volt). Do not adjust R27 past this point as it will cause an error in calibration.

The SM is the second part of the surge calibration. Measure with the DVM the voltage at pin 7 of IC26 when the Intensity control is set at 500 volts. Record that voltage. Then attach DVM to pin 6 of IC32. Turn Intensity to zero and adjust R68 (Surge Modulator Zero Adjust) to a zero reading on DVM. Then turn Intensity control to 500 volts, SURGE ON to 10 seconds and SURGE OFF to 0.5 seconds. During SURGE ON time, as indicated by the SURGE ON and OFF indicator (at maximum voltage of ramp), adjust R32 (Surge Modulator Drive) until pin 6 of IC32 is the same as the voltage previously recorded at pin 7 of IC26.

This completes the adjustment phase of the calibration.

Make the following **Functional Checks**:

1. **Probe Mode -**

Place the unit into the Probe Mode. Check to see that the Probe Intensity control varies the voltage on the intensity meter and the output from zero to 500 volts. Check the operation of the Polarity Switch. With a 500 ohm load, check the operation of the current meter in the probe mode.

2. **Continuous Mode -**

Red and black pad LED's should be lighted. Check power-on-clear function. Check to see that intensity varies from zero to 500 volts. Check operation of current meter on the red and black pads with the 500 ohm load.

3. **Surge Mode -**

Red and black pad LED's should be lighted. Surge On LED should work. with 500 ohm load check to see that current meter follows the surge waveform.

4. **Reciprocate Mode: 2.5, 5 and 10 seconds -**

Visually check period of 2.5, 5 and 10 second reciprocate periods. Red and black pad LED's should light alternately starting with the red LED. Check the operation of the Intensity Balance control. The control should reduce the output into either pad from 500 volts to a maximum of 400 volts. Also check both red and black pad outputs on the current meter using the 500 ohm load.

5. **Reset -**

Check to see that stimulator reset operates properly between all modes except while in reciprocate modes.

6. **Combination Mode -**


Check Combination Mode.

HIGH VOLTAGE CALIBRATION AND TESTING CHECKLIST


Step:	BY	DATE
1. Power supply voltage adjust 15 volts \pm .02 volts. Reference voltage check 9.6 to 10.4 volts.	_____	_____
2. Meter voltage adjustment Collector of Q16 check 0 to 15 volts \pm 1 volt	_____	_____
3. Reciprocate time adjustment 819 Hz \pm 20 Hz	_____	_____
4. Pulse rate calibrate: 120 pps 1 pps Pulse indicator LED	_____ _____ _____	_____ _____ _____
5. Output voltage adjustment 500 volts peak Check waveform of output and 75 μ sec spacing of pulses	_____ _____	_____ _____
6. Current meter calibration: Zero 1000 mA	_____ _____	_____ _____
7. Surge function generator: On Time Off Time Surge Function Zero Surge Modulator: Zero Drive	_____ _____ _____ _____ _____	_____ _____ _____ _____ _____
Functional Checks:		
1. Probe Mode	_____	_____
2. Continuous Mode	_____	_____
3. Surge Mode	_____	_____
4. Reciprocate Mode	_____	_____
5. Reset Check	_____	_____
6. Combination Mode	_____	_____
Safety Tests:		
1. Hi-Pot 2500 volts	_____	_____
2. Ground Resistance Test	_____	_____
3. Leakage Current Test	_____	_____

APPENDIX


Truth Tables for logic used in Intelect 700




A	B	OUT
0	0	0
0	1	0
1	0	0
1	1	1




A	B	OUT
0	0	0
0	1	1
1	0	1
1	1	1




or




A	B	OUT
0	0	1
0	1	1
1	0	1
1	1	0




or



A	B	OUT
0	0	1
0	1	0
1	0	0
1	1	0





A	B	OUT
0	0	0
0	1	1
1	0	1
1	1	0




INPUTS		OUTPUTS			
13	10	15	14	3	2
1	1	1	0	0	0
0	1	0	1	0	0
0	0	0	0	1	0
1	0	0	0	0	1

TIMING DIAGRAM FOR SURGE CIRCUIT






SURGE ON AND OFF INDICATOR (TEST ADAPTER)



**APPLICATOR CURRENT
TRANSFORMER ADAPTER
FOR CALIBRATING INTELECT 200 & 700**



REVIEWS

REV.	E.N. NO.	DATE	DESCRIPTION	ENGINEER	APPROVAL
F	1102	5-5-83	ADD SURGE TO INT. 700	D. BLEY	
G	1110	6-6-83	CORRECTION TO E.N. 1102	D. BLEY	

ORIGIN BLOCK #	CIRCUIT DESIGNATION	DESCRIPTION OR NAME OF SIGNAL	DESTINATION BLOCK #	CIRCUIT DESIGNATION
①	IC 22 PIN 6	SURGE	②	IC 25 PIN 6
①	IC 4 PIN 3	SURGE ENABLE	③	IC 5 PIN 13
①	IC 17 PIN 4	RESET SURGE	③	R150 10K
①	R112 100K	PROBE ENABLE	②	IC 19 PIN 6
①	IC 24 PIN 11	BLACK PAD ENABLE	②	IC 25 PIN 13
①	IC 24 PIN 10	RED PAD ENABLE	⑦	R58 15K
①	IC 16 PIN 11	INTENSITY SQUELCH	⑦	R52 15K
①	IC 15 PIN 4	BLACK + RED SELECT	⑦	R55 15K
①	IC 17 PIN 11	IFF RED SELECT	②	R143 8.2K
①	IC 17 PIN 10	IFF BLACK SELECT	②	IC 19 PIN 13
①	IC 17 PIN 3	RESET CONTINUOUS/RECIPROCAT	②	IC 19 PIN 12
②	IC 26 PIN 7	INTENSITY VOLTAGE SIGNAL	②	IC 19 PIN 5
②	Q16 COLLECTOR	0-15V INTENSITY PROGRAMMING	②	IC 25 PIN 5
③	IC 32 PIN 6	SURGE INTENSITY OUT	③	R63 10K 1%
④	IC 6 PIN 3	CHARGE PULSE	⑥	R69 18K
④	IC 34 PIN 3	TRIGGER PULSE	⑤	Q17 EMITTER
⑤	J5 PIN 13	(+) TWIN PULSE	②	IC 19 PIN 9
⑤	J5 PIN 15	(-) TWIN PULSE	⑤	R16 4.7K
⑤	T4 (●)	CURRENT SIGNAL	⑤	R9 1K
			⑦	P5 PIN 13
			⑦	P5 PIN 15
			⑥	C31 0.1

ORIGIN BLOCK #	CIRCUIT DESIGNATION	DESCRIPTION OR NAME OF SIGNAL	DESTINATION BLOCK #	CIRCUIT DESIGNATION
⑦	J7 PIN 1	PROBE INTENSITY SIGNAL	②	R139 10K
⑧	J3 PIN 1	OVER TEMP. SIGNAL	⑨	P3 PIN 1
⑧	J3 PIN 3	DUTY CYCLE	⑨	P3 PIN 3
⑧	J3 PIN 6	- 12 TO -70V D.C. VARIABLE	⑨	P3 PIN 6
⑩	IC 26 PIN 14	+10V REF. (H.V.)	⑦	J7 PIN 2
⑩	C 34 (-)	STIMULATOR COMMON (V)	②	J4 PIN 17
⑩	VR1 OUT	ULTRASOUND +12VDC	①	- ⑦
⑩	BR1 (-)	-70V UNREGULATED	⑧	- ⑨
⑩	BR1 (+)	ULTRASOUND COMMON	⑧	R14 330Ω
			⑦	- ⑧ - ⑨

NOTES:

1. ALL CAPACITORS IN μ F UNLESS OTHERWISE SPECIFIED
2. ALL RESISTORS IN OHMS, $\frac{1}{4}$ W, 5% UNLESS OTHERWISE SPECIFIED
3. ALL DIODES IN 4148 UNLESS OTHERWISE SPECIFIED

NTE 519



INTELECT 700 SCHEMATIC
SIGNAL LIST

DRAWING NO.

E.S. 72823

SIZE

B

SCALE NONE

SHEET 1

OF 11

UNLESS OTHERWISE SPECIFIED
TOLERANCES ARE:
FRACTIONS DECIMALS ANGLES
XX' XX" XXX'

APPROVALS
ENGINEER D. BLEY
DRAWING
CHECK


DATE
5-11-83
4-4-83

MATERIAL

APPROVAL
D. Bley

FINISH


DO NOT SCALE DRAWING



35

--- PART OF R121 INTENSITY CONTROL
(SEE BLOCK NO. 2 DWG.)


UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE FRACTIONS DECIMALS ANGLES XXX .XXX .XXX		APPROVALS	DATE
ENGINEER		D. BLEY	5-11-83
DRAWING		JSM	3-7-83
CHECK			
MATERIAL			
FINISH			
DO NOT SCALE DRAWING			
CHATTAHOOGA CORPORATION			
INTELECT 700 SCHEMATIC			
STIMULATOR FUNCTION / MODE LOGIC			
SIZE	DRAWING NO.		
B	1		
SCALE NONE	REV. G		
SHEET 2	OF 11		



INTELT 700 SCHEMATIC STIMULATOR INTENSITY CONTROL CIRCUITRY

SIZE	B	DRAWING NO.	E.S. 72823
SCALE	NONE	REV.	G
SHEET	3	OF	11


UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE FRACTIONS DECIMALS ANGLES	APPROVALS	DATE
1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42 43 44 45 46 47 48 49 50 51 52 53 54 55 56 57 58 59 60 61 62 63 64 65 66 67 68 69 70 71 72 73 74 75 76 77 78 79 80 81 82 83 84 85 86 87 88 89 90 91 92 93 94 95 96 97 98 99 100	ENGINEER DRAWING CHECK	5-11-83 3-8-83 5-7-83
MATERIAL		
FINISH		
DO NOT SCALE DRAWING		



APPROVALS		DATE
ENGINEER	<i>D. E. Y</i>	5-11-83
DRAWING	<i>8/8/83</i>	3-10-83
CHECK		
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE FRACTIONS DECIMALS ANGLES XXX .XXX .XXX		
MATERIAL		
FINISH		
DO NOT SCALE DRAWING		
DRAWING NO.		72823
SIZE	B	SCALE NONE
BLOCK	3	REV. G
SHEET	4	OF 11

micke


~500μs - 800μs



UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE FRACTIONS DECIMALS ANGLES 1/16 1/32 1/64 .01 .001 .0005 1/16 1/32 1/64 .01 .001 .0005 1/16 1/32 1/64 .01 .001 .0005		APPROVALS		DATE	
ENGINEER D. B. LEY		DRAWING J. M. M.		5-11-83	
CHECK		APPROVAL		5/1/83	
MATERIAL		FINISH		DO NOT SCALE DRAWING	
INTELECT 700 SCHEMATIC		STIMULATOR PULSE RATE + TIMING GENERATOR		DRAWING NO. E.S. 72823	
SCALE NONE		REV. G		SHEET 5 OF 11	


A Scope GND
B Test Point
C Test Point

500Ω Load Digi to Real (or B/L)
With MAX INT, Check TWIN METER (0-500V)
CONNECT Scope 1 To Common Point Between
C32 & C33 Check wave form & Voltage JLN



IF NOT FAULTY (IE 1 PULSE) CURRENT
NORMALLY CHECK & REPLACE D11 & 12
Q20 & 21 & C32 & 33

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES FRACTIONS DECIMALS ANGLES 1/16 1/32 1/64 1/8 1/4 1/2 3/4 1 1 1/2 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42 43 44 45 46 47 48 49 50 51 52 53 54 55 56 57 58 59 60 61 62 63 64 65 66 67 68 69 70 71 72 73 74 75 76 77 78 79 80 81 82 83 84 85 86 87 88 89 90 91 92 93 94 95 96 97 98 99 100 101 102 103 104 105 106 107 108 109 110 111 112 113 114 115 116 117 118 119 120 121 122 123 124 125 126 127 128 129 130 131 132 133 134 135 136 137 138 139 140 141 142 143 144 145 146 147 148 149 150 151 152 153 154 155 156 157 158 159 160 161 162 163 164 165 166 167 168 169 170 171 172 173 174 175 176 177 178 179 180 181 182 183 184 185 186 187 188 189 190 191 192 193 194 195 196 197 198 199 200 201 202 203 204 205 206 207 208 209 210 211 212 213 214 215 216 217 218 219 220 221 222 223 224 225 226 227 228 229 230 231 232 233 234 235 236 237 238 239 240 241 242 243 244 245 246 247 248 249 250 251 252 253 254 255 256 257 258 259 260 261 262 263 264 265 266 267 268 269 270 271 272 273 274 275 276 277 278 279 280 281 282 283 284 285 286 287 288 289 290 291 292 293 294 295 296 297 298 299 300 301 302 303 304 305 306 307 308 309 310 311 312 313 314 315 316 317 318 319 320 321 322 323 324 325 326 327 328 329 330 331 332 333 334 335 336 337 338 339 340 341 342 343 344 345 346 347 348 349 350 351 352 353 354 355 356 357 358 359 360 361 362 363 364 365 366 367 368 369 370 371 372 373 374 375 376 377 378 379 380 381 382 383 384 385 386 387 388 389 390 391 392 393 394 395 396 397 398 399 400 401 402 403 404 405 406 407 408 409 410 411 412 413 414 415 416 417 418 419 420 421 422 423 424 425 426 427 428 429 430 431 432 433 434 435 436 437 438 439 440 441 442 443 444 445 446 447 448 449 450 451 452 453 454 455 456 457 458 459 460 461 462 463 464 465 466 467 468 469 470 471 472 473 474 475 476 477 478 479 480 481 482 483 484 485 486 487 488 489 490 491 492 493 494 495 496 497 498 499 500 501 502 503 504 505 506 507 508 509 510 511 512 513 514 515 516 517 518 519 520 521 522 523 524 525 526 527 528 529 530 531 532 533 534 535 536 537 538 539 540 541 542 543 544 545 546 547 548 549 550 551 552 553 554 555 556 557 558 559 560 561 562 563 564 565 566 567 568 569 570 571 572 573 574 575 576 577 578 579 580 581 582 583 584 585 586 587 588 589 590 591 592 593 594 595 596 597 598 599 600 601 602 603 604 605 606 607 608 609 610 611 612 613 614 615 616 617 618 619 620 621 622 623 624 625 626 627 628 629 630 631 632 633 634 635 636 637 638 639 640 641 642 643 644 645 646 647 648 649 650 651 652 653 654 655 656 657 658 659 660 661 662 663 664 665 666 667 668 669 670 671 672 673 674 675 676 677 678 679 680 681 682 683 684 685 686 687 688 689 690 691 692 693 694 695 696 697 698 699 700 701 702 703 704 705 706 707 708 709 710 711 712 713 714 715 716 717 718 719 720 721 722 723 724 725 726 727 728 729 730 731 732 733 734 735 736 737 738 739 740 741 742 743 744 745 746 747 748 749 750 751 752 753 754 755 756 757 758 759 760 761 762 763 764 765 766 767 768 769 770 771 772 773 774 775 776 777 778 779 780 781 782 783 784 785 786 787 788 789 790 791 792 793 794 795 796 797 798 799 800 801 802 803 804 805 806 807 808 809 810 811 812 813 814 815 816 817 818 819 820 821 822 823 824 825 826 827 828 829 830 831 832 833 834 835 836 837 838 839 840 841 842 843 844 845 846 847 848 849 850 851 852 853 854 855 856 857 858 859 860 861 862 863 864 865 866 867 868 869 870 871 872 873 874 875 876 877 878 879 880 881 882 883 884 885 886 887 888 889 890 891 892 893 894 895 896 897 898 899 900 901 902 903 904 905 906 907 908 909 910 911 912 913 914 915 916 917 918 919 920 921 922 923 924 925 926 927 928 929 930 931 932 933 934 935 936 937 938 939 940 941 942 943 944 945 946 947 948 949 950 951 952 953 954 955 956 957 958 959 960 961 962 963 964 965 966 967 968 969 970 971 972 973 974 975 976 977 978 979 980 981 982 983 984 985 986 987 988 989 990 991 992 993 994 995 996 997 998 999 1000 1001 1002 1003 1004 1005 1006 1007 1008 1009 1010 1011 1012 1013 1014 1015 1016 1017 1018 1019 1020 1021 1022 1023 1024 1025 1026 1027 1028 1029 1030 1031 1032 1033 1034 1035 1036 1037 1038 1039 1040 1041 1042 1043 1044 1045 1046 1047 1048 1049 1050 1051 1052 1053 1054 1055 1056 1057 1058 1059 1060 1061 1062 1063 1064 1065 1066 1067 1068 1069 1070 1071 1072 1073 1074 1075 1076 1077 1078 1079 1080 1081 1082 1083 1084 1085 1086 1087 1088 1089 1090 1091 1092 1093 1094 1095 1096 1097 1098 1099 1100 1101 1102 1103 1104 1105 1106 1107 1108 1109 1110 1111 1112 1113 1114 1115 1116 1117 1118 1119 1120 1121 1122 1123 1124 1125 1126 1127 1128 1129 1130 1131 1132 1133 1134 1135 1136 1137 1138 1139 1140 1141 1142 1143 1144 1145 1146 1147 1148 1149 1150 1151 1152 1153 1154 1155 1156 1157 1158 1159 1160 1161 1162 1163 1164 1165 1166 1167 1168 1169 1170 1171 1172 1173 1174 1175 1176 1177 1178 1179 1180 1181 1182 1183 1184 1185 1186 1187 1188 1189 1190 1191 1192 1193 1194 1195 1196 1197 1198 1199 1200 1201 1202 1203 1204 1205 1206 1207 1208 1209 1210 1211 1212 1213 1214 1215 1216 1217 1218 1219 1220 1221 1222 1223 1224 1225 1226 1227 1228 1229 1230 1231 1232 1233 1234 1235 1236 1237 1238 1239 1240 1241 1242 1243 1244 1245 1246 1247 1248 1249 1250 1251 1252 1253 1254 1255 1256 1257 1258 1259 1260 1261 1262 1263 1264 1265 1266 1267 1268 1269 1270 1271 1272 1273 1274 1275 1276 1277 1278 1279 1280 1281 1282 1283 1284 1285 1286 1287 1288 1289 1290 1291 1292 1293 1294 1295 1296 1297 1298 1299 1300 1301 1302 1303 1304 1305 1306 1307 1308 1309 1310 1311 1312 1313 1314 1315 1316 1317 1318 1319 1320 1321 1322 1323 1324 1325 1326 1327 1328 1329 1330 1331 1332 1333 1334 1335 1336 1337 1338 1339 1340 1341 1342 1343 1344 1345 1346 1347 1348 1349 1350 1351 1352 1353 1354 1355 1356 1357 1358 1359 1360 1361 1362 1363 1364 1365 1366 1367 1368 1369 1370 1371 1372 1373 1374 1375 1376 1377 1378 1379 1380 1381 1382 1383 1384 1385 1386 1387 1388 1389 1390 1391 1392 1393 1394 1395 1396 1397 1398 1399 1400 1401 1402 1403 1404 1405 1406 1407 1408 1409 1410 1411 1412 1413 1414 1415 1416 1417 1418 1419 1420 1421 1422 1423 1424 1425 1426 1427 1428 1429 1430 1431 1432 1433 1434 1435 1436 1437 1438 1439 1440 1441 1442 1443 1444 1445 1446 1447 1448 1449 1450 1451 1452 1453 1454 1455 1456 1457 1458 1459 1460 1461 1462 1463 1464 1465 1466 1467 1468 1469 1470 1471 1472 1473 1474 1475 1476 1477 1478 1479 1480 1481 1482 1483 1484 1485 1486 1487 1488 1489 1490 1491 1492 1493 1494 1495 1496 1497 1498 1499 1500 1501 1502 1503 1504 1505 1506 1507 1508 1509 1510 1511 1512 1513 1514 1515 1516 1517 1518 1519 1520 1521 1522 1523 1524 1525 1526 1527 1528 1529 1530 1531 1532 1533 1534 1535 1536 1537 1538 1539 1540 1541 1542 1543 1544 1545 1546 1547 1548 1549 1550 1551 1552 1553 1554 1555 1556 1557 1558 1559 1560 1561 1562 1563 1564 1565 1566 1567 1568 1569 1570 1571 1572 1573 1574 1575 1576 1577 1578 1579 1580 1581 1582 1583 1584 1585 1586 1587 1588 1589 1590 1591 1592 1593 1594 1595 1596 1597 1598 1599 1600 1601 1602 1603 1604 1605 1606 1607 1608 1609 1610 1611 1612 1613 1614 1615 1616 1617 1618 1619 1620 1621 1622 1623 1624 1625 1626 1627 1628 1629 1630 1631 1632 1633 1634 1635 1636 1637 1638 1639 1640 1641 1642 1643 1644 1645 1646 1647 1648 1649 1650 1651 1652 1653 1654 1655 1656 1657 1658 1659 1660 1661 1662 1663 1664 1665 1666 1667 1668 1669 1670 1671 1672 1673 1674 1675 1676 1677 1678 1679 1680 1681 1682 1683 1684 1685 1686 1687 1688 1689 1690 1691 1692 1693 1694 1695 1696 1697 1698 1699 1700 1701 1702 1703 1704 1705 1706 1707 1708 1709 1710 1711 1712 1713 1714 1715 1716 1717 1718 1719 1720 1721 1722 1723 1724 1725 1726 1727 1728 1729 1730 1731 1732 1733 1734 1735 1736 1737 1738 1739 1740 1741 1742 1743 1744 1745 1746 1747 1748 1749 1750 1751 1752 1753 1754 1755 1756 1757 1758 1759 1760 1761 1762 1763 1764 1765 1766 1767 1768 1769 1770 1771 1772 1773 1774 1775 1776 1777 1778 1779 1780 1781 1782 1783 1784 1785 1786 1787 1788 1789 1790 1791 1792 1793 1794 1795 1796 1797 1798 1799 1800 1801 1802 1803 1804 1805 1806 1807 1808 1809 1810 1811 1812 1813 1814 1815 1816 1817 1818 1819 1820 1821 1822 1823 1824 1825 1826 1827 1828 1829 1830 1831 1832 1833 1834 1835 1836 1837 1838 1839 1840 1841 1842 1843 1844 1845 1846 1847 1848 1849 1850 1851 1852 1853 1854 1855 1856 1857 1858 1859 1860 1861 1862 1863 1864 1865 1866 1867 1868 1869 1870 1871 1872 1873 1874 1875 1876 1877 1878 1879 1880 1881 1882 1883 1884 1885 1886 1887 1888 1889 1890 1891 1892 1893 1894 1895 1896 1897 1898 1899 1900 1901 1902 1903 1904 1905 1906 1907 1908 1909 1910 1911 1912 1913 1914 1915 1916 1917 1918 1919 1920 1921 1922 1923 1924 1925 1926 1927 1928 1929 1930 1931 1932 1933 1934 1935 1936 1937 1938 1939 1940 1941 1942 1943 1944 1945 1946 1947 1948 1949 1950 1951 1952 1953 1954 1955 1956 1957 1958 1959 1960 1961 1962 1963 1964 1965 1966 1967 1968 1969 1970 1971 1972 1973 1974 1975 1976 1977 1978 1979 1980 1981 1982 1983 1984 1985 1986 1987 1988 1989 1990 1991 1992 1993 1994 1995 1996 1997 1998 1999 2000 2001 2002 2003 2004 2005 2006 2007 2008 2009 2010 2011 2012 2013 2014 2015 2016 2017 2018 2019 2020 2021 2022 2023 2024 2025 2026 2027 2028 2029 2030 2031 2032 2033 2034 2035 2036 2037 2038 2039 2040 2041 2042 2043 2044 2045 2046 2047 2048 2049 2050 2051 2052 2053 2054 2055 2056 2057 2058 2059 2060 2061 2062 2063 2064 2065 2066 2067 2068 2069 2070 2071 2072 2073 2074 2075 2076 2077 2078 2079 2080 2081 2082 2083 2084 2085 2086 2087 2088 2089 2090 2091 2092 2093 2094 2095 2096 2097 2098 2099 2100 2101 2102 2103 2104 2105 2106 2107 2108 2109 2110 2111 2112 2113 2114 2115 2116 2117 2118 2119 2120 2121 2122 2123 2124 2125 2126 2127 2128 2129 2130 2131 2132 2133 2134 2135 2136 2137 2138 2139 2140 2141 2142 2143 2144 2145 2146 2147 2148 2149 2150 2151 2152 2153 2154 2155 2156 2157 2158 2159 2160 2161 2162 2163 2164 2165 2166 2167 2168 2169 2170 2171 2172 2173 2174 2175 2176 2177 2178 2179 2180 2181 2182 2183 2184 2185 2186 2187 2188 2189 2190 2191 2192 2193 2194 2195 2196 2197 2198 2199 2200 2201 2202 2203 2204 2205 2206 2207 2208 2209 2210 2211 2212 2213 2214 2215 2216 2217 2218 2219 2220 2221 2222 2223 2224 2225 2226 2227 2228 2229 2230 2231 2232 2233 2234 2235 2236 2237 2238 2239 2240 2241 2242 2243 2244 2245 2246 2247 2248 2249 2250 2251 2252 2253 2254 2255 2256 2257 2258 2259 2260 2261 2262 2263 2264 2265 2266 2267 2268 2269 2270 2271 2272 2273 2274 2275 2276 2277 2278 2279 2280 2281 2282 2283 2284 2285 2286 2287 2288 2289 2290 2291 2292 2293 2294 2295 2296 2297 2298 2299 2300 2301 2302 2303 2304 2305 2306 2307 2308 2309 2310 2311 2312 2313 2314 2315 2316 2317 2318 2319 2320 2321 2322 2323 2324 2325 2326 2327 2328 2329 2330 2331 2332 2333 2334 2335 2336 2337 2338 2339 2340 2341 2342 2343 2344 2345 2346 2347 2348 2349 2350 2351 2352 2353 2354 2355 2356 2357 2358 2359 2360 2361 2362 2363 2364 2365 2366 2367 2368 2369 2370 2371 2372 2373 2374 2375 2376 2377 2378 2379 2380 2381 2382 2383 2384 2385 2386 2387 2388 2389 2390 2391 2392 2393 2394 2395 2396 2397 2398 2399 2400 2401 2402 2403 2404 2405 2406 2407 2408 2409 2410 2411 2412 2413 2414 2415 2416 2417 2418 2419 2420 2421 2422 2423 2424 2425 2426 2427 2428 2429 2430 2431 2432 2433 2434 2435 2436 2437 2438 2439 2440 2441 2442 2443 2444 2445 2446 2447 2448 2449 2450 2451 2452 2453 2454 2455 2456 2457 2458 2459 2460 2461 2462 2463 2464 2465 2466 2467 2468 2469 2470 2471 2472 2473 2474 2475 2476 2477 2478 2479 2480 2481 2482 2483 2484 2485 2486 2487 2488 2489 2490 2491 2492 2493 2494 2495 2496 2497 2498 2499 2500 2501 2502 2503 2504 2505 2506 2507 2508 2509 2510 2511 2512 2513 2514 2515 2516 2517 2518 2519 2520 2521 2522 2523 2524 2525 2526 2527 2528 2529 2530 2531 2532 2533 2534 2535 2536 2537 2538 2539 2540 2541 2542 2543 2544 2545 2546 2547 2548 2549 2550 2551 2552 2553 2554 2555 2556 2557 2558 2559 2560 2561 2562 2563 2564 2565 2566 2567 2568 2569 2570 2
--




INTELECT 700 SCHEMATIC HIGH VOLTAGE OUTPUT SELECTION

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE FRACTIONS DECIMALS ANGLES XX: ; XXX: ;		APPROVALS		DATE	
MATERIAL	FINISH	ENGINEER D. BLEY		5-11-83	
		DRAWING J & M		3-24-83	
		CHECK			
DO NOT SCALE DRAWING		APPROVAL 5-17-83			

SIZE	DRAWING NO.
B	7
SCALE	REV. 3
NONE	SHEET 8 OF 11

E.S. 72823

NOT checked
14 2-15 210000



UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES FRACTIONS DECIMALS ANGLES XXX XXX XXX		APPROVALS	DATE
ENGINEER <i>D. BLEY</i>		DATE	5-17-83
DRAWING		CHECK	4-13-83
APPROVAL <i>[Signature]</i>		DATE	5-17-83
MATERIAL		FINISH	
DO NOT SCALE DRAWING			


CHATTANOOGA CORPORATION	
INTELECT 700 SCHEMATIC	
ULTRASOUND CONTROL CIRCUIT	
SIZE B	DRAWING NO. E.S. 72823
SCALE NONE	REV. G
SHEET 9 OF 11	

Jumping H.V. Low End
Check at 0.2 R15

NOTE: IF
Q1 2N6027
UPDATE MUST
BE DONE
R32 TO 10K
R32 ADD 6.8K
Twice Legs
ADD R3C 330K

905,000 + 1,050,000

Note: IF 330~ BUINS UP
Check crystal contact!! Q5 2N6545



* SELECTED AT CALIBRATION



INTELECT 700 SCHEMATIC
ULTRASOUND OSCILLATOR

DRAWING NO.

SIZE B

SCALE NONE

REV. G

SHEET 10 OF 11

Block 9

E.S.

72823

UNLESS OTHERWISE SPECIFIED
DIMENSIONS ARE IN INCHES
TOLERANCES ARE
FRACTIONS DECIMALS ANGLES
XXXI

MATERIAL

FINISH

DO NOT SCALE DRAWING

APPROVALS

DATE

ENGINEER

DRAWING

CHECK

DATE

ENGINEER

DRAWING

CHECK

DATE

ENGINEER

DRAWING

CHECK

DATE

ENGINEER

DRAWING

CHECK

DATE

ENGINEER

rear panel designations

infect™ MODEL 700 ULTRASOUND GENERATOR S/N

OSC FREQ 1.0 MHz, CONTINUOUS AND PULSED @ 100 HZ

DUTY CYCLE	TEMPORAL PEAK/AVG INTENSITY RATIO
20%	100 PPS 5
50%	100 PPS 2
100%	CONT 1

THIS DEVICE COMPLIES WITH REQUIREMENTS SET FORTH IN 21CFR 1050.10
FCC TYPE APPROVED

CHATTANOOGA CORPORATION CHATTANOOGA, TN. 37405 70854

GROUNDING RELIABILITY CAN ONLY BE ACHIEVED WHEN THIS EQUIPMENT IS CONNECTED TO AN EQUIVALENT RECEPTACLE MARKED "HOSPITAL GRADE"

DANGER HIGH VOLTAGE REFER SERVICING ONLY TO QUALIFIED PERSONNEL

FCC ID: BWU8UKU700
UNITED STATES OF AMERICA

CHATTANOOGA CORPORATION CHATTANOOGA, TN. 37405

POWER INPUT 207 253V 50/60 HZ 1.25 AMP MAX

WARNING FOR CONTINUED PROTECTION AGAINST FIRE HAZARD REPLACE ONLY WITH SAME TYPE AND RATING OF FUSE

FUSE 1 AMP/250 V. SLO-BLO

CAUTION: CASUALTY (BURN) AND FIRE HAZARD DO NOT USE NEAR CONDUCTIVE MATERIALS SUCH AS METAL BED POSTS, INNER SPRING MATTRESSES, ETC. RENEW ELECTRODE CABLES UPON EVIDENCE OF DETERIORATION

CAUTION: ELECTRIC SHOCK HAZARD DO NOT REMOVE COVER REFER SERVICING TO QUALIFIED SERVICE PERSONNEL


CERTIFIE SELON LES EXIGENCES DU CODE CANADIEN DE L'ÉLECTRICITÉ. L'ACNOR N'A PAS ÉTUDIÉ LES AUTRES EFFETS PHYSIOLOGIQUES POSSIBLES.

CERTIFIED TO THE REQUIREMENTS OF THE CANADIAN ELECTRICAL CODE CSA HAS NOT INVESTIGATED OTHER PHYSIOLOGICAL EFFECTS.

LR-16036

70424

This is the actual serial number of your Intelect 700. ^{T.M.}



CHATTANOOGA CORPORATION
Chattanooga, Tennessee 37405