VWR 12 and 18 litre shaking water bath

INSTRUCTION MANUAL

European Catalogue Numbers

462-0353 18 Litre with UK mains plug

462-0354 12 Litre with UK mains plug

462-0355 18 Litre with EURO mains plug

462-0356 12 Litre with EURO mains plug

Version: 1
Issued: 3/1/2013

EN

FR

DE

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Instruction manual

Introduction

The VWR shaking water baths are a general purpose series of thermostatically controlled shaking bath designed for indoor laboratory use by a professional user. Conforms to DIN 12876 Class 1 for use with non-flammable liquid.

Contents

1 Safety

- 1.1 Warning symbol
- 1.2 Safety certification
- 1.3 Safety features
- 1.4 Before first operating the equipment
- 1.5 Precautions during and after operation

2. Getting started

- 2.1 Unpacking
- 2.2 Optional accessories
- 2.3 Electrical supply
- 2.4 Conditions of use

3. Operation

- 3.1 Water level
- 3.2 Flat bottomed vessels
- 3.3 Operation above 60°C
- 3.4 Fitting the shaking trolley and tray
- 3.5 Setting the temperature and shaking speed of baths

4. Specification

- 4.1 Physical and electrical details
- 4.2 Performance

5. Maintenance

- 5.1 Cleaning
- 5.2 Replacement of fuses
- 5.3 Routine safety tests
- 5.4 Service

6. Technical Tips

- 6.1 Which water should you use in your bath?
- 6.2 How to prevent rust in water baths
- 6.3 How to prevent algae and bacteria

7. Troubleshooting

8. Warranty

9. Disposal

1. Safety

1.1 Warning symbols

The symbols below are marked on the equipment to indicate:

Caution: Surfaces and water can be hot during and after use.

Read this manual before using the bath

1.2 Safety certification

VWR water baths meet the requirements of international safety standard IEC 61010-2-10 and IEC 61010-2-51 and national standards based on them including:

UL 61010A-2-010

UL 61010-2-051-04

CAN/CSA-C22.2 NO. 61010-2-010-04;

CAN/CSA-C22.2 No. 61010-2-051

1.3 Safety features

The shaking trolley is driven by a magnetic drive system. The drive magnet in the tank is magnetically coupled to a motor-driven magnet under the tank. The magnetic coupling cannot exert enough force to cause injury.

To protect samples in case of failure of the primary temperature control system there is a second completely independent user-settable control system.

In addition to the second control system for sample protection, there is a separate fixed-temperature cut-out to protect the bath and surroundings in the very unlikely event of failure of both control systems.

1.4 Before first operating the equipment

Read the whole of these instructions. Safety may be impaired if they are not followed.

If the equipment has been transported or stored in cold or humid conditions, condensation may form inside it. If that could have happened, allow time (at least 2 hours) for the condensation to evaporate before using the equipment.

1.5 Precautions during and after operation

The bath is for use only with water as the bath liquid. Make sure that it cannot become contaminated by other liquids. The use of de-ionized water is recommended to prevent corrosion. A full explanation of this is given in Section 6 of this manual.

CAUTION. PLEASE FILL BATH WITH WATER BEFORE SWITCHING ON.

Change the water regularly and empty when not in use for prolonged periods, to further prevent corrosion and contamination.

Before emptying a bath, allow the water temperature to fall to a safe level. The 18 litre bath, should be emptied moving.

Do not use the equipment in an area where there are aggressive or explosive chemical mixtures or atmosphere.

If potentially hazardous liquid is spilt onto the equipment, disconnect it from the power supply and have it checked by a qualified service technician. It is the user's responsibility to carry out appropriate decontamination if hazardous material is spilled on the equipment.

Do not use the bath to heat any material that could cause a fire or any other kind of hazard.

2. Getting started

2.1 Unpacking

Remove packing materials carefully, and retain for future shipment or storage of the equipment. Standard equipment includes:

- Thermostatic shaking bath
- Mains cord with plug
- Drive magnet
- Shaking trolley
- Universal tray
- Gabled polycarbonate lid
- Operating manual
- Quick start guide

CAUTION. PLEASE FILL BATH WITH WATER BEFORE SWITCHING ON.

2.2 Optional accessories

	VWR European Catalogue number
Universal flask tray - with springs	
For 12 litre baths 462-0354, 462-0356	462-0360
For 18 litre baths 462-0353, 462-0355	462-0359
Plain tray	
For 12 litre baths 462-0354, 462-0356	462-0358
For 18 litre baths 462-0353, 462-0355	462-0357
Test tube tray – to hold test tube racks	
For 12 litre bath, holds 3 test tube racks	462-0361
For 18 litre bath, holds 5 test tube racks	462-0362
Test tube racks	
Rack for 10mm tubes	462-0363
Rack for 13mm tubes	462-0364
Rack for 16mm tubes	462-0365
Rack for 19mm tubes	462-0366
Rack for 25mm tubes	462-0367
Rack for 30mm tubes	462-0368
Rack for 0.5mm micro-centrifuge tubes	462-0369
Rack for 1.5mm micro-centrifuge tubes	462-0370
Replacement non drip Polycarbonate Gabled Lids	
For 12 litre baths 462-0354, 462-0356	462-0373
For 18 litre baths 462-0353, 462-0355	462-0374
Stainless Steel Gabled Lids	
For 12 litre baths 462-0354, 462-0356	462-0371
For 18 litre baths 462-0353, 462-0355	462-0372

2.3 Electrical supply

Check that the supply voltage marked on the serial number plate, and the type of mains plug, are correct for your mains supply outlet (which must have a ground connector).

To disconnect the equipment from the mains supply, remove the mains plug from the mains supply outlet. Make sure that the mains plug is easily accessible.

2.4 Conditions of use

The water baths are for indoor laboratory use only. Check that the environmental conditions of the laboratory are within the following limits:

Temperature	5 to 40 °C
Maximum relative humidity	80 % r.h. in room temperatures up to 31°C decreasing linearly to 50 % r.h. at 40°C
Altitude	Up to 2,000m above sea level

3. Operation

3.1 Water level

Ensure that the water level is at least 5cm above the bottom of the tank, and not higher than 2.5 cm from the top. This applies both without any vessels in the bath and with the maximum contents.

If the bath is allowed to run dry, the safety cutout will trip and disconnect the heater. If this happens, unplug the bath and have the cut-out re-set by a qualified service technician.

Water is the only liquid suitable for use in the baths.

3.2 Operation above 60°C

The supplied gabled lid must always be used above 60°C to optimize the temperature control; for the bath to operate efficiently and to ensure the set temperature is reached and maintained.


Without the lid, excessive evaporation will require the bath to be filled more often.

3.3 Use as a thermostatic bath without shaking

If the shaking trolley and drive magnet (see 3.4 below) are lifted out, a shaking bath can be used as a simple thermostatic bath. When used in this way, if flat bottomed vessels or objects are to be placed in the bath, always use a raised base tray (optional accessory) to avoid damage to the under-tank heater.

3.4 Fitting the shaking trolley and tray

1. Put the drive magnet into the hole in the base of the tank
2. Place the shaking trolley in the tank with the axles parallel to the front panel, and with the peg of the drive arm in the hole in the top of the drive magnet
3. Place the universal tray on the shaking trolley

3.5 Setting the temperature and shaking speed

3.5.1 Setting the control temperature


1. Turn the knob of the sample protection thermostat to maximum
2. Press  control key, the display will show °C, press  control key again. The display will flash once a second with the current set temperature to show that the controller is ready to be re-set
3. Press **up**  or **down**  control keys to set the desired temperature
(If no key is pressed for 4 seconds, the display will revert back to showing the actual temperature and the set point will remain at its original value)
4. Press  control key. This will store the requested value and the display will revert to showing the actual water temperature. The water temperature will change to the new set value
5. When the temperature has stabilized, measure the temperature with a thermometer. If it is different from the set point temperature use the single point calibration method (section 3.5.3) to adjust the display temperature

EN

3.5.2 Setting the sample protection thermostat (if required)

1. Turn the knob of the sample protection thermostat to maximum
2. Set the control temperature 2°C above the desired operating temperature, and wait for the temperature to stabilize
3. Turn the knob of the sample protection thermostat slowly anti-clockwise until a click is heard and **0.0** is displayed
The audible alarm will sound continuously to warn you that the bath is controlling using the sample protection thermostat.
4. Re-set the control temperature to the desired temperature
The audible alarm will sound continuously until the primary control system is re-activated.

In the unlikely event of failure of the primary temperature control system, the sample protection thermostat will maintain the water in the bath at a temperature a few degrees above the set temperature and with greater fluctuations. The display will cycle between **0.0** and the actual bath temperature. If this happens, first check that the sample protection thermostat is not set to a lower temperature than the primary control system. If that is **not** the reason, the primary control system is not functioning correctly. The bath can continue to be used without compromising the safety of persons or the surroundings until current work is completed. However, it is recommended to have the bath checked by a qualified service technician as soon as conveniently possible.


Sample protection thermostat

Control keys

3.5.3 Setting the shaking speed

1. Press control key, the display will flash °C
2. Press **up** control key until the display changes to a moving U symbol:

3. Press control key. The display will show the present shaking speed
4. Press **up** or **down** control keys to select the required speed (off, or between 40 and 400)
5. Press control key to set the speed to the value displayed
(If no key is pressed for 4 seconds, the speed will revert to the previous value)

The maximum speed for a particular application will depend on number and size of vessels, and the depth of the water. The shaking trolley is driven by a drive magnet under the tank, magnetically coupled to the trolley through the tank. As the speed is increased a point is reached where the magnetic coupling de-latches, causing the trolley either to move erratically or to stop.

The shaking speed will revert to previously set value when reconnected after power interruption.

3.5.4 Single point calibration (SPC)

The display of set temperature and actual temperature is accurate to within 1°C at 37°C but may be up to 5° C out at maximum temperature. Using SPC the display can be re-set to be accurate to $\pm 0.1^\circ \text{C}$ at a single point by the following procedure:

1. Set the bath to the calibration temperature and allow the bath temperature to stabilize for 20 minutes and indicated by the display not fluctuating by more than 0.2°C. Measure the actual water temperature with a calibrated thermometer
2. Press  control key. The display will show °C, press the **up**  control key until the display changes to show **SPC**, press  control key to enter SPC mode. The display will flash once a second with the current value and shows that the controller is ready to be re-set
3. Press **up**  or **down**  control keys to set the amount of correction value (positive or negative) as measured by a thermometer. (If no key is pressed for 4 seconds, the display will revert back to showing the actual temperature and the set point will remain at its original value)
4. Press  control key. This will store the new corrected set temperature and the display will revert to showing the new correct actual water temperature

4 Specification

4.1 Physical and electrical details

Mains supply 230V @ 50/60 Hz
Pollution degree: 2
Installation Category: II

Note: Mains supply voltage fluctuations are not to exceed $\pm 10\%$ of the nominal supply voltage

Models	Capacity/Litres	Current Rating/Amps
462-0354 462-0356	12L	2.8A
462-0353 462-0355	18L	5.4A

4.2 Performance

Range	5°C above ambient to 99°C
Display (also used for setting)	10.0 to 99.0°C in steps of 0.1°C
Temperature stability to DIN 12876-3 *	$\pm 0.2^\circ \text{C}$
Shaking speed range	40 to 400 strokes per minute

The maximum shaking speed for a particular application will depend on the number and size of vessels in the shaking tray, and the depth of the water).

* The supplied polycarbonate lid must be used to achieve the quoted performance

5 Maintenance and service

No routine maintenance is required except for cleaning.

5.1 Cleaning

Clean the equipment with a damp cloth, using water only. Do not use chemical cleaning agents. Before using any other cleaning or decontamination method, check with the manufacturer or supplier to make sure that the proposed method will not damage the equipment.

Scale on immersed parts can be removed using chemical de-scaling products designed for use on kitchen equipment which has metal heating elements.

Warning - De-scaling products may be toxic - follow the manufacturer's instructions

5.2 Replacement of fuses

Fuses are Littlefuse 3AB 314 series. Fast-acting, high breaking current (max breaking current at least 750 A). Dimensions; 1.25 inch long, 0.25 inch diameter. Replace fuses only by the same type and rating (250volt).

Models	Capacity	Fuse rating/ Amps
462-0354 462-0356	12L	10
462-0353 462-0355	18L	15

Replace fuses as follows:

1. Disconnect the unit from the power supply
2. Remove the mains input connector from the socket at the back of the bath
3. Press down the fuse drawer catch
4. Pull out the fuse drawer, check the fuse(s) and replace it if necessary, using the fuse type and rating specified above
5. Push back the drawer, and replace the mains input connector

5.3 Routine safety tests

If routine tests are to be made, we recommend a test of the integrity of the protective earth conductor and an insulation test at 500 Vd.c. Routine flash tests are **not** recommended for any electrical equipment, because repeated high voltage tests degrade insulation materials.

5.4 Service

Web Resources

Visit the VWR's website at www.vwr.com for:

- Complete technical service contact information
- Access to VWR's Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or visit. www.vwr.com.

6 Technical Tips

6.1 Which water should you use in your bath?

For the long-term reliability of water baths it is important to use oxygenated water that is free from ions and minerals that can cause corrosion of stainless steel. We recommend the use of distilled water and de-ionized water from modern ion exchange systems that do not use salt back flushing to regenerate the ion-exchange cartridges.

Stainless steel is protected from corrosion by a layer of chromium oxide. If the layer is damaged, oxygen present in water can reform the oxide layer. If the water is still or de-oxygenated, and the oxide layer is damaged, ions can corrode the stainless steel tank. If a water bath has been unused for some time, or water boiled, we recommend changing to fresh distilled water or correct de-ionized water.

Water normally contains calcium or magnesium ions. De-ionized water has most ions removed as indicated by its conductivity level; the purer the water the lower the conductivity. It is important to use only de-ionized water from an ion exchange system with replaceable cartridges. Do not use de-ionized water generated from an ion-exchange system that incorporates a salt back-flush system to regenerate the ion-exchange resin as this can leave sodium ions that are very corrosive to stainless steel.

6.2 How to prevent rust in water baths

Most VWR tanks, as well as immersed parts, are made from type 304 stainless steel, an extremely versatile general purpose grade of stainless steel. It is the excellent forming characteristic that has made this grade dominant in the manufacture of laboratory and industrial water baths, as well as domestic sinks and saucepans. Type 304 stainless steel is highly suitable for applications where hygiene is important; it exhibits good heat resistance and excellent resistance to corrosion.

However, despite resistance to general surface corrosion, stainless steel is susceptible to specific types of corrosion, in particular pitting (small pin hole style corrosion) and stress corrosion cracking. It can also undergo general corrosion in specific environments, such as one containing hydrochloric or sulphuric acids.

Stainless steel is protected by its high content of alloying elements, primarily chromium and nickel. Chromium is the most important with respect to corrosion resistance, although the nickel assists in allowing the chromium to do its job. The chromium forms an oxide layer on the surface of the steel, which inhibits further oxidation. This layer adheres extremely well to the metal substrate, but it is essential that it remains intact, and must be protected from various forms of damage.

If the surface chromium oxide layer becomes damaged, oxygen present in water can partially reform the oxide layer, so it is advisable to ensure that water is always fresh and well oxygenated. Baths that will be out of use for an extended period should be emptied, and all moisture should be wiped from the bottom of the tank.

In some cases a brown layer may appear on the surface of a stainless steel tank. In most of these cases this is not rust, but it may be a surface deposit of minerals from the local water supply, or ferrous particles or salts that have fallen into the tank. These surface deposits can usually be removed by using a household cleaner or metal polish.

6.3 How to prevent algae and bacteria

Water baths provide the ideal environment for the growth of micro-organisms. If left uncontrolled the growth of these organisms can result in a range of serious problems and health risks from pathogenic bacteria.

The growth of algae on the surface of parts will cause biofouling which can reduce performance.

Micro-organisms that produce acidic metabolic by-products can cause bio-corrosion by depolarization of metal surfaces.

There are a number of biocides available on the market.

7. Troubleshooting

No lights or display

Check power source, mains switch and fuses.

Water is not reaching temperature

Check that the set temperature and the sample protection thermostat are set correctly. For high temperatures (>90°C), the calibration procedure 3.5.4 may need to be used to reach the temperature.

Display shows 'OEL' or 'Err'

1. Check that the sample protection thermostat is set correctly. See section 3.5.2
2. Check the bath has water in it. If the it has been run dry the bath has tripped the fixed temperature cut-out, an internal circuit breaker. it will need to be checked by a qualified service technician. Contact your local VWR representative.

Display shows 'Prb'

Internal probe error. Contact your local VWR representative.

Display shows 'B.B.B.' flashing permanently.

System error. Contact your local VWR representative.

Display temperature is different from actual temperature

Use single point calibration (SPE) in section 3.5.4 to correct this.

There is no last digit on the display

Allow the bath to stabilize.

Temperature of water fluctuates

Check that the sample protection thermostat is set correctly.

8. Warranty

VWR International warrants that this product will be free from defects in material and workmanship for a period of three (3) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

9. Disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste. Instead it's your responsibility to correctly dispose of your equipment at lifecycle -end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment. For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

Thank you

9. Compliance

RoHS directive

All the products covered by this manual comply with the requirements of the RoHS Directive (Directive 2002/95/EC).

Electrical safety and electromagnetic compatibility

All the products covered by this manual comply with the requirements of the following European directives.

Low Voltage Directive (2006/95/EC) for electrical safety. Applying standard BS EN 61010 Part 2-010 and BS EN 61010 Part 2-051

EMC directive (2004/108/EC) for electromagnetic compatibility. Applying standard BS EN 61326-1

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VWR International AB
VWR International GmbH
VWR International Ltd