
Alaris® PK Syringe Pump

Sufentanil-Gepts
Model

Contents

	Page
EN - Directions for Use.....	1
DA - Brugsanvisning.....	4
DE - Gebrauchsanweisung	7
ES - Manual del Usuario.....	10
FR - Mode d'emploi	13
IT - Istruzioni per l'uso.....	16
NL - Gebruiksaanwijzing	19
NO - Bruksanvisning.....	22
PT - Instruções de Utilização.....	25
SV - Bruksanvisning	28

Introduction

This addendum should be used in conjunction with the Alaris® PK Syringe Pump Directions For Use. The purpose of this addendum is to provide additional material for operation when using a drug configuration based upon the Sufentanil-Gepts model. To use this drug / model combination a suitable data set must be created, see Alaris® PK Editor Software Directions For Use.

TCI Precautions

When first starting the infusion the pharmacokinetic / pharmacodynamic models within the Alaris® PK Syringe Pump are reset to zero. Therefore, if for any reason the pump is switched off during the surgical procedure, all current pharmacokinetic / pharmacodynamic model information will be lost. Under such circumstances switching the pump on and restarting the infusion whilst the patient contains a significant residual drug dose could result in an over-infusion and, therefore, the pump should not be restarted in TCI mode.

Pharmacokinetic and Pharmacodynamic parameters

Drug : Sufentanil Model: Gepts (not weight adjusted)

Age Limit: 12 years upwards

Unit of Plasma Concentration: ng/ml

Max. Plasma concentration: 2 ng/ml

$V_c = 14.3 \text{ l}$

$k_{10} = 0.0645 \text{ min}^{-1}$

$k_{12} = 0.1086 \text{ min}^{-1}$

$k_{13} = 0.0229 \text{ min}^{-1}$

$k_{21} = 0.0245 \text{ min}^{-1}$

$k_{31} = 0.0013 \text{ min}^{-1}$

Reference from the literature : Gepts et al.: Anesthesiology 1995, 83, 1194-1204

Additional: ke_0 calculated with time to peak effect 5.6 min ($ke_0 = 0.17559 \text{ min}^{-1}$)(reference: Shafer et al Anesthesiology. 1991 Jan;74(1):53-63)

Model Restrictions

An adult three compartmental PK model without any covariate, enlarged with an effect site model for Sufentanil is developed by Gepts et al. The model is developed and clinically validated through publications in a normal adult population. The PK model is also found to be accurate in obese patients. (Barvais et al. J Cardiothorac Vasc Anaest 2000 ; 14 : 402-8 ; Slepchenko et al. Anesthesiology 2003 ; 98:65-73; Pandin et al. Anesthesiology 2000 ; 93: 653-61).

Recommended infusate concentrations from publication are from 1 and 5 µg/ml.

Concerning target concentrations, originally, Bailey et al (Anesth Analgesia 1993; 76 : 247-52) measured very high plasma levels of sufentanil (> 20 ng/ml) when given alone during cardiac anesthesia. As sufentanil should always be combined with an hypnotic agent during anesthesia, much lower concentrations should be given. In more recent studies using TCI sufentanil with the Gepts model in combination of propofol, isoflurane and midazolam used target blood concentration between 1 and 10 ng/ml during cardiac surgery and between 0.1 and 1 during general surgery. (Barvais et al. J Cardiothorac Vasc Anaest 2000 ; 14 : 402-8 ; Slepchenko et al. Anesthesiology 2003 ; 98:65-73; Pandin et al. Anesthesiology 2000 ; 93: 653-61). It has to be noticed that in some countries, recommended concentrations are already included into the sufentanil drug label. In France, target concentrations between 0.4 and 2 ng/ml for cardiac and 0.15 and 0.6 for non-cardiac surgery are included. In the latest version of Miller's Anesthesia, therapeutic concentrations of sufentanil are described as follows : predominant agent (5-10 ng/ml; major surgery (1-3 ng/ml); minor surgery (3-6 ng/ml); spontaneous ventilation (< 0.4 ng/ml) and analgesia (0.2-0.4 ng/ml).

As Gepts et al didn't study the pharmacodynamic of sufentanil in her study, the effect site concentration is calculated using a ke_0 developed by Scott et al using a different study population than Gepts et al and measuring EEG changes as a measure of drug effect, therefore the effect site concentration information should be handled with care (Scott et al, Anesthesiology 1985, 62, 234-241). However, using the method of time to peak effect algorithm whereby t_{peak} is set at 5.6 min after bolus, the stability of the algorithm has been proved, even resulting in beneficial clinical effects when effect compartment TCI has been administered. (Derrodde et al. Brit J Anaesth 2003, 842-9).

In patients with compromised health, (ASA 3-4), targets should be titrated with care. It might be recommended to start with a lower target and increase following onset of drug effect.

The information above is extracted from publications and we recommend the full reading and understanding of the relevant publications on this model, and the prescribing information, prior to use.

At the release of this document, Sufentanil administration mode using the TCI is recommended only in countries where prescribing information is available. For other countries, the "Clinical Trial" option should be set within the Alaris® PK Editor Software until prescribing information is available.

Effect site targeting will not trigger the Soft or Hard limit alarms set for Plasma Concentrations

Profiles from TCI mode

When targeting in TCI Mode the Alaris® PK Syringe Pump will automatically calculate the flow rate profile from the specific pharmacokinetic/pharmacodynamic model for the selected drug. This section of the Directions For Use is intended to help users understand the profiled infusion and the performance accuracy attained from the TCI pump.

Induction Bolus and maintenance rates are displayed before starting the titration. When initially starting the infusion or after increasing the target (plasma or effect) concentration by titration, the pump will first deliver a bolus dose through a typically short, high rate infusion. On completion of this bolus, the pump will immediately switch to a lower maintenance rate (when plasma target mode is used) or will pause for a period of time before switching to a lower maintenance rate (when effect site targeting mode is used). Once the maintenance phase is reached, any reduction made to the target (plasma or effect) concentration will typically result in the infusion rate reducing to zero until the predicted plasma (or effect) concentration reduces the new target value.

The Alaris® PK Syringe Pump updates the pharmacokinetic model driving the plasma (or effect) concentration prediction and the infusion rate every 10 seconds. The infusion rate graph, shown on page 3, were measured in accordance with the protocol described in the IEC60601-2-24¹ Standard, with the data sample period reduced from 30 to 10 seconds.

The pump solves the pharmacokinetic/pharmacodynamic algorithms so that the target (plasma or effect) concentration is attained as rapidly and as accurately as possible. However, the User may need to take into consideration the limitations of the physical system in attaining the target (plasma or effect) concentration; this includes:

- The limit on the flow rate permitted by the infusion pump mechanism;
- The limit on the flow rate permitted by the syringe size;
- The patient / drug dose limitation from the prescribing information to insure the safety of the administration;
- The variation in individual patient response to reach the plasma (or effect) concentration;
- The model specific cap rate.

A true assessment of the performance of the Alaris® PK Syringe Pump can be made if the volumetric error, that is the difference between the actual volume infused and the predicted volume infused, is calculated. For the performance graphs shown on page 3, over a one hour period, the Alaris® PK Syringe Pump has a mean volumetric accuracy in TCI Mode better than $\pm 5\%^2$.

By measuring the volume from the flow rate profile delivered from the Alaris® PK Syringe Pump and then introducing this into a reverse pharmacokinetic model the predicted plasma (or effect) concentration can be calculated from the flow rate. These are illustrated on page 3, showing the typical performance of the system against changes in the target plasma (or effect) concentration for a typical, idealised profile. For the same targeted profile, the deviation of the predicted plasmatic (or effect) concentration (back calculated from the volume collected) from the expected Ideal plasma (or effect) concentration, results from the volumetric inaccuracy of the system (pump and syringe). The Alaris® PK Syringe Pump will track the predicted plasma (or effect) concentration to within $\pm 5\%^2$ of that calculated by pharmacokinetic model over a one hour period. Flow rate inaccuracies and start-up delays may decrease the accuracy of the predicted plasma (or effect) concentration particularly where high syringe drug concentrations are used in conjunction with large sizes of syringes and low target plasma (or effect) concentrations as the syringe plunger motion over time (proportional to the flow rate accuracy) will be significantly reduced.

Note: For a given drug concentration, the volumetric error is proportional to the dose rate error. Knowledge of the system accuracy over different time intervals may be of interest when assessing the impact of administering short-half life drugs. In these circumstances, short-term fluctuation in the infusion rate could have a clinical impact that cannot be determined from the performance profiles shown in Figures below. In general, the volumetric error will increase with small induction and maintenance rates, which may occur with large volume syringes, high syringe concentrations, low patient weights and low target (plasma or effect) concentrations. For applications where system accuracy is important, maintenance rates less than 1.0 ml/h are not recommended; syringe sizes, drug concentrations / dilutions and target (plasma or effect) concentrations should be selected accordingly to ensure the maintenance rate exceeds this lower limit.

Note:

¹ IEC60601-2-24: Particular Requirements for the Safety of Infusion Devices;


² 95% Confidence / 95% Population.

Profiles from TCI mode - BD 50ml Syringe

Target Effect


Infusion Rate vs Target Effect Site Concentration

Drug Concentration: 1.5 µg/ml
Volumetric accuracy: +0.8%


Predicted vs Ideal Effect Site Concentration

Drug Concentration: 1.5 µg/ml
Volumetric accuracy: +1.6%

Predicted Plasma vs Predicted Effect Site Concentration

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
Drog koncentration: 1,5 µg/ml

Alaris is a registered trademark of CareFusion Corporation or one of its subsidiaries. All rights reserved.
All other trademarks are property of their respective owners.

© 2005-2010 CareFusion Corporation or one of its subsidiaries. All rights reserved.

This document contains proprietary information of CareFusion Corporation or one of its subsidiaries, and its receipt or possession does not convey any rights to reproduce its contents, or to manufacture or sell any product described. Reproduction, disclosure, or use other than for the intended purpose without specific written authorization of CareFusion Corporation or one of its subsidiaries is strictly forbidden.

CareFusion Switzerland 317 Sàrl,
CH-1180, Rolle

CareFusion U.K. 305 Ltd., RG22 4BS, UK

1000DF00450 Issue 1