
- (EN) Collection Bottle and Water Seal
- (FR) Flacon de récupération, et joint d'étanchéité
- (NL) Verzamel fles en waterafdichting
- (DE) Sammelflasche und Wasserschloss
- (IT) Vaso di raccolta e tenuta idraulica
- (PT) Garrafa Coletora e Vedante Hidráulico
- (RU) Сборная бутыл и водяной затвор
- (ES) Frasco de recolección y sello de agua
- (SV) Uppsamlingsflaska och vattenlås

Rx Only

(EN) Instructions for use	2
(FR) Mode d'emploi	10
(NL) Gebruiksaanwijzing	18
(DE) Bedienungsanleitung	26
(IT) Istruzioni per l'uso.....	34
(PT) Instruções de utilização	42
(RU) Инструкция по эксплуатации.....	50
(ES) Instrucciones de uso.....	58
(SV) Bruksanvisning	66

Table of Contents

Safety Instructions.....	2
Intended Use.....	2
User Responsibility	2
Definitions	3
Specifications	3
Pre-Use Checkout Procedure	4
Equipment Setup.....	5
Disassembly – Collection Bottle.....	6
Cleaning	6
Parts List	7
Troubleshooting	8

Safety Instructions

This manual provides you with important information about the Collection Bottle and Water Seal. To ensure the safe and proper use of this device, Read and understand all of the safety and operating instructions. If you do not understand these instructions, or have an questions, refer to service manual, contact your supervisor, dealer of the manufacturer before attempting to use the device.

Intended Use

The Ohio Medical Suction Bottle Assembly and Water Seal are designed for the collection of various body fluids.

- The water seal is used with the Thoracic Regulator to prevent air from entering the pleural cavity, and to indicate the presence of air leaks in the lungs and/or the collection system.
- The suction bottle assembly with a float type shut-off is designed to prevent fluids from entering the vacuum regulator or central suction system.
- The collection bottle includes a splash deflector to control frothing.

User Responsibility

This Product will perform as described in this operating manual and accompanying labels and/or inserts, when assembled, operated, maintained and repaired in accordance with the instructions provided. This Product must be checked periodically. A defective Product should not be used. Parts that are broken, missing, worn, distorted or contaminated should be replaced immediately. This Product or any of its parts should not be repaired other than in accordance with written instructions provided by Ohio Medical and by Ohio Medical trained personnel. The Product must not be altered without the prior written approval of Ohio Medical's Quality Assurance Department. The user of this Product shall have the sole responsibility for any malfunction which results from improper use, faulty maintenance, improper repair, damage, or alteration by anyone other than Ohio Medical.

Definitions

WARNING: Possibility of injury to the patient or the operator exists.

CAUTION: Possibility of damage to the equipment exists.

WATER SEAL LINE	= water seal line
UNLOCK LOCK	= unlock, lock
PATIENT	=  patient
VAC	= vacuum

 0086	Symbol indicates the device complies with the requirements of Directive 93/42/EEC concerning medical devices		
	Warnings and Cautions		Manufacturer
REF	Catalog Number	LOT	Manufacturer's Batch Code
	MR Conditional		Consult instructions for use
EC REP	Manufacturer's authorized representative in the European Union		
Rx Only	Caution: Federal (USA) Law restricts this device to sale by or on the order of a licensed healthcare provider.		

Specifications

Vacuum:	0 to 650 mm Hg
Accuracy:	± 10 % full scale
Operating Temperature Range:	39.2°F to 120.2°F (4°C to 49°C)
Storage Temperature Range:	-0.4°F to 159.8°F (-18°C to 71°C)
Operating and Storage Relative Humidity:	5 to 95%
Applied Standards:	ISO 10079-3 : 2014

Usable Volume	3200 ml	1300 ml
Weight	510 g (1.1 lb)	284 g (0.6 lb)
Height	27 cm (10.6 in)	21 cm (8.4 in)
Diameter	16 cm (6.3 in)	12 cm (4.8 in)
Shut Off Volume	3200 ml	1300 ml
Max Volume	3785 ml	1892 ml
Increments	100 ml	50 ml

WARNING

⚠ Do not use this device in the presence of flammable anesthetics. Explosion hazard may exist.

CAUTION

⚠ This device is to be used only by persons who have been adequately instructed in its use.

⚠ Only a trained individual should repair this equipment.

⚠ Not for field or transport use.


Pre-Use Checkout Procedure

WARNING

⚠ Perform the Pre-Use Checkout Procedure on the Collection Bottle before using on each patient. If there is any failure, remove from service and repair.

Note: If leaks occur, check O-rings and or seals.


Collection Bottle

Water Seal


WARNING

⚠ Bubbles in Water Seal indicate leaks in suction system or patient.

Equipment Setup

Use hospital-supplied suction tubing between the end piece and the collection container, and between the patient port and the patient (minimum inside diameter is 6mm).

WARNING:

⚠ This device must be in the vertical position to function.

MR Conditional

MRI testing demonstrated these devices are MR Conditional under the following conditions ONLY:

Static magnetic field of 3-Tesla or less


Maximum spatial gradient magnetic field of 720-Gauss/cm or less

IMPORTANT NOTE: This product is intended for use inside of the MRI environment (e.g., in the MR system room). However, this product should **NOT** be utilized directly inside of the MR system (e.g., inside of the bore of the scanner) during the operation of the scanner. The assessment of magnetic field interactions for this product specifically involved evaluations of magnetic field interactions, **ONLY**, in association with a 3-Tesla MR system

Disassembly

Cap Removal

When the collection bottle shut-off vacuum feature stops vacuum and flow, the device must be emptied to continue the suctioning procedure.

Cleaning

WARNING

- ⚠ After patient use, suction equipment may be contaminated.
- ⚠ To reduce service personnel exposure to hazardous contamination, clean and disinfect all suction equipment before disassembly.

CAUTION

- ⚠ To prevent damage to equipment, completely rinse and wipe dry before autoclaving.
- ⚠ Before cleaning equipment check to be certain the cleaning agent is compatible.
- ⚠ After each cleaning or sterilization, inspect components and replace any that are damaged. Reassemble and perform pre-use check prior to returning to service.

Note: To ensure the longest life of equipment: Clean all parts with a mild soap, rinse and dry. Ohio Medical recommends replacing product after 30 cycles of cleaning and disinfecting or sterilization.

Note: Sparingly lubricate the rubber O-ring on the cap with Vac Kote® 37951 M (6700-0092-200) when reassembling the overflow trap following cleaning and sterilization.

Table 1 – Cleaning Agent Compatibility^β and Sterilization Guidelines

	Collection Bottle	Cap and Float	Water Seal
Recommended Cleaning Agents	Soaps Detergents Ammonia Alcohol†	Soaps Detergents	Soaps Detergents
Steam Sterilization*£	3-5 min. at 269.6°F (132°C)	3-5 min. at 269.6°F (132°C)	3-5 min. at 269.6°F (132°C)
	3 min. at 273.2°F (134°C)	3 min at 273.2°F (134°C)	3 min. at 273.2°F (134°C)
	1 min. at 285.8°F (141°C)	1 min. at 285.8°F (141°C)	1 min. at 285.8°F (141°C)
	15-25 min. at 249.8°F (121°C)	15-25 min. at 249.8°F (121°C)	15-25 min. at 249.8°F (121°C)


* Preferred temperature and duration in bold type.

£ Orientation – Bottle and jar upside down, all other parts upright or flat, parts not touching.

† Alcohol may be used, but may cause cracking or crazing to highly stressed area with repeated use.

β Do not use Sodium Hypochlorite (bleach), Quaternary Ammonium, Glutaraldehyde acetone, methylene chloride or chloroform.

Parts List

1.	0221-6289-300	Bushing, wing	13.	6700-0140-500	Float guide
2.	0221-6265-300	Tube, inlet vacuum bottle	14.	0221-6291-300	Float
3.	0210-0629-300	O-ring 2-010 neoprene	15.	0221-6251-300 0221-6252-300	Cap, underwater seal 1 gal. Cap, underwater seal 1/2 gal.
4.	0221-6237-300	Wing, lock	16.	6700-0001-500	Bracket, wall (int'l, small mouth)
5.	0221-6242-300	Link, conductivity	17.	0221-6267-300	Bracket, bottle wall/stand
6.	0221-6266-300	Cap, vacuum bottle w/gasket	18.	0221-6236-300	Bracket, bottle, floor
7.	0210-1063-300	Gasket, cap vacuum trap bottle	19.	0221-6256-300 0221-6257-300	Tube-underwater seal, 10-5/32L" Tube-underwater seal, 8-15/16L"
8.	6700-0183-400	Nut, tube bushing	20.	0212-0302-600 0212-0301-600	Bottle, plastic (3000ml) Bottle, plastic (1300ml)
9.	0221-6290-300	Deflector, splash			
10.	6700-0092-500	O-ring 2-005 BUNA-N			
11.	6700-0179-400	Stem, stem & seat assembly			
12.	6700-0180-400	Filter screen, cap & float assembly			

Troubleshooting

Problem	Probable Cause	Remedy
Float Stuck	1. Float neck dirty	1. Clean float and float neck
No Suction	1. Leaking Connections 2. Cap Assembly Loose 3. Suction Source not turned "ON" 4. Collection Bottle Full	1. Check connections and reconnect 2. Check Cap Assembly and reconnect to Bottle 3. Turn Suction Source "ON" 4. Empty Collection Bottle
Fluid Exits Bottle	1. Cap Assembly Loose 2. Float doesn't rise, float neck dirty 3. Crack in Bottle	1. Check Cap Assembly and reconnect 2. Clean float and float neck 3. Replace Collection Bottle
Cloudy/Opaque Bottle	1. Improper Cleaning agents used 2. Excessive Autoclaving	1. Clean Bottle with recommended cleaning agents 2. Replace Collection Bottle

© 2018 Ohio Medical, LLC.

This document contains information that is proprietary and confidential to Ohio Medical, LLC. Use of this information is under license from Ohio Medical, LLC. Any use other than that authorized by Ohio Medical, LLC is prohibited.

Ohio Medical and the Ohio Medical Logo are registered trademarks of Ohio Medical, LLC. Loctite is a registered trademark of Henkel Corporation. Cavicide is a registered trademark of Metrex Research Corporation.

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Detta dokument innehåller information som är upphovsrättsskyddad och konfidentiell för Ohio Medical, LLC.
Användning av informationen sker under licens från Ohio Medical, LLC.
All användning som inte godkänns av Ohio Medical, LLC är förbjuden.
Ohio Medical och Ohio Medical-logotypen är registrerade varumärken som tillhör Ohio Medical, LLC.
Loctite är ett registrerat varumärke som tillhör Henkel Corporation.
Cavicide är ett registrerat varumärke som tillhör Metrex Research Corporation.