GENDEX®

765DC

Intraoral x-ray system
USER MANUAL

1-17

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69-85

Table of Contents

General Description	Page 2
Identification Labels	Page 3
Compliance with Applicable Standards	Page 4
Explanation of Symbols on Technical Labels	Page 6
Safety	Page 7
Operating Controls	Page 8
Operating Procedure	Page 9
Configuring Control Settings	Page 10
Default Exposure Tables	Page 12
Error Display and Conditions	Page 13
System Specifications	Page 14
Recommended Maintenance	Page 17
Cleaning and Disinfecting	Page 17


765DC INTRAORAL X-RAY SYSTEM

GENERAL DESCRIPTION

This User Manual is intended to assist the Operator in the safe, and effective operation of the 765DC high frequency Intraoral x-ray system. It contains the information necessary to operate the equipment, the technical specifications, and the procedures to perform its simple maintenance. Before attempting to operate the equipment, read this manual thoroughly. Pay particular attention to all Warnings and Cautions incorporated in it.


✓ The original language of this manual is English.

SYSTEM COMPONENTS

IDENTIFICATION LABELS

The 765DC Tubehead, Master Control, and Cone have identification labels that specify the model number, serial number and applicable product approval listings. On specified components, subject to U.S. Government Radiation Performance Standards 21 CFR, Subchapter J, a certification statement is included with other required information.

COMPLIANCE WITH APPLICABLE STANDARDS

The 765DC complies with the following standards:

34SL

General (electrical/mechanical) safety:

- UL60601-1 File Number: E185414

MEDICAL EQUIPMENT WITH RESPECT TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH UL60601-1, IEC60601-1, CAN/CSA C22.2 NO,601.1-M90, AND TO THE FOLLOWING PARTICULAR STANDARDS, IEC60601-2-32, IEC60601-2-7 IEC60601-2-28, IEC60601-1-2.

- **Radiation protection** IEC60601-1-3 Radiation protection/x-ray equipment.

The certified components of the 765DC system comply with Radiation Performance Standards 21 CFR, Subchapter J, at the time of manufacture.

The equipment must only be installed, and operated in accordance with the safety procedures and operating instructions given in this User Manual and in the Installation/Maintenance Manual for the purposes and applications for which it was designed. Modifications and/or additions to the equipment may only be carried out by Gendex or by third parties expressly authorized by Gendex to do so. Such changes must comply with legal requirements as well as with the generally accepted technical rules. It is the responsibility of the user to ensure that existing legal regulations regarding installation of the equipment with respect to the building are observed.

IMPORTANT! X-RAY PROTECTION

X-ray equipment may cause injury if used improperly.

The instructions contained in this manual must be read and followed when operating the 765DC. Your Gendex dealer will assist you in placing the 765DC in operation.

The 765DC intraoral x-ray system provides a high degree of protection from unnecessary X-radiation. However, no practical design can provide complete protection, nor prevent operators from exposing themselves or others to unnecessary radiation.

Authorized Representative:

GENDEX Dental Systems s.r.l.
Via A, Manzoni 44
20095 Cusano Milanino - Mi
Tel: +39.02.61800.81
Fax: +39.02.61800.809

GENDEX Dental Systems
901 West Oakton Street
Des Plaines, IL 60018-1884
Tel: (847)640-4800
Fax (847)640-4891

EC Declaration of Conformity

A. Name/Description of product:

Gendex 765DC Intraoral X-Ray System

REF: Model 110-0154G1,G2,G3

Class: IIb

B. Reference Standards under which conformity is declared:

IEC60601-1:1988 +A1:1991 +A2:1995

IEC60601-1-1:2000

IEC60601-1-2:2001

IEC60601-2-28:1993

IEC60522:2003

IEC60601-1-3:1994

IEC60601-2-32:1994

IEC60878:2003

IEC60601-1-4:1996

IEC60336:2005

EN 980:1997

IEC60601-2-7:1998






IEC60417:1973

ISO14971:2000

C. GENDEX Dental Systems, Inc. declares that the products described herein meet all the applicable Essential Requirements of the EC Medical Device Directive 93/42/EEC in Annex I. For Class IIb products described herein, the product is manufactured, inspected, tested, and released in accordance with the approved quality assurance system established in accordance with ISO 9001, EN 46001 and Annex II of the EC Medical Device Directive under the supervision of the British Standards Institution, a Notified Body carrying the Notified Body No. 0086.

John R. Miller
Director of QA/RA
GENDEX Dental Systems

EXPLANATION OF SYMBOLS ON TECHNICAL LABELS

	Type b: Protection against electric shock (UL60601-1:2003)
	Consult written Instructions in this Operator's Manual
	X-RAY SOURCE Assembly
	WARNING X-RAY THIS X-RAY UNIT MAY BE DANGEROUS TO PATIENT AND OPERATOR UNLESS SAFE EXPOSURE FACTORS AND OPERATING INSTRUCTIONS ARE OBSERVED.
	X-RAY EMISSION
	Mains HOT WIRE
	Mains NEUTRAL WIRE
	Earth Ground
	Electronic Equipment No recycle into normal trash

Classification

The 765DC is a UL Class I Type b equipment (UL60601-1:2003).

SAFETY

Radiation safety

Only qualified and authorized personnel may operate this equipment observing all laws and regulations concerning Radiation Protection.

The operator at all times must remain 6FT. (2m) from the focal spot and the X-ray beam for operator protection.

Full use must be made of all radiation features on the equipment, and of all radiation protection devices, accessories, and procedures available to protect the Patient and Operator from x-ray radiation.

Electrical safety

Covers on the equipment should only be removed by qualified and authorized service personnel.

This equipment must only be used in rooms or areas that comply with all applicable laws and recommendations concerning electrical safety in rooms used for medical purposes, e.g. IEC, US National Electrical Code, or VDE standards concerning provisions of an additional protective earth (ground) terminal for power supply connection.

This equipment must always be electrically disconnected from the mains electrical supply before cleaning or disinfecting.

ORDINARY MEDICAL EQUIPMENT without protection against ingress of liquids. No water or any other liquid should be allowed to leak inside the equipment, as they may cause a short-circuit and corrosion.

Explosion safety

This equipment must not be used in the presence of flammable or potentially explosive disinfecting gases or vapors, which could ignite causing personal injury and/or damage to the equipment. If such disinfectants are used, the vapor must be allowed to disperse before using the equipment


Mechanical safety

Where complete safeguarding of the equipment is not possible, due care must be taken to ensure that no part of the user's or patient's body or clothing can be trapped or injured by any part of the equipment. In particular, make sure that fingers are not caught or pinched in the articulated arm when closing it.

OPERATING CONTROLS

1. **TIME SELECTION** - Allows the operator to override the anatomical preset time and adjust the exposure time up or down to a specific value on the R10 scale.
2. **ANATOMICAL TIME SELECTION** - Allows the operator to select the exposure time based on the tooth area desired.
3. **IMAGING TYPE SELECTION** - A toggle function that allows the operator to select between the use of film or a digital imaging system. This will modify all of the preset anatomical time selections.
4. **PATIENT SELECTION** - A toggle function that allows the operator to select either an adult or a child. This will modify all of the preset anatomical time selections.
5. **READY INDICATOR LAMP** - Indicates that the "On" switch has been pushed and the system is ready for operation.
6. **COOL-DOWN INDICATOR LAMP** - This lamp will light when the duty cycle of the system requires that the tubehead cool before the next exposure.
7. **RADIATION INDICATOR LAMP** - This lamp will light when the exposure button is pressed indicating that x-ray radiation is being produced.
8. **PUSH-BUTTON EXPOSURE SWITCH** - Allows the operator to initiate the exposure directly from the control panel when it is mounted outside the operator.
9. **COIL-CORD EXPOSURE SWITCH (Optional – Not Shown)** - Allows the operator to exit the operator and then initiate the exposure.
10. **SECONDARY TUBEHEAD CONTROLS** - This allows the operator to verify and adjust if necessary the anatomic exposure times directly at the tubehead. A lamp also indicates whether the adult or child setting has been selected.
11. **POWER SWITCH** - Depressing the switch to the "on" position will supply power to the control. The Ready Indicator Lamp on the control panel will indicate that the system is ready for use.

OPERATING PROCEDURES

1. Switch the power switch to the "On" position. The Ready indicator lamp on the front of the control panel will light.
2. Select or verify that the unit is set for either film or digital imaging. The indicator light should be lit next to the appropriate graphic icon. Press the button between the two icons to toggle to the correct selection.
3. Select or verify that the system is set for the appropriate patient selection (i.e. an adult or child). Press the button between the two icons to toggle to the correct selection.
4. Press the selector button in the center of the control to toggle the indicator lamp until it lights next to the tooth (or bitewing) symbol for the area that you wish to x-ray.

Note: The exposure time defaults (the times displayed when the system is initially turned on) can be modified per the instructions in the Configuring Control Settings section on page 10.

5. If desired, the up or down buttons can be used to select a specific exposure time other than the default setting using the R10 scale. For recommendations for Occlusal film, see the Default and Recommended Exposure Tables on page 12.
6. Position the tubehead to the patient using standard accepted positioning procedures. Note: Care should be taken to not place fingers in areas where they may be potentially pinched during the movement of the arm. Also, be careful to not have the tubehead hit the wall after returning the arm to the storage position.
7. Prior to making the exposure, verify or modify the time settings as necessary using the secondary tubehead controls. Note: The unit is supplied with an 8" (20cm) cone. An optional longer 12" (30cm) cone can be used to sharpen the image as when the paralleling technique is used. A rectangular cone is available as an option to reduce the x-ray field size to that of the film and reduce the radiation to the patient.
8. When using the optional coil-cord exposure switch, it is strongly recommended that the operator exit the operatory. Note: In order to comply with regulations and good safety practices, the technique factors must be visible to the operator from the remote location.
9. Press and hold either the coil-cord or push-button exposure switch until the audible signal and exposure indicator light terminate. **(Note: Releasing the exposure button at any time will immediately terminate the exposure).**

Configuring Control Settings

The 765DC control panel is preset at the factory (called default settings) for operation with a standard configuration. This configuration assumes the use of the standard 8" focal length cone and D-speed film or a DenOptix digital imaging system. This standard configuration presets the exposure times for the anatomical selections. These presets can be changed for use with E-speed or F-speed film or another type of digital imaging system, or with the optional long cone.

- The film speed and digital option can be changed using the control panel buttons and software.
- Compensating for the cone length is done by changing a switch inside of the control box. *Note: It is recommended that this be done only by an authorized Gendex service agent.*

Changing the default settings for the film speed used

The 765DC is set at the factory for the use of D-speed film. This can be verified by turning the control and noting that the time initially displayed is 0.320 seconds. The use of E-speed or F-speed film requires modifying the exposure times used. The system can be programmed to accomplish this automatically.

1. Turn the system off.
2. Press and hold the anatomical time selector switch while turning the system on (see Figure 1).
3. The time display will now read "SP 1". This confirms the system is in the film speed menu.
4. **To change to E-speed.** Press the exposure time 'Up' button. The time display will now read "SP 2". This confirms that the system is now set for E-speed film.
5. Press the anatomical time selector switch again to store the new setting. The time displayed should now indicate 0.200 seconds.
6. **To change to F-speed.** Press the exposure time 'Up' button again until the display reads "SP 3".
7. Press the time selector switch to store the setting. The time displayed should now indicate 0.125 seconds.

Figure 1

Exposure time select button

Anatomical time selector switch

Film or digital selector switch

Changing the Digital Option

The 765DC is set at the factory for use with the Gendex DenOptix digital imaging system. This can be confirmed by turning the system on and pressing the button to select digital. The time displayed should now indicate 0.160 seconds.

1. Turn the system off.
2. Press and hold the anatomical time selector switch while turning the system on (see Figure 1).
3. The time display will now read "SP 1,2, or 3" depending on which film type is configured.
4. Press the film or digital selector switch until the display reads "D1". Press the exposure time 'UP' button. The display will now indicate "D2". This indicates that the system is now configured for a CCD type digital imaging system.
5. Press the anatomical time selector switch again to store the new setting.

Changing the power on selection

The 765DC is set at the factory to have the film option for an adult selected when the unit is first turned on. If desired, the unit can be set so that it automatically selects the digital option or a child instead when the control is turned on.

1. Turn the system off.
2. Press and hold the anatomical time selector switch while turning the system on (see Figure 1).
3. The time display will now read either "SP 1,2, or 3" or "D1" or "D2" depending on whether the film or digital option is selected. Note that the green indicator light next to the selection is flashing.
4. To select for the unit to come on in the digital mode, make certain that the green light is flashing next to the digital option. Press the anatomical selector switch to store the information and turn the system off. The unit should now automatically select the digital option when first turned on.
5. To select for the unit to come on in the child mode, press the adult or child selection button until the green light is flashing next to the child option. Press the anatomical selector switch to store the information and turn the unit off. The unit should now automatically select the child option when the unit is first turned on.

Modifying the system for use with the optional 12" (30 cm) cone

The 765DC is set at the factory for use with the standard supplied 8" (20 cm) cone. The 12" (30 cm) long cone is recommended when using the paralleling film positioning technique along with the Rinn XCP kit. Using the long cone requires the use of longer exposure times, which can be programmed into the system by changing an internal electrical setting.

This modification should be performed only by a qualified Gendex service agent. Contact your Gendex service agent to make the change.

Default and Recommended Exposure Tables

SPEED 1: D-speed Film

	8" (20cm) Cone		12" (30cm) Cone	
Anatomy Selected	Adult	Child	Adult	Child
Bitewing	0.320 second	0.160 second	0.630 second	0.320 second
Lower Molar	0.320 second	0.160 second	0.630 second	0.320 second
Upper Molar	0.400 second	0.200 second	0.800 second	0.400 second
Incisor	0.200 second	0.100 second	0.400 second	0.200 second
Bicuspid	0.250 second	0.125 second	0.500 second	0.250 second
Occlusal	0.500 second	0.250 second	1.00 second	0.500 second

SPEED 2: E-speed Film

	8" (20cm) Cone		12" (30cm) Cone	
Anatomy Selected	Adult	Child	Adult	Child
Bitewing	0.200 second	0.100 second	0.400 second	0.200 second
Lower Molar	0.200 second	0.100 second	0.400 second	0.200 second
Upper Molar	0.250 second	0.125 second	0.500 second	0.250 second
Incisor	0.125 second	0.063 second	0.250 second	0.125 second
Bicuspid	0.160 second	0.080 second	0.320 second	0.160 second
Occlusal	0.320 second	0.160 second	0.630 second	0.320 second

SPEED 3: F-speed Film

	8" (20cm) Cone		12" (30cm) Cone	
Anatomy Selected	Adult	Child	Adult	Child
Bitewing	0.125 second	0.063 second	0.250 second	0.125 second
Lower Molar	0.125 second	0.063 second	0.250 second	0.125 second
Upper Molar	0.160 second	0.080 second	0.320 second	0.160 second
Incisor	0.080 second	0.040 second	0.160 second	0.080 second
Bicuspid	0.100 second	0.050 second	0.200 second	0.100 second
Occlusal	0.200 second	0.100 second	0.400 second	0.200 second

D1: Digital (DenOptix)

	8" (20cm) Cone		12" (30cm) Cone	
Anatomy Selected	Adult	Child	Adult	Child
Bitewing	0.160 second	0.080 second	0.320 second	0.160 second
Lower Molar	0.160 second	0.080 second	0.320 second	0.160 second
Upper Molar	0.200 second	0.100 second	0.400 second	0.200 second
Incisor	0.100 second	0.050 second	0.200 second	0.100 second
Bicuspid	0.125 second	0.063 second	0.250 second	0.125 second
Occlusal	0.250 second	0.125 second	0.500 second	0.250 second

D2: Digital (Sensor)

	8" (20cm) Cone		12" (30cm) Cone	
Anatomy Selected	Adult	Child	Adult	Child
Bitewing	0.080 second	0.040 second	0.160 second	0.080 second
Lower Molar	0.080 second	0.040 second	0.160 second	0.080 second
Upper Molar	0.100 second	0.050 second	0.200 second	0.100 second
Incisor	0.050 second	0.025 second	0.100 second	0.050 second
Bicuspid	0.063 second	0.032 second	0.125 second	0.063 second
Occlusal	0.125 second	0.063 second	0.250 second	0.125 second

Error Display and Conditions

The green 'Ready' lamp, amber 'Cool-down' lamp and 4 digital 7-segment LED displays on the face of the Operator Control panel can provide indications that errors have occurred. The following chart provides an explanation of what the indicators mean and required action.

Error Indication	Explanation	Action Required
The amber 'Cool-down' lamp comes on and exposure cannot be made.	Normal operation by design. 'Cool-down' lamp comes on if too many exposures are made in a short period of time. This feature protects and extends the life of the tubehead.	Wait until lamp goes out indicating that the tube has properly cooled down.
"Err0" flashes on the 7-segment LED displays..	The exposure has been terminated prematurely. This can be caused by an equipment malfunction or the operator releasing the exposure switch prior to the end of the exposure time selected.	Pressing any key except exposure switch on the Operator Control will clear the display and restore normal operation. Care must be taken to press and hold the exposure button until the exposure is finished.
"Err1" flashes on the 7-segment LED.	Power Supply voltage was outside of required range: 108V – 132V 198V – 253V	Wait until line voltage returns to normal (indicated by the lamp being steady) or have a qualified electrician check the power line.
"Err2" is displayed on the 7-segment LED displays.	Power Supply voltage dropped below the minimum requirement (108V or 198V) during the exposure. The exposure is not terminated but output may be below normal.	Pressing any key on the Operator Control will clear the flashing and restore normal operation. If the problem persists, it is strongly recommended that a qualified electrician check the power line.
"Err3" flashes on the 7-segment LED displays.	A hardware fault has occurred.	Reset the unit by Turn OFF the power switch for 5 seconds and Turn ON the power switch, if the problem persists. Contact your authorized Gendex service agent.
"Err4" flashes on the 7-segment LED displays.	A fault has occurred and terminated the exposure.	Pressing any key except exposure switch on the Operator Control may clear the Error and restore normal operation. If the problem persists. Contact your authorized Gendex service agent.

SYSTEM SPECIFICATIONS

Cooling and duty rating specification apply at altitudes up to 12,000 ft. (3,600 m), average relative humidity not exceeding 90 percent, and maximum ambient temperature not exceeding 104 degrees F (40 degrees C).

TUBEHEAD SPECIFICATIONS

Maximum Rated Tube Potential: 65 kV at 7 mA tube current.

Leakage radiation in the loading state: Less than 250 micro Gy/h at 1 meter from the focal spot.

Minimum permissible first Half Value Layer: Greater than 1.5 mm Aluminum equivalent

Nominal Focal spot size: 0.4 mm (IEC 336/1993)

Anode angle: 12.5 degrees

Anode heat storage: 7 kHU

Duty cycle: 1:30

CONE SPECIFICATIONS

X-ray Field Size: Containable in a 6 cm diameter circle, measured at the distal end of the installed cone.

Distance from the focal spot to the distal end of the cone is: 8 or 12 inches (20 or 30 cm).

POWER SUPPLY REQUIREMENTS

<u>Nominal</u>	<u>Tolerance</u>	<u>Useful Range</u>
120 VAC	± 10%	108 - 132 VAC
230 VAC	+ 10%, -14%	198 - 253 VAC

ACCURACY OF TECHNIC FACTORS

Tube Potential	65 kV ± 3 kV
Current	7 mA ± 5%
Time	.02 - .50 Seconds ± .005 Seconds .63 - 2.00 Seconds ± 1%

Intraoral X-Ray Tube Housing

Product Data Sheet

Gendex[®] 765DC

The Gendex[®] 765DC housing is designed to operate with the Gendex[®] GX 70-10DC insert. The housing has been designed to operate on the Gendex[®] 765DC high frequency high voltage control, making up the Gendex[®] 765DC Intraoral Dental X-Ray system.

Specifications

Shockproof Housing: Zinc, hermetically sealed under vacuum with high dielectric insulating oil. Inherent to the housing is an expansion chamber designed to provide adequate compensation for the full temperature range.

Operating Temperature: +50°F to 100°F ≈ (+10°C to +40°C)

Storage Temperature: -40°F to 160°F ≈ (-40°C to +70°C)

Inherent Filtration:

Material	Thickness / mm.		AHE ¹ /mm.	Al. Equivalent	
	Min.	Max.	65 kV	Min.	Max.
Glass	1.27	1.91	0.51	0.65	0.97
Insulating Oil	6.5	7.7	0.08	0.52	0.62
Polycarbonate	16.4	16.7	0.09	1.48	1.50
Total				2.64	3.09

Radiation Output:

Technique 65kV, 7mA, 0.50sec. Source to Detector Distance 230mm (9")
410mR ±125mR (3.55 mGy ± 1.05 mGy)

Typical Half Layer: 1.9mm Aluminum at 65kV

Leakage Radiation:

Complies to 601-1-3 IEC Section 29.204 "Leakage Radiation". Maximum continuous power input equals 1,400 heat units per min. (HU/Min.) tube is the limiting factor. $HU/Min. = (kV) \times (mA) \times (1.4) \times (Exp. Time (sec) / 60)$

$$1,400 HU/Min. = 65 kV \times 0.27mA \times 1.4 \times 60 sec.$$

X-Ray Coverage:

The insert target angle is 12.5°

- Circular cone at 8 inches (200 mm) SSD² 60 mm Dia. +0, -3mm.
- Circular cone at 12 inches (300 mm) SSD² 60 mm Dia. +0, -3mm.
- Rectangular cone at 8 inches (200 mm) SSD² 36mm x 46mm

¹ AHE Aluminum Half value Equivalent correlation

² SSD Source to Skin Distance

Weight:

13.0 lbs.

5.9 kg.

Maximum Head Voltage:

Anode to Cathode

65 kV


Anode to Ground

70 kV

Thermal Characteristics:

- Housing Heat Storage Capacity 250 kHU's (177J)
(see Graph for more detail)
- Maximum Cooling Rate 3.3 kHU's/Min.

$$Heat Units = (kV * mA * time in seconds * 1.4)$$

RECOMMENDED MAINTANCE

In the interest of equipment safety, a maintenance program must be established. It is the owner's responsibility to arrange for this service and to assure that the personnel performing this are fully qualified to service GENDEX dental x-ray equipment.

The calibration and adjustments to be performed are listed in the System and Function Checks section of the 765DC Installation and Maintenance manual supplied with the product.

It is recommended that the following System and Function Checks and Maintenance Procedure be performed at installation and every twelve months thereafter.

- System Function Checks listed in Section E, page 35 of the 765DC Installation / Maintenance Manual.
- Suspension Arm Maintenance as listed in Section F, page 43 of the 765DC Installation / Maintenance Manual.

Note: To avoid any potential hazard to operators or patients, any unusual operation, mechanical functionality, or debris should be reported to your service agent immediately.

CLEANING AND DISINFECTING

All external and cleanable surfaces can be cleaned and/or disinfected using the following list of commonly used solutions:

Isopropyl alcohol
Phenylphenl
Benzyl-P-Chlorophenol
Iodophor
Disobutyl phenoxethoxyethyl Dimethyl Benzyl Ammonium Chloride

Note: Use bleach or glutaraldehydes only occasionally.

Verify the contents of any cleaning or disinfecting solution against the list above applying it the any surfaces of the 765DC.

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Gendex Dental Systems S.r.l.
Ivoclar Vivadent Pty Ltd.