ISO 9001
EN ISO 13485

CE
0499

CFDA

ANVISA

MK_{300L} MK_{400L}

DSM-800S DSM-600S

User Manual

Indications for Safety Use

Read this manual carefully.
This manual is for user's safety
and preventing any property-
loss. Before using our device,
please read this manual
inevitably.

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DS MAREF
DAESUNG MAREF CO.,LTD.

Introduction

Thanks for purchasing our product. This product is an Air Compression Therapy System to improve the blood circulation of the human body. Patient can also purchase or use it personally after receiving doctor's indications definitely.

This product is comprised of intermittent pneumatic controller, cuffs with 6 air chambers, and connectable hose. The operating principle is that the air from the controller will be delivered to the cuff with 6 air chambers and the air pressurizes sequentially the chambers from 1st to 6rd.

This user's manual contains the information how to use, keep and repair this product. Please, be well informed of this manual to use this product correctly and to prevent the malfunction for the best effect.

Please keep this manual with this system always and please carefully read it before using the system properly.

Guidance for indication

 caution	Wound or property-loss can be occurred in case of disregarding the instruction
 warning	Death or serious wound can be occurred if users do not follow the instruction.
 prohibition	Instruction to be prohibited.
 prohibition of disjoining	Do not disjoint the device.

Electrical Safety

This device is for domestic use only. Please use a suitable AC cord for each country. (Fire accidents and electric shocks possible)  prohibition	Stop using the compressor when it makes any kinds of abnormal smells. In this case, top the compressor, and plug it out, and then call customer service of a seller. (It may cause the fire accident or electric shot by immoderate use)  prohibition
Pull out plug when not using or providing against lighting, thunder.  caution	Grip the plug, not a power cord with hands when you plug in or out  prohibition
Do not plug in or out with wet hands (Electric shocks possible)  prohibition	Do not bend or twist the power cord. (Fire accidents & electric shocks possible)  prohibition
Do not touch signal input, signal output or other connectors, and the patient simultaneously.  prohibition	External equipment intended for connection to signal input, signal output or other connectors, shall comply with relevant IEC Standard(e.g., IEC60950 for IT equipment and IEC60601-1 series for medical electrical equipment). In addition, all such combination-system- shall comply with the standard IEC60601-1 and/or IEC60601-1-1 harmonized national standard or the combination. If, in doubt, contact qualified technician or your local representative.  prohibition

Regarding cuff

Do not use these cuffs together with accessories, watch etc. (Physical damages, cuff breakage possible)  caution	Do not use cuffs on undressed body. Use the cuffs with wearing thin under wears. (Cuff's endurance can be weakened due to sweat cosmetics and sebum)  prohibition
Do not keep cuffs near to heating materials such as stove, cigarettes etc. (Air leaking or fire accidents possible)  caution	Do not fold or bend the wire by force and put any heavy object on the wire. (It may cause the fire accident)  caution
Do not use other tools, except our machines, to blow air into cuffs.  prohibition	Do not stain cuffs by oil, benzene, alcohol, gasoline etc. (Anti-endurance possible)  prohibition

Safety for body

Forbid children to use the device. Serious injuries can be occurred when the air blows into cuffs.  warning	Do not use air hose to nose, mouth or ears.(Any injuries possible)  warning
Prohibit to use the device in case of artificial-heart etc.  warning	Use the cuffs with wearing thin cloths. Zip up perfectly the cuff as the unzipped zipper can damage the skin during using.  warning
Prohibit using the device to the patients having any artificial material inserted like artificial-joint, metal-inserting, silicon, etc in their body. (it can be the cause of inflammation can or other side effects)  warning	Use the device after consulting to a doctor' definitely before using the device or do not use it in case of lymph edema, vein diseases vein, heart diseases, vascular cancer, a patient under cancer therapy, variant, gangrene, recent artificial skin graft, or deep vein thrombosis.  warning
Use the device after consulting to a doctor's definitely before using the device or do not use it in case of heart disease, hypertension, malignant tumor, acute inflamed skin disease, high pyrexia, hyposthenia or a patient having any diseases now.  warning	Stop using the device when you feel any abnormal condition during using, and consult to doctor.  warning
Do not use the device in a drunken state.  caution	Use the device according to a doctor's prescription definitely or do not use the device to a patient with any surgical operation, especially, internal organs, nerve or brain surgery or a patient within 1 year after surgery.  prohibition

Use and Maintenance

Use the device indoors of 0~40°C (If not, any defects or fire accidents possible)  prohibition	Do not use the device in high humidity place(Bathroom and sauna etc). (Fire accidents, electric shocks possible)  warning
Do not use or keep the device near to heating equipments or a direct ray of the sun. (Defects, fire accidents, or transformations possible)  prohibition	Forbid children to touch. Disorder can be happened due to children's touch.  prohibition
Keep the device horizontally. (If leaning, defects or malfunctions possible)  caution	

Cleaning and Others

Do not use oil, benzene, alcohol, gasoline, chemicals when cleaning the device. (Use dry towels)  prohibition	Do not use oil, benzene, alcohol, gasoline, chemicals when cleaning cuffs. (Use wet towels.)  prohibition
Do not repair, reproduce or disjoint the device by yourself. (Fire accidents, electric shocks possible)  prohibition of disjointing	Be careful to not put clips, staples, needles, metals or food into the device. (Short circuits, heating generations possible) * If happened, inform it to concerned agent.  caution

Product Description 1) MK300L

(1) Introduction


This product is an Air Compression Therapy System to improve the blood circulation of the human body. There are three different pressurizing modes which the user can select so that it will be very effective to use this product up to patient's state or condition. This product has 2 channels which can use maximum 4 cuffs at the same time. When all four cuffs are connected, each channel works simultaneously and each cuff works sequential mode. If three cuffs are connected to the air socket; the rest socket won't discharge the air on the machine. This is DS Maref's own technique and we have patent for it.

Also there is a safety switch for emergency to make stop working during operation. This principle is to make stop operation by controlling "start/stop" control parts software if user presses the switch during operation. This switch is connected on each channel (two of them). For example, if user presses the CH1's safety switch during operation, the CH1's operation is stopped but CH2 is working normally.

The air chamber of cuff is comprised of 6 chambers and it pressurizes sequentially from 1st chamber to 6th chamber to make blood into the foot circulating from bottom to up.

The pressure setting is that the sensor checks the pressure every time so it will not be pressurized the cuff with the higher pressure than the pressure which the user set up. The time to reach the setting pressure can be different a little according to user's body size but pressure applying to body is kept regularly by the built-in air sensor.



(2) Name of operating function (part)


1	CH1 MODE SELECT BUTTON	6	CH2 MODE SELECT BUTTON
2	CH1 TIME SET UP BUTTON	7	CH2 TIME SET UP BUTTON
3	CH1 START/STOP BUTTON	8	CH2 START/STOP BUTTON
4	CH1 OPERATING PRESSURE SET UP KNOB	9	CH2 OPERATING PRESSURE SET UP KNOB
5	CH1 AIR SOCKET	10	CH2 AIR SOCKET

MODE setting

If the selected mode is touched, blue LED of the selected mode is displayed.

A mode (Wave mode by 1 chamber) 	The selected air chamber(s) inflate(s) and deflates by one chamber sequentially up to thigh from foot. 
B mode (Squeezing mode) 	After all the selected air chamber(s) inflate(s) sequentially up to thigh from foot, the inflated chambers deflate at once. 
C mode (B mode two times + Reversed A mode 1 time) 	It repeats that B mode operates two times and A mode operates reversely one time. 

TIME setting

- If the ▲ button is touched to increase use time, time indicator increases by 1 minutes unit.
- If the ▲ button is touched for a long time, time indicator increases quickly.
- If the ▼ button is touched to decrease use time, time indicator decreases by 1 minutes unit.
- If the ▼ button is touched for a long time, time indicator decreases quickly.
- Time range is 1~90 minutes.

PRESSURE setting

- If Pressure button is turned to "-" direction, pressure is decreased by 10mmHg and to "+" direction, it is increased by 10mmHg.
- Pressure indicator displays the set pressure value.
- Pressure range is 10~200mmHg.

START/STOP

- Individual operation (CH1/CH2)
- In order to operate CH1, the user can touch "Ⓢ" on CH1, then, "Ⓢ" will turn to blue. To stop operating, the user can touch it again, then, it will turn to red.
- In order to operate CH2, the user can touch "Ⓢ" on CH2, then, "Ⓢ" will turn to blue. To stop operating, the user can touch it again, then, it will turn to red.

How to connect hose to the controller

After connecting a connectable hose plug to air socket of the controller, connect the other remained plug of the hose to air socket of cuff.

There's a hole next to the air socket which user can connect the safety switch. Press the safety switch button, stop the operation.

Product standard and specifications

Time control	1~90 minutes
Pressure control	10~200mmHg
Product dimension	415 x 310 x 160 mm
Product weight	6.5kg (Only body)

Rated voltage	AC 220~230V, 50/60Hz
Fuse capacity	AC 250V / 2A
Power consumption	70VA

※ The label is located on the bottom of the device.

Product Description 2) MK400L


(1) Introduction

This product is an Air Compression Therapy System to improve the blood circulation of the human body. There are three pressurizing modes, reverse modes and optional care functions which the user can select so that it will be very effective to use this product up to patient's state or condition.

The air chamber of cuff is comprised of 6 chambers and it pressurizes sequentially from 1st chamber to 6th chamber to make blood into the circulating from bottom to up. The pressure setting is that the sensor checks the pressure every time so it will not be pressurized the cuff with the higher pressure than the pressure which the user set up.

The time to reach the setting pressure can be different a little according to user's body size but pressure applying to body is kept regularly by the built-in air sensor.


(2) Name of operating function (part)

1	OPTIONAL CARE BUTTON	6	OPERATING PRESSURE SET UP KNOB
2	MODE SELECT BUTTON	7	START/STOP BUTTON
3	TIME SET UP BUTTON	8	REVERSE BUTTON
4	SPEED SET UP BUTTON	9	AIR SOCKET
5	REMOTE		

MODE setting

Touch the MODE button until you can select your wanted mode from A mode, B mode, C mode. If you touch "●" of MODE, MODE indicator LED is displayed to the selected mode in order of A mode > B mode > C mode > A mode



A mode (Wave mode by 1 chamber)	The selected air chamber(s) inflate(s) and deflates by one chamber sequentially up to thigh from foot. 
B mode (Squeezing mode)	After all the selected air chamber(s) inflate(s) sequentially up to thigh from foot, the inflated chambers deflates at once. 
C mode (B mode two time + Reversed A mode 1 time)	The selected air chamber(s) inflate(s) and deflates by two chambers sequentially up to thigh from foot. 

REVERSE setting

If "●" button of REVERSE is touched, after LED of "●" is displayed, the selected mode is changed to REVERSE mode automatically. Press "●" button of REVERSE to change one mode of from A mode, B mode, C mode to Reverse mode as the pictures.


A mode (+ reverse mode)	It repeats that A mode is operated reversely after two normal operations of A mode. 
B mode (+ reverse mode)	It repeats that B mode is operated reversely after two normal operations of B mode. 
C mode (+ reverse mode)	It repeats that C mode is operated reversely after two normal operations of C mode. 

■ TIME setting

- If the ▲ button is touched to increase use time, time indicator increases by 1 minutes unit.
- If the ▲ button is touched for a long time, time indicator increases quickly.
- If the ▼ button is touched to decrease use time, time indicator decreases by 1 minutes unit.
- If the ▼ button is touched for a long time, time indicator decreases quickly.
- Time range is 1~90 minutes.

■ SPEED setting

- If the ▲ button is touched to increase the speed to faster level, speed indicator increases by 1 level unit.
- If the ▲ button is touched for a long time, speed indicator increases quickly.
- If the ▼ button is touched to decrease speed, speed indicator decreases by 1 level unit.
- If the ▼ button is touched for a long time, time indicator decreases quickly.
- Speed level range is 1 ~ 6 level

■ OPTIONAL CARE setting


- If you touch the number matched each chamber position, blue light is displayed on LED of the number and the selected chamber(s) is (are) inflated.

■ PRESSURE setting

- If Pressure button is turned to "-" direction, pressure is decreased by 10mmHg and to "+" direction, it is increased by 10mmHg.
- Pressure indicator displays the set pressure value.
- Pressure range is 10 ~ 200mmHg.

■ START/STOP

- If "⏻" of START/STOP is touched and red LED is changed to blue LED, it starts operating. If "⏻" of START/STOP is to stop operating and blue LED is changed to red LED, it stops operating.

■ How to connect hose to the controller

- After connecting a connectable hose plug to air socket of the controller, connect the other remained plug of the hose to air socket of cuff.
- There's a hole next to the air socket which user can connect the safety switch. Press the safety switch button, stop the operation.


■ Product standard and specifications

Time control	1~90 minutes
Pressure control	10~200mmHg
Product dimension	415 x 310 x 160 mm
Product weight	5kg (Only body)





Rated voltage	AC 220~240V, 50/60Hz
Fuse capacity	AC 250V / 2A
Power consumption	50VA
Speed	1~6

※ The label is located on the bottom of the device.

Accessories (MK300L, MK400L)

① One touch leg cuff (L/XL)	② One touch waist cuff	③ One touch arm cuff	④ Leg or arm extension zipper
MK300L-4	MK300L-1	MK300L-1	MK300L-Optional
MK400L-2	MK400L-Optional	MK400L-Optional	MK400L-Optional

⑤ Single hose / Double hose	⑥ Safety switches	⑦ Center body cuff	⑧ Full body cuff
MK300L-4	MK300L-2	MK300L-Optional	MK300L-Optional
MK400L-2	MK400L-1	MK400L-Optional	MK400L-Optional

How to use the device

(1) Notices for the first users

For a user during initial period, we recommend you to use this device from the lowest pressure intensity, and increase the pressure intensity gradually with checking the body conditions under a doctor's advice for more safe and effective use. Please make yourself as a comfortable posture for best effects.

(2) Use and setting

- 1) Set the controller on a desk or flat place horizontally and connect the AC cord (plug) the controller after checking local voltage.
- 2) Insert the plug of a connectable hose into a socket of the controller perfectly.
- 3) Connect a hose to a cuff.
- 4) Apply the selected cuff to a patient. We recommend for a patient to wear thin clothes to protect skin and sweater during operation. (Please zip up perfectly. If unzipped cuff can damage the skin or cuff during operation)
- 5) Set the controllers of MODE, TIME, PRESSURE, SPEED (MK400L) and REVERSE (MK400L).
Or press the START/STOP button without any setting for the previous setting use or factory recommended setting use.
- 6) Press START/STOP button to start operating.

How to use cuff

How to use one cuff (MK300L, MK400L)

- Apply the selected cuff from leg cuff, arm cuff, waist cuff to a patient.
- Connect one cuff and one connectable hose.
- Insert a hose plug into air socket of a machine.
- Put the hose plug into socket of a cuff and a controller perfectly to prevent air leaking during operation.


■ How to use two cuffs (MK300L, MK400L)

- Apply two selected cuffs from leg cuff, arm cuff, waist cuff to a patient.
- Connect two cuffs and two connectable hoses.
- Insert two hose plugs into two air sockets of a machine.
- Put the hose plug into socket of a cuff and a controller perfectly to prevent air leaking during operation.

※ Caution : In the case of two arm cuffs use, as arm cuff is applied near to a heart, please keep a protector near a patient or under a doctor's indications.

■ Three cuffs use (MK300L only)

- Apply three kinds of cuffs from leg cuff, arm cuff, waist cuff to a patient.
- Connect three cuffs and three connectable hoses.
- Insert three hose plugs into three air sockets of a machine.
- Put the hose plug into socket of a cuff and a controller perfectly to prevent air leaking during operation.

※ NOTICE! In the case of waist cuff use, please check whether inflation process from upper abdomen to hip or not.

■ Four cuffs use (MK300L only)


- Apply four kinds of cuffs from leg cuff, arm cuff, waist cuff to a patient
- Connect four cuffs and four connectable hoses.
- Insert four hose plugs into four air sockets of a machine.
- Put the hose plug into socket of a cuff and a controller perfectly to prevent air leaking during operation.

※ NOTICE! In the case of two arm cuffs use simultaneously, it can be dangerous for high blood pressure.


■ How to use extension zipper

- Use an extension zipper if the cuff does not fit you.
- Zip down completely and put the extension zipper between the cuff and zip up.


Leg cuff (L) and extension zipper


Leg cuff (XL) and extension zipper

Arm cuff size (one size)

Full body cuff size (one size)

Clause 5

1. Type of protection against electric shock: Class I equipment.
2. Degree of protection against electric shock: Not classified - BF applied parts
3. Classification according to the degree of protection against ingress of water: IPX0, ordinary equipment.
4. This equipment is not suitable for use in the presence of flammable anesthetics or oxygen.
5. Mode of operation: continuous operation.

Information regarding potential electromagnetic or other interference

This equipment has been tested and found to comply with the limits for medical devices in IEC/EN 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation between the equipment.
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
- Consult the manufacturer or field service technician for help.

WEEE marking

This marking shown on the product or its literature, indicates that it should not be disposed with other household wastes at the end of its working life.

To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate this from other types of waste and recycle it responsibly to promote the sustainable reuse of material resources.

Household users should contact either whether they purchased this product, or their local government office, for details of where and how they can take this item for environmentally safe recycling.

Business users should contact their supplier and check the items and conditions of the purchase contract. This product should no be mixed with other commercial waste for disposal.

Environmental conditions for transport and storage

Operating condition	Temperature (°C)	0-40
	Relative humidity(%)	10-90
	Atmospheric pressure (hPa)	700-1060
Transport and storage condition	Temperature (°C)	-20-60
	Relative humidity(%)	0-90
	Atmospheric pressure (hPa)	500-1060

Contraindications

※ MK300L and MK400L should NOT be used in the following conditions:

1. Severe arteriosclerosis or other ischemic vascular diseases.
2. Known or suspected acute deep vein thrombosis.
3. Severe congestive cardiac failure.
4. Any local condition in which garments would interfere, for example:
 - Gangrene • Untreated, infected leg wounds
 - Recent skin graft • Dermatitis

Caution

※ In case of not following the instructions, the user may suffer from injury or property damage.

Stop the device during operation

※ Stop the operation of a control and follow one's doctor's indications for the below symptoms.

- 1) When feeling difficulty of breathing or giddiness
- 2) When a part of body is felt any so strong compression
- 3) When feeling so serious pain during using a device
- 4) When body is paralyzed temporally or felt any itchy
- 5) When feeling any other strange during using a device

Recommendation

※ As this device is a medical devices for the purpose of the assist therapy, use MK300L, MK400L personally with easy way under their doctor's indication and guidance. What a patient feels weaker the compression of this system as time is goes is due to accustom oneself to the compression. So, do not use it for long time and use it definitely according to a doctor's indications.

Troubleshooting

No.	Condition	Cause	Solution
1	No electrip power	Power connection error	Is the plug correctly inserted into outlet?
2	Power on but not-operating	Power supply error	Turn off and on the power to the controller.
3	Noisy during operating	Setting condition	Horizontally installed?
			Is something laid on the device?
4	No air in the hose	Hose connection error	Is the hose correctly inserted into the device?
		Hose bent	Is there any bending place of the hose?
5	Air is injecting in wrong order.	Connection condition	Is the plug correctly inserted into outlet?

■ Others related to defect

No.	Condition	Cause	Solution
1	Weak air injection	Air hose damage	Plug correctly inserted into outlet?
		Air hose socket damage	
		Inner parts defect	
2	Power on but no operating	Inner parts defect	Contact the seller.

※ We can not be responsible for any defect occurred from user's careless use, even though warranty period.

Service and Maintenance

■ Custody

- 1) Keep the device in the place without the direct rays of the sun.
- 2) Don't keep the device in high humidity place
- 3) Keep the device in clean place to protect from dust.

When the device is not used for long time, be in the custody of it in a box after wrapping with vinyl pack.

■ Clean maintenance

- 1) After for use, clean the device with a dried towel without any moisture.
(Chemical cleaning can be the cause of a color change or defect).
- 2) Keep a cuff or a connectable hose from water or dust.
- 3) Do not disassemble the device at one's option and contact your nearest seller or manufacturer.

Warranty

Much appreciated on using our device. We, Daesung Maref are doing our best to improve the quality of our products.

※ We can not be responsible for any defect occurred from user's careless use or in case of followings, even though warranty period :

1. Disorder happened by strong impact.
2. In case user repair or reproduce internal part arbitrarily.
3. In case of using the device in prohibited place.
4. In case of against our <How To Use>
5. Cuff is articles of consumption.

DESCRIPTION	Air compression therapy system
MODEL NAME	MK300L, MK400L
SERIAL NO.	
WARRANTY	Controller : 1 year
BUYING PLACE	
BUYING DATE	

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DS MAREF
DAESUNG MAREF CO.,LTD.

History & Certificates

- 2013 Commissioned member of the Ministry of Trade Industry Energy planning committee
- 2013 Consultant of National Unification Advisory Council
- 2013 Member of Trade Industry Forum
- 2011 President Award
- 2010 Anvisa in Brazil
- 2007 Authentication of Merit certificate
- 2006 Innovation management awards
- 2006 KOTRA B2B e-Trade awards
- 2006 KFDA awards by Prime minister
- 2006 KGMP registration
- 2004 The 34 precious Technique bronze awards
- 2004 Korean World class products award
- 2004 Vice president award
- 2004 SFDA in China
- 2003 Success Design product
- 2002 ISO9001/EN13485
- 2002 CE marks (DL series, MK series, DL1200, DVT-2600)
- 1986 DS MAREF established

298-24, Gongdan-ro, Gunpo-shi, Gyeonggi-do, 435-862, Korea
Tel : +82-31-459-7211 Fax : +82-31-459-7215
E-mail : ds@dsmaref.com <http://www.dsmaref.com>

EC Representative Supply Medical Service

18, Chemin Des Courses 31100 Toulouse, France
Tel : +33-(0)5-3450-4455 Fax : +33-(0)5-3450-4461

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