GE Healthcare

Optima EMS ECG Management System HIS/EMR Interface Reference Manual

Software Version 1.0
2069773-002 Revision B

Optiam EMS
English
© 2012 General Electric Company.
All Rights Reserved.

Publication Information

The information in this manual applies only to Optima EMS ECG Management System Version 1.0. It does not apply to earlier product versions. Due to continuing product innovation, specifications in this manual are subject to change without notice.

MUSE and MAC are trademarks owned by GE Medical Systems *Information Technologies*, Inc., a General Electric Company going to market as GE Healthcare. All other trademarks contained herein are the property of their respective owners.

The document part number and revision appear at the bottom of each page. The revision identifies the document's update level. The revision history of this document is summarized in the following table.

Revision	Date	Comments
A	1 November 2012	Initial Release
B	10 December 2012	Initial release for customer

To access other GE Healthcare manuals, go to the Common Documentation Library (CDL), located at http://www.gehealthcare.com/usen/service/biomed_tech_selfservice/services_user_doc/products/support.html and click **Cardiology**.

To access Original Equipment Manufacturer (OEM) manuals, go to the device manufacturer's Web site.

Service Manual Language Information

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1

Introduction

This chapter provides general information required for the proper use of the system and this manual. Familiarize yourself with this information before using the system.

Manual Information

This section provides information for the correct use of this manual.

Keep this manual with the equipment at all times and periodically review it. You should request training assistance from GE Healthcare, if needed.

Intended Audience

This manual is intended for use by anyone responsible for integrating the Optima EMS system (also referred to as "the system" throughout this manual) and the HIS/EMR system.

Manual Purpose

This manual details the description of the data flow between the system and HIS/EMR. It describes how to integrate with these systems.

Document Conventions

This manual uses the following conventions.

Typographical Conventions

Convention	Description
Bold Text	Indicates keys on the keyboard, text to enter, or hardware items such as buttons or switches on the equipment.
<i>Italicized-Bold</i> Text	Indicates software terms that identify menu items, buttons or options in various windows.
CTRL+ESC	Indicates a keyboard operation. A plus (+) sign between the names of two keys indicates that while holding the first key, you should press and release the second key. For example, Press CTRL+ESC means to press and hold the CTRL key and then press and release the ESC key.

Convention	Description
<space>	Indicates that you must press the spacebar. When instructions are given for typing a precise text string with one or more spaces, the point where you must press the spacebar is indicated as <space>. This ensures that the correct number of spaces is inserted in the correct positions within the literal text string. The purpose of the < > brackets is to distinguish the command from the literal text within the string.
Enter	Indicates that you must press the Enter or Return key on the keyboard. Do not type Enter .
>	The greater than symbol, or right angle bracket, is a concise method to indicate a sequence of menu selections. For example, the statement "From the main menu, select System > Setup > Options to open the Option Activation window" replaces the following: <ol style="list-style-type: none"> 1. From the main menu, select System to open the System menu. 2. From the System menu, select Setup to open the Setup menu. 3. From the Setup menu, select Options to open the Option Activation window.

Illustrations

All illustrations in the manual are provided as examples only. Depending on system configuration, screens in the manual may differ from the screens on your system.

All patient names and data are fictitious. Any similarity to actual persons is coincidental.

Notes

Notes provide application tips or additional information that, while useful, are not essential to the correct operation of the system. They are called out from the body text through a flag word and indentation, as follows:

NOTE:

The tip or additional information is indented below the **NOTE** flag word.

Additional Assistance

GE Healthcare maintains a trained staff of application and technical experts to answer questions and respond to issues and problems that may arise during the installation, maintenance, and use of this system.


Contact your local GE Healthcare representative to request additional assistance.

2

HL7 Inbound

The system supports HL7 interfaces to accept ADT, OMG, and Order messages from the Hospital Information System and Electronic Medical Record (HIS/EMR). You must enable the appropriate HIS options (ADT, Orders) in the system.

The following diagram depicts a typical production deployment environment for HIS/EMR and the system.

The system receives the inbound ADT, OMG, and Order (ORM) HL7 messages from external systems (HIS/EMR) and routes them to the HL7 service entry of a special site.

The data transfer from HIS/RMR to the system is done using TCP/IP socket-based communication protocol. A pre-defined HL7 message format is used to route data.

The following HL7 event types are supported by the system HL7 Interface for HL7 messages

Supported HL7 event	Sub-event	Supported HL7 version	Supported (Optima EMS) workflow function	Definition
ADT	A01	2.3, 2.3.1 2.4 2.5 2.5.1	Admit patient	Patient is admitted.
ADT	A08	2.3, 2.3.1 2.4 2.5 2.5.1	Update patient information	Patient data has changed.

ADT	A39	2.3, 2.3.1 2.4 2.5 2.5.1	Merge patient	Patient is merged.
ADT	A40	2.3, 2.3.1 2.4 2.5 2.5.1	Merge patient	Patient is merged.
ORM	O01	2.3, 2.3.1 2.4 2.5 2.5.1	Create order	An order is created.
OMG	O19	2.3, 2.3.1 2.4 2.5 2.5.1	Create order	An order is created.

All other HL7 messages and versions are unsupported and will be discarded by the HL7 interface. When an unsupported HL7 message or version is sent to the system HL7 interface, an unsupported message is logged.

Mapping HL7 Messages to the Optimal EMS HL7 Interface


Each centricity hospital in Optimal EMS will configure the HL7 interface to meet the needs of the HIS/EMR systems.

The system provides a configuration file to meet the HL7 interface route and parse HL7 messages. The configuration file name is **hl7_mapping_config.xml** located at the default folder at **C:\Program Files\GE HealthCare\OptimaEMSService**.

Site HL7 configuration

Each centricity hospital has a configuration file with five attributes to configure as indicated in the following table.

Attribute name	Description
Name	The name of the centricity hospital.
ID¹	The site ID of the centricity hospital in the system..
ECG Keyword	The order type description in the HL7 message. If HL7 order (ORM) messages contain this word, The system will mark this order as an ECG order.
Stress Keyword	Similar to the ECG Keyword , the system will mark this order as a stress order.

Holter Keyword	Similar to the ECG Keyword , the system will mark this order as a Holter order.
<p>¹ Obtain the configuration file ID value from the system as follows:</p> <ol style="list-style-type: none"> 1. Open the Site Information window by clicking System Configuration > Site Configuration. 2. The Site ID is listed in the Site ID field of Basic Information. 	

Message Handler Component Interfaces

The system HL7 Interface consists of three message handler component interfaces:

- **ADT**
- **ORM**
- **OMG**

Each HL7 message handler has an associated configuration in the configuration file. It gets configured according to various HIS.

Each HL7 message handler configuration is composed of one or two group sub-sections, and each group has multiple mapping parts. A mapping parts has three attributes.

Key Attribute

The **key** is universally pre-defined and cannot be changed. The keys are mapped are as indicated in the following table.

key definition

Key	Mapped
100	PID of patient
101	Family Name of patient
102	Given Name of patient
103	Gender of patient

key definition (cont'd.)

104	Birth Date of patient
107	Height of patient
108	Weight of patient
200	Admit Date/Time of patient
201	PID of patient visit
202	Room of patient visit
203	Bed of patient visit
204	Diagnosis of patient visit
205	Visit ID of patient visit
300	Order Date of patient order
301	PID of patient order
302	Type of patient order
303	Test Reason of patient order

Path Attribute

The **Path** attribute determines how to get the data item from HL7 messages. It has different values for various HL7 messages.

Explanation Attribute

The **Explanation** attribute contains a description of the mapped parts and can be empty.

A Sample of the Configuration File

A sample of the configuration file is indicated as follows:

The ADT_A01 message will be used as an example. ADT_A01 has a corresponding interface, and has two parameters, **Patient** and **Visit**, each is a **Group** item. Each group indicates what data is needed (the **key** attribute) and where to get it (the **path** attribute).

For the date time type mapping element, a format attribute is needed because different systems may employ different date time formats.

The **explanation** attribute exists only for the user's convenience.

```
<?xml version="1.0" encoding="utf-8" ?>
<HL7Mapping site="gumei" id="1075" ecgkeyword="ecg" stresskeyword="stress"
holterkeyword="holter">
  <Message type="ADT_A01" version="2.3">
    <Group name="Patient">
```

```

<mapping key="100" path=".PID(0)-2" explanation="patient id" />
<mapping key="101" path=".PID(0)-5-1" explanation="family name" />
<mapping key="102" path=".PID(0)-5-2" explanation="given name" />
<mapping key="103" path=".PID(0)-8" explanation="gender" />
<mapping key="104" path=".PID(0)-7" explanation="birthday" format="YYYY-MM-DD"
/>

</Group>
<Group name="Visit">
<mapping key="100" path=".PID(0)-2" explanation="patient id" />
<mapping key="101" path=".PID(0)-5-1" explanation="family name" />
<mapping key="102" path=".PID(0)-5-2" explanation="given name" />
<mapping key="103" path=".PID(0)-8" explanation="gender" />
<mapping key="104" path=".PID(0)-7" explanation="birthday" format="YYYY-MM-DD"
/>
</Group>
</Message>
</HL7Mapping>

```

ADT Message Handler Configuration

A01/A08

A01/A08 has two groups that need to be configured: **patient** and **visit**.

See the following example:

```

<Message type ="ADT_A01" version ="2.3">
  <HL7Mapping site="gumei" id="1075" ecgkeyword="ecg" stresskeyword="stress"
holterkeyword="holter">
    <Group name ="Patient">
      <mapping key="100" path=".PID(0)-2" explanation="patient id" />
      <mapping key="101" path=".PID(0)-5-1" explanation="family name" />
      <mapping key="102" path=".PID(0)-5-2" explanation="given name" />
      <mapping key="103" path=".PID(0)-8" explanation="gender" />
      <mapping key="104" path=".PID(0)-7" explanation="birthday" format="YYYY-MM-DD"
/>
    </Group>
    <Group name="Visit">

```

```

<mapping key="100" path=".PID(0)-2" explanation="patient id" />
<mapping key="101" path=".PID(0)-5-1" explanation="family name" />
<mapping key="102" path=".PID(0)-5-2" explanation="given name" />
<mapping key="103" path=".PID(0)-8" explanation="gender" />
<mapping key="104" path=".PID(0)-7" explanation="birthday" format="YYYY-MM-DD"
/>

</Group>
</Message>

```

A39/A40

A39/A40 has two groups that need to be configured, ***new patient*** and ***old patient***.

See the following example:

```

<Message type="ADT_A39" version="2.3">
  <Groups>
    <Group name="OldPatientID">
      <Mappings>
        <Mapping key ="100" path=".MRG(0)-4" rep ="false" explanation="Old patient ID" />
      </Mappings>
    </Group>
    <Group name="NewPatientID">
      <Mappings>
        <Mapping key ="100" path=".PID(0)-2" rep ="false" explanation="New patient ID" />
      </Mappings>
    </Group>
  </Groups>
</Message>

```

ORM Message Handler Configuration

The ORM message handler has two groups that need to be configured: ***order*** and ***visit***.

See the following example:

```

<Message type="ORM_001" version="2.3">
  <Groups>
    <Group name="Order">
      <Mappings>

```

```

<Mapping key="300" path=".ORC(0)-9" rep ="false" explanation="order date" />
<Mapping key="301" path=".PID(0)-3" rep ="false" explanation="patient id" />
<Mapping key="302" path=".OBR(0)-4-2" rep ="false" explanation="device type" />
<Mapping key="303" path=".OBR(0)-31-2" rep ="false" explanation="test reason" />
</Mappings>
</Group>
<Group name="Visit">
<Mappings>
<Mapping key="200" path=".MSH(0)-7" rep ="false" explanation="visit date time" />
<Mapping key="201" path=".PID(0)-3-1" rep ="false" explanation="patient id" />
<Mapping key="202" path=".PV1(0)-11-2" rep ="false" explanation="room" />
<Mapping key="203" path=".PV1(0)-11-3" rep ="false" explanation="bed" />
<Mapping key="204" path=".DG1(0)-4" rep ="false" explanation="diagnosis"/>
<Mapping key="205" path=".PV1(0)-19-1" rep ="false" explanation="visit id"/>
</Mappings>
</Group>
</Groups>
</Message>

```

OMG Message Handler Configuration

OMG message handler has one group to configure, ***order***.

See the following example:


```

<Message type="OMG_O19" version="2.3">
<Groups>
<Group name="Order">
<Mappings>
<Mapping key="300" path=".ORC(0)-9" rep ="false" explanation="order date" />
<Mapping key="301" path=".PID(0)-3" rep ="false" explanation="patient id" />
<Mapping key="302" path=".OBR(0)-4-2" rep ="false" explanation="device type" />
<Mapping key="303" path=".OBR(0)-31-2" rep ="false" explanation="test reason" />
</Mappings>
</Group>
</Groups>
</Message>

```

Enabling the system HL7 Inbound

Open the **Site Configuration** window by clicking **System Configuration > Site Configuration**.


HL7 Settings are enabled by selecting the following check boxes:

- **Allow HL7 ADT**
- **Allow HL7 Order**

Configuring the System HL7 Service Entry

1. Open the **HL7 Service Entry** window by clicking **System Configuration > HL7 Service Entry**.


- Each site should use an HL7 service entry so that it is mapped.
 - Make sure the port has not been occupied by other program and software. The port is different between each site.
2. In the **HL7 Service Entry** window, click **Start** to start the service entry.


NOTE:

Make sure the status of the service entry is **Running**.

- You can change the configuration of the service entry by clicking **Modify** (Port, remove service entry, and so on).

- If any change has been made to the service entry configuration, you need to restart the service entry by clicking **Restart**.

After the above configurations have been set, the system will accept HL7 ADT and Order messages from the HIS. You need to restart the system service to apply these changes.

3


Non-HL7 Inbound

The system supports non-HL7 interfaces to import patient and order data from the HIS. It provides a web service to import patients and orders for third-party applications. The default service address is <http://Optima EMS Server IP address:8999/MacItWebServiceInBound>.

The following web service functions are provided by the system non-HL7 Interface.

Function Name	Function Prototype
<i>IsPatientExisted</i>	IsPatientExisted (string patientID, string hisID)
<i>AdmitPatient</i>	AdmitPatient (Patient patient, PatientVisit visit, string HisID)
<i>UpdatePatient</i>	UpdatePatient (Patient patient, PatientVisit visit, string HisID)
<i>MergePatient</i>	MergePatient (string oldPatientIDList, string newPatientID, string HisID)
<i>CreateOrder</i>	CreateOrder (Order order, string HisID)

See the following data flow diagram.

IsPatientExisted


This function is used to check whether a patient already exists in the system.

IsPatientExisted has two parameters.

- *PatientID*: the patient identify ID in HIS system
- *HisID*: identifies which site the patient belongs to

The HISID should be same as in the system and cannot be empty. Obtain the value by performing the following procedure:

1. Log into the Optimal EMS system.
2. Open the **Site Configuration** window by clicking **System Configuration > Site Configuration**.
3. At the **Site Configuration** window, set the **HISID** value shown in the left window to the same value as the **HIS Name**.

AdmitPatient

This function is used to create a new patient in the system. It has three parameters.

AdmitPatient has follow sub-data items.

Name	Enable empty value	Description
SocialInsuranceNum	No	Patient social security number. If the patient does not have a social security number, keep this field blank.
SocialInsuranceType	No	Patient social security type. Different areas of China may have different types of social security. If a social security type is not present, keep this field blank.
PatientID	No	Patient ID in the HIS system.
FamilyName	No	Family name of the Patient
GivenName	No	Given name of the Patient

<i>BirthDate</i>	Yes	Date of birth of the Patient
<i>Gender</i>	No	Gender of patient. The value can be one of following: Male, Female, Unknown.
<i>Age</i>	Yes	Age of the Patient
<i>AgeUnit</i>	No	Patient age field. The value can be one of following: Year, Month, Week, Day, Hour, Unknown.
<i>HeightInCm</i>	Yes	Height of the Patient
<i>WeightInKg</i>	Yes	Weight of the Patient

PatientVisit: patient visit information from the HIS has the following sub-data items:

Name	Enable empty value	Description
<i>HisVisitID</i>	No	Patient visit ID in the HIS. It is used to identify multiple patient visits.
<i>PatientID</i>	No	Patient ID in the HIS system.
<i>Room</i>	Yes	Room number of the Patient
<i>Bed</i>	Yes	Bed number of the Patient
<i>Diagnosis</i>	Yes	Diagnosis result of the Patient
<i>AdmitDatetime</i>	No	Patient visit time.

HisID: same as the HisID parameter of “[IsPatientExisted](#)” on page 23.

UpdatePatient

This function is used to update patient information and create a new patient visit information in the system.

UpdatePatient has three parameters.

- **Patient:** has the same patient parameters as the **AdmitPatient** function.
- **PatientVisit:** has the same patient visit parameters as the **AdmitPatient** function
- **HisID:** same as the patient visit parameters of “[IsPatientExisted](#)” on page 23.

MergePatient

This function is used to merge multiple patients into one patient in the system.

MergePatient has three parameters.

- **oldPatientIDList**: the patient ID in the HIS system which will be merged in the system.
- **newPatientID**: the patient ID in the HIS system which will replace other patients in the system.
- **HisID**: same as the HisID parameters of “[IsPatientExisted](#)” on page 23.

CreateOrder

This function is used to create patient order information in the system.

CreateOrder has two parameters.

Order: patient order information from the HIS with the following sub-data items:

Name	Enable empty value	Note
PatientID	No	Patient ID in the HIS system.
HisOrderID	No	Patient order ID in the HIS system.
HisVisitID	No	Patient visit ID in the HIS system.
OrderDate	No	Patient order date in the HIS system.
TestReason	Yes	Reason for the patient test
DeviceType	No	Type of patient order. The value can be one of following: ECG , Stress , Holter , and Unknown .

HisID: same as the HisID parameters of “[IsPatientExisted](#)” on page 23.

4

HIS/EMR Outbound

The system supports exporting of billing and diagnosis result report information to a shared folder. A HIS or EMR system should access this shared folder to get outbound data (billing, diagnosis result report). There are no real-time HL7 outbound messages which communicate between the system and HIS/EMR. The appropriate HIS option (HIS output) should be enabled in the system, to automatically export billing and test reports of a confirmed test. HIS/EMR could access the shared folder on the system server over the network.

Activating the System HIS/EMR Outbound

Enabling the System HIS Outbound


1. Open the **Site Configuration** window by clicking **System Configuration > Site Configuration**.
2. Select the **Allow HL7 Outbound** check box.

The HIS outbound is now enabled.

Configuring the System HIS Outbound Configurations

1. At the **HIS Name** field, type the HIS name, which is the same name as the folder name where all test report are exported to for that particular site. This field should not be left empty.
2. Select the **Default Report Template** by choosing from the list in the drop-down menu.

In the following figure, the **Test report template** was selected.


Configuring the System HIS Outbound Shared Folder

This shared folder is on the system server. The default path is: **C:\Program Files\GE HealthCare\OptimaEMSService\ECG Data\MACIT outbound\ + HisName.**

- Set the permissions of the following folder so that it can be accessed by the HIS/EMR system:

C:\Program Files\GE HealthCare\OptimaEMSService\ECG Data\MACIT outbound\ + HisName.

The system will now automatically export billing and test reports to a shared folder when a user confirms a test.

NOTE:

The system will export billing of an order as a file with the .txt extension. The content is a HL7 DFT_P03 message.

The system will export diagnosis result report of an order as a file with a .pdf extension.

The export file name rule is **patientOrderDatetime_PID_Report.pdf** and **patientOrderDatetime_PID_Billing.txt**

The system supports customized HIS/EMR outbound file name rules when needed.


Customizing HIS/EMR Outbound File Name Rules

Creating A C# Class Library Project


- From your desktop, click **Start > Programs > Microsoft Visual Studio 2010 > Microsoft Visual Studio 2010**.

NOTE:


The **Microsoft Visual Studio 2010** is a software development tool which needs to be installed on your local machine.

The *Microsoft Visual Studio* window opens.


2. Click *File > New > Project....*

The *New Project* window opens.

3. Click the **Class Library** in the **New Project** window.


4. Enter the information displayed at the bottom of the **New Project** window.

- **Name:** the project name
For example, **CustomizeOutBoundFile**.
- **Location:** the project location
It is a local storage file folder path which save this project content.
For example, **D:\CustomizeOutBoundFile folder**.

5. Click **OK**.


A new project is created successfully and will be opened automatically (**CustomizeOutBoundFile** window as example).

Adding MacItHooks.dll Reference

A reference needs to be added to the **MacItHooks.dll** file in a new project.

- Right-click the **Reference** in the **Solution Explorer** panel of the **CustomizeOutBoundFile** window.


The **Add Reference** window opens.

- Click **Browse** to choose the **MacItHooks.dll** from the local file.

NOTE:

The **MacItHooks.dll** is located at the system server installation folder by default: **C:\Program Files\GE HealthCare\OptimaEMSService**.


- Click **OK**.

The reference is added to the **MacItHooks.dll** file successfully.


Adding A New Class Named *CustomizeOutBoundFiler*

1. Right-click the *CustomizeOutBoundFile* in the *Solution Explorer* panel of the *CustomizeOutBoundFile* window.
2. Select *Add > Class....*

The *Add New Item* window opens.

3. Select **Class** in the **Add New Item** window, and enter the **Name** displayed at the bottom of the **Add New Item** window.

4. Click **Add**.
the **Source Code Edit** window will open automatically.

```

1 //using System;
2 //using System.Collections.Generic;
3 //using System.Linq;
4 //using System.Text;
5
6 namespace CustomizeOutBoundFile
7 {
8     public class CustomizeOutBoundFile : GEHC.MAIC.Hooks.IHisOutboundFileNameStrategy
9     {
10        /// <summary>
11        /// Gets the name of his outbound file.
12        /// </summary>
13        /// <param name="siteId">The site id.</param>
14        /// <param name="orderId">The order id.</param>
15        /// <param name="testId">The test id.</param>
16        /// <param name="patientId">The patient id.</param>
17        /// <param name="hisOrderId">His order id.</param>
18        /// <param name="testDate">The test date.</param>
19        /// <returns></returns>
20        public string GetHisOutboundFileName(int siteId, int orderId, int testId, string patientId, string hisOrderId, DateTime testDate)
21        {
22            // add out bound file name code here
23            return string.Empty;
24        }
25    }
26}

```

Implementing IHisOutboundFileNameStrategy Interface in CustomizeOutBoundFiler Class

In the **Source Code Edit** window, you can write your source code to implement the customize function.


1. Mark class *CustomizeOutBoundFiler* as the public prefix.
 2. Inherit the *IHisOutboundFileNameStrategy* interface.
 3. Implement the *GetHisOutboundFileName* function code customized.

The screenshot shows the Microsoft Visual Studio interface with the title bar "CustomizeOutBoundFile - Microsoft Visual Studio". The menu bar includes File, Edit, View, NCrunch, Theme, Project, Build, Debug, Team, ClearCase, Data, Tools, Architecture, Test, ReSharper, Analyze, Window, Help. The toolbar has icons for file operations like Open, Save, Print, and a search bar. The Solution Explorer on the right shows a single solution named "CustomizeOutBoundFile" with one project and one file: "CustomizeOutBoundFile.cs". The code editor window displays the following C# code:

```
1 using System;
2 using System.Collections.Generic;
3 using System.Linq;
4 using System.Text;
5
6 namespace CustomizeOutBoundFile
7 {
8     public class CustomizeOutBoundFile : GEHC.MACIT.Hooks.IHisOutboundFileNameStrategy
9     {
10         ///<summary>
11         /// Gets the name of his outbound file.
12         ///</summary>
13         ///<param name="siteId">The site id.</param>
14         ///<param name="orderId">The order id.</param>
15         ///<param name="testId">The test id.</param>
16         ///<param name="patientId">The patient id.</param>
17         ///<param name="hisOrderId">His order id.</param>
18         ///<param name="testDate">The test date.</param>
19         ///<returns></returns>
20         public string GetHisOutboundFileName(int siteId, int orderId, int testId, string patientId, string hisOrderId, DateTime testDate)
21         {
22             // add out bound file name code here
23             return string.Empty;
24         }
25     }
26 }
27
```


Building C# dll Library Project

- Select **Release** in the **Build Option** field.


- Build the project by clicking **Build > Build Solution**.

This will generate the **Build Log** in the **Output** window.


- Get the **Project Build Result** file.

The **Project Build Result** file is located in **bin\Release** folder under the project location you created in Step 4 of "[Creating A C# Class Library Project](#)" on page [30](#). And the file name is **CustomizeOutBoundFile.dll**.

Releasing Assembly dll File to the System Server

Copy the **Project Build Result** file you got in Step 3 of “Building C# dll Library Project” on page 37 to the Plugin folder under the system server installation folder.

GE Medical Systems
Information Technologies, Inc.
8200 West Tower Avenue
Milwaukee, WI 53223 USA
Tel: +1 414 355 5000
+1 800 558 7044 (US Only)
Fax: +1 414 355 3790

Asia Headquarters

GE Medical Systems
Information Technologies, Inc.
Asia; GE (China) Co., Ltd.
1 Huatuo Road
Zhangjiang Hi-tech Park Pudong
Shanghai, People's Republic of China 201203
Tel: +86 21 3877 7888
Fax: +86 21 3877 7451

GE Medical Systems *Information Technologies, Inc.*, a General Electric Company, going to market as GE Healthcare.

www.gehealthcare.com