stryker*

Operations Manual



Symbols

[ji]	Operating instructions/Consult instructions for use
\wedge	General warning
\triangle	Caution
REF	Catalogue number
SN	Serial number
US Patents	For US Patents see www.stryker.com/patents
***	Manufacturer
<u>^</u>	Safe working load
4	Dangerous voltage
~	Alternating current
===	Direct current
A	Unit provides terminal for connection of a potential equalization conductor. The potential equalization conductor provides direct connection between the unit and potential equalization busbar of the electrical installation.
	Protective earth ground
IPX4	Protection from liquid splash
†	Type B applied part
CUL US	Medical Equipment Classified by Underwriters Laboratories Inc. With Respect to Electric Shock, Fire, and Mechanical Hazards Only in Accordance with ANSI/AAMI ES60601-1: 2005 and CAN/CSA-C22.2 No. 60601-1:08.
X	In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment, this symbol indicates that the product must not be disposed of as unsorted municipal waste, but should be collected separately. Contact your local distributor for disposal information.
(C ₂))	Non-ionizing radiation

www.stryker.com 3005-109-001 REV B

Symbols

	Caution; electrostatic sensitive
	iBed Locator is connected
X	iBed Locator is not connected
(('T'))	Network is connected
1×	Network is not connected
F©	This device complies with Part 15 of the FCC rules

Table of Contents

Summary of safety precautions 4 Introduction 6 Product description 6 Indications for use 6 Indications for use for //Bed Wireless 6 Indications for use for //Bed Wireless 6 Indications for use for //Bed Wireless 6 Indications 7 Contraindications 7 Specifications 7 System requirements and recommendations for ///>///>///>///>//>//>//>//>//>//>//>//
Product description 6 Indications for use 6 Indications for use for /Bed Wireless 6 Expected service life .7 Contraindications .7 Specifications .7 Specifications .7 System requirements and recommendations for /Bed Wireless (option) .8 // Bed Server requirements for /Bed Wireless (option) .8 // Stryker /Bed Wireless Client radio specifications .9 Client device data usage .10 Customer network communication requirements for /Bed Wireless option .11 Product illustration .13 Contact information .13 Serial number location .13 Date of manufacture .14 Setup .15 Equipping the optional nurse call communication .15 Setup .15 Equipping the optional /Bed Wireless .16 Operation .17 Applying or releasing the brakes .17 Transporting the product with steer lock .17 Activating the CPR release
Indications for use 6 Indications for use for /Bed Wireless 6 Expected service life
Indications for use for iBed Wireless 6 Expected service life .7 Contraindications .7 Specifications .7 System requirements and recommendations for iBed Wireless (option) .8 iBed Server requirements for iBed Wireless (option) .8 Stryker iBed Wireless Client radio specifications .9 Client device data usage .10 Customer network communication requirements for iBed Wireless option .11 Product illustration .13 Contact information .13 Serial number location .14 Date of manufacture .14 Setup .15 Equipping the optional nurse call communication .15 Setting up the optional iBed Wireless .16 Operation .17 Applying or releasing the brakes .17 Transporting the product with steer lock .17 Activating the CPR release .18 Raising the lower leg section .20 Attaching a fracture frame .20 Securing a Foley bag to the Foley bag hooks .20 Installing the patient restraint strap tie-ins <t< td=""></t<>
Expected service life.
Contraindications 7 Specifications 7 System requirements and recommendations for iBed Wireless (option) 8 iBed Server requirements for iBed Wireless (option) 8 Stryker iBed Wireless Client radio specifications 9 Client device data usage 10 Customer network communication requirements for iBed Wireless option 11 Product illustration 13 Contact information 13 Serial number location 14 Date of manufacture 14 Setup 15 Equipping the optional nurse call communication 15 Setting up the optional iBed Wireless 16 Operation 17 Applying or releasing the brakes 17 Transporting the product with steer lock 17 Activating the CPR release 18 Raising the lower leg section 20 Attaching a fracture frame 20 Securing a Foley bag to the Foley bag hooks 20 Installing the patient restraint strap tie-ins 21 Raising the siderails 22
Specifications. 7 System requirements and recommendations for /Bed Wireless (option) 8 //Bed Server requirements for /Bed Wireless (option) 8 Stryker /Bed Wireless Client radio specifications 9 Client device data usage 10 Customer network communication requirements for /Bed Wireless option 11 Product illustration 13 Contact information 13 Serial number location 14 Date of manufacture 14 Setup 15 Equipping the optional nurse call communication 15 Setting up the optional /Bed Wireless 16 Operation 17 Applying or releasing the brakes 17 Transporting the product with steer lock 17 Activating the CPR release 18 Raising the lower leg section 20 Attaching a fracture frame 20 Securing a Foley bag to the Foley bag hooks 20 Installing the patient restraint strap tie-ins 21 Raising the siderails 22
System requirements and recommendations for /Bed Wireless (option) 8 //Bed Server requirements for /Bed Wireless (option) 8 Stryker /Bed Wireless Client radio specifications 9 Client device data usage 10 Customer network communication requirements for /Bed Wireless option 11 Product illustration 13 Contact information 13 Serial number location 14 Date of manufacture 14 Setup 15 Equipping the optional nurse call communication 15 Setting up the optional /Bed Wireless 16 Operation 17 Applying or releasing the brakes 17 Transporting the product with steer lock 17 Activating the CPR release 18 Raising the lower leg section 20 Attaching a fracture frame 20 Securing a Foley bag to the Foley bag hooks 20 Installing the patient restraint strap tie-ins 21 Raising the siderails 22
############
Stryker /Bed Wireless Client radio specifications 9 Client device data usage 10 Customer network communication requirements for /Bed Wireless option 11 Product illustration 13 Contact information 13 Serial number location 14 Date of manufacture 14 Setup 15 Equipping the optional nurse call communication 15 Setting up the optional /Bed Wireless 16 Operation 17 Applying or releasing the brakes 17 Transporting the product with steer lock 17 Activating the CPR release 18 Raising the lower leg section 20 Attaching a fracture frame 20 Securing a Foley bag to the Foley bag hooks 20 Installing the patient restraint strap tie-ins 21 Raising the siderails 22
Client device data usage 10 Customer network communication requirements for iBed Wireless option 11 Product illustration 13 Contact information 13 Serial number location 14 Date of manufacture 14 Setup 15 Equipping the optional nurse call communication 15 Setting up the optional iBed Wireless 16 Operation 17 Applying or releasing the brakes 17 Transporting the product with steer lock 17 Activating the CPR release 18 Raising the lower leg section 19 Lowering the lower leg section 20 Attaching a fracture frame 20 Securing a Foley bag to the Foley bag hooks 20 Installing the patient restraint strap tie-ins 21 Raising the siderails 22
Customer network communication requirements for iBed Wireless option 11 Product illustration 13 Contact information 13 Serial number location 14 Date of manufacture 14 Setup 15 Equipping the optional nurse call communication 15 Setting up the optional iBed Wireless 16 Operation 17 Applying or releasing the brakes 17 Transporting the product with steer lock 17 Activating the CPR release 18 Raising the lower leg section 20 Attaching a fracture frame 20 Securing a Foley bag to the Foley bag hooks 20 Installing the patient restraint strap tie-ins 21 Raising the siderails 22
Product illustration 13 Contact information 13 Serial number location 14 Date of manufacture 14 Setup 15 Equipping the optional nurse call communication 15 Setting up the optional /Bed Wireless 16 Operation 17 Applying or releasing the brakes 17 Transporting the product with steer lock 17 Activating the CPR release 18 Raising the lower leg section 19 Lowering the lower leg section 20 Attaching a fracture frame 20 Securing a Foley bag to the Foley bag hooks 20 Installing the patient restraint strap tie-ins 21 Raising the siderails 22
Contact information 13 Serial number location 14 Date of manufacture 14 Setup 15 Equipping the optional nurse call communication 15 Setting up the optional iBed Wireless 16 Operation 17 Applying or releasing the brakes 17 Transporting the product with steer lock 17 Activating the CPR release 18 Raising the lower leg section 19 Lowering the lower leg section 20 Attaching a fracture frame 20 Securing a Foley bag to the Foley bag hooks 20 Installing the patient restraint strap tie-ins 21 Raising the siderails 22
Serial number location
Date of manufacture
Setup
Equipping the optional nurse call communication
Setting up the optional iBed Wireless
Operation17Applying or releasing the brakes17Transporting the product with steer lock17Activating the CPR release18Raising the lower leg section19Lowering the lower leg section20Attaching a fracture frame20Securing a Foley bag to the Foley bag hooks20Installing the patient restraint strap tie-ins21Raising the siderails22
Applying or releasing the brakes
Transporting the product with steer lock
Activating the CPR release
Raising the lower leg section
Lowering the lower leg section
Attaching a fracture frame
Securing a Foley bag to the Foley bag hooks
Installing the patient restraint strap tie-ins
Raising the siderails
•
Lowering the siderails
Positioning the optional two-stage permanently attached IV pole
Positioning the optional removable IV pole
Illuminating the room with the night light
Activating nurse call communication (option)
Replacing the nurse call backup battery (option)
Connecting peripheral equipment to the built-in 110 volt auxiliary power outlet (option)
Operator control panel (outside siderail)
Patient control panel (inside siderail)31
Optional Smart TV control panel (inside siderail)
Footboard control panel - Bed controls
Footboard control panel - Lockouts

Table of Contents

Footboard control panel - Chaperone bed exit (option)	36
Arming or disarming Chaperone bed exit (option)	36
Footboard control panel - Chaperone bed exit with zone control (option)	37
Arming or disarming Chaperone bed exit with zone control (option)	37
Footboard LED indicators	38
Footboard control panel - Scale	40
Weighing a patient	40
Setting the scale to zero	41
Menu display	42
Accessing functions and features with the menu display	43
Menu display with iBed Wireless (option)	44
Viewing the weight log	45
Measuring weight gain or loss	45
Changing equipment	46
Changing the patient weight	47
Changing the scale units	48
Changing the backlight intensity	49
Setting the alarm tones	50
Setting the brake alarm	50
Setting an audible iBed Awareness alarm	50
Setting the iBed Awareness nurse call alarm	51
iBed Awareness (option)	51
Configuring iBed Awareness	52
Acknowledging iBed Awareness status alerts	52
Motion pendant with nurse call (option)	55
Motion and communication pendant (option)	56
Infrared (IR) module (option)	57
iBed Locator (option)	57
Accessories	59
Cleaning	60
Disinfecting	61
Preventive maintenance	62
FCC notification	63
EMC information	64
Warranty	68
Warranty exclusion and damage limitations	68
To obtain parts and service	68
Return authorization	68
Damaged product	68
International warranty clause	68

Warning/Caution/Note Definition

The words WARNING, CAUTION, and NOTE carry special meanings and should be carefully reviewed.



MARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.



CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or

Note: Provides special information to make maintenance easier or important instructions clearer.

Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

★ WARNING

- Always plug the product directly into a grounded, hospital grade wall outlet. You can only achieve grounding reliability when you use a hospital grade wall outlet. This product is equipped with a hospital-grade plug for protection against electric shock hazard.
- Always use a Stryker supplied interface cable. Use of any other cable may cause the bed to function improperly, which may result in patient or user injury.
- Always associate or map the iBed Locator to the room or location to provide accurate location information. Failure to properly associate or map the iBed Locator to the room or location will yield incorrect remote information. If you move an iBed Locator after it has been installed and mapped, you must remap to the new room or location. If the room or location information is changed after initial installation, you must remap the iBed Locator.
- Always apply the brakes when a patient is getting into or out of the product to avoid instability.
- Always apply the brakes when the patient is unattended.
- Do not apply the brakes to slow or stop the product while it is in motion.
- Always lock the siderails in the full up position with the sleep surface horizontal in the lowest position when you transport the product with a patient.
- Always unplug the power cord before you transport the product.
- Always release the brakes before you transport the product. Do not transport the product with the brakes applied.
- Do not transport the product laterally after you apply the steer lock pedal. The product cannot swivel when you transport with steer lock.
- Use only retractable traction or fracture frames. Failure to use a retractable frame may result in injury to the patient and/or damage to the equipment.
- Always set the siderail position to make sure that a patient is safely in the product.
- Always lock the operator control panel and patient control panel when the patient is unattended.
- Always keep the siderails outside of the oxygen tent.
- Only use hospital-grade electric equipment consuming 10A or less with the auxiliary power outlet (option). The use of standard electric equipment may bring the current leakage to a level unacceptable for hospital equipment.
- Do not use the 110V auxiliary power outlet (option) for life sustaining equipment.
- Always lock the control panel when you leave the patient unattended.
- Always lock the patient control panel when the patient's condition requires extra safety measures.
- Do not use bed exit to replace patient monitoring protocol.
- Bed exit is intended only to aid in the detection of a patient exiting the product.
- Do not use bed exit with patients who weigh less than 50 lb (23 kg).
- Do not use the scale system reading as a reference for medical treatment.
- The scale system assists only in the monitoring of the patient's weight variation.
- Do not set the scale to zero or weigh the patient when a support surface therapy is active. Motion from the support surface functions may adversely affect the scale system performance.
- Do not use iBed Awareness to replace your patient monitoring protocol.
- Do not use iBed Awareness as a lock indicator for siderails. iBed Awareness only detects the position of the siderails.

Summary of safety precautions

⚠ CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
- Do not modify the product or any components of the product. Modifying the product can cause unpredictable
 operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
- Always make sure that all persons and equipment are away from the area below and around the Fowler before you
 activate the CPR release. The CPR release is for emergency use only.
- Always determine the proper use of the restraint straps and restraint strap locations. Improperly adjusted restraint straps can cause serious injury to a patient. Stryker is not responsible for the type or use of restraint straps on any Stryker products.
- Do not load the IV pole above the safe working load of 40 lb (18 kg).
- Always raise the siderails when the litter is in its full down position. This prevents the scale system from weighing a
 patient inaccurately.
- Always make sure that the siderails are locked before you arm iBed Awareness.
- · Make sure that you set the desired product parameters before arming iBed Awareness.
- Do not use accessories that cover the footboard LED light bars.
- Do not clean, disinfect, service, or perform maintenance while the product is in use.
- Always unplug the power cord from the wall outlet when large spills occur near the circuit boards, cables, and
 motors. Remove the occupant from the product, clean up the fluid, and inspect the product. Fluids can cause
 unpredictable operation and decreased functionality of any electrical product. Do not return the product to service
 until it is completely dry and you have tested for safe operation.
- Always clean Velcro® after each use. Saturate Velcro with disinfectant and allow disinfectant to evaporate.
 Appropriate disinfectant for nylon Velcro should be determined by the hospital.
- Always wipe down with clean water (or 70% isopropyl alcohol, if using Virex® TB) and dry each product after
 disinfecting. Some disinfectants are corrosive in nature and may cause damage to the product. If you do not rinse
 and dry the product, you may leave a corrosive residue on the surface of the product. This corrosive residue could
 cause premature degradation of critical components. Failure to follow these disinfecting instructions may void your
 warranty.
- · Do not use quaternary disinfectants formulated with glycol ethers.

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

↑ CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
- Do not modify the product or any components of the product. Modifying the product can cause unpredictable
 operation resulting in injury to patient or operator. Modifying the product also voids its warranty.

Notes

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Product description

The Stryker Model 3005 **S3®** MedSurg bed has Fowler, gatch, and lift articulation capabilities that aide in the adjustment of surface contour, angle, and height. The product offers various options outlined in the product operations manual, including but not limited to the *i*Bed® Awareness system, scale, 110 VAC option, IV pole, and defibrillator tray. *i*Bed Awareness allows operators to set various bed parameters to monitor bed positioning. Alerts inform operators to patient movement within a specific zone(s) on a patient surface. The product may be equipped with an integrated scale to weigh the patient in bed. Both the *i*Bed Awareness system and Chaperone® Bed Exit system provide both visual and audible alerts.

Indications for use

The 3005 **S3** MedSurg bed is intended to support and transport patients within the Med/Surg and Critical Care hospital environments. The 3005 **S3** MedSurg bed is typically used in pre-op, post-op and recovery areas of hospital facilities. The intended user for this product is both health care providers (HCP), such as, nurses, nurses' aides, medical doctors, and human patients. Lockout features may limit patient accessible controls. Use this product with a patient sleep surface. The scale output is not intended to be used to determine diagnosis or treatment.

The intended patient population for the 3005 \$3 MedSurg bed includes:

- Patients above 50 lb with a safe working load of 500 lb
- Patients at least 2 years of age
- · Patients less than 84 in. without a bed extender or 96 in. with a bed extender

The product is not intended to support more than one individual at a time.

Indications for use for iBed Wireless

The intended use for *i*Bed Wireless (with *i*Bed Awareness system) is to assist clinical staff to monitor bed parameters on specific Stryker beds. The desired bed parameters are set by operators at the bedside. *i*Bed Wireless is only intended for use with specifically enabled Stryker beds that have been verified and validated with *i*Bed Wireless, and is not intended to provide bed status information for non-Stryker beds. *i*Bed Wireless is not intended to communicate any patient status information, nor to permanently store any type of data. *i*Bed Wireless with *i*Bed Awareness system is not intended to provide automated treatment decisions or as a substitute for professional healthcare judgment. *i*Bed Wireless with *i*Bed Awareness system is not a replacement or substitute for vital signs monitoring or alert equipment. All patient medical diagnosis and treatment are to be performed under direct supervision and oversight of an appropriate healthcare professional.

Expected service life

The 3005 **S3** MedSurg bed has a ten year expected service life under normal use conditions and with appropriate periodic maintenance.

The optional motion/nurse call pendant has a two year expected service life under normal use conditions.

The optional motion/nurse call/SmartTV pendant has a two year expected service life under normal use conditions.

The optional iBed Wireless components have a three year expected service life under normal use conditions.

Contraindications

None known.

Specifications

	Safe working load Note: Safe working load indicates the sum of the occupant and accessory weight.			227 kg
Product weight			570 lb	259 kg
Scale system capacity (optional	equipment), loads wei	ghing up to:	500 lb	227 kg
Scale system accuracy (optional equipment)			± 3 lb for patients weighing 50 to 100 lb ± 3% of the total patient weight for patients weighing 100 to 500 lb	
Overal length and width		Siderails up	95 in. x 41.5 in.	241.3 cm x 105.4 cm
Overal length and width Siderail		Siderails down	95 in. x 39.5 in.	241.3 cm x 100.3 cm
Patient sleep surface			84 in. x 35 in.	213.4 cm x 88.9 cm
Bed height to top of seat litter with 6" casters			16 in. to 30 in.	40.6 cm to 76.2 cm
Litter platform to top of siderail Head end siderail			15 in.	38.1 cm
(full up)	Foot end siderail		15.5 in.	39.37 cm
Space between siderails (full up)			2.5 in.	5.72 cm
Gatch position			0° to 45°	
Fowler position			0° to 60°	
Trendelenburg and Reverse Trendelenburg			+10° to -10° ± 1°	
Electrical requirements - all electrical requirements meet UL 60601 specifications			120 VAC, 60 Hz, 8 A	
Outlet option			110 VAC, 60 Hz, 10 A	
Duty cycle			1 minute 45 seconds ON, 30 minutes OFF	

Specifications (Continued)

Compatible mattress							
Thickness		Width		Length		ILD	
6 in.	15.2 cm	>= 35 in.	>= 88.9 cm	>= 84 in.	>= 213.4 cm	80 lb	36.3

Note: These mattress specifications comply with HBSW and IEC specifications.

Stryker reserves the right to change specifications without notice.

Specifications listed are approximate and may vary slightly from product to product or by power supply fluctuations.

Environmental conditions	Operation	Storage and transportation	
Ambient temperature	50 °F (40 °C) (10 °C)	-22 °F - (-30 °C)	
Relative humidity (non-condendisng)	30%	10% -95%	
Atmospheric pressure	700 hPa	500 hPa 500 hPa	

System requirements and recommendations for iBed Wireless (option)

To implement *i*Bed Wireless, follow these requirements for hardware, software, and communication, product specifications, required settings and recommendations.

Note: If minimum system requirements are not met, system performance will be impacted.

iBed Server requirements for iBed Wireless (option)

The following tables describes the server hardware and software requirements for *iBed* Wireless:

Hardware requirements		
Capacity	Requirement	
1-300 connected devices	2.x GHz processor or better with 4 cores Memory: 8 GB RAM Hard drive: 150 GB	
301-600 connected devices	2.x GHz processor or better with 8 cores Memory: 16 GB RAM Hard drive: 150 GB	
601-800 connected devices	2.x GHz processor or better with 16 cores Memory: 32 GB RAM Hard drive: 150 GB	

System requirements and recommendations for iBed Wireless (option) (Continued)

Hardware requirements		
Capacity	Requirement	
801-1000 connected devices	2.x GHz processor or better with 24 cores Memory: 32 GB RAM Hard drive: 150 GB	
1001-1300 connected devices	2.x GHz processor or better with 32 cores Memory: 64 GB RAM Hard drive: 150 GB	
1300+ connected devices Add additional core per 50 devices		
Two server environments are recommended for the <i>iBed</i> Wireless: TEST and PROD		
The iBed Wireless is supported in either physical or virtual environments		

Software and setup requirements		
	Requirement	
Operating system	Microsoft Windows Server 2008 R2 / 2012 R2	
Server roles	Web server (IIS) Roles services - Application development - ASP.NET - ASP	
Server features	.NET Framework 3.5.1	
Other	All current Microsoft priority updates and optional update .NET Framework 4.5	

Security requirements		
	Requirement	
Credentials	Administrator account on the server machine for installation and configuration	
Remote access	Provide remote access to the <i>i</i> Bed Server machines. For example, VPN, Citrix	

Note:

Backups of the system or a Disaster Recovery Plan is the responsibility of the customer.

Stryker iBed Wireless Client radio specifications

Manufacturer model	Silex SX-SDMAN-2830S
Chipset	AR6233X, Core chip AR6003X CSP (Qualcomm Atheros)

System requirements and recommendations for *i*Bed Wireless (option) (Continued)

IEEE 802.11	a/b/g/n
RF bands	2.4 GHz, 5 GHz
Encryption	AES and TKIP (TKIP is not supported with WPA2)
Authentication	WPA Personal / Enterprise and WPA2 Personal / Enterprise
802.1X	PEAP-MSCHAP - V2
Client certificates	Stryker iBed Wireless Client(s) cannot accept or upload certificates
Supported data rates	802.11b/g: 1-54 Mbps 802.11a: 6-54 Mbps 802.11n: MCS0-6
Channel plan	2.4 GHz: All Channels Supported 5 GHz: All Channels Supported (Recommend against using DFS and ISM Channels)
Other	Stryker <i>i</i> Bed Wireless Client is able to connect to an existing SSID

Client device data usage

- The client uses 10-15 KB per connected device every 40 seconds.
- The client uses an additional 5-21 KB per device for each subscription that is created by a third-party vendor like Connexall, Capsule, Epic, and Cerner.

Note:

Based on network conditions, device messages are typically sent in near real time or in up to five minutes while connected. This depends on device activity like applying the brakes, adjusting the rails, alarms, and how the third-party defines subscription times.

System requirements and recommendations for iBed Wireless (option) (Continued)

Customer network communication requirements for iBed Wireless option

LAN environment		
Client/server communication	IPv4 only	Not applicable
Client device IP allocation	Static	If Static - Unique IP address will be required for each client MAC address
	DHCP	If DHCP and not using a DNS name - Each client MAC address will need a reserved IP address If DHCP and using a DNS name - It is required to create a unique name for each client MAC address for client management Stryker recommends using the Stryker client host name when the Stryker device connects to the wireless network - Example: SYK-00197b12365 so it may look like http://SYK-00197b12365.
Server IP allocation	Static IP required	Not applicable
VLAN	New, existing	Install iBed Wireless on a separate VLAN

IP traffice environment		
Source	Protocol / Port number	Destination
iBed Server	TCP/21	Stryker iBed Wireless Client
iBed Server	TCP/80/443	Third party / Stryker back office
iBed Server	TCP/1639	Stryker iBed Wireless Client
Third party / Stryker iBed Wireless Client	TCP/80/443	iBed Server

Customer WLAN environment		
		Required
Supported wireless vendors	Cisco, Aruba	Yes
Access point (AP) types	Controller-based or autonomous	Yes
Channel width	2.4 GHz: 20 MHz 5 GHz: 20/40 MHz	Yes

System requirements and recommendations for *i*Bed Wireless (option) (Continued)

Customer WLAN environment		
		Required
Channel utilization	Consistently less than 30%	Recommended
Signal strength range (minimum)	2.4 GHz: -67dBm +0/-8dBm 5 GHz: -67dBm +0/-8dBm	Yes
Minimum SNR	Minimum 20dB	Yes
Priority queuing	Prioritized over best effort traffic	Recommended
Client exclusion	Disabled	Recommended
Client load balancing	Disabled	Recommended
Max number of SSIDs	5	Recommended
Authentication timeouts	Add session timeout of at least 24 hours	Recommended

Note: A transmit power asymmetry problem may arise at the edges of virtual cell coverage if an APs transmit power is higher than the Stryker Wireless Client device (~6 mW 2.4 GHz or 12 mW 5 GHz). The received signal strength indicator (RSSI) of the Stryker *iBed* Wireless Client on the AP must be verified. The device should never drop below an RSSI of -75 dBm on the AP.

Product illustration




Figure 1: Model 3005 S3 MedSurg bed

Α	110 VAC outlet (option)	ŀ
В	BackSmart® brake pedal	<u> </u>
С	Caster	
D	CPR release handle	
Е	Foley bag hooks	
F	Footboard	
G	BackSmart footboard control panel	
Н	Headboard	
I	iBed Awareness LED light bar	
J	Integrated pump rack	

K	Isolated Foley bag hooks
L	Mattress retainer
М	Motion interrupt pan
N	Night light
0	BackSmart operator control panel
	5
Р	Patient control panel
P Q	Siderail release handle
	'
Q	Siderail release handle

Contact information

Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.

Contact information (Continued)

Stryker Medical 3800 E. Centre Avenue Portage, MI 49002 USA

To view your operations or maintenance manual online, see https://techweb.stryker.com/.

Have the serial number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

Serial number location

You can find the serial number (A) above the power cord behind the headboard (Figure 2 on page 14).




Figure 2: Serial number location

Date of manufacture

The year of manufacture is the first four digits of the serial number.

Setup



MARNING WARNING

Always plug the product directly into a grounded, hospital grade wall outlet. You can only achieve grounding reliability when you use a hospital grade wall outlet. This product is equipped with a hospital-grade plug for protection against electric shock hazard.

Before you place the product into service, make sure these components are working properly:

- 1. Plug the product into a grounded, hospital grade wall outlet and make sure that the power button LED light at the foot end of the product comes on.
- 2. Make sure that the siderails raise, lower, lock in the up position, lock in the intermediate position when lowered and store smoothly.
- 3. Make sure that all four casters lock when the brake is applied.

Note: Make sure that the Brake LED located on the outside of the head end siderails (Operator control panel (outside siderail) on page 30) and on the footboard control panel (Figure 24 on page 39) blink when the brakes are released.

- 4. Raise the Fowler (head of bed) up to approximately 60°. Squeeze the CPR release handle and make sure that the back will drop with minimal effort.
- 5. Perform each function on the footboard control panel to make sure that each function is working.
- 6. Perform each function on both head end siderails to make sure that each function is working.
- 7. Activate the motion stop system to make sure that it is working. Press the **Bed height down** button to lower the litter. As the litter lowers, push up on the motion interrupt pan under the litter and make sure that the downward motion stops. Release the pan and allow the downward motion to continue.

Note: The product's upward motion or other functions are not disrupted by the motion stop system

Equipping the optional nurse call communication



WARNING

Always use a Stryker supplied interface cable. Use of any other cable may cause the bed to function improperly, which may result in patient or user injury.

If your product is equipped with the optional nurse call communication, follow the setup instructions below before putting the product into service.

1. Plug the interface cable into the 37-pin connector in the litter frame at the head end of the product (B) (Figure 3 on page 16) and into the applicable connection (patient station, head wall, or docker station).

Note: Only connect the 37-pin connector to the head wall output configuration (B) (Figure 3 on page 16) or product Communications Tester (sold separately).

- Test the interface cable to verify connectivity.
- 3. Plug the power cord into the wall outlet.
- 4. Push in the pendant port switch and turn the switch to the on position (90° clockwise) (A) (Figure 3 on page 16).
- 5. Press the Nurse call button (see Operator control panel (outside siderail) on page 30) to verify the connection between the product's nurse call signal and the hospital's nurse call system.

Note: A 9V battery, located in the head end of the product, powers the nurse call signal. If the product is unplugged from a wall outlet while the nurse call pendant port switch is set to the on position, the charge of the 9V battery will begin to drain.

Equipping the optional nurse call communication (Continued)

To activate the optional nurse call communication, see Activating nurse call communication (option) on page 26.




Figure 3: Nurse call pendant port and 37-pin connector

Setting up the optional iBed Wireless



WARNING

Always associate or map the iBed Locator to the room or location to provide accurate location information. Failure to properly associate or map the iBed Locator to the room or location will yield incorrect remote information. If you move an iBed Locator after it has been installed and mapped, you must remap to the new room or location. If the room or location information is changed after initial installation, you must remap the *i*Bed Locator.

To setup your product to receive a wireless connection, you must install the iBed Locator on the wall at the head end of the product. The iBed Locator communicates with the IR module that is installed in your product. For detailed instructions about how to mount the 5212 iBed Locator, see the instructions for use that was included with your optional 5212 iBed Locator installation kit.

Contact Stryker Technical Support at (800) 327-0770 with any installation questions.

Note:

You must load the wireless connection settings before the device will communicate with the iBed Server application. See the iBed Server Installation and Configuration manual.

Applying or releasing the brakes



WARNING

- Always apply the brakes when a patient is getting into or out of the product to avoid instability.
- Always apply the brakes when the patient is unattended.
- Do not apply the brakes to slow or stop the product while it is in motion.

You can find the brake pedals on both the left and right sides of the product.

To apply the brakes, fully depress the brake pedal down toward the floor.

To release the brakes, fully depress the brake pedal down toward the floor again.

Note: The Brake LED on the footboard control panel (Footboard LED indicators on page 38) illuminates when you apply the brakes. You can also set an audible alarm in the advanced menu (Setting the brake alarm on page 50).



Figure 4: Applying or releasing the brakes

Transporting the product with steer lock



WARNING

- Always lock the siderails in the full up position with the sleep surface horizontal in the lowest position when you transport the product with a patient.
- Always unplug the power cord before you transport the product.
- Always release the brakes before you transport the product. Do not transport the product with the brakes applied.
- Do not transport the product laterally after you apply the steer lock pedal. The product cannot swivel when you transport with steer lock.

Steer lock guides the product along a straight line during transport and pivots the product around corners. The steer lock pedal locks the right side caster on the foot end. You can find the steer lock pedal on the head end of the product.

To transport with steer lock:

- Align the wheels to face the direction of transport.
- 2. Push the steer lock pedal to on (Figure 5 on page 18).

To release steer lock, push the steer lock pedal to off (Figure 5 on page 18).

Note: To move the product in any direction, including laterally, release the steer lock pedal.

Transporting the product with steer lock (Continued)



Figure 5: Applying or releasing steer lock

Activating the CPR release



Always make sure that all persons and equipment are away from the area below and around the Fowler before you activate the CPR release. The CPR release is for emergency use only.

When you raise the Fowler and need quick access to the patient, you can position the product to 0° by activating the CPR release.

You can find the two CPR release levers at the head end section on both the left and right sides of the Fowler (A) (Figure 6 on page 19).

To activate the CPR release:

- 1. Grasp and squeeze the lever on either side of the Fowler.
- 2. Guide the Fowler to the flat position.

Activating the CPR release (Continued)




Figure 6: Activating the CPR release

Raising the lower leg section

Prop the footrest to raise the lower leg section manually.

To raise the lower leg section:

- 1. Grasp the lower leg section with both hands (A) (Figure 7 on page 20).
- 2. Swing the foot prop toward the foot end of the bed.
- 3. Release the foot prop when the prop rod is in the gatch bracket.

Raising the lower leg section (Continued)




Figure 7: Raising or lowering the lower leg section

Lowering the lower leg section

Prop the footrest to lower the lower leg section manually.

To lower the lower leg section:

- 1. Lift up the end of the lower leg section to release the prop rod from the gatch bracket.
- 2. Swing the foot prop toward the head end.
- Guide the foot prop down into the litter (B) (Figure 7 on page 20).

Attaching a fracture frame



WARNING

Use only retractable traction or fracture frames. Failure to use a retractable frame may result in injury to the patient and/or damage to the equipment.

You can attach a standard fracture frame to the product using the IV sockets that are located on all four corners of the bed. You can use IV poles with a fracture frame if you purchase the IV pole adaptor sockets (Accessories on page 59).

Securing a Foley bag to the Foley bag hooks

To secure a Foley bag to the Foley bag hook, place the hook of the Foley bag on the Foley bag hook. Make sure that you secure the Foley bag to the Foley bag hook.

You can find two isolated Foley bag hooks under the gatch section on both sides of the product (A) (Figure 8 on page 21). If you weigh the patient with the scale system, the isolated Foley bag weight is not included with the patient weight.

There are four Foley bag hooks under the seat section (B) and foot section (C) on both sides of the product (Figure 9 on page 21). If you weigh the patient with the scale system, the Foley bag weight is included with the patient weight.

Securing a Foley bag to the Foley bag hooks (Continued)





Figure 8: Isolated Foley bag hook

Figure 9: Foley bag hook

Installing the patient restraint strap tie-ins

CAUTION

 Always determine the proper use of the restraint straps and restraint strap locations. Improperly adjusted restraint straps can cause serious injury to a patient. Stryker is not responsible for the type or use of restraint straps on any Stryker products.

Installing the patient restraint strap tie-ins (Continued)

There are eight patient restraint strap tie-in locations on the litter assembly for installing patient restraint straps. Two are located on the Fowler section, two are located on the seat section, and four are located on the foot section (Figure 10 on page 22).




Figure 10: Restraint strap tie-in locations

Raising the siderails



WARNING

- Always set the siderail position to make sure that a patient is safely in the product.
- Always lock the operator control panel and patient control panel when the patient is unattended.
- Always keep the siderails outside of the oxygen tent.

Note: Do not use siderails as a patient restraint device.

When raising the siderails, listen for the click that indicates that the siderail has locked into position. Pull on the siderail to make sure that it is locked.

To raise the siderail to its highest position, grasp the release handle and rotate the siderail toward the head end of the product (Figure 11 on page 23).

Note: The siderail does not lock into the intermediate position when you raise the siderail.

Raising the siderails (Continued)



Figure 11: Siderail highest position

Lowering the siderails



WARNING

- Always set the siderail position to make sure that a patient is safely in the product.
- · Always lock the operator control panel and patient control panel when the patient is unattended.
- · Always keep the siderails outside of the oxygen tent.

Note: Do not use siderails as a patient restraint device.

To lower the siderail to the intermediate position, grasp the release handle and rotate the siderail forward until it locks into the intermediate position (Figure 12 on page 23).



Figure 12: Siderail intermediate position

To lower the siderail to its lowest position, grasp the release handle and rotate the siderail forward (Figure 13 on page 24).

Notes

- · You can stow the siderail under the litter when the siderail is at its lowest position.
- Make sure that the siderail is in the lowest position before you raise the siderail directly to the full up position. If you do not completely lower the siderail, the siderail will lock into the intermediate position when you raise the siderail.

Lowering the siderails (Continued)



Figure 13: Siderail lowest position

Positioning the optional two-stage permanently attached IV pole



Do not load the IV pole above the safe working load of 40 lb (18 kg).

Note: The two-stage permanently attached IV pole is an option that was installed at either the head, foot, or both ends of the bed. This choice was made when the product was purchased.

To position the IV pole:

- 1. Lift and pivot the IV pole from the storage position and push down until the IV pole locks into the receptacle.
- 2. To raise the height of the pole, pull up on the telescoping portion (A) of the pole until it locks into place at its fully raised position (Figure 14 on page 25).
- 3. Rotate the IV hangers (B) to the desired position and hang the IV bags (Figure 14 on page 25).
- 4. To lower the pole, turn the latch (C) clockwise until the telescoping portion (A) lowers into the bottom tube (Figure 14 on page 25).
- 5. Lift up and pivot the pole down into the storage position.

Positioning the optional two-stage permanently attached IV pole (Continued)




Figure 14: Two-stage IV pole

Positioning the optional removable IV pole

↑ CAUTION

- Do not load the IV pole above the safe working load of 40 lb (18 kg).
- 1. Install the pole at any of the four receptacles of the litter, located on all four corners of the frame.
- 2. To raise the height of the pole, turn the knob (A) counterclockwise and pull up on the telescoping portion (B) of the pole (Figure 15 on page 26).
- 3. Turn the knob (A) clockwise to tighten the telescoping portion (B) in place (Figure 15 on page 26).

Positioning the optional removable IV pole (Continued)

4. To lower the height of the pole, turn the knob (A) counterclockwise and lower the telescoping portion (B) of the pole into the bottom tube (Figure 15 on page 26).




Figure 15: Removable IV pole

Illuminating the room with the night light

There are two night lights that illuminate the floor area around the product. You can find the night light switch under the frame on the patient left side (N) (Product illustration on page 13).

To turn on the night light, turn the switch to on (Figure 16 on page 26).

To turn off the night light, turn the switch to off (Figure 16 on page 26).




Figure 16: Night light

Activating nurse call communication (option)

Nurse call allows the patient to send a signal to the nurse station.

To activate nurse call, press the **Nurse call** button (E) (see Operator control panel (outside siderail) on page 30). Communication between the patient and the nurse station is established at the moment when the nursing staff responds to the nurse call signal.

Activating nurse call communication (option) (Continued)

Note: A 9V battery, located in the head end of the product, powers the nurse call signal. If the product is unplugged from a wall outlet while the nurse call pendant port switch is set to the on position (Figure 17 on page 27), the charge of the 9V battery will begin to drain.

If the 9V nurse call battery needs to be replaced, a message will appear on the footboard display. See Replacing the nurse call backup battery (option) on page 27.




Figure 17: Nurse call switch

Replacing the nurse call backup battery (option)

To prevent a low battery condition when the bed is not plugged in, position the cord out switch at the head end of the bed to the off position. If the switch is not positioned as shown below (A) (Figure 18 on page 28) and the bed power cord and pendant cord are unplugged, the life of the backup battery will be significantly reduced.

Replacing the nurse call backup battery (option) (Continued)




Figure 18: Cord out switch position

If the 9V nurse call battery needs to be replaced, a message will appear on the footboard display. The battery is located on the patient's left side at the head end of the bed (B) (Figure 19 on page 28). No tools are required to replace the battery.




Figure 19: Battery housing

To replace the battery:

- 1. Unplug the power cord from the wall outlet.
- 2. Separate the **Velcro** that holds the battery housing (B) to the litter frame and lift the battery housing until it is free from the litter frame (Figure 19 on page 28).
- 3. Discard the battery.

Note: Follow local regulations for disposal or recycling of this part.

4. Reverse steps to reinstall.

Connecting peripheral equipment to the built-in 110 volt auxiliary power outlet (option)

\wedge

WARNING


- Only use hospital-grade electric equipment consuming 10A or less with the auxiliary power outlet (option). The use
 of standard electric equipment may bring the current leakage to a level unacceptable for hospital equipment.
- Do not use the 110V auxiliary power outlet (option) for life sustaining equipment.

The 110 volt auxiliary power outlet is a built-in outlet for peripheral equipment. You can find the outlet under the foot end on the patient right side of the product (Product illustration on page 13).



Figure 20: Optional auxiliary outlet

Operator control panel (outside siderail)



А	Brake LED	Flashes amber when you release the brakes. Brake LED turns off when you apply the brakes.
В	Cardiac chair position	Press and hold to place the product into the cardiac chair position
С	Bed height up	Raises the litter
D	Bed height down	Lowers the litter
Е	Nurse call (optional)	Activates nurse call
F	Gatch up	Raises the gatch
G	Gatch down	Lowers the gatch
Н	Fowler up	Raises the Fowler
I	Fowler down	Lowers the Fowler

Patient control panel (inside siderail)




Figure 21: Patient control panel

Patient control panel (inside siderail) (Continued)





Figure 22: Patient control panel with Smart TV (option)

Α	Room light (option)	Turns the room light on or off
В	Product overhead light (option)	Turns the product overhead light on or off
0	TV/Radio power (option) (Figure 21 on page 31)	Turns on the TV or the radio
С	Smart TV power (option) (Figure 22 on page 32)	Turns on Smart TV
Б	TV/Radio volume up (option) (Figure 21 on page 31)	Increases the volume
D	Smart TV volume up (option) (Figure 22 on page 32)	Increases Smart TV volume
E	TV/Radio volume down (option) (Figure 21 on page 31)	Decreases the volume
	Smart TV volume down (option) (Figure 22 on page 32)	Decreases Smart TV volume
F	Nurse call LED (option)	Illuminates amber when the patient presses the Nurse Call button
G	Nurse call (option)	Activates nurse call
Н	Nurse call answer LED (option)	Illuminates green when a nurse answers a call
I	Fowler up	Raises the Fowler
J	Fowler down	Lowers the Fowler

Patient control panel (inside siderail) (Continued)

K	Gatch up	Raises the gatch
L	Gatch down	Lowers the gatch


Optional Smart TV control panel (inside siderail)



Α	TV channel up	Changes the TV channel up
В	TV channel down	Changes the TV channel down
С	Mute	Turns the volume on and off.
D	Closed caption	Turns closed captions on and off.

www.stryker.com 3005-109-001 REV B 33

Footboard control panel - Bed controls



А	Trendelenburg	Places the product into the Trendelenburg position (head down with foot up)	
В	Reverse Trendelenburg	Places the product into the Reverse Trendelenburg position (head up with foot down)	
С	Bed height up	Raises the litter	
D	Bed height down	Lowers the litter	
Е	Fowler up	Raises the Fowler	
F	Fowler down	Lowers the Fowler	
G	Gatch up	Raises the gatch	
Н	Gatch down	Lowers the gatch	
I	Cardiac chair position	Press and hold to place the product into the cardiac chair position	
J	CPR	Press and hold to flatten the bed and lower it to low height	

Notes

- The CPR button overrides all lockouts.
- The Low Height footboard LED indicator illuminates when you lower the product to low height.


Footboard control panel - Lockouts



WARNING

- Always lock the control panel when you leave the patient unattended.
- · Always lock the patient control panel when the patient's condition requires extra safety measures.

Lockouts can lock out motion control input from the operator control panel, patient control panel, and motion pendants (optional). Bed Exit (option), scale (option), and nurse call (optional) features are still available.



А	Bed motion lock	Locks all motion controls from the operator control panel and patient control panel
В	Fowler 30°+	Press and hold to raise and lock the Fowler at 30°. Note: You can raise the Fowler between 30° and 60° after you lock the Fowler at 30°+.
С	Patient control Fowler lock	Locks or unlocks the Fowler
D	Patient control gatch lock	Locks or unlocks the gatch
Е	Bed height lock	Locks or unlocks the litter


Notes

- The CPR button overrides all lockouts.
- The corresponding lock LED illuminates when you lock a motion control.
- The Bed Motion Lock footboard LED indicator illuminates when you lock bed motion.
- The Fowler 30°+ footboard LED indicator illuminates when you lock the Fowler.
- If the product is held in a specific position when a lock is enabled, the product will be locked in that position.
- · Lock parameters are saved when the product is unplugged or during a power failure.
- Do not lock the control panel functions from the footboard if you must access the control panel functionality when you remove the footboard.

Footboard control panel - Chaperone bed exit (option)

WARNING

- Do not use bed exit to replace patient monitoring protocol.
- Bed exit is intended only to aid in the detection of a patient exiting the product.
- Do not use bed exit with patients who weigh less than 50 lb (23 kg).



А	Zone	The patient can move freely, but the alarm sounds when the patient attempts to leave the product
В	Arm/Disarm	Arms or disarms bed exit

Note: The Bed Exit footboard LED indicator illuminates when you arm Chaperone Bed Exit (Figure 24 on page 39).

Arming or disarming Chaperone bed exit (option)

When armed, Chaperone bed exit monitors the patient's position on the product.

Note: A notification appears if there is not enough weight on the product to arm bed exit.

To arm bed exit:

- 1. Set the scale to zero. See Setting the scale to zero on page 41.
 - Note: If you do not set the scale to zero before placing a patient on the product, bed exit may not operate properly.
- Position the patient on the product.
- Press and hold the Arm/Disarm button (A). See Footboard control panel Chaperone bed exit (option) on page 36.

After arming Chaperone Bed Exit, the LED light bars on the footboard illuminate green, the Bed Exit footboard LED indicator on the footboard LED indicator illuminates, and the selected zone on the footboard control panel illuminates.

If the parameter conditions selected for **Chaperone** bed exit are changed:

- LED light bars on the footboard flash amber
- bed exit indicator LED on the footboard LED indicator flashes
- sound alarm is triggered
- bed exit status alert is displayed on the display screen

Arming or disarming Chaperone bed exit (option) (Continued)


To disarm bed exit, press and hold the **Arm/Disarm** button again (A). See Footboard control panel - **Chaperone** bed exit (option) on page 36.

Footboard control panel - Chaperone bed exit with zone control (option)

Λ

WARNING

- · Do not use bed exit to replace patient monitoring protocol.
- · Bed exit is intended only to aid in the detection of a patient exiting the product.
- · Do not use bed exit with patients who weigh less than 50 lb (23 kg).



А	Zone 1	Allows the patient to move freely on the bed. Alarms when the patient moves 50 percent of their body weight out of the designated zone.	
В	Zone 2	Allows for limited movement. Alarms when the patient approaches th siderail or the foot end of the bed.	
С	Zone 3	Allows minimal movement. Alarms when the patient moves out of the tightly restricted zone.	
D	Arm/Disarm	Arms or disarms bed exit	
E	Zone	Changes the zone	

Note: The Bed Exit footboard LED indicator illuminates when you arm Chaperone Bed Exit (Figure 24 on page 39).

Arming or disarming Chaperone bed exit with zone control (option)

When armed, Chaperone bed exit monitors the patient's position on the product.

Note: A notification appears if there is not enough weight on the product to arm bed exit.

Arming or disarming Chaperone bed exit with zone control (option) (Continued)

To arm bed exit:

- 1. Set the scale to zero. See Setting the scale to zero on page 41.
 - Note: If you do not set the scale to zero before placing a patient on the product, bed exit may not operate properly.
- 2. Position the patient on the product.
- 3. Press and hold the **Arm/Disarm** button (D). See Footboard control panel **Chaperone** bed exit with zone control (option) on page 37.
 - Note: Zone 1 illuminates as the default zone when you arm Chaperone bed exit.
- 4. To change the zone, press the **Zone** button (E). See Footboard control panel **Chaperone** bed exit with zone control (option) on page 37.

After arming **Chaperone** bed exit, the LED light bars on the footboard illuminate green, the bed exit indicator LED on the footboard LED indicator illuminates, and the selected zone on the footboard control panel illuminates.

If the patient moves from the armed zone and changes the bed exit parameter:

- · LED light bars on the footboard flash amber
- bed exit indicator LED on the footboard LED indicator flashes
- sound alarm is triggered
- selected zone on the footboard control panel flashes
- · status alert is displayed on the display screen

To disarm bed exit, press and hold the **Arm/Disarm** button again (D). See Footboard control panel - **Chaperone** bed exit with zone control (option) on page 37.

Footboard LED indicators




Figure 23: Footboard LED indicators

	Indicator	LED illuminates amber
А	Bed motion lock	Bed motion lock is activated or when the patient control (Fowler, gatch, bed up/down) lock buttons are activated
В	Fowler 30°+	Fowler 30°+ is locked. The LED will blink amber if: • iBed Awareness system is on • Fowler 30°+ is being monitored and the Fowler goes below 30 degrees • Fowler 30°+ is turned off
С	Patient control Fowler lock	The patient control Fowler lock is on

Footboard LED indicators (Continued)

	Indicator	LED illuminates amber	
D	Patient control gatch lock	The patient control gatch lock is on	
Е	Patient control bed height lock	The patient control bed height lock is on	
F	Zone 1	Bed exit is on and zone 1 is active. The LED will blink amber if a bed exit event occurs.	
G	Zone 2	Bed exit is on and zone 2 is active. The LED will blink amber if a bed exit event occurs.	
Н	Zone 3	Bed exit is on and zone 3 is active. The LED will blink amber if a bed exit event occurs.	





Figure 24: Footboard LED indicators

	Indicator	LED illuminates amber
А	Bed motion lock	Bed motion lock is activated or when the patient control (Fowler, gatch, bed up/down) lock buttons are activated
В	Fowler 30°+	Fowler 30°+ is locked. The LED will blink amber if: • iBed Awareness system is on • Fowler 30°+ is being monitored and the Fowler goes below 30 degrees • Fowler 30°+ is turned off
С	Low height	Bed is in low height. The LED will blink amber if: • iBed Awareness system is on • Low height is being monitored • Bed is not in low height
D	Brake	Brake is set, and will blink amber if the brake is not set
E	Bed exit (optional)	Bed exit is armed. The LED will blink amber if: Bed exit is turned off while the <i>i</i> Bed Awareness system is turned on Bed exit alarms while monitored by the <i>i</i> Bed Awareness system
F	Bed zero (optional iBed Awareness)	Bed zero is successful
G	Siderail (optional iBed Awareness)	iBed Awareness system is on. The LED will blink amber when siderail state has changed.
Н	Power	Green when bed has power.

Footboard control panel - Scale



Α	Scale	Weighs the patient
В	Zero	Sets the scale to zero

Note: The Bed Zero LED indicator illuminates when you set the scale to zero (Figure 24 on page 39).

Weighing a patient



WARNING

- Do not use the scale system reading as a reference for medical treatment.
- The scale system assists only in the monitoring of the patient's weight variation.
- Do not set the scale to zero or weigh the patient when a support surface therapy is active. Motion from the support surface functions may adversely affect the scale system performance.



/ CAUTION

Always raise the siderails when the litter is in its full down position. This prevents the scale system from weighing a patient inaccurately.

To weigh a patient:

Press and hold the Scale button (A) (see Footboard control panel - Scale on page 40) until Release Button appears (Figure 25 on page 40).



Figure 25: Release button

Footboard control panel - Scale (Continued)

2. Release the Scale button.

Note: Do not touch the product when you weigh the patient (Figure 26 on page 41).



Figure 26: Do not touch bed

A confirmation notification indicates that weighing the patient was successful (Figure 27 on page 41).



Figure 27: Patient weight

Note: To clear the patient weight from the display, press the **Scale** button again. The patient weight is still recorded in the weight log.

Notes

- Always calibrate the bed after adding a support surface or mattress to the bed frame.
- Always set the scale to zero before putting a patient on the product.

Setting the scale to zero

The zero function clears all of the stored values from the weight log, change patient weight, and gain or loss.

Note: Always set the scale to zero before putting a patient on the product.

To set the scale to zero:

1. Press and hold the **Zero** button (B) (see Footboard control panel - Scale on page 40) (Figure 28 on page 41) until the **Release Button** appears (Figure 29 on page 41).



Figure 28: Hold to zero



Figure 29: Release Button

Footboard control panel - Scale (Continued)

2. Release the Zero button.

Note: Do not touch the product when you set the scale to zero (Figure 30 on page 42).



Figure 30: Do not touch bed

A confirmation notification indicates that setting the scale to zero was successful (Figure 31 on page 42).



Figure 31: Zeroing successful

Notes

- · The Bed Zero LED on the footboard LED indicator illuminates.
- · The scale icon and 0.0 appears on the footboard display.



Figure 32: Scale set to zero

Note: If you first receive an **Unable to Zero - Try Again** notification, the scale attempts to set the scale to zero again for 30 seconds. After three attempts, the scale system locks and an **Unable to Zero** notification appears.

Menu display

The S3 footboard control panel has a menu that displays the menus for S3 functions and features.

Menu display (Continued)




Figure 33: Menu

А	Menu	Acesses menu functions
В	Display	Displays menu functions
С	Up arrow	Scroll up through menu functions
D	Down arrow	Scroll down through menu functions
Е	Exit	Exits from menu functions or cancels operation
F	Enter	Selects menu function or saves operation

Menu functions		
Weight Log (Weight log is the default selection)	5. Scale Units (Change scale units)	
2. Gain/Loss	6. Backlight (Backlighting)	
3. Change Equip. (Change equipment)	7. Advanced Options	
4. Change Ptnt Wght (Change patient weight)	8. Exit Menu	

Accessing functions and features with the menu display

To access a menu option, press Menu (A) (see Menu display on page 42).

To scroll through menu options, press Arrow up (C) or Arrow down (D) (see Menu display on page 42).

To select a menu option, press Enter (F) (see Menu display on page 42).

To go back one menu or cancel a request, press Exit (E) (see Menu display on page 42).

Accessing functions and features with the menu display (Continued)

Notes

- · If no control panel or menu activity is detected within 60 seconds, the display backlight dims.
- The status screen shows the current Fowler angle and current Trendelenburg angle values by default.



Figure 34: Status screen

Menu display with iBed Wireless (option)

Information on the display includes the WiFi and *i*Bed Locator connection status, Fowler angle, and Trendelenberg angle values (Figure 35 on page 44).

Notes

- · If no control panel or menu activity is detected within 60 seconds, the display backlight dims.
- · The status screen shows the current Fowler angle and current Trendelenburg angle values by default.




Figure 35: iBed Wireless display screen

Icons	1 ×	T'	('T')	((₁))
Wireless connectivity status	Not connected or trying to connect	Connected		
Signal strength level	None	Low	Good	Excellent
Signal strength, X	X < -90 dB or X = 0 dB	-90 dB ≤ X < -71 dB	-71 dB ≤ X < -57 dB	X ≥ -57 dB

Viewing the weight log

Weight history displays measured and stored weight values. The system stores a maximum of 10 weight measurements. Any measurement taken after the tenth measurement deletes the oldest measurement taken.

Note: If the change in patient weight is less than .2 lb., the weight will not be stored.

To view weight history:

- 1. Press the Menu button (A). See Figure 33 on page 43.
- 2. Scroll to Weight log (Figure 36 on page 45).
- 3. Press the Enter button (F). See Figure 33 on page 43.



Figure 36: Weight log

Measuring weight gain or loss

The gain or loss function compares the initial patient weight to the current patient weight. It then displays the weight the patient has gained or lost since the initial patient weight reading.

To enable gain or loss:

- 1. Press the Menu button (A). See Figure 33 on page 43.
- 2. Scroll to Gain/Loss.
- 3. Press and hold the **Enter** button (F) (seeFigure 33 on page 43) (Figure 37 on page 45) until **Release Button** appears (Figure 38 on page 45).



Figure 37: Hold to enable



Figure 38: Release Button

Measuring weight gain or loss (Continued)

4. Release the Enter button.

Note: Do not touch the product when you measure weight gain or loss (Figure 39 on page 46).



Figure 39: Do not touch bed

A confirmation notification indicates that enabling Gain/Loss was successful.

The base measurement is the initial patient weight registered when you enable **Gain/Loss**. The weight difference between the base weight and the weight gained or lost is displayed on the bottom right hand corner (Figure 40 on page 46).



Figure 40: Base and gain/loss weight measurement

Note: If the weight gained or lost exceeds 99.9 lb., Error -.- appears on the display (Figure 41 on page 46).



Figure 41: Gain/loss error

To disable gain or loss:

- 1. Press the **Menu** button (A). See Figure 33 on page 43.
- 2. Scroll to Gain/Loss.
- 3. Press and hold the Enter button (F). See Figure 33 on page 43.

Changing equipment

Change equipment allows you to add or remove equipment or devices from the product without affecting the patient weight.

To change equipment:

1. Press the Menu button (A). See Figure 33 on page 43.

Changing equipment (Continued)

- 2. Scroll to Change Equip..
- 3. Press and hold the **Enter** button (F) (see Figure 33 on page 43) (Figure 42 on page 47) until **Add / Remove Equipment** appears (Figure 43 on page 47).



Figure 42: Hold to Change Equipment



Figure 43: Add / Remove Equipment

4. Release Enter.

Note: Do not touch the product when you change equipment (Figure 44 on page 47).



Figure 44: Do not touch bed

A confirmation notification indicates when you are able to add/remove equipment (Figure 43 on page 47).

- 5. Press the Enter button (F). See Figure 33 on page 43.
 - Note: Do not touch the product when configuring to change equipment.
- 6. Add or remove equipment or devices from the product.
- 7. After adding or removing equipment or devices from the product, press the **Enter** button (F). See Figure 33 on page 43.

To cancel the request, press the Exit button (E). See Figure 33 on page 43.

Changing the patient weight

To change the patient weight:

- 1. Press the Menu button (A). See Menu display on page 42.
- 2. Scroll to Change Ptnt Wght.

Changing the patient weight (Continued)

3. Press and hold the **Enter** button (F) (see Menu display on page 42) (Figure 45 on page 48) until **Release Button** appears (Figure 46 on page 48).



Figure 45: Hold to change patient weight



Figure 46: Release Button

4. Release the Enter button.

Note: Do not touch the product when configuring patient weight (Figure 47 on page 48).



Figure 47: Do not touch bed

- 5. When the system is ready to change patient weight, press the **Up arrow** button (C) or **Down arrow** button (D) to change the displayed weight (see Menu display on page 42).
- 6. After changing the weight, press the **Enter** button (F) (see Menu display on page 42) and the message **Patient** Weight Changed will display.

To cancel the request, press the **Exit** button (E) (see Menu display on page 42) and the message **Operation Canceled** will display.

Changing the scale units

You can change the measuring unit to pounds (lb) or kilograms (kg) on your display.

Note: The default scale unit is pounds (lb).

To change the displayed scaled units:

- 1. Press the Menu button (A). See Menu display on page 42.
- 2. Scroll to Scale Units.

Changing the scale units (Continued)

- 3. Press the Enter button (F). See Menu display on page 42.
- 4. Select a scale unit (Figure 48 on page 49).




Figure 48: Scale Units

5. Press the Enter button (F). See Menu display on page 42.

To cancel the request, press the Exit button (E). See Menu display on page 42.

Changing the backlight intensity

The backlight changes the LED backlight intensity for all control panels (operator control panel, patient control panel, and footboard control panel).

Note: The default backlight intensity is low.

Five settings are available for the control panel LED intensity.

Setting	LED intensity
1	Off
2	Low
3	Medium
4	High
5	Nurse call only

Note: The nurse call LED backlight on the patient control panel shows the patient which button to press to contact the nurse's station. Turning the nurse call LED backlight light off may compromise this ability in a darkened room.

To change the backlight LED intensity:

- 1. Press the **Menu** button (A). See Menu display on page 42.
- 2. Scroll to Backlight.
- 3. Press the Enter button (F). See Menu display on page 42.
- 4. Select a backlight intensity (Figure 49 on page 49).




Figure 49: Backlight

5. Press the **Enter** button (F). See Menu display on page 42.

Setting the alarm tones

S3 has 10 alarm tone settings.

Note: The alarm tone you choose is the same tone used for all activated alarm options.

To set an alarm tone:

- 1. Press the Menu button (A). See Menu display on page 42.
- 2. Scroll to Advanced options.
- 3. Press the Enter button (F). See Menu display on page 42.
- 4. Scroll to Choose Exit Alarm.
- 5. Press the Enter button (F). See Menu display on page 42.
- 6. Scroll through the alarm tones.

Note: A brief sample of the tone plays as you scroll through the tone options.

7. Press the Enter button (F). See Menu display on page 42.

Note

A confirmation notification indicates that you set the alarm tone.

Setting the brake alarm

You can set an alarm to alert you when the brake is not set and the product is plugged in.

To set the brake alarm:

- 1. Press the Menu button (A). See Menu display on page 42.
- 2. Scroll to Advanced options.
- 3. Press the Enter button (F). See Menu display on page 42.
- 4. Scroll to Brake Alarm.
- 5. Press the Enter button (F). See Menu display on page 42.
- 6. Select On.
- 7. Press the Enter button (F). See Menu display on page 42.

Note:

A confirmation notification indicates that you set the alarm tone.

Setting an audible iBed Awareness alarm

You can set an audible alarm to alert you when an *iBed* Awareness parameter condition is compromised and the product is plugged in.

To set an audible iBed Awareness alarm:

- 1. Press the Menu button (A). See Menu display on page 42.
- 2. Scroll to Advanced options.
- 3. Press the Enter button (F). See Menu display on page 42.
- 4. Scroll to Awareness Alarm.
- 5. Press the **Enter** button (F). See Menu display on page 42.
- 6 Select On
- 7. Press the Enter button (F). See Menu display on page 42.

Note:

A confirmation notification indicates that you set the alarm tone.

Setting the iBed Awareness nurse call alarm


S3 sends a signal through the nurse call system when a parameter condition is compromised. To set up an *iBed* Awareness alarm through the nurse call system:

- 1. Press the Menu button (A). See Menu display on page 42.
- 2. Scroll to Advanced options.
- 3. Press the Enter button (F). See Menu display on page 42.
- 4. Scroll to Status to nurse call.
- 5. Press the Enter button (F). See Menu display on page 42.
- 6. Select On.
- 7. Press the Enter button (F). See Menu display on page 42.

Note

A confirmation notification indicates that you set the alarm tone.

iBed Awareness (option)



Note: For a description of the menu functions, see Menu display on page 42.

А	iBed Arms/disarms iBed Awareness		
В	Menu Accesses menu functions		
С	Display	Displays menu functions	
D	Up arrow	Scroll up through menu functions	
E	Down arrow	Scroll down through menu functions	
F	Exit	Exits from menu functions or cancels operation	
G	Enter	Selects menu functions or saves operation	

iBed Awareness (option) (Continued)

Configuring iBed Awareness



WARNING

- Do not use *i*Bed Awareness to replace your patient monitoring protocol.
- Do not use iBed Awareness as a lock indicator for siderails. iBed Awareness only detects the position of the siderails.

CAUTION

- Always make sure that the siderails are locked before you arm iBed Awareness.
- Make sure that you set the desired product parameters before arming iBed Awareness.
- Do not use accessories that cover the footboard LED light bars.

When enabled, iBed Awareness helps to monitor S3 status and parameter conditions.

To monitor a parameter, place the product to the desired position. You can monitor the low height position, Chaperone bed exit, and Fowler 30°+.

When armed, iBed Awareness automatically monitors all current siderail positions and the brake.

To arm iBed Awareness, press iBed (A) (see iBed Awareness (option) on page 51).

Note: If there is an error in one of the product functions, an error code will appear. iBed Awareness will not arm. For more information about error codes, see the maintenance manual.

After you armibed Awareness, the LED light bars on the footboard illuminate green and the monitored footboard LED indicators on the footboard illuminate.

To set an alarm tone for iBed Awareness, see Setting an audible iBed Awareness alarm on page 50.

To disarm iBed Awareness, press iBed again (A) (see iBed Awareness (option) on page 51).

Note: The settings for lockout controls, scale calibration data, bed exit, and iBed Awareness are preserved when the product is unplugged, or during a power failure.

Acknowledging iBed Awareness status alerts

If the parameter conditions selected for *iBed* Awareness are changed:

- LED light bars on the footboard flash amber
- changed indicator LED on the footboard LED indicator flashes
- sound alarm is triggered
- changed parameter status alert is displayed on the display screen

iBed Awareness (option) (Continued)

If the Low Height position changes (Figure 50 on page 53), return the product to low height (see operator control panel or the footboard control panel).




Figure 50: Low Height status alert

If the brakes are released (Figure 51 on page 53), apply the brake (see Applying or releasing the brakes on page 17).




Figure 51: Brake status alert

If a siderail position changes (Figure 52 on page 53), return the affected siderail to its original position. See Raising the siderails on page 22 and Lowering the siderails on page 23.

Note: The arrow in the status alert points to the affected siderail.




Figure 52: Siderail status alert

If a patient moves from a Chaperone bed exit zone (Figure 53 on page 53), you must:

- disarm bed exit, see Arming or disarming Chaperone bed exit (option) on page 36 or Arming or disarming
 Chaperone bed exit with zone control (option) on page 37
- disarm iBed Awareness, see Configuring iBed Awareness on page 52

Then, return the patient to the product and position them to the monitored zone, arm bed exit, and arm iBed Awareness.




Figure 53: Chaperone Bed Exit status alert

If Fowler 30°+ becomes unlocked (Figure 54 on page 54), lock Fowler 30°+.

iBed Awareness (option) (Continued)




Figure 54: Fowler 30°+ lock status alert

If the Fowler 30°+ position changes (Figure 55 on page 54), return Fowler 30°+ to its original position.






Figure 55: Fowler 30°+ position status alert

Motion pendant with nurse call (option)




Α	Nurse call	Activates nurse call
В	Fowler up	Raises the Fowler
С	Fowler down	Lowers the Fowler
D	Gatch up	Raises the gatch
E	Gatch down	Lowers the gatch

Motion and communication pendant (option)



А	TV	Turns the TV on or off
В	Nurse call	Activates nurse call
С	Volume up	Increases the volume
D	Volume down	Decreases the volume
Е	TV Channel up	Changes the TV channel up
F	TV Channel down	Changes the TV channel down
G	Fowler up	Raises the Fowler
Н	Fowler down	Lowers the Fowler
1	Gatch up	Raises the gatch
J	Gatch down	Lowers the gatch
К	Room light	Turns the room light on or off
L	Product overhead light	Turns the product overhead light on or off

Infrared (IR) module (option)



А	IR (infrared) module cable	Connects to the bed and provides power and signal communications
В	IR (infrared) lens	Provides infrared communications with the <i>i</i> Bed Locator
С	iBed Locator connection LED	Provides connection status for the IR (infrared) communications with <i>i</i> Bed Locator Slow flash - Attempting to connect to <i>i</i> Bed locator Solid LED - <i>i</i> Bed locator connected Rapid flash - Error condition detected Off - <i>i</i> Bed locator is not trying to connect
D	Wireless (WiFi) connection LED	Provides connection status for wireless (WiFi) communications with wireless access point Slow flash - WiFi attempting to connect Solid LED - WiFi connected Rapid flash - WiFi was not connected after six minutes and timed out

iBed Locator (option)

The *i*Bed Locator provides the *i*Bed Locator ID and battery status information to the IR module. See the *i*Bed Locator Instructions For Use manual for installation and operational procedures for the *i*Bed Locator.

The IR (infrared) lens (A) provides infrared communications with the *i*Bed IR module (Figure 56 on page 58).

iBed Locator (option) (Continued)




Figure 56: iBed Locator

Accessories

These accessories may be available for use with your product. Confirm availability for your configuration or region. Call Stryker Customer Service: 1-800-327-0770.

Name	Number
Adapter frame assembly	3006-333-000
Bed extender assembly	3006-346-020S
Bed extender assembly with mattress	3006-700-047
Bed extender mattress	3000-318-050
Defibrillator tray assembly	3006-120-004
IV pole latch assembly	0785-035-103
IV pole, removable	3000-300-080
IV pole, removable assembly	3001-338-010
IV pole, two-stage, dual assembly	2035-113-000
IV pole, two-stage, head left	2035-112-010
IV pole, two-stage, head left assembly	2035-112-000
IV pole, two-stage, head right	2035-113-011
Roller assembly	3006-345-005S
Roller assembly, litter	3006-335-000
Roller assembly, patient helper	3006-345-000
Traction socket adapter, 1/2"	3000-337-050
Traction socket extension, 4" x 1/2"	3000-337-450
Traction socket extension, 4" x 3/4"	3000-337-475
Traction socket extension, 8" x 1/2"	3000-337-850
Traction socket extension, 8" x 3/4"	3000-337-875
Upright oxygen bottle holder assembly	3006-150-000
Wall saver, single	3001-344-835
Wall saver, 10 pack	3001-344-840

Cleaning

! CAUTION

- · Do not clean, disinfect, service, or perform maintenance while the product is in use.
- Always unplug the power cord from the wall outlet when large spills occur near the circuit boards, cables, and
 motors. Remove the occupant from the product, clean up the fluid, and inspect the product. Fluids can cause
 unpredictable operation and decreased functionality of any electrical product. Do not return the product to service
 until it is completely dry and you have tested for safe operation.
- Always clean Velcro® after each use. Saturate Velcro with disinfectant and allow disinfectant to evaporate.
 Appropriate disinfectant for nylon Velcro should be determined by the hospital.

To remove undesirable build-up before disinfecting between uses:

- 1. Hand wash all surfaces on the product with a mild detergent using spray or pre-soaked wipes.
- 2. Clean all exposed surfaces.
- 3. Follow the cleaning solution manufacturer's instructions for appropriate contact time and rinsing requirements.
- 4. Dry the product thoroughly before returning the product to service.

Avoid over saturation. Do not allow the product to remain wet.

Disinfecting

⚠ CAUTION

- · Do not clean, disinfect, service, or perform maintenance while the product is in use.
- Always unplug the power cord from the wall outlet when large spills occur near the circuit boards, cables, and
 motors. Remove the occupant from the product, clean up the fluid, and inspect the product. Fluids can cause
 unpredictable operation and decreased functionality of any electrical product. Do not return the product to service
 until it is completely dry and you have tested for safe operation.
- Always wipe down with clean water (or 70% isopropyl alcohol, if using Virex® TB) and dry each product after
 disinfecting. Some disinfectants are corrosive in nature and may cause damage to the product. If you do not rinse
 and dry the product, you may leave a corrosive residue on the surface of the product. This corrosive residue could
 cause premature degradation of critical components. Failure to follow these disinfecting instructions may void your
 warranty.
- Always clean Velcro® after each use. Saturate Velcro with disinfectant and allow disinfectant to evaporate.
 Appropriate disinfectant for nylon Velcro should be determined by the hospital.
- Do not use quaternary disinfectants formulated with glycol ethers.

The recommended disinfectants for this product's surfaces include the following:

- · Quaternary (active ingredient ammonium chloride)
- Phenolic (active ingredient o-phenylphenol)
- Chlorinated bleach solution (5.25% less than 1 part bleach to 100 parts water)

To wipe down the product with disinfectant between uses:

- 1. Follow the manufacturer's dilution recommendations exactly.
- 2. Apply the recommended disinfectant solution by spray or pre-soaked wipes.
- 3. Hand wash all surfaces of the product with the recommended disinfectant.
- 4. Disinfect all exposed surfaces.
- 5. Follow the disinfecting solution manufacturer's instructions for appropriate contact time and rinsing requirements.
- 6. Dry the product thoroughly before returning the product to service.

Avoid over saturation. Do not allow the product to remain wet.

Follow the manufacturer's dilution recommendations for appropriate contact time and rinsing requirements. Follow the chemical manufacturer's guidelines for proper disinfecting.

Preventive maintenance

Remove the product from service before you perform the preventive maintenance inspection. Make sure that all items listed during annual preventive maintenance for all Stryker Medical products. You may need to perform preventive maintenance checks more often based on your level of product usage. Service only by qualified personnel.

Note: Clean and disinfect the exterior of the support surface before inspection, if applicable.

Inspect	the following items:
	All fasteners are secure
	Engage brake pedal and push on bed to make sure that all casters lock securely
	LED on the footboard and head end siderails blink when brakes are not engaged
	Locking steer caster locks and unlocks
	Siderails move, latch, and stow
	CPR release operable
	Foot prop intact and operable
	IV pole is intact and operable (optional equipment)
	Foley bag hooks intact
	No cracks or splits in headboard, footboard or siderail panels
	No rips or cracks in mattress cover
	All functions on head end siderails operable (includes LEDs)
	All functions on footboard operable (includes LEDs)
	Scale and bed exit system calibrated
	Motion interrupt switches operable
	Night light operable
	Power cord and plug not frayed or damaged
	No cables worn or pinched
	All electrical connections tight
	All grounds secure to the frame
	Ground impedance not more than 200 m Ω (milliohms)
	Current leakage not more than 300 μA (microamps)
	Apply grease to the litter grease points (see maintenance manual for locations)
	Make sure that ground chains are clean, intact, and have at least two links touching the floor
	Make sure that Fowler angle accuracy is 0° - 60°
	Make sure that the Fowler holds position at 30° with patient weight
	Siderail switches operable (iBed Awareness option)
	Center light bar LED and side light LED operable (iBed Awareness option)
	Inspect footboard control labels for signs of degradation
	Inspect siderail gas spring for oil leaks
	Inspect Fowler damper for oil leaks
	Make sure that all motions function
	Make sure that the nurse call functions
	Make sure that the nurse call battery functions (optional equipment)
	Make sure that the iBed Wireless Module and IR Module are intact and footboard icons display (iBed
	Awareness option)
Produ	ct serial number:
Comp	leted by:
Date:	

FCC notification

Notifications

FCC ID: Z7A-SDMANIC NO.: 4919E-SDMAN

Notice

Federal Communication Interference Statement (United States Only)

This equipment has been tested and found to comply with the limits for a class B digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- · Re-orient or relocate the receiving antenna
- · Increase the separation between the equipment and receiver
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected
- Consult the dealer or an experienced radio/TV technician for help

This device and its antenna(s) must not be co-located or operation in conjunction with any other antenna or transmitter.

For product available in the USA/Canada market, only channels 1-11 can be operated. Selection of other channels is not possible.

If this device is to be operated in the 5.15°5.25GHz frequency range, it is restricted to indoor environments only.

Antenna: Proprietary

Antenna gain information: Embedded Antenna: 2.5dBi (2.4 GHz), 3.5dBi (5 GHz)

Frequency Tolerance +/-20ppm

Guidance and manufacturer's declaration - electromagnetic emissions

The 3005 **S3** MedSurg bed is intended for use in the electromagnetic environment specified below. The customer or the user of the 3005 **S3** MedSurg bed should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment
RF Emissions CISPR 11	Group 1	The 3005 S3 MedSurg bed uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The 3005 S3 MedSurg bed is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations Flicker Emissions IEC 61000-3-3	Complies	

Recommended separations distances between portable and mobile RF communication equipment and the 3005 S3 MedSurg bed

The 3005 **S3** MedSurg bed is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the 3005 **S3** MedSurg bed can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the 3005 **S3** MedSurg bed as recommended below, according to the maximum output power of the communications equipment.

	Separation distance according to frequency of transmitter			
Rated maximum output power of transmitter	m			
w	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz	
	D=(1.2) (√ <i>P</i>)	D=(1.2) (√ <i>P</i>)	D=(2.3) (√ <i>P</i>)	
0.01	1.12	0.12	0.23	
0.1	0.38	0.38	0.73	
1	1.2	1.2	2.3	
10	3.8	3.8	7.3	
100	12	12	23	

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration - electromagnetic immunity

The 3005 S3 MedSurg bed is suitable for use in the electromagnetic environment specified below. The customer or the user of the 3005 S3 MedSurg bed should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	<u>+</u> 6 kV contact <u>+</u> 8 kV air	<u>+</u> 6 kV contact <u>+</u> 8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic fast transient/ burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/ output lines	±2 kV for power supply lines ±1 kV for input/ output lines	Main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±8 kV lines to lines ±2 kV lines to earth	±8 kV lines to lines ±2 kV lines to earth	Main power quality should be that of a typical commercial or hospital environment.
Voltage dips, voltage variations and short interruptions on power supply input lines IEC 61000-4-11	$<5\%U_T \ (>95\% \ dip \ in \ U_T)$ for 0,5 cycle $40\%U_T \ (60\% \ dip \ in \ U_T)$ for 5 cycles $70\%U_T \ (30\% \ dip \ in \ U_T)$ for 25 cycles $<5\% \ U_T \ (>95\% \ dip \ in \ U_T)$ for 5 sec.	$ <5\%U_T \ (>95\% \ dip \ in \ U_T) $ for 0,5 cycle $ 40\%U_T \ (60\% \ dip \ in \ U_T) $ for 5 cycles $ 70\%U_T \ (30\% \ dip \ in \ U_T) $ for 25 cycles $ <5\% \ U_T \ (>95\% \ dip \ in \ U_T) $ for 5 sec.	Main power quality should be that of a typical commercial or hospital environment. If the user of the 3005 S3 MedSurg bed requires continued operation during power main interruptions, it is recommended that the device be powered from an uninterrupted power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: U_T is the a.c. mains voltage before applications of the test level.

(Continued)

Conducted RF IEC 61000- 4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 Vrms 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the 3005 S3 MedSurg bed, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended separation distance D=(1.2) (√P) 80 MHz to 800 MHz D=(2.3) (√P) 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site surveya, should be less than the compliance level in each frequency rangeb. Interference may occur in the vicinity of equipment marked with the following symbol:
--	---	-----------------	---

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

(Continued)

^aField strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the 3005 **S3** MedSurg bed is used exceeds the applicable RF compliance level above, the 3005 **S3** MedSurg bed should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the 3005 **S3** MedSurg bed.

bOver the frequency range 150 kHz to 80 MHz, field strengths are less than 3 V/m.

Warranty

Stryker Medical, a division of Stryker Corporation, warrants to the original purchaser the Stryker Model 3005 **S3** MedSurg bed, to be free from defects in material and workmanship for a period of 1 year after date of delivery. Stryker's obligation under this warranty is expressly limited to supplying replacement parts and labor for, or replacing at its option, any product which is, in the sole discretion of Stryker, found to be defective. If requested by Stryker, product or parts for which a warranty claim is made shall be returned prepaid to the factory. Any improper use or any alteration or repair by others in such a manner as in Stryker's judgment affects the product materially and adversely, shall void this warranty. Any repair of Stryker products using parts not provided or authorized by Stryker shall void this warranty. No employee or representative of Stryker is authorized to change this warranty in any way.

Stryker Medical bed products are designed for a 10 year expected service life under normal use, conditions, and with appropriate periodic maintenance as described in the maintenance manual for each device. Stryker warrants to the original purchaser that the welds on its bed products will be free from structural defects for the expected 10 year life of the bed product as long as the original purchaser owns the product.

The above noted warranty periods apply only to the original purchaser of the 3005 **S3** MedSurg bed and begin on the date of delivery to such original purchaser.

Warranty exclusion and damage limitations

The express warranty set forth herein is the only warranty applicable to the product. Any and all other warranties, whether express or implied, including any implied warranty of merchantability or fitness for a particular purpose are expressly excluded by Stryker. In no event shall Stryker be liable for incidental or consequential damages.

To obtain parts and service

Stryker products are supported by a nationwide network of dedicated Stryker Field Service Representatives. These representatives are factory trained, available locally, and carry a substantial spare parts inventory to minimize repair time. Simply call your local representative or call Stryker Customer Service at 1-800-327-0770.

Return authorization

Product cannot be returned without prior approval from the Stryker Customer Service Department. An authorization number will be provided which must be printed on the returned product. Stryker reserves the right to charge shipping and restocking fees on returned product. Special, modified, or discontinued products are not subject to return.

Damaged product

ICC Regulations require that claims for damaged product must be made within fifteen (15) days of receipt of the product. Do not accept damaged shipments unless such damage is noted on the delivery receipt at the time of receipt. Upon prompt notification, Stryker will file a freight claim with the appropriate carrier for damages incurred. Claims will be limited in amount to the actual replacement cost. In the event that this information is not received by Stryker within the fifteen (15) day period following the delivery of the product, or the damage was not noted on the delivery receipt at the time of receipt, the customer will be responsible for payment of the original invoice in full within thirty (30) days of receipt. Claims for any incomplete shipments must be made within thirty (30) days of invoice.

International warranty clause

This warranty reflects U.S. domestic policy. Warranty outside the U.S. may vary by country. Contact your local Stryker Medical representative for additional information.



