## HandStand™ Traction Tower

ENGLISH	1
FRANÇAIS	7
ITALIANO	14
DEUTSCH	21
ESPAÑOL	28
ΕΛΛΗΝΙΚΑ	35
PORTUGUÊS	42
PORTUGUÊS DO BRASIL	49
NEDERLANDS	56
DANSK	63
SVENSKA	70
SUOMI	76
NORSK	83
简体中文	90
日本語	95
POLSKI	101
TÜRKÇE	108
ROMÂNĂ	114
العربية	121
РУССКИЙ	126
HRVATSKI	133
ČEŠTINA	140
MAGYAR	147
LIETUVIŲ K.	154
SLOVENČINA	160
SLOVENŠČINA	167
한국어	173

# CONMED HANDSTAND™ TRACTION TOWER

Catalog Number	Description	Sterilization Method
HS-TT	HandStand™ Traction Tower	Steam

## DESCRIPTION

The **HandStand™ Traction Tower** provides the surgeon with a sterile method of positioning and distracting the patient's arm (elbow to fingertips) for orthopedic hand and wrist surgical procedures. The traction tower features a tension scale which allows for traction maintenance and monitoring via the reference marks and numbers. Another feature of the tower is the adjustable pivot joints which allow for wrist flexion and extension, as well as radial or ulnar deviations, while maintaining traction. Set-up configurations of the traction tower allow for dorsal, as well as volar, approaches. Single-use, sterile Digit Traps will allow for an atraumatic method of securing the patient's fingers. A bead chain at the distal end of each trap allows for quick patient positioning. Single-use, sterile straps also provide an atraumatic method of securing the patient's forearm and biceps to the traction tower. The stand-alone design of the traction tower device provides stability without requiring fixation to an operating table. Between cases, the system is easily disassembled for cleaning and sterilization. An insert tray provides a place for the reusable traction tower to fit into, which then goes into a sterilization container.

## COMPATIBILITY

The HandStand Traction Tower is designed to work with HandStand Digit Traps (HS-DT2) and the HandStand Forearm and Biceps Straps (HS-FS and HS-BS). Additionally, the HandStand Traction Tower is compatible with the ConMed HandStand Sterilization Tray (HS-ST) and ConMed Sterilization Container (INST-C352).

The HandStand Traction Tower is also compatible with our legacy Digit Traps, 9906 and 9903.

Use of the HandStand Traction Tower with other devices not listed in compatibility section is not recommended.

## MATERIAL SPECIFICATIONS

Traction Tower: Stainless Steel, Aluminum, PEEK and Brass

## INTENDED USE

The HandStand Traction Tower System allows a surgeon to position the patient's arm and apply traction to the patient's arm (elbow to fingertips) during orthopedic surgical procedures.

The HandStand Traction Tower System can be used in open, minimally invasive and arthroscopic extremity procedures.

## WARNINGS

- Additional warnings include those applicable to any surgical procedure. In general, careful attention must be paid to asepsis and avoidance of anatomical hazards.
- The device is supplied non-sterile and must be cleaned, disinfected, and sterilized before the first use and after every use.
- When steam autoclave sterilization methods are used before a procedure, adequate cooling time must be allowed before use in surgery in order to prevent injury to patient or user.
- The Scale is not a measurement device, it is to be used as a reference when increasing traction or decreasing traction on patient anatomy, as well an indication that traction is maintained. It is at the discretion of the surgeon to increase or decrease traction on the patient. Careful attention to the Scale should be used when adjusting traction, as well as careful consideration for patient anatomy under traction. Application of excessive traction on patient anatomy may result in patient injury.
- During assembly of the tower, ensure that the sliding pin on the Lower Tower drops and is locked into Base Plate.
- Use caution when adjusting the HandStand Traction Tower while sharp objects are in the patient.
- Always ensure that the Lower Arm is firmly grasped before actuating any levers or buttons, or injury to patient may occur.
- Placing the device on an unstable or uneven surface may lead to the tower moving during use, and possible patient injury. Always ensure the device is placed on a stable surface before use.
- If the Traction Tower is improperly setup or inadvertent adjustment actuation occurs, the Lower Arm and/or the Upper Arm may fall, possibly resulting in patient injury.

## ENGLISH/INSTRUCTIONS FOR USE

10. Any use for specific procedures that require the use of only one Digit Trap should take into consideration the risk to patient anatomy, and never exceed 15 on the Scale.
11. Failure to ensure user and patient are free of pinch points during setup may lead to injury.
12. Overtightening of straps or Digit Traps may lead to patient injury.
13. If patient's fingers do not fit in Digit Trap or remain secure in Digit Trap while under traction, alternative methods for surgery should be pursued, or injury to patient may occur.
14. When maneuvering fully assembled Traction Tower, avoid lifting by the Upper Arm, to prevent a drop hazard.
15. Proper setup of devices should be followed to prevent damage to device and possible injury to patient or user.
16. Traction applied above the 20 lb reference will not be registered on the scale. It is at the discretion of the surgeon whether to induce a greater tension on the patient's limb.

## PRECAUTIONS


1. Federal Law (USA) restricts this device to sale by or on the order of a physician.
2. Prior to use, cleaning or sterilization, remove all protective packaging and tip protector, if applicable.
3. Inspect device prior to use to ensure it is in good physical condition and functions properly. There should be no loose, broken or misaligned parts.
4. Exercise care in the use of the device to minimize side or bending loads.
5. Do not use excessive force on device to avoid damage or breakage during use.
6. Proper setup of device should be followed to prevent damage to device and possible injury to patient or user.
7. Device is designed to work with a large range of patient arm sizes. If patient arm size is too small or large, alternate surgical method must be used.
8. Inspect device after use to ensure it has not been damaged.
9. Device is designed for use by surgeons experienced in the appropriate specialized procedures. It is the responsibility of the surgeon to become familiar with the proper techniques for use.
10. If device needs service, please contact customer service.

## ADVERSE EVENTS

1. Infections.
2. Allergies, tissue irritation/inflammation and other reactions to device materials.

**Diagram for HS-TT**

## Traction Tower Assembly

Components shown in diagram above

1. Base Plate
2. Lower Tower
3. Tower Locking Knob
4. Upper Tower and Rotation Joint
5. Lower Arm
6. Scale Rack
7. Upper Arm and Scale with Traction Knob

## INSTRUCTIONS FOR USE

1. The HandStand Traction Tower is designed to be used in conjunction with either an arm board or handtable. Prior to positioning the patient, all pads must be removed from the arm board or handtable.

**WARNING: Ensure device is placed on stable surface before use.**

2. Assemble the tower by mating components in sequential order of etching.

3. First lock the Lower Tower into the Base Plate by inserting key on bottom of Lower Tower into key slot on Base Plate and rotating until Lower Tower sliding pin locks in place.

**WARNING: Ensure that the sliding pin on the Lower Tower drops and is locked into Base Plate.**

4. Connect the Upper Tower and Rotation Joint into the Lower Tower with the Tower Locking Knob.

The Upper Tower can be adjusted depending on patient arm length.

**WARNING: Ensure the Tower Locking Knob is adequately tightened to hold the Upper Tower to the Lower Tower. Failure to do so could result in the Upper Tower dropping to the lowest adjustment position.**

5. Depress button on the top of the Upper Tower and Rotation Joint and insert Lower Arm into Rotation Joint opening. Release the button and ensure the Lower Arm is locked in the vertical position.

## ENGLISH/INSTRUCTIONS FOR USE

6. Depress lever in Lower Arm and insert Upper Arm. The height of the Upper Arm can be adjusted depending on patient arm length.
7. Hook Scale Rack on threaded rod in Upper Arm. The Traction Knob at the top of the Upper Arm allows for adjustment of traction after the arm is setup in the traction tower.
8. To position the patient's arm for surgery, flex the patient's arm 90° and rest the elbow on Base Plate against the Upper Tower. If a tourniquet is to be used, it should be in place at this point.  
**WARNING: If the Base Plate, Upper Tower, or any other point of contact for the patient or user have not been given time to adequately cool, injury (burn) to patient or user may occur.** 
9. To secure Upper Arm to Traction Tower, place buckles of Biceps Strap in channels on Base Plate until they are next to either side of the patient's arm. Adjust strap length until patient's arm is secured to Base Plate. Ensure foam on Biceps Strap is centered over the patient's arm.
10. Place the patient's forearm against the recess in the Upper Tower. The height of the Upper Tower can be adjusted depending on patient's forearm length.
11. Insert Forearm Strap buckle bracket into Upper Tower and wrap the Forearm Strap around the patient's forearm and lock the buckle into the bracket. Adjust strap length until patient's forearm is secured to Upper Tower.  
**WARNING: Failure to ensure user and patient are free of pinch points during setup may lead to injury.**  
**Do not apply excessive traction on the limb or overtighten straps during surgery, it may result in an injury to the patient.** 
12. Place Digit Traps on appropriate fingers to obtain the desired traction. Digit Traps should extend beyond the proximal interphalangeal joints, to prevent slippage. Adjust cable tie on Digit Trap to lock at base of the finger.  
**WARNING: Overtightening of Digit Traps may lead to patient injury. If overtightening occurs, readjust tightness to prevent injury to patient fingers.** 
13. Place ball chain on Digit Traps in recesses of the Scale Rack to attach Digit Traps to Traction Tower. Ensure ball chain securely rests in scale rack detent.  
**WARNING: Avoid grabbing the threaded rod while attaching Digit Traps to Scale Rack (possible pinch point).** 
14. Traction can now be adjusted to achieve desired distraction. Macro adjustment can be made at the Lower Arm to Upper Arm connection. Micro adjustment can be made by adjusting the knob on the Upper Arm. A reference to the amount of traction will be displayed on the Scale.  
**WARNING: Scale is for reference only, it is to be used as a reference when increasing traction or decreasing traction on patient anatomy, as well an indication that traction is maintained.**  
**Ensure traction is adequately adjusted and not excessive for each use.**  
**WARNING: Traction applied above the 20 lb reference will not be registered on the scale.**   
**WARNING: If patient's fingers do not fit in Digit Trap or remain secure in Digit Trap while under traction, alternative methods for surgery should be pursued, or injury to patient may occur.** 
15. Once traction is applied, ensure patient's fingers remain securely in Digit Traps, before making incisions.  
**WARNING: Ensure the user maintains the position of the Lower Arm prior to depressing the button on the Rotation Joint. Ensure the button lifts into the locked position before releasing the Lower Arm. Failure to do so may result in the Traction Tower Arms falling, causing user/patient injury.** 
16. The wrist may be positioned to allow for wrist flexion, extension, ulnar or volar deviation by depressing the button on the Rotation Joint and adjusting the Upper and Lower Arm to the adequate position.  
**WARNING: Additional, the Upper and Lower Arm may be rotated about the vertical axis of the Rotation Joint to allow for access of instruments into the wrist or to necessitate C-arm imaging.** 

## DISINFECTION,CLEANING, AND STERILIZATION GUIDELINES

### Refer to Instrument Disinfection Guide for Disinfection instructions.

#### **CLEANING**

1. Follow universal precautions for protective apparel when handling and cleaning contaminated equipment.
2. It is recommended that devices are cleaned as soon as reasonably practical following use.

#### **PREPARATION FOR CLEANING**

3. Pre-clean device at point of use to prevent drying of bioburden and organic soil: Open all ports, wipe with moistened towel/sponge, rinse off gross debris using cool water, flush devices with cannulations with cool water then soak in water at 22°C-43°C/ 72°F-110°F to prevent coagulation of blood and proteins on the device. Temperature of pre-rinse water should not exceed 45°C (113°F).  
**NOTE: Saline solution is not recommended as it has a corrosive effect on certain metals.**
4. Maintain moisture during transport from point of use to decontamination area.

**MANUAL CLEANING**

5. Clean under running water to prevent aerosolization.
6. Rinse excess soil from device.
7. Each subassembly should be cleaned using a neutral pH cleaning solution.
8. Follow solution manufacturer's recommendations for preparation and use.
9. Using a soft bristle brush, apply detergent solution to all surfaces ensuring that hinged attachments are cleaned in both open and closed positions.

**NOTE: Clean cannulation and holes using an appropriate brush to ensure that full depth of the feature is reached.**

10. Rinse each subassembly under running tap water (minimum temperature 25°C/77°F) for a minimum of 30 seconds using a minimum of 6 liters of rinse water. Ensure running water passes through cannulation, and that blind holes are repeatedly filled and emptied.
11. Ensure each subassembly is visibly free of detergent or cleaning residues.

**FINAL RINSE**

12. Rinse each subassembly under running deionized water (minimum temperature 25°C/ 77°F) using a minimum volume of 6 liters of rinse water for a minimum of 30 seconds.

**INSPECTION**

13. Inspect for cleanliness. Ensure device is free of visible debris and/or soil. Failure to follow cleaning instructions may compromise the sterility of the device.

**DRYING**

14. Dry the device using a clean, lint free towel until visibly free of water droplets and/or condensation.

**AUTOMATED CLEANING INSTRUCTIONS**

1. Rinse each subassembly under running water (minimum temperature 25°C/77°F) for a minimum volume of 6 liters of water for a minimum of 30 seconds as required to remove all traces of blood and debris.
2. Load all subassembly into washer/sanitizer and operate according to machine manufacturer's instructions. It is recommended to use a neutral pH cleaner. Only use approved solutions for cleaning device according to the manufacturer's recommendations.
3. Run washer/sanitizer according to manufacturer's specifications. A dry cycle is recommended but not required.

**LUBRICANT**

1. If applicable, a commercial water-based lubricant should be used to minimize wear of movable parts.

**STEAM STERILIZATION**

1. Prepare device such that all surfaces have direct contact with steam.
2. The HandStand Traction Tower (HS-TT) may be processed in a pre-vacuum steam sterilization cycle, either wrapped or in the HandStand Sterilization Tray (HS-ST). For sterilization in the tray, consult guidelines provided with the specific sterilization tray.

Method	Cycle	Minimum Temperature	Exposure Time	Dry Time
Steam (individually wrapped or unwrapped within a tray)	Pre-Vacuum	273°F (134°C)	3 minutes	20 minutes
Steam (individually wrapped or unwrapped within a tray)	Pre-Vacuum	270°F (132°C)	3 minutes	20 minutes

**PACKAGING AND LABELING**

1. Upon initial use, this device should only be used if the original packaging and labeling are intact.
2. If packaging has been opened or altered, contact your local CONMED representative, or in the U.S. contact the Customer Service Department (1-866-426-6633).

## SYMBOLS

The following are descriptions of the graphical symbols used to label CONMED Corporation products:

	Caution
	Manufacturer
	Non Sterile
	Do not use if package is damaged
	Keep dry
<b>Rx ONLY</b>	Prescription only, Caution: Federal law restricts this device to sale by or on the order of a physician
<b>EC REP</b>	EC REP
<b>REF</b>	Catalog Number
<b>CE</b>	CE Mark

## CONTACT INFORMATION

For more information or a product demonstration, contact your regional CONMED sales representative, or call 1-866-426-6633 in the U.S.

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©2022 CONMED Corporation. All Rights Reserved.
**MDSS GmbH**
**CONMED UK Ltd**