Surgical Instrument Trays

ENGLISH	1
FRANÇAIS	6
ITALIANO	11
DEUTSCH	16
ESPAÑOL	21
ΕΛΛΗΝΙΚΑ	27
PORTUGUÊS	33
PORTUGUÊS DO BRASIL	39
NEDERLANDS	45
DANSK	50
SVENSKA	55
SUOMI	60
NORSK	65
简体中文	70
日本語	74
POLSKI	79
TÜRKÇE	85
ROMÂNĂ	90
العربية	95
РУССКИЙ	100
HRVATSKI	106
ČEŠTINA	111
MAGYAR	116
LIETUVIŲ K.	121
SLOVENČINA	127
SLOVENŠČINA	132
한국어	137

Surgical Instrument Trays

REF	Description
8538	Zone Specific® II Meniscal Repair System Surgical Instrument Tray
8744B	Cruciate Reconstruction Surgical Instrument Tray
9853	Shutt® Suture Punch System Surgical Instrument Tray
9879A	Linvatec Manual Surgical Instrument Tray
10053	UltraFix® Surgical Instrument Tray
C2607	Soft Tissue Anchoring System Surgical Instrument Tray
C6108	Revo® Surgical Instrument Tray
C6117	Mini-Revo® Surgical Instrument Tray
C6154	BioAnchor® Surgical Instrument Tray
C6178	Bio Instability Surgical Instrument Tray
C9055	SE™ Graft Tensioner System Surgical Instrument Tray
MA0001	Meniscus Arrow™ Surgical Instrument Tray
MAC0001	CrossBow® Meniscus Arrow™ Surgical Inserter Instrument Tray
PS8831	GraFix® Graft Preparation Table Surgical Instrument Tray

DESCRIPTION

The Surgical Instrument Tray is constructed of metal and/or plastic with perforations to facilitate sterilant penetration, evacuation and drying. The tray is designed to fit any standard autoclave and is manufactured from durable, biocompatible materials that are corrosion resistant and compatible with the environment of repeated steam sterilization.

INTENDED USE / INDICATIONS FOR USE

The Surgical Instrument Tray is used for loading surgical instruments in order to conveniently organize, transport and store the instruments between uses. The tray is intended for non-porous loads, e.g. hooks, graspers, forceps, scissors, knot pushers, probes, punches, taps, drivers, rasps, switching sticks, obturators, etc. The Surgical Instrument Tray design does not utilize or incorporate gaskets, filters, valves and/or seals to maintain a sterile barrier. The tray design does not incorporate latches that produce and hold a seal for maintenance of sterile contents. The Surgical Instrument Tray is not intended to maintain sterility. To maintain sterility of the general reusable surgical instruments stored in the Surgical Instrument Tray, the tray shall be enclosed in an FDA cleared sterilization wrap.

CONTRAINDICATIONS

1. The Surgical Instrument Tray design does not utilize or incorporate gaskets, filters, valves and/or seals to maintain a sterile barrier.
2. The tray design does not incorporate latches that produce and hold a seal for maintenance of sterile contents.
3. The Surgical Instrument Tray is not intended to maintain sterility.
4. To maintain sterility of the general reusable surgical instruments stored in the Surgical Instrument Tray, the tray shall be enclosed in an FDA cleared sterilization wrap.

WARNINGS

1. Some instruments loaded into the trays may be extremely sharp. Handle with care to avoid injury.
2. DO NOT use EO or gamma sterilization. Steam sterilize according to the parameters and instructions in this document only.
3. TRAYS ARE NOT DESIGNED FOR IN SITU CLEANING OF CONTENTS.
4. Remove all medical devices contained in the tray for cleaning of instruments and tray prior to sterilization and reuse.
5. DO NOT overload trays. Overloading may inhibit steam flow, cause excessive drying times, and make trays too heavy to safely handle. Follow AORN and hospital guidelines for tray loading.
6. Trays provide no level of sterility protection. An FDA cleared (U.S.) sterile wrap must be used according to the manufacturer's instructions to maintain sterility.

PRECAUTIONS

1. Prior to use, cleaning and sterilization of instruments, remove all protective packaging and tip protectors, when applicable, from instrumentation.
2. Inspect instruments and tray for damage, such as cracks, broken handles or hinges, or broken containment brackets prior to use to ensure they are in good physical condition and functioning properly.
3. Inspect tray for damage, such as a sharp edge, that could cause a tear or rupture in the wrap. Do not use tray if damaged.
4. Inspect instruments after sterilization to ensure there is no damage. Do not use instruments if damaged.
5. Avoid unintended contact with other surgical instruments during use to prevent damage or breakage.
6. DO NOT drop Surgical Instrument Tray or components.
7. The use of nonabsorbent tray liners can cause condensate to pool.
8. DO NOT use contents if the sterilization wrap used to wrap the tray is torn or damaged.
9. DO NOT use strong disinfecting solutions to sterilize the trays or instruments. Trays and instruments should be cleaned using neutral pH-balanced detergent only.
10. DO NOT stack Surgical Instrument Trays on top of one another. Be sure that ventilation holes are not obstructed. For effective sterilization, trays must have adequate steam penetration through all surfaces. They must also be placed upright on shelves for proper ventilation. Condensate can pool on non-absorbent surfaces.
11. DO NOT place trays on their sides or at vertical angles in chamber as it may prevent proper drainage.
12. Only legally marketed, FDA cleared sterilization wrap should be used to wrap the trays.

CLEANING and STERILIZATION GUIDELINES

Inspection and Cleaning of Surgical Instrument Tray

This Surgical Instrument Tray is a reusable medical device, that can be used indefinitely, if not damaged. Surgical Instrument Trays should be cleaned and inspected between uses.

1. Inspection should be conducted by trained personnel.
2. Inspect Surgical Instrument Tray to ensure it is free from any damage such as cracks, broken handles or hinges, or broken containment brackets.
3. Ensure all movable parts are working correctly (i.e. hinges, locking mechanisms). Instrument Tray may be used to hold instruments during sterilization. Trays do not provide any level of sterility protection. An FDA cleared (U.S.) sterile wrap must be used according to the manufacturer's instructions to maintain sterility. Follow manufacturer's instructions for selecting the proper size and thickness of the sterile wrap.

Cleaning the Surgical Instrument Tray

1. Follow universal precautions for protective apparel when handling and cleaning contaminated equipment.
2. Remove all products from the Surgical Instrument Tray prior to cleaning it.
3. Thoroughly clean the Surgical Instrument Tray by wiping with a clean, soft, nonabrasive cloth dampened with a neutral pH-balanced detergent.
4. Rinse tray under running water followed by another rinse with distilled water.
5. Visually inspect the tray to ensure that there is no sign of visible debris left on tray surfaces including but not limited to tray perforations, instrument brackets, and all moving parts such as handles and latches.
6. Repeat steps 3 through 5 as needed until the tray is visually clean.
7. Dry surfaces using a clean, lint free towel.

Cleaning Instruments Prior to Loading in Surgical Instrument Tray

1. Follow universal precautions for protective apparel when handling and cleaning contaminated equipment.
2. Remove instruments from the Surgical Instrument Tray for instrument cleaning.
3. Follow the manufacturer's instructions on how to adequately clean and dry devices prior to using steam sterilization.
4. Complex instruments should be prepared and sterilized according to the manufacturer's written instructions.

Loading the Surgical Instrument Tray

1. Load instruments inside tray. Brackets and/or silicone mats secure instruments in place during transportation and sterilization.
2. DO NOT stack Surgical Instrument Trays on top of one another. Be sure that ventilation holes are not obstructed. For effective sterilization, trays must have adequate steam penetration through all surfaces. They must also be placed upright on shelves for proper ventilation. Condensate can pool on non-absorbent surfaces.
3. DO NOT overload trays. Overloading may inhibit steam flow, cause excessive drying times, and make trays too heavy to safely handle. Follow AORN and hospital guidelines for tray loading.
4. Only one tray placed at the center of the chamber can be processed at a time in an ST8 compliant sterilizer.
5. Small baskets, trays, or other accessories with covers or lids should only be used in trays specifically designed and labeled for the purpose.

STEAM STERILIZATION

Sterilize the Zone Specific® II Meniscal Repair System Surgical Instrument Tray (8538) using the following parameters:

Method	Cycle	Minimum Temperature	Minimum Exposure	Dry Time
Steam (wrapped)	Pre-vacuum	270°F (132°C)	4 minutes	25 minutes
Steam (wrapped)	Gravity	270°F (132°C)	15 minutes	20 minutes

Sterilize the Cruciate Reconstruction Surgical Instrument Tray (8744B) using the following parameters:

Method	Cycle	Minimum Temperature	Minimum Exposure	Dry Time
Steam (wrapped)	Pre-vacuum	270°F (132°C)	3 minutes	30 minutes
Steam (wrapped)	Gravity	270°F (132°C)	10 minutes	30 minutes
Steam (wrapped)	Gravity	250°F (121°C)	50 minutes	30 minutes

Sterilize the Shut® Suture Punch System Surgical Instrument Tray (9853) using the following parameters:

Method	Cycle	Minimum Temperature	Minimum Exposure	Dry Time
Steam (wrapped)	Pre-vacuum	270°F (132°C)	4 minutes	35 minutes*
Steam (wrapped)	Gravity	270°F (132°C)	15 minutes	35 minutes*

Sterilize the Linvatec Manual Surgical Instrument Tray (9879A) using the following parameters:

Method	Cycle	Minimum Temperature	Minimum Exposure	Dry Time
Steam (wrapped)	Pre-vacuum	270°F (132°C)	4 minutes	12 minutes
Steam (wrapped)	Gravity	270°F (132°C)	10 minutes	12 minutes
Steam (wrapped)	Gravity	250°F (121°C)	35 minutes	12 minutes

Sterilize the UltraFix® Surgical Instrument Tray (10053) using the following parameters:

Method	Cycle	Minimum Temperature	Minimum Exposure	Dry Time
Steam (wrapped)	Pre-vacuum	270°F (132°C)	4 minutes	20 minutes
Steam (wrapped)	Gravity	270°F (132°C)	15 minutes	30 minutes

Sterilize the Soft Tissue Anchoring System Surgical Instrument Tray (C2607) using the following parameters:

Method	Cycle	Minimum Temperature	Minimum Exposure	Dry Time
Steam (wrapped)	Pre-vacuum	270°F (132°C)	4 minutes	20 minutes
Steam (wrapped)	Gravity	270°F (132°C)	15 minutes	30 minutes

Sterilize the Revo® Surgical Instrument Tray (C6108) using the following parameters:

Method	Cycle	Minimum Temperature	Minimum Exposure	Dry Time
Steam (wrapped)	Pre-vacuum	270°F (132°C)	4 minutes	35 minutes*
Steam (wrapped)	Gravity	270°F (132°C)	15 minutes	40 minutes

Sterilize the Mini-Revo® Surgical Instrument Tray (C6117) using the following parameters:

Method	Cycle	Minimum Temperature	Minimum Exposure	Dry Time
Steam (wrapped)	Pre-vacuum	270°F (132°C)	4 minutes	35 minutes*
Steam (wrapped)	Gravity	270°F (132°C)	15 minutes	30 minutes

Sterilize the BioAnchor® Surgical Instrument Tray (C6154) using the following parameters:

Method	Cycle	Minimum Temperature	Minimum Exposure	Dry Time
Steam (wrapped)	Pre-vacuum	270°F (132°C)	4 minutes	30 minutes
Steam (wrapped)	Gravity	270°F (132°C)	15 minutes	30 minutes

Sterilize the Bio Instability Surgical Instrument Tray (C6178) using the following parameters:

Method	Cycle	Minimum Temperature	Minimum Exposure	Dry Time
Steam (wrapped)	Pre-vacuum	270°F (132°C)	4 minutes	20 minutes
Steam (wrapped)	Gravity	270°F (132°C)	30 minutes	20 minutes
Steam (wrapped)	Gravity	250°F (121°C)	60 minutes	20 minutes

Sterilize the SE™ Graft Tensioner System Surgical Instrument Tray (C9055) using the following parameters:

Method	Cycle	Minimum Temperature	Minimum Exposure	Dry Time
Steam (wrapped)	Pre-vacuum	270°F (132°C)	4 minutes	20 minutes
Steam (wrapped)	Gravity	270°F (132°C)	15 minutes	20 minutes

Sterilize the Meniscus Arrow™ Surgical Instrument Tray (MA0001) using the following parameters:

Method	Cycle	Minimum Temperature	Minimum Exposure	Dry Time
Steam (wrapped)	Pre-vacuum	270°F (132°C)	4 minutes	20 minutes
Steam (wrapped)	Gravity	270°F (132°C)	15 minutes	25 minutes

Sterilize the CrossBow® Meniscus Arrow™ Surgical Instrument Tray (MAC0001) using the following parameters:

Method	Cycle	Minimum Temperature	Minimum Exposure	Dry Time
Steam (wrapped)	Pre-vacuum	270°F (132°C)	4 minutes	35 minutes*
Steam (wrapped)	Gravity	270°F (132°C)	15 minutes	45 minutes*

Sterilize the GraFix® Graft Preparation Table Surgical Instrument Tray (PS8831) using the following parameters:

Method	Cycle	Minimum Temperature	Minimum Exposure	Dry Time
Steam (wrapped)	Pre-vacuum	270°F (132°C)	4 minutes	30 minutes

*** PRECAUTION: Testing demonstrates that a minimum dry time of 35 minutes is required to prevent wet packs when using the Prevacuum cycle and a minimum 45 minutes is required to prevent wet packs when using the gravity cycle.**

CONTACT INFORMATION

For more information or a product demonstration, contact your regional ConMed Linvatec sales representative, or call 1-866-426-6633 within the U.S.

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All rights reserved.
MDSS GmbH
CONMED UK Ltd
ConMed Corporation
©2022 ConMed Corporation.