GE Healthcare

MobileLink™ Wireless Communication Installation Manual For GEH-BR-4600 Wireless Bridge 2053535-080 C

MobileLink English © 2017-2018 General Electric Company. All Rights Reserved.

Publications Information

NOTE:

The information in this manual applies only to the following MobileLink product configurations:

MAC 2000 Resting ECG Analysis System:

Software version 1.1 or higher

MAC 800 Resting ECG Analysis System:

Software version 1.0.X and 2.0.X

Servers:

Silex Wi-Fi Bridge GEH-BR-4600

MUSE Systems Software Versions:

7.1.1 or higher

8.0.1 or higher

9.0 or higher

The information is provided for reference to support legacy systems with this MobileLink configuration.

MUSE, MAC, and MobileLink are trademarks owned by GE Medical Systems *Information Technologies*, Inc., a General Electric company, going to market as GE Healthcare. All other marks are owned by their respective owners.

The document part number and revision are on each page of the document. The revision identifies the document's update level. The revision history of this document is summarized in the following table.

Revision	Date	Comment
А	8 August 2017	Initial release of GEH-BR-4600 Server only manual
В	6 November 2017	Updated to phase in Silex wireless module.
С	24 May 2018	Update according to the service comments.

To access other GE Healthcare Diagnostic Cardiology documents, go to the Common Documentation Library (CDL), located at www.gehealthcare.com/documents, and click *Cardiology*.

To access Original Equipment Manufacturer (OEM) documents, go to the device manufacturer's website.

NOTE:

All illustrations in this document are provided as examples only. Depending on system configuration, screens in the document may differ from the screens on your system.

All patient names and data are fictitious. Any similarity to actual persons is coincidental.

If you require additional assistance, contact your GE Healthcare representative, or GE Healthcare support at one of the following numbers:

- North America: 1-800-558-7044
- Europe: +49 761 45 43 -0
- Asia: +86 21 3877 7888

Service Manual Language Information

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Contents

1 Introduction

	Manual Information Manual Purpose Additional Assistance Customer Responsibility. On-Line Center Support	 13 13 13 13 13
	Before You Begin Silex GEH-BR-4600 Introduction System Requirements Print the Configuration Capture File Determining Appropriate Installation Procedures	 14 14 14 15 16
	Connectors and Indicators Silex Wireless Bridge GEH-BR-4600WAN2-03-XX	 16
2	Mounting the GEH-BR-4600 Wireless Bridge	
	Mounting the GEH-BR-4600 Wireless Bridge on the MAC 2000 Device	20
	Mounting the GEH-BR-4600 Wireless Bridge on the MAC 800 Device	22
3	Configuring the Silex Wireless Bridge GEH-BR-4600 for the MAC 2000 and MAC 800	2
	Important	25
	Configuring Silex Wireless Bridge GEH-BR-4600 for Wireless LAN	25
	Configuring Silex Wireless Bridge GEH-BR-4600 to Single Client Mode	28
	Configuring Silex Wireless Bridge GEH-BR-4600 to Single Client Mode Configuring Silex Web Page with Password	28 29

4 Troubleshooting

Α

Silex Wireless Bridge GEH-BR-4600	
Using the Serial/IP Tools to Troubleshoot Port Activity	42
Trace Window	
Technical Specifications	
Silex Wireless Bridge GEH-BR-4600	45

1

Introduction

MobileLink is a wireless, integrated ECG communication solution that helps you capture, transmit, and analyze your patient's ECG information quickly and efficiently. MobileLink combines all available security protocols into one superset to guard your facility from the risk of transmission interception.

This chapter provides general information required for the proper use of the system and this manual. Familiarize yourself with this information before using the system.

Manual Information

This section provides information for the correct use of this manual.

Keep this manual with the equipment at all times and periodically review it. You should request training assistance from GE Healthcare, if needed.

Manual Purpose

The purpose of this manual is to provide instruction on how to install MobileLink Wireless Communication to a MAC 2000/MAC 800 using a Silex Wireless Bridge.

Additional Assistance

GE Healthcare maintains a trained staff of application and technical experts to answer questions and respond to issues and problems that may arise during the installation, maintenance, and use of this system.

Contact your local GE Healthcare representative to request additional assistance.

Customer Responsibility

GE Healthcare Service is responsible for the initial installation and configuration of the MobileLink wireless communication for ECG. Any changes made to configuration, security, and so on by the customer after the initial installation are the responsibility of the customer.

On-Line Center Support

Contact the GE Healthcare On-Line Center if you have any questions or problems while performing this installation.

Before You Begin

Before installing the MobileLink option, confirm that your site meets the following system requirements:

- Confirm that your site meets the system requirements for MobileLink. See "System Requirements" on page 14.
- Review the different procedures to install MobileLink and determine which is appropriate for your site. See "Determining Appropriate Installation Procedures" on page 16.

Silex GEH-BR-4600 Introduction

The GEH-BR-4600 bridge is the wireless bridge which allows the MAC 800/MAC 2000 LAN port to be used as a wireless device. A high-performance wireless network adapter which enables the Ethernet LAN port to become an 802.11a/b/g/n dual band wireless network device.

GE Healthcare Part Number	Manufacturing Part Number	Description
2095170-001	GEH-BR-4600WAN2-03-US	SILEX WI-FI BRIDGE KIT MAC2000 US
2098761-001	GEH-BR-4600WAN2-03-WW	MAC 2000 SILEX WIRELESS BRIDGE KIT- GLOBAL VERSION

Key Features:

- Communicate in the 2.4 GHz and 5 GHz frequencies
- Wired LAN port configured for 10base-T-Full
- Uses the MAC 800/2000 MAC address for wireless communications (MAC address transparency feature)
- Configuring the bridge for Single Client Mode disables the wireless bridge functionality when the GEH-BR-4600 bridge is connected to the wired LAN port of another MAC 800/MAC 2000 or another PC (security feature)

System Requirements

Review and confirm that your site meets the following system requirements.

- 1. The MAC system software requirements:
 - MAC 2000 system must be at v1.1 or higher
 - MAC 800 system must be at 1.0.X or 2.0.X
- 2. The MAC system must have the appropriate option activated.
 - a. From the *Main Menu*, press *System Configuration > More > More > Option Setup*.
 - b. Verify that the *LANM* or *LANC* option is activated.

MobileLink[™] Wireless Communication

- 3. A dynamic IP, static IP, or host name must be provided for each MAC system to be configured with the LAN option.
- 4. MUSE system software requirements:
 - v7:7.1.1
 - v8: 8.0.1
 - v9: 9.0.0 (only available for MAC 2000)
 - CardioSoft-CS v6.71 or higher and the connectivity server option

Print the Configuration Capture File

Confirm all the information about this system has been accurately entered into the *Configuration Capture File (2002783-067)*. Prior to this installation, you may have received a copy of this file from Project Management or ND&I with information gathered during pre-installation. If so, use that copy of the *Configuration Capture File*. If you have not received a copy of the *Configuration Capture File*, use the blank version of the file on the MobileLink CD.

Print a copy as a reference to use during the installation. During installation, you will add information to the file and then store the completed file on the MUSE file server.

MobileLink Configuration

The IP Address, Subnet Mask, and Gateway Address If this information cannot be acquired ahead of t After completion, this file should be saved on th

Site Name: IS Contact Name: IS Contact Phone: Notes:

Configuration Information

IP Address: Subnet Mask: Gateway Address: SSID: Station Name: MAC5K001 Port Number: 3001 MAC Address: Client Bridge Serial Number: Encryption Information:

Configuration Capture File Provides data entry for up to 30 MAC systems.

NOTE:

The *Station Name* is a general field that is used to identify a cart. This name can be different from or the same as the System Name, if hostname addressing is used.

NOTE:

Station name and Port Number have been pre-populated in this file with the following suggested convention. The Station Name and Port Numbering system is a suggested convention.

Station Name	Port Number
MAC2K001	3001
MAC2K002	3002

Station Name	Port Number
MAC2K003	3003
etc.	etc.

Determining Appropriate Installation Procedures

These instructions are written for several types of installations, including one type of serial server and two kinds of trolleys. When installing and configuring the MobileLink option, use the section of this manual which is appropriate to the hardware you are installing.

- 1. Use Chapter 3 "Configuring the Silex Wireless Bridge GEH-BR-4600 for the MAC 2000 and MAC 800" on page 25.
- 2. Chapter 2 "Mounting the GEH-BR-4600 Wireless Bridge" on page 19.

See "Mounting the GEH-BR-4600 Wireless Bridge on the MAC 2000 Device" on page 20.

See "Mounting the GEH-BR-4600 Wireless Bridge on the MAC 800 Device" on page 22.

3. Configure the System as described in Chapter 3 "Configuring the Silex Wireless Bridge GEH-BR-4600 for the MAC 2000 and MAC 800" on page 25.

Connectors and Indicators

Silex Wireless Bridge GEH-BR-4600WAN2-03-XX

Silex Wireless Bridge GEH-BR-4600

FRONT

ltem	Name	Description
1	Push switch	Starts this product in Configuration Mode:
		Press and hold the push switch, release the push switch when the WLAN LED and STATUS LED start to blink together. It may take up to 5 seconds for the 2 LEDS to start blinking together.
		Resets this product to factory default settings:
		Press and hold the push switch while turning on this product. Release the switch when the LAN Link LED turns to the Orange and then Green.
2	LAN Status LED (Yellow)	BLINK: Waiting for connection when Link LED is OFF. Transferring data when Link LED is ON.
3	LAN port	Connects to the network cable.
4	LAN Link LED (Green/Orange)	ON (Green): Connected to 10BASE-T/100BASE-TX network
5	Power connector	Connect USB power cable

BOTTOM

Item	Name	Description
6	POWER LED (Orange)	Turns ON when this product is powered on.
7	WLAN LED (Green)	 ON: Running is Infrastructure mode. OFF: Running is Ad hoc mode. * Blink together with STATUS LED when running in Configuration Mode.

Item	Name	Description
8	STATUS LED (Green)	ON : A network connection is established.
		BLINK: Transferring data.
		* Blink together with WLAN LED when running in Configuration Mode .
9	Ethernet address	The Ethernet (MAC) Address of this product is found on the bottom.
		The last 6 digits are the product's serial number.
		Example: If the Ethernet Address is 00:80:92:00:11:22 (noted as 008092001122) then the product's serial number is 001122.
		NOTE: DHCP Reservation by MAC address - Use the MAC 800/MAC 2000 MAC address, not the GEH-BR-4600 MAC address.

2 Mounting the GEH-BR-4600 Wireless Bridge

Use the following instructions for mounting a GEH-BR-4600 wireless bridge.

Mounting the GEH-BR-4600 Wireless Bridge on the MAC 2000 Device

1. Stick a pair of dual locks on the GEH-BR-4600 wireless bridge as shown in the following graphic.

2. Put another pair of dual locks on the dual locks sticked on the Silex wireless bridge and press firmly.

3. Connect the Silex cable to the Silex wireless bridge using the circular connector.

4. Connect the Silex cable to the MAC 2000 device using the USB connector.

- 5. Insert the internet cable to the network connectors of both the Silex wireless bridge and MAC 2000 device.
- 6. Remove the stick cover on the dual locks and place the Silex wireless bridge on the MAC 2000 device as shown in the following graphic.

Mounting the GEH-BR-4600 Wireless Bridge on the MAC 800 Device

1. Select the **Dual Locks** from the **Silex Wireless Bridge Kit** and engage the hook and loop surface of the **Dual Locks** against each other to form a pair.

2. Remove the adhesive liner from one side of the dual lock and adhere each pair to the flat surface of Silex wireless bridge module.

3. Remove the adhesive liner and carefully position the Silex Module on the MAC 800 unit by placing it to the rear of the printer side. The LAN port and power switch should be visible from the rear of the MAC 800 device.

4. Route and connect the Ethernet and USB power cables from the Silex wireless module to the MAC 800 unit.

Mounting the GEH-BR-4600 Wireless Bridge

3 Configuring the Silex Wireless Bridge GEH-BR-4600 for the MAC 2000 and MAC 800

Important

Before you begin this procedure, check for any updates to the procedure that may have been released in a Technical Bulletin.

NOTE:

Before configuring the wireless bridge, be sure you have performed the steps in "Before You Begin" on page 14.

Configuring Silex Wireless Bridge GEH-BR-4600 for Wireless LAN

- 1. Configure PC LAN port for DHCP.
- 2. Connect the Silex wireless bridge to a PC using the standard CAT 5 network cable that came with the product.
- 3. Take the USB power cable and connect the DC connector side to Silex wireless bridge and USB connector side to the same PC.

4. When the Power LED at the top turns **ON**, press and hold the **Push** switch at the front.

5. Release the **Push** switch when **WLAN** LED and **STATUS** LED starts to blink together.

It may take 20 seconds for the 2 LEDs to blink together.

The Silex wireless bridge is now set to **Configuration Mode**.

6. Use the same PC to connect to the bridge for configuring the product. Disable any wireless LAN (if enabled) on the PC.

7. Start a Web Browser on the PC you are using for the setup.

The web page of this product is <u>silex.com</u>. If the webpage is not displayed, enter silex.com in the address bar of the browser and press the **Enter** key.

silex technology	Wireless LAN Configuration			
Select Language English	Easy Configuration Detailed Configuration Smart Wireless Setup The following is the procedure to configure the wireless LAN settings.			R
Configuration – Advanced Configuration – Certificate Configuration – Password Configuration	 Select a wireless network. Wireless Network List Search Wireless Networks 			
Maintenance Reboot Factory Default Configuration Firmware Update Information	Belect SSD df-wrt BLUESSO internet	RSS1 00000 00000 00000	CH Wireless Standard 6 11s 1 11ns 1 11ns	Wireless Mode
- License Information	Enter the Network Key. Network Key			
	3. Save the configuration and reboot this product.			

NOTE:

If the hospital wireless network is not listed on this page, select the **Detailed Configuration** tab.

8. The *Detailed Configuration* window opens.

silex technology	Uireless LAN Configurati	ion
Select Language	Easy Configuration Detailed Configuration	ton Smart Wireless Setup
System Configuration	Operating Mode	Infra. V
- Wireless LAN Configuration - Advanced	Wireless Standard	AUTO
- Certificate Configuration	SSID	test
Configuration	SSID Filter	OFF V
Maintenance Reboot Factory Default Castline maties	Network Authentication	WPA2 V
- Firmware Update	EEE802.1X Authentication	OFF V
 Information License Information 	WPA/WPA2 Configuration	
	Name	Value
	Encryption Mode	AES V
	Pre-Shared Key	•••••
		Submit Reset

- Select *Wireless Standard* and choose 2.4GHZ or 5GHZ if a specific radio frequency is required.
- Enter the **SSID**
- Go to the "Silex wireless Bridge GEH-BR-4600 Encryption/Authentication Modes" on page 30 for the encryption and authentication settings.
- Select Submit after the authentication and encryption settings have been configured.

NOTE:

After *Submit* is selected and the changes activated, the Silex Bridge leaves the configuration mode.

The WLAN and STATUS LEDs do not blink any more.

Configuring Silex Wireless Bridge GEH-BR-4600 to Single Client Mode

CAUTION:

To prevent unauthorized access to the wireless bridge, configure the bridge for single client mode using the MAC 800/MAC 2000 MAC address. Configuring the bridge for Single Client Mode disables the wireless bridge functionality when the GEH-BR-4600 bridge is connected to the wired LAN port of another MAC 800/2000 or another PC (security feature).

- 1. Record the MAC address of the device.
 - a. On the device, go to *Main Menu* > *System Configuration* > *Service Setup*. The password is 7763.
 - b. Select System Diagnostic and press Enter.
 - c. Select Service Report and press Enter.
 - d. The service report lists the MAC address.

- 2. In the Silex configuration page, click *Advanced Configuration*.
- 3. Select 10BASE-T-FULL for LAN Interface field.

Silex technology	Advanced Configurati	on.
Select Language English	Advanced Configuration	
▼ System Configuration - Wireless LAN Configuration - Advanced	Configure the parameters for "Adv Enter the values and click "Submit	anced Configuration". ".
Configuration	Advanced Configuration	
- Certificate	Name	Value
- Password Configuration	LAN Interface	10BASE-T-Full
▼ Maintenance - Reboot - Factory Default	Network Device Address	
Configuration – Firmware Update	Client Mode	Single Client Mode 💙
▼ Information		
- License Information	TCP/IP Configuration	
	Name	Value
	DHCP	DISABLE 🗸
	IP Address	0.0.0.0
	Subnet Mask	0.0.0
	Default Gateway	0.0.0.0

Type the MAC 800/MAC 2000 MAC address in the *Network Device address* field.
 NOTE:

The MAC Address needs to be entered without the hyphen.

- 5. Keep all other fields as default values.
- 6. Click *Submit* to submit the configuration changes.

Configuring Silex Web Page with Password

CAUTION:

To prevent unauthorized change of patient data, enable password while accessing configuration page.

1. In Silex Configuration mode, click *Password configuration*.

silex	II Password Configuration	C
English V	Password Carlieuration	
Contens Configuration Hindex LAN Configuration Enderstand Configuration Configuration	Configure the password used to access the Web page B a password is set, only users with the password will be able to modify the configuration Be predict one to target the password, often will be out into breach be able to configure this protect For any other the new password and citics "Dubries". New Password Confism New Password	Gpictur ^
silex SX-BR-4600WAN		~

2. Type a password in the *New Password* field.

NOTE:

The default password is access. It is recommended that the default password is used.

- 3. Type the same password in the **Confirm New Password** field.
- 4. Click Submit.

The password is set to access the Silex web page. If a password is set, only users with the password will be able to modify the product configuration. Be careful not to forget the password or you will not be able to configure the product.

Silex wireless Bridge GEH-BR-4600 Encryption/Authentication Modes

GEH-BR-4600WAN2-03-XX supports the following security functions:

- WEP (64 Bit/128 Bit)
- WPA/WPA2 Mix Mode (TKIP/AES)
- WPA2-PSK (AES)
- IEEE 802.1X PEAP, EAP-TLS, EAP-TTLS, EAP-FAST, LEAP

Silex Encryption/Authentication Matrix

Authentication	Encryption		
	WEP	WPA/WPA2 Mix Mode	WPA2
Open	supported See "Configuring WEP/Open Authentication" on page 31.	not supported	not supported

Silex Encryption/Authentication Matrix (cont'd.)

PSK	supported See "Configuring WEP/Pre-shared Key" on page 32.	supported See "Configuring WPA/WPA2 Mix Mode Encryption and PSK Authentication" on page 33.	supported See "Configuring WPA2 Encryption and PSK Authentication" on page 33.
LEAP	not supported	supported See "Configuring EAP with LEAP Authentication" on page 35.	supported See "Configuring EAP with LEAP Authentication" on page 35.
PEAP	not supported	supported See "Configuring EAP with PEAP Authentication" on page 34.	supported See "Configuring EAP with PEAP Authentication" on page 34.
EAP-TLS	not supported	supported See "Configuring EAP with TLS Authentication" on page 36.	supported See "Configuring EAP with TLS Authentication" on page 36.
EAP-TTLS	not supported	supported See "Configuring EAP with TTLS Authentication" on page 37.	supported See "Configuring EAP with TTLS Authentication" on page 37.
EAP-FAST	not supported	supported See "Configuring EAP-FAST Authentication" on page 38.	supported See "Configuring EAP-FAST Authentication" on page 38.

Configuring WEP/Open Authentication

1. Type or select the information in the *Detailed Configuration* tab using the information in the following table:

WEP/Open Authentication Selections

Item	Selection
Network	Select Open for WEP.
Autnentication	• Open (Open System) Allow all accesses without using any authentication. This mode uses WEP (64bit/128bit) to encrypt the communication. You can disable this encryption when it is not necessary.
WEP	Select ON or OFF .

WEP/Open Authentication Selections (cont'd.)

Item	Selection
Key Index	Select 1, 2, 3, or 4.
WEP Key	Enter the WEP key of this product.
	If the key format is hexadecimal, enter 10 or 26 hexadecimal digits.
	If the key format is alphanumeric character, enter 5 or 13 character.

NOTE:

Do not change the other settings in this window.

2. Click **Submit**.

Configuring WEP/Pre-shared Key

1. Type or select the information in the *Detailed Configuration* tab using the information in the following table:

WEP/Pre-shared Key Selections

Item	Selection
Network	Select Shared for WEP.
Authentication	• Shared (Pre-Shared Key) Allow access only to network devices with the same WEP Key by using it as an authentication key. This mode uses WEP (64bit/128bit) to encrypt the communication. You can disable this encryption when it is not necessary.
WEP:	Select ON .
Key Index:	Select 1, 2, 3, or 4.
WEP Key:	Enter the WEP key of this product.
	If the key format is hexadecimal, enter 10 or 26 hexadecimal digits.
	If the key format is alphanumberic character, enter 5 or 13 character.

NOTE:

Do not change the other settings in this window.

2. Click **Submit**.

Configuring WPA2 Encryption and PSK Authentication

1. Type or select the information in the *Detailed Configuration* tab using the information in the following table:

Item	Selection
Network Authentication	Select WPA2 .
IEEE802.1x Authentication	Select OFF .
Encryption Mode:	Select AES .
Pre-shared Key	Type the pre-shared key or passphrase.
	NOTE: The key value or passphrase entered here is used to initialize the session with the access point. If a key value is entered, it must be exactly 64 hex characters. A passphrase must be 8 to 63 displayable characters.

WPA2 Encryption and PSK Authentication Selections

NOTE:

Do not change the other settings in this window.

2. Click **Submit**.

Configuring WPA/WPA2 Mix Mode Encryption and PSK Authentication

1. Type or select the information in the *Detailed Configuration* tab using the information in the following table:

Item	Selection
Network Authentication	Select WPA/WPA2 Mix Mode.
IEEE802.1x	Select OFF .

WPA/WPA2 Mix Mode Encryption and PSK Authentication Selections

IEEE802.1x Authentication	Select OFF .
Encryption Mode:	Select AES or Auto .
Pre-shared Key	Type the pre-shared key or passphrase.
	NOTE: The key value or passphrase entered here is used to initialize the session with the access point. If a key value is entered, it must be exactly 64 hex characters. A passphrase must be 8 to 63 displayable characters.

NOTE:

Do not change the other settings in this window.

2. Click **Submit**.

Configuring EAP with PEAP Authentication

NOTE:

This configuration requires Silex firmware version 1.1.1 or higher.

1. Type or select information in the *Detailed Configuration* tab using the information in the following table:

EAP with PEAP Authentication Selections

Item	Selection
Network Authentication	Select WPA2 or WPA/WPA2 Mixed Mode for PEAP. • WPA2
	• WPA/WPA2 Mixed Mode This mode uses network authentication by PSK and encrypts the communication using TKIP/AES. Please note that WPA/WPA2 can be used only when this product is running in Infra mode.
IEEE802.1x Authentication	Select ON .
IEEE802.1x Authentication Type	Select PEAP .
Encryption Mode:	This parameter is determined by the option selected in the detailed Configuration page. If WPA2 is selected the only option is AES , if WPA/WPA2 Mixed mode is selected AES and Auto is selected. The default is AES .
Inner Authentication Method:	This defines the authentication method used for the inner authentication of the bridge. Select MSCHAPv2 .
Server Certificate Verification	Select ON or OFF . This provides an option not to use CA Certificates during the authentication process. Disabling certificate verification weakens the authentication process and is not recommended.
IEEE802.1x User Name	This is the user name defined on the RADIUS server account being used by the target network. This is an ASCII string.
Password	This is the password assigned to the user name defined on the RADIUS server account being used by the target network.
Network Device Address	Type the MAC address of the device. If using IEEE802.1X authentication, register the Ethernet address of the network device that will be connected to a wired LAN.

NOTE:

Do not change the other settings in this window.

If server certificate verification is ON, a CA Certificate is required for PEAP and must be loaded and registered for PEAP to work. See "Configuring Authentication CA Certificate" on page 40 for loading and registering certificates.

2. Click **Submit**.

Configuring EAP with LEAP Authentication

NOTE:

This configuration requires Silex firmware version 1.1.1 or higher.

1. Type or select information in the *Detailed Configuration* tab using the information in the following table:

EAP with LEAP Authentication Selections

Item	Selection
Network Authentication	 Select WPA2 or WPA/WPA2 Mixed Mode for LEAP. WPA2 WPA/WPA2 Mixed Mode This mode uses network authentication by PSK and encrypts the communication using TKIP/AES. Please note that WPA/WPA2 can be used only when this product is running in Infra mode.
IEEE802.1x Authentication	Select ON .
IEEE802.1x Authentication Type	Select <i>LEAP</i> .
Encryption Mode:	This parameter is determined by the option selected in the detailed Configuration page. If WPA2 is selected, the only option is AES , if WPA/WPA2 Mixed mode is selected, AES and Auto are selected. The default is AES .
IEEE802.1x User Name	This is the user name defined on the RADIUS server account being used by the target network. This is an ASCII string.
Password	This is the password assigned to the user name defined on the RADIUS server account being used by the target network.
Network Device Address	Type the MAC address of the device. If using IEEE802.1X authentication, register the Ethernet address of the network device that will be connected to a wired LAN. One device for Single Client Mode and up to 16 devices for Multi-Client Mode can be registered.

NOTE:

Do not change the other settings in this window.

No certificate is required for LEAP.

LEAP is a deprecated authentication protocol, use is discouraged due to security vulnerabilities.

2. Click **Submit**.

Configuring EAP with TLS Authentication

NOTE:

This configuration requires Silex firmware version 1.1.1 or higher.

1. Type or select information in the **Detailed Configuration** tab using the information in the following table:

EAP with TLS Authentication Selections

Item	Selection
Network	Select WPA2 or WPA/WPA2 Mixed Mode for TLS.
Authentication	• WPA2
	• WPA/WPA2 Mixed Mode This mode uses network authentication by PSK and encrypts the communication using TKIP/AES. Please note that WPA/WPA2 can be used only when this product is running in <i>Infra</i> mode.
IEEE802.1x Authentication	Select ON .
IEEE802.1x Authentication Type	Select EAP-TLS .
Encryption Mode:	This parameter is determined by the option selected in the detailed Configuration page. If WPA2 is selected, the only option is AES , if WPA/WPA2 Mixed mode is selected, AES and Auto are selected.
	The default is AES .
User name in IEEE802.1x	This is the user name defined on the RADIUS server account being used by the target network.
	This is an ASCII string.
Network Device	Type the MAC address of the device.
Address	If using IEEE802.1X authentication, register the Ethernet address of the network device that will be connected to a wired LAN.
	One device for Single Client Mode and up to 16 devices for Multi-Client Mode can be registered.

NOTE:

Do not change the other settings in this window.

A CA Certificate and Client Certificate are required for TLS and must be loaded and registered for EAP-TLS to work. See "Configuring Authentication CA Certificate" on page 40 and "Configuring Authentication Client Certificate" on page 40 for loading and registering certificates.

2. Click **Submit**.

Configuring EAP with TTLS Authentication

NOTE:

This configuration requires Silex firmware version 1.1.1 or higher.

1. Type or select information in the *Detailed Configuration* tab using the information in the following table:

EAP with TTLS Authentication Selections

Item	Selection
Network Authentication	Select WPA2 or WPA/WPA2 Mixed Mode for TTLS. • WPA2
	• WPA/WPA2 Mixed Mode This mode uses network authentication by PSK and encrypts the communication using TKIP/AES. Please note that WPA/WPA2 can be used only when this product is running in <i>Infra</i> mode.
IEEE802.1x Authentication	Select ON .
IEEE802.1x Authentication Type	Select EAP-TTLS .
Encryption Mode:	This parameter is determined by the option selected in the detailed Configuration page. If WPA2 is selected the only option is AES , if WPA/WPA2 Mixed mode is selected AES and Auto is selected. The default is AES .
Inner Authentication Method	This defines the authentication method used for the inner authentication of the bridge. The options include PAP , CHAP , MSCHAP and MSCHAPv2 .
Server Certificate Verification	Select ON or OFF . This provides an option not to use CA Certificates during the authentication process. Disabling certificate verification weakens the authentication process and is not recommended.
IEEE802.1x User Name	This is the user name defined on the RADIUS server account being used by the target network. This is an ASCII string.

EAP with TTLS Authentication Selections (cont'd.)

Item	Selection
Password	This is the password assigned to the user name defined on the RADIUS server account being used by the target network.
Network Device Address	Type the MAC address of the device. If using IEEE802.1X authentication, register the Ethernet address of the network device that will be connected to a wired LAN. One device for Single Client Mode and up to 16 devices for Multi-Client Mode can be registered

NOTE:

Do not change the other settings in this window.

NOTE:

If server certificate verification is ON, a CA Certificate is required for TTLS and must be loaded and registered for EAP-TTLS to work. See "Configuring Authentication CA Certificate" on page 40 for loading and registering certificates.

2. Click **Submit**.

Configuring EAP-FAST Authentication

1. Type or select information in the *Detailed Configuration* tab using the information in the following table:

EAP-FAST Authentication Selections

Item	Selection
Network Authentication	Select WPA2 or WPA/WPA2 <i>Mixed Mode</i> for EAP-FAST. • WPA2
	• WPA/WPA2 Mixed Mode This mode uses network authentication by PSK and encrypts the communication using TKIP/AES. Please note that WPA/WPA2 can be used only when this product is running in Infra mode.
IEEE802.1x Authentication	Select ON .
IEEE802.1x Authentication Type	Select EAP-FAST .
Encryption Mode:	This parameter is determined by the option selected in the detailed <i>Configuration</i> page. If <i>WPA2</i> is selected, the only option is <i>AES</i> , if <i>WPA/WPA2 Mixed mode</i> is selected, <i>AES</i> and <i>Auto</i> are selected. The default is <i>AES</i> .

EAP-FAST Authentication Selections (cont'd.)

Item	Selection
IEEE802.1x User Name	This is the user name defined on the RADIUS server account being used by the target network.
	This is an ASCII string.
Password	This is the password assigned to the user name defined on the RADIUS server account being used by the target network.
PAC Auto-distribution	EAP-FAST uses a file called a PAC file to distribute the necessary credentials for the network. EAP-FAST has an option for automatic credential distribution. It is mandatory to enable this option, so that the PAC file can be obtained from the network when initial configuration is done.
Network Device Address	Type the MAC address of the device. If using IEEE802.1X authentication, register the Ethernet address of the network device that will be connected to a wired LAN. One device for Single Client Mode and up to 16 devices for Multi-Client Mode can be registered.

NOTE:

Do not change the other settings in this window.

2. Click **Submit**.

System Configuration – Certificate Configuration

Enterprise security types (EAP types) often require the use of certificates as part of the authentication process. These certificates must be stored on the bridge for them to be used during the outer and inner authentication stages.

The PAC file is a Cisco format file that encapsulates all the necessary client side credentials for the bridge to authentication to an EAP-FAST configured network.

Setting	Description
Client Certificate	This section allows you to choose the client certificate required for the authentication of the bridge to the network. This certificate will be unique to the client. You must locate and select the certificate on the connected laptop/desktop, and then enter the password used to encrypt the certificate.
CA Certificate	This section allows you to choose the server side (Certificate Authority) certificate required for the authentication of the bridge to the network. This certificate will be unique to the authentication server being used by the network. You must locate and select the certificate on the connected laptop/desktop.
Certificate Registration	This section shows the currently loaded and registered certificates on the bridge.

The following settings allow use of EAP-TLS:

The files will not be saved to the bridge until the *Submit* button has been pressed.

Configuring Authentication CA Certificate

NOTE:

This is required for PEAP, TLS or TTLS authentication.

Steps for the Hospital IT Department

- 1. Identify the authorized Server Certificate Authority recognized by the radius server.
- 2. Generate CA certification from radius server with **.PEM** format.
- 3. Provide the **.PEM** file generated in the previous step to the GE Healthcare field engineer.

Steps for the GE Healthcare Field Engineer

- 1. From the *Silex Network Security* window, click on *Configure authentication server certificate*.
- 2. Click *Browse* and locate the exported *.PEM* file exported.
- 3. Click *Submit* and click the link to reset the device.
- 4. Verify that the *root key* was loaded on the *Certificate Root Key* field.

Configuring Authentication Client Certificate

NOTE:

This is required for TLS authentication.

Steps for the Hospital IT Department

- 1. Identify the User/Client Certificate being recognized by the RADIUS server.
- 2. Generate client certification with **.P12** format in radius server.

Make sure to select the option to export the private key. This requires a passphrase to be entered. Record the passphrase.

3. Provide the **.P12** file saved in the previous step to the GE Healthcare field engineer.

Steps for the GE Healthcare Field Engineer

- 1. From the *Silex Network Security* window, click on *Configure Private Key*.
- 2. Type the passphrase entered during the export process and browse for either the private key or client certificate .PEM files and click **Submit**.
- 3. Reset the device and repeat the procedure for the private key .P12 file.

Troubleshooting

Silex Wireless Bridge GEH-BR-4600

Symptom	Possible Cause	Remedy
The MAC device is connected to wireless bridge but unable to communicate	The wireless bridge is running in Configuration mode instead of Infrastructure Mode	Check the status of WLAN LED and STATUS LED. LEDs should blink together, if the wireless bridge is running in Configuration mode . In such case, restart the wireless bridge and check that the LEDS do not blink anymore.
	Problem with the IP address attained by MAC device OR duplicate IP Address is assigned to device as it would cause a conflict if another device on the network has the same address.	Verify if the RJ45 cable is functioning properly OR connected properly to Ethernet port of MAC device. Go to System Configuration > Communication Setup . Check the Wired LAN settings page. Verify if the MAC device has obtained a valid IP.
	The wireless LAN settings on the wireless bridge are not configured properly.	Obtain proper information of wireless settings and ensure that configured settings are correct. If settings are incorrect, disconnect the wireless bridge from MAC device. Connect it to PC and access the configuration web page. Configure the wireless bridge with valid settings again.

Symptom	Possible Cause	Remedy	
	The Settings on the MAC device are not configured properly.	Disconnect the wireless bridge from MAC device. Connect MAC device to wired LAN and check if communication works. If the MAC device is still unable to communicate, check that the communication settings for the MAC device are configured properly.	
	* Call tech support for assistance if problem persists		

Using the Serial/IP Tools to Troubleshoot

NOTE:

This section only applies to MUSE v5.x and MUSE v7.x systems as well CardioSoft V 6.6x and higher that use Tactical software. For MUSE v8.x and v9.x, troubleshooting, see "Testing Connections Using netstat" on page 43.

The Serial/IP Utility can be launched from the MUSE file server by right-clicking the icon in the systray and selecting either **Port Monitor...**, **Trace Window...**, or **Configure...**.

Port Activity

The Port Monitor window shows one COM port configured in Serial/IP. The lights show TR and DR when there is a connection. TD and RD flashes to show the transmission and reception of data.

FACTIC	AL ware	C-Bar	1-2	
ort Activity	Trace Window			
Port	TD RD TR DR CD	Status	IP Address	
COM11		Connected	192 168 0.50	

Trace Window

The **Trace Window** shows the activity of an ECG being sent into the MUSE system. Green text is outgoing data and red text is incoming data. Select the **Enable Trace** check box to display information in this window. By checking the **Auto Scroll** check box, the last entry is displayed automatically at the bottom of the list.

Tracel	utertra 1			
-ort Activity Trabe V	*******]			
08:11:55.726	COM11	:	39	^o
08:11:55.726	COM11	:	29	E^@^@*) ^8^@^@^?^Qf^ç^Î^I^I^E^Í^I^I^D**@
08:11:55.726	COM11	:	304	^D
08:11:55.945	COM11	:	-10	^DE^0^06^00^0^0^7^Q^ï\^î^I^I^D^î^I^I^E"
08:11:55.945	COM11	:	294	^D
08:11:55.945	COM11	:	>>	E^8^8***8^8*8*8*0*?*Qe*c*1*1*1*E*1*1*D****
08:11:55.945	COM11	:	39	^ <u>_</u>
08:11:56.945	COM11	:	-	^DE^8^868P^8^8^?^0^I[^Í^I^I^D^Í^I^I^E"
08:11:56.945	COM11	:	30.	^D
08:11:56.945	COM11	:	20	E^8^8*+^8^8*8*?*Qd^g^1^1^1^E^1^I^I^D****
08:11:56.945	COM11	:	39	<u>^</u>
08:11:57.195	COM11	:	~	*DE*8*8 (*80*8*8*?*0*1X*1*1*D*1*1**

The following **Trace Window** shows how the MUSE system sends data every 7 seconds to keep the CSI connection active. The blue text is informational and, in this case, shows the connection closing after powering down the cart.

e Edit Help	inton 3.0			=10
softwar	-	21	<	
Port Activity Trace \	Window			
08:03:54.226	COM11	:	- 20	^D
08:03:54.226	CON11	:	- 20-	E^8^8(^R^X^8^8^7^Q)^N^1^I^I^E^1^I^I^D"^
08:03:54.226	CON11	:	39	<u>^</u>
08:04:01.226	CON11	:	- 24	⁴ 0
08:04:01.226	COM11	:	29	E^8^8 (^S^X^8^8^7^Q)*N^1^I^I^I^I^I^I^I^N
08:04:01.226	COM11	:	35	^D
08:04:01.241	COM11	:	1	FlushTX
08:04:01.241	COM11	:	1	FlushRX
08:04:01.820	CON11	:	1	FlushTX
08:04:01.820	COM11	:	1	FlushRX
08:04:01.820	CON11		1	Close
08:04:09.820	com11		1	Port close
1				
Clear	Ena	ble	Trac	e F Hex Display V Auto Scroll F Always On Too

Testing Connections Using netstat

With the removal of the Tactical Software in MUSE v8.x and v9.x, you no longer have the ability to open this program and view the connection status of each network connection. To view connection status in MUSE v8 and v9, open a command prompt window on the system where the MUSE Modern service is running. This is typically the MUSE server.

To view the connection status, type the command *netstat -n*.

A listing is displayed showing all open socket connections on the system. If there is a connection established between the MUSE Modem service and the cart, you see an entry that has the "local address" with the IP address of the server, the "foreign address" with the IP address of the cart with the port number following the cart's IP address, and the state of the connection.

Local Address Foreign Address State

Local Address	Foreign Address	State
TCP 3.62.192.148:52379	TCP 3.62.192.159:3001	ESTABLISHED

In the above example, the server is at **3.62.192.148** and the cart is at **3.62.192.159**. The port number configured for this cart is **3001**. The connection state should be established.

If the cart is not connected, you do not see an entry in this list. The MUSE Modem service opens a connection once per the number of seconds equal to the retry interval. If it fails to connect, it closes this connection until the next interval time.

This command only shows established connections. If there is no connection, nothing is listed even though a CSI Network modern may be defined in the MUSE application.

Technical Specifications

Silex Wireless Bridge GEH-BR-4600

Manufacturer/Model	GEH-BR-4600WAN2-03-XX
Physical Requirements	Dimension: 110.5 × 79,0 × 27.6 (mm) Weight: 130 (g)
Interface Requirements	Ethernet Port: 10M/100 Mbps BASE-T, support Auto MDIX Power Connector: 5.5 (outer)/2.1 (Inner) mm Diameter Indication: LED indication for Power , wireless connection, Data communication
Power Requirements	Nominal Input Power: 5V Input Current: 750 mA
Wireless Requirements	Wireless LAN Protocol: IEEE 802.11a/b/g/n Wireless LAN Channel: IEEE 802.11b/g: Ch1~CH13 Encryption: WEP (64/128), WPA-PSK (TKIP/AES), WPA2-PSK (AES) IEEE802.1X enterprise authentication: EAP-PEAP, EAP-TLS, EAP-TTLS, EAP-FAST, EAP-LEAP
Frequency Band	2.4 GHz 5 GHz
Operating frequency range Maximum operating transmit power	Frequency Range: 2400-2483.5 MHz Output Power (including antenna gain): 100 mW (20 dBm)
Operating frequency range Maximum operating transmit power	Frequency Range: 5150-5350 MHz; 5470-5725 MHz Output Power (including antenna gain): 126 mW (21 dBm)

Environmental Requirements	Operating Temperature: 0 ~ 40 °C
	Operating Relative humidity: 20 ~ 80%
	Operating barometric pressures: 700 ~ 1060 hPa (Altitude range: 3010.9 to -381.9 meters)
	Non-Operating Temperature: -10 ~ +50 °C
	Non-Operating Relative humidity: 20% ~ 90%
	Non-Operating barometric pressures: 500 ~ 1060 hPa (Altitude range: 5570 to -380 meters)
Accessories	LAN Cable: Length 250 mm, RJ45 *2 connector
	USB Power Cable: Length 260 mm, Type A Plug, Right angle USB connector, Right angle DC connector, Wire size 24 AWG
Standards Compliance	CE certification
	FCC/IC certification
	 This product complies with the following Regulatory requirements for Australia, New Zealand and Singapore.
	• EMC Directive: EN55032 Class B, EN55024, EN301489-1/-17 v1.8.1
	 RE Directive: EN 300-328 v2.1.1, EN 301-893 v1.8.5 (EN 301-893 v1.8.1 (Adaptivity), EN 301-893 v2.1.0 (Receiver Blocking)), EN 60950-1, EN 62311, EN301-489-1 v2.1.1, EN301-489-17 v3.1.1
	 This product is compliant with the EU's RoHS directive (2011/65/EU or newer)
	• This product is compliant with the EU's WEEE directive (2002/96/EC)

GE Medical Systems Information Technologies, Inc. 8200 West Tower Avenue Milwaukee, WI 53223 USA Tel: +1 414 355 5000 +1 800 558 5120 (US Only)

GE Medical Systems *Information Technologies*, Inc., a General Electric Company, going to market as GE Healthcare.

www.gehealthcare.com