

# 1 Capnography (ETCO<sub>2</sub>)

## Operator's Manual Addendum

### 1.1 Overview

This addendum describes the operation of the Capnography option for the Puritan Bennett™ 980 Series Ventilator. The Capnography option monitors end-tidal carbon dioxide levels (ETCO<sub>2</sub>) only.

The capnography sensor uses infrared absorption spectroscopy to measure mainstream ETCO<sub>2</sub> while avoiding contamination with patient secretions.

### 1.2 Product Description

The Capnography option must be used with either of the following components:

- ETCO<sub>2</sub> airway adapter (for ETCO<sub>2</sub> monitoring only). This adapter is available in neonatal/pediatric and pediatric/adult configurations, and in reusable or disposable versions.
- Flow and CO<sub>2</sub> combination sensor (for use with neonatal patients, only). This component combines both a proximal flow transducer and an ETCO<sub>2</sub> built-in adapter. This sensor can be used when both NeoMode 2.0 and Capnography options are installed. This sensor is for single-use only.

The flow sensing component of the flow and CO<sub>2</sub> combination sensor is not required to be used during ETCO<sub>2</sub> monitoring. If the flow and CO<sub>2</sub> combination sensor is used for ETCO<sub>2</sub> monitoring only, the sensor must be connected to the ventilator's front panel in the location labeled Prox, and the Proximal Flow option disabled. This combination sensor can be used when the Proximal Flow System option is also in use so that flow sensing and ETCO<sub>2</sub> monitoring can both occur simultaneously. More information regarding the Proximal Flow System is available in the *Proximal Flow Option Appendix* and NeoMode 2.0 details are described in the *NeoMode 2.0 Appendix* in this manual.

For general parameter and ventilator setup information Reference *Chapter 4* in this manual.

### 1.3 Intended Use

The Capnography option is used for measuring the partial pressure of carbon dioxide in the exhaled gas of invasively ventilated, neonatal, pediatric, and adult patients at the breathing circuit wye.

## 1.4 Safety Symbol Definitions

This section contains safety information for users who should always exercise appropriate caution while using the ventilator.

**Table 1-1.** Safety Symbol Definitions

Symbol	Definition
	<b>WARNING</b> Warnings alert users to potential serious outcomes (death, injury, or adverse events) to the patient, user, or environment.
	<b>Caution</b> Cautions alert users to exercise appropriate care for safe and effective use of the product.
	<b>Note</b> Notes provide additional guidelines or information.

## 1.5 Safety Information



### **WARNING:**

The Puritan Bennett™ 980 series ventilator contains phthalates. When used as indicated, very limited exposure to trace amounts of phthalates may occur. There is no clear clinical evidence that this degree of exposure increases clinical risk. However, in order to minimize risk of phthalate exposure in children and nursing or pregnant women, this product should only be used as directed.



### **WARNING:**

The ventilator offers a variety of breath delivery options. Throughout the patient's treatment, the clinician should carefully select the ventilation mode and settings to use for that patient, based on clinical judgment, considering the condition and needs of the individual patient, as they change from time to time, and considering the benefits, limitations and operating characteristics of each breath delivery option.



### **WARNING:**

Explosion hazard. Do not use in the presence of flammable anesthetics.



### **WARNING:**

Follow precautions for electromagnetic interference (EMI) to avoid unreliable ventilator readings.

**WARNING:**

Monitor the CO<sub>2</sub> waveform for elevated baseline. An erroneously elevated baseline can be caused by sensor problems.

**WARNING:**

If the CO<sub>2</sub> waveform appears abnormal, inspect the airway adapters and replace if needed.

**WARNING:**

ETCO<sub>2</sub> readings are intended only as an adjunct in patient assessment and must be used in conjunction with clinical signs and symptoms. Do not use ETCO<sub>2</sub> readings as a basis for changes to ventilator parameters without reference to clinical condition and independent monitors such as blood gas.

**WARNING:**

Do not use the flow and CO<sub>2</sub> combination sensor if there are kinks in the tubing.

**WARNING:**

Prior to patient ventilation, run SST with the exact configuration that will be used on the patient. This includes a patient circuit, airway adapter, and all accessories used with the patient circuit. Reference *To run SST* in Chapter 3 of this manual.

**WARNING:**

Changing ventilator accessories can change the system resistance and compliance. Do not add or remove accessories after running SST.

**WARNING:**

Discontinue use if ETCO<sub>2</sub> monitoring fails to respond as described.

**WARNING:**

The flow and CO<sub>2</sub> combination sensor measures gas flow at the patient wye. A system leak, such as that caused by an uncuffed endotracheal tube or a damaged flow and CO<sub>2</sub> combination sensor may significantly affect flow-related readings.

**WARNING:**

Use only Covidien-branded flow and CO<sub>2</sub> combination sensors, capnography sensors, and airway adapters with the Capnography option. Use of other sensors results in an "invalid sensor" message or a "capnography sensor inoperative" alarm and/or erroneous readings.

**WARNING:**

To minimize the potential for condensation or secretions clogging the sensor's pneumatic lines, position the flow and CO<sub>2</sub> combination sensor exactly as described in this addendum.



**WARNING:**

**Do not position the capnography sensor cable or flow and CO<sub>2</sub> combination sensor tubing in any manner that may cause entanglement or strangulation.**



**WARNING:**

**To reduce the risk of extubation or disconnection, do not apply tension to or pull on the capnography sensor or flow and CO<sub>2</sub> combination sensor tubing.**



**WARNING:**

**To reduce the risk of extubation or breathing circuit disconnection, do not rotate the flow and CO<sub>2</sub> combination sensor in the breathing circuit by pulling on the sensor's tubing.**



**WARNING:**

**The cable management clips supplied with each flow and CO<sub>2</sub> combination sensor must be used to mitigate risk of entanglement, kinking, or extubation that could lead to strangulation, hypercarbia, or hypoxemia.**



**WARNING:**

**Do not install the capnography sensor or flow and CO<sub>2</sub> combination sensor in the patient circuit if the sensor is not also connected to the BDU.**



**WARNING:**

**Excessive moisture in the flow and CO<sub>2</sub> combination sensor tubing may affect the accuracy of the measurements. Periodically check the sensor and tubing for excessive moisture or secretion build-up.**



**WARNING:**

**To avoid the risk of bacterial contamination, clean and sterilize reusable airway adapters per the methods described in this addendum.**



**WARNING:**

**Disposable airway adapters are intended for single use only. Do not re-use these items. These sensors and adapters are not compatible with sterilization techniques.**



**WARNING:**

**Nitrous oxide, elevated levels of oxygen, helium, and halogenated hydrocarbons can influence the CO<sub>2</sub> measurement.**



**WARNING:**

**Inspect the airway adapter or flow and CO<sub>2</sub> combination sensor prior to use, and do not use if the adapter, sensor body, tubing, or connector are damaged or broken.**

-  **Caution:**  
**Do not use aerosolized medications when employing CO<sub>2</sub> monitoring with either the airway adapter and capnography sensor or the flow and CO<sub>2</sub> combination sensor. Increased medication viscosity may contaminate the sensor windows and cause the sensor to fail prematurely.**
-  **Caution:**  
**To prevent damage to cables or pneumatic lines, use the included cable management clips.**
-  **Caution:**  
**Insert sensors in the ventilator circuit with the tubes upright to avoid the effects of excessive moisture.**
-  **Caution:**  
**Ensure all connectors are properly connected, fully engaged, and free from moisture.**
-  **Caution:**  
**To avoid possible damage to the ventilator or sensors, follow standard precautions for electrostatic discharge (ESD).**
-  **Caution:**  
**Clean the reusable airway adapters using only the method and cleaning agents described in this addendum.**
-  **Note:**  
The flow and CO<sub>2</sub> combination sensor, capnography sensor, and airway adapters contain no user serviceable parts. Refer service to qualified service personnel.
-  **Note:**  
Dispose of the flow and CO<sub>2</sub> combination sensor and disposable airway adapters in accordance with your institution's protocol.
-  **Note:**  
The white-striped tubing of the flow sensor should always be proximal to the patient.
-  **Note:**  
Position the flow and CO<sub>2</sub> combination sensor with its windows in a vertical, not horizontal position. This helps keep patient secretions from pooling on the windows.

## 1.6 Software Requirements

Purchased software options must be enabled after the option is purchased using an encrypted access code provided to you or the Customer Service Engineer (CSE).

### To install software options

1. Enter Service mode. Reference *To access Service Mode*, in *Chapter 3* of this manual.
2. Touch the options button on the screen that appears.
3. Touch Installed Options.
4. Touch *Update Options*.
5. Enter the option access code on the virtual keyboard and touch *Accept*.
6. Confirm the option is installed by touching *Installed Options*.
7. Attach software option label to the installed software options label located at the back of the ventilator.

## 1.7 Hardware Requirements

The following hardware is required:

- Option Host Card
- Capnography sensor (required for ETCO<sub>2</sub> measurement)
- Airway adapter suitable for the particular patient OR
- Flow and CO<sub>2</sub> combination sensor (neonatal circuit only and supports both proximal flow measurement and/or ETCO<sub>2</sub> measurement)

## 1.8 End Tidal CO<sub>2</sub> Monitoring Description



### Note:

If the current vent type is NIV or vent type changes from INVASIVE to NIV, the Capnography option is automatically disabled.

A capnography sensor is attached to an airway adapter or flow and CO<sub>2</sub> combination sensor which is installed at the patient circuit wye to measure end-tidal CO<sub>2</sub>. The capnography sensor connects to the BDU via a connector located behind the door on the BDU's front panel.

[\*Capnography Sensor Installation to Flow and CO<sub>2</sub> Combination Sensor and Ventilator\*](#), page 1-14.  
Airway adapters are available in pediatric/adult and neonatal sizes.

When the Option Host card is installed in the ventilator, the card is the communication interface for the capnography system. Data measured by the capnography sensor are displayed on the GUI for monitoring purposes, not for ventilator control. When the ventilator has a capnography sensor installed and enabled, end-tidal CO<sub>2</sub> measurements are obtained and displayed on the GUI, and data are updated at the end of each exhalation. CO<sub>2</sub> data can be configured as a patient data parameter. If ETCO<sub>2</sub> is chosen as a patient data parameter, dashes (- -) appear if the ETCO<sub>2</sub> data value is unable to be displayed. The CO<sub>2</sub> waveform can be configured as a waveform layout, if desired. Reference *To configure the patient data displayed on the GUI* and *To configure waveforms and loops* in Chapter 3 of this manual.

The time duration from enabling the ETCO<sub>2</sub> monitoring function to achieving specified accuracy is approximately two minutes due to the capnography sensor warm-up process.

### 1.8.1 ETCO<sub>2</sub> Monitoring Components

ETCO<sub>2</sub> monitoring components include:

**Capnography sensor** — The capnography sensor connects to the airway adapter or the flow and CO<sub>2</sub> combination sensor and contains the optics and electronics for ETCO<sub>2</sub> measurement.

**ETCO<sub>2</sub> Airway Adapter** — (for ETCO<sub>2</sub> monitoring only. This adapter is available in neonatal/pediatric and pediatric/adult configurations, and in reusable or disposable versions.

**Flow and CO<sub>2</sub> Combination Sensor** — This sensor and adapter combination is a single-use device which acts as a combination of a Proximal Flow Sensor and an airway adapter (for use with neonatal patients, only). This component combines both a proximal flow transducer and an ETCO<sub>2</sub> built-in adapter. This sensor can be used for capnography when both NeoMode 2.0 and Capnography options are installed. This sensor is for single-use only.



#### Note:

Either the airway adapter or Flow and CO<sub>2</sub> combination sensor are connected to the capnography sensor for ETCO<sub>2</sub> measurement only, or proximal flow and/or ETCO<sub>2</sub> measurement, respectively.

Use the airway adapters, flow and CO<sub>2</sub> combination sensor, and capnography sensor in the correct combination based upon the patient type and desired monitoring shown in the following tables.

**Table 1-2.** Endotracheal Tube Diameters


Item	Endotracheal Tube Inner Diameter (ID), mm
Neonatal/pediatric airway adapters	≤ 4.0
Pediatric/adult airway adapters	≥ 4.5
Neonatal flow and CO <sub>2</sub> combination sensor	2.5 to 4.0

**Table 1-3.** Items Needed for ETCO<sub>2</sub> or Flow/Volume Monitoring


Desired Monitoring	Items Needed	Patient Type
ETCO <sub>2</sub>	Capnography sensor + neonatal/pediatric airway adapter	Neonatal, pediatric
	Capnography sensor + pediatric/adult airway adapter	Pediatric, adult
Flow/Volume Only	Flow and CO <sub>2</sub> combination sensor	Neonatal
	Proximal flow sensor <sup>1</sup>	Neonatal
ETCO <sub>2</sub> and Flow/Volume	Flow and CO <sub>2</sub> combination sensor and capnography sensor	Neonatal

1. Reference the *Proximal Flow Option Appendix* in this manual.


Reference *Capnography Monitoring System Part Numbers*, page 1-34 to distinguish the components by color.

**Figure 1-1.** Capnography Sensor

VEN\_10323\_B

**Figure 1-2.** ETCO<sub>2</sub> Airway Adapters

VEN\_10325\_B

1 Pediatric/Adult ETCO<sub>2</sub> Airway Adapter      2 Neonatal/Pediatric ETCO<sub>2</sub> Airway Adapter**Figure 1-3.** Flow and CO<sub>2</sub> Combination Sensor

VEN\_10298\_B

## 1.9 On-screen Symbols

When the ETCO<sub>2</sub> monitoring function is enabled, measured ETCO<sub>2</sub> can be configured to display in the patient data banner. Reference *Vital Patient Data* in Chapter 3 of this manual for information on configuring the GUI to display various patient data values. When ETCO<sub>2</sub> data are questionable or invalid, the data are not displayed.

**Figure 1-4.** Sample GUI Screen Showing ETCO<sub>2</sub> Data

VEN\_11537\_B

1 Patient data banner configured with ETCO<sub>2</sub>

## 1.10 SST Requirements

SST is not required for ETCO<sub>2</sub> monitoring only, but in cases where the flow and CO<sub>2</sub> combination sensor is used for proximal flow measurement as well as ETCO<sub>2</sub> measurement, SST is required. SST must be performed with all circuit components in the configuration to be used on the patient in order for the ventilator to calculate the correct compliance and resistance. Reference *To run SST in Chapter 3* of this manual for complete instructions on how to run SST.



### Note:

For SST, the capnography sensor does not need to be connected to the ventilator, but the Flow and CO<sub>2</sub> sensor must connect to the ventilator's port labeled Prox Flow to avoid leaks.

### 1.10.1 Attaching the Flow and CO<sub>2</sub> Combination Sensor for SST

#### To attach the flow and CO<sub>2</sub> combination sensor to the patient circuit

1. Verify the flow and CO<sub>2</sub> combination sensor, pneumatic lines and connector are not damaged.
2. Open the connector panel door and firmly attach the flow and CO<sub>2</sub> combination sensor connector to the receptacle in the BDU's front connector port labeled Prox Flow. Reference *Capnography Sensor Installation to Flow and CO<sub>2</sub> Combination Sensor and Ventilator*, page 1-14.

3. Insert the **smaller**(15 mm diameter) end (opposite end of UP arrow) of the airway adapter portion of the flow and CO<sub>2</sub> combination sensor into the breathing circuit wye, as shown. Reference *Attaching Flow and CO<sub>2</sub> Combination Sensor or Airway Adapter to Breathing Circuit Wye*, page 1-11. Ensure the sensor tubing is in the upward position.
4. To reposition, grasp the sensor plastic body. Do not rotate it by pulling on the tubing.
5. Confirm a tight connection.
6. Space the clips provided evenly to secure the sensor tubing to the breathing circuit.
7. Zero the adapter if it is new or if prompted by a message displayed on the ventilator's GUI. Reference *Zeroing the Sensor*, page 1-15.
8. Run SST. Reference *To run SST* in *Chapter 3* of this manual.



**Note:**

**If using a Heat-Moisture Exchanger (HME) on the endotracheal tube, place the sensor between the HME and the breathing circuit wye.**

**Figure 1-5.** Attaching Flow and CO<sub>2</sub> Combination Sensor or Airway Adapter to Breathing Circuit Wye



VEN\_11540\_A

### 1.10.2 Attaching the airway adapter for SST

#### To attach the airway adapter for SST

1. Attach the capnography sensor and airway adapter so the **small** end of the airway adapter connects to the breathing circuit wye. Reference the figure above. The airway adapter attaches the same way as the Flow and CO<sub>2</sub> Combination sensor.

2. Run SST. Reference *To run SST* in Chapter 3 of this manual.

## 1.11 Using the ETCO<sub>2</sub> Monitoring function

Ensure SST has been run with all accessories installed in the ventilator breathing circuit.

Review and follow all warnings prior to patient ventilation with the ETCO<sub>2</sub> monitoring function. Reference *Safety Information*, page 1-2.


### To connect the capnography sensor to the ventilator

1. Verify the capnography sensor and connector are not damaged in any way.
2. Open the connector panel door and firmly attach the sensor connector to the receptacle in the BDU's front connector port labeled CO<sub>2</sub>.

### To attach the ETCO<sub>2</sub> airway adapter to the capnography sensor

1. Use an appropriate ETCO<sub>2</sub> airway adapter that matches the patient type (neonatal/pediatric or pediatric/adult) and the patient being ventilated.
2. Press the capnography sensor onto the ETCO<sub>2</sub> airway adapter, oriented so that the **smaller** end of the adapter can attach to the breathing circuit wye. The capnography sensor will "click" into place when properly seated.

**Figure 1-6.** Connecting the Capnography Sensor to the ETCO<sub>2</sub> Airway Adapter




VEN\_10324\_B

**To attach the capnography sensor and airway adapter to the patient circuit**

1. Attach the capnography sensor and airway adapter so the **small** end of the airway adapter connects to the breathing circuit wye, and the large end connects to the ET tube. Reference *Connecting Breathing Circuit and the ET Tube*, page 1-13 for an illustration of the breathing circuit and ET tube connection. The capnography sensor and airway adapter connect the same way as the flow and CO<sub>2</sub> combination sensor.


**Figure 1-7.** Connecting Breathing Circuit and the ET Tube



VEN\_10326\_A

**To attach the flow and CO<sub>2</sub> combination sensor to the capnography sensor**

1. Snap the flow and CO<sub>2</sub> combination sensor into the capnography sensor as shown. Reference *Capnography Sensor Installation to Flow and CO<sub>2</sub> Combination Sensor and Ventilator*, page 1-14.

**Figure 1-8.** Capnography Sensor Installation to Flow and CO<sub>2</sub> Combination Sensor and Ventilator

- |   |   |   |                       |
|---|---|---|-----------------------|
| 1 | Capnography sensor (connect to BDU port labeled CO <sub>2</sub> )                   | 5 | Breathing circuit wye |
| 2 | Patient circuit inspiratory and expiratory limbs                                    | 6 | HME                   |
| 3 | Flow and CO <sub>2</sub> combination sensor (connect to BDU port labeled Prox Flow) | 7 | Elbow (optional)      |
| 4 | Cable management clip   | 8 | Endotracheal tube     |

**Note:**

If using a Heat-Moisture Exchanger (HME) on the endotracheal tube, place the flow and CO<sub>2</sub> combination sensor between the HME and the breathing circuit wye.

2. Ensure the sensor tubing is positioned in an upward direction, as shown in the figure above. If the sensor needs repositioning, **DO NOT** rotate it by pulling on the tubing. Reposition as follows:

- a. Grasp the sensor's plastic body with one hand and the breathing circuit wye with the other hand.
  - b. Rotate the sensor body and wye towards each other until the sensor tubing is upright.
  - c. Confirm a tight connection between the sensor and breathing circuit wye.
3. Use cable management clips to attach the sensor tubing to the breathing circuit tubing. Space the clips evenly along the length of the sensor tubing. Twist the ends of each clip to close. Ensure the appropriate clip size is used with the patient circuit in use.
  4. Connect the capnography sensor to the ventilator as described above and shown.
  5. Connect the Flow and CO<sub>2</sub> combination sensor to the BDU's front connector port labeled Prox Flow, as shown. Reference *Capnography Sensor Installation to Flow and CO<sub>2</sub> Combination Sensor and Ventilator*, page 1-14.

**Note:**

When the ventilator is set up for Proximal Flow System operation, the flow and CO<sub>2</sub> combination sensor can be switched as necessary. There is no need to run SST after switching sensors unless the breathing circuit or other ventilator accessories have been changed.

### 1.11.1 Zeroing the Sensor


**WARNING:**

If the sensor/adapter assembly zero process is not executed properly, the CO<sub>2</sub> waveform information may not be representative of actual exhaled CO<sub>2</sub>.

If a problem is detected with the CO<sub>2</sub> measurement the GUI displays a message stating "CO<sub>2</sub> sensor zeroing required" or "Check CO<sub>2</sub> sensor adapter." Zeroing the sensor is typically done when switching between different types of airway adapters. This process corrects for optical differences due to reusable vs. disposable adapters.

**Before zeroing the sensor**

1. Verify the capnography sensor is connected to the port labeled CO<sub>2</sub> on the ventilator.



VEN\_11166\_E

2. Connect the capnography sensor to a clean and dry airway adapter or flow and CO<sub>2</sub> combination sensor. Reference [Connecting the Capnography Sensor to the ETCO<sub>2</sub> Airway Adapter](#), page 1-12 or [Capnography Sensor Installation to Flow and CO<sub>2</sub> Combination Sensor and Ventilator](#), page 1-14. Ensure the capnography sensor and airway adapter or flow and CO<sub>2</sub> combination sensor is motionless, exposed only to room air, and away from all sources of CO<sub>2</sub>, including the ventilator, the patient's breath, and your breath.

### To zero the sensor

1. Use the capnography sensor and airway adapter or flow and CO<sub>2</sub> combination sensor connected to the ventilator from the previous steps.
2. Wait two (2) minutes before proceeding.
3. At the ventilator setup screen, touch the configure icon. A menu containing tabs appears.
4. Touch the Options tab. A screen appears containing Installed Options and CO<sub>2</sub> tabs.
5. Touch the CO<sub>2</sub> tab.
6. Follow zeroing setup instructions on the GUI screen.
7. Touch the Start button to zero the sensor. The length of time to zero the sensor is typically 15 to 20 seconds. During this time, ensure the sensor and adapter is not exposed to any form of CO<sub>2</sub>.
8. Verify the message on the screen indicates "CO<sub>2</sub> zeroing passed."
9. If message reads "Zeroing failed," ensure all zeroing conditions are met and retest or change adapter.
10. Place the successfully zeroed sensor assembly in the patient circuit. Reference [Attaching Flow and CO<sub>2</sub> Combination Sensor or Airway Adapter to Breathing Circuit Wye](#), page 1-11.

### 1.11.2 Disabling/enabling the Capnography Option

The Capnography option has two states: Enabled and Disabled.



#### Note:

If the Capnography option has been disabled or enabled, SST does not have to be re-run unless the breathing circuit or other breathing system accessories have been changed (including the flow and CO<sub>2</sub> combination sensor), removed, or added.



#### To disable or enable the capnography monitoring function

1. At the ventilator setup screen, touch the configure icon. A menu containing tabs appears.
2. Touch the Options tab. A screen appears containing Installed Options and CO<sub>2</sub> tabs.
3. Touch the CO<sub>2</sub> tab.
4. Touch *Enabled* or *Disabled* to enable or disable the Capnography option.

### CO<sub>2</sub> Waveform

The CO<sub>2</sub> waveform appears when the Capnography option is enabled and the waveform's y-axis is configured to show CO<sub>2</sub>. The CO<sub>2</sub> waveform displays the expired CO<sub>2</sub> value as measured and reported by the capnography sensor at the circuit wye.

**Figure 1-9.** CO<sub>2</sub> Waveform



VEN\_11539\_B

## 1.12 Accuracy Check

No user-performed calibration is necessary or possible, except for zeroing the sensor. An accuracy check of the capnography sensor should be performed once per year.

**WARNING:**

**Before performing an accuracy check, ensure no patient is connected to the ventilator.**

Before performing the accuracy check, ensure the Capnography option has been enabled. Reference *Disabling/enabling the Capnography Option*, page 1-17.


Perform an accuracy check using the following items:

**Table 1-4.** Parts Required for Accuracy Check

Manufacturer	Description	Part number
Philips (formerly known as Respi-ronics)	Gas regulator	6081-00
Air Liquide	Calibration gas [5% CO <sub>2</sub> with the balance being Nitrogen] (carton of 4)	T4507NM-4PD (USA)
Air Liquide	Calibration gas [5% CO <sub>2</sub> with the balance being Nitrogen] (sold by the can with a required minimum order of 6)	T4507SRI-PD (Other countries)
External barometer	For use in calculating the expected CO <sub>2</sub> value obtained during the accuracy check.	Local supplier

**To perform an accuracy check**

1. Allow the calibration gas time to reach room temperature (approximately 30 minutes depending upon the ambient storage temperature of the calibration gas).
2. Attach the capnography sensor to the ventilator at the connector port labeled CO<sub>2</sub> on the ventilator's front panel. Ensure the capnography sensor/adapters are not connected to the patient circuit.
3. Connect three identical airway adapters end-to-end to the gas regulator. This configuration is called a stack.

**Figure 1-10.** Airway Adapter Configuration

VEN\_11242\_B

1 Top of stack (smaller diameter)


3 Bottom of stack (larger diameter)

2 Attach capnography sensor to middle  
airway adapter

4. Turn on the ventilator and configure it for an Adult circuit type, if not already configured. To change the circuit type, SST must be performed. Reference *To run SST* in Chapter 3 of this manual.
5. Set the ventilator up for a New Patient. The Quick Start settings may be used.
6. Set the oxygen concentration to 21%.
7. Connect a patient circuit to the ventilator and to a test lung. At the time the circuit and test lung are connected, the ventilator will sense a patient connection, and the waveforms display will appear on the GUI screen.
8. Change the waveform display to show CO<sub>2</sub> in mmHg, by double-tapping the y-axis parameter. A list of buttons appears with parameters.

9. Touch the CO<sub>2</sub> button. The waveform now shows CO<sub>2</sub> in mmHg.
10. Connect the capnography sensor to the middle airway adapter.
11. Zero the sensor. Reference *Zeroing the Sensor*, page 1-15.
12. Turn the gas regulator flow control valve OFF (turned fully clockwise).
13. Attach the gas regulator and hose assembly to the calibration gas cylinder.

**Figure 1-11.** Gas Regulator Attachment



VEN\_11241\_A

1	Regulator output. Attach bottom of airway adapter stack (larger diameter) here	3	Flow control valve
2	Factory preset - do not adjust	4	Pressure gauge

14. Connect the gas regulator output to the bottom of the airway adapter stack described in step 3.
15. Provide a constant flow of the test gas (approximately 2 L/min) through the airway adapters by fully opening the valve, and allow the gas to flow for 30 s.
16. Touch the pause icon for the waveform and use the knob to move the cursor along the waveform. Read the CO<sub>2</sub> measurement appearing on the cursor.
17. Turn the test gas off.

18. Since the gas stabilized at room temperature, use the ambient temperature, pressure displayed when the option was enabled, and CO<sub>2</sub> measurements from the display, and compare with the following table:



**Note:**

To use the table, locate the intersection of the pressure measurement in the Barometric Pressure Reading column, and the ambient temperature in the Gas Temperature row. The intersection of pressure and temperature is the expected CO<sub>2</sub> measurement in mmHg.

**Table 1-5.** Expected CO<sub>2</sub> measurements (mm Hg) vs. Temperature and Pressure

Barometric Pressure Reading (mmHg)	Gas Temperature (°C)						
	20	21	22	23	24	25	26
620	33.1	32.9	32.8	32.7	32.6	32.5	32.4
621	33.1	33.0	32.9	32.8	32.7	32.6	32.5
622	33.2	33.1	32.9	32.8	32.7	32.6	32.5
623	33.2	33.1	33.0	32.9	32.8	32.7	32.6
624	33.3	33.2	33.1	33.0	32.9	32.7	32.6
625	33.3	33.2	33.1	33.0	32.9	32.8	32.7
626	33.4	33.3	33.2	33.1	33.0	32.9	32.8
627	33.4	33.3	33.2	33.1	33.0	32.9	32.8
628	33.5	33.4	33.3	33.2	33.1	33.0	32.9
629	33.5	33.4	33.3	33.2	33.1	33.0	32.9
630	33.6	33.5	33.4	33.3	33.2	33.1	33.0
631	33.6	33.5	33.4	33.3	33.2	33.1	33.0
632	33.7	33.6	33.5	33.4	33.3	33.2	33.1
633	33.7	33.6	33.5	33.4	33.3	33.2	33.1
634	33.8	33.7	33.6	33.5	33.4	33.3	33.2
635	33.8	33.7	33.6	33.5	33.4	33.3	33.2
636	33.9	33.8	33.7	33.6	33.5	33.4	33.3
637	34.0	33.9	33.7	33.6	33.5	33.4	33.3
638	34.0	33.9	33.8	33.7	33.6	33.5	33.4
639	34.1	34.0	33.9	33.7	33.6	33.5	33.4
640	34.1	34.0	33.9	33.8	33.7	33.6	33.5

**Table 1-5.** Expected CO<sub>2</sub> measurements (mm Hg) vs. Temperature and Pressure (Continued)

Barometric Pressure Reading (mmHg)	Gas Temperature (°C)						
	20	21	22	23	24	25	26
641	34.2	34.1	34.0	33.9	33.7	33.6	33.5
642	34.2	34.1	34.0	33.9	33.8	33.7	33.6
643	34.3	34.2	34.1	34.0	33.9	33.7	33.6
644	34.3	34.2	34.1	34.0	33.9	33.8	33.7
645	34.4	34.3	34.2	34.1	34.0	33.9	33.7
646	34.4	34.3	34.2	34.1	34.0	33.9	33.8
647	34.5	34.4	34.3	34.2	34.1	34.0	33.9
648	34.5	34.4	34.3	34.2	34.1	34.0	33.9
649	34.6	34.5	34.4	34.3	34.2	34.1	34.0
650	34.6	34.5	34.4	34.3	34.2	34.1	34.0
651	34.7	34.6	34.5	34.4	34.3	34.2	34.1
652	34.8	34.6	34.5	34.4	34.3	34.2	34.1
653	34.8	34.7	34.6	34.5	34.4	34.3	34.2
654	34.9	34.8	34.6	34.5	34.4	34.3	34.2
655	34.9	34.8	34.7	34.6	34.5	34.4	34.3
656	35.0	34.9	34.8	34.6	34.5	34.4	34.3
657	35.0	34.9	34.8	34.7	34.6	34.5	34.4
658	35.1	35.0	34.9	34.7	34.6	34.5	34.4
659	35.1	35.0	34.9	34.8	34.7	34.6	34.5
660	35.2	35.1	35.0	34.9	34.7	34.6	34.5
661	35.2	35.1	35.0	34.9	34.8	34.7	34.6
662	35.3	35.2	35.2	35.0	34.9	34.7	34.6
663	35.3	35.2	35.1	35.0	34.9	34.8	34.7
664	35.4	35.3	35.2	35.1	35.0	34.8	34.7
665	35.4	35.3	35.2	35.1	35.0	34.9	34.8
666	35.5	35.4	35.3	35.2	35.1	35.0	34.8
667	35.6	35.4	35.3	35.2	35.1	35.0	34.9
668	35.6	35.5	35.4	35.3	35.2	35.1	35.0
669	35.7	35.6	35.4	35.3	35.2	35.1	35.0

**Table 1-5.** Expected CO<sub>2</sub> measurements (mm Hg) vs. Temperature and Pressure (Continued)

Barometric Pressure Reading (mmHg)	Gas Temperature (°C)						
	20	21	22	23	24	25	26
670	35.7	35.6	35.5	35.4	35.3	35.2	35.1
671	35.8	35.7	35.5	35.4	35.3	35.2	35.1
672	35.8	35.7	35.6	35.5	35.4	35.3	35.2
673	35.9	35.8	35.7	35.5	35.4	35.3	35.2
674	35.9	35.8	35.7	35.6	35.5	35.4	35.3
675	36.0	35.9	35.8	35.6	35.5	35.4	35.3
676	36.0	35.9	35.8	35.7	35.6	35.5	35.4
677	36.1	36.0	35.9	35.8	35.6	35.5	35.4
678	36.1	36.0	35.9	35.8	35.7	35.6	35.5
679	36.2	36.1	36.0	35.9	35.7	35.6	35.5
680	36.2	36.1	36.0	35.9	35.8	35.7	35.6
681	36.3	36.2	36.1	36.0	35.9	35.7	35.6
682	36.4	36.2	36.1	36.0	35.9	35.8	35.7
683	36.4	36.3	36.2	36.1	36.0	35.8	35.7
684	36.5	36.3	36.2	36.1	36.0	35.9	35.8
685	36.5	36.4	36.3	36.2	36.1	36.0	35.8
686	36.6	36.5	36.3	36.2	36.1	36.0	35.9
687	36.6	36.5	36.4	36.3	36.2	36.1	35.9
688	36.7	36.6	36.4	36.3	36.2	36.1	36.0
689	36.7	36.6	36.5	36.4	36.3	36.2	36.0
690	36.8	36.7	36.6	36.4	36.3	36.2	36.1
691	36.8	36.7	36.6	36.5	36.4	36.3	36.2
692	36.9	36.8	36.7	36.5	36.4	36.3	36.2
693	36.9	36.8	36.7	36.6	36.5	36.4	36.3
694	37.0	36.9	36.8	36.6	36.5	36.4	36.3
695	37.0	36.9	36.8	36.7	36.6	36.5	36.4
696	37.1	37.0	36.9	36.8	36.6	36.5	36.4
697	37.2	37.0	36.9	36.8	36.7	36.6	36.5
698	37.2	37.1	37.0	36.9	36.7	36.6	36.5

**Table 1-5.** Expected CO<sub>2</sub> measurements (mm Hg) vs. Temperature and Pressure (Continued)

Barometric Pressure Reading (mmHg)	Gas Temperature (°C)						
	20	21	22	23	24	25	26
699	37.3	37.1	37.0	36.9	36.8	36.7	36.6
700	37.3	37.2	37.1	37.0	36.9	36.7	36.6
701	37.4	37.3	37.1	37.0	36.9	36.8	36.7
702	37.4	37.3	37.2	37.1	37.0	36.8	36.7
703	37.5	37.4	37.2	37.1	37.0	36.9	36.8
704	37.5	37.4	37.3	37.2	37.1	36.9	36.8
705	37.6	37.5	37.3	37.2	37.1	37.0	36.9
706	37.6	37.5	37.4	37.3	37.2	37.1	36.9
707	37.7	37.6	37.5	37.3	37.2	37.1	37.0
708	37.7	37.6	37.5	37.4	37.3	37.2	37.0
709	37.8	37.7	37.6	37.4	37.3	37.2	37.1
710	37.8	37.7	37.6	37.5	37.4	37.3	37.1
711	37.9	37.8	37.7	37.5	37.4	37.3	37.2
712	38.0	37.8	37.7	37.6	37.5	37.4	37.3
713	38.0	37.9	37.8	37.7	37.5	37.4	37.3
714	38.1	37.9	37.8	37.7	37.6	37.5	37.4
715	38.1	38.0	37.9	37.8	37.6	37.5	37.4
716	38.2	38.0	37.9	37.8	37.7	37.6	37.5
717	38.2	38.1	38.0	37.9	37.7	37.6	37.5
718	38.3	38.2	38.0	37.9	37.8	37.7	37.6
719	38.3	38.2	38.1	38.0	37.9	37.7	37.6
720	38.4	38.3	38.1	38.0	37.9	37.8	37.7
721	38.4	38.3	38.2	38.1	38.0	37.8	37.7
722	38.5	38.4	38.2	38.1	38.0	37.9	37.8
723	38.5	38.4	38.3	38.2	38.1	37.9	37.8
724	38.6	38.5	38.4	38.2	38.1	38.0	37.9
725	38.6	38.5	38.4	38.3	38.2	38.0	37.9
726	38.7	38.6	38.5	38.3	38.2	38.1	38.0
727	38.8	38.6	38.5	38.4	38.3	38.2	38.0

**Table 1-5.** Expected CO<sub>2</sub> measurements (mm Hg) vs. Temperature and Pressure (Continued)

Barometric Pressure Reading (mmHg)	Gas Temperature (°C)						
	20	21	22	23	24	25	26
728	38.8	38.7	38.6	38.4	38.3	38.2	38.1
729	38.9	38.7	38.6	38.5	38.4	38.3	38.1
730	38.9	38.8	38.7	38.6	38.4	38.3	38.2
731	39.0	38.8	38.7	38.6	38.5	38.4	38.2
732	39.0	38.9	38.8	38.7	38.5	38.4	38.3
733	39.1	39.0	38.8	38.7	38.6	38.5	38.4
734	39.1	39.0	38.9	38.8	38.6	38.5	38.4
735	39.2	39.1	38.9	38.8	38.7	38.6	38.5
736	39.2	39.1	39.0	38.9	38.7	38.6	38.5
737	39.3	39.2	39.0	38.9	38.8	38.7	38.6
738	39.3	39.2	39.1	39.0	38.9	38.7	38.6
739	39.4	39.3	39.1	39.0	38.9	38.8	38.7
740	39.4	39.3	39.2	39.1	39.0	38.8	38.7
741	39.5	39.4	39.3	39.1	39.0	38.9	38.8
742	39.6	39.4	39.3	39.2	39.1	38.9	38.8
743	39.6	39.5	39.4	39.2	39.1	39.0	38.9
744	39.7	39.5	39.4	39.3	39.2	39.0	38.9
745	39.7	39.6	39.5	39.3	39.2	39.1	39.0
746	39.8	39.6	39.5	39.4	39.3	39.2	39.0
747	39.8	39.7	39.6	39.4	39.3	39.2	39.1
748	39.9	39.7	39.6	39.5	39.4	39.3	39.1
749	39.9	39.8	39.7	39.6	39.4	39.3	39.2
750	40.0	39.9	39.7	39.6	39.5	39.4	39.2
751	40.0	39.9	39.8	39.7	39.5	39.4	39.3
752	40.1	40.0	39.8	39.7	39.6	39.5	39.3
753	40.1	40.0	39.9	39.8	39.6	39.5	39.4
754	40.2	40.1	39.9	39.8	39.7	39.6	39.4
755	40.2	40.1	40.0	39.9	39.7	39.6	39.5
756	40.3	40.2	40.0	39.9	39.8	39.7	39.6

**Table 1-5.** Expected CO<sub>2</sub> measurements (mm Hg) vs. Temperature and Pressure (Continued)

Barometric Pressure Reading (mmHg)	Gas Temperature (°C)						
	20	21	22	23	24	25	26
757	40.4	40.2	40.1	40.0	39.9	39.7	39.6
758	40.4	40.3	40.2	40.0	39.9	39.8	39.7
759	40.5	40.3	40.2	40.1	40.0	39.8	39.7
760	40.5	40.4	40.3	40.1	40.0	39.9	39.8
761	40.6	40.4	40.3	40.2	40.1	39.9	39.8
762	40.6	40.5	40.4	40.2	40.1	40.0	39.9
763	40.7	40.5	40.4	40.3	40.2	40.0	39.9
764	40.7	40.6	40.5	40.3	40.2	40.1	40.0
765	40.8	40.7	40.5	40.4	40.3	40.1	40.0
766	40.8	40.7	40.6	40.5	40.3	40.2	40.1
767	40.9	40.8	40.6	40.5	40.4	40.3	40.1
768	40.9	40.8	40.7	40.6	40.4	40.3	40.2
769	41.0	40.9	40.7	40.6	40.5	40.4	40.2
770	41.0	40.9	40.8	40.7	40.5	40.4	40.3
771	41.1	41.0	40.8	40.7	40.6	40.5	40.3
772	41.2	41.0	40.9	40.8	40.6	40.5	40.4
773	41.2	41.1	40.9	40.8	40.7	40.6	40.4
774	41.3	41.1	41.0	40.9	40.7	40.6	40.5
775	41.3	41.2	41.1	40.9	40.8	40.7	40.5
776	41.4	41.2	41.1	41.0	40.9	40.7	40.6
777	41.4	41.3	41.2	41.0	40.9	40.8	40.7
778	41.5	41.3	41.2	41.1	41.0	40.8	40.7
779	41.5	41.4	41.3	41.1	41.0	40.9	40.8

The following example reads an expected CO<sub>2</sub> value from the table above, when the pressure reading is 760 mmHg and temperature is 23°C:

**To find the expected CO<sub>2</sub> value**

1. Find the barometric pressure of 760 mmHg in the table.
2. Find the temperature of 23°C.
3. The expected CO<sub>2</sub> value is the intersection of those two values in the table, or 40.1 mmHg.

If for any reason the temperature and pressure values are not listed in the table above, use the following equation to calculate the expected CO<sub>2</sub> value:

$$CO_2(mmHg) = \frac{CO_2\ gas\ \% \times Barometric\ Pressure \times Gas\ Comp}{1 - (0.003 \times (35 - Temperature))}$$

where Gas Comp =

$$1 + (0.000865 \times (setO_2\ \% - O_2\ gas\ \%))$$

CO<sub>2</sub> gas is the percentage of CO<sub>2</sub> in the specified test gas, which is 5% or 0.05.

Pressure Reading is the pressure in mmHg from the waveform display.

Set O<sub>2</sub>% is the set O<sub>2</sub> concentration on the ventilator which, for Gas Comp calculation purposes is 21.

O<sub>2</sub> gas%, which is zero, as there is no O<sub>2</sub> in the specified test gas.

Temperature is the ambient temperature in °C.

The next example calculates the expected CO<sub>2</sub> value using the equation given above with barometric pressure of 759 mmHg and temperature of 19°C.

$$CO_2(mmHg) = \frac{0.05 \times 759 \times 1.018165}{1 - (0.003 \times (35 - 19))}$$

The expected CO<sub>2</sub> value equals 40.59 mmHg.

Dispose of calibration gas per the institution's protocol.



**Note:**

Gas Comp is always calculated to equal 1.018165 using the equation above, because the set O<sub>2</sub>% should always equal 21 (from step 6 above) and the percentage of O<sub>2</sub> in the specified test gas is 0%.

## 1.13 Cleaning the Capnography Sensor



**Note:**

Allow the capnography sensor to cool to room temperature for 30 minutes before attempting to clean it.



**Caution:**

- **Disconnect the capnography sensor from the BDU before cleaning**
- **Only the outside of the capnography sensor can be cleaned. Do not sterilize or immerse the capnography sensor in liquids, as cable and/or sensor damage could result.**
- **Do not attempt to sterilize the capnography sensor.**
- **Allow the capnography sensor to dry completely before reconnecting to the BDU and breathing circuit.**

**To clean the capnography sensor**

1. Wipe with a cloth dampened with any of the cleaning agents listed in the table on page [1-35](#).
2. Wipe with a clean, water-dampened cloth to rinse and dry before use. Ensure that the sensor windows are clean and dry before reuse.

## 1.14 Cleaning and Sterilizing Reusable Airway Adapters



**WARNING:**

**Do not reuse reusable airway adapters on other patients without first cleaning and sterilizing as described below.**



**Caution:**

**Do not insert any object, such as a brush, into the airway adapter, as damage may occur to the windows.**

**To clean a reusable airway adapter**

1. Wash in a warm, soapy solution.
2. Soak in one of the following liquid disinfectants:
  - 70% isopropyl alcohol

- Aqueous 10% solution of sodium hypochlorite (bleach)
  - Glutaraldehyde 2.4% solution such as Cidex®
  - Steris System 1®
3. Rinse with sterile water and dry.

**To sterilize a reusable airway adapter**

1. Clean the airway adapter before sterilizing.
2. Sterilize using one of the following methods:
  - Ethylene Oxide (ETO) 38°C for three (3) hours
  - Steam autoclave (Pediatric/Adult adapters **only**) 121°C or 134°C for 20 minutes
  - Immerse and soak in Cidex or equivalent 3.4% glutaraldehyde solution 20°C ± 5°C for ten (10) hours
  - Immerse and soak in Perasafe™ or equivalent peracetic acid 0.26% solution 20°C ± 5°C for ten (10) hours
  - Cidex® OPA (follow manufacturer's instructions)
3. Ensure the windows are dry and residue free and that the adapter has not been damaged during handling or the cleaning/disinfecting process.

**WARNING:**

**To avoid the risk of bacterial contamination, do not attempt to clean, sterilize, or re-use single-use components. These sensors and adapters are not compatible with sterilization techniques.**

## 1.15 Alarms

Three alarm events are associated with the ETCO<sub>2</sub> feature:

- High ETCO<sub>2</sub> level ( $\uparrow$ ETCO<sub>2</sub>) where the measured ETCO<sub>2</sub> level is above the operator-set value.
- Low ETCO<sub>2</sub> level ( $\downarrow$ ETCO<sub>2</sub>) where the measured ETCO<sub>2</sub> level is below the operator-set value.
- Capnography sensor inoperative (a malfunction has occurred with the capnography sensor).

Reference *Chapter 4* for information on setting alarms. The  $\uparrow$  and  $\downarrow$  ETCO<sub>2</sub> alarms may be set during ventilator operation, only.

If any of the alarm conditions occur, the GUI displays an alarm message similar the one shown below. Follow the information contained in the remedy message and the prompt area to troubleshoot the alarm.

**Figure 1-12.** Alarm Message for Capnography Sensor Inoperative

VEN\_11538\_A

## 1.16 Messaging

**Messages are displayed either in the Prompt area of the GUI (Reference the figure Areas of the GUI in Chapter 4 of this manual), in the capnography screen, or in the alarm banner.**

**Examples of messages include:**

- CO<sub>2</sub> sensor ready
- CO<sub>2</sub> sensor warming up
- CO<sub>2</sub> sensor zeroing in progress
- CO<sub>2</sub> sensor zeroing passed
- Zeroing attempt rejected, verify procedure
- Capnography is currently disabled

The following messages occur in conjunction with the Capnography inoperative alarm banner:

- CO<sub>2</sub> sensor not connected
- Check CO<sub>2</sub> sensor adapter
- Invalid CO<sub>2</sub> sensor

- CO<sub>2</sub> sensor faulty
- CO<sub>2</sub> sensor not ready
- Zeroing failed. CO<sub>2</sub> sensor zeroing required

## 1.17 CO<sub>2</sub> Specifications

The following tables describe the specifications for the capnography sensor.

**Table 1-6.** Ventilator Patient Data Displayed Range, Resolution, and Accuracy

Parameter	Range, Resolution, Accuracy
CO <sub>2</sub>	Range: 0 mmHg to 150 mmHg Resolution: 0.1 mmHg for values 0 mmHg to 69 mmHg; 0.25 mm Hg for values 70 mm Hg to 150 mmHg Accuracy: 0 mmHg to 40 mmHg ± 2 mmHg 41 mmHg to 70 mmHg ± 5% of reading 71 mmHg to 100 mmHg ± 8% of reading 101 mmHg to 150 mmHg ± 10% of reading Temperature at 35°C (No degradation due to respiratory rate or I:E ratio.) Accuracy (short term drift): ≤ 0.8 mmHg over four hours Accuracy (long term drift): Accuracy specification maintained over 120 hours

**Table 1-7.** Ventilator Alarm settings Range and Resolution

Setting	Description	Range and Resolution
High end tidal CO <sub>2</sub> alarm (↑ETCO <sub>2</sub> )	The ↑ETCO <sub>2</sub> alarm indicates the measured end tidal CO <sub>2</sub> level is ≥ the set alarm limit.	Range: OFF or 10.0 mmHg to 150 mmHg and > ↓ETCO <sub>2</sub> alarm limit Resolution: 1.0 mmHg
Low end tidal CO <sub>2</sub> alarm (↓ETCO <sub>2</sub> )	The ↓ETCO <sub>2</sub> alarm indicates the measured end tidal CO <sub>2</sub> level is ≤ the set alarm limit.	Range: OFF or 5 mmHg to 60 mmHg and < ↑ETCO <sub>2</sub> alarm limit Resolution: 1.0 mmHg

**Table 1-8.** Capnography Sensor Specifications

Parameter	Description	Specification
System response time	Time until 90% of step change of final CO <sub>2</sub> value is displayed	< 60 ms
Sampling rate	Frequency of CO <sub>2</sub> measurements	100 Hz

**Table 1-8.** Capnography Sensor Specifications

Parameter	Description	Specification
Transducer type		Mainstream
Principle of operation	N/A	Non-dispersive infrared (NDIR) single beam optics, dual wavelength, no moving parts.
CO <sub>2</sub> calculation method	Peak of the expired CO <sub>2</sub> waveform	BTPS
Water resistance	Category for liquid ingress	IPX4 (splash proof)
Interfering gases	Accuracy not affected by presence of specified gas concentrations	0.1% ethanol 0.1% isopropanol 0.1% acetone 1% methane
Non-condensing humidity	N/A	Accuracy not affected
Respiratory rate (f)	Number of breaths per minute set by the operator	Range: 0 breaths/min to 150 breaths/min Accuracy: $\pm 1$ breath Accuracy not affected. No limitation on respiratory rate when monitoring CO <sub>2</sub> .
I:E ratio	Ratio of inspiratory time to expiratory time	Accuracy not affected

**Note:**

The capnography sensor is equipped with automatic barometric pressure compensation.

**Table 1-9.** Neonatal/Pediatric Single Use Airway Adapter Specifications

Parameter	Description	Specification
Dead space	The volume of air not included in the ETCO <sub>2</sub> measurement	< 1 mL
Weight	N/A	9.1 g
Pressure drop	The difference in pressure from the inlet to the outlet at a specified flow rate	0.74 cmH <sub>2</sub> O at 10 L/min
Color	N/A	Purple

**Table 1-10.** Neonatal/Pediatric Reusable Airway Adapter Specifications

Parameter	Description	Specification
Dead space	The volume of air not included in the ETCO <sub>2</sub> measurement	< 1mL
Weight	N/A	14.9 g
Pressure drop	The difference in pressure from the inlet to the outlet at a specified flow rate	0.68 cmH <sub>2</sub> O at 10 L/min
Color	N/A	Red

**Table 1-11.** Pediatric/Adult Single Use Airway Adapter Specifications

Parameter	Description	Specification
Dead space	The volume of air not included in the ETCO <sub>2</sub> measurement	<6 mL
Weight	N/A	7.7 g
Pressure drop	The difference in pressure from the inlet to the outlet at a specified flow rate	0.40 cmH <sub>2</sub> O at 60 L/min
Color	N/A	Colorless, transparent

**Table 1-12.** Pediatric/Adult Reusable Airway Adapter Specifications

Parameter	Description	Specification
Dead space	The volume of air not included in the ETCO <sub>2</sub> measurement	<6 mL
Weight	N/A	12.0 g
Pressure drop	The difference in pressure from the inlet to the outlet at a specified flow rate	0.38 cmH <sub>2</sub> O at 60 L/min
Color	N/A	Black

**Table 1-13.** Flow and CO<sub>2</sub> Combination Sensor Specifications

Parameter	Description	Specification
Dead space	The volume of air not included in the ETCO <sub>2</sub> measurement	< 1 mL
Weight	Weight of the Flow and CO <sub>2</sub> combination sensor	9.6 g
Pressure drop	The difference in pressure from the inlet to the outlet at a specified flow rate	3.1 cmH <sub>2</sub> O at 10 L/min

**Table 1-13.** Flow and CO<sub>2</sub> Combination Sensor Specifications (Continued)

Parameter	Description	Specification
Color	N/A	Purple, transparent

**Note:**

The Flow and CO<sub>2</sub> combination sensor is for neonatal use, only.

## 1.18 Part Numbers

The following table lists the part numbers for the ETCO<sub>2</sub> monitoring function individual components.

**Table 1-14.** Capnography Monitoring System Part Numbers

Item	Part Number
Capnography sensor	10087409
Airway Adapter, single-use pediatric/adult (package of 10) (colorless transparent)	10078387
Airway Adapter, single-use neonatal/pediatric (package of 10) (purple transparent in color)	10078386
Airway Adapter, reusable pediatric/adult (black in color)	10083942
Airway Adapter, reusable neonatal/pediatric (red in color)	10083943
Flow and CO <sub>2</sub> Combination Sensor, single-use (package of 10)	10005002
CO <sub>2</sub> sensor cable (internal)	PT00088615

## 1.19 Cleaning Agents

The following table lists cleaners that are compatible with the Capnography sensor.

**Table 1-15.** Capnography Sensor Compatible Cleaning Agents

Cleaning Agent
Isopropyl alcohol 70%
A 10% aqueous solution of 6% chlorine bleach
Steris Coverage® Spray HB
Clinell Wipes®
PDI Sani Cloth Bleach®
PDI Super Sani Cloth AF®
Speedy Clean™
Tuffie™
Tuffie 5™
Accel TB RTU
Bacillol 30 Foam
Bacillol AF
Caltech-Dispatch 5200 Wipes
Hydrogen Peroxide
Meliseptol
Metrex CaviWipes 1
Oxivir TB Wipes
PDI Sani-Cloth Plus Germicidal Wipes
PDI Super Sani Cloth
Revital-Ox Resert™ XL HLD
Sporox™ II Sterilizing & Disinfecting Solution
Viraguard® Isopropanol 70%
Virex Tb
Wipes Plus® Disinfecting Wipes (Item # 74402)

Page Left Intentionally Blank



**Rx**  
ONLY

CE  
0123

 COVIDIEN

Part No. PT00097188 A 2018-12

© 2013 Covidien.

 Covidien IIc  
15 Hampshire Street, Mansfield, MA 02048 USA  
 EC REP Covidien Ireland Limited, IDA Business and Technology Park Tullamore, Ireland.  
[www.covidien.com](http://www.covidien.com)