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
PCR Workstation *Pro*

CONTENTS

INTRODUCTION	2
WARRANTY	2
DELIVERY CONTENTS	3
SAFETY INSTRUCTIONS	4
ASSEMBLING THE PCR WORKSTATION <i>PRO</i>	6
SPECIFICATIONS	10
PEQLAB PCR Workstation <i>Pro</i>	10
Power supply	10
CONNECTIONS	11
AC Power inlet	11
AC power outlets	12
OPERATING THE PCR WORKSTATION <i>PRO</i>	13
Use of PCR equipment inside the workstation	13
Operating the UV and white lights	13
Timer	13
Turning on UV decontamination	14
White light and UV Air Recirculation	14
Changing the decontamination time	14
Setting up a daily decontamination routine	15
Changing date and time	16
Checking the operating hours of the surface UV tubes	16
Resetting the operating hours of the surface UV tubes	17
Checking the operating hours of the UV Air Recirculator UV tube	17
Resetting the operating hours of the UV Air Recirculator UV tube	18
Overhead white light	19
UV Air Recirculator	17
Measuring the UV tube intensity	20
Bacterial destruction chart	21
MAINTENANCE	22
Care and cleaning	22
Replacing the UV germicidal or white light tubes	22
Replacing the UV tube of the UV Air Recirculator	23
Replacing the fuse	23
Replacing the filter of the UV Air Recirculator	23
TECHNICAL SERVICE AND ORDER INFORMATION	24
DISPOSAL INSTRUCTION	24

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INTRODUCTION

The PEQLAB PCR Workstation *Pro* is designed as an ideal environment for the manipulation of DNA and RNA, especially for the setup of PCR assays. Contaminations can lead to false or misleading results which costs time and money. The PEQLAB PCR Workstation *Pro* minimizes the risk of contamination. It provides a 'separated room', e.g. for the setup of PCR reactions.

The high intensity surface UV tubes inactivate DNA as a source of contamination between experiments. In addition, an UV Air Recirculator is integrated into the workstation system which reduces airborne contaminants during experiments. UV light is effectively blocked by solid polycarbonate screens to ensure maximum protection of the user.

The surface of the workstation consists of stainless steel, which is very robust, can easily be cleaned and has an antimicrobial effect.

The PEQLAB PCR Workstation *Pro* offers a controlled environment for PCR and RNA applications that protects your samples and helps to achieve optimal results.

WARRANTY

VWR warrants its products to be free of defects in materials and workmanship for a period of two years from date of purchase. The foregoing warranty of VWR shall be of no force and effect if buyer has modified or damaged the product. Tubes and filters are warranted for 90 days. crazing of the polycarbonate panels may occur within the warranty period and is regarded as normal wear and tear not covered by the warranty.

All warranties or merchantability and fitness for any purpose and all other warranties, expressed or implied, except those expressly set forth herein, are deemed waived and excluded.

VWR's duty under the warranty is limited to replacement and/or repair of the defective part at the option of VWR. VWR shall not be liable for any expenses or damages incurred by the purchaser except as expressly set forth herein, and in no event shall VWR be liable for any special, incidental or consequential damages of any kind. This warranty does not supersede any statutory rights that may be available in certain countries.

For research use only. Not for use in veterinary or human medicine or diagnostic procedures.

DELIVERY CONTENTS


After receipt of the system, please check if the following items are included in the shipment:


① Base (1x)


② Rear panel (1x)

③ Side panel (2x)


④ Front panel (1x)


⑤ Shelf (2x)

⑥ Top assembly (1x)

⑦ Ventilation cover (1x)

⑧ Screw (31x)

⑨ Screwdriver (1x)

⑩ Mains cable (1x)

SAFETY INSTRUCTIONS

The PCR Workstation *Pro* is designed with function, reliability, and safety in mind. Before using the PCR Workstation *Pro* please read this manual entirely and strictly follow the safety instructions!

Caution: The unit includes shortwave UV, which is a powerful source of UV radiation that will cause damage to unprotected eyes and skin if directly exposed to the UV lamps. Before operating any unit, ensure that the PCR Workstation *Pro* is not damaged, and that instructions for use of this equipment are followed to. A safety shut-off switch automatically turns the UV light off when the front panel is open, protecting users from UV exposure. The polycarbonate panels are UV-blocking.

Intended use of the PCR Workstation *Pro*: Working environment for protection of PCR instruments and genetic substances against pollution and contamination.

Safety and protection of the user can not be ensured if the PCR Workstation *Pro* is used in any manner that is not specified in this instruction. Special notice must be taken of the following instructions:

Caution of dangerous voltage

Ensure that the voltage requirements of the PCR Workstation *Pro* exactly match your local AC power supply.
Only use the original AC cable supplied with the PCR Workstation *Pro*.

Caution of liquids

Ensure that no liquids can enter the device. Do not place liquids on top of the PCR Workstation *Pro*.

Caution of dangerous explosive material

Do not use reactive or explosive substances in the PCR Workstation *Pro*.

Caution of damages

Check for transport damages and completeness of the accessories upon arrival. If any damage is noticed, contact your distributor or manufacturer immediately. If you detect any damages of the PCR Workstation *Pro* or one of its accessories, do not connect the PCR Workstation *Pro* to the AC wall outlet and do not use the device in any manner!

Caution of optical radiation


Do not switch on the PCR Workstation *Pro* if not completely assembled. See "ASSEMBLING THE PCR WORKSTATION *PRO*", page 5.
Do not look directly into the UV lamps.

ASSEMBLING THE PCR WORKSTATION *PRO*


The PCR Workstation *Pro* should be assembled by two people, directly at the site of operation. We re-recommend to wear safety gloves for the protection of hands.

Step 1


Put the base ① on an even surface with the rubber feet facing down and the sensor cable outlet facing to the front. Lay the back panel ② onto the base ① with the power outlets facing down, and the power input and the type label facing to the right.


Connect both sensor cables.

Slide the connection assembly into the cavity of the back panel.

Assemble the back panel ② onto the base ①. The connecting link of the base ① must fit in the back panel ②. **Pay attention not to squeeze the sensor cable.**



Connect both parts by screwing 3 screws ⑧ into the back side. Use the screw driver ⑨ to do this. Turn the parts around 90° on the site of operation and incline both screwed parts to the back part of the back panel ②. Screw both panels with 2 screws ⑧ from underneath. Tilt back the panels in starting position with the back panel facing backwards.

Step 2

Remove the protective foil from the side panels.

Insert the side panels ③ into the guiding rails at the base ① and push up against the back panel ② within the guiding rails at both sides. Carefully screw the side panels ③ with 3 screws each in the base ① and the back panel ②, respectively. Use the supplied screws ⑧ and the screwdriver ⑨.

Step 3

Put the shelves ⑤ into the shelf brackets at the back panel ②.

Side panel ③

Remove the protective foil from the front panel.

Put the front panel ④ onto the side panels ③.

Shelves ⑤

Take care that the screw holes of the front panel match the threaded holes in the side panels ③. Hold the front panel ④ in place. Screw the front panel ④ to both side panels ③ using 3 screws ⑧ each. Use the screwdriver ⑨.

Front panel ④


Side panel ③

Step 4


Turn the assembled parts around 180°.

Remove the small UV tube at the back of the top assembly.

Put the top assembly ⑥ onto the panels. Take care not to squeeze the cables protruding from the back panel. Screw the top assembly to the back panel ② at two positions. Use the screws ⑧ and screwdriver ⑨ supplied.


Connect the sensor cables from the back panel ② and the top assembly ⑥. Put the power plug of the back panel ② into the power outlet of the top assembly ⑥.

Step 5

Insert the UV tube back in place.
Cover the service access of the top assembly ⑥ with the ventilation cover ⑦. Screw both parts together with 6 screws ⑧. Use the screw-driver ⑨.
Turn the system around 180° and place it on the site of operation as needed.

Finally connect the system to the power supply by plugging the power plug ⑩ from the back panel ② into the local power outlet.

The system is now ready for use.

SPECIFICATIONS

PEQLAB PCR Workstation *Pro*

Light Sources integrated into the top assembly:

- 2 x 25 W 254 nm UV light tubes – for effective surface decontamination
- 1 x 8 W 254 nm UV light tube – for decontamination of air during work in UV Air Recirculator
- 1 x 15 W white light tube – for a bright work surface

Features:

- Four power outlets
- Two shelves made from stainless steel
- Electronic UV timer
- Working area and back panel made from stainless steel
- Front and side panels made from polycarbonate (Makrolon®) that blocks wavelengths below 400nm

Dimensions (W x H x D):	750 x 780 x 620 mm
Weight:	48 kg
Pollution rate:	2
Environmental temperature range:	+10 °C to +30 °C
Maximum relative humidity:	70 %
Maximum height above sea level:	2000 m
Noise level:	< 51 dB (A)
Circulation rate fan:	56 m ³ /h


Power supply

EU/UK version:	220 - 240 V AC, 50/60 Hz, max. 1650 VA
US version:	120 V, 60 Hz, max. 900 VA
Fuse:	10 AT, 5 x 20 mm

For indoor use only!

CONNECTIONS

AC Power inlet

Before the first use, ensure that the voltage requirements of the PCR Workstation *Pro* exactly match your local AC power supply. Connect the PCR Workstation *Pro* to the AC power supply only by the AC power cord delivered with your PCR Workstation *Pro*. Only connect the PCR Workstation *Pro* to electric power systems with an earth conductor. The AC power inlet and the AC power switch are located at the rear side of the workstation. The AC power switch serves as disconnection of the AC power supply from the device. The accessibility of the AC power switch must be assured and not hindered.

The electric fuse is placed at the back of the device between the power switch and the AC power inlet. After pulling off the AC power cord and pushing upwards a safety catch, the fuse can be changed. Use only a fuse with correct values (see "SPECIFICATIONS", page 8). Indications about the fuse type are also located at the back of the device.


If service is required, contact your local distributor or the manufacturer (see "TECHNICAL SERVICE AND ORDER INFORMATION", page 22).

The serial number of the PCR Workstation *Pro* is located at the rear side of the instrument. The replacement of single components of the device (except polycarbonate panels, UV tubes, fuse and filter of the UV Air Recirculator) must not be done by the user but exclusively by authorized specialists using original replacement parts. Therefore the device must be sent to the manufacturer. Unauthorized work on the device voids the warranty.

The surface of the device should be decontaminated with a lint free cloth soaked with 70 % ethanol before sending. Please fill in and sign the decontamination certificate and send it to the manufacturer together with the device.

For safety reasons the device must be switched off and the mains plug pulled before cleaning is performed!

AC power outlets

The 4 AC power outlets serve as power supplies for devices used inside the PCR chamber.

The sum of the maximum power ratings of devices connected to these power outlets must not exceed the values specified for the PCR Workstation *Pro*. See chapter "OPERATING THE PCR WORKSTATION *PRO*", page 11, for further information.

OPERATING THE PCR WORKSTATION *PRO*

Use of PCR equipment inside the workstation

Shelves are provided for placement of small tools for storage and decontamination.

Power outlets allow use of shakers, rockers, and other equipment for PCR experiments inside the PCR chamber. The equipment can be decontaminated between experiments. To operate equipment within the chamber, plug the equipment into the power plug. Specific power plugs are installed as shown in the illustration "AC power outlets", page 10.

Timing of the UV decontamination can be programmed for time-controlled every day decontamination.

NOTE: When operating equipment plugged into the outlets, combined Amps total is not to exceed 6 Amps.

NOTE: Crazeing is a normal process for polycarbonate panels exposed to UV light. Crazeing will occur over a period of time. Reduce crazeing by keeping exposure to UV to a minimum. Crazeing may occur within the warranty period and is regarded as normal wear and tear not covered by the warranty. Polycarbonate panels are available as spare parts and can be replaced.

Operating the UV and white lights

The PCR Workstation *Pro* provides an automated process for eliminating contamination by utilizing the decontaminating properties of shortwave 254 nm ultraviolet light. The UV light is typically operated when no samples are inside the chamber. The chamber is decontaminated and DNA as possible contamination is inactivated.

Timer

The default setting is 30 minutes adjustable at 5 minutes increments down to a minimum of 5 minutes. For information on setting the timer please see page 13.

Note: A magnetic sensor is built into the bottom of the front panel. When the front panel is opened, the ultraviolet light will automatically shut off!

Note: Do not attempt to perform PCR procedures with the germicidal lamp on. The germicidal lamp is used to decontaminate the chamber between experiments.

Turning on UV decontamination


To turn on UV decontamination of the workspace push the << **UV Light** >> button. It automatically activates the program to decontaminate the workspace. The default time for the decontamination process is 30 min.

The display will show the remaining time of the decontamination program.

White light and UV Air Recirculation

- To turn on white light push the << **White Light** >> button. A red light will light up to indicate that the white lights are turned on.
- To turn off the white light push the << **White Light** >> button again.
- To turn on the UV Air Recirculator that decontaminates the air during work push << **Air Recirculation** >>. A red light will light up to indicate that the UV Air Recirculator is turned on.
- To turn off the UV Air Recirculator push the << **Air Recirculation** >> button again.

Changing the decontamination time

- To run a shorter decontamination press << **Mode** >>.
- The screen will show "Decontaminated". Time can be changed from 5 to 30 min in 5 min steps by pushing the << + >> and << - >> buttons.
- Push << **Set** >> to set the default UV light decontamination time to the indicated value.

Setting up a daily decontamination routine

- To set up a daily start program that starts a decontamination program every day at the same time push << **Mode** >> two times until the display shows "Daily Start".
- Push << **Set** >> to change the settings.
- By pushing << + >> the program can be turned "ON".
- To confirm push << **Set** >>.
- The cursor will automatically switch to set up the minutes.
- Change minutes by pushing << + >> and << - >>.
- Push << **Set** >> to confirm.
- The cursor will jump to hours.
- Change the hour to turn on the PCR Workstation *Pro* by pushing << + >> and << - >>.
- Push << **Set** >> to confirm.

The PCR Workstation *Pro* will now start a decontamination program every day at the set time (e.g. for 30 minutes at 10:30)

Changing date and time

The PEQLAB PCR Workstation *Pro* has an integrated date and time function.


- To change the date and time push the << Mode >> button three times until you see the time and date.
- Push << Set >> to change the seconds using the << + >> and << - >> buttons.
- Push << Set >> to change the minutes using the << + >> and << - >> buttons.
- Push << Set >> to change the hours using the << + >> and << - >> buttons.
- Push << Set >> to change the day using the << + >> and << - >> buttons.
- Push << Set >> to change the month using the << + >> and << - >> buttons.
- Push << Set >> to change the year using the << + >> and << - >> buttons.
- Push << Set >> to confirm.

Checking the operating hours of the surface UV tubes

To see the operating hours of the light tubes push the << **Mode** >> button 4 times. We recommend to change the UV tubes after 1000 – 1500 operating hours as the UV intensity decreases significantly after that time.

Resetting the operating hours of the surface UV tubes

After renewing the surface UV tubes the operating hours should be set back to zero.

To reset the time

- Push << **Mode** >> 4 times until you see the operating hours screen.
- Push << **Set** >>. You will see the Reset Surface Screen.
- To reset the operating hours push << + >>.
- To go back without resetting the operating hours push << - >>.


Checking the operating hours of the UV Air Recirculator UV tube

To see the operating hours of the UV Air Recirculator UV tube push the << **Mode** >> button 5 times.

We recommend to change the UV tube after 1000 – 1500 operating hours as the UV intensity decreases significantly after that time.

Resetting the operating hours of the UV Air Recirculator UV tube

After renewing the UV Air Recirculator UV tube the operating hours should be set back to zero.

To reset the time

- Push << **Mode** >> 5 times until you see the operating hours screen.
- Push << **Set** >>. You will see the Reset Surface Screen.
- To reset the operating hours push << + >>.
- To go back without resetting the operating hours push << - >>.

[Warning] Exposure to UV light is harmful.

The polycarbonate front panel is UV blocking and will not allow UV radiation to pass through the panel. The UV lamp will shut off if the front panel is opened. The PCR Workstation *Pro* contains a powerful source of UV radiation that will cause damage to unprotected eyes and skin. Before operating any unit, be sure all personnel in the area are properly protected. Even though the unit shuts the UV off when the front panel is open, UV blocking eyewear should be worn as well.

Extended periods of time UV exposure can be used to decontaminate the interior chamber by inactivation of DNA. Refer to "Measuring the UV tube intensity", page 18, for information on measuring the intensity of the UV tubes with a UV meter.

Overhead white light


The overhead white light can remain lit at any time throughout experiment procedures or the decontamination process. The fluorescent white light provides a bright light within the PCR Workstation *Pro*.

UV Air Recirculator

The PEQLAB PCR Workstation *Pro* includes a built-in UV Air Recirculator for elimination of air-borne microbes during experiments.

NOTE: The UV Air Recirculator can operate while the front panel is open.

The UV Air Recirculator circulates airflow into and out of the chamber. The air is sucked into the UV Air Recirculator in the top of the chamber, the air outlets are at the bottom of the back panel. The drawing indicates the location of the air vents.


NOTE: Equipment in front of the air vents will block air flow.

Measuring the UV tube intensity

As UV tubes age, the intensity and germicidal destruction rate decreases. We recommend to change the UV tubes after 1000 – 1500 hours of operation (please refer to "Checking the operating hours of the surface tubes/UV Air Recirculator tube", page 14/15). Otherwise it is important to monitor the efficiency of the tubes to ensure the germicidal requirements are met. The germicidal destruction rate is a function of the UV intensity at wavelength 254 nm and exposure time. The lower the tube intensity, the longer the lamps must be on to accomplish the same objective.

An UV Intensity Meter with 254 nm sensor will allow the user to measure the 254 nm emissions from the tubes. When it is time to measure the UV intensity inside the PCR Workstation *Pro*, the sensors can be placed on the floor of the PCR Workstation *Pro*. Close the front panel and turn on the UV tubes. Please refer to the manual of your UV radiometer on how to perform the measurement. It is recommended that the user records an initial value upon receipt and assembly of the PCR Workstation *Pro*.

To measure the UV intensity from the UV Air Recirculator tube, open the cover lid on the back side of the PCR Workstation *Pro* and perform the measurement as described on the manual of the UV radiometer. Wear protective eyewear and shielding when operating the UV Air Recirculator with the ventilation cover open. Replace the cover and screws after the measurement is completed.

Bacterial destruction chart

The germicidal destruction rate calculation:

$$\text{Microwatt seconds/cm}^2 = \text{microwatts/cm}^2 \times \text{seconds of exposure}$$

The bacterial destruction chart below indicates the amount of shortwave (254 nm) UV energy required for complete destruction of various organisms.

Bacteria	Microwatt seconds/cm ²	Other organisms	Microwatt seconds/cm ²
Bacillus anthracis	8700	YEAST	
S. enteritidis	7600	Saccharomyces ellipsoideus	13200
B. Megatherium sp. (veg.)	2500	Saccharomyces sp.	17600
B. Megatherium sp. (spores)	5200	Saccharomyces cerevisiae	13200
B. paratyphosus	6100	Brewer's yeast	6600
B. subtilis	11000	Baker's yeast	8800
B. subtilis spores	22000	Common yeast cake	13200
Clostridium tetani	22000		
Corynebacterium diptheriae	6500	MOLD SPORES	
Eberthella typosa	4100	Penicillium roqueforti	26400
Escherichia coli	6600	Penicillium expansum	22000
Micrococcus cadidus	12300	Penicillium digitatus	88000
Micrococcus sphaeroides	15400	Aspergillus glaucus	88000
Mycobacterium tuberculosis	1000	Aspergillus flavus	99000
Neisseria catarrhalis	8500	Aspergillus niger	330000
Phytomonas tumefaciens	8500	Rhizopus nigricans	220000
Proteus vulgaris	6600	Mucor racemosus A	35200
Pseudomonas aeruginosa	10500	Mucor racemosus B	35200
Pseudomonas fluorescens	6600	Oospora lactis	11000
S. typhimurium	15200		
Salmonella	10000	VIRUS	
Sarcina lutea	26400	Bacteriophage (E. coli)	6600
Serratia marcescens	6160	Tobacco mosaic	44000
Dysentery bacilli	4200	Influenza	6600
Shigella paradysenteriae	3200		
Spirillum rubrum	6160	PROTOZOA	
Staphylococcus albus	5720	Paramecium	200000
Staphylococcus aureus	6600	Nematode eggs	9200
Streptococcus hemolyticus	5500	Chlorella vulgaris (algae)	22000
Streptococcus lactis	8800		
Streptococcus viridans	3800		

MAINTENANCE

NOTE: Polycarbonate is aging in UV light, therefore crazing is considered as a normal process. Crazing will occur with time. Reduce crazing by keeping exposure to UV to a minimum. Crazing may occur within the warranty period and is regarded as normal wear and not covered by the warranty. Polycarbonate panels can be replaced.

Care and cleaning

The front and side panels are made of polycarbonate and are subject to scuffing and scratches if improperly cleaned. Therefore it is strongly recommended to use a soft cloth (preferentially microfiber) for cleaning.

- Wipe excess water from inside the unit and outside the unit with an absorbent soft cloth or sponge.
- Use mild detergent and water with a soft cloth or sponge to clean the exterior and interior of the unit.
- Clean the front and side panels with a mild detergent only.
- Do not use abrasive pads or cleansers.
- A plastic cleaner solution is recommended to clean the front panel and is available from local plastic supply distributors.

Never use organic based compounds like acetone or ammonia containing cleaners for polycarbonate panels!

Due to the specific qualities of polycarbonate panels made of Makrolon® 70 % ethanol can be used for cleaning.

All metal parts consist of stainless steel V2A. Stainless steels have sufficient corrosion protection against normal environmental influences. It is strongly recommended to use a stainless steel cleaning and care spray. This cleans not only, but also leaves a protective layer. Do not use aggressive cleaning agents.


Warning: Unplug the unit before cleaning around the UV tube contacts and internal plug outlets.

Replacing the UV germicidal or white light tubes

Turn the unit power off and unplug the unit from the outlet. Lift up the front panel and locate the tubes. The tube has two prongs on each end that first fit into a socket, then rotate into place.

Carefully hold the tube on the end and insert the tube into the socket. Twist the tube 90° until the tube locks into place.

Drawing shows tube location in top assembly. Outer tubes are the shortwave 254 nm UV and the middle tube is white light.

Replacing the UV tube of the UV Air Recirculator


Locate and remove the six screws on the ventilation cover on the outer back panel of the workstation.

Pull away the ventilation cover from the workstation.

The tube has two prongs on each end that fit into a socket, then rotates into place.

Carefully hold the tube on the ends and insert the tube into the socket. Twist the tube 90° until the tube locks into place.

Replace the cover back into position and reattach the six screws.

Replacing the fuse


The fuse is located at the back right side of the system. To remove the fuse, use a flathead screwdriver to open the lid and pull out the fuse.

Insert the new fuse and using the screwdriver, turn the fuse clockwise to lock.

Replacing the filter of the UV Air Recirculator

The fan filter of the UV Air Recirculator is placed on the inner top of the PCR Workstation *Pro* in front of the light tubes (see picture to the right). To change the filter, take off the plastic cover with the filter attached to it and replace it with a new filter. Plunge the plastic cover back into its fixation.

TECHNICAL SERVICE AND ORDER INFORMATION

Should you have any questions about the PCR Workstation *Pro* or its accessories, please visit www.de.vwr.com to find the respective contact person or send an email to service.peqlab@de.vwr.com.

Please have the unit's serial number (located at the back side of the instrument) available when calling. Should an item require return to VWR for service, a decontamination form must be completed first by the user. Items sent without decontamination form will not be accepted.


All returns must be pre-approved by VWR!

Order information:

Article	Description	Order number
PCR Workstation <i>Pro</i>	PCR Workstation <i>Pro</i> with UV Air Recirculator, DE & AT	732-2840
PCR Workstation <i>Pro</i>	PCR Workstation <i>Pro</i> with UV Air Recirculator, export version	PEQL90-UV/PCR2-EU
PCR Workstation <i>Pro</i>	PCR Workstation <i>Pro</i> with UV Air Recirculator, UK version	PEQL90-UV/PCR2-UK
UV tubes	254 nm UV tubes for PCR Workstation <i>Pro</i> , 2 tubes	PEQL90-UV/PCR-TUV
UV tube	254 nm UV tube for UV Air Recirculator, 1 tube	PEQL90-UV/PCR-TAR
White light tube	White light tube for PCR Workstation <i>Pro</i> , 1 tube	PEQL90-UV/PCR-TWL

DISPOSAL INSTRUCTION

This product is subject to the WEEE disposal regulations and cannot be disposed of with regular waste.

Please contact VWR International GmbH or your local distributor for disposal of the PCR Workstation *Pro*. Old electronic equipment is not useless waste. Environmentally friendly disposal can help to retrieve valuable resources. Thereby you contribute to the protection of the environment and human health.

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