

[RETURN TO THIS MANUAL'S TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

Drägermedical

A Dräger and Siemens Company

Technical Service Manual

Part Number: 4115138-001

Rev: R

Date: 5 May 2004

© 2004 Draeger Medical, Inc.

Narkomed Mobile Anesthesia System

[RETURN TO THIS MANUAL'S TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

Narkomed Mobile Service Manual Table of Contents

What's New in Rev. R

DESCRIPTION	PAGE
Section 1: Introduction	
General Safety Precautions	1-1
1.1 Purpose	1-2
1.2 Recommendations	1-2
1.3 General Troubleshooting Guidelines	1-3
1.4 Related Publications	1-3
1.5 Symbol Definitions	1-3
Section 2: Diagnostics	
2.1 Main Service Screen	2-2
2.2 Service Log	2-5
2.3 PMS Criteria Screen	2-6
2.4 Oxygen Monitor Service Screen	2-7
2.5 Pressure Monitor Service Screen	2-8
2.6 Serial Port Configuration Screen	2-9
Section 3: Troubleshooting	
3.1 Power Supply and Voltage Distribution	3-1
3.2 Battery	3-3
3.3 Troubleshooting Guides	3-3
Section 4: Replacement Procedures	
4.1 Cylinder Yokes and Regulators	4-2
4.2 Cylinder and Pipeline Pressure Gauges	4-5
4.3 Flowmeters	4-9
4.4 Auxiliary O2 Flowmeter	4-12
4.5 Flow Control Valves	4-15
4.6 Oxygen Supply Failure Protection Device	4-18
4.7 Oxygen Supply Pressure Alarm Switch	4-21
4.8 O2 - Air Switch	4-23
4.9 System Power Switch	4-25
4.10 Oxygen Ratio Controller	4-28
4.11 O2 Flush Valve	4-31
4.12 Ventilator Controller (Bezel Assembly)	4-35
4.13 Bellows Valve Assembly	4-37
4.14 Caster Replacement	4-39
4.15 Battery Replacement	4-41
4.16 Power Supply Assembly Replacement	4-43
4.17 Processor Assembly	4-47
4.18 Display Assembly	4-49
4.19 Breathing Pressure Hose Assembly	4-51
4.20 Manual/Automatic Selector Valve	4-53

DESCRIPTION	PAGE
4.21 Ultrasonic Flow Sensor	4-55
4.22 Oxygen Sensor	4-57
Section 5: Adjustment and Calibration Procedures	
5.1 Cylinder Pressure Regulator Adjustment	5-2
5.2 Oxygen Supply Pressure Alarm Switch Adjustment	5-4
5.3 Oxygen Ratio Controller (ORC) Adjustment	5-6
5.4 Oxygen Sensor Calibration	5-8
5.5 Breathing Pressure Monitor Calibration	5-11
Section 6: PMC Procedure, Narkomed Mobile	
6.1 Self-Diagnostics	6-12
6.2 Electrical Safety- One Year Service Interval; Due Date _____	6-12
6.3 Configuration	6-14
6.4 Service Data	6-14
6.5 Calibrations - One Year Service Interval; Due Date _____	6-15
6.6 Absorber Inspection	6-16
6.7 High Pressure Leak	6-17
6.8 Breathing System	6-18
6.9 Oxygen Analyzer	6-21
6.10 Flowmeters & Concentrations	6-22
6.11 High Pressure Regulators - Six Month Service Interval; Due Date _____	6-24
6.12 Low O2 Supply – Six Month Service Interval; Due Date _____	6-25
6.13 Oxygen Supply Failure Protection	6-26
6.14 Pressure Monitor	6-26
6.15 Ventilator	6-27
6.16 Bellows Adult	6-29
6.17 Scavenger- Six Month Service Interval; Due Date _____	6-30
6.18 Final Tests	6-32

NM MOBILE	CONTENTS (continued)
------------------	-----------------------------

DESCRIPTION	PAGE
Section 7: Update Procedure	
7.1 Software Transfer to PC Via Modem	7-1
7.2 Installing Narkomed Mobile Software from a PC	7-1
 Section 8: Spare and Replacement Parts	
Vaporizer & rel. parts, fresh gas outlet, mainswitch, flowmtr sub-asm	8-2, 8-3
Processor Assembly and related parts	8-4, 8-5
Ventilator Controller (Bezel Assembly).	8-6, 8-7
Bellows Valve Assembly.	8-8, 8-9
Bellows Valve Assembly Details	8-9A, 8-9B
Display Assembly and Monitor Support Arm	8-10, 8-11
Pipeline Inlet (Manifold) Assemblies.	8-12, 8-13
Failsafe (OFPD) and ORC Assemblies	8-14, 8-15
Main Switch Assembly.	8-16, 8-17
O2 Supply Pressure Switch, O2 - Air Switch	8-18, 8-19
Flowmeter Shield, Gauges	8-20, 8-21
Flow Tubes, Restrictor Assemblies, Flow Control Valves	8-22, 8-23
Auxiliary O2 Flowmeter.	8-24, 8-25
Cylinder Yokes, Regulators, O2 Flush Valve	8-26, 8-27
Frame, Casters	8-28, 8-29
Power Supply Assembly and related items.	8-30, 8-31
Absorber, Inspiratory Valve, Ultrasonic Flow Sensor	8-32, 8-33
Valve, Man/Auto Selector.	8-34, 8-35
Breathing Pressure Hose Assembly, O2 Sensor, PEEP Valve, Exp. valve . .	8-36, 8-37
Scavenger, Touch-up paint: Euro white, Euro Blue	8-38, 8-39
Miscellaneous Items.	8-40
 Section 9: Specifications	
General	9-1

Contents continues on next page with a List of Figures.

DESCRIPTION	PAGE
-------------	------

List of Figures


Narkomed Mobile Anesthesia System.	1-0
2-1 Power-up Diagnostic Screen	2-1
2-2 Main Service Screen, View Mode	2-2
2-3 Main Service Screen, Service Mode.	2-3
2-4 Main Service Screen, ID Entry	2-4
2-5 Service Log Screen	2-5
2-6 PMS criteria Screen	2-6
2-7 Oxygen Monitor Service Screen.	2-7
2-8 Pressure Monitor Service Screen.	2-8
2-9 Serial Port Configuration Screen	2-9
3-1 Power Supply Voltage Test Points	3-2
3-2 Narkomed Mobile Power Distribution	3-2
4-1 Cylinder Yokes and Regulators	4-3
4-2 Cylinder and Pipeline Pressure Gauges	4-6
4-2A Tubing Connections: Flowmeter Sub-assembly	4-7
4-3 Flow Tube Replacement.	4-10
4-4 Auxiliary O2 Flowmeter.	4-13
4-5 Flow Control Valves	4-16
4-6 OFPD Replacement	4-19
4-7 Oxygen Supply Pressure Alarm Switch	4-22
4-8 O2 - Air Switch	4-24
4-9 System Power Switch.	4-26
4-10 Oxygen Ratio Controller	4-29
4-11 O2 Flush Valve.	4-32
4-12 Ventilator Controller (Bezel Assembly)	4-36
4-13 Bellows Valve Assembly.	4-38
4-14 Caster Replacement	4-40
4-15 Battery Replacement	4-42
4-16 Power Supply Assembly Replacement	4-44
4-16A Removal of Machine from Support Frame	4-45
4-17 Processor Assembly	4-48
4-18 Display Assembly	4-50
4-19 Breathing Pressure Hose Assembly	4-52
4-20 Manual/Automatic Selector Valve.	4-54
4-21 Ultrasonic Flow Sensor	4-56
4-22 Oxygen Sensor	4-58

DESCRIPTION	PAGE
--------------------	-------------

List of Figures (continued)

5-1 Cylinder Pressure Regulator Adjustment	5-3
5-2 Oxygen Supply Pressure Alarm Switch Adjustment	5-5
5-3 Oxygen Ratio Controller (ORC) Adjustment	5-7
5-4 Oxygen Sensor Calibration	5-9
5-5 Pressure Monitor Calibration	5-11
7-1 Narkomed Mobile Serial Port Location	7-2
Vaporizer & rel. parts, fresh gas outlet, mainswitch, flowmtr sub-asm	8-2
Processor Assembly and related parts	8-4
Ventilator Controller (Bezel Assembly).	8-6
Bellows Valve Assembly.	8-8
Display Assembly and Monitor Support Arm	8-10
Pipeline Inlet (Manifold) Assemblies.	8-12
Failsafe (OFPD) and ORC Assemblies	8-14
Main Switch Assembly	8-16
O ₂ Supply Pressure Switch, O ₂ - Air Switch.	8-18
Flowmeter Shield, Gauges	8-20
Flow Tubes, Restrictor Assemblies, Flow Control Valves	8-22
Auxiliary O ₂ Flowmeter	8-24
Cyl. Yokes and Regulators, O ₂ Flush Valve	8-26
Casters,	8-28
Power Supply Assembly and related items	8-30
Absorber, Inspiratory Valve, Ultrasonic Flow Sensor	8-32
Valve, Man/Auto Selector.	8-34
Breathing Pressure Hose Assembly, O ₂ Sensor, PEEP Valve, Expiratory valve .	8-36
Scavenger	8-38

Narkomed Mobile Anesthesia System



General Safety Precautions

The following are general safety precautions that apply during servicing of the equipment covered by this manual. These precautions are repeated elsewhere in this manual where needed.

WARNINGS indicate conditions or practices which if not strictly observed could result in personal injury.

CAUTIONS indicate conditions or practices which if not strictly observed or remedied could result in damage to the equipment.

WARNING: Ensure that AC power is removed from the machine before removing the power supply. Failure to observe this precaution may cause injury by electric shock.

WARNING: Possible explosive hazard if used in the presence of flammable anesthetics.

CAUTION: The flow tube must be properly centered over the guide rings or damage to the flow tube may occur.

CAUTION: Do not over-tighten the retainer. Over-tightening the retainer may break the flowmeter tube.

CAUTION: The controller circuit board contains static sensitive devices. Use ESD protection when handling the controller assembly. Static discharge can damage components on the circuit board.

CAUTION: The processor board contains static sensitive devices. Use ESD protection when handling the processor assembly. Static discharge can damage components on the circuit board.

CAUTION: Observe wiring colors and polarity markings to ensure that the battery is connected correctly. Connecting the battery with reversed polarity may damage the circuitry.

CAUTION: Always operate machine on level surface. Before moving the machine, remove vaporizer, remove items from top shelf and display arm, secure absorber against left side of machine.

1.0 Introduction

1.1 Purpose

This manual provides the information needed to field service and maintain the Narkomed Mobile anesthesia system. The DIAGNOSTICS section describes self-test and service diagnostics for checking the system functions. An understanding of the on-board service capabilities is necessary before any attempt is made to troubleshoot the unit. The TROUBLESHOOTING section shows the electrical distribution scheme and provides troubleshooting guides to assist the TSR in locating the source of a problem. The REPLACEMENT PROCEDURES section contains instructions for removal and replacement of the assemblies that are considered field-replaceable. The ADJUSTMENT AND CALIBRATION PROCEDURES section contains the field procedures needed to restore original system specifications. The Periodic Manufacturer's Certification (PMC) PROCEDURE section outlines the steps required to verify the electrical, mechanical and pneumatic safety of the unit and also identifies components requiring periodic replacement.

1.2 Recommendations

Because of the sophisticated nature of Draeger Medical, Inc. anesthesia equipment and its critical importance in the operating room setting, it is highly recommended that only appropriately trained and experienced professionals be permitted to service and maintain this equipment. Please contact DrägerServiceSM at (800) 543-5047 for service of this equipment.

Draeger Medical, Inc. also recommends that its anesthesia equipment be serviced at three-month intervals. Periodic Manufacturer's Service Agreements are available for equipment manufactured by Draeger Medical, Inc.. For further information concerning these agreements, please contact us at (800) 543-5047.

Draeger Medical, Inc. products/material in need of factory repair shall be sent to:

**DrägerService
3124 Commerce Drive
Telford, PA 18969
(Include RMA Number)**

1.3 General Troubleshooting Guidelines

Troubleshooting the Narkomed Mobile should always begin by communicating with those who observed or experienced a problem with the unit. This may eliminate unnecessary troubleshooting steps. Once a general problem is identified, refer to the troubleshooting flow charts in Section 3 to determine the proper corrective action to be taken.

After a component has been replaced, verify that the unit is operating properly by running the appropriate diagnostic procedure. The PMC PROCEDURE in Section 6 must also be performed after any component has been replaced.

1.4 Related Publications

Narkomed Mobile Operator's Manual, Part Number 4115139-001

1.5 Symbol Definitions



CAUTION: Refer to accompanying documents before operating equipment.



CAUTION: Risk of electric shock. Do not remove cover. Refer servicing to a qualified technical service representative.



Degree of protection against electric shock: Type B.

[RETURN TO THIS MANUAL'S TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

2.0 Diagnostics

The Narkomed Mobile contains a diagnostic system that monitors certain system functions and records their operational status. Following a brief System Startup display at power up, the diagnostics screen shown in Figure 2-1 appears. This display includes one of three messages at the completion of the diagnostics

FUNCTIONAL: This message indicates that the Narkomed Mobile has passed all power-up tests and is fully functional. The machine will proceed to the MACHINE MONITOR screen after a short delay.

CONDITIONALLY FUNCTIONAL: This message indicates that a minor problem has been detected. The screen will retain this display until any key is pressed, then the MACHINE MONITOR screen will be displayed.

NON-FUNCTIONAL: this message indicates that a serious problem has been detected. The machine will not proceed into the MACHINE MONITOR or SYSTEM MONITOR screen.

The PREVENTIVE MAINTENANCE DUE message will appear on the screen if the current date exceeds the Periodic Manufacturer's Service due date stored in the machine.

Further diagnostic functions are available through service screens that can be called up at the display panel. The following paragraphs provide a description of each service screen that can be accessed at the display. If no display is present upon system power-up, refer to Section 3 of this manual for troubleshooting assistance.


COPYRIGHT 1998		NAD INC.
SOFTWARE ID		XXXX
DIAGNOSTIC TESTS		
FIRMWARE	PASS	
RAM	PASS	
VIDEO	PASS	
A/D CONVERTER	PASS	
AUDIO -PRIMARY	PASS	
-BACKUP	PASS	
SERIAL I/O	PASS	
CLOCK	PASS	
NON-VOLATILE MEMORY	PASS	
PREVENTIVE MAINTENANCE DUE		
FUNCTIONAL		


Figure 2-1. Power-Up Diagnostics Screen

2.1 Main Service Screen

2.1.1 View Mode

The Main Service Screen displays the machine serial number, the last service date, hours run since last service and total hours run.

To access the Main Service Screen, press and hold the Oxygen High Limit and Volume Low Limit keys, and press the  key. The View Mode service screen shown in Figure 2-2 will then appear.

Press the  key to proceed to the Service Mode as shown in Figure 2-3, or press the key next to EXIT to return to the monitoring screen.


MAIN SERVICE SCREEN			
SRVC LOG	MACHINE SERIAL NUMBER	:	XXXXXXXX
	LAST SERVICE DATE	:	8-30-1998
	HOURS RUN SINCE LAST SERVICE:		263
	TOTAL HOURS RUN	:	624
PMS SCHED			
MON CAL			SELECT
PORT	SERVICE CODE		EXIT

Figure 2-2. Main Service Screen, View Mode

2.1.2 Service Mode



In this screen, the Service Code changes to SRVC.

Press the key next to **SELECT** to enable the Technical Service ID entry as described on the next page.

MAIN SERVICE SCREEN	
SRVC LOG	MACHINE SERIAL NUMBER : XXXXXXXX
	LAST SERVICE DATE : 8-30-1998
	HOURS RUN SINCE LAST SERVICE: 263
	TOTAL HOURS RUN : 624
PMS SCHED	
MON CAL	SELECT
PORT	EXIT
SERVICE CODE	
	SRVC

Figure 2-3. Main Service Screen, Service Mode

2.1.3 Service Mode: I.D. Entry

The Service Mode screen appears as shown in Figure 2-4. Press the key next to SELECT. Enter the first digit of your service code by using the  and  keys to display the desired character. Press the key next to SELECT to advance to the next digit, and enter the next and remaining I.D. characters in the same manner.

When this screen is entered, an entry is made in the Service Log.

To access any of the other service screens described on the following pages, press the key next to the desired function on the left side of the screen: Service Log, PMS Schedule, Monitor Calibration, or Port communication settings.

Pressing the key next to RESET will reset the HOURS RUN SINCE LAST SERVICE to zero, and the LAST SERVICE DATE to the current date.



If desired, press the key next to EXIT to return to the monitoring screen.

MAIN SERVICE SCREEN	
RESET	
	MACHINE SERIAL NUMBER : XXXXXXXX
	LAST SERVICE DATE : 8-30-1998
SRVC	HOURS RUN SINCE LAST SERVICE: 263
LOG	TOTAL HOURS RUN : 624
PMS SCHED	
MON CAL	<div style="border: 1px solid black; padding: 5px; display: inline-block;"> TECHNICAL SERVICE REPRESENTATIVE I.D. </div>
	SELECT
PORT	<div style="display: flex; justify-content: space-around; width: 100px;"> <div style="border: 1px solid black; padding: 2px 5px;">0</div> <div style="border: 1px solid black; padding: 2px 5px;">0</div> <div style="border: 1px solid black; padding: 2px 5px;">0</div> <div style="border: 1px solid black; padding: 2px 5px;">0</div> </div>
	EXIT

Figure 2-4. Main Service Screen, ID Entry

2.2 Service Log

From the Service Screen (described earlier), press the key next to SRVC LOG.

Figure 2-5 shows an example of the screen that will appear. This screen allows you to view the events recorded in the machine's service log. Use the  and  keys to scroll down or up through the log entries.

Press the key next to EXIT to return to the Main Service Screen.



SERVICE LOG			
DATE	TIME	PARAMETER	CODE
09-11-98	10:26	00000000	0000
SYSTEM POWERUP			
09-11-98	10:30	00000000	E400
AUDIOGEN SPKR CHK			
09-13-98	07:30	00000004	E100
EXIT			

Figure 2-5. Service Log Screen


2.3 PMS Criteria Screen

The PMS Criteria Screen allows you to select the month when the PREVENTIVE MAINTENANCE DUE message appears on the power-up diagnostics screen.

From the Service Screen (described earlier), press the hidden key next to PMS SCHED.

Figure 2-6 shows an example of the screen that will appear. Use the  and  keys to set the desired month.

Press the key next to EXIT to return to the Main Service Screen.



The screenshot shows a screen titled "SELECT CRITERIA FOR 'PREVENTIVE MAINTENANCE' MESSAGE". Below the title, there are two input fields: "MONTH" and "SEP", followed by the year "1998". At the bottom right of the screen, there is an "EXIT" button.

Figure 2-6. PMS Criteria Screen

2.4 Oxygen Monitor Service Screen

The Oxygen Monitor Service Screen shown in Figure 2-7 displays current readings for the O₂ cells, a zero calibration procedure, and the stored calibration values.

From the Service Screen (described earlier), press the key next to MON CAL.

To perform a zero calibration, follow the calibration procedure shown on the screen. Pressing the key next to ZERO stores the current values as the new zero calibration.

To proceed to the Pressure Monitor Service Screen, press the key next to PRES MON. To return to the Main Service Screen, press the key next to EXIT.

OXYGEN MONITOR SERVICE SCREEN	
CURRENT CELL A: 238	ZERO
CURRENT CELL B: 250	
ZERO CALIBRATION PROCEDURE:	
- REMOVE O2 CELL FROM HOUSING	
- LET CURRENT CELL VALUES STABILIZE	
- PRESS "ZERO" KEY TO ENTER	
CALIBRATION VALUES	
- REINSTALL O2 CELL IN SENSOR HOUSING	
	PRES MON
STORED ZERO CELL A: 250	
STORED ZERO CELL B: 250	EXIT

Figure 2-7. Oxygen Monitor Service Screen

2.5 Pressure Monitor Service Screen

The Pressure Monitor Service Screen shown in Figure 2-8 displays the current reading for airway pressure, a procedure for zero and span calibration, and the stored calibration values.

To enter the Pressure Monitor Service Screen from the Oxygen Monitor Service Screen (described earlier), press the key next to PRES MON (ref. Figure 2-7).

To perform a zero calibration, follow the procedure shown on the screen. Pressing the key next to ZERO stores the current value as the new zero calibration.

To perform a span calibration, follow the procedure shown on the screen. Pressing the key next to SPAN stores the current value as the new span calibration.

To return to the Oxygen Monitor Service Screen, press the key next to OXY MON. To return to the Main Service Screen, press the key next to EXIT.



PRESSURE MONITOR SERVICE SCREEN	
CURRENT PRESSURE VALUE: 250	ZERO
ZERO CALIBRATION PROCEDURE:	
- REMOVE PRESSURE SAMPLE LINE FROM ABSORBER, EXPOSE TO AIR.	SPAN
- LET CURRENT PRESSURE VALUE STABILIZE	
- SELECT "ZERO" KEY TO ENTER CALIBRATION VALUES.	
SPAN CALIBRATION PROCEDURE:	
- REMOVE PRESSURE SAMPLE LINE FROM ABSORBER, APPLY 50 CMH ₂ O CONSTANT PRESSURE AT THE SAMPLE LINE, VERIFIED BY A KNOWN, CALIBRATED METER.	OXY MON
- LET PRESSURE VALUE STABILIZE	EXIT
- SELECT THE "SPAN" KEY TO ENTER THE CURRENT VALUE.	

Figure 2-8. Pressure Monitor Service Screen


2.6 Serial Port Configuration Screen

The Serial Port Configuration screen shown in Figure 2-9 allows you to set the machine parameters for communicating with external devices.

From the Service Screen (described earlier), press the key next to PORT.

Use the  and  keys to change the settings; press the key next to SELECT to move to the next setting.

Press the key next to EXIT to return to the Main Service Screen.



SERIAL PORT CONFIGURATION	
BAUD RATE:	<input type="text" value="19200"/>
PARITY :	NONE
STOP BITS:	1
DATA BITS :	8
PROTOCOL:	VLNK
SELECT	
EXIT	

Figure 2-9. Serial Port Configuration Screen

[RETURN TO THIS MANUAL'S TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

3.0 Troubleshooting

This section contains information to assist the Draeger Medical, Inc. qualified Technical Service Representative (TSR) in locating electrical faults affecting the Narkomed Mobile monitoring and display devices. Since most troubleshooting efforts begin with verifying power supply voltages, the following paragraph outlines the voltage distribution scheme within the machine along with test points for each of the voltages.

3.1 Power Supply and Voltage Distribution


In the Narkomed Mobile, +5VDC, +12VDC and -12VDC are supplied to J14 on the processor board; +8VDC is supplied to J2 on the ventilator controller. These voltages can be measured at the connectors shown in Figure 3-1. Output voltage of the primary power supply is measured at J3 on the Condor supply. Table 3-1 lists the acceptable range for each voltage under normal load conditions. Figure 3-2 shows a block diagram of the Narkomed Mobile voltage distribution scheme.

Table 3-1. Test Points and Allowable Ranges

PROCESSOR	VOLTAGE	ACCEPTABLE RANGE
J14-12,14 (Red, Orn)	+ 5 VDC	4.80 to 5.25 VDC
J14-1 (Wht)	+ 12 VDC	11.65 to 12.85 VDC
J14-5 (Gry)	- 12 VDC	-11.50 to -13.00 VDC
J14-7,8,9 (Grn, Blu, Yel)	Common	
VENTILATOR CONTROLLER	VOLTAGE	ACCEPTABLE RANGE
J2-3 (Brn)	+ 8 VDC	7.70 to 8.30 VDC
J2-1 (Yel)	Common	
CONDOR PWR SUPP	VOLTAGE	ACCEPTABLE RANGE
J3-1 (Brn)	+ 15 VDC	14.0 to 16.0 VDC
J3-8 (Wht)	Common	

TROUBLESHOOTING GUIDE (continued)

NM MOBILE



SV00301

Figure 3-1. Power Supply Voltage Test Points




Figure 3-2. Narkomed Mobile Power Distribution

NOTE: The Narkomed Mobile will not turn on, or operate, unless the power cable is connected to J14 on the processor board. Disconnecting this cable breaks a sense connection that automatically powers down +5V, +12V, -12V, and +8V.

3.2 Battery

While the machine is operating from an AC line, the battery voltage at full charge should be within the range of 13.50 to 14.80 VDC. Battery voltage can be measured at the battery terminals. During battery operation, the low battery cutoff voltage should be within the range of 10.5 to 10.0 VDC.


3.3 Troubleshooting Guides

Table 3-2 lists common failure modes and symptoms (excluding simultaneous multiple faults) for the monitoring and display devices in the Narkomed Mobile. Each failure mode or symptom is keyed to a troubleshooting guide flow chart at the back of this section to assist the TSR in locating a problem. These flow charts assume that the machine is plugged into an AC outlet with the correct voltage, and the machine is not running on its backup battery.


Table 3-2. Narkomed Mobile Failure Mode and Symptom List

FAILURE MODE / SYMPTOM	CORRECTIVE ACTION
Loss of O ₂ Monitor	Guide 1
Loss of Breathing Pressure Monitor	Guide 2
Loss of Respiratory Volume Monitor	Guide 3
No Audio Alarms	Guide 4
Serial Port Communication Failure	Guide 5
No Oxygen Supply Pressure Alarms	Guide 6
Display Blank Upon System Power-up	Guide 7
Keypad Inoperative	Guide 8
Ventilator Inoperative	Guide 9


GUIDE 1: Loss of O₂ Monitor




GUIDE 2: Loss of Breathing Pressure Monitor




GUIDE 3: Loss of Respiratory Volume Monitor




GUIDE 4: No Audio Alarms




GUIDE 5: Serial Port Communication Failure




GUIDE 6: No O₂ Supply Pressure Alarms




GUIDE 7: Display Blank Upon System Power-up



GUIDE 8: Keypad Inoperative



GUIDE 9: Ventilator Inoperative



4.0 Replacement Procedures

This section outlines removal and replacement procedures for the field-replaceable assemblies of the Narkomed Anesthesia System.

These procedures are to be performed only by a Draeger Medical, Inc. qualified Technical Service Representative (TSR).

The following are the only procedures authorized by Draeger Medical, Inc. to be performed in the field. All other service procedures shall be referred to Draeger Medical, Inc.'s Technical Service Department.

NOTE: The PMC PROCEDURE given in Section 6 must be performed after any replacement, removal, calibration or adjustment procedure.

4.1 Cylinder Yokes and Regulators

The cylinder yokes and regulators are installed as an assembly consisting of the yoke, check valve, regulator and spacer block. Access to the tubing connections requires removal of the flowmeter housing back cover. Tubing and mounting arrangements are shown in Figure 4-1.

- 4.1.1 Disconnect all pipeline hoses and turn the System Power switch to ON.
- 4.1.2 Close the N₂O cylinder valve; open the oxygen cylinder valve.
- 4.1.3 Set the oxygen flow to 5 liters per min.
- 4.1.4 Open the other gas flow control valves to drain pressure from the system.
- 4.1.5 Close the O₂ cylinder valve, and close the flow control valves. Press the O₂ Flush valve to drain oxygen pressure from the system.
- 4.1.6 Turn the System Power switch to STANDBY.
- 4.1.7 Remove the cylinders from the yokes.
- 4.1.8 Remove the flowmeter housing back cover.
- 4.1.9 Disconnect the tubing from the regulators where connections are accessible.

For regulator connections that are not accessible, disconnect these tubes at their other end.
- 4.1.10 Remove the yoke spacer mounting screws, and lift the assembly from the flowmeter housing.
- 4.1.11 If you are replacing a regulator, record the serial number of the regulator that was removed, and record the serial number of the replacement regulator.
- 4.1.12 Where tubing was removed, transfer the tubing to the corresponding connections on the replacement regulator.




Figure 4-1. Cylinder Yokes and Regulators

REPLACEMENT PROCEDURES (continued)

NM MOBILE

- 4.1.13 Position the replacement yoke and regulator assembly in the flowmeter housing, and secure it with the hardware that was previously removed.
- 4.1.14 Reconnect all tubing that was previously disconnected within the flowmeter housing.
- 4.1.15 Reinstall the cylinders in the yokes.
- 4.1.16 Measure (and adjust if necessary) the regulator output pressure in accordance with the procedure given in Section 5.
- 4.1.17 Reinstall the flowmeter housing back cover.
- 4.1.18 Perform the PMC Procedure given in Section 6.

4.2 Cylinder and Pipeline Pressure Gauges

Replacement of the cylinder and pipeline pressure gauges requires disassembly in the flowmeter sub-assembly area for access to the gauge connections. Figure 4-2 shows gauge mounting and connection details.

- 4.2.1 Disconnect all pipeline hoses and turn the System Power switch to ON.
- 4.2.2 Open the oxygen cylinder valve.
- 4.2.3 Set the oxygen flow to 5 liters per min.
- 4.2.4 Open the other gas flow control valves to drain pressure from the system.
- 4.2.5 Close the O₂ cylinder valve, and close the flow control valves. Press the O₂ Flush valve to drain oxygen pressure from the system.
- 4.2.6 Turn the System Power switch to STANDBY.
- 4.2.7 Remove the oxygen flow control knob.
- 4.2.8 Remove the two screws holding the knob guard in place, and remove the knob guard.
- 4.2.9 Remove the plexiglass flowmeter shield.
- 4.2.10 Remove the flowmeter housing back cover. Be sure to disconnect the ventilator exhaust hose.
- 4.2.11 Disconnect the copper tubing at points **A**, **B**, **C** and **D** as shown in Figure 4-2A.
- 4.2.12 Remove the four screws securing the flowmeter sub-assembly to the flowmeter housing.
- 4.2.13 Pull the flowmeter sub-assembly forward far enough to gain access to the gauge connections.
- 4.2.14 For the cylinder pressure gauges:

Disconnect the 3/16 in. copper tube compression fitting at the back of the gauge.

For the pipeline pressure gauges:

Remove the press-on hose clamp and disconnect the flex tubing from the hose barb at the back of the gauge.





Figure 4-2. Cylinder and Pipeline Pressure Gauges



SV00348

Figure 4-2A. Tubing Connections: Flowmeter Sub-assembly

REPLACEMENT PROCEDURES (continued)

NM MOBILE

- 4.2.15 Remove the gauge mounting nuts, and remove the gauge from the front of the flowmeter housing.
- 4.2.16 Install the replacement gauge in the flowmeter housing and secure it with the hardware that was previously removed.
- 4.2.17 For the cylinder pressure gauges:

Reconnect the 3/16 in. copper tube compression fitting at the back of the gauge.

For the pipeline pressure gauges:

Reconnect the flex tubing from to the hose barb at the back of the gauge, and secure it with the press-on hose clamp.
- 4.2.18 Reinstall the flowmeter sub-assembly, and reconnect all copper tubing.
- 4.2.19 Reinstall the front plexiglass flowmeter shield.
- 4.2.20 Reinstall the knob guard and secure it with the two mounting screws.
- 4.2.21 Reinstall the oxygen flow control knob. Re-set the 'off stop' with the label oriented correctly.
- 4.2.22 Reconnect the ventilator exhaust hose and reinstall the flowmeter housing back cover.
- 4.2.23 Connect the pipeline supplies.
- 4.2.24 Perform the PMC Procedure given in Section 6.

4.3 Flowmeters

The flowmeter tubes are held by compression in gaskets at the top and bottom of each tube. Each upper gasket is seated in an adjustable retainer that allows removal of the tube as shown in Figure 4-3. Access to the flow tubes and their retainers requires removal of the plexiglass flowmeter shield.

- 4.3.1 Disconnect all pipeline hoses and turn the System Power switch to ON.
- 4.3.2 Open the oxygen cylinder valve.
- 4.3.3 Set the oxygen flow to 5 liters per min.
- 4.3.4 Open the other gas flow control valves to drain pressure from the system.
- 4.3.5 Close the O₂ cylinder valve and the O₂ flow control valve. Press the O₂ Flush valve to drain oxygen pressure from the system.
- 4.3.6 Turn the System Power switch to STANDBY.
- 4.3.7 Remove the oxygen flow control knob.
- 4.3.8 Remove the two screws holding the knob guard in place, and remove the knob guard.
- 4.3.9 Remove the plexiglass flowmeter shield.
- 4.3.10 Turn the flow tube retainer as shown in the illustration until you can pull the top of the flow tube outward, and remove the tube.

NOTE: If the bottom of the tube is seated in a restrictor housing, be sure that the arrangement of the restrictor and its gaskets is not disturbed.

- 4.3.11 Make sure that the replacement flow tube bears the correct markings and has a ball.

CAUTION: The flow tube must be properly centered over the guide rings or damage to the flow tube may occur.

- 4.3.12 Place the bottom of the flow tube into the guide ring of the lower gasket, and position the top of the flow tube in the center of the retainer.

CAUTION: Do not over-tighten the retainer. Over-tightening the retainer may break the flowmeter tube.

REPLACEMENT PROCEDURES (continued)

NM MOBILE




Figure 4-3. Flow Tube Replacement

- 4.3.13 Ensure that the markings on the flow tube are facing forward, and turn the retainer as shown in the illustration until the flow tube is firmly held in place.
- 4.3.14 Perform the following leak test on the system:
- Disconnect the absorber hose from the fresh gas outlet. Ensure that all flow control valves are closed.
- Connect a test gauge and B.P. bulb to the fresh gas outlet, and pressurize the system to 50 cm H₂O.
- The pressure should not drop more than 10 cm H₂O in thirty seconds.
- 4.3.15 Disconnect the test gauge and re-connect the absorber hose to the fresh gas outlet.
- 4.3.16 Reinstall the front plexiglass flowmeter shield.
- 4.3.17 Reinstall the knob guard and secure it with the two mounting screws.
- 4.3.18 Reinstall the oxygen flow control knob.
- 4.3.19 Connect the pipeline hoses.
- 4.3.20 Perform the PMC Procedure given in Section 6.

4.4 Auxiliary O₂ Flowmeter

The auxiliary O₂ flowmeter is attached to the side of the flowmeter housing. Access to its attaching hardware requires disassembly in the flowmeter sub-assembly area. Figure 4-4 shows the arrangement of the mounting screws and O₂ supply line.

- 4.4.1 Disconnect all pipeline hoses and turn the System Power switch to ON.
- 4.4.2 Open the oxygen cylinder valve.
- 4.4.3 Set the oxygen flow to 5 liters per min.
- 4.4.4 Open the other gas flow control valves to drain pressure from the system.
- 4.4.5 Close the O₂ cylinder valve and the O₂ flow control valve. Press the O₂ Flush valve to drain oxygen pressure from the system.
- 4.4.6 Turn the System Power switch to STANDBY.
- 4.4.7 Remove the flowmeter housing back cover. Be sure to disconnect the ventilator exhaust hose.
- 4.4.8 Ref. Figure 4-2A: Disconnect the copper tubing at points **A**, **B**, **C** and **D**.
- 4.4.9 Remove the four screws securing the flowmeter sub-assembly to the flowmeter housing.
- 4.4.10 Pull the flowmeter sub-assembly forward far enough to gain access to the auxiliary O₂ flowmeter mounting screws.
- 4.4.11 Remove the press-on hose clamp and disconnect the flexible auxiliary O₂ supply line from the hose barb on the O₂ pipeline inlet assembly.
- 4.4.12 Remove the screws securing the auxiliary O₂ flowmeter to the side of the flowmeter housing, and remove the flowmeter.
- 4.4.13 Position the replacement flowmeter at the side of the flowmeter housing (feed the flex tubing through the clearance hole) and secure the flowmeter with the two screws that were previously removed.




Figure 4-4. Auxiliary O2 Flowmeter

REPLACEMENT PROCEDURES (continued)


NM MOBILE

- 4.4.14 Reconnect the small diameter tubing from the auxiliary O₂ flowmeter to the hose barb on the O₂ pipeline inlet assembly and secure the connection with the press-on hose clamp.
- 4.4.15 Reinstall the flowmeter sub-assembly, and reconnect all copper tubing.
- 4.4.16 Reinstall the front plexiglass flowmeter shield.
- 4.4.17 Reinstall the knob guard and secure it with the two mounting screws.
- 4.4.18 Reinstall the oxygen flow control knob. Re-set the 'off stop' with the label oriented correctly.
- 4.4.19 Reconnect the ventilator exhaust hose and reinstall the flowmeter housing back cover.
- 4.4.20 Connect the pipeline supplies.
- 4.4.21 Perform the PMC Procedure given in Section 6.

4.5 Flow Control Valves

The flow control valves have replaceable elements that are removable from the front of the flowmeter sub-assembly as shown in Figure 4-5. Each flow control knob has a positive stop arrangement that prevents damage to the valve seat. Whenever a valve is replaced the "off stop" must be set as outlined in the following procedure.

- 4.5.1 Disconnect all pipeline hoses and turn the System Power switch to ON.
- 4.5.2 Open the oxygen cylinder valve.
- 4.5.3 Set the oxygen flow to 5 liters per min.
- 4.5.4 Open the other gas flow control valves to drain pressure from the system.
- 4.5.5 Close the O₂ cylinder valve and the O₂ flow control valve. Press the O₂ Flush valve to drain oxygen pressure from the system.
- 4.5.6 Turn the System Power switch to STANDBY.
- 4.5.7 Remove the oxygen flow control knob.
- 4.5.8 Remove the two screws holding the knob guard in place, and remove the knob guard.
- 4.5.9 Remove the plexiglass flowmeter shield.
- 4.5.10 Remove the knob (if not already removed) from the valve that is being replaced, and remove the stop pin nut.
- 4.5.11 Remove the flow control valve by holding it at the wrench flats and turning it counter-clockwise.
- 4.5.12 Install the replacement flow control valve in the flowmeter sub-assembly.
- 4.5.13 Reinstall the stop pin nut.
- 4.5.14 Reinstall the front plexiglass flowmeter shield.
- 4.5.15 Reinstall the knob guard and secure it with the two mounting screws.
- 4.5.16 Connect the pipeline supplies and turn the System Power switch to ON.



SV00306

Figure 4-5. Flow Control Valves

4.5.17 For the O₂ flow control valve:

Turn the flow control valve clockwise until the flow rate will not drop any further. (If the machine has been modified to eliminate the minimum flow feature, turn the valve until the flow rate is zero.)

For the other gas flow control valves:

Set the oxygen flow rate to four liters per minute.

Turn the other gas flow control valve clockwise until the flow rate is zero.

4.5.18 Place the knob on the flow control valve shaft and turn it clockwise until it engages the stop pin. Tighten one of the knob setscrews.

4.5.19 Turn the knob in both directions and ensure that the flow can be controlled over its entire range. When the valve is closed, the knob should be against the clockwise stop. Tighten the remaining set screw.

4.5.20 If the knob label is not horizontal when the valve is closed, remove the label and install a new label in the correct position.

4.5.21 Perform the PMC Procedure given in Section 6.

4.6 Oxygen Supply Failure Protection Device

The air and nitrous oxide supplies within the machine are monitored by oxygen supply failure protection devices (OFPDs) which prevent the flow of these gases if there is insufficient oxygen pressure available. Access to these devices requires disassembly in the flowmeter sub-assembly area. See Figure 4-6.

- 4.6.1 Disconnect all pipeline hoses and turn the System Power switch to ON.
- 4.6.2 Open the oxygen cylinder valve.
- 4.6.3 Set the oxygen flow to 5 liters per min.
- 4.6.4 Open the other gas flow control valves to drain pressure from the system.
- 4.6.5 Close the O₂ cylinder valve and the O₂ flow control valve. Press the O₂ Flush valve to drain oxygen pressure from the system.
- 4.6.6 Turn the System Power switch to STANDBY.
- 4.6.7 Remove the oxygen flow control knob.
- 4.6.8 Remove the two screws holding the knob guard in place, and remove the knob guard.
- 4.6.9 Remove the plexiglass flowmeter shield.
- 4.6.10 Remove the flowmeter housing back cover. Be sure to disconnect the ventilator exhaust hose.
- 4.6.11 Ref. Figure 4-2A: Disconnect the copper tubing at points **A**, **B**, **C** and **D**.
- 4.6.12 Remove the four screws securing the flowmeter sub-assembly to the flowmeter housing.
- 4.6.13 Pull the flowmeter sub-assembly forward far enough to gain access to the OFPDs.
- 4.6.14 Remove the press-on hose clamp and disconnect the small diameter tubing from the hose barb at the top of the OFPD.
- 4.6.15 Air OFPD: Remove the three screws securing the OFPD to the flowmeter sub-assembly, and remove the OFPD.

NOTE: The N₂O OFPD has longer mounting screws, which pass through the oxygen ratio controller (ORC) and into the flowmeter sub-assembly. These screws retain both devices.




Figure 4-6. OFPD Replacement

REPLACEMENT PROCEDURES (continued)

NM MOBILE

N₂O OFPD: Remove the three screws securing the OFPD, and remove the OFPD.

- 4.6.16 Ensure that the O-ring is correctly in place, and install the replacement OFPD with the hardware that was previously removed.
- 4.6.17 Reconnect the small diameter tubing to the hose barb on the OFPD and secure the connection with the press-on hose clamp.
- 4.6.18 Reinstall the flowmeter sub-assembly, and reconnect all copper tubing.
- 4.6.19 Reinstall the front plexiglass flowmeter shield.
- 4.6.20 Reinstall the knob guard and secure it with the two mounting screws.
- 4.6.21 Reinstall the oxygen flow control knob. Re-set the 'off stop' with the label oriented correctly.
- 4.6.22 Reconnect the ventilator exhaust hose and reinstall the flowmeter housing back cover.
- 4.6.23 Connect the pipeline supplies.
- 4.6.24 Perform the PMC Procedure given in Section 6.

4.7 Oxygen Supply Pressure Alarm Switch

The oxygen supply low pressure alarm switch is located inside the flowmeter housing, attached to a bracket on the flowmeter sub-assembly. Access to the switch requires removal of the flowmeter housing back cover. Figure 4-7 shows the pneumatic and electrical connections to the switch.

- 4.7.1 Turn the System Power switch to STANDBY and disconnect the pipeline hoses.
- 4.7.2 Remove the flowmeter housing back cover. Be sure to disconnect the ventilator exhaust hose.
- 4.7.3 Separate the in-line connector on the switch wire harness.
- 4.7.4 Remove the press-on hose clamp and disconnect the tubing from the hose barb on the bottom of the switch.
- 4.7.5 Loosen the four screws holding the switch to the bracket; lift out the switch and screws.
- 4.7.6 Transfer the screws to the replacement switch; position the switch on the bracket, and tighten the four screws to secure the switch to the bracket.
- 4.7.7 Reconnect the tubing to the hose barb on the switch and secure it with the press-on clamp.
- 4.7.8 Join the in-line connector to its corresponding wire harness.
- 4.7.9 Perform the oxygen supply pressure alarm switch adjustment procedure given in Section 5 of this manual.
- 4.7.10 Reinstall the flowmeter housing back cover.
- 4.7.11 Perform the PMC procedure given in Section 6.

REPLACEMENT PROCEDURES (continued)

NM MOBILE




Figure 4-7. Oxygen Supply Pressure Alarm Switch

4.8 O₂ - Air Switch

The ventilator drive gas (O₂ - Air) selector is a manually operated pneumatic switch located at the opening at the side of the flowmeter housing. The switch is attached to a recess housing fixed to the inside wall of the flowmeter housing. Access to the switch requires removal of the flowmeter housing back cover. Figure 4-8 shows the mounting arrangement and connections.

- 4.8.1 Disconnect all pipeline hoses and turn the System Power switch to ON.
- 4.8.2 Open the oxygen cylinder valve.
- 4.8.3 Set the oxygen flow to 5 liters per min.
- 4.8.4 Open the other gas flow control valves to drain pressure from the system.
- 4.8.5 Close the O₂ cylinder valve and the O₂ flow control valve. Press the O₂ Flush valve to drain oxygen pressure from the system.
- 4.8.6 Turn the System Power switch to STANDBY.
- 4.8.7 Disconnect the copper tube at the upper port of the Air pipeline inlet assembly (marked **A** in the illustration).
- 4.8.8 Disconnect the copper tube (O₂) from the O₂ port at the switch.
- 4.8.9 Remove the two screws securing the switch bracket.
- 4.8.10 In the recess housing, remove the panel nut and lock washer from the switch.
- 4.8.11 Pull the switch out far enough to gain access to the output port; disconnect the flex tubing and remove the switch from the flowmeter housing.
- 4.8.12 Remove the copper tube (Air) from the switch, and transfer it to the replacement switch.
- 4.8.13 Transfer the recess housing (around the switch handle) to the replacement switch.
- 4.8.14 Place the switch bracket over the switch and position the replacement switch in the flowmeter housing. Reconnect the flex tubing to the output port on the switch.
- 4.8.15 Reinstall the bracket screws that were previously removed.
- 4.8.16 Connect the O₂ copper tube to the O₂ port on the switch.

REPLACEMENT PROCEDURES (continued)

NM MOBILE

- 4.8.17 Connect the Air copper tube to the upper port of the Air pipeline inlet assembly. Tighten all copper tube compression fittings.
- 4.8.18 Reinstall the flowmeter housing back cover.
- 4.8.19 Perform the PMC Procedure given in Section 6.




Figure 4-8. O2 - Air Switch

4.9 System Power Switch

The system power switch assembly is located in the power supply housing. Access to the switch assembly requires that the machine be separated from its support frame, and removal of the power supply. Figure 4-9 shows the connection and mounting arrangement for the switch assembly.

- 4.9.1 Disconnect all pipeline hoses and turn the System Power switch to ON.
- 4.9.2 Open the oxygen cylinder valve.
- 4.9.3 Set the oxygen flow to 5 liters per min.
- 4.9.4 Open the other gas flow control valves to drain pressure from the system.
- 4.9.5 Close the O₂ cylinder valve and the O₂ flow control valve. Press the O₂ Flush valve to drain oxygen pressure from the system.
- 4.9.6 Turn the System Power switch to STANDBY.
- 4.9.7 Remove the following items from the machine: external monitors, cylinders, vaporizer and display.
- 4.9.8 Disconnect the fresh gas hose.
- 4.9.9 Disconnect the ventilator breathing and scavenger hoses, and the sensor interface connections.
- 4.9.10 Loosen the wingnut on the absorber arm, and remove the absorber assembly.
- 4.9.11 Ref. Figure 4-16A: separate the machine from its support frame.
- 4.9.12 Ref. Paragraph 4.16.8: remove the power supply assembly from its housing.
- 4.9.13 Carefully place the machine upside-down on a suitable surface.
- 4.9.14 Disconnect the tubing at the three compression fittings on the switch assembly (see Figure 4-9).
- 4.9.15 Locate the 2-pin in-line connector that joins the switch wires to the main harness (orange and violet wires).

Separate this in-line connector.
- 4.9.16 Release the two switch wire pins from the connector body, and retrieve the switch wires through the bulkhead bushing.




Figure 4-9. System Power Switch

- 4.9.17 Remove the two screws securing the switch assembly to the power supply housing, and remove the switch assembly from the housing.
- 4.9.18 Install the replacement switch assembly in the power supply housing and secure it with the hardware that was previously removed.
- 4.9.19 Route the switch wires through the bulkhead bushing in the same manner as the original, and install the pins in the connector body: violet wire in Position #1, orange wire in Position #2.
- 4.9.20 Join the in-line connector.
- 4.9.21 Re-connect the tubing to the compression fittings on the replacement switch assembly.
- 4.9.22 Ref. Paragraph 4.16.8: reinstall the power supply assembly.
- 4.9.23 Ref. Figure 4-16A: reattach the machine to its support frame.
- 4.9.24 Reinstall the absorber assembly and reconnect the fresh gas hose.
- 4.9.25 Reinstall the display, vaporizer and cylinders.
- 4.9.26 Reinstall all accessories that were previously removed; restore all breathing, scavenger and sensor interface connections. Refer to the Installation Instructions in the *Narkomed Mobile Setup and Installation Manual*.
- 4.9.27 Reconnect the pipeline supplies.
- 4.9.28 Restore AC power to the machine and ensure that new system power switch is working properly.
- 4.9.29 Perform the PMC Procedure given in Section 6.

4.10 Oxygen Ratio Controller

The Oxygen Ratio Controller (ORC) is part of the N₂O flowmeter sub-assembly and is located within the flowmeter housing. Access to the ORC requires disassembly in the flowmeter sub-assembly area, and removal of the N₂O OFPD assembly. Figure 4-10 shows the ORC location and mounting arrangement.

- 4.10.1 Disconnect all pipeline hoses and turn the System Power switch to ON.
- 4.10.2 Close all cylinder valves except the O₂ cylinder.
- 4.10.3 Set the oxygen flow to 5 liters per min.
- 4.10.4 Open the other gas flow control valves to drain pressure from the system.
- 4.10.5 Close the O₂ cylinder valve and the O₂ flow control valve. Press the O₂ Flush valve to drain oxygen pressure from the system.
- 4.10.6 Turn the System Power switch to STANDBY.
- 4.10.7 Remove the oxygen flow control knob.
- 4.10.8 Remove the two screws holding the knob guard in place, and remove the knob guard.
- 4.10.9 Remove the plexiglass flowmeter shield.
- 4.10.10 Remove the flowmeter housing back cover. Be sure to disconnect the ventilator exhaust hose.
- 4.10.11 Ref. Figure 4-2A: Disconnect the copper tubing at points **A**, **B**, **C** and **D**.
- 4.10.12 Remove the four screws securing the flowmeter sub-assembly to the flowmeter housing.
- 4.10.13 Pull the flowmeter sub-assembly forward far enough to gain access to the OFPDs.
- 4.10.14 Remove the press-on hose clamp and disconnect the small diameter tubing from the hose barb at the top of the OFPD.
- 4.10.15 Remove the press-on hose clamp and disconnect the small diameter tubing from the hose barb on the side of the ORC.

NOTE: The N₂O OFPD has long mounting screws, which pass through the oxygen ratio controller (ORC) and into the flowmeter sub-assembly. These screws retain both devices.




Figure 4-10. Oxygen Ratio Controller

REPLACEMENT PROCEDURES (continued)

NM MOBILE

- 4.10.16 Remove the three screws securing the OFPD and the ORC, and remove the OFPD and ORC.
- 4.10.17 Position the replacement ORC at the back of the N₂O flowmeter sub-assembly; be sure that its O-rings and filter are in place. Reinstall the OFPD, with the mounting screws going through the ORC and into the flowmeter sub-assembly. Tighten the screws.
- 4.10.18 Reconnect the small diameter tubing to the ORC and secure the connection with the press-on hose clamp.
- 4.10.19 Reconnect the small diameter tubing to the hose barb on the OFPD and secure the connection with the press-on hose clamp.
- 4.10.20 Reinstall the flowmeter sub-assembly, and reconnect all copper tubing.
- 4.10.21 Reinstall the front plexiglass flowmeter shield.
- 4.10.22 Reinstall the knob guard and secure it with the two mounting screws.
- 4.10.23 Reinstall the oxygen flow control knob. Re-set the 'off stop' with the label oriented correctly.
- 4.10.24 Reconnect the ventilator exhaust hose and reinstall the flowmeter housing back cover.
- 4.10.25 Connect the pipeline supplies.
- 4.10.26 Perform the PMC Procedure given in Section 6.

4.11 O₂ Flush Valve

The O₂ flush valve is located at the front of the machine next to the fresh gas outlet. Access to the flush valve is through the bottom of the machine, which requires removal of the power supply assembly. You will need to remove the machine from its support frame. Figure 4-11 shows the mounting arrangement of the O₂ flush valve and its tubing connections.

- 4.11.1 Disconnect all pipeline hoses and turn the System Power switch to ON.
- 4.11.2 Close the O₂ cylinder valve and the O₂ flow control valve. Press the O₂ Flush valve to drain oxygen pressure from the system.
- 4.11.3 Turn the System Power switch to STANDBY, and disconnect the AC power cord.
- 4.11.4 Remove the following items from the machine: external monitors, cylinders, vaporizer and display.
- 4.11.5 Disconnect the fresh gas hose.
- 4.11.6 Disconnect the ventilator breathing and scavenger hoses, and the sensor interface connections.
- 4.11.7 Loosen the wingnut on the absorber arm, and remove the absorber assembly.
- 4.11.8 Ref. Figure 4-16A: separate the machine from its support frame.
- 4.11.9 Ref. Paragraph 4.16.8: remove the power supply assembly from its housing.
- 4.11.10 Carefully place the machine upside-down on a suitable surface.
- 4.11.11 Hold the O₂ Flush button in and rotate it until one of its set screws are visible through an access hole in the guard ring, and loosen the set screw.
- 4.11.12 Turn the O₂ Flush button 180 degrees, hold it in and loosen the other set screw.
- 4.11.13 Remove the O₂ Flush button and washer from the valve shaft.
- 4.11.14 Disconnect the compression fittings at the valve. The O₂ Flush valve is retained by the guard ring on the front of the machine. Hold the body of the Clippard valve with an open end wrench; insert a rod or hex wrench through the holes in the guard ring (or use a spanner wrench), and un-screw the guard ring.




Figure 4-11. O₂ Flush Valve

- 4.11.15 Assemble the replacement O₂ Flush valve, spacer, internal tooth lock washer and guard ring through the chassis and tighten the assembly, making sure that the valve is mounted straight. Connect the compression fittings to the valve.
- 4.11.16 Place the washer and the O₂ Flush button on the valve shaft.
- 4.11.17 Hold the O₂ Flush button in and turn it until a set screw is visible through an access hole in the guard ring. Tighten the set screw. Rotate the button 180 degrees until the other set screw is visible, and tighten the set screw.
- 4.11.18 Ref. Paragraph 4.16.8: reinstall the power supply assembly.
- 4.11.19 Ref. Figure 4-16A: reattach the machine to its support frame.
- 4.11.20 Reinstall the absorber assembly and reconnect the fresh gas hose.
- 4.11.21 Reinstall the display, vaporizer and cylinders.
- 4.11.22 Reinstall all accessories that were previously removed; restore all breathing, scavenger and sensor interface connections. Refer to the Installation Instructions in the *Narkomed Mobile Setup and Installation Manual*.
- 4.11.23 Connect a test gauge and B.P. bulb to the fresh gas outlet, and perform the following test:
 - 4.11.23.1 Open the oxygen cylinder valve and allow the pressure to stabilize. (The cylinder pressure must be at least 1000 psi for this test.)
 - 4.11.23.2 Release any pressure that is indicated on the test gauge.
 - 4.11.23.3 Over the next 60 seconds, the test gauge should not show a pressure increase greater than 2 cm H₂O.
 - 4.11.23.4 Increase the pressure to 50 cm H₂O.
 - 4.11.23.5 The pressure should not drop more than 10 cm H₂O in the next 30 seconds.
 - 4.11.23.6 Disconnect the test gauge from the fresh gas outlet.
 - 4.11.23.7 Close the oxygen cylinder valve.
 - 4.11.23.8 The pressure should not drop more than 50 psi in two minutes.

- 4.11.23.9 Connect a volumeter to the fresh gas outlet, and reset the volumeter to zero.
- 4.11.23.10 Press the O₂ Flush button and observe the flow rate. It should be between 45 and 65 liters per minute.
- 4.11.23.11 Disconnect the volumeter from the fresh gas outlet.
- 4.11.24 Connect the absorber fresh gas hose to the fresh gas outlet.
- 4.11.25 Connect the pipeline hoses.
- 4.11.26 Perform the PMC Procedure given in Section 6.

4.12 Ventilator Controller (Bezel Assembly)

The ventilator controller assembly is attached to the inner shelf above the flowmeter housing and bellows box. Figure 4-12 shows the mounting screw locations, pneumatic and electrical connections to the ventilator controller.

4.12.1 Turn the System Power switch to STANDBY and remove AC power from the machine.

4.12.2 Disconnect all pipeline hoses and close the cylinder valves.

CAUTION: The controller circuit board contains static sensitive devices. Use ESD protection when handling the controller assembly. Static discharge can damage components on the circuit board.

4.12.3 Remove the flowmeter housing back cover. Be sure to disconnect the ventilator exhaust hose.

4.12.4 Remove the two screws holding the ventilator controller chassis to the inside shelf; remove the two screws holding the hex standoffs attached to the controller assembly. You will need to open the battery compartment and remove the battery for access to one of the screws.

4.12.5 Pull the assembly out from the front of the machine far enough to gain access to its connections.

4.12.6 Disconnect tube **D** from the fitting at the bottom of the inspiratory flow regulator.

4.12.7 Disconnect the power wiring harness from J2 on the controller circuit board.

4.12.8 Disconnect the switch and solenoid wiring harness from J1 on the controller circuit board.

4.12.9 Disconnect tube **B** from the fitting at the diaphragm valve.

4.12.10 Disconnect small diameter tubes **A** and **C**, and remove the controller assembly from the machine.

4.12.11 Position the replacement controller assembly at the front of the ventilator box and reconnect the four pneumatic lines.

4.12.12 Reconnect wire harnesses to J1 and J2 on the controller circuit board.

4.12.13 Slide the controller into the ventilator box until it is properly seated.

4.12.14 Reinstall the two screws to secure the controller chassis to the shelf.

- 4.12.15 Reinstall the two screws holding the hex standoffs, and reinstall the battery.
- 4.12.16 Reconnect the ventilator exhaust hose and reinstall the flowmeter housing back cover.
- 4.12.17 Perform the PMC Procedure given in Section 6.




Figure 4-12. Ventilator Controller (Bezel Assembly)

4.13 Bellows Valve Assembly

The bellows valve assembly is located in the bellows box directly above the bellows canister. Access to the valve assembly requires removal of the flowmeter housing back cover (which also covers the bellows box). Figure 4-13 shows the pneumatic connections and mounting arrangement of the valve assembly.

- 4.13.1 Turn the System Power switch to STANDBY and remove AC power from the machine.
- 4.13.2 Remove the flowmeter housing back cover.
- 4.13.3 Adjust the TIDAL VOLUME control to raise the volume indicator to its maximum setting.
- 4.13.4 Disconnect the breathing hose and the scavenger hose from the bellows assembly. Loosen the wing nuts and remove the bellows assembly.
- 4.13.5 Remove the canister from the bellows box by pulling it downward. (It is a press-fit.)
- 4.13.6 Disconnect tube **B** from the venturi, and tube **C** from the auto-ranging valve.
- 4.13.7 Remove the four screws securing the valve assembly to the bellows box.
- 4.13.8 Pull the valve assembly toward the back of the machine until the bottom plate of the assembly is able to drop down through the cutouts in the bottom lip on each side of the bellows box.
- 4.13.9 Install the replacement valve assembly - slide the assembly up into the bellows box and forward into position. (It may be necessary to rotate the pressure limit and tidal volume knobs until their drive slots are aligned with the shaft pins, in order to move the valve assembly into its correct position.)
- 4.13.10 Secure the valve assembly to the bellows box with the hardware that was previously removed.
- 4.13.11 Reconnect the tubing to the venturi, and the auto-ranging valve.
- 4.13.12 Reinstall the canister in the bellows box. Ensure that the markings are facing forward.
- 4.13.13 Reinstall the bellows assembly and tighten the wing nuts.
- 4.13.14 Reconnect the breathing and scavenger hoses to the bellows assembly.
- 4.13.15 Reinstall the flowmeter housing back cover.
- 4.13.16 Perform the PMC Procedure given in Section 6.




Figure 4-13. Bellows Valve Assembly

4.14 Caster Replacement

Each caster is retained by a set screw in the side of the lower frame rail as shown in Figure 4-14. Caster replacement requires that the machine be tilted to provide enough clearance for the caster stem to be withdrawn from the bottom of the frame rail.

WARNING: Do not tilt the machine more than 10 degrees or raise the casters more than 3½ inches from the floor. Failure to observe this precaution may result in a tip-over, causing personal injury. Vaporizers containing anesthetic agent may also be damaged.

- 4.14.1 Obtain a brace capable of supporting one side of the machine with its casters two to three inches from the floor.
- 4.14.2 Remove all unsecured equipment and accessories from the machine.
- 4.14.3 Lock the front casters.
- 4.14.4 Using at least two people, tilt the machine until the casters on one side are raised two to three inches from the floor, and position the support brace under the frame rail between the front and back casters.
- 4.14.5 Loosen the set screw until the caster can be removed.
- 4.14.6 Insert the replacement caster into the frame; align the threaded hole in the caster stem with the hole in the frame rail.
- 4.14.7 Tighten the caster stem set screw.
- 4.14.8 Using at least two people, tilt the machine, remove the support brace and carefully lower the machine to the floor.
- 4.14.9 Check for proper operation of the caster and ensure that the front casters lock properly.
- 4.14.10 Perform the PMC Procedure given in Section 6, including a vaporizer calibration verification.
- 4.14.11 Replace any unsecured equipment and accessories that were previously removed.




Figure 4-14. Caster Replacement

4.15 Battery Replacement

The battery is located in a compartment accessible at the back of the bellows box. Figure 4-15 shows the battery compartment door and battery wiring arrangement.

- 4.15.1 Turn the System Power switch to STANDBY and remove AC power from the machine.
- 4.15.2 Unscrew the two captive mounting screws on the battery compartment door, and remove the door.
- 4.15.3 Pull the battery from its compartment; note the wire colors and positions, and disconnect the wires from the tabs on the battery.
- 4.15.4 Connect the wires to the replacement battery in the same manner as the original, and install the battery in the battery compartment.
- 4.15.5 Reinstall the battery compartment door and secure it with the two captive mounting screws.
- 4.15.6 Restore AC power to the machine to allow the battery to charge. Allow 12 hours charging time for a fully discharged battery.

End of life battery disposal:

Dispose of the spent rechargeable, sealed lead-acid battery in conformance with local waste disposal regulations.




Figure 4-15. Battery Replacement

4.16 Power Supply Assembly

The power supply assembly comprises the power inlet connector, circuit breakers, primary power supply and regulator PCB. Later designs do not have the separate regulator PCB. The power supply is field serviced by replacing the complete assembly. You will need to remove the machine from its support frame for access to the power supply assembly. Figure 4-16 shows the mounting arrangement and electrical connections.

4.16.1 Disconnect all pipeline hoses and close all cylinder valves.

WARNING: Ensure that AC power is removed from the machine before removing the power supply. Failure to observe this precaution may cause injury by electric shock.

4.16.2 Turn the System Power switch to STANDBY and disconnect the AC power cord from the power inlet on the back of the machine.

4.16.3 Disable the circuit breakers by pulling out each button with a knife or sharp object.

4.16.4 Remove the following items from the machine: external monitors, cylinders, vaporizer and display.

4.16.5 Disconnect the fresh gas hose.

4.16.6 Disconnect the ventilator breathing and scavenger hoses, and the sensor interface connections.

4.16.7 Loosen the wingnut on the absorber arm, and remove the absorber assembly.

4.16.8 Separate the machine from its support frame as follows (see Figure 4-16A):

4.16.8.1 Remove the two drawer support screws.

4.16.8.2 Loosen (do not remove) the two vertical support arm screws on the inside of the support frame, on each side.


4.16.8.3 Carefully lift the anesthesia machine from the support frame (the vertical support arms remain attached to the machine).

4.16.9 Carefully place the machine upside-down on a suitable surface.

4.16.10 Remove the six screws holding the power supply assembly to the housing, and lift out the assembly far enough for access to its cables.

REPLACEMENT PROCEDURES (continued)

NM MOBILE



SV00317

Figure 4-16. Power Supply Assembly Replacement




Figure 4-16A. Removal of Machine from Support Frame

- 4.16.11 Carefully disconnect the following from the PCB on the assembly:
- Output wire harness from J1
 - 14 cond. ribbon cable from J2
 - Battery wire harness from J6
 - Power applied lamp from J5
- 4.16.12 Reconnect the cables and wire harnesses to J1, J2, J6, and J5 on the replacement power supply assembly.
- 4.16.13 Position the power supply assembly in the housing and secure it with the six mounting screws through the bottom of the power supply.
- 4.16.14 Reattach the machine to its support frame as follows:
- 4.16.14.1 Carefully fit the vertical support arms into the frame; note that the right side arm (facing the front of the machine) must pass through the upper drawer support.
 - 4.16.14.2 Tighten the two support arm screws on each side.
 - 4.16.14.3 Reinstall the two drawer support screws.
- 4.16.15 Reinstall the display, vaporizer, absorber assembly and cylinders.
- 4.16.16 Reconnect the fresh gas hose.
- 4.16.17 Reinstall all accessories that were previously removed; restore all breathing, scavenger and sensor interface connections. Refer to the Installation Instructions in the *Narkomed Mobile Setup and Installation Manual*.
- 4.16.18 Enable the circuit breakers by pushing in each button.
- 4.16.19 Reconnect the AC power cord.
- 4.16.20 Perform the PMC procedure given in Section 6.

4.17 Processor Assembly

Access to the processor assembly requires removal of the top cover from the monitor housing. Figure 4-17 shows the processor mounting arrangement and location of cables that must be disconnected.

- 4.17.1 Turn the System Power switch to STANDBY and remove AC power from the machine.
- 4.17.2 Remove the screws securing the top cover of the monitor housing. Lift the cover and disconnect its ground wire.

CAUTION: The processor board contains static sensitive devices. Use ESD protection when handling the processor assembly. Static discharge can damage components on the circuit board.

- 4.17.3 Disconnect the following items from their connectors on the processor board:

J105: 30 psi switch	J12: Speaker	J4: Keypad cable
J204: Display cable	J102: Spiromed interface	J103: O ₂ Sensor interface
J18: Serial port	J14: Pwr supp wire harness	Pressure sensor tubing

- 4.17.4 Remove the screws securing the processor assembly to the monitor housing, and lift out the processor assembly.
- 4.17.5 Inspect the jumper on JP6 (top PCB) on the replacement processor assembly and ensure it is installed on both pins of JP6 before installing the processor assembly. This jumper is not placed on both pins of JP6 in order to prevent backup battery drain while the processor assembly is in stock.

NOTE: On later design boards this jumper is at JP2 (see illustration).


- 4.17.6 Make sure the configuration jumper is across Pins 1 and 2 of JP101 on the GS personality card. This will configure the software for the ultrasonic flow sensor.

NOTE: On later design boards this jumper is at JP1 (see illustration).

- 4.17.7 Install the replacement processor using the hardware that was removed in the previous step.
- 4.17.8 Reconnect the pneumatic tubing and cables that were previously disconnected.
- 4.17.9 Reconnect the ground wire, and reinstall the monitor housing top cover.
- 4.17.10 Restore power to the machine and observe the Power-Up Diagnostic display (see Section 2) to verify that the replacement processor is working properly.
- 4.17.11 Perform the PMC Procedure given in Section 6.

REPLACEMENT PROCEDURES (continued)

NM MOBILE



SV00318

Figure 4-17. Processor Assembly

4.18 Display Assembly

The display assembly comprises the keypad & display panel, its housing, cable, and mounting rod. The display assembly is attached to the joint assembly (remote bar and monitor support arm) by a screw as shown in Figure 4-18.

4.18.1 Turn the System Power switch to STANDBY.

4.18.2 Unplug the display cable from its port on the side of the machine.

4.18.3 Remove the display assembly retainer screw from the joint assembly (support arm) and lift out the display assembly.

NOTE: If applicable, install a slide-latch retainer kit for the remote display cable. Refer to the Spare and Replacement Parts section for details.

4.18.4 Install the replacement display assembly in the support arm; reinstall and tighten the display assembly retainer screw.

4.18.5 Connect the display cable to its port on the side of the machine.

4.18.6 Power up the machine and observe the display for correct operation; exercise all of the keypad functions to verify their operation.




Figure 4-18. Display Assembly

4.19 Breathing Pressure Hose Assembly

The breathing pressure hose assembly consists of a pressure sensor adapter with a hose attached, with a quick-connect fitting at the outboard end of the hose. The hose connects to the breathing pressure interface panel on the monitor housing. The pressure sensor adapter is inserted between the inspiratory valve dome and the oxygen sensor housing as shown in Figure 4-19.

Removal:

- 4.19.1 Turn the System Power switch to STANDBY.
- 4.19.2 Disconnect the hose fittings from the interface panel and the absorber.
- 4.19.3 Pull the oxygen sensor housing from the pressure sensor adapter (it is a press fit).
- 4.19.4 Pull the pressure sensor adapter from the absorber top dome (it is a press fit).

Installation:

- 4.19.5 Insert the pressure sensor adapter into the absorber top dome.
- 4.19.6 Insert the oxygen sensor housing into the pressure sensor adapter.
- 4.19.7 Connect the short hose to the fitting on the absorber assembly.
- 4.19.8 Connect the longer hose to the breathing pressure interface panel.
- 4.19.9 Restore power to the machine and perform the pressure monitor diagnostic test given in Section 2 to verify operation of the system.




Figure 4-19. Breathing Pressure Hose Assembly

4.20 Manual/Automatic Selector Valve

The manual/automatic selector valve is part of the absorber assembly, and is fixed to the absorber flange with three screws. Figure 4-20 shows the components that are mounted on the valve, and the attaching hardware at the absorber flange.

- 4.20.1 Turn the System Power switch to STANDBY.
- 4.20.2 Disconnect the ventilator hose from the rear 22 mm terminal on the selector valve, and remove the terminal.
- 4.20.3 Disconnect the bag hose from the bottom 22 mm terminal on the selector valve, and remove the terminal.
- 4.20.4 Disconnect the scavenger hose from the 19 mm terminal on the APL valve.
- 4.20.5 Loosen the locking ring on the APL valve with a spanner wrench (P/N S010058), and unscrew (counter-clockwise, viewed from the top) the APL valve from the manual/automatic selector valve.
- 4.20.6 Remove the three screws securing the selector valve to the absorber flange, and remove the valve.
- 4.20.7 Ensure that the O-ring at the absorber flange is in good condition, and attach the replacement manual/automatic selector valve to the absorber flange using the three screws and lock washers that were previously removed.
- 4.20.8 Reinstall the APL valve: Ensure that the fiber washer is in place, and screw the valve in to its desired position. Tighten the locking ring with the spanner wrench.
- 4.20.9 Reinstall the 22 mm hose terminals on the selector valve; ensure that their O-rings are in good condition.
- 4.20.10 Re-connect the ventilator hose to the rear 22 mm terminal, and re-connect the bag hose to the bottom 22 mm terminal.
- 4.20.11 Re-connect the scavenger hose to the 19 mm terminal on the APL valve.
- 4.20.12 Perform the PMC Procedure given in Section 6.




Figure 4-20. Manual/Automatic Selector Valve

4.21 Ultrasonic Flow Sensor

The ultrasonic flow sensor is mounted on a bracket attached to the inspiratory valve mount on the top of the absorber. A short connector hose joins the expiratory valve to the flow sensor. The mounting and connection arrangement is shown in Figure 4-21.

- 4.21.1 Turn the System Power switch to STANDBY.
- 4.21.2 Disconnect the sensor plug from the volume sensor interface panel on the monitor housing.
- 4.21.3 Disconnect the breathing hose from the flow sensor.
- 4.21.4 Remove the connector hose from the other port of the flow sensor by unscrewing the retaining ring on the hose.
- 4.21.5 Pull the the flow sensor up and off the mounting bracket
- 4.21.6 Slide the replacement sensor onto the mounting bracket - oriented in the same manner as the original.
- 4.21.7 Join the connector hose to the threaded port on the flow sensor, and reconnect the breathing hose to the sensor.
- 4.21.8 Connect the sensor plug to the volume sensor interface panel on the monitor housing.
- 4.21.9 Restore power to the machine and perform the respiratory flow monitor calibration procedure given in Section 5.
- 4.21.10 Perform the PMC Procedure given in Section 6.




Figure 4-21. Ultrasonic Flow Sensor

4.22 Oxygen Sensor

The oxygen sensor is installed in the pressure sensor adapter on top of the inspiratory valve. Figure 4-22 shows the arrangement of the sensor capsule and its housing, and also its interface cable connection on the side of the monitor housing.

- 4.22.1 Turn the System Power switch to STANDBY.
- 4.22.2 Pull the oxygen sensor housing from the oxygen sensor adapter. (It is a press fit.)
- 4.22.3 Unscrew the cover from the sensor housing and remove the sensor capsule.
- 4.22.4 Remove the replacement sensor capsule from its shipping container and install it in the housing. Ensure that the copper rings on the capsule mate with the electrical contacts in the sensor housing.
- 4.22.5 Wait 15 minutes to allow the sensor capsule to stabilize.
- 4.22.6 Restore power to the machine and perform the Zero and the 21% calibration procedure for the oxygen monitor given in Section 5.
- 4.22.7 Insert the oxygen sensor assembly into the pressure sensor adapter.
- 4.22.8 Perform the PMC Procedure given in Section 6.




Figure 4-22. Oxygen Sensor

5.0 Adjustment and Calibration Procedures

Equipment Required:

- Test Gauge for setting cylinder pressure regulator, NAD Part No. 4114807
- Oxygen Monitor for adjusting Oxygen Ratio Controller
- Test fixture with breathing pressure line connector, TEE connector, gauge, and inflation device, for breathing pressure monitor calibration

5.1 Cylinder Pressure Regulator Adjustment

5.1.1 Turn the System Power switch to STANDBY.

5.1.2 Remove the flowmeter housing back cover.

NOTE: Minimum cylinder pressures for this adjustment shall be:
N₂O & CO₂: 600 psi; O₂, Air, He, He/O₂, N₂: 1000 psi.

NOTE: Figure 5-1 shows test connections for the O₂ regulator adjustment. If you are adjusting the N₂O regulator, connect the test gauge in series with the N₂O pipeline hose.

5.1.3 Connect test pressure gauge (P/N 4114807) between machine's pipeline inlet connector and the pipeline supply hose.

5.1.4 Open the cylinder valves and turn the System Power switch to ON.

5.1.5 Set the oxygen flow to 4 liters per min. (If you are adjusting the N₂O regulator, also set the nitrous oxide flow to 4 liters per minute.)

5.1.6 Depress the push button on the test device.

5.1.7 Release the push button. After the pressure decay stabilizes, the gauge should indicate 46 psi.

5.1.8 Remove the acorn nut from the regulator to expose the adjusting screw. For N₂O, turn the screw until the test gauge indicates 46 psi. (50 psi for CSA machines). For O₂ and Air, use the compensated regulator output setting based on the cylinder pressure given in the following table.

CAUTION: Based on information supplied by the cylinder regulator manufacturer, when the regulator is used for O₂ or Air, its output pressure will decrease 0.5 psi for every 100 psi increase in cylinder pressure above 1000 psi. Currently, these pressure regulators are calibrated at 47 psi with a cylinder supply of 1000 psi. If a 2000 psi cylinder is then installed, the regulator output will be 42 psi. This change in output must be compensated for to provide accurate performance throughout the cylinder's working range.

NOTE: Cylinder pressure compensation for the N₂O regulator is not required.

Cylinder Pressure (psi)	Compensated Regulator output setting (psi)	Compensated Regulator output tolerances (-4/+2)
2000	42 (*45)	38 - 44 (*41 - 47)
1800	43 (*46)	39 - 45 (*42 - 48)
1600	44 (*47)	40 - 46 (*43 - 49)
1400	45 (*48)	41 - 47 (*44 - 50)
1200	46 (*49)	42 - 48 (*45 - 51)
1000	47 (*50)	43 - 49 (*46 - 52)

* Canada settings

5.1.9 Verify the adjusted regulator maintains its compensated output tolerance.

5.1.10 If the O₂ cylinder regulator output was adjusted, perform the following test:

5.1.10.1 Set the ventilator to Volume Mode ventilation.

5.1.10.2 Set the Auxiliary O₂ and Fresh Gas O₂ flows to 10 L/min.

5.1.10.3 Press the O₂ Flush and verify the Lo O₂ Supply alarm is not active. If the alarm is active refer to Section 5.2: Oxygen Supply Pressure Alarm Switch Adjustment.

5.1.11 Reinstall the acorn nut on the regulator.

5.1.12 Close the cylinder valves and allow pressure to drain from the system.

5.1.13 Close all of the flow control valves and set the System Power switch to STANDBY.

5.1.14 Disconnect the test gauge.

5.1.15 Reinstall the flowmeter housing back cover.

5.1.16 Connect the pipeline hoses.

5.1.17 Perform the PMC Procedure given in Section 6.

This page intentionally left blank




Figure 5-1. Cylinder Pressure Regulator Adjustment

5.2 Oxygen Supply Pressure Alarm Switch Adjustment

- 5.2.1 Disconnect all pipeline hoses and turn the System Power switch to ON.
- 5.2.2 Open the oxygen cylinder valve.
- 5.2.3 Set the oxygen flow to 5 liters per min.
- 5.2.4 Open the other gas flow control valves to drain pressure from the system.
- 5.2.5 Close the O₂ cylinder valve, and close the flow control valves. Press the O₂ Flush valve to drain oxygen pressure from the system.
- 5.2.6 Turn the System Power switch to STANDBY.
- 5.2.7 Remove the rear cover from the flowmeter housing.
- 5.2.8 Connect test pressure gauge (P/N 4114807) between machine's oxygen pipeline inlet connector and the oxygen pipeline supply hose.
- 5.2.9 Open the O₂ cylinder valve and turn the System Power switch to ON.
- 5.2.10 Set the oxygen flow to 200 mL per min.
- 5.2.11 Close the oxygen cylinder valve.
- 5.2.12 As the pressure drops, the O₂ SUPPLY alarm should activate when the pressure is between 40 and 34 psi as shown on the test gauge.
- 5.2.13 If the alarm activates when the pressure is below 34 psi or above 40 psi, turn the adjustment set screw (see illustration); repeat the test and adjust as necessary to bring the set point into the correct range.
- 5.2.14 Turn the System Power switch to STANDBY.
- 5.2.15 Disconnect the test gauge.
- 5.2.16 Reinstall the rear cover and its retaining screws.
- 5.2.17 Connect the pipeline hoses.
- 5.2.18 Perform the PMC Procedure given in Section 6.




Figure 5-2. Oxygen Supply Pressure Alarm Switch Adjustment

5.3 Oxygen Ratio Controller (ORC) Adjustment

- 5.3.1 Remove the rear cover of the flowmeter housing.
- 5.3.2 Connect a calibrated oxygen monitor to the fresh gas outlet.
- 5.3.3 Connect the pipeline hoses.
- 5.3.4 Turn the System Power switch to ON.
- 5.3.5 Set the O₂ flow to 8 l/min.
- 5.3.6 Set the N₂O flow to 8 l/min.
- 5.3.7 Set the O₂ flow to 800 ml/min for one (1) minute. Verify that the O₂ concentration is between 21% and 29% (N₂O flow of 2.7 to 3.0 l/min.). If needed, loosen the locknut on the ORC and turn the adjusting screw (counter-clockwise to decrease N₂O flow, clockwise to increase N₂O flow) to achieve a nominal O₂ concentration of 25%.
- 5.3.8 Repeat the previous three steps until no further adjustment is needed. Tighten the locknut.
- 5.3.9 Adjust the oxygen flow to a point where the nitrous oxide flowmeter indicates 8 l/min.
- 5.3.10 Verify that the O₂ concentration is between 21% and 29% (O₂ flow of 2.1 to 3.3 l/min.).
- 5.3.11 Slowly decrease the oxygen flow to 800 l/min. The nitrous oxide flow should decrease proportionally, and the O₂ concentration should remain between 21% and 29%.
- 5.3.12 Close the O₂ flow control valve, and fully open the N₂O flow control valve. Verify that the O₂ concentration is between 22% and 31%.
- 5.3.13 Close the N₂O flow control valve and turn the System Power switch to STANDBY.
- 5.3.14 Replace the flowmeter housing rear cover.
- 5.3.15 Perform the PMC Procedure given in Section 6.




Figure 5-3. Oxygen Ratio Controller (ORC) Adjustment

5.4 Oxygen Sensor Calibration

5.4.1 Turn the System Power switch to ON.

5.4.2 Enter the Main Service Screen and select the Service Mode (ref. Section 2).

5.4.3 Enter the Oxygen Monitor Service Screen.

5.4.4 Zero Calibration

5.4.4.1 Remove the oxygen sensor capsule from its housing and allow several minutes for the displayed offset readings to stabilize.

NOTE: The difference between the displayed CELL A and CELL B readings should be no greater than 8.

5.4.4.2 Press the key next to ZERO to store the current values as the new zero calibration.

5.4.4.3 Reinstall the sensor capsule in its housing.

5.4.5 21% Calibration

5.4.5.1 Expose the sensor to ambient air only (away from any open part of the breathing system) and allow it to stabilize for several minutes.

5.4.5.2 Press the key next to EXIT to return to the Main Service Screen. Press the key next to EXIT again to return the display to normal operation.

5.4.5.3 Press the CAL key to initiate the 21% O₂ calibration.

During calibration, the LED next to the CAL key lights, and the label CAL appears in the oxygen monitor window. Following successful calibration, the currently sensed oxygen concentration appears in the oxygen monitor window.

5.4.5.4 When calibration is complete, reinstall the sensor assembly in the inspiratory valve dome.

NOTE: If the O₂ sensor will not calibrate properly, refer to the Oxygen Monitoring section of the *Narkomed Mobile OPERATOR'S INSTRUCTION MANUAL* for further information.




Figure 5-4. Oxygen Sensor Calibration

5.5 Breathing Pressure Monitor Calibration

5.5.1 Turn the System Power switch to ON.

5.5.2 Enter the Main Service Screen and select the Service Mode (ref. Section 2).

5.5.3 Proceed to the Pressure Monitor Service Screen.

5.5.4 Zero Calibration

5.5.4.1 Disconnect the breathing pressure hose from the interface panel on the monitor housing, and let the current pressure value stabilize.

5.5.4.2 Press the key next to ZERO to store the current value as the new zero.

5.5.5 Span Calibration

5.5.5.1 With a test fixture connected as shown in Figure 5-5, apply a pressure of 50 cm H₂O to the breathing pressure interface panel.

5.5.5.2 When the displayed current value is stabilized, press the key next to SPAN to store the current value as the new span calibration.

5.5.6 Disconnect the test fixture; reconnect the breathing pressure hose to the interface panel.

5.5.7 Press the key next to EXIT to return to the Main Service Screen.




Figure 5-5. Pressure Monitor Calibration

[RETURN TO THIS MANUAL'S TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

6.0 PMC PROCEDURE, NARKOMED MOBILE

The procedures in this section shall be performed in their entirety each time a component is removed, replaced, calibrated, adjusted and during all scheduled Periodic Manufacturer's Certification (PMC) visits. A PMC Checklist form, P/N 4115597 is available from Draeger Medical, Inc. and shall be completed by the Technical Service Representative each time a PMC is performed. Steps in the procedure marked with (✓) require a response at the corresponding line on the checklist form.

Space is also provided on the PMC checklist form to record the results of a vapor concentration test. Refer to the current Anesthesia Equipment & Monitoring System Service Information CD-ROM Service Procedures section for vapor concentration verification procedures.

NOTE: Test equipment listed below with an asterisk (*) requires calibration at a maximum interval of one year. Verify the dates on test equipment calibration labels. DO NOT USE any test equipment having an expired calibration date. Notify your supervisor immediately if any equipment is found to be out of calibration. In the space provided at the bottom of the PMC checklist form, record the Model and ID number of all calibrated test equipment used.

In the space provided at the bottom of the PMC checklist form, record the Model and ID number of all calibrated test equipment used. Also record the calibration due dates. Examples are: multimeter, digital pressure meter, Riken gas analyzer, safety analyzer, volumeter, trace gas analyzer, simulators.

Test Equipment Required:

- *Electrical Safety Analyzer (Biotek 501 Pro or equivalent)
- *Pressure Gauge with DISS Adapters (P/N 4114807 or equivalent)
- *Flowmeter 0-250 ml min. (P/N S000081 or equivalent)
- *Volume Meter (P/N 2212300 or equivalent)
- *Digital Pressure Manometer (SenSym PDM 200CD or Equivalent)
- *Riken Gas Indicator (Model 18H, or 1802D or equivalent)
- Stop Watch
- Test Lung (P/N 8401892)

Materials Required:


- Breathing Bag 3 liter (P/N 9995330 or equivalent)
- Patient Circuit: Y-piece, elbow, 2x 32" x 22mm hoses
- Hose 22 mm x 32" (P/N 9995132)
- Fresh Gas Outlet Volume Test Device (P/N S010158 or equivalent)
- Fresh Gas Leak Test Adapter (P/N 4115041 or equivalent)
- Volumeter/Fresh Gas Adapter (P/N 4115042)

Materials Required (continued):

- Test Terminal 2x (P/N 4104389 or equivalent)
- Breathing System Leak Test Device (P/N S010159 or equivalent)
- PDM/Suction Adapter (P/N 4115038)
- Scavenger Adapter (P/N 4108114)
- Pressure Monitor Test Adapter (P/N 4115043 or equivalent)

Key test equipment and materials illustrations are shown on following pages.



TEST EQUIPMENT AND FIXTURE ILLUSTRATIONS




SP18201

TEST EQUIPMENT AND FIXTURE ILLUSTRATIONS - Continued


4114807 PRESSURE TEST ASSEMBLY , WITH ADAPTERS




4114830-002




4114830-001




4114830-004



4114830-003




4114830-006




4114830-005


TEST EQUIPMENT AND FIXTURE ILLUSTRATIONS - Continued




TEST EQUIPMENT AND FIXTURE ILLUSTRATIONS - Continued



4115038
PDM TO PATIENT SUCTION ADAPTER



4108114
SCAVENGER ADAPTER




4115043
PDM TO
MONITOR ADAPTER


SP18204

TEST EQUIPMENT AND FIXTURE ILLUSTRATIONS - Continued


S000081
FLOW METER
TEST STAND




2212300
MINUTE
VOLUMETER



8401892
SIEMENS TEST LUNG



4104389
TEST TERMINAL
ADAPTER



SP18205

Periodic Manufacturer's Certification General Instructions

The purpose of these procedures is to provide detailed instructions for performing a Periodic Manufacturer's Certification (PMC) inspection on the Narkomed Mobile anesthesia machine.

A PMC consists of a complete Periodic Manufacturer's Service procedure and a certification level inspection based on Draeger Medical, Inc. Recommendations and equipment performance. Additional inspections are also performed to ensure proper product labeling.

Several additional documents have been created to assist the technician through the process. Following is a brief description of the purpose of each document.

Field Service Procedure:

Periodic Manufacturer's Certification Forms - Part Number SP00175. This procedure illustrates the sample checklists with typical periodic maintenance items filled in, including vapor concentration verification tests, parts replaced, general comments and certification levels. Also included are sample PMC labels marked to show several levels of certifications. An excerpt from Draeger Medical, Inc.'s *Anesthesia System Risk Analysis and Risk Reduction* is included, and also a sample of an Executive Summary to be furnished to the hospital's Risk Manager or Chief of Anesthesia.

Field Service Procedure:

DMI Recommendation Guidelines Index Anesthesia Systems - Part Number S010250. This Guideline was created to provide an assessment of each machine's certification. It contains various comprehensive overviews of possible equipment conditions and their associated certification levels.

The first list in the Recommendation Guidelines is a reference chart for machine certification based on equipment status. The second is an abbreviated summary of all DMI Recommendations and Failure Codes including the Condition Number, Equipment Condition, Recommended Corrections, Certification Code, and Tests Affected when applicable.

There is also a matrix classified as "Failure Codes" which identifies the correct manner in which to document equipment tests that fail, or were unable to be performed due to circumstances beyond the control of the service technician performing the inspection. (Ex: Air cylinder supply is unavailable to perform an Air High Pressure Leak test.) The Failure Codes section also indicates suggested resolution of the situation. Failure Code numbers begin at 34 and use the same certification levels strategy, and carry the same weight as NAD Recommendation equipment condition codes.

The final matrix is the most comprehensive index sorted by machine model and includes Equipment Condition, Certification Code, and DMI Recommendations. It also specifies any suggested upgrade path including ordering information that should be taken such as installing a Bellows with Pressure Limit Control 4109664-S01 Kit, after market modification kit to a machine not equipped with pressure limit control.

The letters A, B, C, D and the Roman Numerals I, II are used as codes in the individual matrix for each model of anesthesia machine. The letters A, B, C, and D are used in descending order to indicate the certification level of the equipment. They are as follows:

- A = Certified
- B = Certified with Recommendations
- C = Conditionally Certified
- D = No Certification

Roman Numerals I and II do not affect the certification level but rather are provided to give further instructions to the end user as follows:

I = The system in its present configuration shall only be used with a CO2 monitor incorporating an apnea warning. The operator of the system is advised to frequently scan the CO2 readings and alarm thresholds.

II = The present configuration of equipment requires that the unit operate at all times with an oxygen analyzer that includes a low oxygen warning. The operator of the system is advised to frequently scan the oxygen readings and alarm limits.

Following is an explanation of machine certification levels:

Certified- No DMI Recommendations or Failure Codes apply to machine being inspected. (Only item number 33 - "No Recommendations" shall apply for this certification level.)

Certified with Recommendations- A numbered DMI Recommendation or Failure Code with a code of B applies to the machine being examined.

Conditionally Certified- A numbered DMI Recommendation or Failure Code with a code of BC I or BC II applies to the machine being examined.

No Certification- A numbered DMI Recommendation or Failure Code with a code of D applies to the machine being examined.

When multiple recommendations apply, "No Certification" would take precedence over "Conditionally Certified" and "Certified with Recommendations". "Conditionally Certified" would take precedence over "Certified with Recommendations".

For example:

A **Narkomed 2B** could have DMI Recommendation number 21 and Failure Code 61.1 that apply. 21 - No ventilator pressure limit control. Code is B. 61.1 - Enflurane agent is unavailable to test. Code is BC. Correct certification for this machine is BC, which means **CONDITIONALLY CERTIFIED WITH RECOMMENDATIONS**.

A **Narkomed 4** could have DMI Recommendation numbers 14 and 21 apply. 14 - CO2/Agent monitor exhaust port is not properly connected to the waste gas scavenger. Code B.

21 - No ventilator pressure limit control. Code B. The correct certification for this machine is B, which means "CERTIFIED WITH RECOMMENDATIONS".

A **Narkomed 2B, 2C or GS** could have DMI Recommendation 30 apply.

30 - Anesthesia machine is equipped with inhalation anesthesia vaporizers without an agent analyzer in the breathing system. Code B. The correct certification for this machine is B, which means "CERTIFIED WITH RECOMMENDATIONS".

A **Narkomed 6000** could have no DMI Recommendations or Failure Codes apply. The correct certification level for this machine is Code A, "CERTIFIED".

The correct certification for this machine is A, which means "CERTIFIED".

Code, D also means "NO CERTIFICATION", also means the machine shall not receive a Periodic Manufacturer's Certification label. The machine shall receive a "WARNING - This System Is Not Certified" label, P/N 4114857. This label shall be placed at a prominent location on the right side of the machine after all other previous PM and "Vigilance Audit® Validation" labels have been removed.

PM Certification Procedure for Narkomed Mobile Anesthesia System

1. Use the PM Certification form for Narkomed Mobile Anesthesia Systems (P/N 4115597).
2. Completely fill in the header information.
3. All Narkomed Mobile machines are equipped with Humphrey valves. No "MJV-2 LUBRICATION" is needed to be performed for this type of anesthesia machine.
4. Replace the VENTILATOR RELIEF VALVE DIAPHRAGM every 12 months in accordance with SP00075. Place a check mark and the replacement date at "VENT VALVE REPLACEMENT" line on the Periodic Manufacturer's Certification form.
5. If machine is equipped with a HALOTHANE Dräger Vapor 19 or 19.1 vaporizer, determine if vaporizer must be inspected for soil condition one. Check the serial number plate located on the rear of the vaporizer for a plus (+) preceding the serial number. A HALOTHANE vaporizer serial number not preceded with a (+) must be tested for soil in accordance with SP00073. If vaporizer does not need to be inspected, indicate so with a plus (+) next to the "Vapor Inspection (H)" line on the Vigilance Audit form. If vaporizer is soil condition 0, indicate so with "SOIL 0" written next to the "Vapor Inspection (H)" line on the Vigilance Audit form. If vaporizer is soil condition one, indicate so with "SOIL 1" written next to the "Vapor Inspection (H)" line on the Vigilance Audit form. Place a "CAUTION DO NOT USE" label (part # 4114327) on the vaporizer, and issue a departmental alert. The TSR shall also seek permission from the equipment operator to remove the failed vaporizer from the machine and apply a replacement vaporizer or an adapter block onto the mount. All "SOIL 1" vaporizers must be removed from service for machine to receive certification.
6. Perform the vapor concentration test on all Dräger vapor vaporizers in accordance with SP00073 at a six month maximum interval. Perform the vaporizer concentration test on all Desflurane vaporizers in accordance with SP00091 for fixed mount vaporizers and SP00189 for user removable D-tec vaporizers at a six month maximum interval. For every vaporizer tested, fill out a "VAPOR VAPORIZER CALIBRATION CHECK" label (part # S010016). Information on this label shall include your signature, type of agent, date tested, a No Agent To Test or the test results @ 1%, 2.5%, 4% for H, E, I, or S vaporizers, or @ 4%, 10%, 12%, 16% for Desflurane vaporizers, and a PASS or FAIL indication. This label shall be attached to the upper right side of the vaporizer. If vaporizer fails the concentration verification, internal leak, or exclusion system tests, check "NO" in the "RECOMMENDED FOR USE" section on the PM Certification form. Place a "CAUTION DO NOT USE" label (part # 4114327) on the vaporizer, and issue a departmental alert. The TSR shall also seek permission from the customer to remove the failed vaporizer from the machine and install a replacement vaporizer or an adapter block onto the mount. All nonfunctional vaporizers must be removed from service for machine to receive certification.

7. Proceed with PM Certification procedure. If any tests fail refer to the "Failure Codes" listing in DMI Recommendations Guidelines Index (P/N S010250) to determine correct certification level starting point. Failure codes shall be documented on the "RECOMMENDATIONS / GENERAL COMMENTS" section of the PM Certification form and on the Executive Summary. If a test fails that has not been identified by the "Failure Codes" list, consult with Draeger Medical, Inc. to assess the proper certification level.
8. Based on the "EQUIPMENT CONDITION" inspect the machine for any "DMI RECOMMENDATIONS" that would apply. Use the Narkomed Mobile section of the "RECOMMENDATION GUIDELINES INDEX" (P/N S010250). Note all applicable DMI recommendations on the Executive Summary. NOTE: If using a carbon form, indicate the Equipment Condition number and to see reverse side under the "RECOMMENDATIONS / GENERAL COMMENTS" section of the form.
9. Determine the correct certification level of the machine based on the combined lowest common denominator of "Equipment Conditions" and "Failure Codes". If the machine is at least conditionally certified fill out the "PM CERTIFICATION" label. Check the box(s) on the validation label where appropriate. Write the month and year, (three months from date of PM Certification) next to "NEXT VISIT DUE:". If certification level is "D", machine shall not receive a "PM CERTIFICATION" label. Any machine not receiving a PM Certification label shall receive a "WARNING NOT CERTIFIED" label, P/N 4114857. This label shall be placed at a prominent location on the left side of the machine after all other previous PMS and Vigilance Audit Validation labels have been removed.
10. In the "CERTIFICATION LEVEL" section of the PM Certification form, record the last visit certification level, the current certification level and the next visit due month and year, (three months from date of PM Certification) in the spaces provided.
11. If applicable, remove the previous PM CERTIFICATION VALIDATION label and attach the new label (P/N S010006 w/phone #, or S010007 w/o phone #) in a prominent location on the rear of the anesthesia machine.
12. Check the appropriate boxes on the "PM CERTIFICATION NOTICE" label, (part # S010011). If the machine is not certified, the last box of this notice label shall be marked. Attach this notice near the flowmeter shield of the anesthesia machine.
13. Have the customer sign each PM Certification form or the Executive Summary, and review any Failure Codes equipment conditions and DMI Recommendations with the customer.
14. Return top copy to Draeger Medical, Inc. Service Department, keep middle copy for service organization records, give bottom copy to customer.

(✓) **6.1 SELF-DIAGNOSTICS**

- 6.1.1 Turn the System Power switch to ON and verify the "ON" LED is lighted?
- 6.1.2 Verify all LED's on the keypad and ventilator displays are lit.
- 6.1.3 Verify that the following is displayed on the Display:

NARKOMED
COPYRIGHT 1997-1998 NAD. INC.
VERSION X.XX NARKOMED
SOFTWARE ID: XXXX

FIRMWARE	PASS
RAM	PASS
VIDEO	PASS
A/D CONVERTER	PASS
AUDIO - PRIMARY	PASS
BACKUP	PASS
SERIAL I/O	PASS
CLOCK	PASS
NON-VOLATILE MEMORY	PASS
FUNCTIONAL	

- (✓) 6.1.4 Record the machine software version on the header of the checklist form.

6.2 ELECTRICAL SAFETY- One Year Service Interval; Due Date _____

- (✓) 6.2.1 Battery Check & Ground Continuity

6.2.1.1 Unplug the AC power cord for all devices mounted to the machine that may provide an alternate path to earth ground, such as a Desflurane vaporizer.

6.2.1.2 Unplug the machine's AC power cord and plug the power cord of the safety analyzer into this AC receptacle.

NOTE: Do not plug the safety analyzer power cord into a line isolation monitor circuit, as inaccurate readings may occur.

NOTE: The BIOTECH 501 PRO will automatically test the source outlet for open ground (or ground resistance of 31 Ohms or higher), reverse polarity, open neutral and open line. (The latter two conditions will prevent the analyzer from powering up.)

6.2.1.3 Is the green "AC POWER" LED off? __ (Y)

6.2.1.4 Does the "AC PWR FAIL" message appear in the Advisory display? __ (Y)


- 6.2.1.5 Press and hold the "BATTERY TEST" button.
- 6.2.1.6 Is green Battery Test LED lighted? __ (Y)
- 6.2.1.7 Release the "BATTERY TEST" button.
- 6.2.1.8 Turn on the safety analyzer and configure it to measure "Resistance" according to the manufacturers documentation. Attach a test lead to the "SINGLE LEAD" connector of the analyzer. Connect the other end of the test lead to the AC receptacle ground socket on the safety analyzer. Verify a displayed resistance of 0.00 ohms or, if necessary, calibrate the unit according to manufacturers documentation.
- 6.2.1.9 Plug the machine's AC power cord into the safety analyzer.
- 6.2.1.10 Apply the analyzer's test lead to a cylinder yoke bolt.
- 6.2.1.11 What is the value displayed on the safety analyzer? __ ohm (0-0.1)

(✓)

6.2.2 Circuit Isolation

- 6.2.2.1 Disconnect the respiratory volume sensor cord from the interface panel.

- 6.2.2.2 Using a multimeter set to its highest resistance range apply the test leads between the yoke bolt and circuit common at the volume interface test pin. Use the position indicated on the illustration for the ultrasonic flow sensor connector. There shall be no continuity between these two points.



- 6.2.2.3 Reconnect the respiratory volume sensor cord to the interface panel.

6.2.3 Chassis Leakage Current

- 6.2.3.1 Apply the analyzer test lead to a cylinder yoke bolt.
- 6.2.3.2 Configure the safety analyzer to the measure "CHASSIS LEAKAGE CURRENT" according to the manufacturers documentation.

6.2.3.3 Record the total leakage current with the Polarity and Ground switches set as follows:

		<u>Ground</u>	<u>Polarity</u>
(✓)	6.2.3.3.1	Normal	Normal
(✓)	6.2.3.3.2	Open	Normal
(✓)	6.2.3.3.3	Open	Reversed
(✓)	6.2.3.3.4	Normal	Reversed

6.2.3.4 Verify that the leakage current is 100* microamps or less in each of the switch positions (110 microamps or less for the 220/240 volt power supply option).

6.2.3.5 300 microamps if external monitors are plugged into convenience receptacles.

6.2.3.6 Shut off and unplug the safety analyzer. Remove the anesthesia machine plug from the analyzer and plug it back into the original AC receptacle.

(✓) **6.3 CONFIGURATION**

6.3.1 Press the CONFIG key.

6.3.2 The CONFIGURE screen is displayed.

6.3.3 Verify the correct Time and Date.

6.3.4 Press the exit key to exit the function.

6.4 SERVICE DATA

6.4.1 Press and hold the Oxygen High Limit key and the Volume Low Limit key, and then press the UP ARROW key.

6.4.2 The Main Service Screen shall appear.

6.4.3 Select and enter the Service Log.

6.4.4 Verify any pertinent information from the Service Log. Contact the Draeger Medical, Inc. Technical Service Department if necessary.

6.4.5 Press EXIT to return to the Main Service screen.

6.4.6 Using the arrow keys, select the SRVC Service Code.

6.4.7 Press the SELECT key and enter your Technical Service Rep. I.D. number.

- (✓) 6.4.8 Press the RESET key. This resets the last service date to the current date and resets the hours run since last service to zero.

6.4.9 Press the PMS SCHED key.

- (✓) 6.4.10 Select and enter the month of the next service due date.

NOTE: The internal clock of the machine limits the amount of date advance to a maximum of six months from the current service date.

6.4.11 Press the EXIT key to return the main service screen. If not performing monitor calibrations, press the EXIT key again to return to normal operation mode.

6.5 CALIBRATIONS - One Year Service Interval; Due Date _____

6.5.1 To bring up the Oxygen Monitor Service Screen, press the Mon Cal key.

6.5.2 Remove the Oxygen sensor from the Breathing Pressure Hose Assembly, and remove the Oxygen Sensor capsule from the Oxygen Sensor Housing.

- (✓) 6.5.3 When the CURRENT CELL A and CURRENT CELL B readings have stabilized, press the ZERO key and verify that the new offset values are stored.

NOTE: The higher the offset, the higher the calculated O₂ concentration appears at high concentrations.

6.5.4 Put the Oxygen Sensor capsule into the Oxygen Sensor Housing.

6.5.5 Press the PRES MON key.

6.5.6 Remove Breathing Pressure Hose Assembly from Inspiratory Valve dome to expose the Breathing Pressure Sensor line to ambient pressure.

6.5.7 Let the Current Pressure Value stabilize and press the ZERO key to store the value.

6.5.8 Remove Breathing Pressure Sensor line from the machine.

6.5.9 Connect a Pressure Monitor Adapter, (P/N 4115043) and calibrated manometer to the Breathing Pressure quik connect on the machine monitor interface.

NOTE: Zero manometer prior to connecting to monitor interface. Slight pressure will be noticed on manometer when connected to monitor interface, DO NOT re-zero manometer.

- 6.5.10 Pressurize the circuit to 50 cm H₂O and allow the [Current Value] to stabilize.
- (✓) 6.5.11 Press the SPAN key and verify that the new span values are stored.
- 6.5.12 Release the pressure, disconnect the manometer and test fixture, reconnect the Breathing Pressure Sensor line to the quik connect on the machine monitor interface and insert the Breathing Pressure Hose Assembly into Inspiratory Valve dome.
- 6.5.13 Press EXIT to return to the Main Service screen.
- 6.5.14 Press EXIT to return to normal operation.

(✓) **6.6 ABSORBER INSPECTION**

Absorber System Inspection and Maintenance

- 6.6.1 Remove the O₂ Sensor Housing and/or Oxygen Sensor Plug and Breathing Pressure Hose Assembly from the inspiratory valve dome and examine the O-rings on each assembly. Replace O-rings as necessary.
- 6.6.2 Remove the Inspiratory and the Expiratory valve domes.
- 6.6.3 Are any pins on the valve crater damaged? Inspiratory ____ (N)
Expiratory ____ (N) Replace complete valves as necessary.
- 6.6.4 Are any pins on the valve domes damaged? Inspiratory ____ (N)
Expiratory ____ (N) Replace complete valves as necessary.
- 6.6.5 Is the valve disc in good condition? Inspiratory ____ (Y) Expiratory ____ (Y)
Replace as necessary.
- 6.6.6 Are the valve dome washers in good condition? ____ (Y) Replace washers as necessary.
- 6.6.7 Reinstall the Inspiratory and Expiratory valve domes.
- (✓) **6.6.8 Ultrasonic Flow Sensor**
 - 6.6.8.1 Remove the Ultrasonic Flow Sensor connector hose.
 - 6.6.8.2 Are the connector hose, connector, and O-ring components in good condition? ____ (Y) Replace components as necessary.
 - 6.6.8.3 Remove the Expiratory valve.
 - 6.6.8.4 Is the washer under the valve in good condition? ____ (Y) Replace washer as necessary.

- 6.6.8.5 Remove the PEEP valve.
- 6.6.8.6 Is the washer under the PEEP valve in good condition? ____ (Y)
Replace washer as necessary.
- 6.6.8.7 Reattach the PEEP valve.
- 6.6.8.8 Reattach the Expiratory valve.
- 6.6.8.9 Remove the Ultrasonic Flow Sensor from the mounting bracket
- 6.6.8.10 Remove the flow housing/transducer assembly from the
electronics housing.
- 6.6.8.11 Remove both transducers from the flow housing; examine each
O-ring and condition of all components, replace O-rings as
necessary, then reassemble the Ultrasonic Flow Sensor.
- 6.6.8.12 Reattach the Ultrasonic Flow Sensor to the mounting bracket.
- 6.6.8.13 Reattach the connector hose between the Ultrasonic Flow
Sensor and Expiratory valve.

(✓) **6.7 HIGH PRESSURE LEAK**

NOTE: Minimum cylinder pressures required for High Pressure Leak tests are:

N₂O: 600 psi;

O₂, Air: 1000 psi.

- 6.7.1 Turn the machine to Standby.
- 6.7.2 Verify the Auxiliary Oxygen flow control valve is closed
- 6.7.3 Disconnect all pipeline supply hoses at the wall outlets.
- 6.7.4 Open then close and remove each cylinder. And if applicable remove the
yoke plug from each additional yoke assembly.
- 6.7.5 Note the reading on each of the cylinder pressure gauges and start a stop
watch.
- 6.7.6 Are there two (2) yoke index pins installed securely in each yoke? ____ (Y)
- 6.7.7 Is the proper gas I.D. label affixed to each yoke? ____ (Y)
- 6.7.8 After two (2) minutes, is the pressure loss equal to or less than 50 PSI?
____ (Y)

- 6.7.9 Verify the presence of only one (1) cylinder washer, then reattach and secure the cylinders to each yoke assembly, then open each cylinder valve.

6.8 BREATHING SYSTEM

- 6.8.1 Absorber & Freshgas Leak/ Vaporizer Leak

6.8.1.1 Close all flow control valves.

6.8.1.2 Set the AUTO/BAG selector to BAG.

6.8.1.3 Close the APL valve.

6.8.1.4 Interconnect the Inspiratory valve and expiratory port on the Ultrasonic Flow Sensor with a 22 mm hose (P/N 9995132).

6.8.1.5 Attach a Freshgas Leak Test Adapter (P/N 4115041), Test Terminal Adapter (P/N 4104389), and Digital Pressure Manometer to the bag mount.

6.8.1.6 Apply 50 cm H₂O test pressure to the absorber system and start a stop watch.

6.8.1.7 Is the pressure on the manometer within 49 to 51 cm H₂O? ____ (Y)

6.8.1.8 Is the pressure displayed on the Breathing Pressure Gauge +/- 3 cm H₂O of the value displayed on the test manometer? ____ (Y)

(✓) 6.8.1.9 After thirty (30) seconds, is the pressure displayed on the manometer equal to or greater than 40 cm H₂O? ____ (Y)

6.8.1.10 If applicable, turn on the vaporizer to the first graduated marking. Reapply 50 cm H₂O of pressure to the system and start a stopwatch.

(✓) 6.8.1.11 After thirty (30) seconds, is the vaporizer test pressure displayed on the manometer equal to or greater than 40 cm H₂O? ____ (Y). Turn off the vaporizer.

(✓) 6.8.2 APL Valve

6.8.2.1 Open the APL valve to its stop.

6.8.2.2 Verify Scavenger is active.

6.8.2.3 Turn the System Power switch to ON.

- 6.8.2.4 Set the Oxygen flow to 8 l/min.
- 6.8.2.5 Is the test pressure within 0 to 3 cm H₂O? ____ (Y)
- 6.8.2.6 Remove Freshgas Leak Test Adapter (P/N 4115041), Test Terminal Adapter (P/N 4104389), and Digital Pressure Manometer (P/N PDM 200CD) from the bag mount.
- (✓) 6.8.3 O₂ Flush
 - 6.8.3.1 Attach a 33 mm x 22 Female Adapter (P/N 4115087) to the top port of the test volumeter.
 - 6.8.3.2 Remove the ventilator hose from the absorber and bellows, and connect to the Bag Mount terminal.
 - 6.8.3.3 Connect test volumeter to the end of ventilator hose.
 - 6.8.3.4 Close the APL valve.
 - 6.8.3.5 Press and hold the O₂ FLUSH button for 6 seconds; multiply the value obtained by 10.
 - 6.8.3.6 Is the calculated Oxygen flush flow rate within 45 to 65 l/min.? ____ (Y)
 - 6.8.3.7 After releasing the flush, does the flow of Oxygen stop immediately? ____ (Y)
 - 6.8.3.8 Reattach ventilator hose between absorber and bellows ventilator terminals.
 - 6.8.3.9 Remove all test equipment. Remove 22mm hose from Inspiratory valve and expiratory port on the Ultrasonic Flow Sensor.
- (✓) 6.8.4 Expiratory Valve Leak
 - 6.8.4.1 Connect a 22 mm hose (P/N 9995132) between the Inspiratory valve and the bag mount.
 - 6.8.4.2 Connect a test terminal (P/N 4104389) to the expiratory port on the Ultrasonic Flow Sensor.
 - 6.8.4.3 Connect a 0-250 ml/min. flowmeter (P/N S000081) to the test terminal.

- 6.8.4.4 Turn up the Oxygen flow until the system pressurizes to 30 cm H₂O. Adjust the APL valve as necessary to maintain 30 cm H₂O.
- 6.8.4.5 Is the value indicated on the flowmeter within 0 to 60 ml/min.? ____ (Y)
- 6.8.4.6 Turn off O₂ flow. Fully close APL valve.
- 6.8.4.7 Remove all test equipment.
- (✓) 6.8.5 Inspiratory Valve Leak
 - 6.8.5.1 Turn the System Power switch to Standby.
 - 6.8.5.2 Connect a test terminal (P/N 4104389) to the Inspiratory valve.
 - 6.8.5.3 Connect a Fresh Gas Leak Adapter (P/N 4115041) and manometer to the test terminal on the Inspiratory valve.
 - 6.8.5.4 Connect another test terminal to the bag mount.
 - 6.8.5.5 Connect a 0-250 ml/min. flowmeter (S000081) to the test terminal on the bag mount.
 - 6.8.5.6 Pressurize the test circuit to 30 cm H₂O.
 - 6.8.5.7 Is the value indicated on the flowmeter within 0 to 60 ml/min.? ____ (Y)
 - 6.8.5.8 Turn the system power switch to ON.
 - 6.8.5.9 Remove all test equipment.
- (✓) 6.8.6 Absorber PEEP Valve
 - 6.8.6.1 Open the APL valve.
 - 6.8.6.2 Disconnect the pressure line going to the Breathing Pressure gauge and connect to a test manometer.
 - 6.8.6.3 Connect a Fresh Gas Outlet Volume Test Device (S010158) to the Fresh gas outlet.
 - 6.8.6.4 Install a test breathing circuit or disposable breathing circuit onto the Inspiratory and Expiratory valves of the absorber.
 - 6.8.6.5 Connect the test terminal, of the Fresh Gas Outlet Volume Test Device, to the Y-piece of the breathing circuit.

- 6.8.6.6 Set the O₂ flow to 5 l/min.
- 6.8.6.7 Adjust the absorber PEEP valve clockwise to the maximum position.
- 6.8.6.8 Does the PEEP valve adjust smoothly? ____ (Y)
- 6.8.6.9 Is the maximum PEEP indicated on the test gauge greater than 15 cm H₂O? ____ (Y)
- 6.8.6.10 Adjust the absorber PEEP valve counterclockwise to its minimum position.
- 6.8.6.11 Is the PEEP indicated on the test gauge less than or equal to 3 cm H₂O? ____ (Y)
- 6.8.6.12 Close the O₂ flow control valve.
- 6.8.6.13 Remove the test equipment and reconnect the pressure line to the Breathing Pressure gauge.

6.9 OXYGEN ANALYZER

- 6.9.1 Press the Cal key.
 - NOTE:** Ensure that the sensor has stabilized in ambient air for several minutes.
- (✓) 6.9.2 After calibration is completed, is the O₂ concentration 21%? ____ (Y)
- 6.9.3 The warning message INSP O₂ LOW shall appear on the central alarm display, and an audible alarm shall sound.
- 6.9.4 Press the ALARM SILENCE key and verify the audio alarm is silenced.
- 6.9.5 Remove the Breathing Pressure Hose Assembly from the Inspiratory valve dome.
- 6.9.6 Place the Oxygen Sensor housing directly into the Inspiratory valve dome.
- 6.9.7 Verify the AUTO/BAG selector is set to BAG.
- 6.9.8 Close the APL valve.
- 6.9.9 Attach a 22 mm hose (P/N 9995132) to the Inspiratory valve.
- 6.9.10 Attach a Breathing System Leak Test Device (P/N S010159) to the bag mount.

- 6.9.11 Press the O₂ Flush.
- 6.9.12 After 20 seconds, is the O₂ concentration within 90 to 100 %? ____ (Y)
- 6.9.13 Release the O₂ Flush, does the flow cease immediately? ____ (Y)
- 6.9.14 Set the Oxygen flow to 8 l/min.
- (✓) 6.9.15 After 1 minute, is the O₂ concentration within 97 to 100 %? ____ (Y)

6.10 FLOWMETERS & CONCENTRATIONS

- (✓) 6.10.1 Oxygen Flowmeter
 - 6.10.1.1 Is it possible to adjust the flow of Oxygen over the full range of the flowmeters? ____ (Y)
 - 6.10.1.2 Set the Oxygen flow to 4 l/min.
 - 6.10.1.3 Is the correct flow control knob and label attached to the Oxygen flow control valve? ____ (Y)
- (✓) 6.10.2 Air Flowmeter
 - 6.10.2.1 If not configured with an Air Cylinder yoke, attach the Air Pipeline hose.
 - 6.10.2.2 Is it possible to adjust the flow of the Air over the full range of the flowmeter? ____ (Y)
 - 6.10.2.3 Set the Air flow to 2 l/min.
 - 6.10.2.4 After the value stabilizes, is the O₂ concentration within 71 to 77 %? ____ (Y)
 - 6.10.2.5 Close the Air flow control valve.
 - 6.10.2.6 Is the correct flow control knob and label attached to the Air flow control valve? ____ (Y)
- 6.10.3 Nitrous Oxide Flowmeter
 - 6.10.3.1 Set the Nitrous Oxide flow to 2 l/min.
 - 6.10.3.2 After the value stabilizes, is the O₂ concentration within 64 to 70 %? ____ (Y)
 - 6.10.3.3 Is the correct flow control knob and label attached to the Nitrous Oxide flow control valve? ____ (Y)

6.10.3.4 Is it possible to adjust the flow of nitrous oxide over the full range of the flowmeter? ____ (Y)

(✓) 6.10.4 Oxygen Ratio Control

6.10.4.1 Open the Nitrous Oxide flow control valve to the stop position.

6.10.4.2 After the value stabilizes, is the O₂ concentration within 21 to 29 %? ____ (Y)

6.10.4.3 Set the Oxygen flow to 2 l/min.

6.10.4.4 After the value stabilizes, is the O₂ concentration within 21 to 29 %? ____ (Y)

6.10.4.5 Set the Oxygen flow to 1 l/min.

6.10.4.6 After the value stabilizes, is the O₂ concentration within 21 to 29 %? ____ (Y)

6.10.4.7 Close the Oxygen flow control valve.

6.10.4.8 Close the Nitrous Oxide flow control valve.

(✓) 6.10.5 Auxiliary Oxygen Flowmeter

6.10.5.1 Connect a test manometer to the Auxiliary Oxygen outlet using a PDM/Suction adapter (P/N 4115038).

6.10.5.2 Slowly open flowmeter valve until a pressure of 50 cm H₂O is achieved.

6.10.5.3 Close flowmeter valve.

6.10.5.4 After 10 seconds, is the pressure within 40 to 60 cm H₂O? ____ (Y)

6.10.5.5 Remove the test manometer and adapter.

6.10.5.6 Is it possible to adjust the Auxiliary Oxygen flowmeter over the full range of the flowmeter? ____ (Y)

6.10.5.7 Set the Auxiliary O₂ flow rate to 5 L/min.

6.10.5.8 Hold the Oxygen sensor at the Auxiliary Oxygen flowmeter outlet.

6.10.5.9 After the value stabilizes, is the O₂ concentration within 80 to 100 %? ____ (Y)

6.10.5.10 Reinstall the Breathing Pressure Hose Assembly.

6.10.5.11 Insert the Oxygen sensor into the Breathing Pressure Hose Assembly.

6.10.5.12 Close the Auxiliary Oxygen flow control valve.

6.11 HIGH PRESSURE REGULATORS – Six Month Service Interval; Due Date _____

(✓) 6.11.1 N2O Regulator- if applicable

NOTE: Minimum cylinder pressure for N2O regulator test is 600 psi.

6.11.1.1 Configure the test gauge (P/N 4114807) using a N2O nut/stem DISS connector (P/N 4114830-004) on the hose, and N2O DISS body connector (P/N 4114830-003) on the valve body side. If the machine is configured with CSA style fittings reverse the position of the connectors.

6.11.1.2 Connect the test fixture hose to the machine's Nitrous Oxide pipeline inlet.

6.11.1.3 Connect the Nitrous Oxide pipeline supply hose to the test fixture.

6.11.1.4 Open the Nitrous Oxide and the Oxygen cylinder valves.

6.11.1.5 Set the Oxygen and Nitrous Oxide flows to 4 l/min.

6.11.1.6 Depress the push button on the test device.

6.11.1.7 Release the push button. After the pressure decay stabilizes, is the regulator output pressure within 40 to 49 psi? ____ (Y)

6.11.1.8 Remove the test fixture.

NOTE: If a pressure decrease does not occur, either the hospital's supply pressure is too low or the regulator pressure is set too high.

(✓) 6.11.2 O2 Regulator

NOTE: Minimum cylinder pressure for O2 regulator test is 1000 psi.

6.11.2.1 Configure a test gauge (P/N 4114807) using an O2 nut/stem DISS connector (P/N 4114830-006) on the hose and an O2 DISS body connector (P/N 4114830-005) on the valve body side. If the machine is configured with CSA style fittings reverse the position of the connectors.

6.11.2.2 Connect the test fixture hose to the machine's Oxygen pipeline inlet.

6.11.2.3 Connect the Oxygen pipeline supply hose to the test fixture.

6.11.2.4 Set the Oxygen flow to 4 l/min.

6.11.2.5 Depress the push button on the test device.

6.11.2.6 Release the push button. After the pressure decay stabilizes, is the regulator output pressure within tolerance given in the following table? ____ (Y)

NOTE: If a pressure decrease does not occur, either the hospital's supply pressure is too low or the regulator pressure is set too high.

Cylinder Pressure psi	USA Compensated Regulator output tolerances	CSA Compensated Regula- tor output tolerances
2000	38 to 44	41 to 47
1800	39 to 45	42 to 48
1600	40 to 46	43 to 49
1400	41 to 47	44 to 50
1200	42 to 48	45 to 51
1000	43 to 49	46 to 52

(✓) **6.12 LOW O₂ SUPPLY – Six Month Service Interval; Due Date _____**

6.12.0.1 Close the Oxygen Cylinder valve and drain Oxygen cylinder pressure.

6.12.0.2 Depress the push button on the test device.

6.12.0.3 Adjust the Oxygen flow to 500 ml/min.

6.12.0.4 Release the test device push button.

6.12.0.5 Is the pressure on the test gauge when the LO O₂ SUPPLY message appears on the monitor within 34 to 40 psi? ____ (Y)

6.12.0.6 Close O₂ Flow control valve.

6.12.0.7 Remove the test equipment.

6.13 OXYGEN SUPPLY FAILURE PROTECTION

- 6.13.1 Connect all pipeline supplies.
- (✓) 6.13.2 *Is the flow of Oxygen 150 to 200 ml/min.? ____ (Y)
- 6.13.3 Open the Nitrous Oxide flow control valve.
- (✓) 6.13.4 *Is the flow of Nitrous Oxide within 375-750 ml/min.? ____ (Y)
- 6.13.5 Adjust the Oxygen, Nitrous Oxide and AIR gases to 4 l/min.
- 6.13.6 Disconnect the Oxygen pipeline supply and close the Oxygen cylinder valve.
- (✓) 6.13.7 Do all fresh gas flows cease when the Oxygen pressure is depleted? ____ (Y)
- 6.13.8 Reconnect the Oxygen pipeline supply.
- 6.13.9 Close all cylinder valves and then disconnect the Nitrous Oxide pipeline supply.
- 6.13.10 Drain the cylinder contents then reconnect the pipeline supplies.
- 6.13.11 Close all flow control valves.

* Nitrous Oxide Bypass flow and Minimum Oxygen flow specifications are given @ 50 psi. Pipeline pressure deviations may affect these tests.

6.14 PRESSURE MONITOR

- 6.14.1 Disconnect the Breathing Pressure Sensor line from the Breathing Pressure interface on the machine.
- 6.14.2 Connect a PDM to Monitor Adapter (P/N 4115043) and test manometer to the Breathing Pressure interface on the machine.
- 6.14.3 Adjust the test pressure to 0 cm H₂O.
- 6.14.4 Simultaneously turn "ON" the ventilator and start a stopwatch.
- (✓) 6.14.5 Does the APNEA PRESSURE alarm appear on the alarm display as a CAUTION within 13 to 17 seconds? ____ (Y)
- (✓) 6.14.6 Increase the test pressure slowly. Does the APNEA PRESSURE alarm deactivate within 10 to 14 cm H₂O? ____ (Y)
- 6.14.7 First decrease the pressure then increase the test pressure above the threshold line shown on the display, and begin timing with a stopwatch.

- (✓) 6.14.8 Does the CONTINUOUS PRES alarm appear as a warning within 13 to 17 seconds? ____ (Y)
- (✓) 6.14.9 Decrease the pressure slowly. Does the CONTINUOUS PRES alarm deactivate within 10 to 14 cm H₂O? ____ (Y)
- (✓) 6.14.10 Increase the test pressure slowly. Does a VENT PRESS HI alarm activate as a warning alarm within 47 to 53 cm H₂O? ____ (Y)
- (✓) 6.14.11 Increase the sub-atmospheric test pressure slowly. Does the SUB ATM PRES warning alarm activate within -7 to -13 cm H₂O? ____ (Y)
- 6.14.12 Turn the ventilator "OFF".
- 6.14.13 Open APL valve.
- 6.14.14 Remove the test equipment and reconnect the Breathing Pressure Sensor line to the Breathing Pressure interface on the machine.

6.15 VENTILATOR

NOTE: Readjustment of inspiratory flow to limit the inspiratory plateau may be required to reduce erratic tidal volumes and breath rates caused by artifact volumes.

- 6.15.1 Remove the bellows vent hose and the scavenger hose at the ventilator relief valve. Remove the bellows sub-assembly and remove bellows.
- 6.15.2 Visually inspect the bellows for deterioration particularly at it's seams and corrugations. Replace as necessary.
- 6.15.3 Verify the presence and condition of it's sealing O-ring. Replace as necessary. Reassemble all components.
- 6.15.4 Set O₂/AIR mode switch to "O₂".
- 6.15.5 Turn the ventilator "ON".
- 6.15.6 Set the FREQUENCY to 10 BPM.
- (✓) 6.15.7 Press and hold the EXTENDED RANGE switch and set the I:E ratio to 2:1. Using a stopwatch, time the extended I:E ratio. Is the inspiratory time within 3.6 to 4.4 seconds and the expiratory time within 1.8 to 2.2 seconds? ____ (Y)
- (✓) 6.15.8 Set the I:E RATIO to 1:2. Using a stopwatch, time the I:E ratio. Is the inspiratory time within 1.8 to 2.2 seconds and the expiratory time within 3.6 to 4.4 seconds? ____ (Y)

- 6.15.9 Adjust the Oxygen flow to 500 ml/min.
- 6.15.10 Set the Tidal Volume to 1200.
- 6.15.11 Attach a patient circuit to the absorber system.
- 6.15.12 Set the pressure limit control to MAX.
- 6.15.13 Adjust the Inspiratory Flow to the bottom of the LOW zone.
- 6.15.14 Occlude the Y-piece.
- 6.15.15 Press the O₂ Flush momentarily to inflate the bellows.
- 6.15.16 Adjust the Inspiratory Flow until a peak pressure of 80 cm H₂O is achieved.
- (✓) 6.15.17 Set the Pressure Limit Control to within the 30 range. Readjust within the band as necessary to achieve proper value, is the peak pressure at the 30 range within 27 to 33 cm H₂O? ____ (Y)
- (✓) 6.15.18 Set the Pressure Limit Control to the MIN position. Is the peak pressure at the MIN range within 0 to 15 cm H₂O? ____ (Y) Return the Pressure Limit control to MAX.
- 6.15.19 Disconnect the hose between the expiratory valve and the ultrasonic flow sensor and blow into it.
- (✓) 6.15.20 Does the Reverse Flow message appear on the display? ____ (Y)
- 6.15.21 Reconnect the hose between the expiratory valve and the flow sensor.
- 6.15.22 Insert a test minute volumeter between the exhalation valve and PEEP valve.
- 6.15.23 Open the Y-piece.
- 6.15.24 Turn the ventilator off and press the ALL STBY key to clear alarms.
- 6.15.25 Turn the ventilator “ON” and start a stop watch.
- 6.15.26 Does APNEA-VOLUME appear as a Caution within 13 to 17 seconds? ____ (Y)
- 6.15.27 Attach a 3 liter breathing bag to the Y-piece.

NOTE: Bag should be placed on a flat horizontal surface to reduce artifact volume.

- 6.15.28 Press the O₂ Flush momentarily to inflate the bellows.
- 6.15.29 Set the Inspiratory Flow to the MED and readjust as necessary to fully collapse the bellows.
- (✓) 6.15.30 Observe the operation of each unidirectional valve disc at eye level and make sure the inspiratory valve disc raises only during the inspiration phase, and the expiratory valve raises only during the exhalation phase.
- (✓) 6.15.31 Is the tidal volume on the volume monitor and on the test volumeter within 20 % of each other? ____ (Y)
- 6.15.32 Does the volume monitor display 10 BPM? ____ (Y)
- 6.15.33 Does the display correctly track the Breathing Pressure waveform? ____ (Y)
- 6.15.34 Verify AIR pipeline hose is connected.
- (✓) 6.15.35 Set O₂/AIR mode switch to "AIR". Does ventilator continue to cycle properly. ____ (Y)
- 6.15.36 Set O₂/AIR mode switch back to "O₂".

6.16 BELLOWS ADULT

- (✓) 6.16.1 What is the tidal volume indicated on the test volumeter within 960 to 1440 ml? ____ (Y)
- (✓) 6.16.2 Does the bellows remain fully inflated during the expiratory pause phase? ____ (Y)
- 6.16.3 Remove the ventilator hose from the bellows hose terminal.
- 6.16.4 Attach a test terminal to the bellows hose terminal.
- 6.16.5 Connect a 0-250 ml/min. flowmeter (P/N S000081) to the test terminal.
- 6.16.6 Set the FREQUENCY to 1 BPM.
- (✓) 6.16.7 Is the drive gas leakage indicated during the inspiratory phase within 0 to 50 ml/min.? ____ (Y)
- 6.16.8 Remove the test equipment from the ventilator hose terminal and reconnect the ventilator hose to the bellows ventilator hose terminal.
- 6.16.9 Set the FREQUENCY to 10 BPM.
- 6.16.10 Adjust the O₂ flow to 8 l/min.

- 6.16.11 Adjust the Tidal Volume to maximum.
- 6.16.12 Press the O₂ Flush momentarily to inflate the bellows.
- 6.16.13 Adjust the INSPIRATORY FLOW to fully compress the bellows.
- (✓) 6.16.14 Is the Tidal Volume on the test volumeter greater than or equal to 1400 ml? ____ (Y)
- (✓) 6.16.15 Is the PEEP value displayed on the monitor within 0 to 3 cm H₂O? ____ (Y)
- 6.16.16 Remove the breathing bag from the Y-piece and replace it with a test lung.
- 6.16.17 Adjust the O₂ flow to 300 ml/min.
- 6.16.18 Adjust the Tidal Volume to 200 ml.
- 6.16.19 Does the bellows stop adjust smoothly and engage properly? ____ (Y)
- 6.16.20 Adjust the INSPIRATORY FLOW to fully compress the bellows.
- (✓) 6.16.21 Is the Tidal Volume on the test volumeter within 125 to 250 ml? ____ (Y)
- 6.16.22 Close the Oxygen flow control valve.
- 6.16.23 Remove the test lung and set the AUTO/BAG selector to BAG.
- 6.16.24 Press the Alarms All Standby key.

6.17 SCAVENGER- Six Month Service Interval; Due Date _____

NOTE: If the ambient air in the local environment contains a significant amount of dust and lint, the cleaning frequency must be increased to compensate for these conditions.

- (✓) 6.17.1 Scavenger - Passive Mode
 - 6.17.1.1 Remove the scavenger hoses and drain all accumulated moisture. Inspect all scavenger hoses for deterioration and replace as necessary. Reinstall hoses and bag.
 - 6.17.1.2 Ensure that the suction needle valve (not used in passive mode) is closed.

6.17.1.3 Positive Pressure Test:

6.17.1.3.1 Connect a 19mm scavenger hose between the Vent Relief valve and the top port of the scavenger. Connect another 19 mm scavenger hose between APL valve and the second scavenger port.

6.17.1.3.2 Interconnect the inspiratory valve and expiratory valve or expiratory port on the ultrasonic flow sensor, if applicable with a 22 mm hose (P/N 9995132).

6.17.1.3.3 Set Man/Auto valve to "Auto".

6.17.1.3.4 Verify PEEP valve is at minimum (counter-clockwise).

6.17.1.3.5 Set O₂ Flow to 8 l/min.

6.17.1.3.6 Occlude the bottom 19 mm connector (EXHAUST).

6.17.1.3.7 After the bellows inflates, the flow of oxygen will exit the system through the positive pressure safety relief valve.

6.17.1.3.8 Is the pressure displayed on absorber cm H₂O gauge is less than or equal to 5 cm H₂O?___(Y)

(✓) 6.17.2 Scavenger - Suction Mode

6.17.2.1 Remove the scavenger hoses and drain all accumulated moisture. Inspect all scavenger hoses and scavenger bag for deterioration and replace as necessary. Reinstall hoses and bag.

6.17.2.2 Remove the filter from the negative pressure relief valve. Clean and reinstall the filter.

6.17.2.3 Negative Pressure Test:

6.17.2.3.1 Interconnect the inspiratory valve and expiratory valve or expiratory port on the ultrasonic flow sensor, if applicable with a 22 mm hose (P/N 9995132).

6.17.2.3.2 Set Man/Auto valve to "BAG" position

6.17.2.3.3 Fully open APL valve.

6.17.2.3.4 Attach a Breathing System Leak Test Device (P/N S010159) to the bag mount (occlude bag mount).

6.17.2.3.5 Verify Suction system is active. Adjust needle valve to allow typical suction through scavenger.

6.17.2.3.6 Verify all flow control valves are closed.

6.17.2.3.7 Install scavenger adapter (P/N 4108114) between either of the 19 mm scavenger connectors and scavenger hose. Connect a manometer to the hose barb on the scavenger adapter and observe reading on manometer. The manometer shall indicate a pressure less than or equal to 0.5 cm H₂O.

6.17.2.4 Positive Pressure Test:

6.17.2.4.1 Fully close scavenger needle valve.

6.17.2.4.2 Push the O₂ flush to inflate scavenger bag. Open oxygen flow control valve to 8 l/min.

6.17.2.4.3 Verify reading on manometer is less than or equal to 5 cm H₂O.

6.17.2.4.4 Remove all test equipment. Adjust scavenger needle valve to allow typical suction through scavenger.

6.18 FINAL TESTS

- (✓) 6.18.1 Is the machine's Operator's Instruction manual in close proximity of the machine? ____ (Y)
- 6.18.2 Verify all cylinder pressure gauges indicate zero.
- 6.18.3 Verify the pipeline hoses are connected to the hospital pipeline.
- 6.18.4 Verify the APL valve knob is turned completely counterclockwise (fully open).
- 6.18.5 Place the AUTO/BAG selector in the BAG position.
- 6.18.6 Verify the ventilator hose is connected between the AUTO/BAG valve and bellows ventilator hose terminal.
- 6.18.7 Verify the Breathing Pressure Sensor line is connected to the Breathing Pressure interface on the machine.
- 6.18.8 Verify the Oxygen sensor is removed from the inspiratory valve dome.
- 6.18.9 Verify that the Inspiratory valve dome is plugged.
- 6.18.10 Verify that the machine is plugged into a live A/C receptacle.

7.0 Software Update Procedure

This section outlines the software installation procedure, including the equipment needed and its connections.

Software updates to the Narkomed Mobile anesthesia system are done through a serial port connection to an external PC using the batch file LOADM.BAT.

7.1 Software Transfer to PC Via Modem

Equipment required:

- Interface Cable, NAD Part No.4109882 P (9-pin to 25-pin)
or 4110328 A (9-pin to 9-pin)
- IBM[®] PC or IBM PC Compatible configured with:
 - PC-DOS or MS-DOS V3.3 or higher
 - RS-232C Serial Port connected to COM 1
 - Hard Drive or Floppy Drive
 - Modem (or external modem)

7.1.1 Download the software to the hard disk or use Drive A (floppy drive) on the PC.

7.2 Installing Narkomed Mobile Software from a PC

7.2.1 Set the System Power switch on the Narkomed Mobile to STANDBY, and the power switch on the PC to OFF.

7.2.2 Connect the appropriate interface cable (9-pin or 25-pin) to COM 1 on the PC, and connect the other end of the cable to the Narkomed Mobile serial interface port as shown in Figure 7-1.

7.2.3 Power up the PC and wait for the DOS prompt to appear on the screen.

7.2.4 Set the PC to read the drive holding the software. For example: if the software was downloaded to drive A, type A: and press ENTER.




Figure 7-1. Narkomed Mobile Serial Port Location

7.2.5 Type LOADNMM and press ENTER.

7.2.6 Turn the System Power switch to ON.

7.2.7 As the software is downloading, the Narkomed Mobile screen will be blank for approximately 5 to 7 minutes, and the incremental number of bytes sent will be displayed on the PC screen. When the download is complete, the PC screen will display

```
\      Bytes sent:    0      Images sent:    1
```

Software installation is complete when the machine resets.


7.2.8 Set the System Power switch on the machine to STANDBY, and the power switch on the PC to OFF. Disconnect the interface cable.

8.0 Spare and Replacement Parts

Part numbers for field-replaceable items on the Narkomed Mobile anesthesia system are listed on the following pages, along with part numbers for related hardware and cables.

The item numbers are keyed to the accompanying illustrations to aid in identifying the item and its location.

ASSEMBLY/PART	PAGE
Vaporizer & rel. parts, fresh gas outlet, mainswitch, flowmtr sub-asm	8-2, 8-3
Processor Assembly and related parts	8-4, 8-5
Ventilator Controller (Bezel Assembly).	8-6, 8-7
Bellows Valve Assembly.	8-8, 8-9
Bellows Valve Assembly Details	8-9A, 8-9B
Display Assembly and Monitor Support Arm.	8-10, 8-11
Pipeline Inlet (Manifold) Assemblies	8-12, 8-13
Failsafe (OFPD) and ORC Assemblies	8-14, 8-15
Main Switch Assembly.	8-16, 8-17
O ₂ Supply Pressure Switch, O ₂ - Air Switch	8-18, 8-19
Flowmeter Shield, Gauges.	8-20, 8-21
Flow Tubes, Restrictor Assemblies, Flow Control Valves	8-22, 8-23
Auxiliary O ₂ Flowmeter.	8-24, 8-25
Cylinder Yokes, Regulators, O ₂ Flush Valve.	8-26, 8-27
Frame, Casters	8-28, 8-29
Power Supply Assembly and related items.	8-30, 8-31
Absorber, Inspiratory Valve, Ultrasonic Flow Sensor	8-32, 8-33
Valve, Man/Auto Selector.	8-34, 8-35
Breathing Pressure Hose Assembly, O ₂ Sensor, PEEP Valve, Expiratory valve . . .	8-36, 8-37
Scavenger, Touch-up paint: Euro white, Euro Blue	8-38, 8-39
Miscellaneous Items.	8-40




SV00344

NM MOBILE	SPARE AND REPLACEMENT PARTS (continued)
------------------	--

ITEM	DESCRIPTION	PART NUMBER
1	Deleted	
2	O-ring, vaporizer post (2x)	4115864
3	Mainswitch Assembly.	4114278
4	O ₂ Flush valve	4103340
5	Fresh Gas outlet assembly.	4108673
6	Vapor block assembly	4115019-004
7	Flowmeter sub-assembly	4114276-001
8	Gasket	1101690-001
9	Rail, remote display	4114136

SPARE AND REPLACEMENT PARTS (continued)


NM MOBILE



NM MOBILE	SPARE AND REPLACEMENT PARTS (continued)
------------------	--

ITEM	DESCRIPTION	PART NUMBER
1	Ventilation Monitor Asm, NMM	4113780-001
	Ventilation Monitor Asm, NMM (Service Exchange.	SE4113780-001
2	Screw, 4-40 x 1.62 in. cap skt hd (4x)	HW01102
3	Lock Washer, #4 int-t (4x)	HW67001
4	Deleted	
5	Screw, 8-32 x 3/8 in. btn hd skt (4x).	HW09018
6	Lock Washer, #8 int-t (4x)	HW67011
7	Hose, 0.13 ID, 8 in.	ML08007
	Fitting, 0.13 hose x 1/8 MPT	4102963
	Nut, Panel 9/16 - 18	4108156
	Quick-disconnect Fitting	4108155
8	O ₂ Interface Cable Assembly	4113909-001
9	Volume Interface Cable Assembly	4113910-001
10	Cable Assembly, Display Interface	4114288
11	Screw, 6-32 x 1/4 in. btn hd skt (2x)	HW09076
12	Lock Washer, #6 split (2x)	HW65001
13	Speaker Mounting Bracket	4113285
	Speaker Assembly (includes wire harness and connector).	4114366
14	Cable Assembly, Serial Port.	4113760-001

Note: Item #1 could be made up of two PCBs, 4113549 and 4113595; or a single PCB, 4114964. They are interchangeable. 4114964 is the newer design.



NM MOBILE	SPARE AND REPLACEMENT PARTS (continued)
------------------	--

ITEM	DESCRIPTION	PART NUMBER
1	Ventilator Controller Assembly (Bezel Assembly)	4114283
	Ventilator Controller Assembly (Bezel Assembly) (Service Exchange)	SE4114283
2	Solenoid	4110906
3	Valve, air piloted	4114047
4	Regulator	4114252
5	Collet (3x)	4112167
6	Knob (3x)	4113281
7	Washer, #10 flat (3x)	HW66003
8	Nut, hex M5 x 0.5 (3x)	4112066
9	Cover (3x)	4113278-002
10	Gauge	4112251-001
	Lens	4112213
	Regulator rebuild kit	4118241

Mounting hardware for controller assembly:


Base plate mounting:

Screw, 8-32 x $\frac{3}{8}$ in. cap skt hd (2x)	HW01012
Lock Washer, #8 split (2x)	HW65011
Flat Washer, #8 (2x)	HW66002

Rear of hex standoffs:

Screw, 6-32 x $\frac{1}{2}$ cap skt hd (2x)	HW01010
Lock washer, #6 int-t (2x)	HW67007
Flat Washer, #6 (2x)	HW66006

SV00334



NM MOBILE


SPARE AND REPLACEMENT PARTS (continued)

ITEM	DESCRIPTION	PART NUMBER
1	Bellows Valve Assembly	4112272-001
	Bellows Valve Assembly (Service Exchange)	SE4112272-001
2	Canister	4106948
3	Collet (3x)	4112167
4	Knob, PLC Adj.	4113279
5	Flat Washer, #10 (3x)	HW66003
6	Hex Nut, M5 x 0.5 (3x)	4112066
7	Knob Cover (2x)	4113278-002
8	Knob, Tidal Vol Adj.	4113280
9	Housing	4114011
10	Housing w/push to turn asm	
11	Knob, Off-On Switch	4113281
12	Bellows Assembly	4114102
13	Relief Valve Assembly	4108050
	Diaphragm Assembly	4110960
14	Urethane (Non-Latex) Bellows Sub-assembly, Adult	4106930-001
	O-ring #217 (neoprene)	4101817
14A	Gasket, Bottom	4105489
15	O-ring #256, canister gasket	4107018
16	Hose barb fitting, 1/16 ID x 10-32 w/seal	4112707-001
17	Knob Cover	4113278-001

SPARE AND REPLACEMENT PARTS (continued)

NM MOBILE

SV00555_a




NM MOBILE	SPARE AND REPLACEMENT PARTS (continued)
------------------	--

ITEM	DESCRIPTION	PART NUMBER
------	-------------	-------------

Bellows Valve Assembly details:

1	Press-on Hose Clamp (2x)	4104161
2	Label, O ₂ Tubing (2x)	4109871
3	Hose, 0.075 I.D..	ML08003
4	Restrictor	4107639
5	Bellows Top Guide	4110735
6	Volume Indicator	4108276
7	Set Screw, 6-32 x ¼ in. cup point (2x)	HW04003
8	Silencer, pneumatic, 3/8 NPT	4109819
9	O-ring, #19 silicone.	4107096
10	O-ring, #22 neoprene	4104595
11	Rod, bellows adjust.	4110727
12	Rod.	4106418
13	Lock washer, #8 split	HW65001
14	Lock washer, #8 int-t	HW67000
15	Screw, 8-32 x 3/8 in. btn hd skt	HW09008




NM MOBILE	SPARE AND REPLACEMENT PARTS (continued)
------------------	--

ITEM	DESCRIPTION	PART NUMBER
1	Display Assembly (early design with thumbscrew retainers on connector)	4114249
1	Later design with slide latch retainer on connector	4115756
	(display assemblies include mounting rod)	
	Kit to convert display asm and machine to slide-latch connector arrangement. . .	4116071
2	Screw, 6-32 x 5/8 in., skt hd	HW01067
3	Joint Assembly and Monitor Support Arm	4114286

SPARE AND REPLACEMENT PARTS (continued)

NM MOBILE



SV00336_a

NM MOBILE	SPARE AND REPLACEMENT PARTS (continued)
------------------	--

ITEM	DESCRIPTION	PART NUMBER
------	-------------	-------------

Pipeline Inlet Assemblies:


1	N ₂ O Inlet Manifold (incl. fittings and filter)	4115047
2	Air Inlet Manifold (incl. fittings and filter)	4115045
3	O ₂ Inlet Manifold (incl. fittings and filter)	4115046
4	Hose Barb Fitting, 1/16 ID x 10-32 w/ seal	4112707-001

Typ., all inlet assemblies:

5	Filter and Connector (assembled, incl. O-ring)	4114345
6	O-ring (neoprene)	4112619-009
7	Plug	4114328


	O-ring for N ₂ O DISS Nipple	4113494
--	---	---------

8	Body, DISS N ₂ O x 1/8 MPT	4111384
9	Body, DISS Air x 1/8 MPT	4102886
10	Body, DISS O ₂ x 1/8 MPT	4102563



NM MOBILE	SPARE AND REPLACEMENT PARTS (continued)
------------------	--


ITEM	DESCRIPTION	PART NUMBER
1	Lock Washer, #8 split (3x)	HW65011
2	Screw, 8-32 x 1 $\frac{3}{8}$ in. cap skt hd (3x).....	HW01110
3	OFPD (Failsafe Assembly) (2x)	4114031
4	Screw, 8-32 x 3 in. cap skt hd (3x)	HW01090
5	Lock Washer, #8 split (3x)	HW65011
6	O-ring, #105 (neoprene)	4111893
7	ORC Assembly	4114277-001
	ORC Assembly (Service Exchange).....	SE4114277-001
	Not shown:	
	Filter	4111805
	O-ring, 0.166 x 0.042 (Buna-N)	4111894



NM MOBILE

SPARE AND REPLACEMENT PARTS (continued)

ITEM	DESCRIPTION	PART NUMBER
1	Main Switch Assembly.	4114278
2	Switch Assembly (incl. leaf switch, wire harness & connector).	4114318
3	Label, Dot, plain white.	4103423
4	Knob, main switch, NMM	4114311
5	Spring, CPRSN, 0.709 OD x 1 1/2 L	4110975-013
6	Housing, main switch, NMM.	4106045-011
7	Cam, main switch actuator	4106043



NM MOBILE	SPARE AND REPLACEMENT PARTS (continued)
------------------	--


ITEM	DESCRIPTION	PART NUMBER
------	-------------	-------------

O₂ Supply Pressure Switch:


1	Screw, 2-56 x 5/8 in. btn hd skt (4x).....	HW09085
2	Lock Washer, #2 int-t (4x)	HW67012
3	Flat Washer, #2 (4x)	HW66009
4	Switch assembly (incl. wire harness)	4114331
5	Fitting, straight 1/16 ID x 10-32m w/seal.....	4112707-001

O₂ - Air Switch:

6	Switch (Valve, 3-way).....	4115136
---	----------------------------	---------



SV00340_a




*R = Loctite #271 (Red)

NM MOBILE	SPARE AND REPLACEMENT PARTS (continued)
------------------	--

ITEM	DESCRIPTION	PART NUMBER
1	Shield, Flowmeter, 3-Gas (Air) NMM	4115124
1A	Flowmeter shield screws	HW09000
2	Gauge, 100 psi NMM (3x)	4114247-001
3	Hose barb ftg, 1/16 ID x 10-32 w/seal (3x) (pipeline gauges).....	4112707-001
4	Kep Nut, (2x per gauge).....	HW55003
5	Straight Fitting (cylinder gauges).....	4109402
6	Gauge, 3000 psi NMM (2x)	4114247-002



SPARE AND REPLACEMENT PARTS (continued)

NM MOBILE



NM MOBILE	SPARE AND REPLACEMENT PARTS (continued)
------------------	--

ITEM	DESCRIPTION	PART NUMBER
1	Flow Control Valve (3x)	4114036
2	Stop Pin Nut (3x)	4111892
3	Knob Guard, 3 Gas	4110574
4	Knob, O ₂	4103156
5	Label, O ₂ Flow Control Knob, Green (USA)	4103178
6	Screw, 6-32 x 7/16 btn hd (2x)	HW09017
7	Label, Air Flow Control Knob, Yellow (USA, Germany)	4103905
8	Label, N ₂ O Flow Control Knob, Blue (USA, UK, Canada)	4103904
9	Knob, Flow (2x)	4103736
10	O-ring, #010 (neoprene) (2x)	4101872
11	Restrictor Housing (2x)	4103440
12	Restrictor, N ₂ O, black.	4110738-005
13	Gasket and Guide Ring (2x each flow tube)	4102724
14	Flow Tube, N ₂ O, 0.1 - 8L	4114263
15	Connector (Flow tube retainer) (3x)	4114017
16	O-ring, #109 (EPDM) (3x)	4112628-001
17	Flow Tube, O ₂ , 0.1 - 8L	4114262
18	Flow Tube, Air, 0.1 - 8L	4114264
19	Restrictor, O ₂ , red.	4110738-003




NM MOBILE	SPARE AND REPLACEMENT PARTS (continued)
------------------	--

ITEM	DESCRIPTION	PART NUMBER
1	Auxiliary O ₂ Flow Meter Assembly, new style	4109310
	Auxiliary O ₂ Flow Meter Assembly, new style (Service Exchange).	SE4109310
	Kit, Flowmeter Auxiliary O ₂ (new installation, includes hardware)	S4109310
2	Screw, 10-32 x 1 1/16 in. btn hd (2x).	HW09043
3	Housing, flowmeter	4111053
4	Screw, set, cup pt, 6-32 x 1/4 in.	HW04003
5	Knob	4111442
6	Nut, kep, 10-32	HW55002
7	Flowmeter, incl. tube & valve	4111460
8	Label, 5/8 w/dot, green & white rings.	4109373
9	Screw, set, cup pt, 10-32 x 7/8 in. (2x) mtg stud	HW04011
10	Label, "Auxiliary Oxygen".	4109381


SPARE AND REPLACEMENT PARTS (continued)

NM MOBILE




NM MOBILE	SPARE AND REPLACEMENT PARTS (continued)
------------------	--

ITEM	DESCRIPTION	PART NUMBER
1	O2 Regulator.	4103590
2	Check valve assembly.	4113932
3	O2 Yoke	1101620
4	N2O Regulator	4103591
5	Check valve assembly.	4113932-002
6	N2O Yoke	1101621
7	Screw (index pin), sltd, 0.157OD x 6-32 x 0.718 L (2x per yoke)	4105929
8	Plug assembly, yoke (2x)	4112755-001
9	Screw, 10-32 x 1 in. cap skt hd (4x)	HW01096
10	O ₂ Flush Valve:	
11	L-fitting, 1/4 tube x 1/8 MPT.	4109410
12	T-fitting, 1/4 tube x 1/4 tube x 1/8 MPT	4109406
13	Spacer-R unthreaded x 0.170	4110792-090
14	Spacer-R unthreaded x 0.049	4110792-021
15	Set screw, 3-48 x 3/16 in. cup pt	HW04020
16	Button	4103249
17	Knob, O2 Flush	4103248
18	Valve, Clippard 2-way	4103340
19	Constrictor, 0.245 OD x 0.046	4115092



NM MOBILE	SPARE AND REPLACEMENT PARTS (continued)
------------------	--

ITEM	DESCRIPTION	PART NUMBER
	Frame Assembly	4115163
1	Caster w/o brake (2x)	4113006-009
2	Caster w/ brake (2x)	4113006-010
3	Foam Pad (two in each vaporizer holder, four supplied w/machine)	4115519



NM MOBILE	SPARE AND REPLACEMENT PARTS (continued)
------------------	--


ITEM	DESCRIPTION	PART NUMBER
1	Power Supply Assembly, NMM (early design)	
2	PCB Assembly	4113579
	PCB Assembly (Service Exchange).	SE4113579
3	Condor Power Supply	4114523
	Condor Power Supply (Service Exchange).	SE4114523
4	Screw, 6-32 x 3/8 in. btn hd skt (6x)	HW09075
5	Lock washer, #6 int-t (6x).	HW67007
6	Power Supply Assembly, NMM (later design)	
7	Condor Power Supply	4115264
	Condor Power Supply (Service Exchange).	SE4115264
8	Battery, 12 V rechargeable.	4114229

Electrical items not shown:

Main Cable Assembly (internal)	4114289
Lamp Assembly (used with early design)	4114330
Lamp Assembly (used with later design)	4115618
Power Cord, 6 ft	4110334


SPARE AND REPLACEMENT PARTS (continued)

NM MOBILE




NM MOBILE	SPARE AND REPLACEMENT PARTS (continued)
------------------	--

ITEM	DESCRIPTION	PART NUMBER
	Absorber Assembly	4114117-001
	Manual/Auto Selector Valve is shown on a subsequent page	
	Expiratory Valve, PEEP Valve and Breathing Pressure hose Assembly are shown on a subsequent page	
1	Inspiratory Valve Assembly w/O ₂ sensor mount	4112773-001
2	Plug Assembly, Oxygen Sensor	4106387
3	Dome & Label, Insp. Valve	4108329
4	Gasket, Valve Dome	2109231
5	Disk	M23225
6	Gasket, Valve Mount	1101690
7	Ultrasonic Flow Sensor	4115754-001
	Service Exchange part number	SE4115754-001
	Flow Housing	4114444
	Transducer (2x)	4114445
	O-ring, set of six	4115147
7A	Flow sensor bracket	4116317
7B	Clamp	4115933
7C	Clamp Screw (2x)	HW01103
7D	Lock washer (2x)	HW65000
8	Connector Hose	4114912
9	Gauge Assembly, (incl. fitting for breathing pressure hose)	4114290
	Replacement Cover	4113387
	Replacement Ring	4113388
10	Gauge Mount Adapter	4114234
11	O-ring, #010 (neoprene)	4101872
12	Canister (Incl lower gasket)	4105851
13	O-ring, #335 (EPDM)	4114167
14	Dust Cup	4114094
15	Fresh Gas Hose	4108577-001
16	Gasket, canister top	4105848
17	Screen, canister	1100022
18	Gasket, canister bottom	4105849
19	Wing nut (2x)	4114087
20	Spring, CPRSN 0.420 OD x 2.5 L	4110975-032




NM MOBILE	SPARE AND REPLACEMENT PARTS (continued)
------------------	--

ITEM	DESCRIPTION	PART NUMBER
	Valve, Man/Auto Selector	4114165
1	APL Valve	4104839
2	Fiber Washer. (Supplied with APL Valve)	
3	O-Ring, #120 EPDM (2x)	4112629-001
4	Connector, 22mm (2x).	4114093
5	O-Ring, #117 Silicone	4105766
6	Lock Washer, #8 split (3x)	HW65011
7	Screw, Selector Valve Mounting, 8-32 x 7/16 in. Skt Hd Cap (3x)	HW01013
8	Auto-Bag Valve Body	Order Valve, Man/Auto Selector



NM MOBILE	SPARE AND REPLACEMENT PARTS (continued)
------------------	--

ITEM	DESCRIPTION	PART NUMBER
1	Breathing Pressure Hose Assembly	4114281
2	Hose barb coupling, 1/8 insert	4114265
3	O-ring (2x)	4106388
4	QDISC w/o valve, 1/4 OD hose	4108137
5	O ₂ Sensor Housing Assembly	4106363
6	O ₂ Sensor Capsule	6850645
	(old number: 6803290)	
7	O-ring, O ₂ sensor housing (2x)	4106388
8	Expiratory Valve Assembly	4112150
9	Dome	2109230
10	Gasket, Valve Dome	2109231
11	Disk	2123249
12	Gasket, Valve Mount (2x).....	1101690
13	PEEP Valve Assembly	4114164
14	Gasket, PEEP Valve to Absorber Mount.....	1101690-001



SV00333

NM MOBILE	SPARE AND REPLACEMENT PARTS (continued)
------------------	--

ITEM	DESCRIPTION	PART NUMBER
1	Scavenger	4114255
2	Mounting Bracket supplied with scavenger	
3	Screw, 6-32 x ¼ in. pan hd (2x)..... supplied with scavenger	
4	Bag, 3L.....	9995430

Not Illustrated:

Lint Filters (10 pcs.).....	1199
Lint Filter Container	14014
Anti-Occlusion Cage.....	4118551

SPARE AND REPLACEMENT PARTS (continued)

NM MOBILE

ITEM	DESCRIPTION	PART NUMBER
------	-------------	-------------

Miscellaneous Items:

	Touch-up paint: Euro white	7900380
	Touch-up paint: Euro blue.....	7901261
	Loctite #271 (Red)	4118558-003
	Loctite #425 (Blue)	4118558-008
	Loctite #222 (Purple).....	4118558-001
	Loctite #747 (Primer)	4118558-014
	Dow Corning High Vacuum Grease	4105908

NM MOBILE

SPECIFICATIONS

General

Anesthesia machine dimensions (approx.) (W x H x D) 23½ x 53½ x 24½ inches
Anesthesia machine weight (approx.) 163 lbs

Environmental

Storage

Temperature -20 to 60° C
Humidity 10–90% relative humidity (noncondensing)
Barometric Pressure 787 to 523 mmHg

Operating

Temperature 10 to 35° C
Humidity 30–70% relative humidity (noncondensing)
Barometric Pressure 787 to 523 mmHg

Electrical

Equipment class UL 2601 Class 1, Type B, continuous operation, IPXO
Leakage current ≤ 300 microamps (UL 2601)
Ground impedance ≤ 0.1 ohm (60 Hz source)
Dielectric withstand ≥ 1500 VAC (UL 2601)
Chassis resistance (between any metallic point
and ground pin on power cord) ≤ 0.1 ohm

Power Supply

Primary input voltage (acceptable range) 100 to 240 VAC @ 50/60 Hz
Primary input current ≤ 2.5 amps @ 50/60 Hz

Backup Battery

Type Sealed Lead Acid, 12 VDC, 3.4 Ah
Charging time ≤ 12 hours
Reserve power time (from full charge) ≥ 90 min

SPECIFICATIONS (continued)

NM MOBILE

Gas Delivery System

Pipeline inlet connections	DISS/male (ANSI B57.1-1977)
Pipeline inlet pressure	50–55 psi (345–380 kPa) (O ₂ , N ₂ O, Air)
Pipeline gauge accuracy	±3 psi (0–25 psi)
	±2 psi (25–75 psi)
	±3 psi (75–100 psi)
Cylinder connection	Pin-indexed hanger yoke (ANSI/CGA V-1-1987)
Over pressure relief valve	95 psi (655 kPa)
Fresh gas common outlet	15 mm female
Fresh gas oxygen concentration (ORC)	25 ±4%
Oxygen flush flow rate	55 (±10) L/min
Minimum oxygen flow (at 50 psi pipeline pressure)	150 ±50 mL/min
Low oxygen supply pressure alarm	34–40 psi
Cylinder gauge accuracy	±90 psi (0–750 psi)
	±60 psi (750–2250 psi)
	±90 psi (2250–3000 psi)

Cylinder Gas Pressure

Oxygen	1900 psi (13100 kPa) E-size cylinder (at 70° F, 21° C)
Nitrous Oxide	745 psi (5130 kPa) E-size cylinder (at 70° F, 21° C)

Flowmeter Accuracy (at 20°C and 760 mmHg)

Oxygen, Nitrous Oxide, Air	
Dual Tapered 0-8 L/min	0.1 to 0.2 L/min ± 50 mL/min
	0.2 to 1.0 L/min ±100 mL/min
	1.0 to 8.0 L/min ±5% FS
Oxygen (Auxiliary Oxygen)	0–10 L/min ±5% FS

SPECIFICATIONS (continued)

NM MOBILE

Oxygen Monitoring

Range	10–100 vol % O ₂
Resolution	1 vol % O ₂
Accuracy	±3 vol % O ₂ (When calibrated within 18 hours, and constant temperature and pressure)
Response time	≤ 25 sec (T90)
Zero drift	≤ 0.1 vol % O ₂ /month
Span drift	≤ 1 vol % O ₂ /8 hours
Temperature error	≤ ± 3% of reading (15° to 40° C)
Sensor service life	≥ 8 months at 25° C, 50% relative humidity, 50% O ₂ gas mixture (or ≥ 5000% hour CO ₂)

Breathing Pressure Monitoring

Numeric display range	-10 to 125 cmH ₂ O
Resolution	1 cmH ₂ O
Accuracy	±3 cmH ₂ O or ±10% of reading, whichever is greater
Waveform display range - full	0–100 cmH ₂ O
Waveform resolution	1 cmH ₂ O
Waveform accuracy	±3 cmH ₂ O or ±10% of reading, whichever is greater
Waveform display scales	0–20, 0–50, 0–100 cmH ₂ O

Respiratory Volume Monitoring

Minute Volume

Display Range 0.2 to 50.0 L
Resolution 0.1 L
Accuracy $\leq +10\%$ or 0.1 L, whichever is greater

Tidal Volume

Display Range 0.05 to 20.1 L
NOTE: The standard bellows will deliver up to 1.5 L
Resolution 0.01 L
Accuracy $\leq \pm 10\%$ or 0.01 L, whichever is greater
Volume Apnea Threshold 0.05 L

Respiratory Rate

Numeric display range 2–99 BPM
Resolution 1 BPM
Accuracy $\leq \pm 10\%$ or ± 1 BPM, whichever is greater

Serial Interface

Serial Port

Type RS-232/422
Baud Rate 300 to 38400
Parity Odd, Even, or None
Data Bits 7 or 8
Stop Bits 1 or 2
Protocol Vitalink

[RETURN TO THIS MANUAL'S TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

[RETURN TO THIS MANUAL'S TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

[RETURN TO THIS MANUAL'S TABLE OF CONTENTS](#)
[RETURN TO CD-ROM TABLE OF CONTENTS](#)

Drägermedical

A Dräger and Siemens Company

DrägerService is a division of
Dräger Medical, Inc.
3122 Commerce Drive
Telford, PA 18969
Tel: (215) 721-5402
(800) 543-5047
Fax: (215) 721-5784
Web: www.draegermedical.com
Printed in the U.S.A.

Narkomed Mobile Service Manual

Rev. R summary of changes

Page	Description
Cover.....	Changed Rev., and Date
6-1.....	Changed P/N of PMC Form
8-33.....	Changed P/N of item #5

Rev. P summary of changes

Page	Description
Cover.....	Changed Rev., and Date
6-1.....	Changed P/N of Test Lung
6-7.....	Changed P/N of Test Lung

Rev. N summary of changes

Page	Description
Cover.....	Changed Rev., and Date
TOC.....	Added new Miscellaneous Items to Section 8
8-9A.....	Added Loctite Reference on illustration
8-10.....	Added Loctite Reference on illustration
8-12.....	Added Loctite Reference on illustration
8-16.....	Added Loctite Reference on illustration
8-20.....	Added Loctite Reference on illustration
8-26.....	Added Loctite Reference on illustration
8-28.....	Added Loctite Reference on illustration
8-39.....	Added Items not Shown
8-40.....	Added Miscellaneous Items section

Rev. M summary of changes

Page	Description
Cover.....	Changed Logo, Rev., and Date

TOCChanged Logo
8-7Change P/N of regulator rebuild kit
Back CoverChanged Logo

Rev. L summary of changes

PageDescription

Section 6 PMCREvised entire section for PMC Reduction Effort

Rev. K summary of changes

PageDescription

[8-33](#)Revised ultrasonic flow sensor part numbers