
Alaris® PK Syringe Pump

Propofol-Paedfusor
Model

Contents

EN - Directions for Use.....	1
DE - Gebrauchsanweisung	4
ES - Manual del Usuario.....	7
FR - Mode d'emploi	10
IT - Istruzioni per l'uso	13
NL - Gebruiksaanwijzing	16
NO - Bruksanvisning.....	19
PT - Instruções de utilização.....	22
SE - Bruksanvisning	25

Introduction

This addendum should be used in conjunction with the Alaris® PK Syringe Pump Directions For Use. The purpose of this addendum is to provide additional material for operation when using a drug configuration based upon the Propofol - Paedfusor model. To use this drug / model combination a suitable data set must be created, see Alaris® PK Editor Software Directions For Use.

TCI Precautions

When first starting the infusion the pharmacokinetic / pharmacodynamic models and memory registers within the Alaris® PK Syringe Pump are reset to zero. Therefore, if for any reason the pump is switched off during the surgical procedure, all current pharmacokinetic / pharmacodynamic model information, including the volume of propofol already administered, will be lost. Under such circumstances switching the pump on and restarting the infusion whilst the patient contains a significant residual drug dose could result in an over-infusion and significant cardiovascular instability. Therefore, in these circumstances, the pump should not be restarted in TCI mode.

Pharmacokinetic and Pharmacodynamic parameters

Drug: Propofol

Model: Paedfusor

Ages 1 – 12 years:

$V1 = 458.4 \times \text{weight}$

$V2 = V1 \times k_{12} / k_{21}$

$V3 = V1 \times k_{13} / k_{31}$

$k_{10} = 0.1527 \times \text{weight}^{-0.3}$

$k_{12} = 0.114$

$k_{13} = 0.0419$

$k_{21} = 0.055$

$k_{31} = 0.0033$

Age ≥ 13 :

$V1 = 400.0 \times \text{weight}$

$k_{10} = 0.0678$

Other micro-constants as above

Age ≥ 14 :

$V1 = 342.0 \times \text{weight}$

$k_{10} = 0.0792$

Other micro-constants as above

Age ≥ 15 :

$V1 = 284.0 \times \text{weight}$

$k_{10} = 0.0954$

Other micro-constants as above

Age ≥ 16 :

$V1 = 228.57 \times \text{weight}$

$k_{10} = 0.119$

Other micro-constants as above

Model Restrictions

The Paedfusor model is a three-compartment PK model developed for use in children. It uses age and weight as covariates, and has been validated in one clinical study involving children undergoing cardiac surgery or cardiac catheterisation¹. We recommend that users read the publication carefully, in order to understand the limitations of the model and the publication.

Although the age range in the study was between 1 and 15 years, there were only 5 patients younger than 3 years, of which only 1 was younger than 2 years. The weight of the children in the study varied from 5 to 53kg, but only 3 children had a body weight of <10kg. Thus the model should not be used in children younger than 1 year, or who weigh less than 5kg. In children younger than 3 years, or who weigh less than 10kg, the model should only be used by consultant paediatric anaesthetists with extensive experience and a special interest in total intravenous anaesthesia in children.

It should be remembered that model accuracy in children with normal cardiac function may be different to that found in the study population. Also, children in the study received fentanyl or alfentanil, both of which exhibit pharmacokinetic interactions with propofol (causing greater propofol concentrations than would have been found if propofol alone had been administered).

As in adults, no single target concentration is likely to be sufficient for all patients. Thus it remains necessary to observe the patient responses closely and titrate the target concentration according to the clinical effect. In children who are unwell (ASA status 3 – 4), particularly those with hypotension or impaired hepatic function, target concentrations should be selected with care. To reduce the risk of cardio-vascular instability, lower initial target concentrations should be selected, and the target only be increased after careful assessment of the patient.

The recommended infusate concentration is 10 mg/ml (1%) or 20 mg/ml (2%). Note that the administration of 2% propofol from large syringes in small children, or at low target concentrations, may result in sub-optimal system accuracy (see below). At the release of this document, no prescribing information is recommended for TCI propofol administration in children. It is recommended that the 'Clinical Trial' option is set within the Alaris® PK Editor software until prescribing information for the Paedfusor model is available.

References:

- 1 Absalom A, Amutike D, Lal A, White M, Kenny GNC: Accuracy of the 'Paedfusor' in children undergoing cardiac surgery or catheterisation *British Journal of Anaesthesia* **91** (4): 507-13 (2003)
- 2 Absalom et al: 'Paedfusor' pharmacokinetic dataset *British Journal of Anaesthesia* **95** (1): 110-13 (2005)

Profiles from TCI mode

When targeting in TCI Mode the Alaris® PK Syringe Pump will automatically calculate the flow rate profile from the specific pharmacokinetic/pharmacodynamic model for the selected drug. This section of the Directions For Use is intended to help users understand the profiled infusion and the performance accuracy attained from the TCI pump.

Induction Bolus and maintenance rates are displayed before starting the titration. When initially starting the infusion or after increasing the target (plasma or effect) concentration by titration, the pump will first deliver a bolus dose through a typically short, high rate infusion. On completion of this bolus, the pump will immediately switch to a lower maintenance rate (when plasma target mode is used) or will pause for a period of time before switching to a lower maintenance rate (when effect site targeting mode is used). Once the maintenance phase is reached, any reduction made to the target (plasma or effect) concentration will typically result in the infusion rate reducing to zero until the predicted plasma (or effect) concentration reduces the new target value.

The Alaris® PK Syringe Pump updates the pharmacokinetic model driving the plasma (or effect) concentration prediction and the infusion rate every 10 seconds. The infusion rate graph, shown on page 3, were measured in accordance with the protocol described in the IEC60601-2-24¹ Standard, with the data sample period reduced from 30 to 10 seconds.

The pump solves the pharmacokinetic/pharmacodynamic algorithms so that the target (plasma or effect) concentration is attained as rapidly and as accurately as possible. However, the User may need to take into consideration the limitations of the physical system in attaining the target (plasma or effect) concentration; this includes:

- The limit on the flow rate permitted by the infusion pump mechanism;
- The limit on the flow rate permitted by the syringe size;
- The patient / drug dose limitation from the prescribing information to insure the safety of the administration;
- The variation in individual patient response to reach the plasma (or effect) concentration;
- The model specific cap rate.

A true assessment of the performance of the Alaris® PK Syringe Pump can be made if the volumetric error, that is the difference between the actual volume infused and the predicted volume infused, is calculated. For the performance graphs shown on page 3, over a one hour period, the Alaris® PK Syringe Pump has a mean volumetric accuracy in TCI Mode better than $\pm 5\%^2$.

By measuring the volume from the flow rate profile delivered from the Alaris® PK Syringe Pump and then introducing this into a reverse pharmacokinetic model the predicted plasma (or effect) concentration can be calculated from the flow rate. These are illustrated on page 3, showing the typical performance of the system against changes in the target plasma (or effect) concentration for a typical, idealised profile. For the same targeted profile, the deviation of the predicted plasma (or effect) concentration (back calculated from the volume collected) from the expected Ideal plasma (or effect) concentration, results from the volumetric inaccuracy of the system (pump and syringe). The Alaris® PK Syringe Pump will track the predicted plasma (or effect) concentration to within $\pm 5\%^2$ of that calculated by pharmacokinetic model over a one hour period. Flow rate inaccuracies and start-up delays may decrease the accuracy of the predicted plasma (or effect) concentration particularly where high syringe drug concentrations are used in conjunction with large sizes of syringes and low target plasma (or effect) concentrations as the syringe plunger motion over time (proportional to the flow rate accuracy) will be significantly reduced.

Note: For a given drug concentration, the volumetric error is proportional to the dose rate error. Knowledge of the system accuracy over different time intervals may be of interest when assessing the impact of administering short-half life drugs. In these circumstances, short-term fluctuation in the infusion rate could have a clinical impact that cannot be determined from the performance profiles shown in Figures below. In general, the volumetric error will increase with small induction and maintenance rates, which may occur with large volume syringes, high syringe concentrations, low patient weights and low target (plasma or effect) concentrations. For applications where system accuracy is important, maintenance rates less than 1.0 ml/h are not recommended; syringe sizes, drug concentrations / dilutions and target (plasma or effect) concentrations should be selected accordingly to ensure the maintenance rate exceeds this lower limit.

Note:

¹ IEC60601-2-24: Particular Requirements for the Safety of Infusion Devices;

² 95% Confidence / 95% Population.

Profiles from TCI mode


The Paedfusor model should only be used with Alaris® PK Syringe Pumps with software revision v3.2.12 and above. It should be noted, with revision v3.2.12, the minimum age and weight available is 3 years and 15kg respectively. Later versions allow settings down to 1 year and 5kg.

Profiles from TCI mode - BD 50ml Syringe

Target Plasma, Propofol 1%


Infusion Rate vs Target Plasma Concentration

Patient age: 9 years
Patient Weight: 38 kg
Drug Concentration: 10 mg/ml
Volumetric accuracy: -0.50%

Predicted vs Ideal Plasma Concentration


Patient age: 9 years
Patient Weight: 38 kg
Drug Concentration: 10 mg/ml
Plasma conc. accuracy: -0.48%

Target Plasma, Propofol 2%


Infusion Rate vs Target Plasma Concentration

Patient age: 9 years
Patient Weight: 38 kg
Drug Concentration: 20 mg/ml
Volumetric accuracy: -0.61%

Predicted vs Ideal Plasma Concentration

Patient age: 9 years
Patient Weight: 38 kg
Drug Concentration: 20 mg/ml
Plasma conc. accuracy: -0.56%

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
Alaris is a registered trademark of CareFusion Corporation or one of its subsidiaries. All rights reserved.
All other trademarks are property of their respective owners.

© 2005-2010 CareFusion Corporation or one of its subsidiaries. All rights reserved.

This document contains proprietary information of CareFusion Corporation or one of its subsidiaries, and its receipt or possession does not convey any rights to reproduce its contents, or to manufacture or sell any product described. Reproduction, disclosure, or use other than for the intended purpose without specific written authorization of CareFusion Corporation or one of its subsidiaries is strictly forbidden.

CareFusion Switzerland 317 Sàrl,
CH-1180, Rolle

CareFusion U.K. 305 Ltd., RG22 4BS, UK

1000DF00447 Issue 1