

Typhoon



The model name and serial number should always be provided when ordering spare parts or during queries per telephone or written correspondence:

Serial number of the machine:

.....

Machine model:

.....

WARNING

To avoid injury, always read this *Instructions for Use* and accompanied documents before using the product.



Mandatory to read the *Instructions for Use*.

Design Policy and Copyright

® and ™ are trademarks belonging to the Arjo group of companies.

© Arjo 2019.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of Arjo.

Contents

PREFACE.....	5
1 SAFETY PRECAUTIONS.....	6
1.1 Important.....	6
1.2 Isolator switch.....	7
1.3 In an emergency.....	7
1.4 Product liability.....	7
1.5 Attention symbols.....	7
2 PRESENTATION	8
2.1 Intended field of application.....	8
2.2 Components of flusher-disinfector	8
2.3 Cleaning principles	8
2.4 Container	8
2.5 Symbols on the machine	9
2.6 Control panel	9
2.6.1 Components.....	9
2.6.2 Functions	10
2.7 Display options	11
2.7.1 Displayed while a program is running.....	11
2.7.2 Display after program end.....	11
2.8 Door function	12
2.8.1 Opening the door	12
2.8.2 Closing the door.....	12
2.8.3 Aborting door closing	12
3 OPERATING INSTRUCTIONS.....	13
3.1 Program overview.....	13
3.2 Run a program.....	14
3.2.1 Start a program.....	14
3.2.2 During ongoing program	14
3.2.3 After the program ends	14
3.3 Examples of goods and program selection.....	15
4 MAINTENANCE.....	16
4.1 Flushing chamber	16
4.2 The outside	16
4.3 Preventive maintenance	16
4.3.1 Periodic maintenance	16
4.3.2 Service table	17
4.4 Machines that are not in use.....	19
5 ALARM.....	20
5.1 Alarm.....	20
5.2 Acknowledgement of error message	20
5.3 During a power failure.....	20

6	TROUBLESHOOTING.....	21
6.1	List of codes.....	21
6.2	Types of codes.....	21
7	ENVIRONMENTAL PRODUCT DECLARATION	22
7.1	Comments	22
7.2	Distribution.....	22
7.3	Use	22
7.4	Waste disposal	22
8	DETERGENT.....	23
8.1	Flushing detergent.....	23
9	WATER QUALITY	24
10	END OF LIFE DISPOSAL.....	26

PREFACE

The user manual describes the design and operation of the machine and the maintenance for which the user is responsible. The installation manual describes the machine's installation. The technical manual describes necessary information for operators and maintenance personnel. The information in the manual must be complied with to ensure safe and error-free operation.



Read the manual before use.

Users must read through the manual before using the machine for the first time, as well as familiarizing themselves with the operation of the machine and its safety instructions. Operators and maintenance personnel must have completed training through the marketing organization of Arjo.

The information in this manual describes the machine as dispatched from Arjo. There may be differences due to customization for customers or countries.

The machine is accompanied by the following documentation:

- User manual
- Installation manual
- Quick guide

The following documentation shall be provided by the sales company to the end-user on request:

- Technical manual
- Spare parts list

The accompanying Quick guide must be hung up and clearly visible in the installation area.

Arjo reserves the right to change the specification and design without prior notice. The information in this manual was up-to-date on the date that the manual was issued.

1 SAFETY PRECAUTIONS

The machine has a number of built-in safety mechanisms. To avoid injury, it is highly important not to bypass these safety mechanisms.

1.1 Important

- Read the manual thoroughly before use.
- The machine should only be handled by authorized personnel. Personnel must also receive regular training.
- If any part is missing or damaged - do NOT use the product.
- Exercise caution when using the machine, as it uses hot water and possibly steam. Because the machine uses hot water, there may be hot surfaces in the facility.
- Exercise caution when handling chemical agents used in the machine. Follow the instructions for the packaging the agent comes in:
 - if the agent is swallowed, comes into contact with the eyes or skin or if vapors are inhaled.
 - for instructions on dosage amount and washing temperatures.
 - for instructions regarding storage and sorting/disposal of packing.
- Make sure that the machine has access to detergent before the process starts. This is important both for the operation of the machine and for keeping the washing chamber free from staining.
- The machine must be kept clean if it is to work properly.
- Do not wash down or hose down the machine with water.
- Installation and service work must be performed by trained personnel.
- Never bypass the machine's door switch under any circumstances.
- Leakage in the system, due to a worn door seal for example, must be repaired without delay.
- Spare parts must only be obtained from Arjo sales companies.
- If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

1.2 Isolator switch

The machine must have a lockable switch for the electric power supply. The switch must be easily accessible on a wall close to the machine. The installation must conform to and be marked in accordance with local provisions. A switch is not required on machines fitted with an electric plug.

1.3 In an emergency

- Switch off the main switch.
- Close shut off valves in the water supply line.

1.4 Product liability

CE marking indicating conformity with European Community harmonised legislation.

Any modification or incorrect use of the equipment without the approval of Arjo negates Arjo's product liability.



CE marking indicating conformity with European Community harmonised legislation
Figures indicate Notified Body supervision.



Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745

1.5 Attention symbols

Some of the warnings, instructions and advice in this manual require extra attention. The symbols and designs used are:



Risk of injuries or machine damage.

2 PRESENTATION

2.1 Intended field of application

Flusher-disinfector for emptying, cleaning and disinfecting, bedpans, urine bottles and other containers for human waste, including suction bottles.

The customer is responsible for ensuring that an Installation Qualification, Operating Qualification and Performance Qualification according to ISO 15883 are carried out before the product is put into use.

2.2 Components of flusher-disinfector

Typhoon is one of Arjo's series of flusher-disinfectors.




Figure 1. Components of flusher-disinfector

1. Control panel
2. Door to machine cabinet
3. Door
4. Flushing chamber
5. Container for descaler and/or detergent

2.3 Cleaning principles

The goods are cleaned by flushing with cold and hot water and then disinfecting with steam.

2.4 Container

The holder in the flushing chamber, in which the goods are placed, is available in several models. All models come with separate instructions for positioning the goods.

2.5 Symbols on the machine



Warning. See manual



Hot surface. Risk of injury.

2.6 Control panel

2.6.1 Components















Figure 2. Control panel

1. Yellow: Process running
2. Green: Process complete
3. Red: Disinfection not achieved or other error code F
4. Display
5. Rim flushing
6. Option (detergent)
7. Open/close door.
8. Sensor to open/close door and to start preselected program
9. Economy program
10. Normal program
11. Intensive program

2.6.2 Functions

Note! Disinfection is not included in the rim flushing program.

Button/lamp	Function
	When the yellow light comes on: The machine is running a program. The temperature/A0-value/time and cleaning phase can be seen on the display.
	When the green light comes on: The program is complete. The green lamp goes out when the door is opened and the maximum temperature of the process shows on the display.
	When the red light comes on: The machine has not reached the disinfection temperature or another error has occurred (Error code F).
	Button to start the economy program for urine bottles.
	Button to start the normal program for normally soiled goods.
	Button to start the intensive program for heavily soiled goods.
	The button combination for selecting programs with cleaning agents for heavily soiled goods.
	The button combination for selecting the program with cleaning agents for general circulation goods.
	Button to start rim flushing Without disinfection!
	Sensor to open/close door and to start preselected program
	Open/close door

2.7 Display options

A service technician can set the information to be shown in the display.

2.7.1 Displayed while a program is running

The display can be set to show one of the following alternatives while a program is running:

- Temperature
- A_0 value
- A_0 value / Temperature
- Remaining Time
- Temperature / Remaining time
- A_0 value / Remaining time
- A_0 value / Remaining time / Temperature

If an alternative with more than one parameter is selected, the parameters will be displayed alternately. The A_0 value always remains at 0 until the thermal disinfection begins. When the temperature reaches 85°C, the A_0 value will increase in line with the temperature and time.

2.7.2 Display after program end

The display can be set to show one of the following alternatives once a program has finished:

- Max. temperature
- A_0 value
- A_0 value / Max temperature

If an alternative with more than one parameter is selected, the parameters will be displayed alternately.

2.8 Door function

2.8.1 Opening the door

- Move your hand toward the sensor, stopping about 3-4 cm from the panel.
- Keep your hand still for one 1 second.
- The door opens automatically.




Figure 3. Open door using sensor




Figure 4. Door opens

2.8.2 Closing the door



WARNING!
Keep hands and arms away when an automatic door has been activated in order to avoid the risk of getting stuck.




- Activate the sensor.
- Alternatively, the door can be closed by choosing a program.




2.8.3 Aborting door closing

- Activate the sensor again when the door is coming down. The closing sequence stops and the door goes up again.

3 OPERATING INSTRUCTIONS

3.1 Program overview

Program	Economy	Normal	Intensive
Degree of soiling	Lightly soiled	Normal soiling	Heavily soiled
Press button			
Program sequence	Cleaning Disinfection Cooling End of program	Cleaning Disinfection Cooling End of program	Cleaning Disinfection Cooling End of program
Program time	6-8 min	7-9 min	9-10 min
Energy	0.25 kWh	0.25 kWh	0.25 kWh
Water (including cooling)	25 l ± 10%	31 l ± 10%	39 l ± 10%






The flushing-disinfector's three standard programs are described above. Other programs may be available using the buttons, , , and/or , see Technical Manual.

3.2 Run a program

3.2.1 Start a program



If the machine has not been used for more than 72 hours, it is necessary to run a program with an empty chamber before using the machine to clean goods.

1. Open the door. (Use the door sensor on the machine.)
2. Place the goods to be cleaned in the holder according to the attached instructions.
3. If the holder has become stuck in the wrong position, press  (open door) to reset the holder to the starting position.
4. The door closes when a program is selected, or when the door sensor is activated.
5. Choose a program from the table below.
 - Economy program: Press .
 - Normal program: Press .
 - Intensive program: Press  if the goods are badly soiled.
 - To add detergent, press .
6. The yellow light under will blink for 6 seconds and during the time it is possible to abort the program.

3.2.2 During ongoing program




If the red indicator lamp lights up an error has occurred. The process is cancelled and must be repeated.

3.2.3 After the program ends



The goods may be hot at the end of the program. If the goods are warmer than 60°C, U7 appears in the display.


1. The green light comes on  when the program is complete. The door then opens automatically.
2. When the door has opened, take out the goods.

3.3 Examples of goods and program selection


Economy program



Normal/Intensive program



Normal/Intensive program



4 MAINTENANCE

4.1 Flushing chamber

If limescale deposits occur, they must be removed with ordinary acid-based cleaning agent. There is no need to rinse it off. This happens automatically during the next program.




Figure 5. The wash chamber is cleaned using acid-based detergent

4.2 The outside



Do not wash the machine down with water.

The front and sides and the control panel can be cleaned and disinfected with an alcohol-based cleaner, medical spirits or a common cleaning agent. Remove stains with an ordinary cleaner for stainless steel.

4.3 Preventive maintenance

For safety purposes, never modify equipment or use incompatible parts.

4.3.1 Periodic maintenance

Periodic maintenance and system testing must be performed to ensure safety and the machine's proper operation.

The amount of maintenance required depends largely on the quality of the incoming water and how often the machine is used. The maintenance interval will have to be determined in each individual case. Arjo recommends that the stated maintenance operations be performed at intervals according to the service table.

4.3.2 Service table

NOTE! The service table is only shown for information purposes.



Maintenance may only be performed by authorized service technicians.

The service table below shows recommended maintenance intervals.

In addition to system testing, we recommend that appropriate cleaning tests and temperature validation are performed according to EN ISO 15883.

	Measure	Every year / 10,000 cycles	Every other year / 20,000 cycles	Time* (minutes)
1	General			
1.1	Check the cabling and connection points.	•		10
1.2	Check that the panel decal is firmly affixed, is undamaged and legible.	•		2
1.3	Check that an goods placement sign and quick help instructions have been put up.	•		1
1.4	Check error codes and the number of processes.	•		2
2	The chamber			
2.1	Check that the rotating nozzles are moving freely. Clean as required.	•		20
2.2	Check that the fixed nozzles are not clogged with grime and deposits. Clean as required.	•		20
2.3	Check for leakage in the nozzle attachments and the hose connections to the chamber.		•	-
2.4	Check for proper operation of the temperature sensor.	•		10
2.5	Check that the door seal and chamber seal are tight. Clean or replace as necessary.	•		10
2.6	Check the mounting for the goods holder.	•		2
2.7	Check that the overflow drain connection to the chamber is tight.		•	-
2.8	Clean the inside of the flushing chamber.	•		10
3	Process tank			
3.1	Clean as required.	•		-

	Measure	Every year / 10,000 cycles	Every other year / 20,000 cycles	Time* (minutes)
3.2	Check and clean the level sensor.	•		5
3.3	Check that the valves are working and that the connections are not leaking. Clean as required.	•		10
4	Steam Generators			
4.1	Check the connections to the steam generator for leakage and make sure the surrounding insulation is intact and no hot surfaces are exposed.	•		5
4.2	Check that the steam generator is working properly.	•		10
5	Main pump			
5.1	Check that the main pump is firmly attached and that there are no leaks in the connections.	•		5
6	Dosing system			
6.1	Check the hoses between the dosing pump and the detergent container. Replace if necessary.	•		5-10
6.2	Check the dosage amount. Adjust if necessary.	•		10-30
6.3	Replace hose in the dosing pump.	•		5
6.4	Check that the empty container alarm works. Clean or replace as necessary.	•		5
6.5	Check the function of the check valve in the steam generator	•		5
7	Drainage connection			
7.1	Check that the drain connection is firmly attached and that it is not leaking.	•		5
7.2	Check the condensation cooling nozzle head (option).	•		5
8	Steam inspection (Option)			
8.1	Check and adjust the needle valve.	•		10
9	Automatic door			
9.1	Check the automatic door for proper operation.	•		5
10	Function check (at each service)			
10.1	Run a complete process and be sure there are no deviations according to the periodic checkpoints.			
10.2	Check that the panel's lamps and alarms are functioning.			


*The time required is estimated with the top off and may vary depending on the installation settings, operating conditions and equipment level.

4.4 Machines that are not in use

A machine that is not being used should be drained. Contact a service technician.



5 ALARM

5.1 Alarm


If the red lamp at  lights up, disinfection was aborted because of low temperature or some other error. The display shows an error code.

5.2 Acknowledgement of error message

Some error messages need to be acknowledged before the machine can be returned to run mode. When the cause of the error has been found and the error has been put right, the error message is acknowledged as follows:

First press the  button (1), then the  button (2). Keep both buttons pressed at the same time for 10 seconds.

The machine is now ready to operate. Re-run the aborted program if the goods are still in the machine. If the error persists, call service personnel.



5.3 During a power failure

If the power fails during a process, the machine remains locked. When the power returns, the machine restarts the process automatically.

6 TROUBLESHOOTING

6.1 List of codes

Error code	Description of error	Measure
U1	Lime binding agent low	Top up lime binding agent
U4	Process agent low level	Top up process agent
U6	Inspection request	Contact authorized service personnel
U7	Warning of high temperature for goodss	Let the goods cool down
U8	External communication	Contact authorized service personnel
H7	Battery error	Contact authorized service personnel
F1	Faulty temp sensor	Contact authorized service personnel
F2	Level sensor error, tank	Contact authorized service personnel
F3	Level sensor error, detergent	Contact authorized service personnel
F5	Logic error, machine configuration	Contact authorized service personnel
F6	Inspection request	Contact authorized service personnel
F7	Disinfection temperature not reached	Contact authorized service personnel
F8	Door open/unlocked during process	Contact authorized service personnel
F9	Settings incorrect	Contact authorized service personnel
F10	Could not fill/empty tank	Contact authorized service personnel
F11	Leakage	Contact authorized service personnel
F14	Operating time, door	Contact authorized service personnel
F15	Operating time, turner	Contact authorized service personnel
F16	Log. error, turner	Contact authorized service personnel
F18	Wrong machine type	Contact authorized service personnel

6.2 Types of codes

- Ux is an information code.
- Hx is a handling code.
- Fx is an error code.

7 ENVIRONMENTAL PRODUCT DECLARATION

Component	Weight in kg (total 135 kg)	Weight as a percentage
Mechanical parts		
Stainless steel	80	59
Steel	20	15
Glass	2	1.5
Brass	2	1.5
Plastic	7	5
Rubber	3	2
Other mech. parts	4	3
Electrical / Electronics		
Wiring	1.5	1
Circuit boards	0.5	0.5
Other components	15	11

7.1 Comments

The “Other mech. parts” and “Other components” fractions comprise pumps, door motors, etc., which cannot be weighed or divided up in the same way as other fractions.

7.2 Distribution

The machines are packed on wooden pallets and clad with wood and Masonite, or wrapped in plastic. Users must make sure that the packaging material is sorted and recycled.

7.3 Use

The flusher-disinfector affects the environment during its entire life cycle. Our assessment is that electric power and water consumption have the greatest impact on the environment. Based on a service life of ten years and 7,000 processes per year, the machines consume about 1,750 kWh of electricity and about 215 cubic meters of water. About 35 liters of agent are discharged into the drainage system, depending on the programs/doses used in the machine. To reduce the environmental impact the user may use green electric power.

8 DETERGENT

WARNING!



To avoid eye and skin damage, always use protective glasses and gloves. If contact occurs rinse with plenty of water. If eyes or skin becomes irritated, contact a physician. Always read the material safety data sheet of the detergent.



Using agents other than those tested and approved by Arjo may damage the machine (dosing pumps and hoses).



Follow the instructions for the container for additional instructions on using the agent.

The machine has been proven and tested with Arjo's detergents. We recommend that you use Arjo detergents.

8.1 Flushing detergent

Following detergents can be ordered from Arjo sales companies.

Agent	Description
■ ARJO LIQUIDS FLUSHER DETERGENT	A mild alkaline detergent to be used for automated cleaning of human waste containers.
■ ARJO LIQUIDS FLUSHER RINSE	A mild alkaline solution to be used for automated rinsing of human waste containers and system piping to prevent hard water scale build up.

9 WATER QUALITY



It is the customer's responsibility to supply the washer-disinfector with water of the right quality.

The quality of the water used in all stages of cleaning is important for good cleaning results. The water used in each step should be compatible with:

- The material that the washer-disinfector is made of.
- The chemicals used in the process.
- Process requirements for the various stages of the process

The main factors for good water quality are:

Hardness	High hardness will cause limescale deposits in the washer-disinfector, leading to poor cleaning results.
Ionic contaminants	A high concentration of ionic contaminants may cause corrosion and pitting on stainless steel. Heavy metals such as iron, manganese or copper cause instruments to tarnish.
Microbial contaminants	Goods to be washed are treated for microorganisms and their residual products, which may cause fever-like symptoms if they get into the human body. Use water that does not increase the biological load.
Sanitary chemicals	High concentrations of, and high exposure to, sanitary chemicals may cause corrosion and pitting on stainless steel.

Arjo recommends that water used in the pre-rinsing, washing and final rinsing phases should be of drinking quality in accordance with the guidelines. Recommended water quality is drinking water with max 5°dH. Detailed information about acceptable water quality can be found in "Guidelines for Drinking Water Quality, published by WHO.

Follow local standards if they are stricter than Arjo's recommendations. Especially when treated water is used for the final washing/disinfection phase. A typical specification for treated water is:

pH	5.5 to 8
Conductivity	<30 µs.cm-1
TDS	<40 mg/l
Maximum hardness	<50 mg/l
Chlorine	<10 mg/l
Heavy metals	<10 mg/l
Phosphates	<0.2 mg/l as P2O5
Silicates	<0.2 mg/l as SiO2
Endotoxins	<0.25 EU/ml
Total number of micro-organisms	<100 per 100 ml

Additional advice should also be obtained from the manufacturers of chemical and medical equipment.

10 END OF LIFE DISPOSAL

Equipment that has electrical and electronic components should be disassembled and recycled per Waste of Electrical and Electronic Equipment (WEEE) or in accordance with local or national regulation.

Disinfection liquids – if you have extra liquid, it's safe to pour it down the drain with running water. These products are water soluble and usually mixed with water when used for cleaning, your municipality will treat these cleaners the same way they would if you used them for cleaning.

The packaging containers should be recycled according to the national or local regulations.



WARNING!

The product may be contaminated and has to be disinfected before recycling.




Figure 6. Recycling symbol

000052

AUSTRALIA

Arjo Australia Pty Ltd
78, Forsyth Street
O'Connor
AU-6163 Western Australia
Tel: +61 89337 4111
Free: +1 800 072 040
Fax: + 61 89337 9077

BELGIQUE / BELGIË

Arjo NV/SA
Evenbroekveld 16
BE-9420 ERPE-MERE
Tél/Tel: +32 (0) 53 60 73 80
Fax: +32 (0) 53 60 73 81
E-mail: info.belgium@arjo.be

BRASIL

Arjo Brasil Equipamentos Médicos Ltda
Rua Marina Ciufuli Zanfelicce, 329 PB02 Galpão
- Lapa
São Paulo – SP – Brasil
CEP: 05040-000
Phone: 55-11-3588-5088
E-mail: vendas.latam@arjo.com
E-mail: servicios.latam@arjo.com

CANADA

Arjo Canada Inc.
90 Matheson Boulevard West
Suite 300
CA-MISSISSAUGA, ON, L5R 3R3
Tel/Tél: +1 905 238 7880
Free: +1 800 665 4831 Institutional
Free: +1 800 868 0441 Home Care
Fax: +1 905 238 7881
E-mail: info.canada@arjo.com

ČESKÁ REPUBLIKA

Arjo Czech Republic s.r.o.
Na Strži 1702/65
140 00 Praha
Czech Republic
Phone No: +420225092307
e-mail: info.cz@arjo.com

DANMARK

Arjo A/S
Vassingerødvej 52
DK-3540 LYNGE
Tel: +45 49 13 84 86
Fax: +45 49 13 84 87
E-mail:
dk_kundeservice@arjo.com

DEUTSCHLAND

Arjo GmbH
Peter-Sander-Strasse 10
DE-55252 MAINZ-KASTEL
Tel: +49 (0) 6134 186 0
Fax: +49 (0) 6134 186 160
E-mail: info-de@arjo.com

ESPAÑA

Arjo Ibérica S.L.
Parque Empresarial Rivas Futura, C/Marie
Curie 5
Edificio Alfa Planta 6 oficina 6.1-.62
ES-28521 Rivas Vacia, MADRID
Tel: +34 93 583 11 20
Fax: +34 93 583 11 22
E-mail: info.es@arjo.com

FRANCE

Arjo SAS
2 Avenue Alcide de Gasperi
CS 70133
FR-59436 RONCQ CEDEX
Tél: +33 (0) 3 20 28 13 13
Fax: +33 (0) 3 20 28 13 14
E-mail: info.france@arjo.com

HONG KONG

Arjo Hong Kong Limited
Room 411-414, 4/F, Manhattan Centre,
8 Kwai Cheong Road, Kwai Chung, N.T.,
HONG KONG
Tel: +852 2960 7600
Fax: +852 2960 1711

ITALIA

Arjo Italia S.p.A.
Via Giacomo Peroni 400-402
IT-00131 ROMA
Tel: +39 (0) 6 87426211
Fax: +39 (0) 6 87426222
E-mail: Italy.promo@arjo.com

MIDDLE EAST

Arjo Middle East FZ-LLC
Office 908, 9th Floor,
HQ Building, North Tower,
Dubai Science Park,
Al Barsha South
P.O Box 11488, Dubai,
United Arab Emirates
Direct +971 487 48053
Fax +971 487 48072
Email: Info.ME@arjo.com

NEDERLAND

Arjo BV
Biezenwei 21
4004 MB TIEL
Postbus 6116
4000 HC TIEL
Tel: +31 (0) 344 64 08 00
Fax: +31 (0) 344 64 08 85
E-mail: info.nl@arjo.com

NEW ZEALAND

Arjo Ltd
34 Vestey Drive
Mount Wellington
NZ-AUCKLAND 1060
Tel: +64 (0) 9 573 5344
Free Call: 0800 000 151
Fax: +64 (0) 9 573 5384
E-mail: nz.info@Arjo.com

NORGE

Arjo Norway AS
Olaf Helsets vei 5
N-0694 OSLO
Tel: +47 22 08 00 50
Faks: +47 22 08 00 51
E-mail: no.kundeservice@arjo.com

ÖSTERREICH

Arjo GmbH
Lemböckgasse 49 / Stiege A / 4.OG
A-1230 Wien
Tel: +43 1 8 66 56
Fax: +43 1 866 56 7000

POLSKA

Arjo Polska Sp. z o.o.
ul. Ks Piotra Wawrzyniaka 2
PL-62-052 KOMORNIKI (Poznań)
Tel: +48 61 662 15 50
Fax: +48 61 662 15 90
E-mail: arjo@arjo.com

PORTUGAL

Arjo em Portugal
MAQUET Portugal, Lda.
(Distribuidor Exclusivo)
Rua Poeta Bocage n.º 2 - 2G
PT-1600-233 Lisboa
Tel: +351 214 189 815
Fax: +351 214 177 413
E-mail: Portugal@arjo.com

SUISSE / SCHWEIZ

Arjo AG
Fabrikstrasse 8
Postfach
CH-4614 HÄGENDORF
Tél/Tel: +41 (0) 61 337 97 77
Fax: +41 (0) 61 311 97 42

SUOMI

Arjo Scandinavia AB
Riihitontuntie 7 C
02200 Espoo
Finland
Puh: +358 9 6824 1260
E-mail: Asiakaspalvelu.finland@arjo.com

SVERIGE

Arjo International HQ
Hans Michelsensgatan 10
SE-211 20 MALMÖ
Tel: +46 (0) 10 494 7760
Fax: +46 (0) 10 494 7761
E-mail: kundservice@arjo.com

UNITED KINGDOM

Arjo UK and Ireland
Houghton Hall Park
Houghton Regis
UK-DUNSTABLE LU5 5XF
Tel: +44 (0) 1582 745 700
Fax: +44 (0) 1582 745 745
E-mail: sales.admin@arjo.com

USA

Arjo Inc.
2349 W Lake Street Suite 250
US-Addison, IL 60101
Tel: +1 630 307 2756
Free: +1 800 323 1245 Institutional
Free: +1 800 868 0441 Home Care
Fax: +1 630 307 6195
E-mail: us.info@arjo.com

JAPAN

Arjo Japan K.K.
東京都港区虎ノ門三丁目7番8号 ランディック
第2虎ノ門ビル9階
電話 : +81 (0)3-6435-6401

Address page - REV 24: 04/2019

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.



ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden
www.arjo.com

arjo

CE
2797